
**SIEMENS**

# Sirona M1

Schaltunterlagen  
Wiring References  
Schémas de circuits  
Colección de esquemas

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## Blatt

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## List of contents

## Sheet

<b>AE</b> Electrical circuit diagram dentist element	1/2
<b>ST</b> Electrical circuit diagram patient's chair	3
<b>HE</b> Electrical circuit diagram assistant element with amalgam-separation	4
<b>HE</b> Electrical circuit diagram assistant element without amalgam-separation	4 rear
Explanation of symbols for circuit diagrams	5
Hydraulic-pneumatic diagram	5 rear
Explanation of symbols for hydraulic-pneumatic diagram	5 rear
Adjusting instruction	6

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









Legend of electric circuit diagram

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- Legend for hydropneumatic diagram**
- AE Dentist element**  
AA1 Fast coupling Sprayvit 4000  
AA2 Hose coupling for burr drives and ultrasound handpiece  
AF Water filter  
AO Sound absorber / oil separator  
M Motor with spray regulation  
MV22 Solenoid valve for water  
T Turbine with spray regulation  
US Ultrasound handpiece
  - BS Burr drive assembly**  
AC Double non-return valve  
AR4 Return suction device  
MV21.1-4 Solenoid valve for driving air  
PV1 Pneumatic valve for driving air  
PV2 Pneumatic valve for spray  
PV3 Pneumatic valve for chip-blower  
PVS Pneumatic valve for manometer air  
RD3 Control choke for driving air
  - BT Connection block**  
MV11 Driving air valve  
MV12 Chip-blower valve  
MV13 Spray air valve  
RD1 Spray air adjusting choke  
DA Connection for manometer
  - GA1 Unit adapter**  
H Waterheater  
R1 Heater cartridge
  - H R1 Waterheater**  
Heater cartridge
  - SB Sprayvit-block**  
Water filter  
AF Flow switch for air  
DA Solenoid valve for mixer  
TM Thermostatic mixer  
MV23 MS

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- Assistant element**  
HE Fast coupling Sprayvit 4000  
AA1 Sprayvit 4000  
AZ Distributor  
AP Pump for automatic separator  
AC Check valve  
M Pump motor  
AW Amalgam separator  
DZ Switch for solid matter separator  
FA Solid matter separator  
BC Water jet injector  
MV37 Solenoid valve for water jet injector  
BK Solid matter collector (Syphon)  
AC Check valve  
GA3 Connection flange  
H Water heater  
AA5 Fast-coupling for round-flushing and hydrocolloid  
AA6 Fast-coupling for tumbler-filling tube  
AE Nozzle for heater overflow  
AH Tumbler-filling tube  
AT Round-flushing on cuspidor bowl  
MV34 Solenoid valve for tumbler filling  
MV35 Solenoid valve for round-flushing  
R1 Heater cartridge  
BN Turbine for round-flushing and hydrocolloid  
RD35 Flow limiter 1 l/min.  
SD34 Flow limiter 2 l/min.

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- Automatic separator**  
AA3 Fast-coupling saliva ejector  
AA4 Fast-coupling suction hose  
AB5 Mechanical valve for suction hose (only for assistant element without amalgam separator)  
AF4 Coarse filter  
AF7 Filter for suction hose  
AJ1 Suction cannula  
AL Saliva ejector  
BP Separating vessel  
BQ Selector valve  
DS Level switch in separating vessel  
MV33 Suction control valve  
SV30 Control valve for suction hose (without function)  
SV31 Control valve for saliva ejector  
SV32 Control valve for suction hose

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- Connection box**  
AB1 Mechanical valve for air  
AB2 Mechanical valve for water  
AD1 Pressure control air  
AD2 Pressure control water  
AD5 Diffence pressure control air  
AF1 Filter for air  
AF2 Filter for water  
MA1 Manometer for air  
MA2 Manometer for water  
MV1 Solenoid valve for air  
MV2 Solenoid valve for water
- Purification system**  
AB3 Safety valve for venting  
AC1 Non-return valve  
AC2 Non-return valve  
AD3 Safety valve  
AF3 Filter for vaccine  
AN Level switch in tank  
MV3 Solenoid valve for vaccination block  
T Tank for purification agent

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## Adjusting instructions for control electronics ST on PC-board MD

Solder jumper from BR2 to BR1.

Use Multimeter Ri  $\geq 20\text{ k}\Omega$

Remove motor from receptacle and unscrew from the supply hose.

Select 100 % on foil button U.

Connect meter to MP7- / MP8 + . Operate footswitch. Adjust R3 to obtain a reference voltage

$U_{ref} = \text{minus } 5\text{ V}$ .

Connect meter to MP10 / MP11. Operate footswitch. Adjust R4 to obtain a voltage

$U_{max} = 24\text{ V DC}$ .

Select 1 % on foil button U. Operate footswitch. Adjust R5 to obtain a voltage

$U_{min} = 1.5\text{ V DC}$ .

Solder jumper from BR1 to BR2.

## Adjusting instructions for Siroson electronics US on PC-board MD

Use Multimeter Ri  $\geq 20\text{ k}\Omega$ .

Remove fuse F1 on board MD. Take US -handpiece out of receptacle. Select 1 % on foil button U.

Connect meter to MP10 / MP11. Operate footswitch. Adjust R7 to obtain a voltage

$U_{US} = 14\text{ V DC}$ .

Re-insert fuse F1.

## Adjusting instructions for Sirotom electronics HF on PC-board MD

Use Multimeter Ri  $\geq 20\text{ k}\Omega$

Remove fuse F1 on board MD. Take HF -handpiece out of receptacle. Select 1 % on foil button U.

Connect meter to MP10 / MP11. Operate footswitch. Adjust R6 to obtain a voltage

$U_{HF} = 6\text{ V DC}$ .

Re-insert fuse F1.

## Adjusting instructions for illumination on drilling instrument on PC-Board MD

The dental instrument lighting can be adjusted to remain lit after drilling is stopped : Set this post-lit period between approx. 5 sec. and 30 sec. with R2.


## Adjusting instructions for minimum drilling speed on PC-board MS

Select 1 % on foil button U.


Set potentiometer R105 to min.

Actuate footswitch FS.

Turn potentiometer R105 to the right, until the motor begins to run.

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