
Neu ab:
New since: 09.2008
Nouveau depuis:
Nuevo desde:

sirona
The Dental Company

C3⁺, C4⁺

Schaltunterlagen
Wiring References
Schémas de circuits
Colección de esquemas

! ACHTUNG

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! ATTENTION

Interference with electromedical devices by radio telephones:
To guarantee the operational safety of electromedical devices, the operation of mobile radio telephones in the medical practice or hospital area is prohibited.

! ATTENTION

When opening the equipment:
Please observe the safety measures for handling PC boards. Touch a ground point to remove any personal electrostatic charge before touching the components.

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Modification compared with last edition:
09.2007

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List of Contents Sheet

Electrical circuit diagram:
 Patient's Chair 1/1A
 Dentist Element C3+ 2
 Dentist Element C4+ 2A
 Assistant Element 3
 Water Unit 4
 Water Unit, option tumbler filling 4A
 Video Dentist Element C3+ 5/5A
 Video Dentist Element C4+ 5B/5C
 Video Assistant Element C3+, C4+ 6/6A
 CEREC Chairline C4+ 6B
 Legend 7A

Hydraulic-pneumatic diagram:
 C3+, C4+ Cart/Cabinet with disinfection system 8
 C3+, C4+ Cart/Cabinet without disinfection system 8A
 C4+ with disinfection system 9
 C4+ without disinfection system 9A
 Legend 10A

CAN BUS, data communication 11
 RS232 Interface 12A

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V107 LED +5V Ref.

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L1, L2, .Leitungen**1, 2, . . Ferritkerne**

Legend for electrical circuit diagram**ST Patient chair**

F1	Power fuse
S1	Power switch
M1	Motor for lift
M2	Motor for backrest
M3	Motor for headrest tilt
M5	Motor for headrest up/down
MV1	Solenoid valve for water
MV35	Solenoid valve for cuspidor rinsing
MV P1	Valve for driving air
MV P2	Valve for chip blower
MV P3	Valve for spray air
R1	Transmitter pot. for lift
R2	Transmitter pot. for backrest
R3	Transmitter pot. for headrest tilt
S1	Limit switch lift up
S2	Limit switch lift down
S3	Limit switch backrest up
S4	Limit switch backrest down
S5	Limit switch headrest up
S6	Limit switch headrest down
SS1	Safety switch for backrest
SS2	Safety switch for base
SS3	Safety switch for headrest

CF Fuse board**CJ Control board ST**

V401	LED shifting track/backward
V403	LED shifting track/forward
V405	LED coupling/shifting track
V439	LED headrest up
V440	LED headrest down
V441	LED headrest in
V442	LED headrest out
V443	LED M2 backrest up
V444	LED M2 backrest down
V445	LED M1 lift up
V446	LED M1 lift down
V447	LED limit switch
V448	LED programming mode
V449	LED switch activated
V450	LED emergency shutdown active
V451	LED on if a safety switch is active
V452	LED D-Error, a motor does not react
V453	LED +32V
V454	LED +5V
V455	LED +24V
V473	LED relay 115V
V474	LED Reset

CK Switch board for tilting part

S1-S8 Chair movements

CSF Footswitch

S2-S4, S6	Cursor
S5	Switch activate instruments
S7	Switch Chip blower
S8	Switch Spray

FSB Bluetooth Basisboard

V1 Footswitch active

FT Foot button for Sirocam C in assistant element**GA1 Connection flange AE for C4⁺**

X1 Plug

KS Foot control

S9	Switch chair up ↑
S10	Switch chair down ↓
S11	Switch back or headrest up ↑
S12	Switch back or headrest down ↓

MD Monitor**NM Power pack/Monitor****SL Wiring board****SR Relay board**

K1	Relay for key Δ
K2	Relay for key #
K3	Relay for suction machine

SVI**SVC Video adapter**

V11	LED +5V
V97	LED +12V
V98	LED display off

T1 Transformer**AE Dentist's element**

MV21.1-4	Solenoid valve for bur drive module
MV22	Solenoid valve f. ultrasonic handpiece water
MV23	Solenoid valve f. ultrasonic handpiece air
MV24	Solenoid valve for SPRAYVIT water
MV25	Solenoid valve for SPRAYVIT air
MV6	Solenoid valve for brake C3 ⁺
S _{TA} 1/2	Switch for brake C3 ⁺

A7C Camera module C3⁺

V1	LED +5V
V11	LED +5V
V14	LED +12V
S1	Selector switch: Instrument holder or external holder for camera

AC7 Camera module C4⁺

V1	LED +5V
V11	LED +5V
V14	LED +12V

ADC Connection board to control module (ADX)

V3	LED +5V
V4	LED +5V

ADX Control module AE

V6	LED salt solution pump ON
V12	LED fuse salt solution pump defective
V27	LED Sprayvit water heater relay on
V30	LED Sprayvit air heater relay on
V41	LED reset
V49	LED MV22 US water
V57	LED claw-type relay K4 for motor voltage
V60	LED claw-type relay K3 for motor voltage
V63	LED claw-type relay K2 for motor voltage
V69	LED MV21.1 bur drive module 1
V74	LED MV21.2 bur drive module 2
V79	LED MV21.3 bur drive module 3
V84	LED MV23 US air
V89	LED MV24 Sprayvit water
V98	LED MV25 Sprayvit air
V110	LED MV21.4 bur drive module
V115	LED +24V
V116	LED +5V
V117	LED +35V
V134	LED +15V
V135	LED +8V
V136	LED +15V
V141	LED heater ON
V145/146	LED without function

AL Module for Electromotor

V30	LED reset
V31	LED motor runs (run)
V32	LED error (error)
V34	LED +32V
V35	LED +8V

APC3⁺/**APC4⁺ Control panel**

V50	LED +5V
V51	LED reset
V52	LED service

AU**Module US - Handpiece SIROSON**

V1 LED flashes 5x: standby ok

GA1 Connection flange ST for C4⁺

X1 Plug

H**Water heating**

F1	Thermal fuse
R1	Heater cartridge
TS1	Temperature sensor

IR4**Light barrier board AE for C4⁺**

V1/11	Light barrier claw 1 SPRAYVIT
V2/12	Light barrier claw 2
V3/13	Light barrier claw 3
V4/14	Light barrier claw 4
V5/15	Light barrier claw 5 (on right)

KSP**Salt solution pump****LCD****Liquid crystal display****MI****SATELEC MiniLED DC/DC module**

PI	Light barrier board AE for C3 ⁺
V1	Light barrier claw 1 SPRAYVIT
V2	Light barrier claw 2
V3	Light barrier claw 3
V4	Light barrier claw 4
V5	Light barrier claw 5 (on right)

RF**Board for X-ray image viewer****ROEBI****X-ray image viewer****SP**

AR2	Air heater cartridge
WR2	Air heater cartridge
AR	Signal generator air
WR	Signal generator water

HE**Assistant's element**

MV31	Solenoid valve saliva ejector precontrol
MV32	Solenoid valve vacuum hose precontrol
MV36	Solenoid valve surgical cannula precontrol
MV38	Solenoid valve for SPRAYVIT water
MV39	Solenoid valve for SPRAYVIT air

HC7**Camera module HE**

V1	LED +24V
V11	LED +5V
V14	LED +12V
V42	LED Camera on
V46	LED image stop

HC**Control board HE**

V7	LED +5V
V14	LED +24V
V22	LED reset
V26	LED Sprayvit light

HEP**Control panel****HP****Supply board HE****MI****SATELEC MiniLED DC/DC module****PL****POLYLUX**

V	Fan
H	Lamp
V1	LED +5V supply voltage CPU
V39	LED +12V
V40	LED +5V operating voltage

RK**Light barrier board HE**

V1/11	Light barrier claw 4
V2/12	Light barrier claw 3
V3/13	Light barrier claw 2 SPRAYVIT
V4/14	Light barrier claw 1

SP**SPRAYVIT**

AR2	Air heater cartridge
WR2	Air heater cartridge
AR	Signal generator air
WR	Signal generator water

WE**Water unit**

MV43	Solenoid valve for cuspidor movement
SS1	Safety switch for cuspidor
SS2/SS3	Safety switch for HE

H1/2**X-ray image viewer**

Lamps

AP**Pump for automatic separator**

M Motor for pump and centrifuge

AW**Amalgam separator**

S3 Safety switch for amalgam separator

DNA**Dürr Wet Suction**

MSBV Mouth-rinsing basin valve

DS Level switch

EAD**Disinfection system with DENTOSEPT**

DS31 Level switch in mixing tank

DS32 Level switch in tank

MV3 Solenoid valve for DENTOSEPT

H**Water heating**

F2	Thermal fuse
MV34	Solenoid valve for tumbler filling
R2	Heater cartridge
TS2	Temperature sensor

P**Pump**

DS33/34 Protective gas contact for pump control

MV40 Solenoid valve for pump

SA**Automatic separator**

DS30	Level switch in separator tank
V1	Speed transmitter light barrier
MV33	Precontrol solenoid valve f. unit selector valve

SC**Flushing vessel for automatic separator**

MV41 Solenoid valve for flushing system

WM**Water alarm**

DN30 Pressure switch

MV37 Solenoid valve for drain flap

WE**Control board WE**

V40 LED for MV33

V41 LED for water heater for tumbler filling

V69 LED on=amalgam separator vessel empty, level switch DS30

V70 LED on=amalgam rotor closed, safety switch S3

V71 LED on=no water, pressure switch DN30

V72/73 LED for protective gas contact DS33/34 water pump limit switches

V74 LED on=tank filled, for level switch DS32

V75 LED without function

V76 LED on=tank filled, for level switch DS31

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Legend for hydraulic-pneumatic circuit diagram**AK Connection box****ANB Connection block AK**

AB	Mechanical valve for air and water
AD1/2	Pressure controller for air
AD3	Pressure controller for water
AF1	Filter for air
AF2	Filter for water
MA1/2	Pressure gauge for air
MV1	Solenoid valve for water
MV35	Solenoid valve for centrifugal rinsing

VB Valve block

MV P1	Valve for driving air
MV P2	Valve for chipblower
MV P3	Valve for spray air

SDI Plug connection**ZG Assembly for additional units****WE Water unit**

AC37	Check valve
AO33	Acoustic damper
RD35	Flow reglator
V1-V3	Three-way/four-way manifold

AP Pump for automatic separator

AC31	Check valve
M	Motor for pump and centrifuge
SN30	Siphon

AW Amalgam separator

FA	Solids separator with switch S3
SN31	Siphon

BK Solid matter collector

SN31	Siphon
------	--------

DNA Dürr Wet Suction

MSBV	Mouth-rinsing basin valve
DS	Level switch

EAD Disinf. system with DENTOSEPT

AF3	Filter for DENTOSEPT
AO34	Acoustic damper
DS31	Level switch in mixing tank
DS32	Level switch in tank
E	Funnel with vent for DENTOSEPT
MB	Mixing tank
MV3	Solenoid valve for DENTOSEPT
SD36	Flow limiter
T	Tank for DENTOSEPT

H Water heater

MV34	Solenoid valve for tumbler filling
R2	Heater with temperature sensor TS2 and overheating protection

SWS Swivellable cuspidor

AC43	Check valve
AE	Nozzle for heater overflow
AF4	Rough filter
AH	Cuspidor
AT	Tumbler filling pipe
AO43	Acoustic damper
BN	Cuspidor flushing
D43	Regulator
MV43	Solenoid valve for cuspidor movement
PZ43	Pneumatic cylinder
STA	Plug adapter for tumbler filling/hydrocolloid

P Pump

AC32/33	Check valve
DS33/34	Protective gas contact for pump control
EV30	Dewatering valve
MV40	Solenoid valve for pump

SA Automatic separator

AF33	Filter
BP	Separator tank with light barrier V1
BQ	Unit selector valve
DS30	Level switch in separator tank
MV33	Precontrol solenoid valve f. unit selec. valve
RD33	Flow reglator

SC Flushing vessel for automatic separator

MV41	Solenoid valve for flushing system
SD41	Flow limiter

WM Water alarm

AO31	Acoustic damper
DN30	Pressure switch
MV37	Solenoid valve for drain flap
PK	Pneumatically operated drain flap

AE Dentist's element

M	Motor SL
T	T1 CONTROL
US	SIROSON L ultrasonic handpiece
AA1	Hose coupling for SPRAYVIT
AA2	Hose coupling for bur drives
AA10	Hose coupling for ultrasonic handpiece
AD 4	Pressure controller for spray air
DM21	Pressure alarm, measuring air

AB Connection block AE

MV24	Solenoid valve for SPRAYVIT water
MV25	Solenoid valve for SPRAYVIT air
AO26-28	Acoustic damper

BS Bur drive module

AC	Double check valve
AC1	Check valve for bur drive
AR4	Back suction device
MV21.1-4	Solenoid valve bur drive module
PV1	Pneumatic valve for driving air
PV2	Pneumatic valve for spray
PV3	Pneumatic valve for chip-blower
PV5	Pneumatic valve for pressure gauge air
RD3/5	Control valve for driving air

H Water heater

AF	Water filter
R1	Heater with temperature sensor TS1 and overheating protection

US Ultrasonic block

MV22	Solenoid valve for water US
NV23	solenoid valve for air US

VA Supply connection**HE Assistant's element**

AA1	Fast-action coupling for SPRAYVIT
AA3	Fast-action coupling for saliva ejector
AA4	Fast-action coupling for vacuum hose
AC34	Check valve for vacuum hose
AC35	Check valve for saliva ejector
AF5/7/8	Filter
AJ1	Suction cannula
AL1	Saliva ejector
MV31	Precontrol valve for saliva ejector
MV32	Precontrol valve for vacuum hose
MV38	Solenoid valve for water SPRAYVIT
MV39	Solenoid valve for air SPRAYVIT
SV31	Control valve for saliva ejector
SV32	Control valve for vacuum hose

CK Surgical cannula

AA7	Fast-action coupling for surgical cannula
AC36	Check valve for surgical cannula
AF5/9	Filter
MV36	Precontrol valve for surgical cannula
SV36	Control valve for surgical cannula

TA Support arm C3⁺

AO6	Acoustic damper
MV6	Solenoid valve for brake
PR	Brake

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RS232 Interface**Dental workstation → PC**

An extended Mouse System Protocol is used to transmit defined conditions or commands from the dental workstation to a connected PC. For the camera and foot switch conditions in byte 1, the bits 0-2 are used to control the video application, to transmit additional commands, the bits 3-6 are used. The bytes 2-5 are 0 and are not used.

7	6	5	4	3	2	1	0
1	User command			Video function			

Byte 1 of the Mouse System

In each case only one user command or information to the video function can be transmitted in one data package. The bits of the other block are set to 0.

The following events and interface commands are defined:

The commands are generated in each case when the event is triggered.

Event	Remark	RS232 command 32 (1st byte of the 5- byte protocol)
Rhomb key activated		88H
Rhomb deactivated	only if switching mode	90H
Call key		98H
X-ray image viewer key		A0H
CCW rotation key (with instruments deposited or video camera removed)		A8H
Instrument light key (with instruments deposited or video camera removed)		B0H
Quick setting key 100 (with instruments deposited or video camera removed)	also be operated via FS	B8H
Quick setting key 75 (with instruments deposited or video camera removed)	cannot be operated via FS	C0H
Quick setting key 50 (with instruments deposited or video camera removed)	cannot be operated via FS	C8H
Quick setting key 25 (with instruments deposited or video camera removed)	cannot be operated via FS	D0H
Quick setting key 1 (with instruments deposited or video camera removed)	also be operated via FS	D8H
Forward cursor from chair key block	only be operated via FS	E0H
Forward cursor from instrument key block	only be operated via FS	E8H
Spray key Left discrete C+ FS button active with camera removed		F0H
Chip blower key Right discrete C+ FS button active with camera removed		F8H
Camera removed		81H
Camera deposited		83H
Foot switch active with camera removed		85H
Foot switch inactive with camera removed		81H
When video camera is assembled on HE side and when in program SIUCOM.exe mode: "Sirocam C direct" is selected:		
Press foot switch		E8H
Press button on camera		A0H

PC → Dental workstation

The RS232 configuration and service information can be output via PC with an 1-byte command.

Command	Action	Remark
0DH	Output of configuration and service information	

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Property of Sirona. All rights reserved.
Sirona Dental Systems GmbH