



Service Guide

IntelliVue Patient Monitor

MX600/MX700

Release J.x

Patient Monitoring

PHILIPS

Table of Contents

| | |
|--|-----|
| 1 Introduction | 7 |
| <hr/> | |
| Who Should Use This Guide | 7 |
| How to Use This Guide | 7 |
| Abbreviations | 8 |
| Responsibility of the Manufacturer | 8 |
| Passwords | 9 |
| Warnings and Cautions | 9 |
| 2 Theory of Operation | 11 |
| <hr/> | |
| Integrated Monitor Theory of Operation | 11 |
| 3 Testing and Maintenance | 29 |
| <hr/> | |
| Introduction | 29 |
| Terminology and Definitions | 29 |
| Recommended Frequency | 30 |
| When to Perform Tests | 31 |
| Testing Sequence | 35 |
| Visual Inspection | 35 |
| Safety Tests | 36 |
| System Test | 57 |
| Preventive Maintenance Procedures | 69 |
| Performance Assurance Tests | 69 |
| Reporting of Test Results | 97 |
| Other Regular Tests | 101 |
| Touchscreen Calibration | 101 |
| Disabling/Enabling Touch Operation | 102 |
| Printer Test Report | 103 |
| After Installation, Testing or Repair | 103 |
| 4 Troubleshooting | 105 |
| <hr/> | |
| Introduction | 105 |
| How To Use This Section | 105 |
| Who Should Perform Repairs | 105 |
| Replacement Level Supported | 105 |
| Software Revision Check | 106 |
| Software Compatibility Matrix | 106 |
| Obtaining Replacement Parts | 108 |
| Troubleshooting Guide | 108 |

| | |
|---|-----|
| 5 Repair and Disassembly | 143 |
| <hr/> | |
| Tools Required | 143 |
| Monitor Disassembly | 143 |
| iPC Disassembly | 172 |
| 8-Slot Flexible Module Rack (FMS-8) Disassembly | 175 |
| 4-Slot Flexible Module Rack (FMS-4) Disassembly (without MMS Mount) | 181 |
| 4-Slot Flexible Module Rack (FMS-4) Disassembly (with MMS Mount) | 185 |
| Plug-in Modules | 192 |
| Multi-Measurement Module (MMS) Disassembly | 195 |
| MMS Extensions - Exchanging the Top Cover, MSL Flex Cable and the Dual Link Bar | 211 |
| Disassembly Procedures for the M3015A MMS Extension (HW Rev. A) | 222 |
| 6 Parts | 227 |
| <hr/> | |
| MX600/MX700 Parts | 228 |
| Remote Control Parts | 232 |
| 8-Slot Flexible Module Rack (FMS-8) Parts | 233 |
| 4-Slot Flexible Module Rack (FMS-4) Parts | 235 |
| Multi-Measurement Module (MMS) Parts | 237 |
| MMS Extension Parts (M3012A, M3014A, M3015A/B) | 244 |
| IntelliVue X2 Part Numbers | 246 |
| Plug-in Modules Part Numbers | 247 |
| External Display Part Numbers | 259 |
| Test and Service Tools | 260 |
| 7 Installation Instructions | 263 |
| <hr/> | |
| Electromagnetic Emissions | 263 |
| Electromagnetic Interference (SRR) | 264 |
| Installation Checklist | 264 |
| Unpacking the Equipment | 264 |
| Initial Inspection | 265 |
| Installing the MX600/700 Monitor | 266 |
| Connecting the Monitor to AC Mains | 282 |
| 8 Site Preparation | 315 |
| <hr/> | |
| Introduction | 315 |
| Monitor MX600/700 Site Requirements | 318 |
| Electrical and Safety Requirements (Customer or Philips) | 319 |
| Remote Device Site Requirements | 319 |
| RS232/MIB/LAN Interface | 327 |
| Nurse Call Paging Cable | 328 |
| ECG Out Interface | 329 |
| 9 Gas Analyzers | 331 |
| <hr/> | |
| 10 Specifications | 333 |
| <hr/> | |
| Essential Performance | 333 |



| | |
|--|-----|
| 11 Appendix - Philips 15210B Calibration Unit | 337 |
| <hr/> | |
| Unpacking the Instrument | 337 |
| Routine Maintenance | 339 |
| Theory of Operation | 340 |
| Gas Flow Performance Check | 340 |
| Disassembly | 342 |
| Parts List | 343 |
| 12 IntelliVue MX600/700 Product Structure | 345 |
| <hr/> | |
| IntelliVue MX600/700 Upgrade Structure | 350 |



Introduction

This Service Guide contains technical details for the IntelliVue MX600/700 Patient Monitor, the measurement modules, the Multi-Measurement Module (MMS), the IntelliVue X2, the Flexible Module Rack (FMS) and the Measurement Server Extensions.

This guide provides a technical foundation to support effective troubleshooting and repair. It is not a comprehensive, in-depth explanation of the product architecture or technical implementation. It offers enough information on the functions and operations of the monitoring systems so that engineers who repair them are better able to understand how they work.

It covers the physiological measurements that the products provide, the Measurement Server that acquires those measurements, and the monitoring system that displays them.

Who Should Use This Guide

This guide is for biomedical engineers or technicians responsible for troubleshooting, repairing, and maintaining Philips' patient monitoring systems.

How to Use This Guide

This guide is divided into ten sections. Navigate through the table of contents at the left of the screen to select the desired topic. Links to other relevant sections are also provided within the individual topics. In addition, scrolling through the topics with the page up and page down keys is also possible.

Abbreviations

Abbreviations used throughout this guide are:

| Name | Abbreviation |
|--------------------------------------|---------------------|
| IntelliVue MX600/700 Patient Monitor | the monitor |
| Flexible Module Rack | FMS |
| Multi-Measurement Module | MMS |
| Measurement Link | MSL |
| Medical Information Bus | MIB |
| Anesthetic Gas Module | AGM |

Responsibility of the Manufacturer

Philips only considers itself responsible for any effects on safety, EMC, reliability and performance of the equipment if:

- assembly operations, extensions, re-adjustments, modifications or repairs are carried out by persons authorized by Philips, and
- the electrical installation of the relevant room complies with national standards, and
- the instrument is used in accordance with the instructions for use.

To ensure safety and EMC, use only those Philips parts and accessories specified for use with the monitor. If non-Philips parts are used, Philips is not liable for any damage that these parts may cause to the equipment.

This document contains proprietary information which is protected by copyright. All Rights Reserved. Reproduction, adaptation, or translation without prior written permission is prohibited, except as allowed under the copyright laws.

Philips Medizin Systeme Böblingen GmbH

Hewlett-Packard Str. 2

71034 Böblingen, Germany

The information contained in this document is subject to change without notice.

Philips makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties or merchantability and fitness for a particular purpose.

Philips shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

Passwords

In order to access different modes within the monitor a password may be required. The passwords are listed below.

Monitoring Mode: No password required

Configuration Mode: 71034

Demo Mode: 14432

Service Mode: 1345

Consult the configuration guide before making any changes to the monitor configuration.

iPC Passwords:

The passwords for the pre-installed user accounts should be changed to ensure appropriate security.

The default passwords and settings are shown in the following table.

| User Name | ClinicalUser | ServiceUser |
|------------------|----------------------|------------------------|
| Long Name | Clinical User | Service User |
| Password | BedPC | Supp0rtuzr (0 is zero) |
| Password expires | Never | Never |
| Group | Users | Administrators |
| Comment | For normal operation | For service purposes |

Warnings and Cautions

In this guide:

- A **warning** alerts you to a potential serious outcome, adverse event or safety hazard. Failure to observe a warning may result in death or serious injury to the user or patient.
- A **caution** alerts you where special care is necessary for the safe and effective use of the product. Failure to observe a caution may result in minor or moderate personal injury or damage to the product or other property, and possibly in a remote risk of more serious injury.

Theory of Operation

Integrated Monitor Theory of Operation


The IntelliVue MX600/700 Patient Monitor:

- displays real-time data
- controls the attached measurement servers
- alarms in the case of patient or equipment problems
- offers limited data storage and retrieval (trending)
- interfaces to the Philips Clinical Network and other equipment

A monitor with just a single integrated measurement server can be connected to additional building blocks to form a monitoring system with a large number of measurements, additional interface capabilities and multiple slave displays. These elements cooperate as one single integrated real-time measurement system.

System Boundaries

The following diagram discusses specific boundaries within the overall system with respect to their openness and real-time requirements:



System Boundaries

| | |
|---|---|
| ● | Measurement connections Built-in measurement block |
| ○ | Philips Clinical Network (wired LAN) connects multiple patient monitors, information centers, application servers; closed system, only Philips qualified products (tested and with regulatory approval) are connected, Philips is responsible for guaranteed real-time functionality and performance |
| ⋯ | Philips Clinical Network (wireless) like Philips Clinical Network (wired) LAN, however due to current wireless technologies available it has reduced bandwidth, longer latencies, reduced functionality |
| ○ | Hospital LAN, Internet Standard Network, not under Philips control, no guaranteed service, no real-time requirements |


Hardware Building Blocks

The following hardware building blocks make up the monitoring system:

IntelliVue MX600/700

The IntelliVue MX600/700 Monitor:

- integrates the display and processing unit into a single package
- uses a 15" TFT WSXGA+ Color display
- The MX700 uses the touch screen as primary input device; a remote control and computer devices such as mice, trackball, and keyboard can be added optionally. The MX600 uses the trimknob as primary input device. It does not have a touch screen; a remote control and computer devices such as mice, trackball, and keyboard can be added optionally.
- supports the Flexible Module Rack (FMS)



NOTE

The 802.11 Bedside Adapter (WLAN) and IIT are mutually exclusive.


Compatible Devices



M8048A 8-Slot Flexible Module Rack (FMS-8)




865243 4-Slot Flexible Module Rack (FMS-4) (without and with MMS Mount)



M3001A Multi-Measurement Module (MMS)



M3002A IntelliVue X2




M3012A, M3014A, M3015A, M3016A MMS Extensions



865244 Remote Control

Power Supply



The AC/DC converter transforms the AC power coming from the plug into 24V/ 150W DC source and isolates the monitoring system from the AC power mains. The 24V is distributed via power bus either directly or over additional converters to all components of the system:

The 56V DC power needed for the MSL is created by an isolating DC/DC converter.

The iPC and the backlights converter are supplied with 24V. The power needed for the panel is converted to 12V by the panel adapter DC/DC converter.


The isolated interfaces are supplied with 2.5V AC. The main board is supplied with 5V, 3.3V and 1.8V.

Additionally, for some infrastructural functions 3.6V is provided to the main board

The mixer board requires a power of 3.3V, 2.5V and 1.2V.

CPU Boards

The CPU boards have an MPC860 50 MHz or MPC86x 100 MHz processor that provides a number of on-chip, configurable interfaces. An array of 12 fast UARTS with configurable protocol options are implemented in an ASIC (along with other system functions such as independent watchdogs etc.), providing interfacing capabilities to measurement modules and I/O boards. The serial interfaces can easily be electrically isolated. The main board contains additional video hardware.



The CPUs provide two LAN interfaces to interconnect CPUs (via the MSL) and to connect to the Philips Clinical Network.

The CPU capabilities are identical. Different loading options are coded on serial EEPROMs to support the automatic configuration of the operating system at boot time.

iPC


The iPC is a full standard PC that is built into the MX600/700 patient monitors as a hardware option

I/O Boards

A single MIB/RS232 or Flexible Nurse Call Relay board can be added optionally.

Data Flow

The following diagram shows how data is passed through the monitoring system. The individual stages of data flow are explained below.



Data Flow

Data Acquisition

Monitoring data (for example patient measurement data in the form of waves, numerics and alerts) is acquired from a variety of sources:

- **Measurement Servers**
The Measurement Servers connected to the internal LAN convert patient signals to digital data and apply measurement algorithms to analyze the signals.
- **External measurement devices**
Data can be also acquired from devices connected to interface boards of the monitor. Software modules dedicated to such specific devices convert the data received from an external device to the format used internally. This applies to parameter modules and the Anesthetic Gas Module.
- **Server systems on the Philips Clinical Network**
To enable networked applications such as the other bed overview, data can be acquired from server systems attached to the Philips Clinical Network, for example a Philips Information Center

Data Provider System Service

All data that is acquired from measurement servers or external measurement devices is temporarily stored by a dedicated data provider system service. All monitor applications use this central service to access the data in a consistent and synchronized way rather than talking to the interfaces directly.

This service makes the applications independent of the actual type of data acquisition device.

The amount of data stored in the data provider system service varies for the different data types. For example several seconds of wave forms and the full set of current numerical values are temporarily stored in RAM.

Persistent Data Storage System Service

Some applications require storage of data over longer periods of time. They can use the persistent data storage system service. Dependent on the application requirements, this service can store data either in battery backed-up (buffered) memory or in flash memory. The buffered memory will lose its contents if the monitor is without power (not connected to mains) for an extended period of time. The flash memory does not lose its contents.

The trend application for example stores vital signs data in a combination of flash memory and buffered memory, while the system configuration information (profiles) is kept purely in flash memory.

Display and User Interface Service

Applications can use high level commands to display monitoring data or status and command windows on the internal LCD panel. These commands are interpreted by the display manager application. This application controls the dedicated video hardware which includes video memory and a special ASIC.

User input is acquired from a variety of input devices, for example the SpeedPoint, the touchscreen or other standard input devices (keyboard, mouse) which may be attached to I/O boards. The system software makes sure that the user input is directed to the application which has the operating focus.

Data Output

The monitoring system is very flexible and customizable regarding its data output devices. Built-in devices (for example LAN, alarm lamps, speaker, video) provide the basic output capabilities.

These capabilities can be enhanced by adding additional I/O boards, as required in the specific end-user setup. The additional I/O boards typically provide data to externally attached devices, for example to RS232 based data collection devices.

The monitor can identify I/O boards by means of a serial EEPROM device that stores type and version information. The operating system detects the I/O boards and automatically connects them with the associated (interface driver) application. For some multi-purpose cards it is necessary to configure the card for a particular purpose first (for example the MIB/RS232 card can support external touch display , data import, data export).

Monitor Applications

The monitor applications provide additional system functionality over the basic measurement and monitoring capabilities. This includes for example trending, report generating, event storage or derived measurements.

In general, the monitor applications use the data provider system service to access the measurement data. Application interfaces to the other system services allow the application to visualize data, to store data over extended periods of time or to output data to other devices.

Internal LAN (Measurement Link)

All components of the monitoring system (including measurement servers and CPUs in the monitor) communicate using an IEEE802.3/ Ethernet LAN in the Measurement Link (MSL). This network is used to distribute data between the components, for example:


- Digitized patient signals including wave data, numerical data and status information (typically from the measurement server to a display unit)
- Control data representing user interactions (typically from the display unit to a measurement server)
- Shared data structures, for example representing patient demographical data and global configuration items

The internal LAN allows plug and play configuration of the monitoring system. The system automatically detects plugging or unplugging of measurement servers and configures the system accordingly.

The components on the internal LAN are time-synchronized to keep signal data consistent in the system. Dedicated hardware support for synchronization eliminates any latency of the network driver software.

The integrated LAN provides deterministic bandwidth allocation/reservation mechanisms so that the real-time characteristic of signal data and control data exchange is guaranteed. This applies to the data flow from the measurement server to the monitor (for example measurement signal data) and the data flow from the monitor to a measurement server (for example to feed data to a recorder module).

Integrated communication hubs in the monitor and the FMS allow flexible cabling options (star topology, daisy chaining of servers).



Philips Clinical Network

The monitoring system may be connected to the Philips Clinical Network, for example to provide central monitoring capabilities or other network services. This connection may be through a normal wired connection or through a wireless connection.

The monitor supports the connection of an internal wireless adapter (#J35). Switching between wired and wireless networks is automatically triggered by the plugging or unplugging of the network cable.

The Philips Clinical Network protocols function very similarly to the protocols used on the internal LAN.

After configuration, the monitoring system sends the digitized patient signals including wave data, numerical data and status information onto the network. Control data representing user interactions can be exchanged between the monitoring system and a central station bi-directionally.

Additional protocols are supported for networked applications, for example for the other bed overview function, which allows viewing of monitoring data from other patients on the network.

For plug and play operation, the monitoring system uses the standard BootP protocol to automatically acquire a network address.

How does the Support Tool Work with the Monitor

The support tool is a Windows application typically installed on the laptop of a customer engineer or a biomedical engineer working in the customer's own service department.


The purpose of the support tool is to upgrade, configure and diagnose all monitoring components (modules, measurement servers, and monitors) in the system over the network.

The service protocol developed for this purpose uses a raw access to the devices without the need for IP addresses etc. over a standard customer network installation, so that even defective devices can be upgraded as long as the few kBytes of initial boot code are working. The boot code itself can also be upgraded using the same protocol.

The tool allows access to internal service information and to serial numbers. It can be remote-controlled, for example via a dial-up connection from a response center, provided the proper infrastructure is in place.

For details see the Instructions for Use for the Support Tool.

Monitor Software Block Diagram



Block Diagram Legend

| Functional Block | Description |
|---------------------------|--|
| Services | |
| Operating System | The Operating System (OS) provides a layer of isolation between the specific hardware implementation and the application software. The OS performs system checks and allocates resources to ensure safe operation when the system is first started. This includes internal self-tests on several hardware modules and configuration checks for validity of configuration with the operating software. During normal operation, the OS continues to run checks on system integrity. If error conditions are detected the OS will halt monitoring operations and inform the operator about the error condition. |
| System Services | The System Services provide generic common system services. In particular: They use a real-time clock component to track time. They synchronize to network time sources and verify the accuracy of the system time information. They are also responsible for managing persistent user configuration data for all Measurement Servers, Flexible Module Racks and IntelliVue Patient Monitoring System software modules. User configuration data is stored in a non-volatile read/write storage device |
| Applications | |
| Application Server Client | The Application Server Client provides the Citrix1 thin client functionality. |
| Reports | <p>The Reports Service retrieves current and stored physiological data and status data to format reports for printing paper documentation. The following reports are supported:</p> <ul style="list-style-type: none"> • Vital Signs Report • Graphical Trend Report • Event Review Report • Event Episode Report • ECG Report (12 Lead/Multi-Lead) • Cardiac Output Report • Calculations Report (Hemodynamic/Oxygenation/Ventilation) • Calculations Review Report • Wedge Report • Test Report • Other reports (e.g. Loops, Review Applications, Drug report) <p>The Reports service generates report data which can be printed on a local or a central printer.</p> |

| Functional Block | Description |
|------------------|--|
| Record | The Record Service retrieves current and stored physiological data and status data to format a continuous strip recording. A recording can be triggered manually by the operator or automatically by an alarm condition. The Record Service uses the services of the Recorder Interface to control an M1116B Recorder in the FMS. The Record Service can also send data to a central recorder. |
| Alarm | The Alarm Service contains logic that prioritizes alarm conditions that are generated either by the Measurement Servers, Flexible Module Rack, or by IntelliVue Patient Monitoring System software modules. Visual alarm signals (messages) are displayed at the top of the IntelliVue Patient Monitoring System display and alarm sounds are generated by a loudspeaker. Alarm conditions may be generated when a physiological parameter exceeds preselected alarm limits or when a physiological parameter or any other software module reports an inoperative status (technical alarm, for example, the ECG leads may have fallen off the patient). The Alarm service manages the alarm inactivation states, for example suspension of alarms, silencing of alarms, and alarm reminder. Alarm signals may also be configured as latching (alarm signals are issued until they are acknowledged by the operator, even when the alarm condition is no longer true). The Alarm service controls the visual alarm signals (alarm lamps). |
| Trend | The Trend service stores the sample values of physiological data and status data with a resolution of 12 seconds, 1 minute or 5 minutes for a period of up to 48 hours. The data is kept in battery buffered read/write storage and flash memory devices to be preserved across power failures. The stored data is protected via consistency checks and checksums. When a new patient is admitted, the trend database erases all data of the previous patient. |
| OxyCRG | The OxyCRG (Oxygen CardioRespiroGram) service derives a high-resolution trend graph from the Beat-to-Beat Heart Rate, SpO2 or tcpO2, and Respiration physiological data. The OxyCRG is specialized for neonatal applications, allowing the operator to identify sudden drops in Heart Rate (Bradycardia) and SpO2 or tcpO2 (Desaturations), and supporting the operator in visualizing Apnea situations. |
| ADT | The ADT (Admit/Discharge/Transmit) service maintains the patient demographics information. The operator may admit a new patient, discharge the old patient and enter or modify the patient demographics. The ADT service also supports the transport of a patient (trend database) with the M3001A Multi-Measurement Module. The ADT service controls the deletion of old patient data, the upload of trend data from the M3001A and the switching back of all settings to user defaults. It also synchronizes patient information with a central station on the network. |

| Functional Block | Description |
|------------------|---|
| Events | The Events Application captures physiological data from episodes for later review and documentation purposes. Events can be triggered automatically by an alarm condition, by user-defined conditions or manually by the operator. |
| Protocol Watch | ProtocolWatch allows the execution of pre-defined clinical protocols in the IntelliVue patient monitor by combining events such as automatically triggered events, time and manually triggered events with textbook knowledge thus aiding the operator to follow clinical guidelines. ProtocolWatch notifies the operator when certain combinations of clinical conditions occur and it documents the developments and clinician actions in a log which can be reviewed on the monitor and documented on a printer. |
| Calc Param | The Calc Param (Calculated Parameters) service accesses current, stored and manually entered physiological data as input to calculation formulas. With these formulas, derived hemodynamic, oxygenation and ventilation variables are computed. The calculation results, including the input parameters, are stored for later review using the Trend service. |
| Heart Mgr. | The Heart Manager Application allows the selection of the alarming source to be either heart rate (from ECG) or the system pulse rate. The system pulse rate can be chosen from any of the possible pulse rate sources (e.g., SpO2 and invasive pressures). The module implements automatic fall-backs when selected signal sources are not available. |
| Drug Calc | The Drug Calc application aids in calculating drug dosages for patients. |
| AGM/EGM | AGM (Anesthesia Gas Module) and EGM (extensible Gas Module) interface ane-asthesis gas measurement devices. The AGM/EGM Module (Anesthesia Gas Module) interfaces the M1013A, M1019A or M1026A/B Gas Analyzer devices. The AGM Module retrieves the measurement data and controls the external device. It provides numerical data, wave form data and alarm data for the gas parameters measured by the attached analyzers. |
| PV Loops | The PV Loops application compares graphic representations of airway waves to help detect changes in the patient airway condition. |

| Functional Block | Description |
|------------------------------|--|
| Interface Managers | |
| MDSE | <p>The MDSE (Medical Data Service Element) Interface Manager is responsible for the exchange of real-time data between the IntelliVue Patient Monitoring System display unit and the Measurement Servers and Flexible Module Rack as well as between the IntelliVue Patient Monitoring System display unit and other devices attached to the network. MDSE establishes and maintains a data communication link between the devices. It provides configuration information about the remote device to applications in the local device and it allows the exchange of measurement data and status information between the devices.</p> |
| Printer | <p>The Printer Interface Manager provides a high level interface to a printer. It provides means to:</p> <ul style="list-style-type: none"> • establish a connection to the printer • transfer data to the printer • get status of the printer • close connection to the printer <p>The Printer Interface Manager also supervises the connection to the printer and whether the printer accepts data (for example paper out). The Printer Interface Manager notifies the operator in such cases.</p> |
| Display & Operator Interface | <p>The Display and Operator Interface Manager performs the following tasks:</p> <ul style="list-style-type: none"> • Screen presentation of real-time and stored physiological measurement data, alarm condition data and status information received from the MDSE interface manager, the Alarm service or other IntelliVue Patient Monitoring System modules • Screen presentation of operating controls (control windows) • Processing of operating control commands received from HIF Control interface. The module verifies and interprets the received commands and forwards them to other software modules of the IntelliVue Patient Monitoring System display unit, Measurement Servers or Flexible Module Rack • Sound generation (issues audible alarm signals and generates audible information signals, for example QRS and SpO2 tones, operator audible feedback) |
| LabData/Manual Data | <p>The Laboratory Data/ Manual Data Entry Interface Manager allows acquisition of laboratory data (e.g. acquired by the central station from a laboratory information system). It also allows to manually enter measurement data to make additional, manually acquired measurements available to internal applications and to the system.</p> |

| Functional Block | Description |
|------------------------------------|---|
| Wireless Measurement Manager (WMM) | The WMM Interface Manager provides connectivity to the SRR interface. It establishes communication between SRR enabled devices and the ASW module that manages the data provided by the device |
| Interfaces | |
| LAN | The LAN interface implements the physical layer of IEEE 802.3. The LAN interface performs Manchester encoding/decoding, receive clock recovery, transmit pulse shaping, jabber, link integrity testing, reverse polarity detection/correction, electrical isolation, and ESD protection. Electronically separated interfaces are used for communication to the Measurement Servers or Flexible Module Rack and to the network. |
| WLAN | The WLAN Interface is a network interface that provides access to an IEEE 802.11 wireless Local Area Network. The configuration of this interface is done by an OS Service. |
| Display Controller | The Display Controller Interface consists of a video controller chip, video RAM and the controlling software. The Display Controller interface processes the high level display commands (character and graphic generation, wave drawing) and translates them into pixels, which are written into the video RAM where the video controller chip generates the video synchronization signals and the pixel stream for the Color LCD Display. |
| HIF Control | The HIF (Human Interface Control) interface scans the Human Interface devices for operator controls (Touch Screen, and USB devices), formats the collected data and sends it to the display and Operating Interface. |
| ECG-Out | The ECG Out interface receives the ECG waveform directly from the ECG/Resp Arrhythmia ST-Segment physiological algorithm via an RS-422 serial interface and converts the digital ECG signal to an analog ECG signal. |
| Sync Out (ECG) | A pulse signal is provided on the Sync Out connector to allow synchronization with other medical devices. |
| RS-232 | The RS-232 component represents a generic serial communication interface to connect external devices as shown in the diagram, also providing power in some cases. |
| RS-422 | The serial link RS-422 interface communicates the ECG signal to the ECG Output of the IntelliVue Patient Monitoring System display unit. The interface is a serial, differential, full-duplex link. The interface is ESD protected. |
| Nurse Call | The Nurse Call board has a modular jack 6P6C connector. The connector has an open and close contact on alarm. |
| Flexible Nurse Call | The Flexible Nurse Call board has a multi-port connector with three alarm relays which are configurable to be open or closed on alarm. In addition, this interface has an audible alert capability for loss of AC power. |

| Functional Block | Description |
|------------------|---|
| MIB | The MIB interface allows full-duplex, short-haul asynchronous binary communication between the monitor and an arbitrary (medical/non-medical) device using an eight-pin RJ45 modular connector. Switching between MIB and RS232 protocol is possible. |
| IIT Interface | The IIT interface allows operation of the monitors with IntelliVue Instrument Telemetry. |
| SRR | The built-in SRR interface allows wireless communication of the monitor with an IntelliVue Instrument Telemetry Transceiver. |
| USB Interface | The USB interface allows connection of USB devices (Mouse, Keyboard, Barcode Scanner, Printer) to the monitor. |
| iPC | The iPC is a full standard PC that is built into the MX600/700 patient monitors as a hardware option. |
| Remote Control | The remote control allows remote operation of the monitor either wireless or via a USB cable connection. |

Testing and Maintenance

Introduction

This chapter provides a checklist of the testing and maintenance procedures to ensure the performance and safety of the monitor, the Multi-Measurement Module (MMS), the MMS Extensions and the Flexible Module Rack (FMS) associated modules.

These tests must be performed only by qualified personnel certified by the responsible organization. Qualifications required are: training on the subject, knowledge, experience and acquaintance with the relevant technologies, standards and local regulations. The personnel assessing safety must be able to recognize possible consequences and risks arising from non-conforming equipment.

All recurring safety and performance assurance tests must be performed under equal environmental conditions to be comparable.

Preventive Maintenance refers specifically to the series of tests required to make sure the measurement results are accurate. The accuracy and performance procedures are designed to be completed as specified in the following sections or when readings are in question.

For detailed instructions on the maintenance and cleaning of the monitor and its accessories, see *Care and Cleaning*, *Using Batteries* and *Maintenance and Troubleshooting* in the monitor's *Instructions for Use*.

Terminology and Definitions

The following terms and definitions are used throughout this chapter and taken from the international standards IEC 60601-1, IEC 60601-1-1 and IEC 62353.

- **Medical System:** a medical electrical system is a combination of at least one medical electrical device and other electrical equipment, interconnected by functional connection or use of a multiple portable socket-outlet.
- **Patient Environment:** any area in which intentional or unintentional contact can occur between the patient and parts of the medical system or between the patient and other persons who have had contact with parts of the medical system. The patient environment is defined anywhere within 1.5m (5 feet) of the perimeter of the patient's bed and 2.5m (8.2 feet) from the floor.
- **Separation Device/Transformer:** a component or arrangement of components with input parts and output parts that, for safety reasons, prevent a transfer of unwanted voltage or current between parts of a medical system.
- **Multiple Portable Socket-Outlet:** a combination of two or more socket-outlets intended to be connected to or integrated with flexible cables or cords, which can easily be moved from one place to another while connected to the power mains.

3 Testing and Maintenance

- **Functional Connection:** an electrical connection for transfer of signals and/or power.
- **Tests:** Safety or Performance Assurance test procedures which may consist of several steps.

Recommended Frequency

Perform the procedures as indicated in the suggested testing timetable. These timetable recommendations do not supersede local requirements.

Table 1 Table 1: Suggested Testing Timetable

| Tests | | Frequency |
|--|--|--|
| Preventive Maintenance* | NBP Performance | Once every two years, or more often if specified by local laws. |
| | Microstream CO ₂ Calibration | Once a year or after 4000 hours of continuous use and following any instrument repairs or the replacement of any instrument parts. |
| Other Regular Tests | Visual Inspection | Before each use. |
| | Power On Test | |
| Performance Assurance Tests | ECG/Resp Performance | Once every two years, or if you suspect the measurement is incorrect, except Mainstream CO ₂ Accuracy Check, Sidestream CO ₂ Accuracy Check and Flow Check - required once a year. |
| | ECG Out Performance | |
| | SpO ₂ Performance | |
| | NBP Performance | |
| | Invasive Pressure Performance | |
| | Temperature Accuracy | |
| | M3014A Capnography Extension Performance Tests | |
| | Microstream CO ₂ Performance Test | |
| | Spirometry Accuracy Test | |
| | C.O. Performance | |
| | BIS Performance | |
| | NMT Performance | |
| | VueLink Performance | |
| | IntelliBridge Performance Test | |
| | Nurse Call Relay Performance | |
| | MSL Assurance Test | |
| Power Loss Alarm Buzzer Performance Test | | |
| Mounting Integrity Test | | |

| Tests | | | Frequency |
|--------------|-------------|------------------------------|---|
| Safety Tests | Visual | Visual Inspection | After each service event |
| | Electrical | Protective Earth | Once every two years and after repairs where the power supply has been removed or replaced or the monitor has been damaged by impact. |
| | | Equipment Leakage Current | |
| | | Applied Part Leakage Current | |
| | System Test | Once every two years | |

*M3015A with the old hardware Rev. A (i.e. Serial No. DE020xxxxx) also require the CO₂ pump/CO₂ scrubber replacement procedure. This is required every three years or after 15000 operating hours.

NOTE

The EEG, SvO₂ (SO₂) and tcGas parameters do not require performance testing. See “EEG, SvO₂ (SO₂) and tcGas Performance Tests” on page 91 for details.

When to Perform Tests

This table tells you when to perform specific tests. The corresponding test procedures are described in the following sections **All tests listed below must be performed on the monitor itself, any attached MMS/X2 and FMS incl. parameter modules.**

Table 2 When to perform tests

| Service Event (When performing... | Tests Required ...Complete these tests) |
|---|--|
| Installation | |
| Installation of a monitor in combination with a medical or non-medical device connected to the same multiple socket outlet. | Perform Visual Inspection, Power On and System Tests |
| Installation of a monitor (with or without iPC) with no display connected to the video output | Perform Visual Inspection and Power On Test |
| Installation of a monitor (with or without iPC) with a medical display specified by Philips | Perform Visual Inspection and Power On Test |
| Installation of a monitor without iPC with an off-the-shelf display (non-compliant with IEC 60601-1) | Perform Visual Inspection, Power On and System Test (per each affected video port) |
| Installation of a monitor with iPC with off-the-shelf displays (non-compliant with IEC 60601-1) connected to both the monitor and the iPC | Perform Visual Inspection, Power On and System Test for both the monitor-display connection and the iPC-display connection |
| Installation of a monitor with AGM or IntelliVue G1/G5, connected to separate mains sockets. | Perform Visual Inspection and Power On Tests |

3 Testing and Maintenance

| Service Event (When performing... | Tests Required ...Complete these tests) |
|--|--|
| Installation of a monitor with a Vuelink connection to another medical device (compliant with IEC 60601-1), connected to separate mains sockets. | Perform Visual Inspection and Power On Tests |
| Installation of a monitor with an IntelliBridge connection to another medical device (compliant with IEC 60601-1), connected to separate mains sockets. | Perform Visual Inspection and Power On Tests |
| Installation of a monitor (with or without iPC) with IT equipment e.g. printer, PC connected via a functional connection USB. | Perform Visual Inspection, Power On and System Tests |
| Installation of monitor with IntelliVue 802.11 Bedside Adapter | Perform Visual Inspection, Power On and IntelliVue 802.11 Bedside Adapter Communication Test |
| Installation of monitor with IntelliVue Instrument Telemetry | Perform Visual Inspection, Power On and IIT Communication Test |
| Installation of monitor with Short Range Radio (SRR) | Perform Visual Inspection, Power On and SRR Communication Test |
| Installation of networked monitor (LAN) without iPC | Perform Visual Inspection and Power On Test |
| Installation of networked monitor (LAN) with iPC (iPC LAN connected to non-Philips network) | Perform Visual Inspection, Power On Test and System Tests for the iPC LAN connected to non-Philips network |
| Preventive Maintenance | |
| Preventive Maintenance* | Perform preventive maintenance tests and procedures: NBP calibration Microstream CO2 calibration |
| Other Regular Tests and Tasks | |
| Visual Inspection | Perform Visual Inspection |
| Power On Test | Perform Power On test |
| Repairs | |
| Repairs where the monitor, FMS, parameter modules, MMS or X2 or iPC have been damaged by impact, liquid ingress, fire, short circuit or electrical surge. | Perform Visual Inspection, Power On, all Safety Tests and Full Performance Assurance Tests |
| Repairs where the power supply, the mains socket or an interface board of the monitor or iPC is removed or replaced or the protective earth ground connection is disrupted. | Perform Visual Inspection, Power On, all Safety Tests and Basic Performance Assurance Test |

| Service Event (When performing... | Tests Required ...Complete these tests) |
|---|--|
| Repairs of IntelliVue 802.11 Bedside Adapter | Perform Visual Inspection, Power On and IntelliVue 802.11 Bedside Adapter Communication Test |
| Repairs of IntelliVue Instrument Telemetry (IIT) Module | Perform Visual Inspection, Power On and IIT Communication Test |
| Repairs of Short Range Radio (SRR) Interface | Perform Visual Inspection, Power On and SRR Communication Test |
| Repairs of the parameter modules, FMS, MMS or X2 (all service events where the parameter modules, FMS, MMS or X2 have been opened) | Perform Visual Inspection, Power On, all Safety Tests and Basic Performance Assurance Test. If a certain parameter seems suspicious, perform Full Performance Assurance Test for this parameter. |
| Repairs where the NBP pump of the MMS or X2 has been replaced | Perform Visual Inspection, Power On, all Safety Tests, Basic Performance Assurance Test and NBP Performance Test and Calibration |
| Repairs where the parameter module, MMS or X2 has been replaced. | Perform Visual Inspection, Power On and Basic Performance Assurance |
| Repairs of the AGM or IntelliVue G1/G5 | Perform Basic Performance Assurance Test. For further testing requirements, see AGM or IntelliVue G1/G5 Service Guide |
| Repairs where the printer connected to the monitor or iPC via connector board has been replaced. | Perform Visual Inspection, Power On, System Test and Printer Test. |
| All other IntelliVue Monitoring System repairs (except when power supply is removed) | Perform Visual Inspection, Power On Test and Basic Performance Assurance Test |
| Performance Assurance | |
| Basic Performance Assurance | Perform basic performance assurance tests for the respective monitoring system component. |
| Full Performance Assurance | Perform all accuracy and performance test procedures listed in the following sections. If a particular measurement is in question, perform the measurement performance test only. |
| Upgrades | |
| Software Upgrades | Perform Visual Inspection, Power On Test and Basic Performance Assurance Test unless otherwise specified in the Upgrade Installation Notes shipped with the upgrade. |
| Hardware Upgrades | Perform Visual Inspection, Power On Test and Basic Performance Assurance Test unless otherwise specified in the Upgrade Installation Notes shipped with the upgrade. |

3 Testing and Maintenance

| Service Event (When performing... | Tests Required ...Complete these tests) |
|--|---|
| Hardware Upgrades where IntelliVue 802.11 Bedside Adapter is installed | Perform Visual Inspection, Power On Test, Basic Performance Assurance Test and IntelliVue 802.11 Bedside Adapter Communication Test |
| Hardware Upgrades where IntelliVue Instrument Telemetry (IIT) is installed | Perform Visual Inspection, Power On Test, Basic Performance Assurance Test and IIT Communication Test |
| Hardware Upgrades where Short Range Radio (SRR) is installed | Perform Visual Inspection, Power On Test, Basic Performance Assurance Test and SRR Communication Test |
| Installation of Interfaces or Hardware Upgrades where the power supply of the monitor or iPC or interface boards of the monitor need to be removed. | Perform Visual Inspection, Power On Test, Basic Performance Tests and all Safety Tests |
| Combining or Exchanging System Components (non-medical equipment connected to an IntelliVue monitor or medical system equipment operated on a multiple socket outlet) | Perform the System Test for the respective system components |


*M3015A with the old hardware Rev. A (i.e. Serial No. DE020xxxxx) also require the pump and scrubber replacement procedures.

NOTE

It is the responsibility of the facility operator or their designee to obtain reference values for recurring safety and system tests. These reference values are the results of the first test cycles after an installation. You may also purchase this service from Philips.

Testing Sequence

Summary of the recommended sequence of testing:



NOTE

If any single test fails, testing must be discontinued immediately and the device under test must be repaired or labeled as defective.

Visual Inspection

Before Each Use

Check all exterior housings for cracks and damage. Check the condition of all external cables, especially for splits or cracks and signs of twisting. If serious damage is evident, the cable should be replaced immediately. Check that all mountings are correctly installed and secure. Refer to the instructions that accompany the relevant mounting solution.

After Each Service, Maintenance or Repair Event

Ensure all fuses accessible from the outside comply with the manufacturer's specification.

Check:

- the integrity of mechanical parts, internally and externally.
- any damage or contamination, internally and externally

- that no loose parts or foreign bodies remain in the device after servicing or repair.
- the integrity of all relevant accessories.

Power On Test

- 1 Connect the monitoring system to mains and switch it on. This includes connected displays, MMS, MMS Extensions, X2, FMS and FMS associated modules, gas analyzers and Vuelink devices.
- 2 Make sure that all steps listed in the table *Initial Instrument Boot Phase* in the Troubleshooting section are completed successfully and that an ECG wave appears on the screen.

The expected test result is pass: the monitor boots up and displays an ECG wave. The wave might be a flat line if no simulator is attached.

Safety Tests

Safety tests are comprised of the following tests performed on the monitoring system:

- protective earth resistance
- equipment leakage current
- applied part leakage current
- system test (if applicable)

Safety test requirements are set according to international standards, their national deviations and specific local requirements. The safety tests detailed in this Service Guide are derived from international standards but may not be sufficient to meet local requirements. We recommend that you file the results of safety tests. This may help to identify a problem early particularly if the test results deteriorate over a period of time.

Each individual piece of equipment which has its own connection to mains or which can be connected or disconnected from mains without the use of a tool must be tested individually. The monitoring system as a whole must be tested according to the procedure described in “System Test” on page 57.

Accessories which can affect the safety of the equipment under test or the results of the safety test must be included in the tests and documented.

Warnings, Cautions, and Safety Precautions

- These tests are well established procedures of detecting abnormalities that, if undetected, could result in danger to either the patient or the operator.
- Disconnect the device under test from the patient before performing safety tests.
- Disconnect the device under test from mains before performing safety tests. If this is not possible, ensure that the performance of these tests does not result in danger to the safety analyzer operator, patients or other individuals.
- Test equipment (for example, a *Safety Analyzer*) is required to perform the safety tests. Please refer to Annex C of IEC/EN 62353 for exact requirements for the measurement equipment and for measurement circuits for protective earth resistance and leakage currents. Refer to the documentation that accompanies the test equipment. Only certified technicians should perform safety testing.
- The consistent use of a *Safety Analyzer* as a routine step in closing a repair or upgrade is emphasized as a mandatory step to maintain user and patient safety. You can also use the *Safety*

Analyzer as a troubleshooting tool to detect abnormalities of line voltage and grounding plus total current loads.

- During safety testing, mains voltage and electrical currents are applied to the device under test. Ensure that there are no open electrical conductive parts during the performance of these tests. Avoid that users, patients or other individuals come into contact with touch voltage.
- For Europe and Asia/Pacific, the monitor complies with:
IEC 60601-1:1988 + A1:1991 + A2:1995; EN60601-1:1990 + A1:1993 + A2:1995; IEC 60601-1-1:2001; EN 60601-1-1:2001; IEC 60601-1-2:2001; EN 60601-1-2:2001.
For USA, the monitor complies with:
UL60601-1
For Canada, CAN/CSA C22.2#601.1-M90
- Local regulations supersede the testing requirements listed in this chapter.
- If a non-medical electrical device is connected to a medical electrical device, the resulting medical electrical system must comply with IEC/EN 60601-1-1.
- Perform safety tests as described on the following pages.

Safety Test Procedures

Use the test procedures outlined here **only** for verifying and recording the initial values prior to or at installation, safe installation or service of the product, and for periodic recurrent testing. The setups used for these tests and the acceptable ranges of values are derived from local and international standards but may not be equivalent. These tests are not a substitute for local safety testing where it is required for an installation or a service event. If using an approved safety tester, perform the tests in accordance with the information provided by the manufacturer of the tester and in accordance with your local regulations, for example IEC/EN 60601-1, UL60601-1 (US), IEC/EN 62353, and IEC/EN 60601-1-1. The safety tester should print results as detailed in this chapter, together with other data.

Please refer to Annex C of IEC/EN 62353 for requirements for the measurement equipment and for measurement circuits for protective earth resistance and leakage currents.

3 Testing and Maintenance

The following symbols are used in the diagrams illustrating the safety tests:



Supply mains




Protective earth

L, N


Supply mains terminals

PE


Protective earth terminal




Mains part



Applied part



F-type applied part



Measuring device



Resistance measuring device



Connection to accessible conductive parts

.....

Optional connection

CAUTION

After each service, maintenance or repair event:

Ensure all fuses accessible from the outside comply with the manufacturer's specification.

Check:

- the integrity of mechanical parts, internally and externally.
 - any damage or contamination, internally and externally.
 - that no loose parts or foreign bodies remain in the device after servicing or repair.
 - the integrity of all relevant accessories.
-

Hints for Correct Performance of Safety Tests


- Perform a visual inspection on all detachable power cords used with the monitoring system and include these in all safety test procedures.

- Connection lines such as data lines or functional earth conductors may appear to act like protective earth connections. These may lead to incorrect measurements and need to be considered during testing. If necessary, unplug these connections.
- During measurements, the device under test shall be isolated from earth (e.g. test on an insulated work bench), except the protective earth conductor in the power supply cord.
- Position all cables and cords in such a manner that they do not influence the safety tests.
- Measurement of insulation resistance is not required.
- When testing a medical electrical system, where possible, test it such that potential ground voltage variations are present as they may be during actual use.

Guideline for Performance of Safety Tests

This section introduces the general principle of performing recurrent safety tests. Product specific test descriptions are described in the following sections.


Connect the detachable power cord of the device under test to the safety analyzer's test mains port. Connect the enclosure test lead of the safety analyzer to the enclosure of the device under test, e.g. to the equipotential connector or unearthed conductive accessible parts where applicable during Equipment Leakage Current Tests and Applied Part Leakage Current Tests. For testing the applied part leakage current, connect all applied parts to the safety analyzer using the appropriate patient lead or adapter cable. For the ECG parameter all ten ECG-leads need to be connected to the safety analyzer. If necessary, use an adapter cable to connect all ten ECG-leads. If necessary, repeat the safety test procedure until all available applied parts have been tested. Refer to the documentation that accompanies the safety analyzer for further details on how to set up and perform the test.



Protective Earth Resistance Test - Setup Example

NOTE


The test lead needs to go to parts that require protective earthing. This may be a single connection or several tested after each other



Equipment Leakage Current Test - Setup Example

NOTE

The test lead needs to go to the grounded enclosure parts, the ungrounded enclosure parts and all of the applied parts connected together.



Applied Part Current Test - Setup Example

NOTE

The above graphics resemble the Metron QA-90 setup and are protected by copyright. Copyright owned by Fluke (Metron).


Safety Test Adapter Cable - Schematics

The following graphics provide schematics of safety test (patient lead) adapter cables which can be used for electrical safety testing. These schematics can also be used as a guideline for making your own safety test adapter cables. Alternatively, other methods to make safety test adapter cables can be used, e.g. using a modified accessory cable.


NOTE

You may not need all of the cables displayed below for electrical safety testing of your respective monitor.


ECG




Wiring Schematic




SpO2 (MP2/X2, MP5, M3001A & M1020B #A01, #A02, #A03, #A04)




Wiring Schematic




SpO2 (M1020A)




Wiring Schematic



Invasive Pressure



Wiring Schematic



Invasive Pressure (M1006B #C01)

Wiring Schematic



Temperature

Wiring Schematic



Plug 2-Contact Male


Conductor Coloring

Plug Banana



CO2 (MP5, M3014A)

Wiring Schematic




PLUG 8-PIN MALE

Conductor Coloring

Plug Banana



CO2 (M1016A, M3016A)




Wiring Schematic

PLUG 12-PIN
MALE


Plug Banana

Conductor Coloring




4 = all resistors 120 KOhm

Cardiac Output




Wiring Schematic




BIS


Use Clamp Adapter Cable and BIS sensor simulator (P/N: M1034-61650, 453563233731).




Wiring Schematic



VueLink




Wiring Schematic



4 = 220 Ohm

IntelliBridge



Wiring Schematic


Plug RJ-45

Plug Banana

Conductor Coloring



EEG



Wiring Schematic


Plug 12-Contact

Plug Banana


Male Conductor Coloring





SvO2 (M1021A)




Wiring Schematic




ScVO2 (M1011A)





tcpO2/tcpCO2





1
Wiring Schematic





MP5 Predictive Temperature



MP5 TAAP




TcG10



Electrical Safety Testing

S(1): Protective Earth Resistance Test

Test to perform:



Measuring circuit for the measurement of Protective Earth Resistance in medical electrical equipment that is disconnected from the supply mains.

This measures the impedance of the Protective Earth (PE) terminal to all exposed metal parts of the Device under Test (DUT), which are for safety reasons connected to the Protective Earth (PE).

You can find metal parts of the device at the equipotential connector.

Measurements shall be performed using a measuring device capable to deliver a current of at least 200 mA into 500 mOhms with maximum open circuit voltage of 24V

This safety test is based on IEC/EN 62353.

Report the highest value (X1).


| Test | Expected test results |
|---|-----------------------|
| Protective Earth Resistance Test (with mains cable) | X1 \leq 300mOhms |

NOTE

- If the protective earth resistance test fails, testing must be discontinued immediately and the device under test must be repaired or labeled as defective.
- All values for current and voltage are the root mean square (r.m.s.) values, unless otherwise stated.
- Flex the power cord during the protective earth resistance test to evaluate its integrity. If it does not pass the test, exchange the power cord. Then repeat the test. If it still does not pass, follow the instructions in the first bullet point of this note above.

S(2): Equipment Leakage Current Test - Normal Condition

Test to perform:



Measuring circuit for the measurement of Equipment Leakage Current - Direct method according to IEC/EN 62353.

This test measures leakage current of accessible conductive and non-conductive metal parts of the monitor and the functional earth leakage current. It tests normal and reversed polarity. Perform the test with S1 closed (Normal Condition).

There are no parts of the equipment that are not protectively earthed. Disconnect any data cables and any connections that may provide an extraneous earth path. Test the device under test (DUT) on an insulated surface. Do not touch the DUT during testing.

This safety test is based on IEC/EN 62353.

3 Testing and Maintenance

For measurement limits, refer to Safety (2) test, Test and Inspection Matrix.

Report the highest value (X1).

| Test | Expected test results |
|---|--------------------------|
| Equipment Leakage Current Test (Normal Condition - with mains cable) | $X1 \leq 100\mu\text{A}$ |


NOTE

All values for current and voltage are the root mean square (r.m.s.) values, unless otherwise stated.

In case of an IT-power system, this safety test measurement requires a special measuring circuit, for example with its own integrated TN-system or use of an external isolation transformer attached to the safety test device.

S(3): Equipment Leakage Current Test - Single Fault Condition

Test to perform:



Measuring circuit for the measurement of Equipment Leakage Current - Direct method according to IEC/EN 62353.

This test measures leakage current of accessible conductive and non-conductive metal parts of the monitor and the functional earth leakage current. It tests normal and reversed polarity. Perform the test with S1 open (Single Fault Condition).

There are no parts of the equipment that are not protectively earthed. Disconnect any data cables and any connections that may provide an extraneous earth path. Test the device under test (DUT) on an insulated surface. Do not touch the DUT during testing.

This safety test is based on IEC/EN 62353.

For measurement limits, refer to Safety (3) test, Test and Inspection Matrix.

Report the highest value (X2).

| Test | Expected test results |
|--|--------------------------|
| Equipment Leakage Current Test (Single Fault Condition - with mains cable) | $X2 \leq 300\mu\text{A}$ |

NOTE

All values for current and voltage are the root mean square (r.m.s.) values, unless otherwise stated.


In case of an IT-power system, this safety test measurement requires a special measuring circuit, for example with its own integrated TN-system or use of an external isolation transformer attached to the safety test device.

S(4): Applied Part Leakage Current - Mains on Applied Part**NOTE**

During measurement of the Applied Part Leakage Current it is possible that the measured current can exceed the allowed limit (per IEC/EN 60601-1 or IEC/EN 62353).

This can occur when the safety tester is connected to the invasive blood pressure and temperature connectors at the same time during the applied leakage current measurement.


The connectors for the invasive blood pressure and temperature are independently functioning connectors.



Although there are individual connectors on the front end, internally those parameters use the same electrical insulation interface and are hardwired to each other. This results in an electrical short of those connectors during measurement if a test current is applied simultaneously. Therefore this should be avoided.

Due to the combined insulation interface, it is sufficient to connect to only one parameter interface (that is, Invasive Blood Pressure or Temperature) of the invasive blood pressure/temperature measurement block. This avoids a short and the potential of exceeding the limit for the current.

Test to perform:



Measuring circuit for the measurement of Applied Part Leakage Current - Direct method according to IEC/EN 62353.

This test measures applied part leakage current from applied part to earth caused by external main voltage on the applied part. Each polarity combination possible shall be tested. This test is applicable to each Applied Part tested and results recorded in turn with all other Applied Parts left floating. Applied Parts with multiple connections (e.g. ECG) are tested with the connections short-circuited.

There are no parts of the equipment that are not protectively earthed.

This safety test is based on IEC/EN 62353.

For measurement limits and test voltage, refer to Safety (4) test, Test and Inspection Matrix.

Report the highest value. (X1).

| Test | Expected test results |
|--|-----------------------|
| Applied Part Leakage Current Test (Single Fault Condition - mains on applied part) | S4 ≤ 50µA (CF) |

NOTE

All values for current and voltage are the root mean square (r.m.s.) values, unless otherwise stated.

In case of an IT-power system, this safety test measurement requires a special measuring circuit, for example with its own integrated TN-system or use of an external isolation transformer attached to the safety test device.

Reference: Allowable Values for IEC 60601-1:1998 and UL 60601-1 Measurements

Protective Earth resistance (between the PROTECTIVE EARTH TERMINAL and any ACCESSIBLE METAL PART which is PROTECTIVELY EARTHED, w/o power cord): 100mOhms

Protective Earth resistance of power cord: 100mOhms

Enclosure leakage current (IEC 60601-1 and UL60601-1): 100 µA (N.C.)

Enclosure leakage current:(IEC 60601-1): 500 µA (S.F.C)

Enclosure leakage current (UL 60601-1): 300 µA (S.F.C)

Patient leakage current: (IEC 60601-1 and UL60601-1): 100 µA (N.C.) for BF

Patient leakage current: (IEC 60601-1 and UL60601-1): 500 µA (S.F.C.) for BF

Patient leakage current: (IEC 60601-1 and UL60601-1): 10 µA (N.C.) for CF

Patient leakage current: (IEC 60601-1 and UL60601-1): 50 µA (S.F.C.) for CF

All values for current and voltage are the root mean square (r.m.s.) values, unless otherwise stated

Insulation Resistance

It is not recommended to perform measurements of the insulation resistance. Refer to IEC 62353 for details about methods of the insulation resistance measurement.

System Test

After mounting and setting up a system, perform system safety tests according to IEC/EN 60601-1-1.

What is a Medical Electrical System?

A medical electrical system is a combination of at least one medical electrical piece of equipment and other electrical equipment, interconnected by functional connection or use of a multiple portable socket-outlet.

- Devices forming a medical electrical system must comply either with IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16.
- Any electrical device such as IT equipment that is connected to the medical electrical equipment must comply either with IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16 and be tested accordingly.
- Non-medical electrical equipment may require connection through a separating device (e.g. an isolation transformer).

General Requirements for a System


After installation or subsequent modification, a system must comply with the requirements of the system standard IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16. Compliance is checked by inspection, testing or analysis, as specified in the IEC/EN 60601-1-1 or in this book.

Medical electrical equipment must comply with the requirements of the general standard IEC/EN 60601-1, its relevant particular standards and specific national deviations. Non-medical electrical equipment shall comply with IEC safety standards that are relevant to that equipment.

Relevant standards for some non-medical electrical equipment may have limits for equipment leakage currents higher than required by the standard IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16. These higher limits are acceptable only outside the patient environment. It is essential to reduce equipment leakage currents to values specified in IEC/EN 60601-1 when non-medical electrical equipment is to be used within the patient environment.

System Example

This illustration shows a system where both the medical electrical equipment and the non-medical electrical equipment are situated at the patient's bedside.



WARNING

- Do not use additional AC mains extension cords or multiple portable socket-outlets. If a multiple portable socket-outlet is used, the resulting system must be compliant with IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16. Do not place multiple socket-outlets on the floor. Do not exceed the maximum permitted load for multiple socket-outlets used with the system. Do not plug additional multiple socket outlets or extension cords into multiple socket outlets or extension cords used within the medical electrical system.
 - Do not connect any devices that are not supported as part of a system.
 - Do not use a device in the patient vicinity if it does not comply with IEC/EN 60601-1 or IEC 60601-1 edition 3 clause 16. The whole installation, including devices outside of the patient vicinity, must comply with IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16. Any non-medical device placed and operated in the patient's vicinity must be powered via a separating transformer (compliant with IEC/EN 60601-1-1 or IEC 60601-1 edition 3 clause 16) that ensures mechanical fixing of the power cords and covering of any unused power outlets.
-


System Installation Requirements

- Ensure that the medical electrical system is installed in a way that the user achieves optimal use.
- Make sure the user is informed about the required cleaning, adjustment, sterilization and disinfection procedures listed in the Instructions for Use.
- The medical electrical system must be installed in such a way that the user is able to carry out the necessary cleaning, adjustment, sterilization and disinfection procedures listed in the Instructions for Use.

- Ensure that the medical electrical system is installed in a way that an interruption and restoration of power to any part of the medical electrical system does not result in a safety hazard.
- We recommend using fixed mains socket outlets to power the medical system or parts thereof. Avoid using multiple portable socket-outlets.
- Any multiple portable socket outlets used must be compliant with IEC 60884-1 and IEC 60601-1-1 or IEC 60601-1 edition 3 clause 16.
- Ensure that any part of the system connected to multiple portable socket-outlets is only removable with a tool, i.e. the multiple portable socket-outlet provides a locking mechanism to prevent power cords from being plugged or unplugged unintentionally. Otherwise, the multiple portable socket-outlet must be connected to a separation device. Multiple Socket Outlets used within the medical electrical system must only be used for powering medical electrical equipment which is part of the system.
- Ensure that any functional connections between parts of the medical electrical system are isolated by a separation device according to IEC 60601-1-1 or IEC 60601-1 edition 3 clause 16 to limit increased equipment leakage currents caused by current flow through the signal connections where necessary (e.g. leakage current coming from non-medical electrical equipment into medical electrical equipment or building ground voltage differences providing leakage current through grounded data cables). This only works if the equipment leakage current of the respective medical electrical system parts is not exceeded under normal conditions. This isolation is especially important where the non-medical electrical equipment leakage currents can pass to the medical electrical equipment in the system or building ground voltage differences can pass to the medical electrical equipment via ground in a data cable connection in the system
- Avoid increase of equipment leakage currents when non-medical electrical equipment within the medical electrical system is used. This only applies when if the equipment leakage current of the respective medical electrical system parts is not exceeded under normal conditions. Use of an additional protective earth connection, separation device or additional non-conductive enclosures are options that can prevent a problem.
- Within the patient environment it is important to limit electrical potential differences between different parts of a system. If necessary, use potential equalization equipment (equipotential cable) or additional protective earth connections.
- Medical electrical equipment used in medical rooms must be connected to potential equalization equipment (equipotential cable) to avoid electrical potential differences. Check your local requirements for details.

Required Protective Measures at System Installation

For any IT equipment (IEC60950-1) operated in the patient environment ensure that the equipment leakage current does not exceed the limits described in IEC 60601-1. Use a separation device to ensure compliance. After installation of IT equipment in the patient environment, an equipment leakage current test is required.




Case 1: Medical Device Combined with Medical Device


If you combine a medical device with another medical device (incl. Philips specified displays) to form a medical electrical system according to IEC60601-1-1 or IEC/EN 60601-1 edition 3 clause 16, no additional protective measures are required. The medical electrical devices may be located in or outside the patient vicinity in a medically used room. This is valid as long as the medical devices are connected to separate mains outlets. No system test is required.

NOTE


The pictures below and in the following chapters show the MX800 as an example. All cases apply to the MX600/700 as well.



3 Testing and Maintenance




If the combined medical devices are connected to the same multiple portable socket outlet an enclosure leakage current test of the entire device combination on the multiple portable socket outlet is required to ensure that the resulting protective earth leakage current and equipment leakage current does not exceed the limits of IEC 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16. Avoid using multiple portable socket outlets. The medical electrical devices may be located in or outside the patient vicinity in a medically used room. If the limits are exceeded, additional protective measures are required, e.g. a separation device or the connection of each device to separate mains.





Case 2: Medical Device Combined with a Non-Medical Device


If you combine a medical device with a non-medical device to form a medical electrical system according to IEC60601-1-1 or IEC/EN 60601-1 edition 3 clause 16, additional protective measures are required, e.g. usage of a separation device. The medical electrical devices or the IT equipment may be located in or outside the patient vicinity in a medically used room. After system installation incl. protective measures, a system test is required to ensure that the resulting equipment leakage current and applied part leakage current does not exceed the limits of IEC 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16.




For any IT equipment (IEC60950-1) operated in patient vicinity ensure that the equipment leakage current does not exceed the limits described in IEC 60601-1 or IEC/EN 60601-1 edition 3 clause 16. Use a separation device to ensure compliance. After installation of IT equipment in patient vicinity, an equipment leakage current test is required.



If the combined devices forming the medical electrical system are connected to the same multiple portable socket outlet, ensure that the resulting protective earth leakage current **and** equipment leakage current do not exceed the limits of IEC 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16. The medical electrical devices or IT equipment may be located in or outside the patient vicinity in a medically used room. Avoid using multiple portable socket outlets. If the limits of IEC 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16 are exceeded, additional protective measures are required, e.g. a separation device or the connection of each device to separate mains.





For any IT equipment (IEC60950-1) operated in patient vicinity ensure that the equipment leakage current does not exceed the limits described in IEC 60601-1 or IEC/EN 60601-1 edition 3 clause 16. Use a separation device to ensure compliance. After installation of IT equipment in patient vicinity, an equipment leakage current test is required.




Case 3: Medical Device Combined with a Medical or Non-Medical Device with one Device in a Non-Medically-Used Room


If you combine a medical device with a medical or non-medical device to form a medical electrical system according to IEC60601-1-1 or IEC/EN 60601-1 edition 3 clause 16 using a common protective earth connection and one of the devices is located in a non-medically used room, additional protective measures are required, e.g. usage of a separation device or additional protective earth connection. The medical electrical devices or IT equipment may be located in or outside the patient vicinity. After system installation incl. protective measures, a system test is required to ensure that the resulting equipment leakage current does not exceed the limits of IEC 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16.





If you combine a medical device with a medical or non-medical device to form a medical electrical system according to IEC60601-1-1 or IEC/EN 60601-1 edition 3 clause 16 using two separate protective earth connections and one of the devices is located in a non-medically used room creating a potential voltage difference, additional protective measures are required, e.g. usage of a separation device or additional protective earth connection. The medical electrical devices or IT equipment may be located in or outside the patient vicinity. After system installation incl. protective measures, a system test is required to ensure that the resulting equipment leakage current does not exceed the limits of IEC 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16.





System Test Procedure

If the medical electrical device has already been tested as a standalone device e.g. during factory safety testing, an equipment leakage current test must only be performed once the device is connected to another electrical device/system. If the medical electrical system has not been tested as a standalone device, the device has to be tested as a standalone device (without connection to the system) and as part of the system (with connection to the system).

Connect the detachable power cord of the device under test to the safety analyzer's test mains port. Connect the enclosure test lead of the safety analyzer to the enclosure of the device under test as described in the "Equipment Leakage Test" section. Refer to the documentation that accompanies the safety analyzer for further details on how to set up the test.

| Test | Expected test results |
|---|-------------------------|
| Equipment Leakage Current Test (Normal Condition) | Sys1 \leq 100 μ A |
| Equipment Leakage Current Test (Single Fault Condition) | Sys2 \leq 300 μ A |

After the testing of the device as a standalone device and as part of the system, check that the resulting values (without connection and with connection to the system) do not differ by more than +/- 10% from each other.

If the devices in the medical electrical system are connected to a multiple portable socket outlet the resulting protective earth leakage current needs to be determined. All system components must be connected to the multiple portable socket outlet and be switched on during this measurement.

| Test | Expected test results |
|---|-------------------------|
| Protective Earth Leakage Current of Multiple Socket Outlets | Sys3 \leq 300 μ A |

Refer to the documentation that accompanies the safety analyzer for further details on how to set up the test.

Preventive Maintenance Procedures

Noninvasive Blood Pressure Measurement Calibration

Carry out the noninvasive blood pressure measurement performance tests at least every two years , or as specified by local laws (whichever comes first).

Microstream CO₂ Calibration

Carry out the Microstream CO₂ calibration once a year or after 4000 hours of continuous use and following any instrument repairs or the replacement of any instrument parts.

Performance Assurance Tests

Some of the following test procedures must be performed in service mode. To enter service mode select **Operating Modes** in the main menu. Then select **Service Mode** and enter the password.

If required, open the screen menu in the monitor info line at the top of the screen and select **Service** to access the service screen. This is required particularly for Anesthetic Gas Module testing procedures.

Basic Performance Assurance Test

This section describes the basic performance test procedure. Please refer to the section for detailed information on when which test procedure is required.

Procedure:

Power on the monitoring system and go into demo mode. Check that each connected parameter (module, MMS, Gas Analyzer, Vuelink connected device) displays values.

Full Performance Assurance Test

The following sections describe the full performance testing procedures i.e. detailed testing of each parameter with a patient simulator or specified tools. Please refer to the section for information on when which testing procedure is required.

ECG/Resp Performance Test

This test checks the performance of the ECG and respiration measurements.

Tools required: Patient simulator.

ECG Performance

- 1 Connect the patient simulator to the ECG/Resp connector on the MMS/IntelliVue X2.
- 2 Configure the patient simulator as follows:
 - ECG sinus rhythm.
 - HR = 100 bpm or 120 bpm (depending on your patient simulator).
- 3 Check the displayed ECG wave and HR value against the simulator configuration.
- 4 The value should be 100bpm or 120 bpm +/- 2 bpm.

Respiration Performance

- 1 Change the Patient Simulator configuration to:
 - Base impedance line 1500 Ohm.
 - Delta impedance 0.5 Ohm.
 - Respiration rate 40 rpm or 45 rpm.
- 2 The value should be 40 rpm +/- 2 rpm or 45 rpm +/- 2 rpm.

| Test | Expected test results |
|------------------------------|---|
| ECG Performance Test | 100 bpm +/- 2 bpm or 120 bpm +/- 2 bpm |
| Respiration Performance Test | 40 rpm +/- 2 rpm or 45 rpm +/- 2 rpm |

ECG Out Performance Test (not available via SRR)

This test checks the performance of ECG synchronization between the monitor and a defibrillator. It only needs to be performed when this feature is in use as a protocol at the customer site.

Tools required:

- Defibrillator with ECG Input.
 - Patient simulator.
- 1 Connect the patient simulator to the ECG connector of the MMS and the defibrillator to the ECG Output on the monitor with the ECG Sync cable.
 - 2 Set the patient simulator to the following configuration:
 - HR = 100 bpm or 120 bpm (depending on your patient simulator).
 - ECG sinus rhythm.
 - 3 Switch the defibrillator to simulation mode.
 - 4 Check that the ECG signal is displayed.

| Test | Expected test results |
|--------------------------|-------------------------------------|
| ECG Out Performance Test | ECG signal is displayed (pass/fail) |

SpO2 Performance Test

This test checks the performance of the SpO₂ measurement.

Tools required: none

- 1 Connect an adult SpO₂ transducer to the SpO₂ cradle.
- 2 Measure the SpO₂ value on your finger (this assumes that you are healthy).
- 3 The value should be between 95% and 100%.

| Test | Expected test results |
|-----------------------------------|-----------------------|
| SpO ₂ Performance Test | 95% and 100% |

Measurement Validation

The SpO₂ accuracy has been validated in human studies against arterial blood sample reference measured with a CO-oximeter. In a controlled desaturation study, healthy adult volunteers with

saturation levels between 70% and 100% SaO₂ were studied. The population characteristics for those studies were:

- about 50% female and 50% male subjects
- age range: 19 to 39
- skin tone: from light to dark brown

NOTE

A functional tester cannot be used to assess the accuracy of a pulse oximeter monitor. However, it can be used to demonstrate that a particular pulse oximeter monitor reproduces a calibration curve that has been independently demonstrated to fulfill a particular accuracy specification.


NBP Performance Test

This section describes NBP test procedures. The monitor must be in service mode and the screen “Service A” must be selected to perform these tests. The NBP Performance Test consists of:

- NBP Accuracy Test
- NBP Leakage Test
- NBP Linearity Test
- Valve Test

NBP Accuracy Test

This test checks the performance of the non-invasive blood pressure measurement. Connect the equipment as shown:



Tools required:

- Reference manometer (includes hand pump and valve), accuracy 0.2% of reading.
- Expansion chamber (volume 250 ml +/- 10%)
- Appropriate tubing.

In service mode, the systolic and diastolic readings indicate the noise of NBP channels 1 and 2 respectively. When static pressure is applied, the reading in NBP channel 1 should be below 50. The value in parentheses indicates the actual pressure applied to the system.

- 1 Connect the manometer and the pump with tubing to the NBP connector on the MMS and to the expansion chamber.
- 2 In service mode, select the **Setup NBP** menu.
- 3 Select **Close Valves: On**

3 Testing and Maintenance

- 4 Raise the pressure to 280 mmHg with the manometer pump.
- 5 Wait 10 seconds for the measurement to stabilize.
- 6 Compare the manometer values with the displayed values.
- 7 Document the value displayed by the monitor (**x1**).
- 8 If the difference between the manometer and displayed values is greater than 3 mmHg, calibrate the MMS. If not, proceed to the leakage test.
- 9 To calibrate the MMS, select **Close Valves off** then **Calibrate NBP** and wait for the instrument to pump up the expansion chamber. Wait a few seconds after pumping stops until **EnterPrVal** is highlighted and then move the cursor to the value shown on the manometer. If one of the following prompt messages appears during this step, check whether there is leakage in the setup:
 - NBP unable to calibrate—cannot adjust pressure
 - NBP unable to calibrate—unstable signal
- 10 Press **Confirm**.

If the INOP NBP Equipment Malfunction message occurs in monitoring mode, go back to service mode and repeat the calibration procedure.

NBP Leakage Test

The NBP leakage test checks the integrity of the system and of the valve. It is required once every two years and when you repair the MMS or X2 or replace parts.

- 1 If you have calibrated, repeat steps 2 to 6 from the accuracy test procedure so that you have 280 mmHg pressure on the expansion chamber.
- 2 Watch the pressure value for 60 seconds.
- 3 Calculate and document the leakage test value (**x2**).
 $x_2 = P_1 - P_2$
where P1 is the pressure at the beginning of the leakage test and P2 is the pressure displayed after 60 seconds.
The leakage test value should be less than 6 mmHg.

NBP Linearity Test

- 1 Reduce the manometer pressure to 150 mmHg.
- 2 Wait 10 seconds for the measurement to stabilize.
- 3 After these 10 seconds, compare the manometer value with the displayed value.
- 4 Document the value displayed by the monitor (**x3**).
- 5 If the difference is greater than 3 mmHg, calibrate the MMS or X2 (see steps 9 to 10 in the accuracy test procedure).

Valve Test

- 1 Raise the pressure again to 280 mmHg.
- 2 Select **Close valves: Off**.
- 3 Wait five seconds and then document the value displayed. The value should be less than 10 mmHg.

- 4 Document the value displayed by the monitor (x4).

| Test | Expected test results |
|----------------|---|
| Accuracy test | x1 = 280 ± 3 mmHg Difference ≤ 3 mmHg |
| Leakage test | x2 = leakage test value x2 < 6 mmHg |
| Linearity test | x3 = 150 ± 3 mmHg Difference ≤ 3 mmHg |
| Valve Test | x4 = value < 10 mmHg |

Invasive Pressure Performance Test

This test checks the performance of the invasive pressure measurement.

Tools required: Patient simulator.

- 1 Connect the patient simulator to the pressure connector.
- 2 Set the patient simulator to 0 pressure.
- 3 Make a zero calibration.
- 4 Configure the patient simulator as P(static) = 200 mmHg.
- 5 Wait for the display.
- 6 The value should be $200 \text{ mmHg} \pm 5 \text{ mmHg}$. If the value is outside these tolerances, calibrate the Invasive Pressure measurement. If the measurement was calibrated with a dedicated reusable catheter, check the calibration together with this catheter.

| Test | Expected test results |
|------------------------------------|---------------------------------------|
| Invasive Pressure Performance Test | $200 \text{ mmHg} \pm 5 \text{ mmHg}$ |

Temperature Performance Test

This test checks the performance of the temperature measurement.

Tools required: Patient simulator (with 0.1°C or 0.2°F).

- 1 Connect the patient simulator to the temperature connector.
- 2 Configure the patient simulator to 40°C or 100°F .
- 3 The value should be $40^\circ\text{C} \pm 0.2^\circ\text{C}$ or $100^\circ\text{F} \pm 0.4^\circ\text{F}$.

| Test | Expected test results |
|------------------------------|---|
| Temperature Performance Test | $40^\circ\text{C} \pm 0.2^\circ\text{C}$ or $100^\circ\text{F} \pm 0.4^\circ\text{F}$ |

M3014A Capnography Extension Performance Tests

The procedures below describe the mainstream and sidestream CO₂ performance tests for the M3014A Capnography Extension.

Mainstream CO₂ Accuracy Check

Tools Required:

- three airway adapters
- Verification Gas M2506A
- Gas cylinder regulator M2505A

You also need a local barometric pressure rating received from a reliable local source (airport, regional weather station or hospital weather station) which is located at the same altitude as the hospital.

Procedure:

- 1 Attach the M2501A CO₂ sensor to the patient monitor. Attach an airway adapter to the sensor. Make sure that the sensor is disconnected from the patient circuit.
- 2 Switch on the patient monitor.
- 3 Enter the monitor's Service Mode.
- 4 Using the sensor status provided in the M2501A Serial protocol, wait for the M2501A sensor to warm up to its operating temperature.
- 5 The default setting for gas temperature is 22°C. If the gas temperature is significantly above or below this value, correct the gas temperature setting.
- 6 Zero the sensor on the airway adapter being used in this test. Ensure Zero Gas is set to Room Air
- 7 Attach a regulated flowing gas mixture of 5% CO₂, balance N₂ to the airway adapter.
- 8 Set the gas correction to off.
- 9 Allow a few seconds for the gas mixture to stabilize and observe the CO₂ value. The expected value is 5% of the ambient pressure ± 2 mmHg

NOTE

Make sure that you follow the above steps correctly. If the sensor fails this check it must be exchanged. The sensor cannot be calibrated.

Example for an expected test result:

The expected test result for an altitude of 0 m (sea level) at approximately 760 mmHg ambient pressure is:

| Test | Expected test results (x1) | Acceptance Range |
|--|--------------------------------------|-------------------|
| Mainstream CO ₂ Accuracy Test | 5% of 760 mmHg pressure ± 2 mmHg | 36 mmHg - 40 mmHg |

NOTE

The expected test results will differ depending on the conditions (i.e. altitude or ambient pressure).

Sidestream CO₂ Accuracy Check

Tools Required:

- Cal gas flow regulator M2267A
- Cal tube 13907A
- Verification Gas M2506A
- Straight Sample Line M2776A

You also need a local barometric pressure rating received from a reliable local source (airport, regional weather station or hospital weather station) which is located at the same altitude as the hospital.

Procedure:

- 1 Attach the M2741A CO₂ sensor to the patient monitor. Attach the sample line and the cal tube to the sensor. Make sure that the sensor is disconnected from the patient circuit.
- 2 Switch on the patient monitor.
- 3 Enter the monitor's Service Mode.
- 4 Using the sensor status provided in the M2741A Serial protocol, wait for the M2741A sensor to warm up to its operating temperature.
- 5 Zero the sensor. Ensure Zero Gas is set to Room Air
- 6 Attach a regulated flowing gas mixture of 5% CO₂, balance N₂ to the cal tube.
- 7 Set the gas correction to off.
- 8 Allow a few seconds for the gas mixture to stabilize and observe the CO₂ value. The expected value is 5% of the ambient pressure ± 2 mmHg

NOTE

Make sure that you follow the above steps correctly. If the sensor fails this check it must be exchanged. The sensor cannot be calibrated

Example for an expected test result:

The expected test result for an altitude of 0 m (sea level) at approximately 760 mmHg ambient pressure is:

| Test | Expected test results (x2) | Acceptance Range |
|--|--------------------------------------|-------------------|
| Sidestream CO ₂ Accuracy Test | 5% of 760 mmHg pressure ± 2 mmHg | 36 mmHg - 40 mmHg |

NOTE

The expected test results will differ depending on the conditions (i.e. altitude or ambient pressure).

Sidestream CO₂ Flow Check

Check the flow rate in the Sidestream CO₂ extension as follows:

- 1 Connect the flowmeter to the sample line
- 2 Check on the flowmeter the flow that the Sidestream CO₂ extension pump draws. It should be 50 ml/min ± 10 ml/min. If the value is not within tolerance check your setup again and perform another flow check. If it fails again, the sensor must be replaced. The sensor cannot be calibrated.

Example for an expected test result:

The expected test result for an altitude of 0 m (sea level) at approximately 760 mmHg ambient pressure is:

| Test | Expected test results (x3) | Acceptance Range |
|---------------------------|----------------------------|-----------------------|
| Sidestream CO2 Flow Check | 50 ml/min \pm 10 ml/min | 40 ml/min - 60 ml/min |

NOTE

The expected test results will differ depending on the conditions (i.e. altitude or ambient pressure).

Microstream CO2 Performance Test

Allow five seconds between individual service procedures to ensure stable equipment conditions. When certain monitor procedures are running, service procedures are not possible and trying to start them will result in a message **Service Operation Failed** in the monitor's status line. Wait until the monitor completes the current operation, then restart the service procedure.

This test checks the performance of the Microstream CO2 measurement. The Microstream CO2 measurement can either be integrated into the IntelliVue MP5 monitor or, for other IntelliVue monitors, into the M3015A/B MMS Extensions. The Microstream CO2 performance test is required once per year or after 4000 hours of continuous use and when the instrument is repaired or when parts are replaced.

This test uses calibration equipment that you can order (see the *Parts* section for the part number). The procedure is summarized in the following steps. Refer to the documentation accompanying the equipment for detailed instructions.

Tools Required:

- Standard tools, such as screwdriver, tweezers
- Electronic flowmeter, M1026-60144 or Mass Flowmeter 453564178121
- Digital Barometer \pm 2mbar or better
- Gas calibration equipment:
- Cal 1 gas 15210-64010 (5% CO₂)
- Cal 2 gas 15210-64020 (10% CO₂)
- Cal gas flow regulator M2267A
- Cal tube 13907A
- Calibration Line M3015-47301
- Leakage Test Kit M1013-64002 (451261014851) (only required for leakage test without M1026-60144 Flowmeter)
- Flexible Connecting Tube

You also need a local barometric pressure rating received from a reliable local source (airport, regional weather station or hospital weather station) which is located at the same altitude as the hospital.

The CO2 calibration for the Microstream extension consists of the following steps:

- Leakage check, **either** with M1026-60144 Flowmeter **or** with 453564178121 Mass Flowmeter*
- Barometric pressure check and calibration, if required.*
- Pump check*

- Flow check and calibration, if required
- Noise check
- CO₂ Cal check and calibration, if required
- CO₂ Cal verification

Perform all checks in the same session.

* Not applicable for all HW Revisions. See individual test sections for details.

NOTE

The M3015A/B HW Rev C is indicated as HW Rev. Q.xx.xx in the IntelliVue Revision Screen.

Leakage Check with M1026-60144 Flowmeter (only for M3015A with HW Rev. A and B and Firmware Revision < P.01.32)

The leakage check consists of checking the tubing between:

- the pump outlet and the mCO₂ outlet and
- the pump inlet and calibration line inlet.

Check the user's guide of the flowmeter for details on how to make a correct flow reading.

Part 1

- 1 Go into service mode and select **Setup CO₂** menu.
- 2 Connect a calibration line to the Microstream CO₂ input to start the pump running.
- 3 Check the ambient pressure and the cell pressure shown in the monitor's status line. The cell pressure should be approximately 20 mmHg lower than ambient pressure. (This test is only to check that the pump starts and is running, which is also indicated by the noise generated by the running pump.)
- 4 Connect the flowmeter outlet to the calibration line inlet using a flexible connecting tube.
- 5 Block the mCO₂ outlet using your fingertip and observe the flowmeter display. The value on the flowmeter (**x1**) should decrease to between 0 and 4 ml/min, accompanied by an audible increase in pump noise. If the value is within the tolerance limits, continue with part 2 of the leakage check.
- 6 If the value is outside the tolerance limits, there is a leakage between the pump outlet and the mCO₂ outlet.
- 7 Open the MMS Extension or MP5 and check the tubing connections at the pump outlet and the extension gas outlet. If the connections are good, then there is a leakage in the tubing and you must exchange the MMS Extension or the mCO₂ Assembly of the MP5 respectively.

Part 2

- 1 Disconnect the flowmeter from the Part 1 setup and connect the flowmeter inlet to the M3015A gas outlet or the MP5 mCO₂ gas outlet.
- 2 Leave the calibration line connected to the M3015A inlet or the MP5 mCO₂ inlet..
- 3 Block the inlet of the calibration line using your fingertip and observe the flowmeter display. The value on the flowmeter (**x2**) should decrease to between 0 and 4 ml/min, accompanied by an audible increase in pump noise. The cell pressure shown in the status line on the display should decrease to between 300 and 500 mmHg. Do not block the inlet for longer than 25 seconds as this will lead to an "Occlusion" INOP. If the value is within the tolerance limits, there are no leakages and the leakage check is completed; proceed to the pump check.

3 Testing and Maintenance

- 4 If the value is not within the tolerance limits, there is a leakage between the calibration line inlet and the pump inlet.
- 5 Check the calibration line connections and open the M3015A or MP5 to check the tubing connections at the pump inlet and the M3015A or MP5 mCO₂ gas inlet. If the connections are good, try replacing the calibration line and repeating the leakage check. If the situation remains, there is a leakage in the tubing and the M3015A or the mCO₂ assembly of the MP5 must be exchanged.

| Test | Expected test results |
|-----------------------------|--|
| Leakage Check Parts 1 and 2 | x1 = value of part 1 leakage check on flowmeter (x1 < 4.0 ml/min) x2 = value of part 2 leakage check on flowmeter (x2 < 4.0 ml/min) |


Leakage Check without for M3015B and M3015A with HW Rev C or M3015A with HW Rev. A/B without M1026-60144 Flowmeter

Preparation of Leakage Test Kit:

Remove two Luer connectors from the Leakage Test Kit, as shown in the following picture.

NOTE

These Luer connectors are not required for the actual Leakage Check. However, you should keep them, as they are required for other tests (e.g. for the kit leak test as documented later in this section).



Test Setup:

- 1 Connect the Calibration Line (M3015-47301) to the inlet of the M8105A/M3015A/B (the M8105A/M3015A/B must be switched off, either by disconnecting from the host monitor or by switching off the monitor).

- 2 Connect the leakage test tubing to the outlet of the M8105A/M3015A/B, to the digital barometer, to the calibration line, and the (empty) syringe as shown below. Make sure all connections have a tight fit!




Test Procedure:

- 1 Open the 3-way stopcock for all three limbs.



- 2 Switch on the digital barometer (the digital barometer should now display the actual ambient pressure).
- 3 Now slowly draw at the syringe, as if filling the syringe, until the pressure (as displayed on the digital barometer) drops to approximately 350 mbar below ambient pressure. Then close the line to the syringe at the 3-way stopcock to syringe (circled in picture below).
- 4 Let the reading on the digital barometer stabilize for a moment and then perform the leakage check: for 30 seconds the change of the pressure reading should be less than 20 mbar.



3 Testing and Maintenance

5 If the leakage test is NOT passed, check all connections once more and repeat the test.

| Test | Expected test results |
|---------------|---|
| Leakage Check | Reading on the digital barometer change is less than 20 mbar for 30 seconds (pass/fail) |

NOTE

To ensure the integrity of the Leakage Test Kit (M1013-64002, 451261014851) the following Kit Leak Test Procedure must be performed:

- a. Form a loop with the leakage test kit as shown in the picture below.



- b. Connect the syringe to the 3-way stopcock and the digital barometer to the open tubing.
- c. Draw at the syringe until the digital barometer shows approximately 350 mbar below ambient pressure.
- d. Close the 3-way stopcock to the syringe and wait 5 - 10 seconds. In this time, the overall pressure should stabilize.
- e. After 1 minute, check the pressure. The pressure should not increase more than 8 mbar in 1 minute for the test to pass.
- f. If this test fails, exchange the leakage test kit.

Barometric Pressure Check and Calibration

NOTE

The M3015A with HW Rev C and the M3015B do not require calibration of the barometric pressure. Therefore you will not be able to activate a barometric pressure calibration. If you are using a HW Rev C M3015A or M3015B, perform the barometric pressure check as described below, making sure that only a sample line is connected to the MMS Extension. If the pressure check fails, the M3015A/B needs to be exchanged.

Check the barometric pressure value in the M3015A/B MMS Extension or the MP5 as follows:

- 1 Go into service mode and select **Setup CO₂** menu.
- 2 Connect a calibration line to the Microstream CO₂ input. This activates the pump in the M3015A/B MMS Extension or the MP5.

- 3 The status line at the bottom of the screen displays “CO₂ pressure reading (ambient/cell) xxx/yyy” where xxx is the ambient pressure and yyy is the measured cell pressure. Check whether the ambient pressure value (**x3**) matches (within the acceptable tolerance of ± 12 mm Hg) the reference value you have received. If so, proceed to the leakage check.
If the value is not correct, calibrate as follows.
- Select **CO₂** then select **Barom.Press** to activate a table of values.
 - Select the value in the table which matches the reference value received from a reliable local source (airport, regional weather station or hospital weather station). (The values are displayed with a resolution of 2 mmHg up to 500 mmHg and a resolution of 1 mmHg from 500 mmHg to 825 mmHg.) Note: the selected value must be within $\pm 10\%$ of the current measured ambient pressure, otherwise an error message will occur at restarting the monitor.
 - Confirm the barometric pressure setting.
 - Check that the ambient pressure displayed in the status line at the bottom of the screen is the same as the value which you selected from the list in step b.

| Test | Expected test results |
|---------------------------|--|
| Barometric Pressure Check | x3 = difference between the reference pressure and the measured ambient pressure displayed on the monitor (x3 < 12 mmHg) |

Pump Check (M3015A HW A only)

- 1 Connect the flowmeter inlet to the mCO₂ gas outlet using a flexible connecting tube.
- 2 Connect the calibration line to the mCO₂ inlet.
- 3 Block the inlet of the calibration line using your fingertip and observe the cell pressure on the monitor display. The cell pressure (**x4**) should be more than 120 mmHg below the ambient pressure shown. If the pressure difference is less than 120 mmHg, the pump is not strong enough and you should replace it, irrespective of the Pump OpTime.

| Test | Expected test results |
|------------|---|
| Pump Check | x4 = difference in pressure between cell pressure and ambient pressure displayed on the monitor during occlusion (x4 > 120 mmHg) |

Flow Rate Check and Calibration

Check the flow rate in the M3015A/B MMS Extension or the MP5 as follows:

- 1 Connect the calibration line to the mCO₂ inlet and the flowmeter outlet to the calibration line.
 - 2 Check on the flowmeter the flow that the M3015A/B MMS Extension or MP5 mCO₂ pump draws (**x5**). It should be 50 +15/-7.5 ml/min. If the value is within tolerance, proceed to the CO₂ Gas calibration check.
If the value is not within tolerance, calibrate as follows.
- a. Adjust the flow in the instrument by selecting **Increase Flow** or **Decrease Flow** until it is as close as possible to 50 ml per minute as indicated on the flowmeter gauge.


3 Testing and Maintenance

- b. When you are satisfied that the flow is set as close as possible to 50 ml per minute, select **Store Flow** and confirm the setting. If you do not store the adjusted flow within 60 seconds of the adjustment, the old flow setting is restored.
- c. If you cannot adjust the flow to within tolerance, replace the pump . If you still cannot make the flow adjustment, this indicates a fault in the measurement extension, which must be replaced.
Note that the pump can only be replaced on M3015A with the old hardware Rev. A (i.e. Serial No. DE020xxxxx)

| Test | Expected test results |
|-----------------|---------------------------------|
| Flow Rate Check | Flow rate is 50 +15/-7.5 ml/min |

Noise Check

- 1 With the monitor in service mode, select **Setup CO₂** menu.
- 2 Connect the calibration line, the cal tube, the flow regulator and the 5% calibration gas to the mCO₂ inlet.




- 3 Open the valve to apply the 5% calibration gas and wait until the value is stable.
- 4 Check the noise index (**x6**) displayed next to the CO₂ value on the display (this indicates the level of noise on the CO₂ wave). If the value exceeds 3 mmHg, replace the measurement extension.

| Test | Expected test results |
|-------------|---|
| Noise Check | x6 = noise index displayed on monitor (x6 <3.0) |

CO₂ Cal Check and Calibration

After switching the measurement extension on, wait at least 20 minutes before checking the calibration. Check the calibration of the CO₂ gas measurement as follows:

- 1 Connect the calibration line, the cal tube, the flow regulator and the 5% calibration gas to the mCO₂ inlet.



- 2 Calculate the expected measurement value in mmHg as follows:
 $0.05 \times (\text{ambient pressure}) = \text{value mmHg}$
 for example $0.05 \times 736 = 36.8 \text{ mmHg}$ (with an ambient pressure of 736 mmHg)
- 3 Open the valve on the flow regulator to allow 5% CO₂ gas to flow into the extension. Allow the value to stabilize.
- 4 Check that the value on the instrument (measurement value on the main screen, **x7**) matches the calculated mmHg value $\pm 2.6 \text{ mmHg}$.
 If the value is outside the tolerance, calibrate as described in step 8a to 8e below.
- 5 Disconnect the 5% calibration gas and connect the 10% calibration gas.
- 6 Calculate the expected measurement value and tolerance in mmHg as follows:
 $0.1 \times (\text{ambient pressure}) = \text{value mmHg}$
 $\pm 0.07 \times (\text{value mmHg}) = \text{tolerance}$
 for example $0.1 \times 737 \text{ mmHg} = 73.7 \text{ mmHg}$ (with an ambient pressure of 737 mmHg)
 $\pm 0.07 \times 73.7 \text{ mmHg} = \pm 5.16 \text{ mmHg tolerance}$
- 7 Open the valve on the flow regulator to allow 10% CO₂ gas to flow into the extension. Allow the value to stabilize.
- 8 Check that the value on the instrument (**x8**) matches the calculated mmHg value within the calculated tolerance. If so, the measurement extension is correctly calibrated.
 If the value is outside the tolerance, calibrate as follows.
 - a. Keep the same setup and connect the 5% calibration gas.
 - b. Select **Cal. CO₂**.
 - c. Select the value for the calibration gas. (The default value is 5.0%.)

3 Testing and Maintenance

- d. Open the valve on the calibration gas to allow CO₂ gas to flow into the extension. Allow the value to stabilize before the start of the calibration. Leave the valve open until the instrument gives a prompt that gas can be removed.
- e. The extension calibrates and prompts when calibration is successful.

| Test | Expected test results |
|---------------------------|---|
| CO ₂ Cal Check | x7 = calculated mmHg value ± 2.6 mmHg x8 = calculated mmHg value within calculated tolerance |

Calibration Verification

- 1 Keep the same setup as described in “CO₂ Cal Check and Calibration” on page 83.
- 2 Reopen the 5% gas valve and allow the value to stabilize.
- 3 Check that the value displayed on the monitor is correct within the tolerance (see step above).
- 4 Disconnect the 5% calibration gas and connect the 10% calibration gas.
- 5 Open the valve on the flow regulator to allow 10% CO₂ gas to flow into the extension. Allow the value to stabilize.
- 6 Check that the value displayed on the monitor is correct within the tolerance (see step above).

If one or both values are not within tolerances, you must exchange the M3015A/B MMS Extension or the MP5 mCO₂ Assembly.

| Test | Expected Test Results |
|---------------------------------------|---|
| Leakage Check parts 1 and 2* | x1 = value of part 1 leakage check on flowmeter (x1 < 4.0 ml/min) x2 = value of part 2 leakage check on flowmeter (x2 < 4.0 ml/min) |
| Leakage Check without Flowmeter | reading on the digital barometer change is less than 20 mbar for 30 seconds |
| Barometric Pressure Check | x3 = difference between the reference pressure and the measured ambient pressure displayed on the monitor (x3 < 12 mmHg) |
| Pump Check** | x4 = difference in pressure between cell pressure and ambient pressure displayed on the monitor during occlusion (x4 > 120 mmHg) |
| Flow Check | x5 = difference between measured value and 50.0 ml/min (x5 = 50 + 15 / -7.5 ml/min) |
| Noise Check | x6 = noise index displayed on monitor (x6 < 3.0) |
| CO ₂ Gas Calibration Check | x7 = difference between measured CO ₂ value and calculated value, based on 5% CO ₂ cal. gas. (x7 < 2.6 mmHg) |
| CO ₂ Cal Verification | x8 = difference between measured CO ₂ value and calculated value, based on 10% CO ₂ cal. gas. (x8 < $\pm \{0.07 \times \text{value calculated}\}$) |

* M3015A HW Rev. A and B and FW Revision < P.01.32 only

** M3015A HW Rev. A only

Reset Time Counters

NOTE

This procedure only applies to M3015A with the old hardware Rev. A (i.e. Serial No. DE020xxxxx)

You must check the time counters on the Microstream CO₂ extension before calibrating the instrument. As well, when parts are replaced, the appropriate counters must be reset to zero.

The counters for CO₂ pump, IR Src and Last Cal are displayed in the status line. The values are updated when entering the **Setup CO2** menu.

Observe the following guidelines:

- When calibrating the CO₂ extension, if no parts have been replaced, check the displayed values of **Reset PumpOpTime** and **Reset IRSourceTime** selections to make sure that they are within suggested guidelines for use (15, 000 hours of continuous use). If the counter time is greater than 15, 000 hours, replace the appropriate part. See *Repair and Disassembly* for details.
- When calibrating the CO₂ extension, if parts have been replaced, reset the appropriate values using the **Reset PumpOpTime** and **Reset IRSourceTime** selections. See *Repair and Disassembly* for details.

Resetting the PumpOpTime generates the INOP: “CO₂ OCCLUSION”. To clear this INOP you must perform a flow check and store the flow in service mode (select **Store Flow**).

CO2 Pump / CO2 Scrubber Replacement

NOTE

This procedure only applies to M3015A with the old hardware Rev. A (i.e. Serial No. DE020xxxxx)

Refer to the Repair and Disassembly section for the replacement procedures.

Spirometry Performance Tests

These tests verify the performance accuracy of the M1014A Spirometry module.

Equipment Required

- Leak test kit (Part number: M1014-64100)
- calibrated barometer
- M2785A Pediatric/Adult Flow Sensor
- 500ml calibration syringe, Hans Rudolph model 5550 or equivalent

Flow Test

- 1 Connect the M1014A Spirometry Module to the host monitor and go into service mode.
- 2 Connect the flow sensor to the module.
- 3 Connect the 500ml calibration syringe to the flow sensor. Make sure the syringe is set to the “empty” position.
- 4 Press the **Setup** key on the module and select **Show all Values** in the **Setup Spirometry** menu.

3 Testing and Maintenance

- 5 Pump the calibration syringe back and forth with a steady motion at a rate of 20 cycles and verify that the readings for TV_{exp} and TV_{in} are 500 ± 25 ml.

If the readings are not within the specified range, try another flow sensor. Ensure that the syringe is calibrated correctly and that the procedure is performed exactly as described above. If the test fails again, replace the module.

| Test | Expected test results |
|-----------|--|
| Flow Test | TV _{exp} and TV _{in} are 500 ± 25 ml |

Leakage Test

- 1 Connect the M1014A Spirometry Module to the host monitor and go into service mode.
- 2 Connect the leak test adapter to the module.
- 3 Press the **Setup** key on the module and then select **Show all Values** in the **Setup Spirometry** menu.
- 4 Press the Purge key on the module and start a purge cycle. At the end of the purge cycle, the values for Paw and P_{peak} should both be above 100 cmH₂O.
- 5 Verify that the pressure difference between P_{peak} and Paw remains less than 10 cmH₂O after 30 seconds.

If the readings are not within the specified range or if an INOP (e.g. SPIRO PURGE FAILED) is issued, check the leak test adapter for any leaks. Disconnect the adapter from the module and start the test procedure from the beginning. If the test fails again, replace the module.

| Test | Expected test results |
|--------------|---|
| Leakage Test | Paw and P _{peak} >100 cmH ₂ O |

Barometer Check

- 1 Connect the M1014A Spirometry Module to the host monitor and go into service mode.
- 2 Attach any airway adapter to the module.
- 3 Press the Setup key on the module and then select **Show all Values** in the **Setup Spirometry** menu.
- 4 Check that the barometric reading (PB) is within ± 5 mmHg of a reference barometer.
- 5 If the readings are not within the specified range, check the accuracy of the barometric pressure reference again. If the test fails again, replace the module.

| Test | Expected test results |
|-----------------|--|
| Barometer Check | PB is within ± 5 mmHg of a reference barometer |

NOTE

The built-in barometer cannot be recalibrated.

Cardiac Output (C.O.) Performance Test

These tests check the performance of the cardiac output measurement.

- 1 Connect the patient simulator to the C.O. module using the patient cable.
- 2 Configure the patient simulator as follows:
Injection temperature: 2 °C
Computation Const: 0.542
(Edward's Catheter)
Flow: 5 l/min
- 3 Check displayed value against the simulator configuration.
- 4 Expected test result: C.O. = 5 +/- 1 l/min.

| Test | Expected test results |
|---------------------------------|-----------------------|
| Cardiac Output Performance Test | C.O. = 5 +/- 1 l/min |

Service Tool Procedure, Version 1

This procedure applies for Service Tool M1012-61601 in combination with C.O. modules without option C10 and M3012A MMS extensions with option C05.

- 1 In monitoring mode, connect the C.O. interface cable to the module.
- 2 Connect one side of the service tool to the injectate receptacle of C.O. interface cable and the other side to catheter cable receptacle.
- 3 Enter the **C.O. Procedure** window and check the results. The expected test result is:

$$T_{\text{blood}} = 37.0^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$$

| Test | Expected test results |
|---|---|
| Cardiac Output Service Tool Procedure Version 1 | $T_{\text{blood}} = 37.0^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$ |

Service Tool Procedure, Version 2

This procedure applies only for Service Tool M1012-61601 in combination with C.O. modules with option C10 and for the M3012A MMS Extension with option C10.

- 1 In monitoring mode, connect the C.O. interface cable to the module.
- 2 Connect one side of the service tool to the injectate receptacle of the C.O. interface cable and the other side to the catheter cable receptacle.
- 3 Enter **C.O. Procedure** window and check results for:
 - Method of measurement
 - Arterial Catheter constant
 - T_{blood}
 The expected results are:
 - Transpulmonary
 - 341
 - $T_{\text{blood}} = 37.0^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$
- 4 Make sure the main alarms are switched on.
- 5 Disconnect the Catheter cable receptacle from the service tool

3 Testing and Maintenance

- 6 Enter the Setup C.O Window and change the method of measurement to “Right Heart”
- 7 Enter the C.O. Procedure window and check the T_{inj} value. The expected result is:
 $T_{inj} = 0.0^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$

| Test | Expected test results |
|--|---|
| Cardiac Output Service Tool Procedure Version 2 | $T_{inj} = 0.0^{\circ}\text{C} \pm 0.1^{\circ}\text{C}$ |

BIS Performance Test

These tests check the performance of the BIS measurement.

BISx/BIS Engine Test


- 1 In monitoring mode connect the BIS sensor simulator (P/N: M1034-61650, 453563233731) (for maximum usage please refer to the documentation delivered with the sensor simulator) to the patient interface cable.
- 2 Enter the BIS menu and select **Show Sensor**.
- 3 Start impedance check by pressing **StartCyclicCheck**. Check the displayed results. Expected results are:

| Test | Expected test results |
|----------------------|---|
| BIS Performance Test | Electrode 1 (+): 4-6 k Ω Electrode 2 (Ref): 8-17 k Ω Electrode 3 (1-): 2-4 k Ω Electrode 4 (2-): 3-5 k Ω |

NMT Performance Test

NMT Stimulation Output Test

- 1 Short circuit the stimulation cables by connecting the two cable clamps to each other as shown below.




- 2 In service mode, select the **Setup NMT** menu.
- 3 Select **Start Test**.

- 4 Select **Confirm**.

| Test | Expected test results |
|-----------------------------|--|
| NMT Stimulation Output Test | NMT Stimulation Output Test passed is displayed on the monitor. |


NMT Transducer Test

- 1 Go into Service Mode. In Service Mode the NMT Bar Graph only contains three bars instead of four.
- 2 Place the NMT Transducer on a flat surface with the flat side facing downwards. Two of the three bars in the NMT bar graph should be at the same level and the third one should be higher than the other two.



3 Testing and Maintenance

- 3 Turn the NMT Transducer 180° and place it on a flat surface with the flat side facing upwards. The bar that was higher than the other two before should now be lower than the other two by approximately the same amount.



| Test | Expected test results |
|---------------------|---|
| NMT Transducer Test | First two bars in the NMT bar graph are at the same level, third bar is higher when the flat side of the transducer is facing downwards and lower by the same amount when the transducer is facing upwards. |

Vuelink Performance Test

This test checks the performance of the Vuelink modules.

Tools required: none / external device (i.e. ventilator) and the required Vuelink cable

- 1 Plug the VueLink module into the Philips patient monitor.
- 2 Switch to Configuration Mode of your monitor.
- 3 Depending on your external device, configure the VueLink module as described in the Philips M1032 VueLink Module Handbook "Configuring the VueLink Module (CMS or V24/26)" on page 13 or "Configuring the VueLink Module (IntelliVue Patient Monitor)" on page 45. (Ensure that you have stored the configuration settings before continuing.)
- 4 Change the operating mode of the monitor to Monitoring Mode.
- 5 Press the Setup key on the front of the VueLink module.
- 6 Press the Setup VueLink pop-up key, if setup menu is not already shown.
- 7 In the Setup VueLink menu select Device, and then select the required Device driver.
- 8 Select Confirm to store the selection and wait for the message "Switched to new device"
- 9 Connect the module by plugging one end of the cable connector into the VueLink Module, and the other end into the connector of the external device. Make sure that you use the correct cable option for that device.
- 10 Select the wave segment on the screen, in which you want the waves to be displayed. In the pop-up menu, select Change Wave, and then select WAVE.

- 11 Switch on the external device. After communication is established, information from the external device will be available on the Philips patient monitor.
- 12 We recommend that you confirm with the user that waves and numerics required from the external device are being accurately received.

| Test | Expected test results |
|--------------------------|---|
| Vuelink Performance Test | Information from the external device is available at the Philips patient monitor. |

IntelliBridge Performance Test

This test checks the performance of the IntelliBridge EC10 & EC5 modules.

Tools required: none / external device (i.e. ventilator) and the required IntelliBridge EC5 Module

- 1 Plug the IntelliBridge EC10 module into the Philips patient monitor.
- 2 Connect the Service PC to the IntelliBridge EC10 module and make sure the correct drivers for the external devices are installed. (See the chapter for details).
- 3 Depending on your external device, connect the appropriate EC5 ID module (indicated on the EC5 label) to the external device.
- 4 Connect the EC5 to the EC10 module using the supplied cable.
- 5 Switch the external device on. The connection status LED will flash green until it has correctly identified the external device and started communication. Check that the connection status LED then lights green continuously indicating that communication has been established. Information from the external device should now be available on the Philips patient monitor.
- 6 Select **Main Setup** -> **Measurements** -> **<External Device Name>** to enter the setup menu for the connected device.
- 7 Select **Setup Waves** or **Setup Numerics** and make any required changes.
- 8 Close the setup menu.
- 9 Select the wave segment on the screen, in which you want the waves to be displayed. In the pop-up menu, select Change Wave, and then select WAVE.
- 10 We recommend that you confirm with the user that waves and numerics required from the external device are being accurately received. If the external device has a demo mode, use this.

| Test | Expected test results |
|--------------------------------|--|
| IntelliBridge Performance Test | Numerics are visible on screen (pass/fail) |

EEG, SvO₂ (SO₂) and tcGas Performance Tests

The EEG and SvO₂ (SO₂) parameters do not require performance tests because the modules perform internal self-tests regularly. These tests suffice for performance testing of these three parameters.

Since the tcGas Module is calibrated regularly it also does not require a separate performance test.

Nurse Call Relay Performance Test


The nurse call relay performance test can be performed at the Modular Jack 6P6C connector.

This test checks the operation of the Nurse Call Relay. The Nurse Call Relay test is recommended for customer sites where the nurse call is in use. The Nurse Call relay functions as follows:

- Standard Operation—connector contact 1-2 open; 1-5 closed.
- Alarm Condition—connector contact 1-2 closed; 1-5 open.

Tools required: Ohmmeter.

- 1 Plug a 6P6C Modular Plug into the Nurse Call Relay connector.
- 2 Connect the ohmmeter.
- 3 When no alarm occurs, connector contacts 1-2 are open and connector contacts 1-5 are closed. When an alarm occurs, connector contacts 1-2 are closed and connector contacts 1-5 are open.



- 4 The expected test result is: Alarm condition - Connector contacts 1-2 are closed and connector contacts 1-5 are open.

| Test | Expected test results |
|-----------------------------------|--|
| Nurse Call Relay Performance Test | Alarm Condition— Connector contacts 1-2 are closed and Connector contacts 1-5 are open |

Multi-Port Nurse Call Connector Test (Flexible Nurse Call)

This test checks the operation of the Flexible Nurse Call Relay. The Nurse Call Relay test is recommended for customer sites where the nurse call is in use. The following diagram and table show the pins and relay identifiers of the connector:

| Pin | Cable Color Coding | Relay |
|-----|--------------------|------------|
| 1 | black | R2-closure |
| 2 | brown | R2-middle |
| 3 | red | R2-opener |
| 4 | orange | R3-closure |
| 5 | yellow | R3-middle |
| 6 | green | R3-opener |
| 7 | blue | n/a |
| 8 | purple | n/a |
| 9 | gray | n/a |

| Pin | Cable Color Coding | Relay |
|-----|--------------------|-------------------|
| 10 | white | n/a |
| 11 | pink | R1-closure |
| 12 | light green | R1-middle |
| 13 | black/white | R1-opener |
| 14 | brown/white | n/a |
| 15 | red/white | n/a |
| 16 | orange/white | n/a |
| 17 | blue/white | R_failure_closure |
| 18 | purple/white | R_failure_middle |
| 19 | green/white | R_failure_opener |
| 20 | red/black | n/a |

The Nurse Call relay functions as follows:

- During standard operation R1,R2,R3 _opener are closed; R1,R2,R3_closure are open.
- During alarm condition—R1,R2,R3_opener are open; R1,R2,R3_closure are closed.

Tools required: Ohmmeter.


- 1 Plug an M8087-61001 cable into the Nurse Call Relay connector.
- 2 Connect the ohmmeter and measure the pins as indicated in the diagram and table.
- 3 The relay contacts should behave as described above. The behavior may vary depending on configuration choices. See the Configuration Guide for details on Alarm Relay settings.
- 4 The expected test results depend on the relay contact used. Please check that the correct relay activity is initiated during alarm condition.

| Test | Expected test results |
|--------------------------------------|--|
| Multi-Port Nurse Call Connector Test | Correct relay activity is initiated during alarm condition (pass/fail) |

MSL Assurance Test

Visually inspect all MSL connector sockets (cable/monitor/FMS/MMS).

- 1 Make sure that the pins of the connectors are not jolted.
- 2 Make sure that no pin is bent inwards or outwards.
- 3 Exchange connectors that show any evidence of damage or breakage



Examples of damaged connectors

| Test | Expected test results |
|--------------------|---|
| MSL Assurance Test | Pins of connector not jolted/bent (pass/fail) |

Power Loss Alarm Buzzer Performance Test

- 1 Switch on the monitor.
- 2 Disconnect the monitor from AC power.
- 3 The Power Loss Alarm Buzzer should beep for about one minute.
- 4 To switch off the alarm sound, either press the power button or connect the monitor to AC power

| Test | Expected test results |
|--|-----------------------|
| Power Loss Alarm Buzzer Performance Test | Beep for one minute |

IntelliVue 802.11 Bedside Adapter Communication Test (Rev H.0 or higher)

- 1 Make sure the LAN cable is disconnected from the rear of the monitor, then switch on the monitor.
- 2 Go into Service Mode and select **Main Setup -> Network -> Setup WLAN**. In the **Setup WLAN** menu:
 - set **Mode** to either **802.11Ah**, **802.11G**, **802.11Bg** (not recommended), **Auto** (not recommended) or **None** (this setting disables the wireless LAN functionality permanently), to match your wireless infrastructure installation.
 - set **SSID** to match your installation.

- set the **Country** code to “1000”. Setting the country code to this value will automatically adjust the regulatory domain to match the configuration of the infrastructure. Do not set the country code to values other than “1000” unless otherwise instructed.
- set the **Security Mode** to either **WEP**, **WPA(PSK)** or **WPA2(PSK)** to match your installation.
- Enter the required keys/passwords for **WEP**, **WPA(PSK)** or **WPA2(PSK)**.

NOTE

WPA2(PSK) was introduced with monitor software revision G.0.(Required WLAN adapter firmware for MP5-90: A.02.24, for MP2/X2: A.00.20)

- 3 Select **Main Setup** -> **Network** -> **WLAN Diagnostic** to access the WLAN Diagnostic window.
- 4 Proper installation of the IntelliVue 802.11 Bedside Adapter is assured by connecting to an access point over the wireless link. Place the monitor with the IntelliVue 802.11 Bedside Adapter installed in close proximity to the access point (e.g. if the access point is mounted on the ceiling, place the monitor directly below). Wait until the **Conn.Status** field in the service window shows *Authenticated* (for Rel. C.0 monitors) or *Connected* (for Rel D.0 or higher). Take the monitor approximately 5 m away from the access point. There should be no walls or other obstacles between the monitor and the access point. The following should apply:
 - Observe the **RSSI** (Received Signal Strength Indicator) value for at least 5 - 10 seconds. The **RSSI** value will fluctuate but should stay above 30 in a 5 m distance from the access point used. The wireless link should be active, i.e. the **Conn.Status** field should be *Authenticated* (for Rel. C.0 monitors) or *Connected* (for Rel D.0 or higher), and the other fields should contain values. If the **RSSI** value is significantly lower, check the distance to the access point and the antenna orientation at the monitor. The antenna orientation should be vertical, but the physical placement of the monitor or other equipment within its vicinity as well as walls or other obstacles may influence the antenna orientation required to receive the best RSSI value.
- 5 If this test fails, retry in a different physical area with a different access point.
- 6 Perform the Wireless Switch test blocks as described in the Philips IntelliVue 802.11 a/g Infrastructure Installation and Configuration Guide.

| Test | Expected test results |
|--|-----------------------|
| IntelliVue 802.11 Bedside Adapter Performance Test | RSSI value above 30 |

IIT Communication Test (Rev H.0 or higher)


- 1 Make sure the LAN cable is disconnected from the rear of the monitor, then switch on the monitor.
- 2 Go into Service mode and, select **Main Setup** -> **Network** -> **Setup IIT**. In the Setup IIT menu, set the **RF Access Code** in each profile to match your installation.
- 3 Go into Service Mode. Select **Main Setup** -> **Network** -> **IIT Diagnostic** to access the Instrument Telemetry Diagnostic window.
- 4 Proper installation of the IIT module is assured by connecting to an access point over the wireless link. Place the monitor with the IIT module installed in close proximity to the access point (e.g. if the access point is mounted on the ceiling, place the monitor directly below). Wait until the **Conn.Status** field in the Instrument Telemetry Service window shows *Active*. *Take the monitor approximately 5 m away from the access point. There should be no walls or other obstacles between the monitor and the access point. The following should apply:*

3 Testing and Maintenance

- Observe the **RSSI** (Received Signal Strength Indicator) value for at least 5 - 10 seconds. The **RSSI** value should be around -50 ± 10 in a 5 m distance from the access point used and the IIT link should be active, i.e. the **Conn.Status** field should be *Active* and the other fields should contain values. If the **RSSI** value is significantly lower, check the distance to the access point and the antenna orientation at both the monitor and the access point (both should be vertical).
 - Remove the antenna. The **RSSI** value should be around -90 ± 10 . The IIT link may be active but the connection could be unreliable. The **Conn. Status** field may toggle between *Inactive* and *Seeking*. If the difference between the **RSSI** values measured with and without antenna is significantly lower, check the antenna and the antenna connector for damage and verify that the cable from the IIT adapter to the antenna connector plate is connected properly.
- 5 If this test fails, retry in a different physical area with a different access point.
- Error Conditions:**
- The field **MAC IIT** should show a value **unequal to** 0000 0000 0000. If it does not, there is a communication problem between the monitor and the IIT adapter.
 - With an incorrect RF Access Code or an incorrect or defective antenna installation, the fields **IP Address, Server IP, Subnet Mask**, and **RSSI** in the Instrument Telemetry Service window will stay blank. The field **Conn. Status** will slowly toggle between *Inactive* and *Seeking*.
- 6 Perform the Access Point Controller (APC) test blocks as described in the Philips IntelliVue Wireless Network Installation and Configuration Guide.

| Test | Expected test results |
|------------------------|--|
| IIT Communication Test | IIT Communication without interference |

Short Range Radio (SRR) Performance Test

- 1 Make sure that the short range radio interface is configured as follows: SRR On and appropriate channel selected.
- 2 Assign a telemetry transceiver or IntelliVue CL Cableless Measurement device to the IntelliVue Monitor according to the procedure described in the Instructions for Use of the patient monitor.
- 3 Check that the following conditions are fulfilled:
 - a. Place the telemetry transceiver or CL device close to the monitor.
 - b. The telemetry transceiver or CL device status is displayed on the monitor in the measurement selection window.
 - c. Waves or numerics from the telemetry transceiver or CL device are displayed on the monitor. There are no dropouts or gaps in waves or numeric transmission.
 - d. The battery status of the telemetry transceiver or CL device is displayed in the measurement selection window.
 - e. The Signal Quality Indicator shows at least 
- 4 Check that the data from the telemetry transceiver or CL device is transmitted to the monitor within a 1m radius and that there are no dropouts or gaps in waves or numerics.
- 5 Check whether the connection remains stable within a 5m radius from the monitor.
- 6 Switch on all telemetry transceivers or CL devices used on the site and check that there are no interferences between the transceivers and their assigned monitors.

- 7 Check and record the coverage area of the telemetry transceivers or CL devices and inform the customer about this coverage area.

| Test | Expected test results |
|--|--|
| Short Range Radio (SRR) Performance Test | SRR Communication without interference within the coverage |

Mounting Integrity Test

Perform the Mounting Integrity Test

- whenever you have removed and reassembled a quick mount
- if one or both of the quick mount screws are loose
- if there is a clearance between the quick mount and the monitor bottom housing
- if the monitor mounting is unstable

Remove the monitor from the mount and disassemble the quick mount. Ensure that the that the threading of the MX600/700 is not damaged or separated from the chassis.

If the quick mount is damaged, exchange the quick mount.

Ensure that all quick mount screws are tight (3.5 Nm). Test the quick mount by pressing the quick release button. If it comes back out gradually and regularly, the quick mount is inserted correctly. If it gets stuck, the quick mount is not centered and must be reinserted correctly.

If you notice any damage to the threading of the MX600/700 chassis, send the MX600/700 in for bench repair.

| Test | Expected test results |
|-------------------------|---|
| Mounting Integrity Test | All quick mount screws are tight. No damage to quick mount. No damage to threading of MX600/700. Quick release button comes back out gradually and regularly. |

Reporting of Test Results

Philips recommends all test results are documented in accordance with local laws. Authorized Philips personnel report the test result back to Philips. While hospital personnel (biomedical engineers or technicians) do not need to report results to Philips, Philips recommends that they record and store the test results in accordance with local laws.

The following table lists what to record after completing the tests in this chapter. Record the results in the empty column in the Test and Inspection Matrix.

The following is a guide as to what your documentation should include:

- Identification of the testing body (for example, which company or department carried out the tests).
- Name of the person(s) who performed the tests and the concluding evaluation.
- Identification of the device(s) and accessories being tested (serial number, etc.).
- The actual tests (incl. visual inspections, performance tests, safety and system tests) and measurements required
- Date of testing and of the concluding evaluation.

3 Testing and Maintenance

- A record of the actual values of the test results, and whether these values passed or failed the tests.
- Date and confirmation of the person who performed the tests and evaluation.

The device under test should be marked according to the test result: passed or failed.

Carrying Out and Reporting Tests

Test Report

| | |
|---|---|
| Testing Organization: | (Check one of the following three options) Test before putting into service (reference value) Recurrent Test Test after Repair |
| Name of testing person: | |
| Date: | |
| Responsible Organization: | |
| Device Under Test: | ID-Number |
| Product Number: | Serial No.: |
| Accessories: | |
| Measurement Equipment (Manufacturer, Type, Serial No., Calibration Date): | |
| Safety Test Method used | |
| Functional Test (parameters tested): | |
| Mains voltage and frequency used during safety testing: | |

Test and Inspection Matrix

| Test | Test or Inspection to be Performed | Expected Test Results | Record the Results (mandatory for Philips Personnel only) | |
|--|---|--|---|----------------|
| | | | What to record | Actual Results |
| Visual Inspection | Perform Visual Inspection | Pass or Fail | V:P or V:F | |
| Power On | Power on the unit. Does the self-test complete successfully | If Yes, Power On test is passed | PO:P or PO:F | |
| Noninvasive Blood Pressure Performance Tests | Perform the Accuracy Test | X1 = value displayed by monitor Difference \leq 3mmHg | PN:P/X1 or PN:F/X1 | |
| | Performance Leakage Test | X2 = leakage test value X2 < 6 mmHg | PN:P/X2 or PN:F/X2 | |
| | Performance Linearity Test | X3 = value displayed by monitor Difference \leq 3mmHg | PN:P/X3 or PN:F/X3 | |
| | Performance Valve Test | X4 = value < 10 mmHg | PN:P/X4 or PN:F/X4 | |
| Temperature Performance Test | Perform the Temperature Performance Test | X1 = 40°C \pm 0.2°C or 100°F \pm 0.4°F | PT: P/X1 or PT: F/X1 | |

3 Testing and Maintenance

| Test | Test or Inspection to be Performed | Expected Test Results | Record the Results (mandatory for Philips Personnel only) | |
|-----------------------------|---|--|---|----------------|
| | | | What to record | Actual Results |
| All other performance tests | Perform the remaining parameter performance tests, if applicable | See expected results in test procedures | P: P or P: F | |
| Safety (1) | Perform Safety Test (1): Protective Earth Resistance | With mains cable: Maximum impedance (X1): ≤300 mOhms | S(1):P/X1 or S(1):F/X1 | * |
| Safety (2) | Perform Safety Test (2): Equipment Leakage Current - Normal Condition. | With mains cable: Maximum leakage current (X1):≤ 100 μA | S(2): P/X1 or S(2): F/X1 | * |
| Safety (3) | Perform Safety Test (3): Equipment Leakage Current - Single Fault Condition (Open Earth) | With mains cable: Maximum leakage current (X2):≤ 300 μA | S(3): P/X2 or S(3): F/X2 | * |
| Safety (4) | Perform Safety Test (4): Applied Part Leakage Current - Single Fault Condition, mains on applied part. | Maximum leakage current (X1): ≤50 μA (CF) | S(4): P/X1 or S(4): F/X1 | * |
| System (Sys 1-2) | Perform the system test according to subclause 19.201 of IEC/EN 60601-1-1 or IEC/EN 60601-1 edition 3, clause 16, if applicable, after forming a system | Equipment Leakage Current: Sys1 ≤ 100 μA (Normal Condition) Sys2 ≤ 300μA (Single Fault Condition) | Sys: PSys1/PSys2 or Sys: FSys1/FSys2 | * |
| System (Sys 3) | Perform the system test according to subclause 19.201 of IEC/EN 60601-1-1 or IEC/EN 60601-1 edition 3 clause 16, if applicable, after forming a system | Protective Earth Leakage Current if medical electrical system components are connected to the same Multiple Portable Socket Outlet: Sys3 ≤ 300 μA | Sys: PSys3 or Sys: FSys3 | * |

Key: P = Pass, F = Fail, X or Sys = test value to be recorded, * = Record the worst-case results and the associated switch positions (e.g. normal/reverse polarity)

NOTE

All values for current and voltage are the root mean square (r.m.s.) values, unless otherwise stated.

Evaluation

| | Yes | No |
|---|-----|----|
| Safety and Functional Test passed | | |
| Repair required at a later date, safety and functional test passed | | |
| Device must be taken out of operation until repair and passed tests | | |
| Device failed and must be taken out of operation. | | |
| Notes: | | |
| | | |
| Next Recurrent Test: | | |
| Name: _____ | | |
| Date/Signature: _____ | | |

Evaluation of Test Results

The evaluation of the test results must be performed by appropriately trained personnel with sufficient product, safety testing and application knowledge.

If any test results are between 90% and 100% of the respective expected result, the previously measured reference values must be taken into consideration for the assessment of the electrical safety of the device under test. If no reference values are available, you should consider shorter intervals between upcoming recurrent tests.

NOTE

If any single test fails, testing must be discontinued immediately and the device under test must be repaired or labeled as defective. Be sure to inform the user about the test failure in writing.

Other Regular Tests


The care and cleaning requirements that apply to the monitor and its accessories are described in the Instructions for Use. This section details periodic maintenance procedures recommended for the monitor and its accessories.

Touchscreen Calibration

To access the touchscreen calibration screen:

- 1 Enter service mode
- 2 Select **Main Setup**
- 3 Select **Hardware**

4 Select **Touch Calibration**



Touchscreen Calibration Screen

Make sure you complete the calibration procedure without powering off the monitor mid-way. If the monitor is powered off after the first point is touched, the touch panel will be deactivated until the touch calibration is performed again.

If the touchscreen is accidentally mis-calibrated by selecting the wrong spot, you must use another input device to re-enter calibration mode. If you have the support tool, you can select **Reset Touch Calibration to Default** and it will create a rough calibration which will allow you to access the calibration menu again via the touchscreen.

Please refer to the documentation shipped with your selected display for further details on touchscreen calibration procedures.

NOTE

If a touchscreen calibration is started on a multiple display system, the calibration is started for all displays at the same time.

Disabling/Enabling Touch Operation

There are two ways to disable/enable touchscreen operation:


- 1 To *temporarily* disable touchscreen operation of the monitor, press and hold the **Main Screen** key. A padlock symbol will appear on the key. Press and hold the **Main Screen** key again to re-enable touchscreen operation.
- 2 To *permanently* disable touchscreen operation:
 - a. Enter Service Mode.
 - b. Select **Main Setup**
 - c. Select **User Interface**
 - d. Change the **Touch Enable** selection to **no**.
To re-enable touchscreen functionality change the **Touch Enable** selection to **yes**.

Printer Test Report

To verify your printer configuration you may want to print a test report.

To print a test report select **Main Setup -> Reports -> Setup Printers -> Print Test Rep.**

Your test report should look like this:



After Installation, Testing or Repair

Before handing the patient monitor over to the end-user, make sure it is configured appropriately and that it is in monitoring mode. Ensure that the user receives the current revision of the monitor documentation.

Troubleshooting

Introduction

This section explains how to troubleshoot the monitor if problems arise. Links to tables that list possible monitor difficulties are supplied, along with probable causes, and recommended actions to correct the difficulty.

How To Use This Section

Use this section in conjunction with the sections *Testing and Maintenance* and *Parts*. To remove and replace a part you suspect is defective, follow the instructions in the section *Repair and Disassembly*. The *Theory of Operation* section offers information on how the monitor functions.

Who Should Perform Repairs

Only qualified service personnel should open the monitor housing, remove and replace components, or make adjustments. If your medical facility does not have qualified service personnel, contact Philips' Response Center or your local Philips representative.

WARNING

High Voltage - Voltages dangerous to life are present in the instrument when it is connected to the mains power supply. Do not perform any disassembly procedures (other than server and extension removal) with power applied to the instrument. Failure to adhere to this warning could cause serious injury or death.

Replacement Level Supported

The replacement level supported for this product is to the printed circuit board (PCB) and major subassembly level. Once you isolate a suspected PCB, follow the procedures in the Repair and Disassembly section, to replace the PCB with a known good PCB. Check to see if the symptom disappears and that the monitor passes all performance tests. If the symptom persists, swap back the replacement PCB with the suspected malfunctioning PCB (the original PCB that was installed when you started troubleshooting) and continue troubleshooting as directed in this section.

Software Revision Check

Some troubleshooting tasks may require that you identify the Software Revision of your monitor. You can find the software revision along with other information, such as the system serial number, in the monitor revision screen. To access the monitor revision screen:

- 1 Enter the Main Setup menu and select **Revision**
- 2 Select **Product**
- 3 Select **Software Revision**
- 4 Select the pop-up key for the device you want to check (e.g. **M8004A** or **M3001A**)

NOTE

The part numbers listed in the monitor revision screen do not necessarily reflect the part numbers required for ordering parts. Please refer to the *Parts* section for the ordering numbers.

NOTE

The system serial number can also be found on the lower right corner on the front of the monitor.

Software Compatibility Matrix

For a detailed software compatibility matrix, please refer to the IntelliVue Compatibility Matrix on your Documentation DVD.

For further information on M3001A HW/SW compatibility, please refer to the *Parts* section.

Compatibility with MMS

The following table shows the compatibility between the monitor and MMS software revisions.

| Monitor Software | MMS Software | | | | | | | | | | |
|------------------|--------------|-----|-----|-----|-----|-----|-----|-----|-----|------|-----|
| | A.2 | B.0 | B.1 | C.0 | D.0 | E.0 | F.0 | G.0 | H.0 | H.1x | J.x |
| H.1x | No | No | No | No | No | No | Yes | Yes | Yes | Yes | Yes |
| J.x | No | No | No | No | No | No | No | Yes | Yes | Yes | Yes |

Compatibility with FMS-8

The following table shows the compatibility between the monitor and FMS-8 software revisions.

| Monitor Software | FMS-8 Software | | | | | | | | | | |
|------------------|----------------|-----|-----|-----|-----|-----|-----|-----|-----|------|-----|
| | A.2 | B.0 | B.1 | C.0 | D.0 | E.0 | F.0 | G.0 | H.0 | H.1x | J.x |
| H.1x | No | No | No | No | No | No | No | No | No | Yes | No |
| J.x | No | No | No | No | No | No | No | No | No | No | Yes |

Compatibility with FMS-4

The following table shows the compatibility between the monitor and FMS-4 software revisions.

| Monitor Software | FMS-4 Software | |
|------------------|----------------|-----|
| | H.1x | J.x |
| H.1x | Yes | No |
| J.x | No | Yes |

Compatibility with Information Center

The following table shows the compatibility between the monitor and Information Center software revisions.

| Monitor Software | Information Center Software | | | | | | | | | | | |
|------------------|-----------------------------|-----|------|-----|-----|-----|-----|-----|-----|-----|-----|-------------|
| | D.01 | E.0 | E.01 | F.0 | G.0 | H.0 | J.0 | K.0 | L.0 | M.0 | N.0 | PIIC iX A.0 |
| H.1x | No | No | No | No | No | No | Yes | Yes | Yes | Yes | Yes | No |
| J.x | No | No | No | No | No | No | No | No | Yes | Yes | Yes | Yes |

The following table shows the compatibility between the MP2/X2/MP5 and the Information Center software revisions if the X2/MP5 are used in companion mode i.e. as monitor and measurement module.

Table 3 Compatibility of X2/MP5 with the IntelliVue Information Center (companion mode use model)

| MP5/X2 Software | Information Center Software | | | | | | | | | | | |
|-----------------|-----------------------------|-----|------|-----|-----|-----|-----|-----|-----|-----|-----|-------------|
| | D.01 | E.0 | E.01 | F.0 | G.0 | H.0 | J.0 | K.0 | L.0 | M.0 | N.0 | PIIC iX A.0 |
| H.1x | No | No | No | No | No | No | No | Yes | Yes | Yes | Yes | No |
| J.x | No | No | No | No | No | No | No | No | Yes | Yes | Yes | Yes |

Number of Supported Parameter Modules

| Module Number | Module Description | Number of Supported Modules | Comments |
|---------------|---------------------------|-----------------------------|----------|
| M1006B | Pressure | 4 | |
| M1006B #C01 | Pressure w/ Analog Out | 4 | |
| M1011A | ScvO2 | 1 | |
| M1012A | Cardiac Output | 1 | |
| M1012A #C10 | Cardiac Output with PiCCO | 1 | |

| Module Number | Module Description | Number of Supported Modules | Comments |
|---|----------------------------------|-----------------------------|----------|
| M1014A | Spirometry | 1 | |
| M1018A | tcpO2 | 1 | |
| M1020B #A01 | SpO2 Philips FAST for IntelliVue | 2 | |
| M1020B #A02 | SpO2 Nellcor OxiMax compatible | 2 | |
| Masimo SET IVM Module distributed by Philips (internal Order No. M1020B #A03) | SpO2 Masimo SET | 2 | |
| M1021A | SvO2 | 1 | |
| M1027A | EEG | 1 | |
| M1029A | Temp | 4 | |
| M1032A #A01 | VueLink Auxiliary | 2 | |
| M1032A #A02 | VueLink Ventilator | 2 | |
| M1032A #A03 | VueLink Gas Analyzer | 2 | |
| M1032A #A04 | VueLink Anesthesia | 2 | |
| M1032A #A05 | VueLink Auxiliary+ | 2 | |
| M1034A | BIS | 1 | |
| M1116B | Recorder | 1 | |
| 865115 | EC10 IntelliBridge | 2 | |

Obtaining Replacement Parts

See *Parts* section for details on part replacements.

Troubleshooting Guide

Problems with the monitor are separated into the categories indicated in the following sections and tables. Check for obvious problems first. If further troubleshooting instructions are required refer to the Troubleshooting Tables.

Taking the recommended actions discussed in this section will correct the majority of problems you may encounter. However, problems not covered here can be resolved by calling Philips Response Center or your local representative.

Checks for Obvious Problems

When first troubleshooting the monitor, check for obvious problems by answering basic questions such as the following:

- 1 Is the power switch turned on?
- 2 Is the AC power cord connected to the instrument and plugged into an AC outlet?
- 3 Is the MSL cable connected correctly?
- 4 Are the MMS and, if present, the MMS Extension connected correctly?
- 5 Are the cables connected properly to the FMS?
- 6 Are the parameter modules plugged into the FMS correctly?

Checks Before Opening the Instrument


You can isolate many problems by observing indicators on the instrument before it is necessary to open the instrument.

Checks with the Instrument switched Off

- AC connected:
 - AC Power LED is on (green).
- No AC connected:
 - All LEDs are off.

Checks with the Instrument Switched On, AC connected

When the monitor is first switched on, all the front-panel LEDs the Power on/Error LED and the AC Power LED light up momentarily. The location of the front-panel LEDs is shown in the following photograph:



Initial Instrument Boot Phase

The following tables describe the regular initial boot phase of the monitor and its components. If the boot phase does not proceed as described below go to [for Troubleshooting information](#).

Monitor Boot Phase:

4 Troubleshooting

For these steps it is assumed that the Monitor is powered correctly. This is indicated by the green Power On LED.

| Time (sec.) after Power On | Event |
|----------------------------|---|
| 0 | AC Power LED is always on when monitor is connected to AC Power. When the Power On/Off button is pressed the red error LED switches on immediately. |
| 3 | The alarm LEDs are switched on with low intensity. Colors: Left LED: cyan; Middle LED: red; Alarm Suspend LED (right): red. The red error LED is switched to green On/Standby LED. |
| 4 | Boot Screen with the Philips Logo appears on the display. Test Sound is issued. |
| 5 | All Alarm LEDs are switched off. |
| 6 | Alarm LEDs are tested in the following sequence: Cyan on-off (left LED only) Yellow on-off (left & middle LED) Red on-off (all LEDs) |
| 8 | Boot Screen with the Philips Logo disappears Fixed screen elements (for example smart keys, alarm fields) appear on the screen. |
| 15-30 | First measurement information appears on the screen, user input devices (for example Navigation Point, Mouse, Touch) are functional |

Flexible Module Rack Boot Phase

For these steps it is assumed that the Flexible Module Rack is connected via MSL-cable to the monitor.

| Time (sec.) after Monitor Power On | Event |
|------------------------------------|---|
| 0 | Red Error LED switches on immediately |
| 1 | Green "Ready" LED switches on |
| 3 | Red Error LED is switched off |
| 5 | Module Power is switched on |
| 5-8 | Module Status LEDs blink once or twice (Module dependent) |

NOTE

The boot phase times may vary depending on the hardware and software revision of your monitor.

Troubleshooting Tables

The following tables list troubleshooting activities sorted according to symptoms. The possible causes of failure and the remedies listed in the troubleshooting tables should be checked and performed in the order they appear in the tables. Always move on to the next symptom until the problem is solved.

How to use the Troubleshooting tables

The possible causes of failure and the remedies listed in the troubleshooting tables should be checked and performed in the order they appear in the tables. Always move on to the next symptom until the problem is solved.

- “Boot Phase Failures” on page 112
- “Integrated Display is blank” on page 115
- “Integrated Touch Screen not functioning” on page 116
- “External Display is blank (Slave Display)” on page 116
- “External Touch Display not functioning” on page 118
- “Remote Control (wired)” on page 119
- “Remote Control (wireless)” on page 120
- “iPC” on page 120
- “Keyboard/Mouse not functioning” on page 122
- “Bedside Network Status Icons” on page 122
- “Network related problems” on page 124
- “IntelliVue 802.11 Bedside Adapter Problems” on page 126
- “IIT-related Problems” on page 128
- “Short Range Radio Interface Problems” on page 129
- “Multi-Measurement Module” on page 130
- “Alarm Lamps” on page 132
- “Alarm Tones” on page 133
- “Individual Parameter INOPS” on page 133
- “Flexible Module Rack” on page 134
- “Printer” on page 135
- “Recorder” on page 136
- “MIB / RS232” on page 137
- “USB” on page 138
- “Nurse Call Relay” on page 139
- “Troubleshooting the ECG OUT” on page 139

Boot Phase Failures

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--------------------------|--|---|
| AC LED does not light up | AC Connection not ok | Check that the AC-Mains are powered and the power cord is ok and connected |
| | LED defective | Try to switch on the monitor. If it operates normally , the LED is defective => exchange Power Switch board |
| | Power Switch board not connected to the main board | Check if power switch board is connected correctly to the Main Board |
| | Power supply defective | Remove power supply and check if output voltage is within the specifications (23.5V - 24.5V).Exchange power supply if defective |
| | iPC defective | Disconnect Power Cable to iPC and check again |
| | I/O board defective | Remove I/O board and check again |
| | Connector board defective | Remove connector board and check again |
| | Panel adapter board defective | Remove panel adapter board and check again |
| | Mixer board defective | Remove mixer board and check again |
| | Main Board defective | Exchange Main Board |

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|--|--|
| Green Power On LED and Red Error LED remain off after pressing power on button: | Power Switch Micro Controller hung | Unplug AC Mains and replug after 10 seconds. Try to switch on the monitor again. |
| | Power switch board not connected to the main board | Check if power switch board is connected correctly to the main board. |
| | Power Switch Board defective | Exchange Power Switch Board and try to switch the monitor on again. |
| | I/O Board defective | Remove all I/O boards and try to switch the monitor on again |
| | Connector board defective | Remove connector board and check again |
| | Mixer board defective IIT defective WLAN board defective Panel adapter board defective Alarm LED board defective | Disconnect all cables and boards (except Power Cable to the main board and Power Switch cable): - mixer board - IIT - WLAN - panel adapter - Alarm LED then try to switch the monitor on again |
| | Main board defective | Exchange main board. Add boards in reverse order and try again with each board. |
| Red Error LED stays on continuously | External connected device defective | disconnect all external cables (except AC) and switch the monitor on again |
| | I/O Board defective | Remove all I/O boards and switch the monitor on again. |
| | Connector board defective | Remove connector board and check again. |
| | Mixer board defective IIT defective WLAN board defective Panel adapter board defective Alarm LED board defective | Disconnect all cables and boards (except Power Cable to the main board and Power Switch cable): - mixer board - IIT - WLAN - panel adapter - Alarm LED then try to switch the monitor on again |
| | Main board defective | Exchange Main board |

4 Troubleshooting

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|------------------------------|---|
| Red Error LED blinks (indicating cyclic reboots) | Hardware Failure | connect Support Tool directly to monitor with crossover cable and start “search for defective devices” If no device is detected, proceed as described above in section “Red error LED stays on continuously” |
| | Software Fault | If the Support Tool can detect the device and it indicates the Operating Mode is ‘Boot’, download and store the status log. Reload software and re-clone the monitor. If this fixes the problem e-mail the status log to your local response center |
| | Hardware Failure | If this does not rectify the problem follow instructions under “Red Error LED stays on continuously”. |
| Alarm LEDs remain off: | Alarm LED board is defective | Check for INOPS and follow instructions Exchange Alarm LED board |
| | Main board defective | Exchange Main board |
| No Test Sound issued | | check for INOPs and follow instructions |
| | Speaker defective | exchange speaker |
| | Main board defective | exchange main board |

Integrated Display is blank

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|--|---|
| Integrated display is blank or brightness is reduced (The information listed in this table is only valid if the boot phase has completed without error. See table for a description of the Boot phase.) | Display brightness is reduced when room temperature, or instruments placed near patient monitor, causes the monitor display to overheat. | Instrument should be placed in an environment that does not exceed 40 degrees C or below 5 degrees C. |
| | | If you have an external display, connect it to the video port. If the external display works, you can eliminate the connector board and the main board as the cause of failure. |
| | Backlight Inverter Cable not connected | Check cable connection of Panel Adapter Board to Backlight Inverter Board |
| | Backlight tubes defective | Replace backlight tubes |
| | Backlight Inverter board defective | If backlight tubes have already been replaced, replace backlight inverter board. |
| | Panel Adapterboard defective | Replace panel adapter board |
| | LCD Flat panel defective | Replace LCD Flat panel |
| Main board defective | Replace main board | |

Integrated Touch Screen not functioning

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|------------------------------|---|--|
| Touch Screen not functioning | Touchscreen functionality has been temporarily disabled | Check if touchscreen functionality has been temporarily disabled (padlock symbol on Main Screen key). If yes, press and hold the Main Screen key to re-enable touchscreen operation. |
| | Touchscreen functionality has been permanently disabled | In service mode, select Main Setup -> User Interface and change the “Touch Enable” selection to “yes”. |
| | Touch screen not connected | Check connection from touch screen to panel adapter board |
| | Panel adapter board defective | Replace panel adapter board |
| | Touch screen defective | Replace touch screen assembly |
| | Main board defective | Replace main board |
| Touch Position invalid | Touch not calibrated | Perform touch calibration: 1. Enter the Main Setup Menu 2. Select Hardware 3. Select Touch Calibration |

External Display is blank (Slave Display)

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---------------------------|---|--|
| External Display is blank | | If integrated display is also blank proceed as described under “Integrated Display is blank” |
| | Video cable to external display not connected | Check video cable connection to external display |
| | External display has no power | Check electricity supply of external display |
| | External display is defective | Check external display and video cable on another monitor or PC |
| | Connector Board defective | Replace Connector board |
| | Main board defective | Replace main board |

External Display Connected to iPC

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---------------------------|---|--|
| External Display is blank | | If integrated display is also blank proceed as described under “Integrated Display is blank” |
| | Video cable to external display not connected | Check video cable connection to external display |
| | External display has no power | Check electricity supply of external display |
| | External display is defective | Check external display and video cable on another monitor or PC |
| | iPC defective | Replace iPC |

External Touch Display not functioning

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|------------------------------|---|---|
| Touch Screen not functioning | Touchscreen functionality has been temporarily disabled | Check if touchscreen functionality has been temporarily disabled (padlock symbol on Main Screen key). If yes, press and hold the Main Screen key to re-enable touchscreen operation. |
| | Touchscreen functionality has been permanently disabled | In service mode, select Main Setup -> User Interface and change the "Touch Enable" selection to "yes". |
| | External Touch cable not connected | Check cable connection from external touch to connector board |
| | External Touch driver configuration | Check connector configuration: 1. Enter Main Setup menu 2. Select Monitor 3. Select Hardware 4. Reconfigure RS232/MIB drivers 5. if problem persists, proceed to the next step |
| | Connector board defective | Replace connector board |
| | External touch defective | Replace external touch |
| | Main board defective | Replace Main board |
| Touch position invalid | Touch not calibrated | Perform touch calibration: 1. Enter Main Setup menu 2. Select Hardware 3. Select Touch Calibration |

General Monitor INOP Messages

| INOP Message | Possible Causes of Failure | Failure Isolation and Remedy |
|--|--|---|
| Checkinternvoltage Check Monitor Func | Problem with too low voltages (5V, 12V) in the monitor. Alarm lamps, display or interfaces may not function correctly. | Remove all I/O boards, connector board, and iPC, and put them back in and reconnect cables one at a time to isolate any defective board. If this does not resolve the problem, replace the main board |
| Check Monitor Temp | The temperature inside the monitor is too high | Check the environment for possible causes |
| | Monitor ventilation obstructed | Clean the monitor ventilation internally and then cool monitor down for 8 hours |
| | Main Board defective | replace Main Board |
| Settings Malfunction | Problem during cloning process. | Reclone configuration file |
| | Memory space in which the settings are stored has been corrupted | Reclone configuration file. This will reload the memory space. |
| | Main board defective | Replace Main board |
| Internal Comm.Malf. | Problem with the I2C Bus communication in the monitor | Disconnect the external display and try another one |
| | Connector board defective | Replace connector board |
| | Main board defective | Replace Main board |
| MCC Unsupported | An MSL coupling cable has been connected to a device which does not support MSL coupling. | Use the MSL coupling cable only when connecting Dual CPU MP90 monitors to a D80 Intelligent display. |

Remote Control (wired)

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|--|--|
| Remote Control not recognized by monitor | Connector Board defective | Replace connector board |
| | Remote Control not plugged according to USB connection rules | See "Connection of USB Devices" in this service guide for the correct connection of USB devices. |
| | Remote Control defective | Exchange Remote Control |

Remote Control (wireless)

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|--|--|
| Remote Control not functioning | Low battery | Exchange battery |
| | Remote Control defective | Exchange Remote Control |
| Remote Control cannot be assigned to a monitor | Another remote control is assigned to the monitor. | Remove any other remote control before assigning a new remote control to a monitor |

Navigation Point

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|--|---|
| Navigation Point not functioning | Navigation Point not connected properly | Check cabling |
| | Navigation Point defective | Replace Navigation Point |
| Navigation Point Knob Rotation or Selection control not functioning | Navigation Point defective | Replace Navigation Point |
| INOP Message Check Input Devices is issued | Navigation Point or other input device defective | Perform a visual and functional check of all the monitor input devices. Replace input devices if necessary. |

iPC

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|--|---|
| iPC menu not available | iPC is not connected correctly | Check all cable connections to the iPC |
| Prompt Message "PC HW Malfunction" is issued | Battery empty | Exchange the battery of the iPC |
| | PC temperature is too high | Let the PC cool off and then try again |
| | PC defective | Exchange iPC |
| PC does not start automatically and/or cannot be started manually | HW setting incorrect | Check Global Settings and set to Autostart PC if Auto Start is desired. |
| | Flat ribbon cable to iPC not connected correctly | Check flat ribbon cable connection to iPC |
| | iPC defective | Exchange iPC |
| PC Operating System does not start up | Cable connection to hard disk loose | Check cable connections to hard disk |
| | Hard disk defective | Exchange hard disk |





| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|---|---|
| Operating System crashes during startup | Windows Installation malfunction | Perform a windows recovery by pressing F8 while the Philips Logo appears on the monitor screen. This will open the boot menu. Use regular windows recovery procedures |
| No sounds from PC | | Verify that the patient monitor sounds are functioning correctly |
| | PC sounds disabled | Check Global Settings and set to PC Audio On |
| | Flat ribbon cable to iPC not connected correctly | Check flat ribbon cable connection to iPC |
| | No sounds issued from PC | Check windows control panel settings. If you are using Windows XP, install the correct Audio drivers. Make sure the PC volume is not set to "0". |
| Undesired PC sounds from internal speaker | Incorrect HW setting | Set PC Audio to "Off" |
| DVI, VGA or USB connections not functioning | iPC defective | Exchange iPC |
| PC is unstable or shows colors incorrectly | RAM faulty | Exchange RAM |
| LAN connection not functioning | Flat ribbon cable to iPC not connected correctly | Check flat ribbon cable connection to iPC |
| | iPC defective | Exchange iPC |
| Date and Time not displaying correctly | Windows settings not synchronized with patient monitor. | Check windows settings |
| Display Resolution incorrect | Incorrect settings | Change the display resolution settings to match the patient monitor display resolution |







Keyboard/Mouse not functioning

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|---------------------------------------|------------------------------|
| Keyboard/Mouse attached directly to the monitor not functioning | Keyboard/Mouse not connected properly | Check cabling |
| | Keyboard/Mouse defective | Replace Keyboard/Mouse |
| | Connector board defective | replace connector board |
| Keyboard/Mouse connected to the iPC not functioning | Keyboard/Mouse not connected properly | Check cabling |
| | Keyboard/Mouse defective | Replace Keyboard/Mouse |
| | iPC defective | replace iPC |

Bedside Network Status Icons

The following table shows the icons displayed on the monitor when network related issues occur.

| Wireless Icon | Wired Icon | Inverse Video | Blinks | Icon Comments | INOP Message | What does it mean? |
|---|---|---------------|--------|------------------------|------------------------------------|---|
| No Icon | No Icon | - | - | - | - | MONITOR does not have a LAN connection (Wireless Monitor cannot find an access point to talk to, wired Monitor cannot hear anything on its LAN connection) |
|  |  | Yes | Yes | Central - outline only | "Unsupported LAN" (after 1 minute) | MONITOR has a LAN connection but does not have an IP address assignment (Wireless MONITOR has found an access point to talk to, wired MONITOR hears traffic on the LAN) |
|  |  | No | No | Central - outline only | "No Central Monitoring" | MONITOR is connected to the LAN and has an IP address assignment, but the bed is not being monitored at the central 1. MONITOR is not assigned to a sector 2. There is another monitor on the network with the same "Equipment Label" |

| Wireless Icon | Wired Icon | Inverse Video | Blinks | Icon Comments | INOP Message | What does it mean? |
|---|---|---------------|--------|---|-------------------------|--|
|  |  | No | No | Central - solid box | - | Normal Operation - MONITOR assigned to a sector and is being monitored by a central |
| - |  | No | No | Central - solid box, network line extended | - | Normal Operation MONITOR assigned to a sector and is being monitored by a central. This monitor also has OVERVIEW functionality on other beds. |
|  | - | No | Yes | Central - solid box | "Wireless Out Of Range" | Wireless MONITOR that currently is being monitored by a central is losing contact with the access point and cannot find another to talk to. |
|  |  | Yes | Yes | Central - outline only, line for broken connection to central | "No Central Monitoring" | Monitor lost connection to the Information Center: <ol style="list-style-type: none"> 1. LAN cable was disconnected 2. Information Center was disconnected 3. Network infrastructure failure (switch, etc.) 4. Out of range (wireless MONITOR) |

Network related problems

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|---|--|
| Prompt Message “no central assigned to this bed” is issued | The monitor label is not set in the monitor (if the beds are “monitor labeled” in the Philips Information Center) | Set Monitor Label in Config Mode |
| | Problem with the Philips Information Center to Switch communication (if the beds are “port mapped” in the Philips Information center) | Check PIC to Switch communication, Switch configuration and Firmware status |
| INOP “Unsupported LAN” is issued | Network failure | Check if switches, Philips Information Center and Database Server are all running and connected to the network |
| | Monitor connected to wrong network | Check if monitor has been connected for example to a different hospital network instead of the Philips Clinical Network |
| | IP address conflict after infrastructure re-installation | Reboot Database Server and Philips Information Center |
| | IIT is enabled but no IIT infrastructure can be found | Move the monitor into the range of the IIT infrastructure or disable IIT in the Setup IIT menu if no IIT infrastructure is available. |
| No connectivity to PIC, no prompt or error message on monitor | Hardware Defect | Check LAN cable connection |
| | | Check NGN Connector board in Monitor |
| | | Check Switch |
| | Configuration problem | Check switch configuration and firmware revision |
| Status Message “Incompatible SW Revision versions” is issued | Monitor and PIC software are not compatible | Check Software compatibility and upgrade to compatible software |

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--------------------------------------|--|--|
| Other Bed Overview not available | Configuration Problem | Check configuration in PIC regarding other bed overview (care group assignment) |
| | | Verify configuration of switch (setting of multicast filters) |
| | This function is not available for IntelliVue Instrument Telemetry and, in combination with earlier IntelliVue Information Center (IIC) revisions, for WLAN (IntelliVue 802.11 Bedside Adapter). | If you are using an IntelliVue 802.11 Bedside Adapter, check the software revision of the IntelliVue Information Center (IIC) to make sure it is compatible. If the software revision of the IIC is incompatible or you are using IIT, switch to a wired configuration |
| “Other Bed” Alarms are not appearing | Configuration problem | Verify configuration in PIC, in Monitor (Config Mode) and check that the feature is not temporarily disabled by the user (Bed Info Window) |

IntelliVue 802.11 Bedside Adapter Problems

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|--|---|---|
| No Network icon or Network icon flashes. No association to central station. | Communication problem between the monitor and the IntelliVue 802.11 Bedside Adapter or RSSI value below 30. | <p>Ensure that the network infrastructure is functioning properly. See Troubleshooting tables in the IntelliVue 802.11 a/g Infrastructure Installation and Configuration Guide for details.</p> <p>Check the antenna cable connection on the IntelliVue 802.11 Bedside Adapter.</p> <p>Check that the IntelliVue 802.11 Bedside Adapter is correctly connected to the panel adapter board.</p> <p>Check that the indicator behind the RSSI value (Main Setup -> Network -> WLAN Diagnostic -> RSSI) is rotating. If it is not, check IntelliVue 802.11 Bedside Adapter hardware.</p> <p>Replace antenna or IntelliVue 802.11 Bedside Adapter if necessary.</p> |
| | IntelliVue 802.11 Bedside Adapter not yet operational | <p>Check menu line 'Wireless LAN' (Main Setup -> Network -> WLAN Diagnostic -> Wireless LAN). If it shows 'Off', the wireless adapter is not yet operational. This does NOT indicate that WLAN has been disabled by a setting.</p> <p>If problem persists, check for an installed wired LAN cable.</p> |
| | Configuration problem using WEP, WPA(PSK), WPA2(PSK). | <p>Make sure that the Mode, SSID, Country and Security settings in the Setup WLAN menu match your installation</p> |
| | Configuration problem using WPA Enterprise or WPA2 Enterprise | <p>1. Check the connection status.(Main Setup -> Network -> WLAN Diagnostics -> Conn.Status)</p> <p>If the state only shows 'Scanning', make sure that the Mode, SSID, Country and Security settings in the Setup WLAN menu match your installation.</p> |

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|----------|------------------|---|
| | | <p>2. Check the connection status. If the device shows the state 'Authenticating', your SSID, Mode, Country and Security settings are correct.</p> <p>You already have a WLAN connection to your Access Point, but the device fails to authenticate, check your authentication server and WLAN controller error log.</p> <hr/> <p>3. As an investigation step, disable the CertificateCheck.</p> <p>(Main Menu -> Network -> WLAN Setup -> CertificateCheck)</p> <p>If authentication is now possible, proceed with step 4.</p> <p>Otherwise double check your authentication server configuration, WLAN controller configuration and the user credentials (User Name, Password, Anonymous Identity).</p> <p>Note:</p> <p>If the previously used credential settings were wrong, the device is perhaps on the exclude list of your WLAN Controller. Resolve this issue on your WLAN controller.</p> <p>Note:</p> <p>Do not forget to re-enable the certificate check.</p> |

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|----------|------------------|--|
| | | 4. Check the time setting of the device (Main Setup -> Date, Time). If not correctly set, the used certificates are detected as invalid. Adjust to the correct time. |
| | | 5. Check the installed CA certificate using the support tool. - Task -> Clone from Medical Device - Open the cloned file using Configuration -> Configuration Editor - In Configuration Editor check Configuration -> Hardware -> Network -> Certificate 1 for validity(Valid from, Valid until) |
| | | 6. Make sure that the installed CA certificate is the root certificate of your authentication server certificate chain. |

IIT-related Problems

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|---|--|--|
| No Network icon or Network icon flashes. No association to central station. | Communication problem between the monitor and the IIT adapter. MAC Instr. Tele. field in Instrument Telemetry Service Window is 0000 0000 0000 | Check that RF Access Code is set correctly and the network is correctly set up. Check the cable connection to the IIT module. Check the antenna cable connection between the IIT module and the antenna. Replace cable, antenna or IIT module if necessary. |
| | Incorrect RF Access Code. No IP Address. | Check that RF Access Code is set correctly. Make sure that network is set up correctly. |

Short Range Radio Interface Problems

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|---|---|---|
| Measurement selection icon does not change to SRR. | Assignment of SRR device to monitor not possible | Check SRR Configuration Settings. Replace defective SRR interface or cable, if necessary. Make sure SRR interface is installed. |
| | SRR interface of telemetry transceiver defective or incompatible | Make sure the telemetry transceiver SRR interface is compatible and functional. |
| Measurement selection icon changes to SRR but Assignment of SRR device to monitor fails. SRR Interference INOP is issued | RF Interferences | Check location for RF interferences and free frequencies by performing a site survey (e.g. with air magnet tool). |
| Communication Dropouts or gaps in parameter waves. SRR Interference INOP may be issued | RF Interferences | Check location for RF interferences and free frequencies by performing a site survey (e.g. with air magnet tool). |
| | Too many SRR devices allocated to one SRR channel | Up to two SRR connections can be established per channel. Check SRR Configuration Settings. |
| SRR communication aborted. SRR Interference or SRR Invalid Chan INOP may be issued. | RF Interferences | Check location for RF interferences and free frequencies by performing a site survey (e.g. with air magnet tool). |
| | Too many SRR devices allocated to one SRR channel | Up to two SRR connections can be established per channel. Check SRR Configuration Settings. |
| | SRR device out of range (either monitor or Telemetry Transceiver) | Position the SRR devices closer to each other. Check SRR signal quality indicator for signal strength. |
| Telemetry Device using SRR not recognized by the monitor. | Telemetry Device not supported by the SRR adapter | Make sure you use a telemetry device which is compatible with SRR. |

Multi-Measurement Module

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|--|---|
| Prompt message “Measurement Server Configuration not supported” is issued | An unsupported MMS Extension has been connected | Disconnect the MMS Extension |
| | MMS Extension is defective | Replace MMS Extension |
| | Measurement Server defective | Replace Measurement Server |
| INOP Message “MsmtSrv not Supp” is issued | Wrong Software Revision | Upgrade monitor and/or measurement server to a matching software version. Refer to for a list of compatible measurement servers. |
| | Too many measurement servers connected | Disconnect unsupported measurement servers for proper operation. |
| | Unsupported type of measurement server (for example M3000A) connected. | Disconnect the unsupported measurement server. Refer to for a list of compatible measurement servers. |
| Prompt message “Measurement Server not supported, unplug device, switch monitor off/on” and INOP “Bad Measurement Server are issued | M3000A Measurement Server Revision A is plugged. This Measurement Server is not compatible with the IntelliVue patient monitors. Parameter board defective. | Disconnect the measurement server and cycle power. Check if all measurements are displayed in the measurement selection window. Exchange MMS/Repair parameter board, if necessary. |

MSL-related problems

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|--|--|
| Measurement Server does not start up (no LEDs active), no INOP or prompt displayed | No Power | Check MSL cable and replace if necessary |
| | | Replace the Connector Board |
| Measurement Server does not start but LEDs are normal | Communication lines in MSL cable or MSL connector broken | Check MSL cable and MSL connectors |
| | Connector board or main board defective | Check connector board and replace if necessary. If problem persists, replace main board. |

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|--|---|
| MSL Power High INOP is issued Note: if this condition persists for longer than 15 minutes, the INOP MSL Power Off will appear (see below) | Attached devices drawing too much power from the monitor. Too many FMS and MMS connected to the monitor | Reduce to a limit of 2 FMS and 1 MMS connected to the monitor |
| MSL Power Off INOP is issued | Attached devices drawing too much power from the monitor. Too many FMS and MMS connected to the monitor | Disconnect all FMS and MMS from the monitor Cycle power to restore power to the MSL devices. If the message disappears, reconnect FMS and MMS one at a time, waiting 15 minutes between each device to see if message reoccurs. If yes, the respective MMS or FMS is faulty. See “Multi-Measurement Module” on page 130 for troubleshooting tasks. If no, add front-end modules one at a time, waiting 15 minutes between each module to see if message reappears, Replace module if faulty. Note: If an individual defective device is connected the MSL Power High or MSL Power Overload INOPs will appear initially. The MSL Power Off INOP will not occur for at least 15 minutes. |
| MSL Power Overload INOP is issued | Short Circuit within MSL system | Disconnect all MSL connections and reconnect devices one at a time. If message persists, replace connector board. |
| INOP Bad Server Link is issued | Unexpected data detected on MSL | Check cable and power cycle the monitor |
| | An FMS or MMS with an incompatible software revision is connected to the monitor. | Connect FMS or MMS with compatible software revision |
| | Communication between the components not functioning | Check software versions and model number of devices for compatibility |

4 Troubleshooting

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|---|--|
| INOP Message Serverlink Malf is displayed, audible indicator: a beep every two seconds | The hardware for communicating with the Multi-Measurement Server is faulty. | Check MSL cable, replace if necessary. |
| | | Check connector board. Replace if necessary |
| A measurement supported by a server does not come up on the monitor | Label conflict | A parameter label from this measurement is already in use in the monitor. Check the conflict window to select the measurement. |
| Prompt message “Too many <label> modules connected” is issued | There are more modules of the type <label> connected than supported by the software | Remove the unsupported module or use the label manager application in the monitor to disable the module. |
| The ECG Out function does not function | Hardware problem | Check MSL cable |
| | | Check ECG Out Hardware in the monitor (Power Switch Board) |
| | | Check the MSL connector in the measurement server |

Alarm Lamps

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|--|------------------------------------|--|
| INOP Message Check Alarm Lamps is issued | Alarm LED board cable disconnected | reconnect Alarm LED board to mainboard |
| | Alarm LED board defective | replace Alarm LED board |
| | Main board defective | replace Main board |
| Alarm occurs, but no LED lights up | Environmental lighting too bright | Place monitor in a darker environment |
| | Alarm LED board cable disconnected | reconnect Alarm LED board to mainboard |
| | Alarm LED board defective | Replace Alarm LED board |
| | Main Board defective | Main board |

Alarm Tones

| Symptoms | Possible Causes of Failure | Failure Isolation and Remedy |
|---|---|---|
| INOP Message Speaker Malfunction is displayed | Speaker cable disconnected | Reconnect speaker cable |
| | Speaker defective | Replace speaker |
| | Sound amplifier on main board defective | Replace main board |
| Alarm occurs but no alarm sound is issued | Audible alarm indicators have been switched off | Switch audible alarm indicators back on |
| | Volume set to 0 | Increase volume |
| | Speaker defective | Replace speaker |
| | Sound amplifier on main board defective | Replace main board |
| Alarm occurs on device connected to VueLink but no alarm sound is issued on the monitor | Configuration of VueLink is incorrect | Check VueLink configuration |

Alarm Behavior

If your monitor did not alarm in the way in which the end user expected, please consult the Instructions for Use for possible setup issues or configuration settings which could affect alarm behavior.

Individual Parameter INOPS

If any of the following parameter INOP messages are issued try the respective parameter in another device. If the INOP message persists replace the parameter module, the MMS or other indicated device.

- **CO2 Equip Malf**
- **ECG Equip Malf**
- **NBP Equip Malf**
- **<Pressure Label> Equip Malf**
- **RESP Equip Malf**
- **SpO2 Equip Malf**
- **SpO2 Transduc Malf**
- **SvO2 Equip Malf**
- **tcpO2 (or tcpCO2) Equip Malf**
- **<Temp Label> Equip Malf**
- **VueLnk Equ. Malf**
- **BIS Equip Malfunc**
- **BISx Malfunction**

- NMT Equip Malfunc
- NMT Incompatible
- NMT Cal Failed

Flexible Module Rack

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|---|---|---|
| Prompt Message “Unrecognized Measurement Module in slot s” is issued | An unsupported module has been plugged into the Flexible Module Rack | Unplug the unsupported module |
| Prompt message “Measurement Module in slot n is currently ignored” is issued | Too many modules of the same kind have been plugged into the Flexible Module Rack | Unplug module in slot n |
| Red Error LED stays on | Unrecoverable hardware selftest error: | Try to attach the MMS directly to the MSL cable. If the measurements show up on the screen, the fault is in the FMS |
| | MSL cable defective | If the measurements do not show up when the MMS is connected directly to the MSL cable, then replace MSL cable |
| | Flex connector from main board to MSL defective | Replace connector |
| | CPU module defective | Replace CPU module |
| Red Error LED flashes | Hardware selftest error | If system comes up, check status log. Otherwise see above |
| Flexible Module Rack LEDs ok, Front End Measurement Module not recognized (no prompt or INOP) | Measurement Module or Measurement Module Connector defective | Replace Measurement Module |
| | No Front-End power because MSL voltage from the monitor is too high or too low | Try a new MSL cable. Replace if failure is rectified. Replace mother board |
| | Mother board or connector on Mother Board defective | If the voltage is in range, or there is obvious damage to a connector, replace mother board |

Printer

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|---|---|--|
| Prompt message “Print job could not be queued” is issued. No print device is found. | Printer is disabled in the Setup Printers menu Paper size of printer does not match paper size of report | Enable the correct printer in the Setup Printers menu Change paper size of the printer in the Setup Printers menu or change paper size of the report in the Setup Reports menu. |
| Status message “Print device Local 1 (Local 2) unavailable” is issued. Printer job is stalled. | Printer not switched on Printer paper tray empty Cabling not connected correctly Connector board defective | Switch on printer power fill printer paper tray Check cabling replace connector board |
| Status message “Print device Remote 1 (Remote 2, Remote 3) unavailable” is issued. Printer job is stalled | Print error on Philips Information Center Network Connection to Philips Information Center not functioning | Print a test report on the Philips Information center. If this fails, refer to Philips Information Center documentation Check that the network connection between the monitor and the Philips Information Center is working |
| Status message “Printing on device Remote 1... (Remote 2, Remote 3)” is issued but no report is printed | Print queue on Philips Information Center is full. Reasons for this may be: - Printer is not switched on - Printer paper tray is empty | Switch on printer power Fill printer paper tray |
| Printouts are not as expected | Printer paper size is not correctly configured Printer resolution is not correctly configured Printer color support is configured to “On” although the printer does not support color Printer not compatible | Configure the paper size according to the inserted print media Configure the printer resolution according to the printer capabilities Configure the printer color support to “Off” Check specifications |

Recorder

| Symptom | Possible Cause | Corrective Action |
|---|---|---|
| System thinks that door is open when it is not. | Defective door switch. | Replace door switch. Exchange module. |
| System thinks that the recorder is out of paper when it is not. | Paper-out sensor dirty. | Clean paper-out sensor. |
| Recorder not communicating with System. | Poor connection to the front-end FMS. | Unplug the module. Plug it back in and try it again in a few seconds. (Watch for the LED to flash.) |
| | Only one recorder module may be used with each monitor. | Remove one of the recorder modules. |
| | System not configured properly. | Check the configuration of the connected monitor. |
| | Too many modules connected. | Check and remove the extra modules. |
| Recorder won't run. | Recorder interface not working correctly. | Unplug the module. Plug it back in and try it again in a few seconds. (Watch for the LED to flash.) |
| Poor print quality. | Printhead dirty. | Clean the Printhead. |
| | Printhead failure. | Exchange the module. |
| Paper not feeding properly. | Paper roll off center. | Center paper roll on roller guides. |
| | Dirty roller. | Clean roller. |
| Module does not lock into FMS. | Locking plates defective. | Remove and exchange the locking plates. |

MIB / RS232

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|--|--|--|
| AGM connected to an RS232 port not functioning | The MIB/RS232 port is not configured for AGM | Check configuration of the MIB/RS232 ports in configuration mode |
| | The cable between AGM and monitor is not connected correctly or defective | Check cable connection, replace cable if necessary |
| | The MIB/RS232 board is in a wrong slot (slot has been changed after software configuration or an additional board has been plugged in) | Verify correct placement of the I/O boards |
| | The MIB/RS232 board or the connector board (depending on which RS232 port is used) is defective | Check board and replace if necessary |

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|------------------------------------|--|--|
| External device not receiving data | The MIB/RS232 port is not configured for data export | Check configuration of the MIB/RS232 ports in configuration mode |
| | The wrong data export protocol driver is configured in the monitor | Check the export protocol required by the attached device and configure the monitor accordingly |
| | The cable between the external device and the monitor is not connected correctly or defective | Check cable and replace if necessary |
| | The external device does not support the version of the data export protocol used in the monitor | Check if the device supports the version of the data export protocol. Upgrade device or monitor if necessary (if matching versions exist). |
| | A terminal concentrator is used in between the device and the monitor and a protocol with dynamic speed negotiation is used | Some terminal concentrators do not support changing the transmission speed (baud rate) dynamically. Check if the connection works without the concentrator |
| | The MIB/RS232 board is in a wrong slot (slot has been changed after software configuration or an additional board has been plugged in) | Verify correct placement of the I/O boards |
| | The MIB/RS232 board or the connector board (depending on which RS232 port is used) is defective | Check board and replace if necessary |
| Detailed Protocol Problem | | Consult the Data Export Protocol document. |

USB

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|--|--|--|
| None of the connected devices are functioning. | Invalid combination of connected devices or connected devices defective. | Make sure the combination of connected devices is valid. Replace defective devices if necessary. |

Nurse Call Relay

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|---|---|---|
| INOP message Check Nurse Call Relay is issued | Connector board defective | Replace connector board |
| Monitor alarmed, Nurse Call did not activate | Incorrect configuration (Relay latency, Relay trigger) | Check monitor configuration (see configuration guide) |
| | Connection of cable to monitor or nurse call system not correct | Check cable connection |
| | Connector board is defective | Replace connector board |

Flexible Nurse Call Relay

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|---|---|---|
| INOP message Check Nurse Call Relay is issued | Flexible Nurse Call Relay board defective | Replace Flexible Nurse Call Relay I/O board. |
| Monitor alarmed, Nurse Call did not activate | Incorrect configuration (Relay latency, Relay trigger) | Check monitor configuration (see configuration guide) |
| | Connection of cable to monitor or nurse call system not correct | Check cable connections |
| | The Flexible Nurse Call Relay board is defective | Replace Flexible Nurse Call Relay board |

Troubleshooting the ECG OUT

| Symptoms | Cause of Failure | Failure Isolation and Remedy |
|----------------------------------|--|--|
| INOP EcgOut Equip Malf is issued | Communication Problem or Power Switch/ECG OUT board defective. | Check that the ECG OUT cable is securely connected and that all MSL connections are properly made. Check that the MSL cable and the MSL connectors are not defective. If the problem persists, replace the Power Switch/ECG OUT board. |
| No ECG-OUT signal to the Defib | | Check the Defib cable and the cable connection from the Power Switch/ECG Out board to the main board. Exchange Power Switch/ECG Out Board if necessary. If problem persists exchange main board. |


Image Sticking

If a static image is displayed for a long time on an LCD display, image sticking, i.e. a temporarily retained image, may occur. To eliminate image sticking, switch off the display and switch it back on again. It is also recommended to use the moving image in standby mode.

Testing the Functionality of the Power Supply

If you are unsure whether the power supply or the main board is defective, perform the following test, to check whether the power supply is faulty.

Place one multimeter probe on the power supply housing and the other at the test point of the fuse (see picture below). If the measured voltage is 24V ($\pm 5\%$) the power supply is functioning correctly.



WARNING

High Voltage - Voltages dangerous to life are present in the instrument when it is connected to the mains power supply.

Status Log

Many events that occur during start-up or regular monitoring are logged in the Status Log. The Status Log can be printed and cleared. Not all entries in the Status Log are errors.

| Monitor | | | | |
|---------|------|-------|---|----------------|
| H | 1720 | 20050 | 1 | 4 Apr 02 16:37 |
| C | 1721 | 21050 | 1 | 4 Apr 02 15:37 |

The window title is either **Monitor** or **MeasServ**, dependent on which system component's status log is currently displayed.

The Status Log window shows logged events which caused a reboot of the system component (monitor or measurement server).

The first column in the log identifies the event class ("C": caused a cold start, "H": caused a hot start, "N": no restart, for information only). Column 3 and 4 identify the event source and event code. Column 4 counts the number of occurrences of the event. The last column shows the time and date of the last occurrence of the event.

The following pop-up keys overlay the SmartKeys:

| | | | | | |
|--------------------------|-----------------|--|---------------|---------------|---------------|
| Clear StatLog | Revision | | 865241 | M8048A | M3001A |
|--------------------------|-----------------|--|---------------|---------------|---------------|

Clear StatLog

This key clears the currently displayed Status Log

Revision

This key switches to the Revision Screen of the currently displayed system component

865241 (MX700) or **865242** (MX600)

This key switches to the Monitor Revision Window

M8048A (FMS-8) or **865243** (FMS-4)

This key switches to the Flexible Module Rack (FMS) Revision Window

M3001A

NOTE

This key switches to the Multi Measurement Server (MMS) Revision Window

If an event occurs repeatedly, contact your Philips Service Representative.

NOTE

It is possible, using the support tool, to download the status log and send it to your Philips Service Representative as a file (for example via e-mail).

List of Error Codes

There are no error codes at this point.

Troubleshooting with the Support Tool

Using the support tool you can:

- access the full status log which can be saved as a file
- reload software
- identify defective devices
- reset touch screen calibration

For details on how to perform these tasks see the Support Tool User Manual.

Troubleshooting the Individual Measurements or Applications

For problems isolated to an individual parameter or application such as event review, please consult the Instructions for Use and configuration information.

If the instructions for use did not resolve an individual parameter problem, then another module or measurement server should be tried.

If you are getting questionable readings for individual measurements you may want to do the Performance Verification tests in the *Testing and Maintenance* section.

4 Troubleshooting

The performance of the individual applications (event review, arrhythmia, trending) are affected by the configuration of the monitor. When contacting Philips support you may be asked about the configuration of the monitor to aid in troubleshooting.

Repair and Disassembly

The following section describes the disassembly and reassembly procedures for the monitor and its components to the extent required to remove and replace faulty assemblies. Do not further disassemble the product past the point described in these procedures.

WARNING

High Voltage - Voltages dangerous to life are present in the instrument. Do not perform any disassembly or reassembly procedures (other than MMS, MMS extension or parameter module removal) with power applied to the instrument. Failure to adhere to this warning could cause serious injury or death.

Before doing any disassembly, turn power off, disconnect the Local Distribution Cable, **AC power cable**, MSL cable, Defib sync' cable and RS232 cable (where appropriate), disconnect the MMS and FMS.

Tools Required

- Torx screwdrivers (T10, T20)
- 1 small flat head screwdriver
- Needle Nose Pliers
- ESD mat and wrist strap
- 1 small Pozi or Philips head screwdriver (PH0 x 60)

Monitor Disassembly

NOTE


- The reassembly procedures are the reverse procedures of the disassembly procedures unless otherwise noted.
- Your monitor may look slightly different than on the pictures in this chapter, depending on the options ordered.
- When disassembling, position the monitor such that the weight of the monitor is not on the Navigation Point.

Removing the Cable Cover

- 1 Push in the release lever of the cable cover housing and remove the cable cover.





- 2 Remove the four screws and pull off the cable cover housing.




Removing the Bottom Housing

- 1 Remove the five screws and pull off the bottom housing.



Removing the Power Button

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145
- 2 Remove the two self-cutting screws inside the bottom housing and pull off the tappet guide.



- 3 Remove the power button.




Removing the Power Switch/ECG Sync Out Board


- 1 Remove the bottom housing as described in “Removing the Bottom Housing” on page 145
- 2 Unplug the flat ribbon cable from the Power Switch/ECG Sync Out Board.




- 3 Remove the three screws from the Power Switch/ECG Sync Out Board.



- 4 Remove the light guide from the Power Switch/ECG Sync Out Board.




- 5 Pull the Power Switch/ECG Sync Out board off its guidance pin and remove the board.



Removing the Quick Mount


- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the four screws from the quick mount.




- 3 Remove the quick mount.

Removing the Housing Rear

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Pull out the two white plastic pins.




- 3 Remove the housing rear.




Removing the I/O Cards

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as as described in the section “Removing the Housing Rear” on page 148.
- 3 Remove all I/O cards and blank I/O slot covers by pushing the release lever and pulling them out at the same time.



Removing the Connector Board


- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.
- 3 Remove the I/O cards as described in the section “Removing the I/O Cards” on page 148.
- 4 Loosen the two captive screws on the connector board assembly.



- 5 Remove the two screws holding the connector board.



- 6 Remove the connector board.



Removing the AC/DC Power Supply

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.

5 Repair and Disassembly

- 3 Unplug the AC/DC power supply cable by pressing the two latches on the side of the connector.



- 4 Remove the four screws.



- 5 Remove the AC/DC power supply.



Removing the iPC


This procedure only applies if the iPC is installed.

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.

- 3 Unplug the three iPC cables from the iPC main board (Video, USB, Power).



- 4 Remove the four screws from the iPC.



- 5 Remove the side USB cable from the iPC housing.




- 6 Remove the iPC.



Separating the Front and Back Half of the Monitor

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Unplug the flat ribbon cable from the Power Switch/ECG Sync Out board.




- 3 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.
- 4 Remove the I/O Boards as described in the section “Removing the I/O Cards” on page 148.
- 5 Remove the connector board as described in the section “Removing the Connector Board” on page 149.
- 6 Unplug the AC/DC power supply cable.




- 7 If installed, unplug the three iPC cables as described in the section “Removing the iPC” on page 150.


- 8 With the display facing downwards, remove the bottom row of screws (four screws).



- 9 Remove the four pins.




- 10 Remove the arched element of the housing containing the AC/DC power supply and, if installed, the iPC.



Removing the Rear Display Housing

- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Pull out the nine release clasps until you hear a click and the padlock symbol on the release clasp is visible.



- 3 Remove the rear display housing.




Reassembling the Rear Display Housing

- 1 Push the nine release clasps back in before repositioning the rear display housing.
- 2 Reinsert the rear display housing and press it down carefully.


Removing the Loudspeaker

- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.

- 3 Unplug the loudspeaker cable..




- 4 Remove the three screws from the loudspeaker.




- 5 Remove the loudspeaker.

Removing the Alarm LED Board


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Unplug the alarm LED cable.



- 4 Remove the two self-cutting screws.



- 5 Remove the alarm LED board.



Removing the Light Pipes


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Remove the Alarm LED board as described in the section “Removing the Alarm LED Board” on page 155.
- 4 Remove the light pipes.




Removing the Video Mixer Board

This procedure only applies if an iPC is installed.


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the Rear Display Housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Unplug the three cables from the video mixer board.



- 4 Remove the four screws from the video mixer board.




- 5 Remove the video mixer board.




Removing the SRR Board

- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 If an iPC is installed, unscrew the screws on the video mixer board, unplug the video mixer board from the main board and lift up the video mixer board from the housing, leaving all other cable connections to the video mixer board plugged. Be careful not to damage the cable connectors on the video mixer board.
- 4 Unplug the SRR cable from the main board.



- 5 Release and remove the SRR board by pushing the latch in the monitor housing outwards.




Removing the WLAN/IIT Antenna

This procedure only applies if WLAN or IIT are installed.

- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Unplug the WLAN/IIT antenna from the WLAN board or IIT module.


- 4 Remove the antenna by pushing the latch in the housing outwards as shown below.




Removing the WLAN board

This procedure only applies if WLAN is installed.


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Remove the WLAN antenna by pushing the latch in the housing outwards as shown below.



- 4 Unplug the WLAN board including its holder from the panel adapter board.




- 5 Reassemble the WLAN Board by performing the above steps in reverse order. When reinserting the WLAN board, make sure to place the WLAN antenna cable above the WLAN board.




Removing the IIT Module

This procedure only applies if IIT is installed.


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Remove the IIT/WLAN antenna incl. cable as described in the section “Removing the WLAN/IIT Antenna” on page 158.
- 4 Remove the IIT cable from the main board.



- 5 Remove the two screws from the IIT module.




- 6 Remove the IIT module by pulling it towards you as shown below.




Removing the Panel Adapter Board


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 If installed, remove the WLAN board and WLAN board holder as described in the section “Removing the WLAN board” on page 159.
- 4 Unplug the flat ribbon cable from the panel adapter board.




- 5 Unplug the touch cable from the panel adapter board by pressing the connector latch down and pulling out the cable. Note that non-touch monitors will not have this cable.



- 6 Remove the two screws from the panel adapter board.




- 7 Remove the panel adapter board.




Removing the Backlight Inverter Board

- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 Unplug the cable between the panel adapter board and the backlight inverter board.




- 4 If IIT is installed, remove the IIT antenna as described in the section “Removing the WLAN/IIT Antenna” on page 158.


- 5 Unplug the four backlight cables.



- 6 Remove the four flat-head screws from the backlight inverter board.



- 7 Remove the backlight inverter board.
- 8 Reassemble the backlight inverter board by performing the above steps in reverse order. Make sure that the backlight cables are inserted into the cable holder as shown below.




Removing the Main Board


- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 If an iPC is installed, remove the video mixer board as described in the section “Removing the Video Mixer Board” on page 157.
- 4 If IIT is installed, unplug the IIT cable from the main board.

5 Repair and Disassembly

- 5 Unplug the multi-colored display cable, the loudspeaker cable, the SRR cable, the alarm led board cable and the HIF board cable from the main board.



- 6 Remove the panel adapter board as described in the section “Removing the Panel Adapter Board” on page 161.
- 7 If installed, remove the video mixer board incl. cables as described in the section “Removing the Video Mixer Board” on page 157.
- 8 Remove the 11 screws from the main board.




- 9 Remove the main board.


Removing the Touch Bezel Assembly

- 1 Separate the front and back half of the monitor as described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 If the IIT module is installed, remove the IIT module incl. antenna as described in the section “Removing the IIT Module” on page 160.
- 4 Unplug the touch cable from the panel adapter board and remove it from the guidance latches in the display chassis assembly.
- 5 Remove the Alarm LED board as described in the section “Removing the Alarm LED Board” on page 155.
- 6 Remove the three light pipes as described in the section “Removing the Light Pipes” on page 156.
- 7 Remove the SRR board as described in the section “Removing the SRR Board” on page 158.

- 8 If WLAN is installed, remove the WLAN board as described in the section “Removing the WLAN board” on page 159.
- 9 Unplug the HIF board connector from the main board.
- 10 Remove the ten self-cutting screws from the touch bezel assembly.




- 11 Remove the display chassis assembly by pulling it up out of the touch bezel assembly.



NOTE


Place the display chassis assembly on a flat, smooth surface after removal to avoid scratching the LCD panel.

- 12 Reassemble the touch bezel assembly by performing the above steps in reverse order. When reassembling, make sure that the antenna contact sheet is inside the bezel.




Removing the HIF Board, Navigation Point and Encoder


- 1 Remove the touch bezel assembly as described in the section “Removing the Touch Bezel Assembly” on page 164.
- 2 Remove the three screws from the HIF board.



- 3 Unplug both cables from the HIF board.



- 4 Push out the silicone pad from the back.



- 5 Pull off the Navigation Point using a screwdriver.




- 6 Remove the nut from the Navigation Point bolt using pliers.



- 7 Pull out the Navigation Point bolt from the bottom.



- 8 Remove the encoder as shown below.




- 9 Remove the serial number plate and the feature plate.




Removing the LCD Panel


- 1 Separate the front and back half of the monitors described in the section “Separating the Front and Back Half of the Monitor” on page 152.
- 2 Remove the rear display housing as described in the section “Removing the Rear Display Housing” on page 154.
- 3 If the IIT module is installed, remove the IIT module incl. antenna as described in the section “Removing the IIT Module” on page 160.
- 4 Unplug the touch cable from the panel adapter board.
- 5 Unplug the alarm LED board cable from the alarm LED board.
- 6 Unplug the SRR cable and remove the SRR board as described in the section “Removing the SRR Board” on page 158.
- 7 If WLAN is installed, unplug the WLAN antenna cable.
- 8 Unplug the backlight cables from the backlight inverter board.
- 9 Unplug the multi-colored display cable from the main board.
- 10 Unplug the cable between the backlight inverter board and panel adapter board from the backlight inverter board.
- 11 If an iPC is installed, remove the video mixer board as described in the section “Removing the Video Mixer Board” on page 157.
- 12 Remove the fourteen screws from the display chassis assembly.



- 13 Remove the display chassis assembly and remove the display.



- 14 Remove the four screws on the sides of the display chassis assembly (two screws on each side) and remove the angled metal sheets from the sides of the LCD panel.




Removing the Backlights

- 1 Remove the LCD Panel as described in the section “Removing the LCD Panel” on page 168.
- 2 Remove the screws from the individual backlights (PH0 x 60).





- 3 Pull out the backlights.



Removing the Blank Side Cover or USB Side Cover

If you have iPC installed you will have a USB side cover. If not, you will have a blank side cover.


- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.
- 3 Remove the side cover or the USB side cover from the bottom housing using a screwdriver. The covers are attached to a hook. Press the hook with the screwdriver and pull out the cover.



Removing the ECG Sync Side Plate

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.
- 3 Remove the Power Switch/ECG Sync Out Board as described in the section “Removing the Power Switch/ECG Sync Out Board” on page 146.



- 4 Remove the ECG Sync side plate from the bottom housing using a screwdriver. The ECG Sync side plate is attached to a hook. Press the hook with the screwdriver and pull out the side plate.




Removing the Side USB Cable

This procedure applies only if iPC is installed.

- 1 Remove the bottom housing as described in the section “Removing the Bottom Housing” on page 145.
- 2 Remove the housing rear as described in the section “Removing the Housing Rear” on page 148.
- 3 Remove the side cover or the USB side cover from the bottom housing as described in the section “Removing the Blank Side Cover or USB Side Cover” on page 170.
- 4 Remove the two cable holders by pulling on the cable.



- 5 Unplug the USB Cable.



NOTE




When reinserting the USB cable, make sure it is reassembled in exactly the same way it was before, otherwise the monitor cannot be reassembled.

iPC Disassembly

The procedures in this section only apply if an iPC is installed.



Removing the Battery

- 1 Remove the iPC from the monitor as described in the section .
- 2 Pull back the latch holding the battery and remove the battery.



Removing the iPC Main Board Assembly

- 1 Remove the iPC from the monitor as described in the section .
- 2 Remove the eight screws from the iPC (six on the side, two on the back).




- 3 Unplug the two cables coming from the hard disk and separate the main board assembly from the hard disk assembly.




Removing the RAM

- 1 Remove the iPC main board assembly as described in the section “Removing the iPC Main Board Assembly” on page 172.
- 2 Press the latches holding the RAM boards to the side.




- 3 Pull out the RAM boards carefully.




Removing the Hard Drive


- 1 Remove the iPC main board assembly as described in the section “Removing the iPC Main Board Assembly” on page 172.
- 2 Pull the two cables out of the metal latches.



- 3 Turn the hard drive assembly around so the hard drive is facing downwards. Then remove the four screws from the chassis.



- 4 Remove the hard drive.




8-Slot Flexible Module Rack (FMS-8) Disassembly



Removing the Handle and the Measurement Server Mount

Please note that any combination of handles and mounts is possible.

- 1 Remove the two screws on the bottom with a T20 screwdriver.




- 2 Slide the handle up and pull it out.
- 3 Remove the MMS if connected. Slide the MMS mount up and remove it by pulling directly perpendicular to the FMS.



NOTE

There is a connector located on the MMS mount. If you tilt the MMS mount as you remove it, you may damage the connector.





5 Repair and Disassembly



- 4 Remove the connector housings on each side of the FMS by compressing the cover slightly using two screwdrivers.





- 5 Remove the two white pins on each side with a small screwdriver.




- 6 Take off the rear housing.



- 7 Remove the four screws on the CPU board, pull it gently off the mother board, unplugging the connector at the same time.



- 8 Remove the four remaining screws on the mainboard.




- 9 Pull off side connector brackets by pulling them gently away from the housing on each side and lifting carefully.



5 Repair and Disassembly


- 10 Lift up and pull on the tabs to remove the connector holders.

Reassembly Note: The connector holders are side specific.



- 11 Unsnap the lightpipe.


Reassembly Note: You must snap the lightpipe bottom into place before inserting the top into the tab.



- 12 Pull off the main board.



Reassembly Note: Make sure that the rubber seal around the module connectors is inserted properly into the front housing.




Please make sure to set the exchange part data (serial number) with the support tool after reassembly. For details please refer to the support tool Instructions for Use.

4-Slot Flexible Module Rack (FMS-4) Disassembly (without MMS Mount)


Separating the Front and Back of the 4-Slot FMS

- 1 If the Universal Mounting Clamp is installed, remove the two screws (T20) from the Universal Mounting Clamp and remove the clamp.



If any other mounting bracket is installed, remove it.

- 2 Remove the three screws (T10) from the back of the FMS.




- 3 Push the front housing off with your hands.




5 Repair and Disassembly

- 4 Flip the back of the FMS down, resting the front of the FMS on the rear housing as shown below. Then unplug the two cables from the main board.




- 5 Reassemble the FMS by performing the above steps in reverse order. When reassembling the front and back of the FMS, make sure the cables are inserted into the main board as shown below.





Make sure that the screws in the rear cover of the FMS are tight enough so the FMS is properly sealed.

Removing the Main Board


- 1 If installed, remove the Universal Mounting Clamp.
- 2 Separate the front and back of the FMS as described in “Separating the Front and Back of the 4-Slot FMS” on page 181.
- 3 Remove the two screws (T10) from the main board.




- 4 Pull out the two connector holders with a flathead screwdriver.



- 5 Remove the main board.




- 6 Reassemble the FMS by performing the above steps in reverse order. Make sure that the light guide is in its correct position (see picture below) and reinsert the connector holders by pushing them in until you hear a click.





Removing the MSL Boards


- 1 Separate the front and back of the FMS as described in “Separating the Front and Back of the 4-Slot FMS” on page 181
- 2 Remove the two screws from the MSL holder.




- 3 Push the MSL Holder forwards to remove it.




- 4 Push the MSL UP board forwards to remove it.



- Remove the MSL DOWN board accordingly by performing the above steps for the other MSL board.



- Reassemble the MSL boards by performing the above steps in reverse order. When reassembling, push down MSL holder before tightening screws.



4-Slot Flexible Module Rack (FMS-4) Disassembly (with MMS Mount)


Removing the Universal Mounting Clamp

- If the Universal Mounting Clamp is installed, remove the two screws (T20) from the Universal Mounting Clamp and remove the clamp. If any other mounting bracket is installed, remove it.



Removing the MSL DOWN board in MMS Mount


- 1 Remove the three screws (T20) from the bottom cover.



- 2 Remove the bottom cover. If cable management hooks are installed, remove them.




- 3 Remove the two pins as shown below.




- 4 Position FMS with the MMS mount pointing upwards.


- 5 Remove the four silicone screw covers as shown below.



- 6 Remove the four screws (T10) from the MMS mount.




- 7 Flip FMS so it is positioned right side up.
- 8 Pull off the MMS mount.




5 Repair and Disassembly


- 9 Remove the two screws (T10) from the MSL holder.




- 10 Remove MSL holder as shown below.



- 11 Remove the MSL Down board.




- 12 Remove the MSL cable from the board.




- 13 Reassemble the FMS by performing the above steps in reverse order. When reassembling insert the silicone screw covers with the holes pointing towards the bottom.

Separating the Front and Back of the 4-Slot FMS with MMS Mount


- 1 Remove the Universal Mounting Clamp or any other mounting bracket as described in “Removing the Universal Mounting Clamp” on page 185.
- 2 Remove the MMS Mount as described in steps 1 to 11 of the section “Removing the MSL DOWN board in MMS Mount” on page 186
- 3 Remove the three screws (T10) from the back of the FMS.



- 4 Push the front housing off with your hands as shown below.




- 5 Unplug the cables from the main board.



5 Repair and Disassembly

- 6 Reassemble the FMS by performing the above steps in reverse order. When reassembling the front and back of the FMS, make sure the cables are inserted into the main board as shown below.





Make sure that the screws in the rear cover are tight enough so that the FMS is properly sealed and that the silicone screw covers in the MMS mount are inserted with the hole facing downwards.

Removing the Main Board


- 1 Separate the front and back of the FMS-4 as described in “Separating the Front and Back of the 4-Slot FMS with MMS Mount” on page 189.
- 2 Remove the two screws (T10) from the main board.




- 3 Pull out the two connector holders with a flathead screwdriver.



- 4 Remove the main board.




- 5 Reassemble the FMS by performing the above steps in reverse order. Make sure that the light guide is in its correct position (see picture below) and reinsert the connector holders by pushing them in until you hear a click..





Removing the MSL UP Board


- 1 Separate the front and back of the FMS as described in “Separating the Front and Back of the 4-Slot FMS with MMS Mount” on page 189.
- 2 Remove the two screws from the MSL holder.



- 3 Push the MSL Holder forwards to remove it.





- 4 Push the MSL board forwards to remove it.



- 5 Reassemble the MSL UP board by performing the above steps in reverse order. When reassembling, push down MSL holder before tightening screws.

Plug-in Modules

The snap lock holds the plug-in module in the FMS.



To remove the snap lock:


- 1 Grip the module firmly in one hand and using your thumb, pull the front edge of the snap lock away from the plug-in module so that the lug on the snap lock clears the retaining edge of the module.
- 2 Push on the rear edge of the snap lock to move the snap lock through the slot toward the front of the module until it is clear.

To replace the snap lock:

- 1 Locate the snap lock into the slot on the bottom of the module.
- 2 Slide the snap lock toward the rear of the module until the lock snaps into position.

Plug-In Module Disassembly

Disassembly of the parameter module enables replacement of the front assembly.



Removing the Module Front Housing

WARNING

When you disassemble/assemble a plug-in module an applied part leakage current test must be performed before it is used again for monitoring.

To disassemble a plug-in module:


Remove the front housing.

- Place the module on a flat surface and insert a card (similar to a credit or cheque type card) into one side of the module to disengage the 2 tabs securing the front housing to the module housing.
- Pull the edge of the front housing away from the module housing.
- Carefully turn the module over so the free edge does not reengage and repeat the first two steps on the other side of the module. The front housing should now be free of the module housing.

To reassemble a plug-in module:

Snap-fit the front housing onto the front of the module case so the openings in the front housing match the LEDs and keys.

tcpO2/tcpCO2 Calibration Chamber Kit



M1018A New Style Module

NOTE

You must order a new front housing AND a new calibration chamber kit when repairing a traditional CMS-Style M1018A module. The calibration chamber must be replaced first for the new style housing to fit properly

To remove the calibration chamber


- 1 Using a flat-tipped screwdriver, remove the screw holding the calibration chamber in place on the front of the plug-in module.
- 2 Lift the chamber off the plug-in module. Ensure that the white plastic switch tip located in the module is not lost.

To replace the calibration chamber

- 1 Ensure the white plastic switch tip is in place in the plug-in module.
- 2 Place the calibration chamber in the allocated position on the plug-in module.
- 3 Insert and tighten the screw into the calibration chamber, securing it to the plug-in module.

Recorder Module Paper

The recorder will not run when the door is open or when the recorder is out of paper. To prevent damage to the recorder module, use only Philips approved paper (Philips re-order number 40477A/B)



To load paper into the recorder module:

- 1 Remove the empty core from the previous roll of paper.
- 2 Cut off and discard the first few inches of paper to eliminate any traces of adhesive.
- 3 Pull out several inches of paper from the new roll, holding the roll with the loose end hanging over the top toward you.
- 4 Open the door and push the paper roll into the holders in the recorder.
- 5 Thread the paper under the roller and over the plastic shelf far enough so it goes around the roller and comes out above it.
- 6 Drape the paper over the end of the door and close the door. The paper should be visible and draped down in front of the door.

Multi-Measurement Module (MMS) Disassembly

Please follow the disassembly and reassembly steps below closely. Do not disassemble the MMS past the point described in the procedures below.

Tools required

- thin-bladed screwdriver
- ESD mat and wrist strap

WARNING

- Do not open the MMS while it is connected to a monitor.
 - Parts inside the instrument may be contaminated with bacteria. Protect yourself from possible infection by wearing examination gloves during this procedure.
-

Removing the Front Cover

- 1 Position the thin-bladed screwdriver in the small slot provided for this purpose. Remove the front cover by pulling it away from the MMS until it snaps off. There may be a slight resistance when removing the front cover.



Removing the Mounting Pin

- 1 Position the MMS with the connectors facing towards you. There are four long mounting pins threaded into the MMS in each of the four corners under the cover. Locate the heads of the two long mounting pins on the top cover and only remove these.
- 2 Use the thin-bladed screwdriver to lift the pins gently out, far enough that they can be removed manually.



- 3 Remove the two pins and set them aside for refitting.



NOTE

Without these long mounting pins, the MMS will not function properly

Removing the Top Cover

Begin by gently pulling the top cover away from the MMS. The top cover is press-latched at the MMS connector. There might be a resistance due to the rubber sealing. Remove the cover slowly, without hitting or touching the inside of the MMS.



Exchanging / Removing the DC/DC Board

NOTE


The HW Rev C MMS (S/N prefix DE610xxxxx) does not have a separate DC/DC board anymore.

The DC/DC board is connected to the main board. Loosen the pin connection to the main board and remove the DC/DC board by gently lifting it up. Avoid touching the surface of the board. Set it aside where it is ESD protected.




Removing the MSL Flex Assembly


- 1 After the DC/DC board is removed, lift up the MSL frame connector to which the MSL Flex is attached.



At the beginning there might be resistance due to the special fixing mechanism shown in the picture below.




- 2 Lift up the flex connector carefully. Do not bend the connector pins on the main board.




Reassembling the MSL Flex Assembly


- 1 Insert the MSL Flex layer into the frame connector as shown below by moving it into the appropriate dove tail.



- 2 To insert the MSL Flex into the MMS, it has to be bent carefully. Bend the MSL Flex in a 180 degree angle as shown below. Do not crease the flex.




- 3 The second bend has to be done as shown below. To be able to connect the MSL flex to the main board afterwards, the flex has to be bent in a 90 degree angle as shown in the picture. Do not crease the flex.




5 Repair and Disassembly

- 4 Insert the frame connector with the attached and bent MSL Flex. Be careful not to damage the MSL flex when pushing the frame connector downwards.



- 5 Position the connector correctly and push it into place.




Removing the NBP pump

- 1 Remove the pump by lifting it up. Set the pump aside. Also remove the old silicon tubes.




- 2 Remove the connector of the NBP pump assembly. The connector may sit tightly. Gently loosen the connector.




Removing the Keypad

- 1 Remove the keypad by lifting it straight up.




NOTE

The rubber sealing and the guides (marked with red circles) may hold the keypad firmly in place. Therefore it may be necessary to loosen the keypad first. Do this carefully to avoid any damage.




Removing the Main Board

- 1 Lift up the main board as shown below. Then turn it over and continue with step 2.



- 2 Remove the connectors to the main board starting from the right side.




Removing the Measurement Board

- 1 Position the MMS with the connectors facing towards you. There are four long mounting pins threaded into the MMS in each of the four corners under the cover. Locate the heads of the two long mounting pins on the bottom cover. Only these need to be removed.
- 2 Use a thin-bladed screwdriver to gently lift the pins out far enough so they can be removed with pliers.



- 3 Remove the two pins and set them aside for refitting.




NOTE

Without these long mounting pins, the MMS will not function properly.

- 4 Loosen the bottom cover gently. Use the screwdriver and position it in the gap between bottom cover and measurement block, then twist the screwdriver. **Do not push the screwdriver into the device as you may damage electronic components inside.**




- 5 Remove the bottom cover. There may be a slight resistance when opening the cover.




5 Repair and Disassembly


- Loosen the measurement block by pushing the it block forward while holding plastic chassis. Then remove the measurement block.



- The spacer keeps the measurement boards in place. Remove the spacer by lifting it up.




- Remove the measurement ribbon cable. Twist the cable slightly in order to loosen it.




Reassembling the Measurement Block

- 1 Insert the spacer to keep the measurement block in place.
- 2 Make sure that you insert the measurement cable correctly. Please refer to the picture below for the appropriate orientation. A cable inserted incorrectly may cause a <Measurement>Malfunction INOP.




- 3 Make sure the measurement block is inserted as shown below. Make sure that there is no gap between the chassis and the bottom.




Refitting the Main Board

- 1 Insert the new main board. Make sure the main board is seated correctly.




- 2 Re-establish the connection to the measurements. Make sure the connection is tight. If it is not connected correctly, a corresponding measurement malfunction INOP may occur.



Refitting the Keypad

- 1 Reinsert the key pad. Make sure it is positioned correctly.



Refitting the new NBP Pump

- 1 Insert new silicon tubes. Make sure they are seated correctly by pressing them into their position.




- 2 Insert the new pump assembly. Lift up the back and press the airways onto the silicon tubes.



- 3 Make sure the airways have a tight connection to the silicon tubes.



- 4 Insert the connector of the NBP assembly into the connector on the main board. Do not crease the flex cable. M3001A HW A/B and M3000A have a post connector with long pins. Press down the connector until there is no gap between the connectors.




Refitting the DC/DC board

NOTE

This step only has to be done on HW A/B

Position the DC/DC board and press it down gently. Make sure it is connected properly to both connectors indicated in the picture.



Refitting the Cover


NOTE

To change the top cover of a HW Rev C MMS (S/N prefix DE610xxxxx) you have to attach the two cushions which are part of the Top Cover Assembly. These two cushions secure the connection of the MSL Flex and the NBP Flex.


NOTE

Perform the following two steps only on an MMS HW Rev C

- 1 Stick the two cushions onto the marked positions inside the top cover.



- 2 Position the top cover, then press it back into place until you hear a click or there is no longer a gap between the two covers.
- 3 The cover has a rubber seal, press the cover firmly together.



- 4 Holding the bottom cover firmly in place, slide the two long mounting pins completely back into the MMS. Make sure there is no gap between the top and bottom cover.



Refitting the Front Cover

To refit the front cover, press it back into place over the measurement connector hardware until you hear a click.



Final Inspection

Perform a final inspection to ensure that:

- The MSL connector is positioned correctly
- There are no gaps between the MSL connector and the cover
- there is no gap between the top and bottom cover

Testing

To ensure that the MMS is functioning correctly, you must perform safety tests and a performance check on it. Please refer to the "Testing and Maintenance" chapter of this service guide.

WARNING

When you disassemble/assemble an MMS, an applied part leakage current test must be performed before it is used again for monitoring.

MMS Extensions - Exchanging the Top Cover, MSL Flex Cable and the Dual Link Bar

This section describes the exchange procedures for:

- The Top Cover with new release mechanism
- The Dual Link Bar incl. the MSL Flex Cable.

for all MMS Extension (MSE) types (M3012A, M3014A, M3015A, M3016A).



Exchange Procedures

NOTE

Please follow the disassembly and reassembly steps closely.

Tools Required:

A thin-bladed screwdriver and a thick-bladed screwdriver, ESD mat and wrist strap

WARNING

- Do not open the MSE while it is connected to a monitor.
 - Parts inside the instrument may be contaminated with bacteria. Protect yourself from possible infection by wearing examination gloves during this procedure.
-

NOTE

Once you have reassembled the MSE, you must perform a performance check on it. Please refer to the "Testing and Maintenance" chapter of this service guide .

Removing the Front Cover

- 1 Position the thin-bladed screwdriver in the small slot provided for this purpose. The front cover (Bezel) then clicks away from the Extension. Remove the front cover



NOTE

There might be a slight resistance when you remove the front cover.




Removing the Mounting Pin

- 1 Position the MSE on the dual link bar with the measurement connector hardware facing upwards and the arm of the dual link bar away from you. There are four long mounting pins threaded into the MSE in each of the four corners under the cover. Locate the heads of the two long mounting pins on the top housing and only remove these.
- 2 Use the thin-bladed screwdriver to lift the pins gently out far enough so they can be removed manually.



- 3 Remove the two pins and set them aside for refitting.




NOTE

Without these long mounting pins the MSE will not function properly.

Removing the Dual Link Bar

The Dual Link Bar consists of three parts as shown below. Follow the specific steps carefully to remove the Link Bar.



CAUTION

Do not try to remove the link bar with force as this can damage the MSL Flex Cable


- 1 Position the MSE with the measurement connector hardware facing towards you.

5 Repair and Disassembly


- 2 Hold the link bar as shown below. While pressing gently on part B, insert a thick-bladed screwdriver between the MSL connector and part A. Twist the screwdriver to the left and at the same time slide part B to the right, so it is released at the top.




- 3 Repeat Step 2 at the bottom.



- 4 Slide part B to the right. If part B fails to move to the side, please repeat steps 2 and 3.




- 5 Now the MSL Flex connector can be moved to the right.




NOTE

Make sure that the movement of the screwdriver does not pinch the MSL flex cable.

- 6 Insert the thin-blade screwdriver behind the release mechanism of part C. Carefully twist the screwdriver, then press gently so that part C drops down.



- 7 Lift part A upwards. It is fixed in a dovetail. Be careful with the MSL flex.



Removing the Top Cover

Begin by gently pulling away the top cover from the MSE. The top cover is press-latched at the link bar end. Remove it slowly, without hitting or touching the inside of the MSE.




Replacing the Flex Cable Assembly


- 1 Hold the Extension firmly and push upwards against the connector. Then slide connector (together with the connector holder) out of the dovetail connection.

NOTE


You will probably need to apply some more force at first until the holder slides out of its mechanical lock.




- 2 Slide the connector out of its holder.




- 3 Remove the flex cable connector on the MSE board. Be careful not to bend any pins on the female part of the MSE connector.

**NOTE**

Some units may have a foam pad on the connector of the inner flex cable of the MSEs (as shown below) and some units may not. This has no impact on the functionality of these units.





- 4 Stick the correct foam pad on the rear side of the inner connector. Use the thick pad for : M3012A, M3014A, M3016A. Use the thin pad for: M3015A. You can also check the old flex cable for the correct pad.



5 Repair and Disassembly


- 5 Insert the flex cable connector into the female receptacle on the MSE board. **Check from the side and the front that the connector is inserted correctly (there is no mechanical guidance) and that no pins are bent, otherwise you may damage the MSE when powering it on.**



WARNING

A misplaced connector might damage the MSE or the monitor.

- 6 Slide the connector into the holder as shown below. Arrange the flex cable in the space beside and underneath the board (be careful not to bend the cable) while positioning the holder for insertion.



- 7 Insert the holder with the connector into the dovetail connection and slide it down until you hear a click.




Refitting the Top Cover


NOTE

Be careful with the MSL Flex cable. Make sure it does not get stuck between the covers.

- 1 Position top cover, then press the bottom cover back into place until a click is heard.
- 2 The cover has a rubber seal. Press the covers firmly together and make sure there is no gap between the top and bottom cover.



- 3 Holding the bottom cover firmly in place, slide the two long mounting pins completely back into the MSE.



Assembling the dual Link Bar


CAUTION

Do not try to assemble any part of the link bar with force as this can damage the MSL Flex Cable.

- 1 Position part A into the dovetail and slide it down.




- 2 Make sure the MSL Flex connector is positioned in the correct slot (See indicated slots below). Then push it gently into part A.



- 3 Making sure the MSL flex cable lies flat in part A of the assembly, place part B into the dovetail and close the open link bar.



- 4 Turn the MSE around and insert part C into the bottom part of the link bar. When you hear a click, part C is correctly inserted.



Refitting the Front Cover

To refit the front cover, press it back into place over the measurement connector hardware until you hear a click.



Final Inspection

Perform a final inspection to ensure that:

- The link bar is positioned correctly
- There are no gaps between the link bar parts
- There is no gap between the top and bottom cover



Testing

To ensure that the MSE is functioning correctly, you must perform a performance check on it. Please refer to the "Testing and Maintenance" chapter of this service guide.

WARNING

When you disassemble/assemble an MMS Extension, an applied part leakage current test must be performed before it is used again for monitoring.

Disassembly Procedures for the M3015A MMS Extension (HW Rev. A)

NOTE

These procedures apply only to M3015A MMS with Serial Numbers DE020xxxxx.

It is recommended that you replace all the replaceable parts in the Extension (CO2 Scrubber and Pump) after 15 000 hours (approximately 3 years) of continuous use.

Tools Required:

- A thin-bladed screwdriver.
- A pair of large tweezers.
- In addition, for removing the pump, you will need a large-bladed screwdriver.

WARNING

There is high voltage inside the Instrument (800V). Do not connect the MMS Extension to a Monitor while the Extension housing is open.

As well, parts inside the Instrument may be contaminated with bacteria. Protect yourself from possible infection by wearing examination gloves during these procedures.

Removing the Front Cover

To remove the front cover, do the following:

- 1 Remove the server and the monitor from the extension.
- 2 Use a thin-bladed screwdriver to prise the grey front cover (the console covering the measurement connector hardware) gently from the bottom of the extension. Position the screwdriver in the small slits provided for this purpose. The front cover then clicks away from the extension.
- 3 Remove the front cover.

Removing the Extension Bottom Cover

To remove the Extension bottom cover, do the following:

- 1 Position the extension on the dual link bar with the measurement connector hardware facing upwards and the arm of the dual link bar towards you. There are four long mounting pins threaded into the extension in each of the four corners under the cover. Locate the heads of the two long mounting pins on the side away from you
- 2 Use tweezers to prise the pins gently out enough to be removed by hand.
- 3 Remove the two pins and set them aside for refitting.

NOTE


Do not lose these long mounting pins since the Extension will not function unless they are in place.

- 4 Using your hands, gently pry the bottom cover away from the Extension at the link bar end first. The bottom cover is press-latched at the link bar end. Remove it gently making sure not to bang or touch the inside of the Extension.

NOTE

If you accidentally try to remove the wrong side of the bottom cover, you will notice that it is attached to the inside of the Extension with a ribbon connector and that the dual link bar prevents you from removing it completely. **Do not try to forcibly remove the wrong side of the M3015A cover; you cannot access replaceable parts from this side.**

The following illustration shows the location of the replaceable parts in the M3015A Measurement Server Extension.



Removing the CO2 Scrubber

To remove the CO2 Scrubber, do the following:

- 1 Locate the CO2 Scrubber in the Extension.
- 2 Being careful not to touch anything else in the Extension, use tweezers to pull the body of the CO2 Scrubber out of the bracket.
- 3 Holding the body of the CO2 Scrubber with your fingers, carefully disconnect the Extension intake tube from the scrubber end and remove the CO2 Scrubber from the Extension.
- 4 Dispose of the CO2 Scrubber according to local legal requirements for low volume chemical waste.


NOTE

Now that it is exposed, do **NOT** allow anything to fall into the Infrared Lamp assembly.

Removing the Pump

To remove the Pump, do the following:

- 1 Locate the Pump in the Extension.
- 2 Being careful not to touch anything else in the Extension, unscrew the screw holding the pump bracket in position. Lift the top part of the bracket away and lift out the pump.



- 3 Gently disconnect the flow tubing attached to the Extension from the Pump.

NOTE

Be sure to note which tube attaches to the inlet and which tube attaches to the outlet.

- 4 Gently disconnect the power lead which attaches the Pump to the Extension.
- 5 Remove the Pump.

NOTE

After replacing the Pump, reset the displayed value displayed using the Reset PumpOpTime selection (Service Mode>CO2 Setup). When the PumpOpTime has been reset an INOP will be generated: "CO₂ OCCLUSION". To clear this INOP you must perform a flow check and store the flow in Service Mode (select "Store Flow")

Refit Procedures for the MMS Extension

Tools Required:

- A thin-bladed screwdriver.
- A pair of large tweezers.
- In addition, for refitting the Pump, you will need a large-bladed screwdriver.

WARNING

There is high voltage inside the Instrument (800V). Do not connect the MMS Extension to a Monitor while the Extension housing is open.

As well, parts inside the instrument may be contaminated with bacteria. protect yourself from possible infection by wearing examination gloves during these procedures.

Refitting the CO2 Scrubber

WARNING

The CO2 Scrubber contains lithium hydroxide monohydrate. This is a strong base. Do not open or damage the CO2 Scrubber. If you come into contact with the CO2 Scrubber material, flush the area immediately with water and consult a doctor.

To refit the CO2 Scrubber, do the following:

- 1 O2 Scrubber through the bracket to meet the Extension intake tube.
- 2 Push the intake tube firmly into the scrubber end to connect it.
- 3 Holding the body of the CO2 Scrubber with tweezers, feed the CO2 Scrubber fresh air intake under the second bracket and position it.

Refitting the Pump

To refit the Pump, do the following:

- 1 Gently connect the power lead to the Extension.

NOTE

The power lead can only be connected one way. Do not try to force the power lead into position. Instead, align it correctly and connect it gently.

- 2 Connect the flow tubing to the Pump.

NOTE

Be sure to reconnect the inlet tube to the inlet valve and the outlet tube to the outlet valve.

- 3 Being careful not to touch anything else in the Extension, insert the pump into the bracket on the PC board. Make sure that the pump is horizontal and does not touch the PC board. (Vibration from the pump in operation will damage the Extension if the pump touches the PC board.)
- 4 Replace the top part of the bracket and screw firmly into position.


NOTE

After replacing the Pump, reset the displayed value using the Reset PumpOpTime selection (Service Mode>CO2 Setup). When the PumpOpTime has been reset an INOP will be generated: “CO₂ OCCLUSION”. To clear this INOP you must perform a flow check and store the flow in Service Mode (select “Store Flow”).

Refitting the Extension Bottom Cover

To refit the Extension bottom cover, do the following:

- 1 Latch the link bar end into place then press-click the bottom cover back into place covering the interior of the Extension.
- 2 Holding the bottom cover firmly in place, thread the two long mounting pins back into the Extension making sure to thread them all the way to the end.



Refitting the Front Cover

To refit the front cover, press-click it back into place over the measurement connector hardware.

General Reassembly/Refitting Comments

- Ribbon Connections—Make sure male-female ribbon connections are correctly lined-up.
- Open Component—Do not allow anything to fall into the open component.

Testing

To ensure that the MMS Extension is functioning correctly, you must perform a performance check on it. Please refer to the *Testing and Maintenance* chapter of this service guide.


When you disassemble/assemble an MMS Extension, an applied part leakage current test must be performed before it is used again for monitoring.

Parts

This section lists the replacement and exchange parts for the following Philips IntelliVue Patient Monitoring System components:

- “MX600/MX700 Parts” on page 228
- “Remote Control Parts” on page 232
- “8-Slot Flexible Module Rack (FMS-8) Parts” on page 233
- “4-Slot Flexible Module Rack (FMS-4) Parts” on page 235
- “Multi-Measurement Module (MMS) Parts” on page 237
- “MMS Extension Parts (M3012A, M3014A, M3015A/B)” on page 244
- “Plug-in Modules Part Numbers” on page 247
- “External Display Part Numbers” on page 259

MX600/MX700 Parts




Exchange Parts

| Part number | 12NC | Description | No in diagram |
|--------------|--------------|--------------------------------------|---------------|
| 453564204481 | 453564204481 | IV2-STAT I/F; Connector Board | 1 |
| 453564212421 | 453564212421 | IV2-STAT PCA Video Mixer board | 2 |
| 453564204431 | 453564204431 | IV2-STAT Main board | 3 |
| 453564172631 | 453564172631 | IV2-STAT ASSY-PWR AC/DC Power Supply | 4 |
| 453564204591 | 453564204591 | IV2-STAT iPC w/o disk incl. Cable | 5 |


Replacement Parts

| Part number | 12NC | Description | No in diagram |
|--------------|--------------|------------------------------------|---------------|
| 453564260811 | 453564260811 | IV2-STAT I/F; Flexible Nurse Call | 6 |
| 453564204361 | 453564204361 | IV2-STAT I/F; Single MIB/ RS232 | 6 |
| 453564204381 | 453564204381 | IV2-STAT Loudspeaker | 7 |
| 453564260821 | 453564260821 | IV2-STAT Panel adapter board MX700 | 8 |
| 453564260831 | 453564260831 | IV2-STAT Panel adapter board MX600 | 8 |
| 453564204411 | 453564204411 | IV2-STAT Alarm LED board | 9 |

| Part number | 12NC | Description | No in diagram |
|--------------|--------------|--|---------------|
| 453564053561 | 453564053561 | IV IIT Module 2.4 GHz | 10 |
| M4840-68708 | 451261009041 | IV IIT Module 1.4 GHz | 10 |
| 453564258611 | 453564258611 | IV SRR brd ver2 (for MP5/T/SC+MX6/7/800) | 11 |
| 453564204471 | 453564204471 | IV2-STAT Power switch/ ECG sync out bd | 12 |
| 453564211691 | 453564211691 | IV2-STAT Triband antenna for IIT/ WLAN | 13 |
| 453564207391 | 453564207391 | IV2-STAT iPC RAM DDR3 2 GB | n/a |
| 453564207371 | 453564207371 | IV2-STAT iPC SSD with pre-installed OS | 14 |
| 453563464231 | 453563464231 | BAT 3V Lithium CR2032 | n/a |
| 453564260851 | 453564260851 | IV2-STAT DSPL Backlights (2x) MX600/700 | n/a |
| 453564260861 | 453564260861 | IV2-STAT DSPL LCD panel MX600/700 | 15 |
| 453564260871 | 453564260871 | IV2-STAT Touch bezel assembly MX700 | 16 |
| 453564260881 | 453564260881 | IV2-STAT Non Touch bezel assembly MX600 | 16 |
| 453564204641 | 453564204641 | IV2-STAT Quick mount | 17 |
| 453564204651 | 453564204651 | IV2-STAT Power button | 18 |
| 453564260891 | 453564260891 | IV2-STAT Rear dspl housing MX600/700 | 19 |
| 453564204691 | 453564204691 | IV2-STAT Rear housing | 20 |
| 453564204701 | 453564204701 | IV2-STAT Bottom housing | 21 |
| 453564260911 | 453564260911 | IV2-STAT Key Silicone Pad | 22 |
| 453564260951 | 453564260951 | IV2-STAT Trim Knob | 23 |
| 453564260991 | 453564260991 | IV2-STAT HIF board | 24 |
| 453564261021 | 453564261021 | IV2-STAT Encoder | 25 |
| 453564204491 | 453564204491 | IV2-STAT WLAN Assy | 26 |
| 453564204511 | 453564204511 | IV2-STAT iPC cable kit | n/a |
| 453564204521 | 453564204521 | IV2-STAT cable management kit | 27 |
| 453564204541 | 453564204541 | IV2-STAT Small parts kit | n/a |
| 453564261261 | 453564261261 | IV2-STAT Basic cable kit MX600/700 | n/a |
| 453564239731 | 453564239731 | IV2 MECHASY Table Top Mount Kit | n/a |
| 453564260841 | 453564260841 | IV2-STAT DSPL Backlt inver. bdMX600/700 | 28 |

WLAN Kit (453564204491) Contents:


| Item | Quantity |
|--|----------|
| PCA MMS+ WLAN Radio Board | 1 |
| IV2-STAT PLAST HOLDER WLAN | 1 |
| IV2-STAT PCA Antenna triband | 1 |
| CBL ASSY AMC RIGHT ANGL PLG JMPER 14MM | 1 |

iPC Cable Kit (453564204511) Contents:

| Item | Quantity |
|----------------------------------|----------|
| IV2-STAT CBL MB-Mixer | 1 |
| CBL Ay IV2-STAT SAM-Mixer | 1 |
| IV2-STAT MECHASY Assy Holder USB | 1 |
| Cable clip | 2 |

Cable Management Kit (453564204521) Contents:

| Item | Quantity |
|-------------------------------------|----------|
| SCRW TPG M3.0 X 1.34 8MM-LG WSHR-HD | 4 |
| IV2-STAT PLAST Cablehouse | 1 |
| IV2-STAT PLAST Cablehouse Lid | 1 |

Small Parts Kit (453564204541) Contents:

| Item | Quantity |
|------------------------------------|----------|
| IV2-STAT PLAST Lightguide P/S | 1 |
| IV2-STAT PLAST IO-Slot Cover Blank | 2 |
| IV2-STAT PLAST Cover USB Side | 1 |
| IV2-STAT PLAST Sync Side Plate | 1 |
| IV2-STAT PLAST LIGHTPIPE LARGE | 1 |
| IV2-STAT PLAST LIGHTPIPE SMALL | 2 |
| IV2-STAT PLAST PIN 265MM | 2 |
| Cable clip | 2 |

6 Parts

| Item | Quantity |
|-------|----------|
| K30x8 | 4 |
| M3x5 | 4 |
| M3x8 | 8 |
| M4x25 | 4 |

Basic Cable Kit (453564261261) Contents:




| Item | Quantity |
|------------------------------|----------|
| IV2-STAT CBL MB-ALARM | 1 |
| IV2-STAT CBL MB-IIT | 1 |
| IV2-STAT CBL MB-SRR | 1 |
| IV2-STAT CBL MB-AC/DC | 1 |
| IV2-STAT CBL MB-Power Switch | 1 |
| IV2-STAT CBL PA-Backlight | 1 |
| IV2-STAT CBL MB-19" Display | 1 |
| CBL MB-HIF | |

Remote Control Parts

| Part number | 12NC | Description |
|--------------|--------------|-------------------------------------|
| 453564212481 | 453564212481 | IV2-RC Remote Control |
| 453564212401 | 453564212401 | IV2-RC Holder for Remote Control |
| 453564212411 | 453564212411 | IV2-RC USB cable for Remote Control |
| 453564262641 | 453564262641 | IV2-RC CBL Tethering cable |

8-Slot Flexible Module Rack (FMS-8) Parts



Flexible Module Rack (FMS-8) Parts

Exchange and Replacement Parts


| Part number | 12NC | Description | No. in Diagram |
|-------------|--------------|--|----------------|
| M4041-68401 | 453563459411 | IV-FMS Mother board assembly | 1 |
| M8055-68401 | 453563459441 | IV CPU board assembly | 2 |
| M8048-64002 | 453563456901 | IV-FMS Small Parts kit | 3 |
| M8048-64001 | 453563456891 | IV-FMS Housing kit | 4 |
| M4041-60005 | 453563477961 | IV-FMS Mounting Plate Assy | 5 |
| M4041-42303 | 453563494101 | IV-FMS PLAST Cable Management | 6 |
| M4041-22302 | 451261011861 | IV-FMS BRKT Mounting Plate Adapter Clamp | |
| M8040-60100 | 451261030081 | IV SWITCH Universal Clamp | |




M8048-64002 Small Parts Kit Contents

| Description | Quantity | Comments |
|------------------------|----------|---|
| Torx M3 x 6 mm screws | 8 | used for securing mainboard to FMS housing |
| Torx M3 x 20 mm screws | 8 | used to connect the CPU board w/ spacers to the mainboard |
| Torx M4 x 8 mm screws | 4 | used to secure the handle or the MMS mount to the FMS housing |
| Torx M4 x 12 mm screws | 4 | used to secure the universal clamp to the FMS housing |
| Connector Holder (FE) | 2 | |
| Connector Holder (SRL) | 2 | |
| Lightpipe | 1 | |
| Cover Connector SRL | 2 | |
| Cover Side | 1 | |
| Pin | 4 | |
| Cover Seal Connector | 2 | |
| Handle Assembly | 1 | |

4-Slot Flexible Module Rack (FMS-4) Parts




4-Slot Flexible Module Rack (FMS-4) without MMS Mount



4-Slot Flexible Module Rack (FMS-4) with MMS Mount

Exchange and Replacement Parts

| Part number | 12NC | Description | No in diagram |
|--------------------------|--------------|---------------------------------------|---------------|
| Exchange Parts | | | |
| 453564261731 | 453564261731 | IV2-FMS4 Main board | 1 |
| Replacement Parts | | | |
| 453564262101 | 453564262101 | IV2-FMS4 MSL up board | 2 |
| 453564262111 | 453564262111 | IV2-FMS4 MSL down board | 3 |
| 453564262121 | 453564262121 | IV2-FMS4 Cable management hooks | 4 |
| 453564262131 | 453564262131 | IV2-FMS4 MMS mount | 5 |
| 453564262141 | 453564262141 | IV2-FMS4 Assy Front MMS Option | 6 |
| 453564262151 | 453564262151 | IV2-FMS4 MECHASY Front w/o MMS Option | 7 |
| 453564262161 | 453564262161 | IV2-FMS4 MECHASY Rear w/o MMS Option | |
| 453564262171 | 453564262171 | IV2-FMS4 MECHASY Rear MMS Option | |
| 453564262181 | 453564262181 | IV2-FMS4 Side MSL Holder | 8 |
| 453564262191 | 453564262191 | IV2-FMS4 Small Parts Kit | n/a |
| M8040-60100 | M8040-60100 | IV SWITCH Universal Clamp | n/a |




453564262191 Small Parts Kit Contents

| Description | Quantity |
|-------------------------|----------|
| Light guide | 1 |
| Pin | 2 |
| M4 x 12 mm screws (T20) | 5 |
| K30 x 8 mm screws (T10) | 5 |
| M3 x 8 mm screws (T10) | 10 |
| Connector Holder | 2 |
| MSL Holder | 2 |

| Description | Quantity |
|-----------------------------|----------|
| Cable MSL Up | 1 |
| Cable MSL down | 1 |
| IV2-FMS4 RUBBER Screw Cover | 12 |
| IV2-FMS4 PLAST Cover Bottom | 1 |

Multi-Measurement Module (MMS) Parts

The primary support strategy for the Multi-Measurement Module is a unit exchange. However, some exchange parts are available. In order to determine which exchange parts need to be ordered check the serial number and the option string of the MMS as described below.



M3001A Multi-Measurement Module

MMS Part Number Overview and Identification

Identify the correct MMS part number by checking the serial number prefix and the option string on the label on the rear of the MMS housing. The following picture shows the label of an M3001A. Use the table below to determine the hardware revision and the required exchange part.



| Option | M3001A Option Description |
|--------|--|
| #A01 | Philips FAST SpO ² |
| #A02 | Nellcor OxiMax compatible |
| #A03 | Masimo SET Technology |
| #A04 | Nellcor OxiMax Technology |
| (#C00) | Standard |
| #C06 | Add Pressure/Temp |
| #C12 | Add 12 Lead ECG (only older revisions - see table below) |
| (#C18) | (Add Pressure/Temp and 12 Lead ECG) old |

With Rel. G.0 option #C18 is split into #C06#C12

| HW Revision | Possible MMS Software Revision | Serial Number Prefix | SW of monitor the MMS is connected to | Option String | Exchange MMS (for 12NC information please refer to tables in the following sections) |
|-------------|--------------------------------|----------------------|---------------------------------------|--|--|
| HW A | A.0 to H.0 | DE227 | A.0, A.1 | --, C06, C12, C18, C06C12 | M3001-68x10 |
| | | DE441 | A.0, A.1 | A01, A01C06, A01C12, A01C18, A01C06C12 | |
| | | DE227 | >=A.2 | --, C06, C12, C18, A01C06C12 | M3001-68x02 |
| | | DE441 | >=A.2 | A01, A01C06, A01C12, A01C18, A01C06C12 | |
| HW B | B.1 to H.0 | DE441 | >=A.2 | A02, A02C06, A02C18, A01C06C12 | M3001-68113 M3001-68x03 |
| | | DE512 | >=A.2 | A02, A02C06, A02C18, A01C06C12 | |
| | | | >=A.2 | A01,A01C06, A01C12, A01C18, A01C06C12 | M3001-68114 M3001-68x04 |

| HW Revision | Possible MMS Software Revision | Serial Number Prefix | SW of monitor the MMS is connected to | Option String | Exchange MMS (for 12NC information please refer to tables in the following sections) |
|-------------|--------------------------------|----------------------|---------------------------------------|--------------------------------|--|
| HW C | D.0 to H.0 | >=DE610 | >=A.2 | A01, A01C06, A01C18, A01C06C12 | M3001-68x05 |
| | | | >=A.2 | A02, A02C06, A02C18, A01C06C12 | M3001-68x08 |
| | | >=DE632 | >=A.2 | A03, A03C06, A03C18, A01C06C12 | M3001-68x07 |
| | | >=DE907 | >=A.2 | A04, A04C06, A04C06C12 | M3001-68x06 |

For further compatibility information please refer to the Software Compatibility Matrix in the Troubleshooting section.

Exchange Multi-Measurement Modules are shipped with English front bezels only. If you require a bezel in another language (compare the part numbers of your language to the English ones to check this) the front bezel has to be ordered additionally. Attach the appropriate bezel before putting the MMS into operation.

MMS Firmware

NOTE

Multi-Measurement Modules (MMS) with HW Rev. A have a fixed firmware that cannot be upgraded. HW B and HW C MMSs allow upgrading of the SpO2 and ECG firmware.

To perform a FW upgrade, the MMS must have SW Revision C.0 or higher and be connected to an IntelliVue patient monitor.

SpO2

| Option | HW Revision | SpO2 FW Rev | Comment |
|--------------|----------------|----------------------|---|
| #A01 | HW A | A.01.04 | The HW and its interface do not allow a FW upgrade |
| #A01 #A02 | HW B | (A.01.41) -> A.01.46 | Must be upgraded to A.01.46 - see internal Field Notification |
| | | (A.01.42) -> A.01.47 | |
| #A01 #A02 | HW C | A.01.46 | |
| | HW C (B.00.03) | A.04.12 | Currently shipped FW |
| #A03 | HW C | Masimo SET FW | Not upgradeable by Philips |
| #A04 | HW C | Nellcor OxiMax FW | Philips upgradeable |

ECG

| Option | HW Revision | ECG FW Rev | Comment |
|------------------------------|-------------|------------|---|
| #A01 | HW A | C.00.13 | The HW and its interface do not allow a FW upgrade |
| #A01 | HW A | C.01.19 | |
| #A01 #A02 #A03 #A04 | HW B/C | D.01.70 | Upgradeable - see internal Field Notifications |
| | | D.01.76 | |
| | | D.01.77 | |
| | | D.01.78 | |
| | | D.01.89 | |
| | | D.02.02 | |
| | | D.02.05 | Currently shipped FW (requires MMS SW Revision F.0) |

MMS Part Numbers - Front Bezels

| Part number | 12NC | Description | Options |
|--------------|--------------|---------------------------------------|------------|
| 451261024391 | 451261024391 | MS_X1 Bezel w/o P/T Eng. Text | #A01, #A03 |
| 451261024401 | 451261024401 | MS_X1 Bezel w P/T Eng. Text | #A01, #A03 |
| 451261024411 | 451261024411 | MS_X1 Bezel w/o P/T Symbols | #A01, #A03 |
| 451261024421 | 451261024421 | MS_X1 Bezel w P/T Symbols | #A01, #A03 |
| 451261024431 | 451261024431 | MS_X1 Bezel Nellcor w/o P/T Eng. Text | #A02, #A04 |
| 451261024441 | 451261024441 | MS_X1 Bezel Nellcor w P/T Eng. Text | #A02, #A04 |
| 451261024451 | 451261024451 | MS_X1 Bezel Nellcor w/o P/T Symbols | #A02, #A04 |
| 451261024461 | 451261024461 | MS_X1 Bezel Nellcor w P/T Symbols | #A02, #A04 |

MMS Support Parts

| Part number | 12NC | Description | Comments |
|--------------|--------------|---|-----------|
| M3000-66541 | 453564107971 | MS_X1 PCA DC/DC Board for HW A/B | |
| M3001-64050 | 451261016391 | MS_X1 MMS MSL Connector Assy | |
| M3001-64500 | 451261020561 | MS_X1 NBP Assembly for X1/X2/MP2 | |
| M3001-68557 | 451261020751 | MS_X1 PCA ECG-5ld/Fast SpO2 | HW C only |
| 453564186021 | 453564186021 | MS_X1 PCA ECG-5ld/Fast SpO2/w PressTemp | HW C only |
| 453564186041 | 453564186041 | MS_X1 PCA ECG-5ld/OxiMax comp SpO2 | HW C only |
| 453564186031 | 453564186031 | MS_X1 PCA ECG-5ld/OxiMax comp SpO2/w PT | HW C only |

| Part number | 12NC | Description | Comments |
|--------------|--------------|--|-----------|
| M3001-68553 | 451261020771 | MS_X1 PCA ECG-5ld/Masimo SpO2 | HW C only |
| 453564186051 | 453564186051 | MS_X1 PCA ECG-5ld/Masimo SpO2/w PT | HW C only |
| M3001-68555 | 451261020791 | MS_X1 PCA ECG-5ld / Nellcor SpO2 | HW C only |
| 453564186081 | 453564186081 | MS_X1 PCA ECG-5ld/Nellcor SpO2/w PT | HW C only |
| M3001-68425 | 453564177921 | MS_X1 Main board for HW C | HW C only |
| 5041-8114 | 453563100081 | MS_X1 Mounting Pin | |

MMS Part Numbers - Top Covers

| Part number | 12NC | Description |
|-------------|--------------|--|
| M3001-68010 | 451261016401 | MS_X1 Top Cover 5ld w/o P/T Text FAST |
| M3001-68011 | 451261016411 | MS_X1 Top Cover 5ld w P/T Text FAST |
| M3001-68012 | 451261016421 | MS_X1 Top Cover 12ld w/o P/T Text FAST |
| M3001-68013 | 451261016431 | MS_X1 Top Cover 12ld w P/T Text FAST |
| M3001-68014 | 451261016441 | MS_X1 Top Cover 5ld w/o P/T Sym FAST |
| M3001-68015 | 451261016451 | MS_X1 Top Cover 5ld w P/T Sym FAST |
| M3001-68016 | 451261016461 | MS_X1 Top Cover 12ld w/o P/T Sym FAST |
| M3001-68017 | 451261016471 | MS_X1 Top Cover 12ld w P/T Sym FAST |
| M3001-68018 | 451261016481 | MS_X1 Top Cover 5ld w/o P/T Text NELLCOR |
| M3001-68019 | 451261016491 | MS_X1 Top Cover 5ld w P/T Text NELLCOR |
| M3001-68020 | 451261016501 | MS_X1 Top Cover 12ld w P/T Text NELLCOR |
| M3001-68021 | 451261016511 | MS_X1 Top Cover 5ld w/o P/T Sym NELLCOR |
| M3001-68022 | 451261016521 | MS_X1 Top Cover 5ld w P/T Sym NELLCOR |
| M3001-68023 | 451261016531 | MS_X1 Top Cover 12ld w P/T Sym NELLCOR |
| M3001-68030 | 453564147051 | MS_X1 TopCover 12ld w/o P/T Text NELLCOR |
| M3001-68031 | 453564147061 | MS_X1 Top Cover 12ld w/o P/T Sym NELLCOR |
| M3001-68024 | 451261016541 | MS_X1 Top Cover 5ld w/o P/T Text MASIMO |
| M3001-68025 | 451261016551 | MS_X1 Top Cover 5ld w P/T Text MASIMO |
| M3001-68026 | 451261016561 | MS_X1 Top Cover 12ld w P/T Text MASIMO |
| M3001-68027 | 451261016571 | MS_X1 Top Cover 5ld w/o P/T Sym MASIMO |
| M3001-68028 | 451261016581 | MS_X1 Top Cover 5ld w P/T Sym MASIMO |
| M3001-68029 | 451261016591 | MS_X1 Top Cover 12ld w P/T Sym MASIMO |
| M3001-68032 | 453564147071 | MS_X1 Top Cover 12ld w/o P/T Text MASIMO |
| M3001-68033 | 453564147081 | MS_X1 Top Cover 12ld w/o P/T Sym MASIMO |

MMS Exchange Part Numbers

NOTE

The MMS always ships with the latest Software Revision. In order to make it compatible with the respective monitor the MMS may need to be upgraded or downgraded. From Support Tool version E.03.01 onwards the MMS can be up- or downgraded with every support tool license key (except general).

| Part number | 12NC | Description |
|--|--------------|--|
| X1 HW Rev: A; Option: A01; SN Prefix: DE227, DE441; SW Rev: Latest SW | | |
| M3001-68102 | 453563462911 | MS_X1 Exch Std 5ld w/o P/T-Eng latest SW |
| M3001-68202 | 453563486921 | MS_X1 Std 5ld w P/T #C06-Eng latest SW |
| M3001-68302 | 453563486931 | MS_X1 12ld w/o P/T #C12-Eng latest SW |
| M3001-68402 | 453563486941 | MS_X1 Std 12ld w P/T #C18-Eng latest SW |
| M3001-68502 | 453563486951 | MS_X1 Exch Std 5ld w/o P/T-Sym latest SW |
| M3001-68602 | 453563486961 | MS_X1 Std 5ld w P/T #C06-Sym latest SW |
| M3001-68702 | 453563486971 | MS_X1 12ld w/o P/T #C12-Sym latest SW |
| M3001-68802 | 453563486981 | MS_X1 Std 12ld w P/T #C18-Sym latest SW |
| X1 HW Rev: B; Option: A02; SN Prefix: DE441, DE512; SW Rev: Latest SW | | |
| M3001-68113 | 451261005361 | MS_X1 5ld w/o P/T-Eng - OxiMax |
| M3001-68203 | 451261005381 | MS_X1 5ld w P/T-Eng - OxiMax HWB |
| M3001-68403 | 451261005401 | MS_X1 12ld w P/T-Eng - OxiMax HWB |
| M3001-68503 | 451261005421 | MS_X1 Std 5ld w/o P/T-Sym - OxiMax HWB |
| M3001-68603 | 451261005441 | MS_X1 5ld w P/T-Sym - OxiMax HWB |
| M3001-68803 | 451261005461 | MS_X1 12ld w P/T-Sym - OxiMax HWB |
| X1 HW Rev: B; Option: A01; SN Prefix: DE512; SW Rev: Latest SW | | |
| M3001-68114 | 451261006041 | MS_X1 Std 5ld w/o P/T-Eng - HW B |
| M3001-68204 | 451261006061 | MS_X1 5ld w P/T-Eng - HW B |
| M3001-68304 | 451261006081 | MS_X1 12ld w/o P/T #C12-Eng - HW B |
| M3001-68404 | 451261006101 | MS_X1 12ld w P/T-Eng - HW B |
| M3001-68504 | 451261006121 | MS_X1 5ld w/o P/T-Sym - HW B |
| M3001-68604 | 451261006141 | MS_X1 5ld w P/T-Sym - HW B |
| M3001-68704 | 451261006161 | MS_X1 12ld w/o P/T #C12-Sym - HW B |
| M3001-68804 | 451261006181 | MS_X1 12ld w P/T-Sym - HW B |
| X1 HW Rev: C; Option: A01; SN Prefix: DE610 <= ; SW Rev: Latest SW | | |
| M3001-68105 | 451261013041 | MS_X1 5ld w/o P/T-Eng - HW C |
| M3001-68205 | 451261013061 | MS_X1 5ld w P/T - Eng - HW C |
| M3001-68305 | 453564146221 | MS_X1 12ld wo P/T-Eng - HW C |
| M3001-68405 | 451261013081 | MS_X1 12ld w P/T -Eng - HW C |

| Part number | 12NC | Description |
|--|--------------|-------------------------------------|
| M3001-68505 | 451261013101 | MS_X1 5ld w/o P/T-Sym - HW C |
| M3001-68605 | 451261013121 | MS_X1 5ld w P/T -Sym - HW C |
| M3001-68705 | 453564146241 | MS_X1 12ld w/o P/T-Eng - HW C |
| M3001-68805 | 451261013141 | MS_X1 12ld w P/T - Sym - HW C |
| X1 HW Rev: C; Option: A04; SN Prefix: DE907 <= ; SW Rev: Latest SW | | |
| M3001-68106 | 451261013161 | MS_X1 5ld w/o P/T-Eng - Nellcor |
| M3001-68206 | 451261013181 | MS_X1 5ld w P/T -Eng - Nellcor |
| M3001-68306 | 453564146171 | MS_X1 12ld wo P/T -Eng- Nellcor |
| M3001-68406 | 451261013201 | MS_X1 12ld w P/T -Eng - Nellcor |
| M3001-68506 | 451261013221 | MS_X1 5ld w/o P/T-Sym - Nellcor |
| M3001-68606 | 451261013241 | MS_X1 5ld w P/T -Sym - Nellcor |
| M3001-68706 | 453564146181 | MS_X1 12ld wo P/T -Sym- Nellcor |
| M3001-68806 | 451261013261 | MS_X1 12ld w P/T -Sym - Nellcor |
| X1 HW Rev: C; Option: A03; SN Prefix: DE632 <= ; SW Rev: Latest SW | | |
| M3001-68107 | 451261013281 | MS_X1 5ld w/o P/T-Eng-Masimo-HWC |
| M3001-68207 | 451261013301 | MS_X1 5ld w P/T-Eng-Masimo-HWC |
| M3001-68307 | 453564146301 | MS_X1 12ld wo P/T-Eng - Masimo HW C |
| M3001-68407 | 451261013321 | MS_X1 12ld w P/T -Eng-Masimo-HWC |
| M3001-68507 | 451261013341 | MS_X1 5ld w/o P/T-Sym-Masimo-HWC |
| M3001-68607 | 451261013361 | MS_X1 5ld w P/T -Sym-Masimo-HWC |
| M3001-68707 | 453564146321 | MS_X1 12ld w/o P/T-Eng- Masimo-HWC |
| M3001-68807 | 451261013381 | MS_X1 12ld w P/T - Sym-Masimo-HWC |
| X1 HW Rev: C; Option: A02; SN Prefix: DE610 <= ; SW Rev: Latest SW | | |
| M3001-68108 | 451261015171 | MS_X1 5ld w/o P/T Eng-OxiMax-HWC |
| M3001-68208 | 451261015191 | MS_X1 5ld w P/T-Eng-OxiMax-HWC |
| M3001-68308 | 453564146351 | MS_X1 12ld wo P/T-Eng- OxiMax HW C |
| M3001-68408 | 451261015211 | MS_X1 12ld w P/T -Eng-OxiMax-HWC |
| M3001-68508 | 451261015231 | MS_X1 5ld w/o P/T-Sym-OxiMax-HWC |
| M3001-68608 | 451261015251 | MS_X1 5ld w P/T -Sym-OxiMax-HWC |
| M3001-68708 | 453564146361 | MS_X1 12ld w/o P/T Eng-OxiMax-HWC |
| M3001-68808 | 451261015271 | MS_X1 12ld w P/T-Sym-OxiMax-HWC |

MMS Extension Parts (M3012A, M3014A, M3015A/B)

Exchange MMS Extensions are shipped with English front bezels only. If you require a bezel in another language (compare the part numbers of your language to the English ones to check this) the front bezel has to be ordered additionally. Attach the appropriate bezel before putting the MMS extension into operation.

The part numbers in the following parts table below, are used to order parts from your Philips representative.

MMS Extension General Support Parts

| Part number | 12NC | Description | Comments |
|--------------|--------------|-------------------------------------|--|
| M3014-64200 | 451261012731 | MS_RCO2 MMS Ext Clips+Springs, 10ea | Old top housing |
| M3012-64600 | 451261012721 | MS_HMS MMS Ext. Lever Locks (5ea) | Grey lever release mechanism |
| M3012-64620 | 451261016601 | MS_HMS MSE Top Cover Assy | |
| M3012-64621 | 451261016611 | MS_HMS MSE Link Bar Assy | Without Flex, must be ordered separately |
| 453564088851 | 453564088851 | MS_MCO2 Flex Cable/Connector Assy | |
| 5041-8114 | 453563100081 | MS_X1 Mounting Pin | |



MMS Extension Front Bezels and Accessories

| Part number | 12NC | Description |
|-----------------------------|--------------|---|
| M3012A BEZELS | | |
| 451261024471 | 451261024471 | MS_HMS Bezel P, T, P/T Eng. Text |
| 451261024481 | 451261024481 | MS_HMS Bezel P, T, P/T Symbols |
| 451261024491 | 451261024491 | MS_HMS Bezel C.O.,P, T, P/T Eng. Text |
| 451261024501 | 451261024501 | MS_HMS Bezel C.O.,P, T, P/T Symbols |
| 451261024511 | 451261024511 | MS_HMS Bezel CCO,C.O.,P,T,P/T Eng Text |
| 451261024521 | 451261024521 | MS_HMS Bezel CCO, C.O.,P,T,P/T Sym |
| M3012A Misc | | |
| 14454A | 453563057651 | I.V POLE MOUNT |
| M3014A Bezels | | |
| 451261024531 | 451261024531 | MS_RCO2 Bezel CO2 only Eng. Text/Sym |
| 451261024541 | 451261024541 | MS_RCO2 Bezel C.O., CO2,P,P/T Eng. Text |
| 451261024551 | 451261024551 | MS_RCO2 Bezel C.O., CO2, P, P/T Symbols |
| 451261024561 | 451261024561 | MS_RCO2 Bezel CO2, P, P/T Eng. Text |
| 451261024571 | 451261024571 | MS_RCO2 Bezel CO2, P, P/T Symbols |
| 451261024581 | 451261024581 | MS_RCO2 Bezel C.O., CCO,CO2,P,P/T Text |
| 451261024591 | 451261024591 | MS_RCO2 Bezel C.O.,CCO, CO2, P, P/T Sym |
| Canpnograhpy Sensors | | |
| M2741-68000 | 451261011291 | SNSR M2741A Sidestream Sensor |
| M2501-68000 | 451261006391 | MS_MCO2 M2501A MainstreamSensor |
| M2741-60000 | 451261011731 | Mounting Bracket Replacement |
| M3015A Bezel | | |
| 451261024601 | 451261024601 | MS_SCO2 Bezel CO2 w P/T Eng. Text |
| 451261024611 | 451261024611 | MS_SCO2 Bezel CO2 w P/T Symbols |
| 451261024621 | 451261024621 | MS_SCO2 Bezel CO2 w/o P/T Eng. Text |
| 451261024631 | 451261024631 | MS_SCO2 Bezel CO2 w/o P/T Symbols |
| M3015B Bezel | | |
| 453564270051 | 453564270051 | MS_SCO2 Bezel Symbols M3015B CO2. P, P, T |
| 453564270061 | 453564270061 | MS_SCO2 Bezel Text M3015B CO2. P, P, T |
| Repair Kit M3015A | | |
| M3015-29303 | 453563332261 | MS_SCO2 OEM M3015A pump kit +CO2scrub |
| M3015-29314 | 453563332281 | MS_SCO2 LAMP M3015A Infra Red lamp |

Exchange Parts List

Exchange parts are parts that have been returned to Philips and reconditioned for further use. Parts offered as exchange parts are in excellent service order according to rigorous Philips standards but offer you a considerable price advantage.

A front bezel with symbols instead of English text is provided with each exchange MMS Extension.

| Part number | 12NC | Description |
|----------------|---------------|---|
| M3012A | 862111 | Hemodynamic Measurement Server Extension |
| M3012-6801A | 451261000201 | MS_HMS OEM EXCH.UNIT ENGLISH (#C00) |
| M3012-6831A | 451261000341 | MS_HMS OEM EXCH.UNIT ENGLISH (#C05) |
| M3012-6861A | 451261000491 | MS_HMS OEM EXCH.UNIT ENGLISH (#C10) |
| M3014A | 862187 | Capnography extension |
| M3014-6801A | 451261009281 | MS_RCO2 OEM EXCH.UNIT ENGLISH (#C00) |
| M3014-6831A | 451261009311 | MS_RCO2 OEM EXCH.UNIT ENGLISH (#C05) |
| M3014-6861A | 451261009601 | MS_RCO2 OEM EXCH.UNIT ENGLISH (#C10) |
| M3014-6891A | 451261009461 | MS_RCO2 OEM EXCH.UNIT ENGLISH (#C07) |
| M3015A | 862393 | Microstream CO2 Extension |
| M3015-6801A | 453563332431 | MS_SCO2 OEM Exch. Unit with P/T (HW Rev. A) |
| M3015-6831A | 453563477871 | MS_SCO2 Exch. Unit without P/T (HW Rev. A) |
| M3015-6802A | 451261005311 | MS_SCO2 M3015A Unit with P/T (HW Rev. B) |
| M3015-6832A | 451261005331 | MS_SCO2 M3015A Unit wo P/T (HW Rev. B) |
| 453564293881 | 453564293881 | MS_SCO2 OEM Exch. Unit with P/T (HW Rev. C) |
| 453564293891 | 453564293891 | MS_SCO2 Exch. Unit without P/T (HW Rev. C) |
| M3015B* | 865377 | Microstream CO2 Extension |
| 453564270041 | 453564270041 | MS_SCO2 M3015B with CO2,P, P and T |

* Requires SW Rev. H or higher

IntelliVue X2 Part Numbers

Please refer to the IntelliVue X2 Service Guide for IntelliVue X2 part numbers.

Plug-in Modules Part Numbers

For inspection procedures; preventive maintenance procedures; cleaning procedures; maintenance, and good practices used to maintain the instrument in good working order, see *Testing and Maintenance*.

| Part number | 12NC | Description |
|-------------|--------------|--------------------------------------|
| M1001-45011 | 453563490691 | M_ECG PLAST SNAP LOCK SINGLE |
| 5040-4249 | 453563099411 | PLAST Housing Module Rear new color |
| 5040-4250 | 453563099421 | PLAST Housing Module Front new color |
| 5040-4247 | 453563099401 | PLAST Snap Lock (CV color) |

M1006B Invasive Blood Pressure Module

| Part number | 12NC | Description |
|-----------------------|--------------|---------------------------------|
| M1006-69601 | 453563463061 | M_PRS IP Mod multilanguage |
| M1006-69603 | 453563462811 | M_PRS IP Mod German |
| M1006-69604 | 453563461771 | M_PRS IP Mod Dutch |
| M1006-69605 | 453563461781 | M_PRS IP Mod Spanish |
| M1006-69607 | 453563461791 | M_PRS IP Mod Norwegian |
| M1006-69608 | 453563461731 | M_PRS IP Mod Swedish |
| M1006-69609 | 453563461741 | M_PRS IP Mod Finnish |
| M1006-69610 | 453563461751 | M_PRS IP Mod Japanese |
| M1006-69613 | 453563461761 | M_PRS IP Mod Chinese |
| M1006-69614 | 453563461271 | M_PRS IP Mod Portugese |
| M1006-69615 | 453563461281 | M_PRS IP Mod Greek |
| M1006-69619 | 453563461291 | M_PRS IP Mod Czech |
| M1006-69620 | 453563461301 | M_PRS IP Mod Polish |
| Front Housings | | |
| M1006-60201 | 453563462101 | M_PRS Frt Housg Assy Press, ENG |
| M1006-60203 | 453563461551 | M_PRS Frt Housg Assy Press, GER |
| M1006-60204 | 453563461581 | M_PRS Frt Housg Assy Press, DUT |
| M1006-60205 | 453563461591 | M_PRS Frt Housg Assy Press, SPA |
| M1006-60207 | 453563461601 | M_PRS Frt Housg Assy Press, NOR |
| M1006-60208 | 453563461611 | M_PRS Frt Housg Assy Press, SWE |
| M1006-60209 | 453563461621 | M_PRS Frt Housg Assy Press, FIN |
| M1006-60210 | 453563461631 | M_PRS Frt Housg Assy Press, JAP |
| M1006-60213 | 453563461641 | M_PRS Frt Housg Assy Press, PRC |
| M1006-60214 | 453563461651 | M_PRS Frt Housg Assy Press, POR |
| M1006-60215 | 453563461661 | M_PRS Frt Housg Assy Press, GRK |
| M1006-60219 | 453563461671 | M_PRS Frt Housg Assy Press, CZE |

| Part number | 12NC | Description |
|---------------------------------------|--------------|---|
| M1006-60220 | 453563461681 | M_PRS Frt Housg Assy Press, POL |
| Overlay old style modules | | |
| M1006-44501 | 453563225431 | M_PRS KBD M1006B OVERLAY English |
| M1006-44502 | 453563225441 | M_PRS KBD M1006B OVERLAY French |
| M1006-44503 | 453563225451 | M_PRS KBD M1006B OVERLAY German |
| M1006-44504 | 453563225461 | M_PRS KBD M1006B OVERLAY Dutch |
| M1006-44505 | 453563225471 | M_PRS KBD M1006B OVERLAY Spanish |
| M1006-44507 | 453563225481 | M_PRS KBD M1006B OVERLAY Swedish |
| M1006-44508 | 453563225491 | M_PRS KBD M1006B OVERLAY Finnish |
| M1006-44509 | 453563225501 | M_PRS KBD M1006B OVERLAY Norwegian |
| M1006-44511 | 453563225511 | M_PRS KBD M1006B OVERLAY Japanese |
| M1006-44512 | 453563225521 | M_PRS KBD M1006B OVERLAY Chinese |
| Modules with analog out | | |
| M1006-69651 | 453563463071 | M_PRS IP Mod multi, w. Analog-out |
| M1006-69653 | 453563461241 | M_PRS IP Mod German, w. Analog-out |
| M1006-69654 | 453563461251 | M_PRS IP Mod Dutch., w. Analog-out |
| M1006-69655 | 453563461261 | M_PRS IP Mod Spanish, w. Analog-out |
| M1006-69657 | 453563460031 | M_PRS IP Mod Norwegian, w. Analog-out |
| M1006-69658 | 453563460041 | M_PRS IP Mod Swedish, w. Analog-out |
| M1006-69659 | 453563460051 | M_PRS IP Mod Finnish, w. Analog-out |
| M1006-69660 | 453563460061 | M_PRS IP Mod Japanese, w. Analog-out |
| M1006-69663 | 453563459011 | M_PRS IP Mod Simpl. Chin, w. Analog-out |
| M1006-69664 | 453563458321 | M_PRS IP Mod Portugese, w. Analog-out |
| M1006-69665 | 453563459021 | M_PRS IP Mod Greek, w. Analog-out |
| M1006-69669 | 453563459031 | M_PRS IP Mod Czech, w. Analog-out |
| M1006-69670 | 453563459041 | M_PRS IP Mod Polish, w. Analog-out |
| Front Housings with analog out | | |
| M1006-60251 | 453563462091 | M_PRS Frt Housg Assy Press C01, ENG |
| M1006-60253 | 453563461691 | M_PRS Frt Housg Assy Press C01, GER |
| M1006-60254 | 453563461701 | M_PRS Frt Housg Assy Press C01, DUT |
| M1006-60255 | 453563461711 | M_PRS Frt Housg Assy Press C01, SPA |
| M1006-60257 | 453563461721 | M_PRS Frt Housg Assy Press C01, NOR |
| M1006-60258 | 453563461821 | M_PRS Frt Housg Assy Press C01, SWE |
| M1006-60259 | 453563461831 | M_PRS Frt Housg Assy Press C01, FIN |
| M1006-60260 | 453563461841 | M_PRS Frt Housg Assy Press C01, JAP |
| M1006-60263 | 453563461851 | M_PRS Frt Housg Assy Press C01, PRC |
| M1006-60264 | 453563461861 | M_PRS Frt Housg Assy Press C01, POR |

| Part number | 12NC | Description |
|----------------------------------|--------------|--|
| M1006-60265 | 453563461871 | M_PRS Frt Housg Assy Press C01, GRK |
| M1006-60269 | 453563461881 | M_PRS Frt Housg Assy Press C01, CZE |
| M1006-60270 | 453563461891 | M_PRS Frt Housg Assy Press C01, POL |
| Overlay old style modules | | |
| M1006-44521 | 453563225531 | M_PRS KBD M1006B #C01 OVERLAY English |
| M1006-44522 | 453563225541 | M_PRS KBD M1006B #C01 OVERLAY French |
| M1006-44523 | 453563225551 | M_PRS KBD M1006B #C01 OVERLAY German |
| M1006-44524 | 453563225561 | M_PRS KBD M1006B #C01 OVERLAY Dutch |
| M1006-44525 | 453563225571 | M_PRS KBD M1006B #C01 OVERLAY Spanish |
| M1006-44526 | 453563225581 | M_PRS KBD M1006B #C01 OVERLAY Italian |
| M1006-44527 | 453563225591 | M_PRS KBD M1006B #C01 OVERLAY Swedish |
| M1006-44528 | 453563225601 | M_PRS KBD M1006B #C01 OVERLAY Finnish |
| M1006-44529 | 453563225611 | M_PRS KBD M1006B #C01 OVERLAY Norwegian |
| M1006-44531 | 453563225621 | M_PRS KBD M1006B #C01 OVERLAY Japanese |
| M1006-44532 | 453563225631 | M_PRS KBD M1006B #C01 OVERLAY Chinese |
| Bezel | | |
| M1006-42202 | 453563463211 | M_PRS PLAST Bezel Press-Round |
| M1006-42201 | 453563469771 | M_PRS PLAST Bezel Press-Square |
| Cables | | |
| M1006-61689 | 453563225911 | M_PRS CBL ADPTR M1006B#C01 1/4" phone JK |

M1011A Intravascular Oxygen Saturation Module

| Part number | 12NC | Description |
|--------------|--------------|---------------------------------------|
| 453564120301 | 453564120301 | M_SO2 Module Exchange, ENG |
| 453564120311 | 453564120311 | M_SO2 Module Exchange, SYMBOLS |
| 453564105801 | 453564105801 | M_SO2 Front Housing Assy SO2 ENGLISH |
| 453564105811 | 453564105811 | M_SO2 Front Housing Assy SO2 SYMBOL |
| M1011-42201 | 453564120341 | M_SO2 Bezel |
| 989803151591 | 989803151591 | SO2 Optical Module |
| 453564097201 | 453564097201 | M_SO2 OEM Exchange SO2 Optical Module |

M1012A Cardiac Output Parameter Module

| Part number | 12NC | Description |
|---------------------------------|--------------|-------------------------------------|
| M1012-69601 | 453563458801 | M_CO Module C.O., ENG |
| M1012-69602 | 453563458761 | M_CO Module C.O., FRC |
| M1012-69603 | 453563458771 | M_CO Module C.O., GER |
| M1012-69605 | 453563458781 | M_CO Module C.O., SPA |
| M1012-69610 | 453563458791 | M_CO Module C.O., JAP |
| M1012-69613 | 453563458831 | M_CO Module C.O., PRC |
| M1012-69615 | 453563460931 | M_CO Module C.O., GRK |
| Front housing | | |
| M1012-60201 | 453563462021 | M_CO Frt Housg Assy C.O., ENG |
| M1012-60202 | 453563460161 | M_CO Frt Housg Assy C.O., FRC |
| M1012-60203 | 453563460171 | M_CO Frt Housg Assy C.O., GER |
| M1012-60205 | 453563460181 | M_CO Frt Housg Assy C.O., SPA |
| M1012-60210 | 453563460191 | M_CO Frt Housg Assy C.O., JAP |
| M1012-60213 | 453563460201 | M_CO Frt Housg Assy C.O., PRC |
| M1012-60215 | 453563460211 | M_CO Frt Housg Assy C.O., GRK |
| Module with Picco | | |
| M1012-69651 | 453563463011 | M_CO Module C10 PiCCO, ENG |
| M1012-69652 | 453563460941 | M_CO Module C10 PiCCO, FRC |
| M1012-69653 | 453563464731 | M_CO Module C10 PiCCO, GER |
| M1012-69655 | 453563460311 | M_CO Module C10 PiCCO, SPA |
| M1012-69660 | 453563460321 | M_CO Module C10 PiCCO, JAP |
| M1012-69663 | 453563460331 | M_CO Module C10 PiCCO, PRC |
| M1012-69665 | 453563460341 | M_CO Module C10 PiCCO, GRK |
| Front housing | | |
| M1012-60251 | 453563462031 | M_CO Frt Housg Assy C.O. PiCCO, ENG |
| M1012-60252 | 453563460221 | M_CO Frt Housg Assy C.O. PiCCO, FRC |
| M1012-60253 | 453563460141 | M_CO Frt Housg Assy C.O. PiCCO, GER |
| M1012-60255 | 453563460151 | M_CO Frt Housg Assy C.O. PiCCO, SPA |
| M1012-60260 | 453563462051 | M_CO Frt Housg Assy C.O. PiCCO, JAP |
| M1012-60263 | 453563462061 | M_CO Frt Housg Assy C.O. PiCCO, PRC |
| M1012-60265 | 453563462041 | M_CO Frt Housg Assy C.O. PiCCO, GRK |
| Overlay old style module | | |
| M1012-44471 | 453563227501 | M_CO KBD M1012A Overlay ENG |
| M1012-44472 | 453563227511 | M_CO KBD M1012A Overlay FRC |
| M1012-44473 | 453563227521 | M_CO KBD M1012A Overlay GER |

| Part number | 12NC | Description |
|--------------|--------------|-----------------------------|
| M1012-44475 | 453563227531 | M_CO KBD M1012A Overlay SPA |
| M1012-44476 | 453563227541 | M_CO KBD M1012A Overlay ITA |
| M1012-44481 | 453563227551 | M_CO KBD M1012A Overlay JAP |
| M1012-44482 | 453563227561 | M_CO KBD M1012A Overlay PRC |
| Bezel | | |
| M1012-42201 | 453563463241 | M_CO PLAST Bezel CO Module |

M1014A Spirometry Module

| Part number | 12NC | Description |
|----------------------|--------------|------------------------------------|
| M1014-69601 | 451261014451 | M_SPR Module Spiro, Eng |
| M1014-69602 | 451261014461 | M_SPR Module Spiro, Symbol |
| Front housing | | |
| M1014-60201 | 451261014491 | M_SPR Frt Housg Assy Spiro, ENG |
| M1014-60202 | 451261014501 | M_SPR Frt Housg Assy Spiro, SYMBOL |
| Bezel | | |
| M1014-42201 | 451261014511 | M_SPR PLAST Bezel Spiro |
| Misc | | |
| M1014-64100 | 451261014521 | M_SPR TUBING Leak Test Kit Spiro |

M1018A Transcutaneous pO₂/pCO₂ Parameter Module

| Part number | 12NC | Description |
|----------------------|--------------|---|
| M1018-69601 | 453563459211 | M_TCPO2 M1018A Exch,Eng,Ger,Dut,... |
| M1018-69602 | 453563460491 | M_TCPO2 M1018A tcpO ₂ /CO ₂ Exch, Fre |
| M1018-69610 | 453563458711 | M_TCPO2 M1018A tcpO ₂ /CO ₂ Exch, Jap |
| M1018-69613 | 453563460541 | M_TCPO2 M1018A tcpO ₂ /CO ₂ Exch, CHN |
| M1018-69614 | 453563460551 | M_TCPO2 M1018A tcpO ₂ /CO ₂ Exch, Por |
| Front housing | | |
| M1018-60201 | 453563461441 | M_TCPO2 Frt Housg Assy ENG |
| M1018-60202 | 453563462211 | M_TCPO2 Frt Housg Assy FRC |
| M1018-60210 | 453563462221 | M_TCPO2 Frt Housg Assy JAP |
| M1018-60213 | 453563462231 | M_TCPO2 Frt Housg Assy PRC |
| M1018-60214 | 453563462241 | M_TCPO2 Frt Housg Assy POR |
| Bezel | | |
| M1018-42201 | 453563463651 | M_TC PLAST Plast Bezel |

| Part number | 12NC | Description |
|-------------|--------------|--------------------------------------|
| M1018-60602 | 453563460501 | M_TC PLAST Calibration Chamber Assy |
| Misc | | |
| M2205A | 989803105991 | Calibration Tube, 5/pouch |
| M1918A | 989803105521 | TCP02/C02 Xducer, solid state design |

M1020B Pulse Oximetry Module

| Part number | 12NC | Description |
|---------------------------------|--------------|--|
| M1020-69651 | 451261000061 | M_SpO2 M1020B Philips/FAST |
| M1020-60251 | 451261000081 | M_SpO2 Frt Housg Assy M1020B PhilipsFAST |
| M1020-66511 | 453563480451 | M_SpO2 SpO2 PCA |
| M1020-42201 | 453563463641 | M_SpO2 PLAST Bezel SpO2 Type 1 Module |
| OxiMax Compatible Module | | |
| M1020-69652 | 451261000101 | M_SpO2 M1020B NELLCOR Oxi |
| M1020-60252 | 451261000121 | M_SpO2 Frt Housg Assy M1020B NELLCOR |
| Masimo SET Module | | |
| M1020-69653 | 451261000131 | M_SpO2 OEM M1020B Masimo |

M1021A Mixed Venous Oxygen Saturation Module

| Part number | 12NC | Description |
|----------------------|--------------|---|
| M1021-69601 | 453563462881 | M_SvO2 M1021A SVO2 Exch,all Lang ex.sCh |
| M1021-69613 | 453563460411 | M_SvO2 M1021A SVO2 Exch,CHN Simpl. |
| Front housing | | |
| M1021-60201 | 453563461461 | M_SvO2 Frt Housg Assy SvO2, ENG |
| M1021-60213 | 453563462151 | M_SvO2 Frt Housg Assy SvO2, PRC |
| Bezel | | |
| M1021-42201 | 453563463621 | M_SvO2 PLAST Plast Bezel |

M1027A EEG Module

| Part number | 12NC | Description |
|-------------|--------------|--|
| M1027-69601 | 453563459151 | M_EEG M1027A Exch, all Lang ex Jap&Gre |
| M1027-69610 | 453563459161 | M_EEG M1027A EEG Exch, JAP |
| M1027-69615 | 453563460481 | M_EEG M1027A EEG Exch, Greek |
| M1027-60201 | 453563461471 | M_EEG Frt Housg Assy EEG, ENG |

| Part number | 12NC | Description |
|-------------|--------------|----------------------------------|
| M1027-60210 | 453563462131 | M_EEG Frt Housg Assy EEG, JAP |
| M1027-60215 | 453563462141 | M_EEG Frt Housg Assy EEG, GRK |
| M1027-42201 | 453563463611 | M_EEG PLAST Plast Bezel |
| M1027-61601 | 453563231141 | M_EEG CBL EEG MODULE TEST DEVICE |

M1029A Temperature Module

| Part number | 12NC | Description |
|---------------------------------|--------------|--|
| M1029-69601 | 453563459291 | M_TMP Temp Module multilanguage |
| M1029-69609 | 453563460581 | M_TMP Temp Module Finish |
| M1029-69610 | 453563460561 | M_TMP Temp Module Japanese |
| M1029-69613 | 453563460571 | M_TMP Temp ModuleChinese Simpl. |
| M1029-69615 | 453563459891 | M_TMP Temp Module Greek |
| M1029-69619 | 453563459901 | M_TMP Temp Module Czech |
| Front housing | | |
| M1029-60201 | 453563461451 | M_TMP Frt Housg Assy Temp, ENG |
| M1029-60209 | 453563462161 | M_TMP Frt Housg Assy Temp, FIN |
| M1029-60210 | 453563462171 | M_TMP Frt Housg Assy Temp, JAP |
| M1029-60213 | 453563462181 | M_TMP Frt Housg Assy Temp, PRC |
| M1029-60215 | 453563462191 | M_TMP Frt Housg Assy Temp, GRK |
| M1029-60219 | 453563462201 | M_TMP Frt Housg Assy Temp, CZE |
| Overlay old style module | | |
| M1029-44471 | 453563231261 | M_TMP KBD M1029A Overlay Multilanguage |
| M1029-44478 | 453563231281 | M_TMP KBD M1029A Overlay Finnish |
| M1029-44481 | 453563231291 | M_TMP KBD M1029A Overlay Chinese |
| M1029-44482 | 453563231301 | M_TMP KBD M1029A Overlay Japanese |
| Bezel | | |
| M1029-42201 | 453563463631 | M_TMP PLAST Temp Module Bezel |

M1032A External Device Interface Module

| Part number | 12NC | Description |
|-------------|--------------|-----------------------------------|
| M1032-69801 | 453563458381 | M_LNK A01 Auxiliary exch module |
| M1032-69802 | 453563458391 | M_LNK A02 Ventilator exch module |
| M1032-69803 | 453563458401 | M_LNK A03 Gas Module exch module |
| M1032-69804 | 453563458411 | M_LNK A04 Anesth mach exch module |
| M1032-69805 | 453563458421 | M_LNK A05 AUX.Plus exch module |


| Part number | 12NC | Description |
|---------------|--------------|---|
| M1032-70806 | 451261002041 | M_LNK Tested PHILIPS SPO2 #A06 |
| M1032-60201 | 453563461401 | M_LNK Frt Housg Kit, all Languages |
| M1032-44401 | 453563231691 | M_LNK KBD Overlay all Languages |
| M1032-42201 | 453563463671 | M_LNK PLAST Vuelink Module Bezel |
| Cables | | |
| 453564086631 | 453564086631 | M_LNK CBL MAQUET JOSTRA HL-20 |
| M1032-61601 | 453563231931 | M_LNK CBL Draeger Cicero (4m) |
| M1032-61602 | 453563231941 | M_LNK CBL Draeger Cato (1.4m) |
| M1032-61603 | 453563231951 | M_LNK CBL Ohmeda Modulus CD (2m) |
| M1032-61604 | 453564044241 | M_LNK CBL Draeger Graph Scr Savina 4m |
| M1032-61605 | 453563231961 | M_LNK CBL NAD Narkomed 2B/2C 2m |
| M1032-61606 | 453564089711 | M_LNK CBL GE DINAMAP PRO 100-400 |
| M1032-61607 | 453564089701 | M_LNK CBL GE Engstroem Carestation |
| M1032-61611 | 453563231971 | M_LNK CBL Free Analog (2m) |
| M1032-61612 | 453563231981 | M_LNK CBL Free Analog (4m) |
| M1032-61613 | 453563231991 | M_LNK CBL Critikon 1846/1846SX (2m) |
| M1032-61614 | 453563232001 | M_LNK CBL Critikon 1846/1846SX (4m) |
| M1032-61615 | 453563232011 | M_LNK CBL Nellcor N-100 (2m) |
| M1032-61616 | 453563232021 | M_LNK CBL Nellcor N-100 (4m) |
| M1032-61617 | 453563498811 | M_LNK CBL Nellcor N-200 (2m) |
| M1032-61618 | 453563232031 | M_LNK CBL Nellcor N-200 (4m) |
| M1032-61619 | 453563232041 | M_LNK CBL ABBOTT OXI 3 |
| M1032-61620 | 453563232051 | M_LNK CBL Abbott Oximetrix 3 (4m) |
| M1032-61621 | 453563232061 | M_LNK CBL Puritan Bennett 7200A/AE (2m) |
| M1032-61622 | 453563232071 | M_LNK CBL Puritan Bennett 7200A/AE (4m) |
| M1032-61623 | 453563232081 | M_LNK CBL Siemens 900 C/D/E (2m) |
| M1032-61624 | 453563232091 | M_LNK CBL Siemens 900 C/D/E (4m) |
| M1032-61625 | 453563498821 | M_LNK CBL Siemens SCM 990 (2m) |
| M1032-61626 | 453563232101 | M_LNK CBL Siemens SCM 990 (4m) |
| M1032-61629 | 453563232111 | M_LNK CBL Draeger Evita (2m) |
| M1032-61630 | 453563232121 | M_LNK CBL Draeger Evita (4m) |
| M1032-61631 | 453563232131 | M_LNK CBL M1025A/B (2m) |
| M1032-61633 | 453563232141 | M_LNK CBL Datex Cap./Ulti.(1.4m) |
| M1032-61635 | 453563232151 | M_LNK CBL Ohmeda RGM 5250 (1.4m) |
| M1032-61636 | 453563232161 | M_LNK CBL Ohmeda RGM 5250 Rev 6 (2m) |
| M1032-61642 | 453563232171 | M_LNK CBL Siemens Servo 300/300A (4m) |
| M1032-61643 | 453563232181 | M_LNK CBL Draeger Evita 2 (4m) |
| M1032-61644 | 453563232191 | M_LNK CBL Draeger Babylog 8000 (4m) |

| Part number | 12NC | Description |
|-------------|--------------|--|
| M1032-61645 | 453563232201 | M_LNK CBL Infrasonic Infant Star/ISV 4m |
| M1032-61648 | 453563232221 | M_LNK CBL Braun FM Sys SW rev.<3.0 (4m) |
| M1032-61649 | 453564113041 | M_LNK CBL BBraun Space System (4m) |
| M1032-61651 | 453563232231 | M_LNK CBL Baxter Explorer (2m) |
| M1032-61652 | 453563232241 | M_LNK CBL Baxter Vigilance(rev.< 6.00)2m |
| M1032-61653 | 453563232251 | M_LNK CBL Mortara ELI 100STM(4m) |
| M1032-61654 | 453563232261 | M_LNK CBL Open Interface 25-pin (4m) |
| M1032-61657 | 453563232271 | M_LNK CBL Bear 1000 (4m) |
| M1032-61658 | 453563232281 | M_LNK CBL Ohmeda 7800 (2m) |
| M1032-61659 | 453563232291 | M_LNK CBL Ohmeda 7810 (2m) |
| M1032-61663 | 453563232311 | M_LNK CBL Support Vuelink TO PC (2m) |
| M1032-61664 | 453563232321 | M_LNK CBL Ohmeda Rascal II (2m) |
| M1032-61665 | 453563232331 | M_LNK CBL NPB 740/760/840 (4m) |
| M1032-61666 | 451261009001 | M_LNK CBL Draeger Zeus (4m) |
| M1032-61667 | 451261011051 | M_LNK CBL I-Stat 1 Analyzer with RIBS 2m |
| M1032-61673 | 453563232351 | M_LNK CBL Hamilton Amadeus/Veolar 4m |
| M1032-61674 | 453563232361 | M_LNK CBL Taema Alis (2m) |
| M1032-61675 | 453563232371 | M_LNK CBL Draeger Cicero EM mono (2m) |
| M1032-61676 | 453563232381 | M_LNK CBL Draeger PM 8050 (2m) |
| M1032-61678 | 453563232391 | M_LNK CBL Ohmeda 7900 |
| M1032-61680 | 453563232401 | M_LNK CBL Draeger Evita XL./4/2dura (2m) |
| M1032-61681 | 453563232411 | M_LNK CBL Draeger Julian (2m) |
| M1032-61682 | 453563232421 | M_LNK CBL Fresenius Vial (4m) |
| M1032-61684 | 453563232431 | M_LNK CBL GE Aesp./Aest./Avance/Aisys 2m |
| M1032-61685 | 453563232441 | M_LNK CBL Draeger Cicero EM Color (2m) |
| M1032-61687 | 453563498831 | M_LNK CBL Aspect BIS A2000 (4m) |
| M1032-61688 | 453563498841 | M_LNK CBL Diametrics IRMA SL (4m) |
| M1032-61689 | 453563498851 | M_LNK CBL Diametrics TrendCare (4m) |
| M1032-61690 | 453563498861 | M_LNK CBL KIT FOR TOF WATCH (2m) |
| M1032-61691 | 453563232451 | M_LNK CBL BRAUN FM System SW rev. >= 3.1 |
| M1032-61692 | 453563498871 | M_LNK CBL Edwards Vigilance/2/Vigileo 2m |
| M1032-61693 | 453563458271 | M_LNK CBL Danm. AAI (SSI)/AEP Monitor 2m |
| M1032-61694 | 453563458281 | M_LNK CBL DraegerPrimus/Apollo/Pallas 2m |
| M1032-61699 | 451261003891 | M_LNK CBL Open Interf 9-pin dig. only 4m |
| M1032-61695 | 453563458261 | M_LNK CBL FOR DRAEGER VAMOS |
| M1032-61700 | 451261013521 | M_LNK CBL Draeger Fabius GS / Tiro (2m) |
| M1032-61696 | 453563498881 | M_LNK CBL Maquet (Siemens) Servo-i/-s 4m |
| M1032-9100L | 451261001271 | M_LNK LP Inst. Note TOF-Watch SX Mnt-eng |

M1034A BIS Module


| Part number | 12NC | Description |
|-------------|--------------|--------------------------------------|
| M1034-69601 | 453563462841 | M_BIS BIS Interface Module, exchange |
| M1034-61630 | 453563233721 | M_BIS CBL Patient Interface Cable |
| M1034-68520 | 453563233761 | M_BIS BIS Engine, exchange |
| M1034-68102 | 453563233741 | M_BIS BIS DSC-XP, exchange |
| M1034-61610 | 453563233701 | M_BIS CBL Engine Cable |
| M1034-61620 | 453563233711 | M_BIS CBL Engine Cable (long) |
| M1034-68521 | 451261003621 | M_BIS BISx Power Link, exchange |
| M1034-61660 | 451261005261 | M_BIS CBL BISx Host Cable |
| M1034-47600 | 451261005271 | M_BIS CBL BISx Bulkhead Connector |
| M1034-41200 | 453563233641 | M_BIS Rack Mount |
| M1034-60104 | 453563490591 | M_BIS MECHASY BIS FMS MOUNT |
| M1034-60201 | 453563461411 | M_BIS Frt Housg Assy BIS, ENG |
| M1034-61650 | 453563233731 | M_BIS BIS Sensor Simulator |
| M1034-42201 | 453563463661 | M_BIS PLAST PLAST Bezel |
| M1034AX | 862375 | BISx |

BIS Solution Components



- 1 BIS Engine Cable
- 2 BIS Module
- 3 BIS Engine
- 4 Digital Signal Converter (DSC)
- 5 Patient Interface Cable
- 6 BIS Sensor

BISx Solution Components



- 1 BIS Module
- 2 BISx
- 3 Patient Interface Cable
- 4 BIS Sensor
- 5 BIS Sensor
- 6 BIS Interface Board (MP20/30 only)
- 7 BISx

M1116B Thermal Array Recorder Module

| Part number | 12NC | Description |
|----------------------|--------------|---|
| M1116-68603 | 453563466701 | 2" RECDR ECMS EXCH ENGLISH |
| M1116-68604 | 453563466711 | 2" RECDR ECMS EXCH GERMAN |
| M1116-68605 | 453563466721 | 2" RECDR ECMS EXCH SPANISH |
| M1116-68606 | 453563466731 | 2" RECDR ECMS EXCH FRENCH |
| M1116-68607 | 453563466741 | 2" RECDR ECMS EXCH SWEDISH |
| M1116-68608 | 453563466751 | 2" RECDR ECMS EXCH CHINESE |
| M1116-68609 | 453563466761 | 2" RECDR ECMS EXCH JAPANESE |
| M1116-68610 | 453563466771 | 2" RECDR ECMS EXCH ITALIAN |
| M1116-68620 | 453563466781 | 2" RECORDER ECMS EXCH FINNISH |
| M1116-68621 | 453563466791 | 2" RECORDER ECMS EXCH PORT |
| Front housing | | |
| M1116-60203 | 453563462301 | PLAST Front Housing English |
| M1116-60204 | 453563489221 | PLAST Front Housing German |
| M1116-60205 | 453563462311 | PLAST Front Housing Spanish |
| M1116-60206 | 453563462321 | PLAST Front Housing French |
| M1116-60207 | 453563489231 | PLAST Front Housing Swedish |
| M1116-60208 | 453563462331 | PLAST Front Housing Chinese |
| M1116-60209 | 453563462341 | PLAST Front Housing Japanese |
| M1116-60210 | 453563462351 | PLAST Front Housing Italian |
| M1116-60220 | 453563489241 | PLAST Front Housing Finnish |
| M1116-60221 | 453563489251 | PLAST Front Housing Assembly Portuguese |
| M1116-60301 | 453563243911 | HSG Rear DBL (CV COLOR |
| M1116-80040 | 453563244171 | PLAST Left Half Chassis |
| M1116-80230 | 453564017111 | PLAST Right Half Chassis |
| 5040-4255 | 453563099431 | PLAST Housing Double Module |
| 5040-4256 | 453563099441 | PLAST Housing Double Module |
| M1116-40041 | 453563243811 | Housing Top (CV color) |
| M1116-40230 | 453563243851 | PLAST Front Housing English |
| M1116-40240 | 453563243861 | Front Door (CV color) |
| 1500-0802 | 453563059491 | MECH Flat Belt |
| 1810-1339 | 453563064911 | ASSY Thermal Print Head |
| M1116-00030 | 453563243781 | MET Spring Leaf |
| M1116-40060 | 453563243821 | MACH Timing Pulley |
| M1116-40070 | 453563243831 | MECHASY Drive Roller |
| M1116-40095 | 453563243841 | SWITCH CAP EXTENDER |

| Part number | 12NC | Description |
|-------------|--------------|---|
| M1116-41050 | 453563243871 | STAMPING Buff Cam |
| M1116-60200 | 453563243881 | EMCH Drive Motor |
| M1116-60201 | 453563243891 | Printhead Cleaning Kit |
| M1116-83002 | 453563244191 | LBL Inside Door |
| 0515-0890 | 453563480701 | SCRW MACH M3 X 0.5 6MM-LG 90-DEG-FLH-HD |
| M1001-60620 | 453563224011 | M_ECG PLAST RECORDER FRONT ASSY, BLANK |

865115 IntelliBridge EC10 Module

| Part number | 12NC | Description |
|--------------|--------------|---|
| 453564116661 | 453564116661 | IB-EC10 Module RS232/LAN |
| 453564116931 | 453564116931 | IB-EC10 Frt Housg Assy Module RS232/LAN |
| M1031-42201 | 453564116891 | IB-EC10 Bezel RJ45 |

865383 NMT Module

| Part number | 12NC | Description |
|--------------|--------------|--------------------------------|
| 453564279141 | 453564279141 | M_NMT Module NMT, US |
| 453564279161 | 453564279161 | M_NMT Module NMT, Intl |
| 453564257851 | 453564257851 | M_NMT PLAST Bezel |
| 453564279181 | 453564279181 | M_NMT Frt Housg Assy NMT, US |
| 453564279191 | 453564279191 | M_NMT Frt Housg Assy NMT, Intl |
| 989803174581 | 989803174581 | CBL NMT Patient Cable* |

* Not orderable via SPS. Must be ordered through Supplies.

External Display Part Numbers

| Part number | 12NC | Description |
|---------------|---------------|--|
| M8031B | 862137 | 15in TFT Touch XGA Display (FIMI) |
| M8031-68001 | 451261001941 | IV DSPL 15" Medical Grade w Touch |
| M8031-60005 | 451261001921 | IV ASSY-PWR 12V for M8031B Display |
| M8031-64001 | 451261001931 | IV ASSY Pwr Supply Holder 15"/19" displ |
| M8033-04701 | 453563480981 | IV DSPL Desk Stand for M8033A/M8033B |
| 2090-0985 | 451261014381 | DSPL Backlights (2x) (Version 2) |
| 2090-0860 | 453563463201 | Backlights (2x) (Version 1) |

| Part number | 12NC | Description |
|---------------|---------------|---|
| M8033C | M8033C | 17 in. TFT Touch (S)XGA Display |
| M8033-68071 | 451261009161 | IV DSPL 17" Medical Grade with Touch |
| M8031-04701 | 451261001901 | IV MECHASY Mon Desk Stand M8031B/M8033C |
| M8033-64603 | 451920880311 | IV DSPL Backlights for M8033C |
| 453564263911 | 453564263911 | IV MECHASY 17" display replacement kit ¹ |
| 865299 | 865299 | 19 in. TFT Touch (S)XGA Display |
| 453564116741 | 453564116741 | ASSY-PWR - E539821 Power Brick - 12V |
| 453564192181 | 453564192181 | IV DSPL 19" Medical Grade w Touch 865299 |
| 451261001931 | 451261001931 | IV ASSY Pwr Supply Holder 15"/19" displ |

¹ Must be ordered together with 453564192181 IV DSPL 19" Medical Grade w Touch 865299

Test and Service Tools

Table 4 UTP LAN Crossover

| Part number | 12NC | Description |
|-------------|--------------|------------------------------|
| M3199-60101 | 453563337371 | PIC CBL NI 3FT CROSSOVER UTP |
| M3199-60102 | 453563337381 | NI CBL UTP Crossover 12ft |

Table 5 Grounding

| Part number | 12NC | Description |
|-------------|--------------|--------------------------------------|
| 8120-4808 | 453563199211 | CMS CBL EXT GND ASSY (gnd lug con) |
| 8120-2961 | 453563198651 | CMS CBL EXT GND ASSY(Crocodile clip) |

Table 6 Test Cables

| Part number | 12NC | Description |
|--------------|--------------|--|
| 451261026081 | 451261026081 | CBL Safety Test ECG |
| 451261026141 | 451261026141 | CBL Safety Test IBP |
| 451261026071 | 451261026071 | CBL Safety Test M1006A/B #C01 Phone Jack |
| 451261026041 | 451261026041 | CBL Safety Test SpO2 (MMS/M1020B) |
| 451261026171 | 451261026171 | CBL Safety Test C.O. |
| 451261026091 | 451261026091 | CBL Safety Test EEG |
| 451261026161 | 451261026161 | CBL Safety Test CO2 |
| 451261026131 | 451261026131 | CBL Safety Test Temp |
| 451261026031 | 451261026031 | CBL Safety Test CO2 (M3014A) |

| Part number | 12NC | Description |
|--------------|--------------|---------------------------------------|
| 453564127781 | 453564127781 | CBL Safety Test IntelliBridge |
| 451261026181 | 451261026181 | CBL Safety Test VueLink I |
| 451261026191 | 451261026191 | CBL Safety Test Defi paddles |
| 453564127771 | 453564127771 | CBL Cable AY. ScvO2 |
| M1012-61601 | 453563227731 | M_CO CONN Test Adapter Cardiac Output |

Table 7 Capnography - Respironics Mainstream

| Part number | 12NC | Description |
|-------------|------|------------------------|
| M2506A | | GAS Verification Gas |
| M2505A | | GAS CYLINDER REGULATOR |

Table 8 Capnography - Respironics Sidestream

| Part number | 12NC | Description |
|--------------|--------------|---------------------------------|
| M2267A | 989803106081 | Calibration Regulator |
| 13907A | 989803100361 | Calibration Tube Assembly |
| M1026-60144 | 453563230731 | M_AGM Electronic Mass FlowMeter |
| 453564178121 | 453564178121 | M_AGM Flowmeter TSI |
| M2506A | | GAS Verification Gas |
| M2776A | | Straight Sample Line |

Table 9 Capnography - Microstream Oridion

| Part number | 12NC | Description |
|--------------|--------------|---|
| M1026-60144 | 453563230731 | M_AGM Electronic Mass FlowMeter |
| 453564178121 | 453564178121 | M_AGM Flowmeter TSI |
| 15210-64010 | 989803100841 | GAS Cal 1 cylinders for TCPC02, 6/bx. |
| 15210-64020 | 989803100851 | GAS Cal 2 Cylinders for TCPC02, 6/bx. |
| M2267A | 989803106081 | Calibration Regulator |
| M3015-47301 | 989803143081 | MS_SCO2 LBSPLY Calibration Line |
| 13907A | 989803100361 | Calibration Tube Assembly |
| M1013-64002 | 451261014851 | M_G1 IV G1/G5 Leakage Test Kit (only required for leakage test without M1026-60144 flowmeter) |

Installation Instructions

The information contained in this chapter should enable the MX600/MX700 to be installed ready for use (the preparation and planning should be adhered to as specified in the *Site Preparation* section). Configuration of the system is explained in the Configuration Guide.

Installation should be carried out by qualified service personnel, either by the hospital's biomedical department, or by Philips Support.

The monitor is suitable for use in all medically used rooms which fulfill the requirements regarding electrical installation according to IEC60364-7-710 "Requirements for special installations or locations- Medical locations, or corresponding local regulations.

For mechanical and electrical installation, you need technically qualified personnel with a knowledge of english. Additionally, for monitor configuration, you need clinically qualified personnel with a knowledge of the use environment.

As the first step in preparing the monitor for use, follow the installation instructions given in this chapter.

Electromagnetic Emissions

The monitor is suitable for use in the electromagnetic environment specified in the table below. You must ensure that it is used in such an environment.

| Emissions test | Compliance | Avoiding Electromagnetic Interference |
|------------------------------------|------------|--|
| Radio Frequency (RF) emissions | Group 1 | The monitor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment |
| RF emissions CISPR 11 | Class A | The monitor is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes |
| Harmonic emissions IEC 61000-3-2 | n/a | |
| Voltage fluctuations IEC 61000-3-3 | n/a | |

WARNING

The monitor should not be used next to or stacked with other equipment. If you must stack the monitor, you must check that normal operation is possible in the necessary configuration before you start monitoring patients.

Electromagnetic Interference (SRR)

Commercially available Short Range Radio 802.15.4 transceivers operate at very low RF power levels to transmit data and need to have high sensitivity receivers to achieve a good link budget. Due to technological limitations the selectivity of the receiver is limited. Consequently, the SRR link is susceptible to other strong RF transmitters not only in the operating frequency band and 5% around it, but also to non-transient RF disturbances stronger than 1V/m at frequencies close to the operating frequency band (2.0 to 2.3 GHz)

Installation Checklist

Use this checklist to document your installation. Please file this installation record

| Step | Task | Check Box when Task Done |
|------|---|--------------------------|
| 1 | Perform initial inspection of delivery, unpack and check the shipment | <input type="checkbox"/> |
| 2 | Mount the monitor as appropriate for your installation | <input type="checkbox"/> |
| 3 | Connect the monitor to AC mains using the supplied power cord | <input type="checkbox"/> |
| 4 | Perform Visual, Power On and Functional test blocks | <input type="checkbox"/> |
| 5 | Perform Safety Tests, if required by local laws and regulations | <input type="checkbox"/> |
| 6 | Load paper into the recorder, if present | <input type="checkbox"/> |
| 7 | Check/set the time and date | <input type="checkbox"/> |
| 8 | Check that the country-specific default settings are appropriate | <input type="checkbox"/> |
| 9 | Perform System Test as necessary | <input type="checkbox"/> |

Unpacking the Equipment

Your equipment will arrive in a carton similar to the ones pictured below. All components of the monitoring system are consolidated into a single packing crate. The user documentation is provided in a separate package. The contents of the monitoring system crate depend on the options you have purchased. In addition to the monitor it can contain the following:

- MMS and user manuals
- FMS
- Parameter modules
- MMS Extensions and accessories



Accessory Packaging (Remove upper box to reveal monitor)



Monitor Packaging

In the unlikely event of a defect on arrival, please keep the packing materials until you have completed the initial inspection.

Initial Inspection

Mechanical Inspection

Open the shipping container(s) and examine each part of the instrument for visible damage, such as broken connectors or controls, or scratches on the equipment surfaces. If the shipping carton/ container is undamaged, check the cushioning material and note any signs of severe stress as an indication of rough handling in transit. This may be necessary to support claims for hidden damage that may only become apparent during subsequent testing.

Electrical Inspection

The instrument has undergone extensive testing prior to shipment. Safety testing at installation is not required (except in situations where devices are interconnected forming a system, see *Connecting Non-Medical Devices* in the Site Preparation Chapter of this manual). An extensive self check may be performed. This recommendation does not supersede local requirements.

All tests are described in the *Testing and Maintenance* section of this manual.

Claims for Damage and Repackaging

Claims for Damage

When the equipment is received, if physical damage is evident or if the monitor does not meet the specified operational requirements of the patient safety checks or the extended self check, notify the carrier and the nearest Philips Sales/Support Office at once. Philips will arrange for immediate repair or replacement of the instrument without waiting for the claim settlement by the carrier.

Repackaging for Shipment or Storage

If the instrument is to be shipped to a Philips Sales/Support Office, securely attach a label showing the name and address of the owner, the instrument model and serial numbers, and the repair required (or symptoms of the fault). If available and reusable, the original Philips packaging should be used to provide adequate protection during transit. If the original Philips packaging is not available or reusable please contact the Philips Sales/Support Office who will provide information about adequate packaging materials and methods.

Installing the MX600/700 Monitor

NOTE

There are different mounting options available for the monitor. This section covers the general concepts of safe mount installations and specific steps for the mounting options sold by Philips. Instructions which ship with a mounting solution should always take precedence over the instructions described in this chapter.

You **MUST** follow the instructions that ship with the mounting solution, regardless of manufacturer.

Mounting Instructions

Assembling Mounts

- The table mount ships with the monitor. Every type of compatible mounting solution is delivered with a complete set of mounting hardware and instructions. Refer to the documentation delivered with the mounting hardware for instructions on assembling mounts. Compatible table mounts are:
 - M8000-64100 (table mount with gray tops and marked with date code 10/31 or higher), or
 - 453564239731 (completely made of stainless steel).


WARNING

- It is the customer's responsibility to have the attachment of the mounting hardware to the ceiling, wall, or mounting rail and the construction of the ceiling, wall, or mounting rail evaluated for structural integrity and compliance with all local, state and any other required codes by a registered, professional, structural and/or mechanical engineer.
 - Ensure that this commitment has been met before assembling mounts.
 - Incorrect mounting and use of inappropriate mounting material may lead to injury. It is the customer's responsibility to ensure that the mounting procedures have been performed correctly and the appropriate mounting devices have been used.
-

Please mount the monitor using either the Philips Quick Mount or Fix Mount solution or another approved mounting solution. The mounting shall be done in a manner that no patient, operator or other person can be harmed by a monitor removed intentionally or released accidentally from the mount. When using the Quick Mount, be aware of the danger of accidental activation of the Quick Mount release button when lifting or moving items located under the monitor, such as pole mounts, etc. If in doubt, use the Philips Fix Mount solution to avoid such situations.

For instructions on how to mount the monitor using the Quick Mount table mount refer to the Assembly Instructions delivered with the mounting kit M8000-64100 or 453564239731.

An alternative mounting solution is the Fix Mount, which should be used for those installations where the Quick Mount might not be appropriate. The Fix Mount is already integrated into the monitor bottom housing.



7 Installation Instructions





Table Mount (M4046-64100, 12NC: 451261001381)

| | A | (B) | C1 | C2 |
|-----------|--------|------|--------|---------|
| MX600/700 | 29.5mm | 57mm | 93.1mm | 164.7mm |


Mounting:

The Quick Mount and Fix Mount solutions are both shipped as standard.



Installing the FMS-4 Cable Management Hooks (without MMS Mount)

- 1 Remove the three screws (T20) from the bottom of the FMS.



- 2 Remove the bottom cover.



7 Installation Instructions

- 3 Insert the cable management hooks in the desired position




- 4 Reinsert the bottom cover at the back first and then flip it down.



- 5 Reinsert and tighten the three screws on the bottom cover.




- 6 You can install up to three cable management hooks.



Installing the FMS-4 Cable Management Hooks (with MMS Mount)


- 1 Remove the three screws (T20) from the bottom cover of the FMS.



- 2 Remove the bottom cover.




- 3 Install the cable management hooks in the desired position.




7 Installation Instructions

- 4 Reinsert the bottom cover at the back first and then flip it down.



- 5 Reinsert and tighten the three screws in the bottom cover.



- 6 You can install as many cable management hooks as you wish.



Installing the Wired Remote Control


NOTE

Connection of a remote control requires SW Rev H.0 or higher
The remote control comes with a cradle and USB cable.




To install the remote control:

- 1 Remove the two screws on the back of the remote control.




- 2 There are two ways to install the USB cable - at the top of the remote control or at the bottom. If you want to install the cable at the bottom of the remote control, proceed to step 3. To install the cable at the top of the remote control:
 - a. Plug the USB cable and lead it through the remote control housing as shown below.





7 Installation Instructions


- b. Remove the plastic cover from the top of the remote control rear housing and reinsert it in the corner of the housing as shown below.



- c. Break open the rubber seal for the USB cable at the perforation as shown below.



- d. Insert the USB cable into the rubber seal.



- e. Insert the rubber seal into the cable slot of the remote control housing.



NOTE

Be sure to always insert the rubber seal when installing the remote control, as it prevents liquid from running into the remote control.

- f. Put the housing of the remote control back together and reinsert the screws.





- g. Perform a functional test to make sure the remote control is functioning correctly.

3 To install the cable at the bottom of the remote control:

- a. Plug the USB cable and lead it through the remote control housing as shown below.





- b. Remove the plastic cover from the bottom of the remote control rear housing and reinsert it in the corner of the housing as shown below.




7 Installation Instructions


- c. Break open the rubber seal for the USB cable at the perforation as shown below.



- d. Insert the USB cable into the rubber seal.



- e. Insert the rubber seal into the cable slot of the remote control housing.



NOTE

Be sure to always insert the rubber seal when installing the remote control, as it prevents liquid from running into the remote control.

- f. Put the housing of the remote control back together and reinsert the screws.



- g. Perform a functional test to make sure the remote control is functioning correctly.

Mounting the Wired Remote Control


The Remote Control can either be mounted to a wall or with the Universal Mounting Clamp.

To mount the remote control to a wall:

- 1 Remove the screw from the remote control holder.



- 2 Separate the back plate from the remote control holder.



7 Installation Instructions


- 3 Insert the USB cable into the holder on the back plate at the desired length.



- 4 Screw the holder back plate to a wall and then reattach the holder to the back plate.

To mount the remote control using the Universal Mounting Clamp:

- 1 Make sure the back plate is attached to the cradle
- 2 Screw the Universal Mounting Clamp to the back of the cradle in the desired direction.



- 3 Insert the USB cable into the holder on the back plate at the desired length.

Installing the Wireless Remote Control

In order to use the wireless remote control with the monitor, the remote control must be assigned to the monitor first. To assign a remote control to the monitor:

- 1 Go into Service Mode or Configuration Mode.
- 2 Press the back key for more than two seconds. This initiates the remote control discovery procedure.



The **Add Cableless** key appears instead of the measurement selection key.

- 3 Select the **Add Cableless** key to open the **Add Cableless** window. The Remote Control is shown in the window with a symbol and its label.



- 4 Select the remote control in the window. The monitor displays the assignment prompt message: **cl RC added <clRC serial number>** in the Status Line of the monitor.

NOTE

- While the remote control discovery procedure is active, the functionality of the remote control is disabled.
- If there is already a remote control assigned to the monitor, this remote control must be unassigned before a new remote control is assigned.

To remove a remote control:

- 1 Select **Main Setup** → **Hardware** → **cl Remote Control**
- 2 Select **Remove RemCtrl**
- 3 Select **Confirm** to unassign the remote control

CAUTION

When using a remote control without a cable, it is important that the user knows which remote control is assigned to which monitor. Use the tethering cable delivered with the remote control to attach it to a bed rail or IV pole, or label the remote control with the bed or monitor ID.


Wireless Remote Control Assignment Rules

- Only one remote control can be assigned to a monitor at a time. If you want to assign another remote control, you must unassign the remote control which is already assigned first.
- If you switch from wireless remote control to USB connected remote control operation with the same remote control, the assignment to the monitor will be lost. If you switch back to wireless remote control operation you will have to reassign the remote control.
- If you assign a wireless remote control to a second monitor without unassigning it from the first monitor and then come back to the monitor it was originally assigned to, press the back key on the remote control to start discovery mode. Once the remote control has been recognized by the monitor you can reassign it immediately. It does not have to be removed first.
- An assignment can only be initiated from the remote control itself.
- An unassignment can only be initiated from the monitor.
- If you exchange the battery of a remote control the assignment to a monitor is kept. You do not need to reassign the remote control.

Installing the Tethering Cable

To attach the tethering cable to a remote control:

- 1 Pull apart the ring on one end of the tethering cable.



- 2 Insert the end of the metal ring spiral into the hole at the top left corner of the remote control.




Exchanging the Remote Control Battery

- 1 Loosen the two screws on the back of the remote control and open the remote control.



7 Installation Instructions

- 2 Remove the old batteries and insert the new batteries as shown below.



- 3 Close the remote control and tighten the screws on the back cover.
- 4 Perform a functional test to make sure the remote control is functioning correctly.

Connecting the Monitor to AC Mains

The monitor has a wide-range power supply that allows you to operate the monitor from an AC (alternating current) power source of 100 V to 240 V ($\pm 10\%$) and 50/60 Hz ($\pm 5\%$).

WARNING

- Always use the supplied power cord with the earthed mains plug to connect the monitor to an earthed AC mains socket. Never adapt the mains plug from the power supply to fit an unearthed AC mains socket.
 - Do not use AC mains extension cords or multiple portable socket-outlets. If a multiple portable socket-outlet without an approved isolation transformer is used, the interruption of its protective earthing may result in enclosure leakage currents equal to the sum of the individual earth leakage currents, so exceeding allowable limits.
 - Do not connect any devices that are not supported as part of a system.
 - Any non-medical device placed and operated in the patient's vicinity must be powered via an approved isolation transformer that ensures mechanical fixing of the power cords and covering of any unused power outlets.
-

Connections

The following figure shows the cable and interface board connections:



| | |
|--|---|
| 1 | Serial/MIB (RS232) connectors, type RJ45, Flexible Nurse Call |
| 2 | AC power input |
| 3 | Protective earth screw hole |
| 4 | Equipotential ground connector |
| 5 | Measurement Server Link (MSL) |
| 6 | USB rear connectors (for remote control, keyboard, pointing devices, printer) |
| 7 | Serial RS232 connector |
| 8 | Nurse Call |
| 9 | Wired network connector |
| 10 | Video out connector (digital/analog) |
| The following connectors are only present with the iPC | |
| 11 | USB rear connectors (for keyboard, pointing devices, printer) |
| 12 | Audio in/out |
| 13 | Local Area Network |
| 14 | Video out connector (digital/analog) |

NOTE

For installation of software on the iPC, refer to the documentation provided with the software you want to install.



- 1 ECG Sync Output/Analog ECG output connector



- 1 USB side connector (only present with the iPC)

Video Interface on the Connector Board

NOTE

The DVI video interface on the connector board has slave display capability only. Slave displays must have the same resolution as the MX600/700's built-in display. If you connect a slave display with a different resolution, you may see distortion or black bars on the edge of your screen.

VESA-compatible displays can be connected to the DVI video interface on the connector board. In rare cases, a specific VESA-compatible display may not be supported, because the video timing of the DVI video interface differs slightly from the VESA standard. We recommend testing the functionality of any external display before installation.

Audio Interface (for iPC only)

The audio interface does not provide an electrical separation. When connecting an audio device which uses an additional power supply e.g. an active speaker, a separation device according to EN/IEC 606010101 or IEC 60601-1 edition 3 clause 16 is required. After installation, a system test is required.

to ensure that the resulting equipment leakage current does not exceed the limits of EN/IEC 60601-1-1 or IEC 60601-1 edition 3 clause 16. For detailed information see the Testing and Maintenance chapter.

Connection of Devices via the RS232 Connector on the Connector Board

The following devices can be connected to the RS232 connector on the connector board:

- IntelliVue G1 Anesthetic Gas Module
- IntelliVue G5 Anesthetic Gas Module
- M1026B Anesthetic Gas Module
- Barcode Reader
- Touchscreen

NOTE

The RS232 connector on the connector board has no data export capability

Connection of Devices via the MIB/RS232 I/O Board

The configuration of a specific MIB/RS232 port can be viewed in config mode and altered in service mode. This is required, for example, when a slave display with touchscreen is installed. To alter the configuration of an MIB/RS232 port select **Main Setup** then **Hardware** then **Interfaces**. You can configure **GM**, **Touch 1** and **Touch 2** to the MIB/RS232 port.

Data Out can be configured up to two times (on two or more MIB/RS232 boards). Note that only the first MIB/RS232 port configured to **Data Out** (i.e. the first one to receive a request) provides wave export. A second MIB/RS232 port configured to **Data Out** will only export numerics.

NOTE

Be aware that if you change a port assignment this assignment is not reset upon boot up. If the MIB/RS232 board is removed and replaced with a different type of board the settings are deleted. If the MIB/RS232 board is then refitted, you must reconfigure the MIB/RS232 port. The configuration of MIB/RS232 is not cloned between monitors. **GM** can only be configured **once** to an MIB/RS232 port.

Connection of USB Devices

The USB ports on the connector board support the following USB devices:

- printer
- keyboard
- barcode scanner
- computer mouse or trackball
- remote control


NOTE

- Connect only the above mentioned devices to the USB Interface. Other devices are not supported.
- Connection of a remote control requires SW Rev H.0 or higher
- A remote control is treated as a keyboard. Either a remote control or a keyboard can be connected to one group.


CAUTION

Do NOT connect a remote control to a USB connector on the iPC (including the side USB Connector) as it will not function properly. To make sure you have connected the remote control to the connector board and not the iPC, press the SmartKeys button on the remote control. The Smartkeys window should appear. If it does not, check the USB connection of the remote control.

As the patient monitor software only supports two input devices, only two input devices can be connected to the USB ports on the connector board. For this purpose, the USB ports are divided into two groups, "A" and "B". Only one input device per group is allowed. In the graphic below, a mouse is connected to a port in each group. Therefore the keyboard is not recognized.



It is, however, possible to connect a mouse, a keyboard and a barcode scanner. In this case, the keyboard and barcode scanner are treated as one input device and must be connected to two ports of the same group.



NOTE


Other USB devices, e.g. USB sticks, iPods etc. are not supported by the USB IF board. **Do not** use USB adapters to connect PS/2 or other devices to the USB board.

Possible USB Device Combinations at the Connector Board


The following three groups of USB device combinations are supported:

- mouse/mouse combination
- keyboard/mouse combination
- keyboard/keyboard combination.


The figures below show the possible device combinations for each of the three groups.



Possible Mouse/Mouse combinations



Possible Mouse/Keyboard combinations



Possible Keyboard/Keyboard combinations

NOTE

Every time a new device is connected to the USB ports, all connected USB devices are stopped and the ports are scanned. Depending on whether the combination is allowed or not, the devices will function again after the scan.

Combined input devices such as a keyboard with an integrated trackball are also supported. However, no additional mouse can be connected in this case. Multiple combined devices are also not allowed.

CAUTION

The USB ports do not provide an electrical separation. When connecting a USB device which uses an additional power supply e.g. printer, a separation device according to EN/IEC 60601-1-1 is required. After installation a system test is required to ensure that the resulting equipment leakage current does not exceed the limits of EN/IEC 60601-1-1. For detailed information see the Testing and Maintenance chapter.

Possible USB Device Combinations at the iPC

The iPC is treated as a regular PC device. There are no specific rules for the connection of USB devices.

You can connect input devices to the rear connectors of the iPC or to the front USB connector and use them for the iPC. Input devices connected to the iPC can also be shared with the patient monitor. Sharing input devices and touchscreen input requires a specific software application; refer to the *IntelliVue XDS Application Instructions for Use* for details.

If the iPC is switched off when input devices are shared, any input devices connected to it will not be available for use with the monitor.

Configuring the USB Interface for Use with a Programmable Barcode Scanner

Requires SW Rev. H.0 or Higher If you want to use a programmable barcode scanner, the USB Interface group you want to connect the barcode scanner to must be configured to Barcode.

- 1 Go into Service Mode.
- 2 Select **Hardware -> Interfaces->USB**
- 3 Select the **Change Setting** key
- 4 Change the setting to **Barcode** and select the **Done** key.

NOTE

A programmable barcode scanner must be pre-configured to provide the codes in the correct form to the monitor. Refer to the Installation Note "Installing and Testing the 4600g Barcode Scanner" on the IntelliVue Documentation DVD for details.

Installing Remote Devices

This section provides instructions for Philips products. Installation instructions for devices not sold by Philips must be provided by the device manufacturer.

Mounting the 15" Remote Display (M8031B)


Mounting solutions for the M8031B must be purchased separately. Please refer to the installation instructions which ship with the mounting solution purchased.

Connections

Connect the cables to the display as shown in the photographs below.



Video Cable Wiring Schematics




| CIRCUIT DIAGRAM | | |
|-----------------|-------------------|-------|
| P1 | WIRE | P2 |
| 1 | RED COAX, CENTER | 1 |
| 2 | GREY COAX, CENTER | 2 |
| 3 | BLUE COAX, CENTER | 3 |
| 5 | BLACK | 5 |
| 10 | | 10 |
| 6 | RED COAX, SHIELD | 6 |
| 7 | GREY COAX, SHIELD | 7 |
| 8 | BLUE COAX, SHIELD | 8 |
| 11 | BROWN | 11 |
| 13 | YELLOW | 13 |
| 14 | WHITE | 14 |
| SHELL | OVER SHIELD | SHELL |

Analog Video Cable Wiring Schematic

Hardware Settings

- This section lists settings grouped in the Hardware Settings Block which are available in Service Mode. These settings are set once per monitor and are the same in every profile. Any changes you make to the hardware settings configuration are automatically stored, there is no need to save them in an extra step. Hardware settings must be entered for each monitor individually, they are stored in the monitor, and they are not cloned. To enter the hardware settings menu, select **Main Setup** - > **Hardware**.
- Keyboard - this setting allows you to select the language of the keyboard connected via USB. See USB Keyboard/Mouse *section* in this chapter for details.
- Setup Video - this setting allows you to set the correct display resolution for the external displays and to choose between the basic (black background) or classic (photo of nurse and child in the background) standby/boot screen.




- MIB/RS232 - see “Connection of Devices via the MIB/RS232 I/O Board” on page 285 section in this chapter for details.

Flexible Module Rack and/or Multi-Measurement Module

Attaching the MMS to a Mount

- 1 Make sure the Measurement Server is oriented correctly relative to the mount (see the picture below).
- 2 Place the Measurement Server on the back mount. If it is not tight against the mount, slip it in the direction of the measurement connectors until it is.
- 3 Slip the Measurement Server forward until it clicks into place.




Detaching the Measurement Server from a Mount

- 1 Press and hold the latch (in the middle at the top of the mount) away from the Measurement Server.
- 2 Slide the Measurement Server off the mount in the direction of the measurement connectors.

Positioning the Measurement Server on a Clamp Mount

If you have your Measurement Server on the clamp mount, you can have it in one of four positions. You can reposition it as follows:

- 1 Press and hold the mount latch toward the clamp screw.



Rotate the Measurement Server and mount until you get it to the position you want.



- 2 Release the mount latch, and make sure it is clicked into one of the four slots on the back of the mount.

Mounting the MMS Mount to the FMS (M8048A)

- 1 Connect the MMS Mount to the FMS and snap it into place.



- 2 Insert and tighten the screw at the bottom of the FMS




Mounting the BIS Engine to the FMS


- 1 Remove the existing handle for the FMS



- 2 Attach the bracket to the FMS using a M4 x 8mm PHMS






Mounting Bracket



M4 x 8mm PHMS



- 3 Attach the BIS Engine onto the mount



Mounting the FMS

A universal clamp for vertical rail or pole solutions ships with each FMS.



Connections

The cable specifications and part numbers for through wall solutions of the M8048A and M3001A are described in the *Site Preparation* section of this manual.

MSL Cable Termination

The following installation procedure describes how to install the wall installation cable kit when the patient monitor and the measurement server are not located at the same site. The kit consists of two connector boxes and a cable (15m or 25m).



For this procedure you need the insertion tool (M3086-43801) and a small screwdriver.

- 1 Draw the MSL cable through the wall from the site of the monitor to the site of the measurement server.
Each MSL face plate kit contains two connector boxes; one in-going and one out-going. (The US version contains an additional rectangular wall-mounting plate).


NOTE


The installation procedure is the same for both connector boxes. This means you must perform steps 3 to 8 of this procedure twice.

The connectors on each box are different, so you must ensure that the correct box is placed at the correct location. The dots on the plastic angled cover indicates at which site you should install the box:



If there are no dots on the cover, symbols are used:

Symbol:  is connector box (in) and must be placed at the monitor site.



Symbol:  is connector box (out) and must be placed at the measurement server site.

The correct connector cable (M3081-61601, M3081-61602 or M3081-61603) has the opposite symbol:

- 2 Detach the PCB assembly (in/out) from the metallic mounting flange.
- 3 Use the Insertion Tool (M3086-43801) to position each wire on the PCB according the wiring schematic below, where each color corresponds to a number.

NOTE

The Insertion Tool should be set to cutting mode on.

*Wiring Schematic*

- 4 Use a small screwdriver to connect the two drain wires to the PCB, see the wiring schematic in .
- 5 Slide the PCB back on to the metallic mounting flange.
- 6 Use screws to fasten the mounting flange to the wall.

NOTE

US version only: Fasten the rectangular wall-mounting plate to the wall. Attach the mounting flange to the wall-mounting plate.

- 7 Mount the plastic cover. The plastic cover consists of two pieces:
 - Frame
 - Angled cover
 Put the frame over the mounting insert and the PCB. Place the angled cover on top and fasten with two screws.
- 8 Connect the monitor and the measurement server to the wall installation.
- 9 Perform the following tests as described in the Test and Maintenance section of this manual:
 - Power-on test blocks
 - Safety test blocks
 - ECG Sync Performance Test

Philips Clinical Network (Wired)

Installation of the Philips Clinical Network should be performed by Philips service personnel. Use unshielded twisted pair (UTP) cables for installation of the clinical network. Refer to the installation instructions in the M3185A Installation Manual for further details.

WARNING

In order to maintain the galvanic isolation of the IntelliVue monitor, it is essential that UTP (Unshielded Twisted Pair) LAN cables is used to connect the IntelliVue monitor to other devices.

Philips Clinical Network (Wireless)

Refer to the installation instructions in the M3185A Philips Clinical Network Installation Manual for network installation instructions when using the wireless ethernet adapter. For instructions on connecting the wireless ethernet adapter, please refer to the instruction sheet shipped with the mounting device for the adapter.

NOTE

The wireless ethernet adapter is not waterproof and therefore should not be installed anywhere where liquid could spill onto it.

Refer to the IntelliVue 802.11 a/g Infrastructure Installation and Configuration Guide for network installation instructions when using the IntelliVue 802.11 Bedside Adapter. For instructions on connecting the IntelliVue 802.11 Bedside Adapter, please refer to the Hardware Upgrade Guide for your bedside monitor.

IntelliVue Instrument Telemetry (IIT)

Frequency Coordination (USA only):

Frequency coordination is a registration and coordination process for wireless medical telemetry devices used in the U.S.A. which operate in the FCC-allocated Wireless Medical Telemetry Service (WMTS) bands (608-614 MHz, 1395-1400 MHz, 1427-1432 MHz). The 865241/865242 #J45 operates in both of the 1395-1400 and 1427-1432 MHz bands.

Under U.S. Federal Communications Commission (FCC) rules, authorized healthcare providers must register their WMTS devices with an authorized Frequency Coordinator designated by the FCC. The American Society for Healthcare Engineering (ASHE) is the current designated Frequency Coordinator.

Registration/Coordination is a two-step process.

Step 1: Registration: Register the healthcare facility on-line, from the ASHE website). Click on the link for Wireless Medical Telemetry Service and come to the registration page. Fill out the details, and pay the associated fee as per the instructions provided. You will receive confirmation of this registration. Confirmation must be received before proceeding to the next step.

Step 2: Frequency Coordination: Along with confirmation of registration, you will receive access information necessary to perform this second step, frequency coordination. This step involves logging the equipment and frequencies used into the FCC's database, so as to identify any existing potential interference and to help prevent potential future interference. Coordination is accomplished via the ASHE website. Click on the links for Wireless Medical Telemetry Service and then Frequency Coordination. The way the coordination process is executed as of today, it will need to be repeated twice for the M4840A system; once for the 1395-1400 MHz band, and then again for the 1427-1432

MHz band, both of which are used concurrently by the Philips product. There is a separate fee for each coordination request, which varies between \$250 and \$2000, depending upon the number of transmitting devices used and the band/s of operation. Coordination is executed by a company named Comsearch, on behalf of ASHE.

To fill in the frequency coordination forms, you'll need to know the following:

- The county.
- Latitude and longitude that represents the center of the area where the transmitting devices will be deployed. Comsearch can help provide this information; www.comsearch.com provides contact information.
- The name/s of the Clinical Unit/s using the devices (e.g. ICU4, CCU-West, ER1, Step-Down North, etc).
- The radius of deployment, expressed in meters. Imagine drawing a circle around the center of the clinical unit, that encloses/encompasses the unit. What is its radius?
- The number of the highest floor on which a transmitting device will operate.
- How many transmitting devices will be used, i.e. the total number of IntelliVue Instrument Telemetry adapter devices combined.
- The Effective Radiating Power: 6.3 mW.
- The Equipment Manufacturer: Philips Medical Systems.
- The Model numbers: 865241/865242 #J45 IntelliVue Instrument Telemetry adapter used with 865241 (MX700) or 865242 (MX600)
- The Frequency Range to be used: Two separate coordinations are required: For the first one, click on the range of 1395.0 through 1400.0 MHz. For the second one, click on all the frequency ranges listed in the range of 1427.0 through 1432.0 MHz.

When both Registration and Frequency Coordination have been successfully completed, the IntelliVue Instrument Telemetry System can be activated. Note that this process is the responsibility of the customer, as the final “operator” of the transmitting equipment.

Philips IntelliVue Information Center

Please refer to the installation instructions and Instructions for Use of the IntelliVue Information Center Rev. System J or higher.

Short Range Radio

Installation of the Short Range Radio interface should be performed by Philips service personnel.

Before installing an SRR infrastructure it might be necessary to perform a site survey to determine available channels. This should be performed by Philips telemetry installation experts.

Configuring SRR Channels

Hardware Setting: **Main Setup -> Hardware -> SRR Channel**

SRR channel settings only apply for monitors that have a short range radio interface installed. They must be set to match the hospital's wireless infrastructure. SRR channel settings are hardware settings and will typically be set by service personnel at installation.

Refer to your configuration guide for details.

SRR Channel Settings Configuration Implications

Channel Use this setting to configure the SRR channel the monitor should use. SRR provides a total of 16 channels in the ISM (2.4 GHz) band. The channels are labeled 11 to 26. Up to two SRR connections can be established per channel. The ISM band is not exclusively reserved for SRR applications. It is also used by, for example, Wireless LAN (WLAN) and the IntelliVue Telemetry network (except for the US). For this reason, depending on the hospital’s existing wireless infrastructure, a number of SRR channels might already be occupied by other wireless applications.

To achieve the best SRR performance possible, follow these recommendations:

- Usage of WLAN together with SRR may cause interferences. Each WLAN network uses at least four of the 16 SRR channels. If the use of WLAN cannot be avoided, limit the number of channels used for the WLAN infrastructure to a minimum.
- Usage of Bluetooth devices together with SRR may cause interferences. Bluetooth devices automatically change channels regardless of whether a channel is already used by another component of the wireless infrastructure and therefore interfere with SRR connections.
- Usage of cordless phones using the ISM band in the vicinity of SRR devices may cause interferences.
- Usage of wireless PC keyboards or mice using the ISM band in the vicinity of SRR devices may cause interferences.


To assign SRR channels to all monitors in a unit that should be used with SRR connections,

- 1 Identify unused SRR channels. This can be done by using commercially available tools, such as AirMagnet.
- 2 Obtain a floor plan of the unit and identify where the monitors with SRR interface are located.
- 3 Determine SRR groups. An SRR group may contain a maximum of two monitors which share the same SRR channel. Monitors belonging to an SRR group should be located close to each other.
- 4 For each SRR group, assign the same SRR channel to all monitors belonging to a group.

SRR Channel Restrictions with WLAN, IIT, and DECT Devices

The following table and graphic show the restrictions of WLAN, IIT, or DECT Device usage together with SRR.

| | | | | | | | | | | | | | | | | | |
|--|-----------------------------------|-----------------------------------|-----------------------------------|----|----|---------|----|----|----|---------|----|----|----|----|----|-----------------------------------|--|
| US WLAN (802.11) | 1 | | | | 6 | | | | 11 | | | | | | | | |
| Europe WLAN (802.11) | 1 | | | | 7 | | | | 13 | | | | | | | | |
| IIT 2.4 GHz (Smart Hopping Channels) | 0 - 13 | | | | 14 | 15 - 30 | | | | 29 - 47 | | | | | | | |
| SRR (802.15.4) | 11 <small>2405 MHz</small> | 12 <small>2410 MHz</small> | 13 <small>2415 MHz</small> | 14 | 15 | 16 | 17 | 18 | 19 | 20 | 21 | 22 | 23 | 24 | 25 | 26 <small>2480 MHz</small> | |



For a successful SRR deployment, the SRR channels must be located in RF spectra where they are least likely to be interfered with. Choosing appropriate channels after reviewing the Spectrum Analyzer data is critical. In hospitals, 802.11 systems are most the likely source of interference with SRR channels. The figures above show the relationship between 802.11, IIT, and DECT Devices. For example, if the site uses European 802.11 channel 1 for WLAN and has no IIT or DECT devices in the SRR channels 15 or 16, these channels can be used for SRR. Philips telemetry experts will identify available SRR channels by performing a site survey.

When using the Philips IntelliVue 802.11 Bedside Adapter we recommend that you use the 5 GHz band to free the 2.4 GHz band for SRR usage.

NOTE

Short range radio signals are low power signals and therefore have a relatively short range. You can use this fact if the number of unused channels is low, and you run out of channels. Provided the distance between two SRR groups is large enough, i.e. none of the short range radio signals transmitted by the one group can interfere with signals of the other group, you may attempt to assign the same SRR channel to both groups. Take into consideration that portable components (such as Telemetry transceiver, MP5/MP5T or an X2) belonging to one group may be temporarily used within the range of another group.

The range of SRR signals cannot be clearly defined as it depends on external factors such as the components and structure of walls, ceilings, etc.

Electromagnetic Interference (SRR)

Commercially available Short Range Radio 802.15.4 transceivers operate at very low RF power levels to transmit data and need to have high sensitivity receivers to achieve a good link budget. Due to technological limitations the selectivity of the receiver is limited. Consequently, the SRR link is susceptible to other strong RF transmitters not only in the operating frequency band and 5% around it, but also to non-transient RF disturbances stronger than 1V/m at frequencies close to the operating frequency band (2.0 to 2.3 GHz)

ECG Out Functionality

Connections




ECG Out Connector

The cables 8120-1022 and M1181-61625 have both ends terminated. The photograph above shows the monitor side connection.

If using a non-terminated cable:

- 1 Strip 5 mm (3/16") insulation from leads and twist conductor strands tightly.
- 2 Solder leads to the connector as shown in the following diagram.



WARNING

According to AAMI specifications the peak of the synchronized defibrillator discharge should be delivered within 60 ms of the peak of the R wave. The signal at the ECG output on the IntelliVue patient monitors is delayed by a maximum of 30 ms. Your biomedical engineer should verify that your ECG/Defibrillator combination does not exceed the recommended maximum delay of 60 ms.

Computer Client or Hospital Network connected to MX600/MX700 Internal LAN or iPC LAN Interface

You must use Unshielded Twisted Pair (UTP) LAN cables to connect the IntelliVue monitor to other devices. The Computer Client and network infrastructure devices typically are not classified as medical devices and must be located outside the patient vicinity. The patient vicinity is defined as an area within 6ft (1.85m) of the perimeter of the patient's bed or within 7.5ft (2.3m) of the floor.

If the Computer Client is installed in the patient vicinity and connected to the monitoring device, it must be correctly isolated from the mains power supply by an isolation transformer.

If the Computer Client is installed in the patient vicinity and a network switch or hub is used to connect it to a monitoring device, it must be correctly isolated from the mains power supply by an isolation transformer.

WARNING

All external devices in the patient vicinity must comply with IEC 60601-1:1988/A1:1991A2:1995 or EN 60601-1:1990/A1:1993/A2:1995. This applies also to all signal connections, entering the patient vicinity. Additional safety equipment, e.g. isolation transformers might be used.

The installation procedures e.g. for electrical connections as documented in the Instructions for Use must be strictly followed.

Configuration Tasks

You must configure these settings during installation in configuration mode.

- Line Frequency
- Printer
- Altitude
- Equipment Label (for wireless networked monitors, or when the Information center is in flexible monitoring mode).
- ECG cable colors
- Height and Weight units
- Setup Network
- Setup WLAN
- Setup IIT

Checking Country-Specific Default Settings

Some settings are made in the factory to match the typical requirements in a specific country. Line frequency, units for weight and height, and ECG cable colors (AAMI or IEC) have been set to appropriate values. If you suspect that these settings may not match your institution's requirements, check the settings and change them if necessary as described in the *Configuration Guide*.

WARNING

Before starting monitoring, check that the configuration meets your requirements, especially patient category, alarm limits and paced setting.

If you need to enter configuration mode:

- 1 In the **Main Setup** menu, select **Operating Modes**.
- 2 Select **Config** and enter the passcode.
The passcode for configuration mode is given in the monitor's Service Guide.

The monitor displays **Config** at the right hand side of the status line and in the center of the Screen while you are in configuration mode.

Before you leave configuration mode, always be sure to store any changes you made. You must store changes made to each Settings Block and to each Profile, individually. As it may be difficult to remember whether the settings you changed belong to a Monitor Settings block or a Measurement Settings block, we recommend that you store each block before you leave configuration mode.

To leave configuration mode:

- ◆ In the **Main Setup** menu, select **Operating Modes** and then select **Monitoring**.

Setting Altitude, Line Frequency, ECG Cable Colors and Height & Weight Units

You require a local barometric pressure rating from a reliable source (such as airport, regional weather station, or hospital weather station) that is located at the same altitude as the institution.

- 1 From the **Main Setup** menu, select **Global Setting**. Select **Altitude** and enter the altitude.
- 2 From the **Main Setup** menu, select **Global Setting**. Select **Line Frequency** and choose the Line Frequency.
- 3 From the **Main Setup** menu, select **Global Setting**. Select **ECG Cable Color** and choose the Cable Color.
- 4 From the **Main Setup** menu, select **Global Setting**. Select **Height Unit** and choose the Height unit.
- 5 From the **Main Setup** menu, select **Global Setting**. Select **Weight Unit** and choose the Weight unit.

Setting Altitude and Line Frequency

You require a local barometric pressure rating from a reliable source (such as airport, regional weather station, or hospital weather station) that is located at the same altitude as the institution.

- 1 From the **Main Setup** menu, select **Global Setting**. Select **Altitude** and enter the altitude.
- 2 From the **Main Setup** menu, select **Global Setting**. Select **Line Frequency** and choose the Line Frequency.

Configuring the Equipment Label

If the Information Center is in fixed monitoring mode, it controls the equipment label. You do not need to follow this procedure.

However, if you are on a wireless network, or your Information Center is configured for flexible monitoring mode, you must set the equipment label. This associates the monitor with a central monitoring sector. An identical monitor label must also be configured in the Information Center.

- 1 Select **Main Setup** -> **Bed Information** to call up the **Bed Information** menu.
- 2 Select **Equipment Label** to call up the onscreen keyboard.
- 3 Enter the system identifier. This needs to be set up in either the monitor or the Information Center. If the Information Center is in flexible monitoring mode, the monitor must be setup to match the Information Center's monitor label.

Configuring the Printer

- 1 From the **Main Setup** menu select **Reports**.
- 2 Select **Printer Settings** and configure **Local** to enabled if the printer is connected directly to the monitor. See configuration guide for further details.

Setup Network

Main Setup => Network => Setup Network

| Network Setup | |
|-----------------|-----------|
| IP Config | BOOTP |
| Dynamic IP | Mandatory |
| IP Address | 0.0.0.0 |
| Subnet Mask | 0.0.0.0 |
| Default Gateway | 0.0.0.0 |
| IGMP | Off |
| CI Mode | Broadcast |
| CI Address | 0.0.0.0 |
| CI TTL | 1 |
| QoS State | Off |
| QoS Level | 0 |
| Directory Mode | Broadcast |
| Directory Addr | 0.0.0.0 |

The Setup Network Menu allows you to configure the following items:

- IP Config** Allows you to change the mode to acquire the IP address. Choices are **BOOTP**, **DHCP**, **DHCP restricted** and **Manual**. In **DHCP restricted** mode any true BOOTP responses are ignored.
- Dynamic IP** Can be either **Mandatory** (default) or **Optional**. If **Mandatory** is set the bedside will announce an INOP if an IPv4 address could not be acquired.
- IP Address** If IP Config is set to **Manual**, the IP Address should be a valid IP address, e.g. not 0.0.0.0. The configured value is ignored when IP Config is not set to **Manual**.
- Subnet Mask** The Subnet Mask must be provided for manual IP addresses. The Subnet Mask must consist of a single consecutive series of “1” bits; e.g. 255.255.248.0. The configured value is ignored when IP Config is not set to **Manual**.
- Default Gateway** The IP Address of the Default Gateway can be optionally configured. The configured value must be within the range of the Subnet Mask. The configured value is ignored when IP Config is not set to **Manual**.
- IGMP** Enables or disables IGMP support (On/ Off).
- CI Mode** The mode in which CI messages (Connect Indication messages) are sent (**Broadcast**, **Multicast**, **Manual**).
- CI Address** Multicast Address for Device Connect Indication (only used if CI Mode is set to **Manual**). If **CI Mode** is set to **Broadcast**, the CI Address is implicitly the subnet broadcast address. If **CI Mode** is set to **Multicast** the CI Address is implicitly 224.0.23.63
- CI TTL** Sets the TTL (Time To Live) of the CI message. Defaults to **1**.
- QoS State** QoS (Quality of Service) can be switched **On** or **Off**. Defaults to **Off**.
- QoS Level** The QoS Level (DSCP CoS value) can be entered between 0 and 7. Defaults to 0.

7 Installation Instructions

Directory Mode The mode in which Directory Service requests are received (Broadcast, Multicast)

Directory Addr Mode Multicast address for Directory Service requests (if Directory Mode set to Multicast).

Store The entered values are verified when the “Store” softkey is pressed. A valid configuration is immediately activated and the window closed. If any field is invalid (e.g. invalid subnet mask or gateway not in subnet), an error popup window is shown and the configuration is not applied.

Network Status Window

The network status window provides network status information which can be used for troubleshooting purposes.

Main Setup => Network => Network Status

| Network Status | |
|-----------------|---------------------------|
| Network Type | LAN |
| MAC Address | 0030 D301 376B |
| MAC IIT | 0097 3498 87EC |
| RF Access Code | 2 |
| IP Config | BOOTP |
| Dynamic IP | Mandatory |
| IP Address | 172.31.10.03 |
| Subnet Mask | 255.255.0.0 |
| Default Gateway | 172.31.20.10 |
| BOOTP Server | 172.31.10.165 |
| Server IP | 172.31.10.200 |
| IPv6 Address | FE80::0209:FBFF:FE6E:D163 |
| IGMP | Off |
| CI Mode | Broadcast |
| CI Address | 172.31.255.255 |
| CI TTL | 1 |
| QoS State | Off |
| QoS Level | 0 |
| Directory Mode | Broadcast |
| Directory Addr. | 172.31.255.255 |

Network Type Displays the current active network type (**LAN, IIT, WLAN**)

IP Config Shows the IP Configuration mode. Choices are **BOOTP, DHCP, DHCP restricted** and **Manual**. In **DHCP restricted** mode any true BOOTP responses are ignored.

Dynamic IP Can be either **Mandatory** (default) or **Optional**. If **Mandatory** the bedside monitor will issue an INOP if an IPv4 address could not be acquired.

The **MAC Address**, current IP address configuration (**IP Address, Subnet Mask, Default Gateway, Server IP address and CI mode/ CI Address/ CI TTL**) and **BOOTP/ DHCP Server** is displayed. The IP address configuration can be changed in Service Mode only. Manually entered IP addresses are marked by the suffix “Manual”, invalid fields are marked by the suffix “Invalid”.

If the Instrument Telemetry interface is available, its **MAC Address** and **RF Access Code** are displayed. In addition, the bedside monitor's IPv6 link local address is shown.

Setup WLAN

Main Setup => Network => Setup WLAN

| Setup Wireless LAN | |
|----------------------|--------------|
| Country | 1000 |
| IP Address | 172.31.10.04 |
| WMM Mode | disabled |
| Mode | 802.11bg |
| SSID | Philips |
| Security Mode | WPA2Enterpr. |
| WEP Key Size | 40 bit |
| WEP Key | * |
| WEP Key Index | 2 |
| WPA Password | * |
| Authentication | PEAP |
| Inner Authentication | MSCHAPv2 |
| PEAP Version | Default |
| PEAP Label | Default |
| User Name | username |
| Password | secret |
| Anonymous Ident. | anonymous |
| CertificateCheck | Enabled |
| CA Certificate | Installed |

The Setup Wireless LAN Menu allows you to configure the following items:

Country Country setting for the WLAN card to adapt to the local frequency and transmit power regulations. The default value 1000 should be left to support adapting to country provided by the AP.

7 Installation Instructions

IP Address IP Address of the WLAN card. Typically the automatic configuration via the BOOTP/ DHCP Server of the central station is used. In this case the field is set to 0.0.0.0. For special requirements, it is possible to switch to a manual/fix IP address configuration.

Manually entered IP addresses are marked by the suffix “Manual”

NOTE

- The subnet mask and the gateway address of the WLAN card will be taken from the corresponding fields in the **Bed Information** or **Setup Network** menu.
- Only limited checks of the manual values are possible. Therefore it is mandatory that a manual configuration is only performed by experienced service personnel to prevent problems such as duplicate IP addresses, non matching subnet mask, etc.

WMM Mode Changes to WMM (Wireless Multimedia Mode) mode for the WLAN adapter supporting this feature.

Mode Defines the WLAN operating mode/ IEEE 802.11 Standard. Valid values are either 802.11ah, 802.11bg, 802.11g, Auto or None. None indicates that no wireless connection should be established although a wireless adapter is available and operational.

SSID Set Identifier: Logical WLAN Network Name.

Security Mode **WEP, WPA(PSK) or WPA2(PSK), WPA-Enterprise and WPA2-Enterprise** with either Protected EAP (PEAP) or Tunneled TLS (TTLS) as authentication methods.

WEP Key Size 40 bit or 104 bit.

WEP Key The number of hex chars for the WEP key depends on the WEP key size chosen. For a 40 bit WEP key size the WEP key must be 10 hex chars long, for a 104 bit key the WEP key must be 26 hex chars long. Will be shown as “*” after the user entered the key.

WEP Key Index Defines the transmit WEP Key Index. This entry must match the WEP Key Index configured at the infrastructure device, i.e. on a WLAN Access Point, and ranges from 1 to 4.

WPA Password In WPA(PSK) or WPA2(PSK) mode this entry defines the Pre-Shared-Secret or Password with 8 to 63 alpha-numeric characters. Will be shown as “*” after the password is entered.

In WPA-Enterprise or WPA2-Enterprise mode the following read only fields are available. The configuration can only be performed via the IntelliVue Support Tool:

Authentication Authentication method can be either Protected EAP (PEAP) or Tunneled TLS (TTLS).

Inner Authentication PEAP and TTLS can be used with several different Inner Authentication methods. PEAP with MSCHAPv2 and TTLS with PAP, CHAP, MSCHAP or MSCHAPv2

PEAP Version This setting describes the PEAP protocol version to be used while authenticating against the authentication server. Valid values are either Default, Version 0 and Version 1. If set to Default the decision is up to the wireless adapter. Version 0 or 1 forces the wireless adapter to use the protocol version required for a certain authentication server. This setting is intended for experts only.

PEAP Label The PEAP label setting defines the string to be use to signal EAP-PEAP encryption to the authentication server. Valid values are Default, EAP or PEAP. Default lets the decision up to wireless

adapter. Both EAP and PEAP forces the wireless adapter to use this setting. This setting is intended for experts only.

- Username** The username used in the encrypted tunnel with 1-63 alpha-numeric characters. It is also used as outer identity as long as the Anonymous Identity is not set.
- Password** The password used in the encrypted tunnel with 8-63 alpha-numeric characters. Will be shown as four stars “****” after the user entered the password.
- Anonymous Identity** The identity used for the outer PEAP or TTLS authentication, which may be “unprotected”. Thus, the identity should be different to the Username for enhanced security. The Anonymous Identity contains 1-63 characters. It can be set to NotConfigured by clearing it.
- CA Certificate** Indicates that the certificate of the Certification Authority (CA) has been installed on the device. The CA certificate is the root certificate of the server delivered certificate chain. The certificate chain ends with the server certificate, which contains the key material used to build up the encrypted tunnel for PEAP or TTLS. Values are either Installed or NotInstalled. This field is not user configurable. The CA certificate can only be installed with the Intellivue Support Tool.
- Certificate Check** As long the Certificate Check is set to Enabled, the CA Certificate is used to verify the authenticity of the certificate chain delivered by the authentication server. The verification involves also the system time to check the validity period of every certificate in the chain. This item can only be set to Enabled, if an CA Certificate has been installed. Valid values are Disabled or Enabled.

WLAN Diagnostic Window

The WLAN Diagnostic window provides WLAN status information which can be used for troubleshooting purposes.

Main Setup => Network => WLAN Diagnostic

| WLAN Diagnostic | |
|------------------|----------------|
| Wireless LAN | On |
| Check WLAN | Security Mode |
| MAC WLAN | 0097 3498 87EC |
| IP Address WLAN | 172.31.11.15 |
| IP Address | 172.31.10.03 |
| Server IP | 172.31.10.165 |
| Subnet Mask | 255.255.0.0 |
| Country | 1000 / 1276 |
| WMM Mode | Disabled |
| Mode | 802.11bg |
| Security Mode | WPA2Enterpr. |
| Authentication | PEAP |
| Inner Authentic. | MSCHAPv2 |
| PEAP Version | Default |
| PEAP Label | Default |
| User Name | username |
| Anonymous Ident. | anonymous |
| CertificateCheck | Enabled |
| CA Certificate | Installed |
| SSID | Philips |
| Rate | 54,0 Mb/s |
| Active Channel | 9 |
| RSSI | 65 - |
| Conn. Status | Connected |
| MAC AP WLAN | 00A0 F8CE D231 |

In the WLAN Diagnostic window, the wireless LAN address information is shown (**MAC address, IP addresses** and **Subnet Mask**) as well as the active security settings, WLAN **Mode**, **SSID** as well as **Country** setting, which are all not editable in this window. In addition, the following status information is shown:

Wireless LAN State of the wireless adapter. Off indicates inactive or just starting, On indicates operational. Does not take the wireless state into account.

| | |
|-----------------------|---|
| Check WLAN | Indicates that a configuration issue has been detected. The message text depends on the operating mode. |
| Country | The configured country code and the dynamically chosen country code based on access point setting. |
| WMM Mode | The WMM mode being used. |
| RSSI | Received Signal Strength Indicator. |
| Conn. Status | Current wireless LAN connection status (None, Scanning, Authenticating, Authenticated Shared, Associating, Connected, Link Problem). |
| Active Channel | Current radio channel. |
| Rate | Currently selected transmission rate (adapts dynamically based on wireless signal propagation behavior). |
| MAC AP WLAN | The MAC address of the access point to which a connection has been established. |

WLAN Configuration Check

Before a wireless connection is established, the configuration is checked for basic issues or incompatibilities. If issues are found during this check, a status message will be shown in the status prompt area of the monitor's resting display. Additionally, the Check WLAN line of the WLAN Diagnostic window will display the status message. In Service Mode, a detailed status message is shown, while a simplified message is shown in all other operating modes.

The table below lists the status messages displayed in service mode, depending on the issue found.

| Issue | Status Message |
|----------------------|--------------------------------------|
| Security Mode | Check WLAN Security Mode |
| Country Setting | Check WLAN Country Setting |
| IEEE Wireless Mode | Check WLAN IEEE Mode |
| WMM Mode | Check WLAN WMM Mode |
| WEP Encryption | Check WLAN WEP Settings |
| WPA Encryption | Check WLAN WPA/WPA2 Key |
| Authentication | Check WLAN Authentication Mode |
| Inner Authentication | Check WLAN Inner Authentication Mode |
| PEAP settings | Check WLAN PEAP Properties |
| Certificate Check | Check WLAN Certificate Check Setting |
| User Name,Password | Check WLAN EAP Credentials |

The simplified message in all other operating modes is "Check WLAN Settings"

Setup IIT

Main Setup => Network => Setup IIT

| IIT Setup | |
|----------------|---------|
| IIT | Enabled |
| RF Access Code | 1 |

The Setup IIT Menu allows you to configure the following items:

- Instr. Telemetry** The Instrument Telemetry network interface can be disabled by a Global Setting to suppress network related technical INOPS if the IIT-capable device is operated in a non-IIT environment.
- RF Access Code** Instrument telemetry specific partitioning of the radio domain. Displays “Disabled” instead of RF Access Code if Instrument Telemetry is disabled by Global Setting.

Instrument Telemetry Diagnostic Window

The Instrument Telemetry Diagnostic window provides IIT status information which can be used for troubleshooting purposes.

Main Setup => Network => IIT Diagnostic

| Instrument Telemetry Diagnostic | |
|---------------------------------|----------------|
| MAC IIT | 0097 3498 87EC |
| RF Access Code | 2 |
| IP Address | 172.31.10.03 |
| Server IP | 172.31.10.165 |
| Subnet Mask | 255.255.0.0 |
| RSSI | -50 |
| Conn. Status | Active |

Network addresses To support troubleshooting of the network connection, the following fields are displayed: Instrument Telemetry **MAC Address**, **IP Address**, **BOOTP/ DHCP Server IP address** and **Subnet Mask**.

RF Access Code Instrument telemetry specific partitioning of the radio domain. Displays “Disabled” instead of RF Access Code if Instrument Telemetry is disabled by Global Setting.

RSSI Received Signal Strength Indicator.

Conn. Status Current connection status of the Instrument Telemetry Subsystem.

IntelliBridge EC10

The IntelliBridge EC10 web based service interface allows you to:

- upgrade the EC10 Firmware
- upload or remove device drivers


- generate and upload clone files
- view system information

Accessing the IntelliBridge EC10 Service Interface

- 1 Connect the Service PC to the IntelliBridge EC10 module using a standard LAN cable.
- 2 Make sure that the TCP/IP settings of your Local Area Connection Properties on the Service PC are set to "**obtain IP Address automatically**".
- 3 Open your Internet Browser. Internet Explorer 7 and Firefox 3.0 are supported.
- 4 Enter **192.168.100.2** in the navigation field of the internet browser.
- 5 The EC10 Configuration screen will open up. Enter the following login information:


Username: service

Password: IBEC10



Firmware Upgrade


- 1 In the EC10 Configuration Screen select the tab **Upgrade**.
- 2 The Firmware upgrade window will open up.



- 3 Click on the **Browse** button and select the correct FW file.
- 4 Click **Go**.


Uploading and Removing Device Drivers

- 1 In the EC10 Configuration Screen, select the **Drivers** Tab.
- 2 A list of the available drivers on the EC10 module will appear.
- 3 To upload a new driver, click the **Add** button, then click the **Browse** button and select the driver file. To remove one or more drivers, select the driver(s) from the list of available drivers and click the **Remove** button.




Generating and Uploading Clone Files

- 1 In the EC10 Configuration Screen, select the **Cloning** tab.
- 2 To generate a clone file, click the **Go** button in the *Save clone file from IntelliBridge EC10* window. A hyperlink will appear. Right click this link and select "**Save Target as...**" to save it to the desired location.
- 3 To upload a clone file click the **Browse** button in the *Load Clone File to IntelliBridge EC10* window and select the clone file. Then click **Go**.



Viewing System Information

- 1 In the EC10 Configuration Screen, select the **System** Tab. This screen allows you to view Log files and other system information.



Handing Over the Monitor

If you are handing over the monitor to the end-users directly after configuration, make sure that it is in Monitoring mode.

Ensure that the users have access to the following documentation delivered with the monitor:

- Training Program - M8000-9461x - for self-training on the monitor before use
- Instructions for Use - M8000-9001x - for more detailed questions during use

WARNING

All users must complete the training program and read the Instructions for Use before working with the monitor.

These training materials (in combination with this service guide) can also be used to train service personnel on how to use and service monitor.

Please refer to the Testing and Maintenance section of this service guide to determine the respective safety and performance tests to be performed before handing over the monitor to the customer.

Site Preparation

Introduction

This section describes the procedures you should follow to plan and prepare a site for an IntelliVue monitor installation. It describes:

- Site planning.
- Roles and responsibilities for local and Philips personnel.
- Remote installation planning.

Site Planning

The careful planning of the site for the IntelliVue monitor is essential for its safe and efficient operation. *A consulting schedule should be established between the Customer and Philips Sales and Support Representatives, to ensure that all preparations are completed when the system is delivered.*

The site planning phases prior to equipment installation are:

Location: Planning the location of the various system components.

Environment: Confirming and correcting, as necessary, the environment of the proposed installation site(s).

System Capabilities: Explaining the possibilities for system expansion.

Mounting: Referencing the mounting hardware information website for the listing of suitable mounting hardware recommended for use with the various system components, and all details on the available mounts and accessories.

Cabling: Identifying the requirements for the cabling, conduiting and faceplates for connecting the various system components.

Roles & Responsibilities

This section describes the procedures necessary to prepare a site for a system installation. The procedures are grouped into two parts: procedures that local staff or contractors are responsible for, and procedures that Philips personnel are responsible for.

Site Preparation Responsibilities

Local Staff

- Ensure that all safety, environmental and power requirements are met.
- Provide power outlets.
- Prepare mounts.
- Pull cables, install conduit, install wallboxes.
- Terminate network cables if a Philips Clinical Network is in use.
- It may be necessary to certify the network cable plant, see Philips Clinical Network Installation Manual for details.

Philips Personnel

- Provide the customer with the safety, environmental and power requirements.
- Assemble mounts.
- Prepare monitor remote cabling.

Procedures for Local Staff

The following tasks must be completed **before** the procedures for Philips personnel may be started.

- **Providing Power Outlets**
One power outlet for each display and for any peripheral device (for example, a printer or slave display) is required by the system. Provide a power outlet in the vicinity (1 m or 3 ft) of each component that requires power.

WARNING

Only the power cables provided with the system may be used. For reasons of safety, power (mains) extension cables or adapters shall not be used.

- **Preparing Mounts**
Where ceiling, wall, or shelf mounts are required for mounting the equipment, the customer is responsible for the following:
 - Providing and installing all hardware which is required to install the mounting hardware supplied by Philips as detailed in the installation notes.
 - Making sure that all ceilings, walls, and mounting rails that supports mounting hardware are suitable for their proposed load.


WARNING

It is the customer's responsibility to have the attachment of the mounting hardware to the ceiling, wall, or mounting rail and the construction of the ceiling, wall, or mounting rail evaluated for structural integrity and compliance with all local, state and any other required codes by a registered, professional, structural and/or mechanical engineer.

Although considerable effort has been made to ensure the safety of the ceiling mount installation and or mounting guidelines, it is to be understood that the installation itself is beyond the control of Philips

Medical Systems. Accordingly, Philips Medical Systems will not be responsible for the failure of any such installation.

The following figures show the dimensions required for the M1180A #C53 wall and the table mounting bracket which ships with the monitor.



Wall Mounting Bracket Dimensions




Table Mounting Bracket Dimensions

- Providing Conduit

Where a remote installation is required, for example the installation of a remote display, the customer is responsible for the following hardware installations:

- Providing conduit and/or trunking of a sufficient cross-sectional area for the planned cables and possible future expansion (for additional components or systems). See “Cabling Options and Conduit Size Requirements” on page 323 for cable specifications for remote installations.
- Providing and/or installing suitable wall boxes to accommodate the faceplates.

- Pulling Cables

WARNING

NEVER run power cables through the same conduit or trunking used for system cables.

- Installing Wall Boxes

It is the customer's responsibility to provide and install wallboxes to house faceplates. The customer must notify the Philips installation coordinator of which size is to be used.

Procedures for Philips Personnel

Before you begin the procedures in the installation sections, ensure that the customer has completed all necessary preparations outlined in the previous section, "Procedures for Local Staff."

Monitor MX600/700 Site Requirements

Space Requirements

The siting of the monitor should be planned such that the nursing staff are able to monitor the patient with relative ease, with all patient connectors and controls readily available and the displays clearly visible. The location should also allow access to service personnel without excessive disruption and should have sufficient clearance all round to allow air circulation.

Maximum dimensions and weight:

Size (W x H x D)

392mm x 321mm x 163mm (15.43" x 12.64" x 6.42") (without cable management)

392mm x 321mm x 235mm (15.43" x 12.64" x 9.25") (with cable management)

Weight

9.5kg (20.94lb) without options

Environmental Requirements

The environment where the MX600/700 monitor will be used should be reasonably free from vibration, dust and corrosive or explosive gases. The ambient operating and storage conditions for the MX600/700 monitor must be observed. If these conditions are not met, the accuracy of the system will be affected and damage can occur.

Temperature

Operating: 0 to 40°C (32 to 104°F) without iPC, 0 to 35°C (32 to 95°F) with iPC

Storage: -20 to 60°C (-4 to 140°F)

Humidity

Operating: 15% to 95% Relative Humidity (RH) (non-condensing)

Storage: 5% to 95% Relative Humidity (RH)

Altitude

Operating: -500m to 3000m (10000 ft.)

Storage: -500m to 4600m (15000 ft.)

Electrical and Safety Requirements (Customer or Philips)

Safety Requirements

If the monitor is to be used in internal examinations on the heart or brain ensure that the monitor is connected to an equipotential grounding system.

Grounding

The monitor **MUST** be grounded during operation (Class I equipment according to IEC 60601-1). If a three-wire receptacle is not available then the hospital electrician must be consulted to ensure that proper grounding is available on installation. **NEVER** attempt to use a three-wire to two-wire adapter with the monitor.

WARNING

Each component must be individually grounded for safety and interference suppression purposes.

Electrical Requirements

Line Voltage Connection

The MX600/700 monitor uses < 200 W.

Line Voltage

The MX600/700 monitor may be operated on ac line voltage ranges of 100 to 240V (50/60 Hz).

Remote Device Site Requirements

The system can be installed with one or more combinations of the following remote devices.

Flexible Module Rack or Multi-Measurement Module

Remote Display

Remote Control

Where more than one site is used for locating equipment (a remote installation), the following sections should be considered for EACH device:


- Space Requirements
- Environmental Requirements
- Mounting
- Electrical and Safety Requirements
- Cabling Options and Conduit Size Requirements

Connecting Non-Medical Devices

The standard IEC-60601-1-1 applies to any combination of medical and non-medical electrical devices, where at least one is a medical electrical device. Therefore IEC-60601-1-1 must still be met after all devices are connected.

WARNING

Do not use a device in the patient vicinity if it does not comply with IEC-60601-1. The whole installation, including devices outside of the patient vicinity, must comply with IEC-60601-1-1; one reasonable solution may be the use of an isolation transformer.



Equipment Location in the Patient Vicinity

NOTE

The site planning requirements, with the exception of the cabling, must be provided by the device manufacturer, if the remote device is not purchased from Philips.

Multi-Measurement Module (MMS) M3001A, IntelliVue X2 M3002A, 8-Slot Flexible Module Rack (FMS-8) M8048A, 4-Slot Flexible Module Rack (FMS-4) 865243 or Remote Control 865244

Space Requirements Multi-Measurement Module M3001A

Size (W x D x H)

188.0mm x 96.5 mm x 51.5 mm

(7.40" x 3.80" x 2.03")

Weight

650g (1.4 lb)

Space Requirements IntelliVue X2 M3002A

Size (W x D x H)

188.0mm x 99 mm x 86 mm

(7.4" x 3.9" x 3.4")

Weight

1.5 kg (3.3 lb)

Space Requirements 8-Slot Flexible Module Rack (FMS-8) M8048A

Size (W x D x H)

320 mm x 120mm x 135mm (12.6" x 4.72" x 5.3")

Weight

< 3500g (7.7lb)

Space Requirements 4-Slot Flexible Module Rack (FMS-4) 865243

Size with MMS Mount (W x H x D)

232 mm x 139mm x 188mm (9.13" x 5.47" x 7.4")

Size without MMS Mount (W x H x D)

194 mm x 139 mm x 110 mm (7.64" x 5.47" x 4.33")

Weight

< 1.1kg (2.43lb)

Space Requirements Remote Control 865244

Size (W x H x D)

53 mm x 165 mm x 23 mm

(2.1" x 6.5" x 0.9")

Weight

<250 g (0.55 lb)

Environmental Requirements Multi-Measurement Module M3001A

Temperature

Operating: 0 to 45°C (32 to 113°F)

Storage: -40 to 70°C (-40 to 158°F)

Humidity

Operating: 95% relative humidity (RH) max. @ 40°C (104°F)

Storage: 90% relative humidity (RH) max. @ 65°C (150°F)

Altitude

Operating: -500m to 4600m (-1600 to 15000 ft.)

Storage: -500m to 15300m (-1600 to 50000 ft.)

Environmental Requirements IntelliVue X2 M3002A

Temperature

Operating: 0 to 40°C (32 to 100°F)

0 to 35°C (32 to 95°F) (while charging batteries or when equipped with IIT or IntelliVue 802.11 Bedside Adapter)

Storage: -20 to 60°C (-4 to 140°F)

Humidity

Operating: 15 to 95% relative humidity (RH) max. @ 40°C (104°F)

Storage: 5 to 95% relative humidity (RH) max. @ 65°C (150°F)

Altitude

Operating: -500m to 3000m (-1600 to 10000 ft.)

Storage: -500m to 4600m (-1600 to 15000 ft.)

Environmental Requirements 8-Slot Flexible Module Rack M8048A

Temperature

Operating: 0 to 45°C (32 to 113°F)

Storage: -40 to 70°C (-40 to 158°F)

Humidity

Operating: 95% relative humidity (RH) max. @ 40°C (104°F)

Storage: 90% relative humidity (RH) max. @ 65°C (150°F)

Altitude

Operating: -500m to 4600m (-1600 to 15000 ft.)

Storage: -500m to 15300m (-1600 to 50000 ft.)

Environmental Requirements 4-Slot Flexible Module Rack 865243

Temperature

Operating: 0 to 40°C (32 to 100°F)

Storage: -20 to 60°C (-4 to 140°F)

Humidity

Operating: 15 to 95% relative humidity (RH)

Storage: 5 to 90% relative humidity (RH)

Altitude

Operating: -500m to 3000m (-1600 to 10000 ft.)

Storage: -500m to 4600m (-1600 to 15000 ft.)

Environmental Requirements Remote Control 865244

Temperature

Operating: 0 to 40°C (32 to 100°F)

Storage: -20 to 60°C (-4 to 140°F)

Humidity

Operating: 15% to 95% relative humidity (RH) (non-condensing)

Storage: 5% to 90% relative humidity (RH)

Altitude

Operating: -500m to 3000m (-1600 to 10000 ft.)

Storage: -500m to 4600m (-1600 to 15000 ft.)

Cabling Options and Conduit Size Requirements

The following table describes the cabling options for the FMS and the MMS.

Table 10 M8048A, M3001A and M3002A Cables

| Part number | 12NC | Description | Conduit sizes | Bend Radius | Connector Size (L x W) |
|---|--------------|--------------------------------|-------------------|-------------|------------------------|
| Both ends are terminated with MSL connectors | | | | | |
| M3081-61626 | 453563474781 | CBL MSL 0.75m | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61601 | 453563402721 | CBL MSL 1m | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61602 | 453563377851 | CBL MSL 2m | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61603 | 453563402731 | CBL MSL 4m | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61627 | 453563484501 | CBL MSL 10m terminated | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61628 | 453563484511 | CBL MSL 15m terminated | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61629 | 453563484521 | CBL 25m MSL Install terminated | 72mm ² | 40 mm | 40 mm x 17 mm |

8 Site Preparation

| Part number | 12NC | Description | Conduit sizes | Bend Radius | Connector Size (L x W) |
|---|--------------|---|-------------------|-------------|------------------------|
| Both ends are unterminated, use faceplates | | | | | |
| M3081-61615 | 453563484481 | CBL MSL Instal.15m unterminated | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-61625 | 453563484491 | CBL MSL Instal.25m unterminated | 72mm ² | 40 mm | 40 mm x 17 mm |
| Faceplates | | | | | |
| M3081-68707 | 453563484531 | CBL MSLFace Plate pair of boxes NON US | 72mm ² | 40 mm | 40 mm x 17 mm |
| M3081-68708 | 453563484541 | CBL MSL Face Plate pair of boxes US | 72mm ² | 40 mm | 40 mm x 17 mm |
| Insertion Tool | | | | | |
| M3086-43801 | 453563484771 | MSL Insertion Tool for unterminated CBL | 72mm ² | 40 mm | 40 mm x 17 mm |

Mounting

Table 11 M3001A Multi-Measurement Module (MMS) and M3002A IntelliVue X2 Mounts

| Product Option Number | Part Number 12NC Part No. | Description |
|-----------------------|------------------------------|---------------------|
| M3080A #A01 | n/a | Wall Mounting Plate |

Table 12 M8048A Flexible Module Rack Mounts

| Product Option Number | Part Number | 12NC Part No. | Description |
|-----------------------|-------------|---------------|--|
| M8048A #E15 | M4041-42303 | 453563494101 | Cable Management |
| M8048A #E20 | M4041-60005 | 453563477961 | MMS Mount |
| | M8048-64002 | 453563456901 | Handle (part of small parts kit) |
| M1180A #C33 | M1034-60104 | 453563490591 | BIS Mount |
| n/a | M4041-22302 | 451261011861 | Mounting Plate Adapter clamp (for wall mounting) |

Remote Displays (M8031B)

Space Requirements

Size (W x D x H)

372mm x 308mm x 74.1mm (14.65" x 12.13" x 2.92")

Weight

Without deskstand: 5200g (11.5lb)

With desk stand: 9000g (19.8lb)

Environmental Requirements

Temperature

Operating: 0 to 40°C (32 to 104°F)

Storage: -20 to 60°C (-4 to 140°F)

Humidity

Operating: 20 to 85% RH (Non-condensing)

Storage: 5 to 85% RH (Non-condensing)

Altitude

Operating: Up to 4000m (13123.36 ft.)

Storage: Up to 12000m (39370,08 ft.)

Electrical and Safety Requirements

Voltage ranges:

90V to 264V

Voltage selection:

Wide range input, no voltage selection required

Power consumption: ~30W

Remote Displays (865299)

Space Requirements

Size (W x H x D)

429mm x 359mm x 61mm (16.89" x 14.13" x 2.4")

Weight

6.1 kg (13.45 lbs.)

Environmental Requirements

Temperature

Operating: 0 to 40°C (32 to 104°F)

Storage: -20 to 60°C (-4 to 140°F)

8 Site Preparation

Humidity

Operating: 30 to 70% RH (Non-condensing)

Storage: 10 to 90% RH (Non-condensing)

Altitude

Operating: Up to 4000m (13123.36 ft.)

Storage: Up to 12000m (39370,08 ft.)

Electrical and Safety Requirements

Voltage ranges:

100V to 240V

Voltage selection:

Wide range input, no voltage selection required

Power consumption: ~48W

Cabling Options and Conduit Size Requirements

The following table describes the cabling options for the M8031A/B 15", the M8033A/B 17" and the 865299 19" TFT Medical Grade Touch Displays.

Table 13 Digital Video cables

| Part number | 12NC | Description |
|-------------|--------------|-----------------------|
| M8071-61001 | 453563484551 | IV CBL 1.5m Cable Kit |
| M8071-61002 | 453563484561 | IV CBL 3m Cable Kit |
| M8071-61003 | 453563484571 | IV CBL 10m Cable Kit. |

Touch Cable

| Part number | 12NC | Description | Conduit sizes | Bend Radius | Connector Size (L x W) |
|-------------|--------------|--------------------|-------------------|-------------|------------------------|
| M8081-61010 | 451261006551 | IV CBL Touch, 1.5m | 30mm ² | 25 mm | 35 x 16 mm |
| M8081-61011 | 451261006561 | IV CBL Touch, 3m | 30mm ² | 25 mm | 35 x 16 mm |
| M8081-61012 | 451261006571 | IV CBL Touch, 10m | 30mm ² | 25 mm | 35 x 16 mm |
| M8081-61013 | 451261006581 | IV CBL Touch, 15m | 30mm ² | 25 mm | 35 x 16 mm |
| M8081-61014 | 451261006591 | IV CBL Touch, 25m | 30mm ² | 25 mm | 35 x 16 mm |

IntelliBridge

Cabling Options

| Part number | 12NC | Description | Conduit Sizes | Max. Bend Angle | Connector Size (L x W) |
|-------------|--------------|------------------------------------|--------------------|-----------------|------------------------|
| M8081-61001 | 453563484591 | IV CBL 1.5m LAN/ RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |
| M8081-61002 | 453563484601 | IV CBL 3m LAN/RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |
| M8081-61003 | 453563484611 | IV CBL 10m LAN/ RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |

Local Printer

See printer documentation

RS232/MIB/LAN Interface

Table 14 MIB Cable and Serial Cable

| Part number | 12NC | Description | Conduit Sizes | Max. Bend Angle | Connector Size (L x W) |
|-------------|--------------|------------------------------------|--------------------|-----------------|------------------------|
| M8081-61001 | 453563484591 | IV CBL 1.5m LAN/ RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |
| M8081-61002 | 453563484601 | IV CBL 3m LAN/RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |
| M8081-61003 | 453563484611 | IV CBL 10m LAN/ RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |
| M8081-61004 | 453563484621 | IV CBL 15m LAN/ RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |
| M8081-61005 | 453563484631 | IV CBL 25m LAN/ RS232 straight | 30 mm ² | 25 mm | 15 x 15 mm |

Cable and Adapter Set

Nurse Call Paging Cable

| Part number | 12NC | Description | Conduit Sizes | Bend Radius | Connector Size |
|--------------|--------------|---|-----------------------|-------------|---|
| 453564214531 | 453564214531 | CBL Nurse paging cable 3m | 0.125 mm ² | 12 mm | 9 x 12 mm One end terminated with 6P6C connector. one end without connector. |
| 453564220531 | 453564220531 | CBL Nurse paging cable 10m | 0.125 mm ² | 12 mm | 9 x 12 mm One end terminated with 6P6C connector. one end without connector. |
| M8087-61001 | 453563484741 | IV CBL 10m flex nurse paging MDR&open con | 40 mm ² | 30 mm | 35 x 16 mm One end terminated with straight MDR connector, one end without connector |

ECG Out Interface

Table 15 ECG Out Cable

| Part number | 12NC | Description | Conduit Sizes | Bend Radius | Diameter | Comments |
|-------------|--------------|---|--------------------|-------------|----------|--|
| 8120-1022 | 453563198151 | CBL ECG Out 3m Phone plug terminated | 40 mm ² | 30 mm | 13 mm | 3m cable (Both ends are terminated with .25" phone plugs) |
| M1181-61625 | 453563255091 | CMS CBL ECG OUT CABLE KIT 25m | | | | cable kit consisting of: 25 m raw cable, 2 x 1/4" socket, 1 x 1/4" plug |
| M1783A | 989803105251 | Sync Cable 12- pin 2.5M. (8ft.) | | | | one end is terminated with ECG 12 pin and one end is terminated with phone plug |
| M5526A | 989803129001 | 24' Sync Cable | | | | one end is terminated with ECG 12 pin and one end is terminated with phone plug |

Gas Analyzers

For details on the M1026A/B Anesthetic Gas Module, the M1013A IntelliVue G1 and the M1019A IntelliVue G5, please refer to the respective Service Guides on your documentation CD.

For details on configuring the RS232/MIB port for the gas analyzers, see *Connection of Devices via the MIB/RS232 Interface (Rev. D.00.xx or higher)* in the Installation Instructions chapter

Specifications

Essential Performance

The following tables list the essential performance for IntelliVue Patient Monitors, in combination with Multi-Measurement Modules (MMS), MMS Extensions, the Flexible Module Rack (FMS) and Measurement Modules. The essential performance is maintained under the electromagnetic conditions specified in chapter "Safety and Performance Test" of the Instructions for Use. The essential performance is specified separately for non-transient and for transient electromagnetic interference.

Essential Performance for Non-Transient EMC Phenomena

- Radiated electromagnetic fields
- Conducted disturbances induced by RF fields
- Conducted disturbances induced by magnetic fields
- Voltage dips / voltage variations

| Parameter | Essential Performance |
|---------------------|--|
| General | No interruption or cessation of current operating mode (e.g. no reboot, display ok, ...) No spontaneous operation of controls (e.g. no activation of touch screen controls without user interaction) |
| BIS | Measurement of the bispectral index (BIS). |
| C.O. / CCO | Measurement of the blood temperature within the specified error limits Alarming on blood temperature limit violations. |
| ECG (Monitoring) | Measurement of the heart rate within the specified error measurement ranges (adult, pediatric and neonatal) Alarming on asystolie, or heart rate limit violations as specified in IEC 60601-2-27 :2005. Detection of VFIB and alarming. ECG-Output Signal delay |

10 Specifications

| Parameter | Essential Performance |
|---|--|
| ECG (Diagnostic) | Diagnostic bandwidth |
| | Measurement of the spectral edge frequency, mean dominant frequency, peak power frequency, total power as well as the percentage of power in bands alpha, beta, theta and delta. |
| Invasive Pressure | Measurement of systolic, diastolic and mean blood pressure and pulse rate within the specified error limits. |
| | Alarming on systolic, diastolic and mean blood pressure and pulse rate limit violations. |
| NBP | Measurement of systolic, diastolic and mean blood pressure and pulse rate within the specified error limits. |
| | Automatic cycling |
| | Alarming on systolic, diastolic and mean blood pressure limit violations. |
| Respiration | Measurement of the respiration rate within the specified error limits. |
| | Alarming on apnea and respiration rate limit violations. |
| Respiratory Gases (CO ₂ , O ₂ , N ₂ O, Anesthetic Agents) | Measurement of the respiratory gases within the error limits of EN ISO 21647:2004 clause 51.101.1 |
| | Measurement of the airway respiration rate within the specified error limits. |
| | Alarming on apnea, airway respiration rate and respiratory gas limit violations. |
| Spirometry | Measurement of spirometry values within the specified error limits. |
| | Alarming on apnea, and PEEP, PIP, MV and respiration rate limit violations. |
| SpO ₂ | Measurement of oxygen saturation and pulse rate within the specified error limits |
| | Alarming on oxygen saturation and pulse rate limit violations. |
| Intravascular Oxygen Saturation | Measurement of intravascular oxygen saturation within the specified error limits. |
| | Alarming on oxygen saturation limit violations. |
| tcpO ₂ / tcpCO ₂ | Measurement of transcutaneous partial pressures O ₂ and CO ₂ within the specified error limits. |
| | Alarming on transcutaneous partial pressure O ₂ and CO ₂ limit violations. |
| Temperature, SureTemp+ | Measurement of temperature within the specified error limits. |
| | Alarming on temperature limit violations. |

Essential Performance Transient EMC Phenomena

- Electrostatic Discharge (ESD)
- Electrical Fast Transients / Bursts
- Surges
- Voltage interruptions
- Electrosurgery (ESU)

| Parameter | Essential Performance |
|-----------|--|
| All | After electrostatic discharge, fast transients / bursts, surges and electro surgery interference the equipment will return to previous operation mode within 10 seconds (SpO2 30 seconds) without loss of any stored data. |
| | After voltage interruptions the equipment returns to previous state without operator intervention and loss of data |

MDD Classification

According to the Council Directive 93/42/EEC (Medical Devices Directive) the device classification is 2B, Rule 10.


Safety and Regulatory Information

For details about safety and regulatory information refer to the Instruction for Use.

Appendix - Philips 15210B Calibration Unit

The Philips 15210B Calibration Unit consists of a gas cylinder connected to a gas outlet via a time controlled valve. The valve is normally closed, ensuring that no gas is lost when the unit is not in use. When the unit is set up for use and the timer control knob is turned, gas is directed to the calibration chamber on the tcpO₂/tcpCO₂ module (Philips M1018A) for a period of up to 20 minutes. After this time the valve automatically closes.

This section provides the necessary information for you to install and service the Philips 15210B.



Unpacking the Instrument

If external damage to the shipping carton is evident, ask the carrier's agent to be present when the unit is unpacked.

Initial Inspection

Check the instrument for any external damage such as dents and scratches on panel surfaces. If the shipping carton is not damaged, check the cushioning material and note any signs of severe stress as an indication of rough handling in transit. Retain the packaging material for possible repacking.

Claims for Damage

If physical damage is evident when the Calibration Unit is received or the unit does not meet the specified operational requirements, please notify the carrier and the nearest Philips Sales/Service office

immediately. The Sales/Service office will arrange for repair or replacement without waiting for settlement of the claim against the carrier.

Repacking for Shipment or Storage

If the Unit is to be shipped to a Philips Sales/Service office, securely attach a tag showing the name and address of the owner, the model and serial number, and the repair required or symptoms of the fault. If available and reusable, the original shipping carton and packaging material should be used to provide adequate protection during shipping. The Philips Sales/Service office will provide information and recommendations on materials to be used if the original material is not available or reusable.

Instrument Identification

Philips uses a nine character sequence for instrument identification. The serial number is located on a plate attached to the rear panel of the instrument.

Specification

| | |
|--------------------|---|
| Gas Supply: | 1 low pressure cylinder. |
| Gas Flow: | 8 ml +4/-2 ml per minute for 15210-64010, 12 ml +4/-2 ml per minute for 15210-60010. |
| Cylinder Pressure: | Indicated by an integral pressure manometer. |
| Timer Period: | 20 minutes. |
| Dimensions: | 90mm (35.4in) high x 220mm (86.6in) wide x 235mm (92.5in) deep, (without cylinder). |
| Weight: | 2.4 kg (5.3lbs), (without cylinder). |

NOTE

The 15210B is intended for use with Philips “CAL 1” gas cylinders (part number 15210-60010 or 15210-64010 for Europe and Japan).

Operating Environment

The environment where the Philips 15210B will be used should be reasonably free from vibration, dust, corrosive or explosive gases, extremes of temperature, humidity, etc. The Philips 15210B operates within specifications at ambient temperatures between 0°C and 55°C. The maximum operating relative humidity is 95% at 40°C. Ambient temperatures or humidities which exceed these limits could affect the accuracy of the calibration unit and cause damage to components.

Operating Information

Each Philips 15210B is delivered with a multilanguage collection of stick-on operating labels. Each label summarizes day-to-day operating procedures using the Calibration Unit. It is intended to be stuck to the top surface of the Unit; however, it may be attached to any flat, grease-free surface.

To attach label: Clean the surface where the label is to be placed with soapy water to remove any dirt or grease. Dry the surface thoroughly. Peel off the paper backing and carefully place the label in the required position. Press down firmly with a clean dry cloth, paying particular attention to the edges.

Fitting the Gas Cylinders

When the Calibration Unit is delivered, no gas cylinder is fitted. Before putting the unit into service, screw the cylinder into the opening in the rear panel (See “Routine Maintenance” on page 339” on page , next section).

When new, the calibration unit will contain a small amount of normal air. To expel this air before use and thus prevent inaccurate calibration, turn the timer control fully clockwise after fitting the gas cylinders and allow it to run for the full period. The calibration unit is now ready for use.

Storage of Gas Cylinders

New gas cylinders should be stored in a cool place and not exposed to direct sunlight.

Disposal of Used Gas Cylinders

Do not crush or incinerate used gas cylinders. They may be disposed of as scrap metal.

Routine Maintenance

Changing the Gas Cylinders

- 1 Before each calibration the gas pressure indicator on the Philips 15210B front panel should be read. If the indicator is in the “black” zone, change the gas cylinder as follows:
- 2 From the rear of the unit turn the empty gas cylinder anti-clockwise until the cylinder is free (3-6 turns).
- 3 Withdraw the empty cylinder.
- 4 Take a full gas cylinder and insert it squarely into the rear of the unit. Turn clockwise until hand tight.
- 5 Check that the pressure indicator is no longer in the “black” zone.

Care and Cleaning


Keep the surfaces of the calibration unit clean and free of dust and dirt. Clean regularly with a lint-free cloth or sponge dampened in soapy water. Avoid using alcohol or ammonia based cleaners which may damage the Calibration Unit. Other strong cleaners such as Povidine RR, Lysol R and Mikroklene R are not recommended since they may stain the unit. Do not pour any liquid on the instrument while cleaning. Never use an abrasive material such as steel wool or metal polish. Cleaning agents and disinfectants should only be used in cases of stubborn dirt. If used, carefully remove any remaining traces of cleaning agent or disinfectant with clean water.

NOTE

Do not allow water to enter the gas outlet.

To clean the gas outlet: Use cotton wool soaked in soapy water to remove any deposits which may collect in the outlet. Dry the outlet thoroughly after cleaning. In the case of severe blockages, a thin length of wire may be used to free the outlet pipe.

Theory of Operation



Block diagram - Internal Components

The gas cylinder is screwed directly into a pressure *regulator block*. This block ensures that, in combination with the *restriction*, the gas flow remains constant as the pressure in the cylinders falls with use. From the regulation block the gas is channelled to the *control block*. The gas passes into the control block via an opening in the side sealed with an “O” ring and filter. The control block acts as a switch.


A restriction piece is fitted in the tubing connecting the control block to the gas outlet. The restriction helps to regulate the gas flow.

Gas Flow Performance Check

Philips recommends that the following gas flow check is conducted once a year.

Test Procedure

- 1 Check that the pressure indicator is not in the black zone (i.e. that there is an adequate supply of gas in the cylinder).
- 2 Fit gas tubing to the gas outlet, then take the free end and fit it to a water-filled syringe in a glass of water.
- 3 Turn the timer control fully clockwise and note the volume of water displaced after 60 seconds.



Test Procedure

Action if outside specification


The volume of water displaced in 60 seconds should be 8ml (4/-2ml) for the 15210-64010 gas cylinder or 12ml (+4/-2ml) for the 15210-60010 gas cylinder. If the displacement is within the appropriate one of these ranges, the supply of gas is within specification.

If the gas flow is less than the permitted minimum, remove the Calibration Unit cover (see *Cover Removal* in the next section) and look for an occlusion or leakage.

If the gas flow is greater than the permitted maximum, follow the procedure below:

- 1 Replace the gas cylinder with a new cylinder.

- 2 Turn the *Gas flow adjuster screw* on the underside of the unit to reduce the gas flow (see next figure to locate the gas flow adjuster).



Calibration Unit viewed from underneath


- 3 Repeat the gas flow performance check described above. If the gas flow is still significantly greater than the permitted maximum, remove the cover following the procedure (*Cover Removal*) below.
- 4 Remove the flow restriction by pulling the tubing off (see Figure 3-4), select a new restriction from the set of restriction pieces (part number 15210-68703) and fit in the unit.
- 5 Reassemble the unit and repeat the gas flow performance check.
- 6 Turn the gas flow adjuster screw as necessary.
- 7 If the gas flow is still greater than the specified rate, repeat the above steps, inserting a longer restriction.

Disassembly

Tools Required: Pozidrive screwdriver, size GN1, Normal screwdriver, size 1/7, Hex-key (Allen-key), size SW 3mm.

- 1 Cover Removal
 - a. Remove the gas cylinder from the calibration unit.
 - b. Remove the four screws on the base of the unit (labeled A – see figure).
 - c. Slide the cover off towards the rear of the unit.
- 1 Timer Control Knob
 - a. The timer control knob is secured with a “grub-screw” located in the side of the knob. Loosen this screw approximately 2 turns. The knob can now be pulled off.
 - b. Regulator / Control Block Removal
 - c. Complete operations 1 and 2 above.
 - d. Remove the connection pipe from the rear of the Calibration Chamber.

- e. Unscrew the four remaining screws on the unit base (labeled B in the figure) to release the Regulator / Control Block.
- f. The two screws on the regulator block side can now be removed to separate the regulator block from the control block. Be careful not to misplace the “O” ring and filter which are fitted between the two blocks.



Exploded view - Regulator and Control Blocks

Parts List

Table 16 Replaceable Parts for the 15210B

| Part Number | Description |
|-------------|--|
| 15210-47101 | Flat sealing ring (to seal gas bottle) |
| 15210-47106 | Membrane foil |
| 15210-62401 | Clock |
| 15201-67711 | Valve control block |
| 15210-67202 | Regulator block left |
| 15210-23711 | Flow regulator restriction |
| 0905-0678 | 8mm ring - between valve control block and regulator block |
| 15210-27401 | Timer control knob |
| 0515-0777 | Screw M6x8 (for timer control knob) |
| 15210-04111 | Cover - bottom |
| 15210-04102 | Cover - top |
| 15210-24702 | Spacer - hexagonal nut for mounting regulator bloc |
| 15210-62302 | Gas outlet block |
| M2205A | Calibration tubing (set of 5) |

IntelliVue MX600/700 Product Structure

The following table shows the product option structure for the IntelliVue MX600/700 patient monitors:

| IntelliVue MX600/700 Product Structure Rev. H.1 | MX600 865242 | MX700 865241 |
|--|-------------------------|-------------------------|
| User Interface | | |
| Navigation Point | Basic | Basic |
| Touch | n/a | Basic |
| Remote Control | Optional | Optional |
| Multiple Profiles | Basic | Basic |
| Multiple Screen Layouts | Basic | Basic |
| Full Customization Capability | Basic | Basic |
| Patient Data Management and documentation | | |
| Tabular Trends | Basic | Basic |
| Graphical Trends | Basic | Basic |
| Short Trends (Split Screen Trends) | Basic | Basic |
| High Resolution Trends | Basic | Basic |
| Recorder Support | Basic | Basic |
| Print Support inc. Reports | Basic | Basic |
| Clinical Decision Support | | |
| OxyCRG | Basic | Basic |
| Neonatal Event Review | Optional | Optional |
| Car Seat Test | Optional | Optional |
| Basic Event Surveillance | Optional | Optional |
| Parameter Histograms | Optional | Optional |
| Horizon Trends | Basic | Basic |
| PV Loops | Basic | Basic |
| Information Portal | Basic | Basic |

| IntelliVue MX600/700 Product Structure Rev. H.1 | MX600 865242 | MX700 865241 |
|--|-------------------------|-------------------------|
| Cardiac Applications | | |
| Arrhythmia | Basic | Basic |
| Full Arrhythmia | Basic | Basic |
| ST Segment Analysis | Basic | Basic |
| ST Map | Basic | Basic |
| QT/QTc | Basic | Basic |
| 12 Lead Display Application | Basic | Basic |
| Clinical Applications | | |
| Alarm Review | Basic | Basic |
| Hemodynamic Calculations | Basic | Basic |
| Ventilation Calculations | Basic | Basic |
| Oxygenation Calculations | Basic | Basic |
| Drug Calculator | Optional | Optional |
| Mathematic Calculator | Basic | Basic |
| Interfaces | | |
| Networking (Fast Ethernet) | Basic | Basic |
| Basic Nurse Call | Basic | Basic |
| Remote Control Interface | Optional | Optional |
| Serial Interface ports (RS232) (Standard/Max.) | 1/1 | 1/1 |
| Serial Interface ports (RS232/MIB) for data export (Standard/Max.) | 0/1 | 0/1 |
| USB 2.0 ports Full Speed (Standard/Max.) | 4/4 | 4/4 |
| VGA/DVI Video Out | Basic | Basic |
| ECG Sync Out | 2 ports | 2 ports |
| MSL Interface (Standard/Max.) | 1/1 | 1/1 |
| Measurement Support | | |
| Number of racks supported (Standard/Max.) | 1/1 | 1/1 |
| Number of modules supported (Standard/Max) | 4/4 | 4/4 |
| Number of IBP measurements (Standard/Max.) | 2/4 | 2/4 |
| Number of Temperature Measurements (Standard/Max.) | 2/2 | 2/2 |
| Number of SpO2 Measurements (Standard/Max.) | 1/2 | 1/2 |
| Number of VueLink/IntelliBridge Connections (Standard/Max.) | 2/2 | 2/2 |
| Add-Ons | | |
| Quick Release Mount | Basic | Basic |
| Required Purchases | MX600 | MX700 |
| Application Areas | | |

| IntelliVue MX600/700 Product Structure Rev. H.1 | | | MX600 865242 | MX700 865241 |
|---|--------------------------|---|-----------------|-----------------|
| H10 | General / ICU Config. | Option H10 - General and Intensive Care software. Adds application-specific functionality and configurations | Optional | Optional |
| H20 | Neonatal Config | Option H20 - Neonatal software. Adds application-specific functionality and configurations | Optional | Optional |
| H30 | OR / Anesthesia Config. | Option H30 - OR/Anesthesia software. Adds application-specific functionality and configurations | Optional | Optional |
| H40 | Cardiac Config | Option H40 - Cardiac software. Adds application-specific functionality and configurations | Optional | Optional |
| Waves | | | | |
| A04 | 4-Wave Capability | Option A04 - The addition of 4-wave capability which can be used for the display of realtime waveforms and/or high resolution trends | Optional | Optional |
| A06 | 6-Wave Capability | Option A06 - The addition of 6-wave capability, which can be used for the display of realtime waveforms and/or high resolution trends | Optional | Optional |
| Add-On Options | | | MX600 | MX700 |
| Clinical Applications | | | | |
| C04 | Neonatal CDS Package | Option C04 - The Neonatal CDS package consists of Car Seat Testing, Parameter Histograms and Neonatal Event Review which provides automatic detection, display, storage and documentation of up to 50 neonatal events over 24 hours. | Optional | Optional |
| C05 | Drug Calculator | Option C05 - Drug Calculator | Optional | Optional |
| C06 | Basic Event Surveillance | Option C06 - Basic Event Surveillance allows the user to define events based on the behavior of up to three parameters. When the event triggers, it can be programmed to capture either a waveform snippet (15 seconds) or standard trend (20 minute trend, 12 second resolution) | Optional | Optional |
| C09 | Parameter Histograms | Option C09 - Parameter Histograms provide a time-distribution bar graph of the numerical values during a given time period. X-axis represents parameter range, Y-axis represents % of time. | Optional | Optional |
| Protocol Watch | | | | |
| P01 | Severe Sepsis Screening | Option P01 - The Severe Sepsis Screening option provides sepsis screening capability according to the SSC guidelines. | Optional | Optional |

| IntelliVue MX600/700 Product Structure Rev. H.1 | | | | MX600 865242 | MX700 865241 |
|--|----------------------------------|---------------------|---|-------------------------|-------------------------|
| P02 | | SSC Sepsis Protocol | Option P02 - Includes the Sepsis Screening and SSC Sepsis Bundles (resuscitation and management). | Optional | Optional |
| XDS external display solution | | | | | |
| X00 | XDS Connectivity | | Option X00 - Provides connectivity to a PC running the XDS application. Allows the PC to display independent waveforms, up to the number of waves defined by the Axx option of the host monitor. Includes remote control of the monitor. | Optional | Optional |
| X30 | XDS Clinical Workstation | | Option X30 - Provides ability to launch external IT applications on a PC running the XDS application. Transfers patient context information from the patient monitor to the IT application. | Optional | Optional |
| Measurement Support | | | | | |
| M06 | Support 2 adtl IBPs | | Option M06 extends the capability to measure invasive pressure by two additional lines | Optional | Optional |
| M20 | Support 1 adtl. SpO ₂ | | Option M20 provides the capability to measure two SpO ₂ values | Optional | Optional |
| Hardware Add-Ons | | | | | |
| E00 | Remote Control | | Option E00 - the IntelliVue Remote Control allows single handed control of the IntelliVue patient monitors. It can be connected either wired (USB) or wireless (Short Range Radio). It comes with a cradle, 3.0 meter USB cable, a wall installation kit and a protection for unintended removal. | Optional | Optional |
| E04 | 4-slot rack with MMS mount | | Option E04 - The 8-slot rack with the built-in MMS mounting plate provides 4 module slots for use with individual measurement modules and the recorder. | Optional | Optional |
| iPC Performance Options | | | | | |
| PC0 | Integrated PC (iPC) | | Option PC0 - the integrated PC with passive cooling system is optimized for use in a clinical environment. Comes with operating system and basic applications software pre-installed | Optional | Optional |
| Wired Interfaces | | | | | |
| J13 | | RS232/MIB interface | Option J13 - Provides Serial/MIB interface (1 port). Allows connectivity to RS232 and MIB-ready devices | Optional | Optional |

| IntelliVue MX600/700 Product Structure Rev. H.1 | | | MX600 865242 | MX700 865241 |
|---|---------------------------|--|-----------------|-----------------|
| J30 | Flexible Nurse Call IF | Option J30 - Flexible Nurse Call Relay provides three alarm relays and one power-fail alarm. Both polarities are supported | Optional | Optional |
| Wireless Interfaces | | | | |
| J35 | 802.11 Wireless IF | Option J35 - IntelliVue 802.11 Bedside Adapter consists of the internal radio interface board and the connecting cable. | Optional | Optional |
| J45 | Smart Hopping IF 1.4 GHz | Option J45 - The IntelliVue 1.4 GHz Smart Hopping Interface consists of the internal radio interface board and the connecting cable. | Optional | Optional |
| J46 | Short Range Radio | Option J46 - The Short Range Radio interface consists of the internal radio board and proved connectivity to the remote control (option E00) | Optional | Optional |
| J47 | Smart Hopping IF 2.4 GHz | Option J47 - The IntelliVue 2.4 GHz Smart Hopping Interface consists of the internal radio interface board and the connecting cable. | Optional | Optional |
| Cables | | | | |
| SC1 | MSL system cable - 0.75 m | Option SC1 - The 0.75 m long MSL cable connects IntelliVue patient monitors with the MMS (M3001A), X2 (M3002A), MP5 (M8105A) and the FMS (M8048A). | Optional | Optional |
| SC2 | MSL system cable - 2.0 m | Option SC2 - The 2.0 m long MSL cable connects IntelliVue patient monitors with the MMS (M3001A), X2 (M3002A), MP5 (M8105A) and the FMS (M8048A). | Optional | Optional |
| SC4 | MSL system cable - 4.0 m | Option SC4 - The 4.0 m long MSL cable connects IntelliVue patient monitors with the MMS (M3001A), X2 (M3002A), MP5 (M8105A) and the FMS (M8048A). | Optional | Optional |
| SC6 | MSL system cable - 10.0 m | Option SC6 - The 10.0 m long MSL cable connects IntelliVue patient monitors with the MMS (M3001A), X2 (M3002A), MP5 (M8105A) and the FMS (M8048A). | Optional | Optional |

IntelliVue MX600/700 Upgrade Structure

The following table shows the available upgrade options for the IntelliVue MX600/700 patient monitors:

| IntelliVue MX600/700 Upgrade Structure 86305(MX600) /86304 (MX700) | |
|---|---------------------------|
| Waves | |
| A06 | 4- to 6-Wave Upgrade |
| Clinical Applications | |
| C04 | Neonatal CDS Package |
| C06 | Basic Event Surveillance |
| C09 | Parameter Histograms |
| Protocol Watch | |
| P01 | Severe Sepsis Screening |
| P02 | SSC Sepsis Protocol |
| XDS external display solution | |
| X00 | XDS Connectivity |
| X30 | XDS Clinical Workstation |
| Measurement Support | |
| M06 | Support 2 additional IBPs |
| M20 | Support 1 additional SpO2 |
| iPC Performance Options | |
| PC0 | Integrated PC (iPC) |
| Wired Interfaces | |
| J13 | RS232/MIB interface |
| J30 | Flexible Nurse Call IF |
| Wireless Interfaces | |
| J35 | 802.11 Wireless IF |
| J45 | Smart Hopping IF 1.4 GHz |
| J46 | Short Range Radio |
| J47 | Smart Hopping IF 2.4 GHz |
| Software | |
| SU0 | Current Software Revision |

Part Number 453564322791
Published in Germany 05/12



PHILIPS