Curapuls 670

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

Copyright:

**Enraf
Nonius**

Enraf-Nonius B.V.

P.O. Box 810

NL-2600 AV DELFT

Netherlands

Tel.: +31 (0)15 - 26 98 400

Fax: +31 (0)15 - 25 61 686


Article number 1403.751-42

December 1999

Curapuls 670

Bedieningshandleiding
Operating Instructions
Gebrauchsanweisung
Mode d'emploi
Instrucciones de manejo

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The Curapuls 670 is a microprocessor-controlled unit for pulsed shortwave therapy, applied by means of one or two inductive electrodes. The equipment is only meant to be used by competent personnel in physiotherapy, rehabilitation or adjacent areas.

Simple and safe operation

The various parameters are selected by means of push-button controls, and are set using one central control. The microprocessor controls and monitors all functions during operation. Directly after switching on, the microprocessor automatically tests the electronics of the Curapuls 670 for safety and correct functioning.

Optimum energy transfer

Considerable attention has been paid to the way in which the high-frequency energy is transferred to the tissue. Continuous microprocessor-controlled fine tuning ensures optimum adaptation to the tissue characteristics.

Special high-output inductive electrodes have been developed for use with the Curapuls 670. These are of the "Circuplode" type, and are available in various diameters and shapes. They are provided with screening for the electrical field so that only a magnetic field is created. The type of electrode (large/small/ellipse-shaped) connected to each channel is recognized automatically, eliminating the chance of incorrect settings.

If the load on the electrode is insufficient, no power will be supplied. This is indicated by a warning lamp on the electrode and by the sounding of a buzzer. This form of power transfer indication prevents inefficient treatment, thus contributing to optimum treatment results.

General

The Curapuls 670 is a modern, safe and easy-to-handle apparatus, which has been developed, produced, calibrated and packed with the greatest care.

To operate the Curapuls in the correct and proper manner, it is important that you read this operating manual before use. Also, make sure this manual is available to all personnel at all times. For use of this apparatus, we refer to the therapy books "Shortwave therapy" and "Pulsed short wave therapy with the Curapuls 670".

Enraf-Nonius cannot be held responsible for any consequences to the operator or the patient due to faulty diagnosis, misuse or mishandling of the equipment or its accessories, misunderstanding of the operating instructions, faulty connection(s) with combined equipment or incorrect interpretation of these operating instructions. The consequences of incorrect maintenance are also excluded.

The Curapuls 670 meets the requirements of the international standard for the safety of medical electrical equipment IEC 60601-1 (General standard), IEC 60601-2-3 (Standard for shortwave therapy equipment) and IEC 60601-1-2 (Standard for EMC). According to the IEC 60601-1-2 (1993) directive, clause 36.201.1.8, class A equipment and/or systems (CISPR 11 Classification) are allowed in domestic establishments when used under the jurisdiction of a health care professional.

Service and maintenance

We advise to have the Curapuls 670 checked per service manual test routines and safety tested once each year at your local authorized service office. It is also recommended that a record be kept of all activities relating to service and maintenance. In some countries this is even obligatory. The check up and/or technical maintenance must be carried out conform the procedure described in the service manual of the unit. Opening of the equipment by an unauthorized agency is not allowed and will cancel the right of guarantee. The manufacturer will not be held responsible for damage to equipment or injuries to persons operating or being treated with the equipment that result from service or maintenance performed by companies or persons not authorized by the manufacturer.

USA ONLY

Safety precautions

This unit operates with high voltages. Servicing of the Curapuls 670 should be referred to qualified service technicians or returned directly to the distributor. To assure continued compliance with FDA, 21 C.F.R. 1050.10 standards, the Curapuls 670 should be calibrated and safety tested at least once each

year. Enraf-Nonius recommends that all service work be performed by the distributor. However, a service manual can be obtained by writing the company for those who are qualified to service the device and are fully aware of their responsibilities.

Caution :

Use of controls, adjustments or procedures other than those specified herein may result in hazardous exposure to shortwave energy.

Installation

Good air circulation is essential to prevent internal heat build-up.

- Do not install the unit in a location near heat sources such as radiators or air ducts, and prevent exposure to direct sunlight.
- The unit should be installed in such a way, that it's warmth can be passed on to the surroundings. Don't cover the ventilation openings.
- Furthermore, it is recommended that the Curapuls 670 should not be used after extreme temperature fluctuations.
- We advise against using wireless telephones in the vicinity of the unit.
- When moving or transporting the unit, never pull at the electrodes or the electrode arm, to avoid the risk of tilting the Curapuls 670.
- For safety reasons, the brakes on the wheels of the unit should be locked before the beginning of a treatment.
- The Curapuls 670 is not intended for use in so-called "wet rooms" (hydrotherapy rooms). Should any liquid or foreign object fall into the cabinet, turn off the power, unplug the unit and have it checked by your nearest authorized service dealer before any further operation.

Connection

The Curapuls 670 is a class I, type BF apparatus according to IEC 60601-1. This means that the Curapuls 670 must be connected to a grounded power supply. Use only the power cable supplied with the apparatus.

Prior to connection to the mains, check that the voltage and frequency stated on the Curapuls 670 correspond with the available power supply rating. Check that the available power supply complies with local requirements regarding medical rooms (electrical codes).

General

We strongly advise that you regularly consult the therapy manual "Pulsed short wave therapy with the Curapuls 670". The chapters regarding indications and contra-indications are especially important with respect to safe and effective treatment.

Also pay attention to the following:

- From considerations of safety, only original, undamaged electrode cables and electrodes should be used.
- Electronic hearing-aids must be taken off if damage or malfunctions are to be avoided.
- Patients under treatment must not be able to touch surrounding metal objects such as metal furniture, radiators or window-frames.
- Patients with a reduced thermal sensibility in the area to be treated need special caution. Only very low dosage and very short sessions are recommended.
- Take special notice of the patient's sensations and check the objective reactions (e.g. on the skin).
- It is strongly advised against treatment of locally anaesthetized areas.
- The patient must be in sight of the therapist at all times.

Furniture and clothing

Metal parts in furniture might cause intensity concentrations in the high-frequency field. Therefore only use metal-free wooden furniture.

For all treatments with shortwave it is advisable to undress the patient, and to apply a towel between the skin and the electrode (against perspiration and to avoid them making contact). There are two advantages in treating patients this way: firstly to ensure no metal objects are in the high-frequency field (metal parts in the clothing or in ornaments under the clothing), secondly to prevent undesired concentrations of energy due to perspiration or for example the wearing of nylon or leather clothing.

Implanted electronic devices

It is not allowed to treat a patient with an implanted electronic device (for example a cardiac pacemaker) with shortwave equipment; they are not even allowed in the vicinity of an operating shortwave unit.

Implanted metals

Metal concentrates high-frequency energy. To prevent any undesirable concentrations around a metal implant and the resulting dangers (burns), short-wave therapy with a thermal effect should basically not be used.

Consult the therapy manual "Pulsed short wave therapy with the Curapuls 670", the chapter on 'Contra-Indications'.

Combination therapy

Since the Curapuls 670 radiates a high-frequency field that can affect electronic circuits inside nearby equipment, you are strongly advised against any sort of combination treatment using a second (electro medical) apparatus.

High-frequency interference

The Curapuls 670 complies with the present CE-demands regarding HF-emission and immunity.

The Curapuls 670 produces an high-frequency field that is designed to generate thermal and biological effects in tissues and joints of the human body. In free air, the strength of this high-frequency field generated from the electrodes, cables and treated patient is reduced quadratically with the distance. This high-frequency field can affect electronic circuits of equipment positioned in the vicinity.

This high-frequency interference is in almost all cases the cause of so-called direct radiation on the statutory permitted frequency of 27.12 MHz. Despite the fact that the affected equipment is not tuned to this frequency, detection appears in several circuits due to the strong radiation, causing the modulation of the high-frequency wave to become audible and/or visible (through the neon checklight). With low- and medium-frequency physiotherapy equipment this modulation might become perceptible to the patient as well.

The susceptibility to high-frequency interference of equipment from various manufacturers can be very different. This is due to the fact that some manufacturers are still shielding their equipment insufficiently against external high-frequency fields. Should interference arise, we advise consulting the manufacturer concerned.

We recommend paying attention to the following items in order to minimize the chance of high-frequency interference.

- The larger the distance between the shortwave unit and other equipment, the smaller the chance of high-frequency interference. Place your shortwave unit, including the patient under treatment, at least 2 metres from surrounding equipment and/or oblong metal objects.
- Use a separate mains group (phase) for the shortwave apparatus. Keep mains leads wide apart.

- Take care that the mains lead of the shortwave apparatus does not come into the vicinity of the patient and/or electrode leads, or of mains leads from other (electro-medical) equipment.
- Use only a completely metal-free treatment chair or couch.

The metal fixation plate at the broad end of the electrode arm(s) is to be connected to the back of the Curapuls by means of the supplied hollow-head screw wrench. You can connect the arm with the knobs either towards the outside, or towards the unit.

The fixation knobs can be adjusted by hand. We advise you to adjust so tight that the arm can still be adjusted, without having to untighten the knobs.

In the top and bottom of the broad side of the electrode arm are holes in which you can click the (supplied) cable guides. If desired, you can also apply the blind plugs.

Should the purple-ish shell be removed from the small end of the arm, replace it and fix it by means of the fixation knob at the small end of the arm, before mounting a Circuplode. You cannot mount a Circuplode, when the purple-ish shell is not in its place.

(See fold out page at front)

1 Illuminated LCD display

2 Start/Stop key

With this key you start/stop the treatment.

3 Enter key

Herewith you confirm the choice you made.

4 Replacement keys ("arrow keys")

With these keys you can move the selection bar in the display to the option desired.

5 Central control

With this regulator you can change the parameters in the display.

6 On/off switch

7 Ventilation openings

The ventilation openings in the housing provide adequate air circulation for dissipating heat generated by the electronics. It is important to ensure that these openings are never covered.

8 Connection for Circuplodes

The electrode is connected to this connection via the HF cable. The length and impedance of this cable are of critical importance. For this reason, only use the cable supplied by Enraf-Nonius.

Connection of any cable other than the one supplied for the purpose by Enraf-Nonius can adversely affect the safety of the patient and the correct functioning of the equipment, and is therefore not permitted.

9 Attachment points for the electrode arms (4x)

10 Mains input with fuses

The mains cable supplied is connected to this socket, which is fitted with two mains fuses. For replacement of defective fuses, see chapter "Fault conditions".

11 Connection for potential equalization cable

This connection is used in areas where local regulations require connection of a potential equalization cable. A special cable can be supplied for this purpose. Please refer to our catalogue.


12 Type number plate

This plate provides all information for identifying the equipment, such as type, version and serial numbers (for service, guarantee etc.) together with connection data such as mains voltage and current consumption.

Pulse power and mean power

Pulse power (P_p)

The pulse power is the peak power of the pulse, i.e. the height of the pulse as shown in the diagram at the right.

Mean power (P)

The mean power is determined by the pulse power, the pulse duration and the pulse frequency. The mean power resulting from the parameters set can be displayed by switching from the pulse power selector to the mean power selector.

The figure below indicates the mean power at different pulse frequency and pulse width at a pulse power of 200 Watt. The grey columns indicate the values with use of Circuplode-E. The grey columns are not applicable for the Circuplodes 90 and 140 mm. These columns indicate additional values for the Circuplode-E.

Hz μs	26	35	46	62	82	110	150	200	300	400	600	800
65	0,34	0,46	0,60	0,80	1,06	1,44	1,96	2,60	3,90	5,20	7,80	10,40
82	0,42	0,58	0,76	1,02	1,34	1,80	2,46	3,28	4,92	5,56	9,84	13,10
110	0,58	0,78	1,02	1,36	1,80	2,42	3,30	4,40	6,60	8,80	13,20	17,60
150	0,78	1,06	1,38	1,86	2,46	3,30	4,50	6,00	9,00	12,00	18,00	24,00
200	1,04	3,40	1,84	2,48	3,28	4,40	6,00	8,00	12,00	16,00	24,00	32,00
300	1,56	2,10	2,76	3,72	4,92	6,60	9,00	12,00	18,00	24,00	36,00	48,00
400	2,08	2,80	3,68	5,12	6,56	8,80	12,00	16,00	24,00	32,00	48,00	64,00

Three different types of electrodes are available for the Curapuls 670.

With the **Circuplode diam. 140 mm** a high in-depth radiation in small regions can be achieved. It is the most universal and therefore the most used emitter.

The **Circuplode diam. 90 mm** is perfectly suited for treatment of small joints (e.g. the wrist joint) and in children's physiotherapy.

The **Circuplode-E (ellipse-shaped Circuplode)** is especially for for treatment of larger parts of the body, such as the hipjoint and large oblong muscles. The Circuplode-E can generate a frequency of 800 Hz. Ideal for thermal treatment.

With most pre-programmed treatment protocols an advice is given regarding the best suited Circuplode for the selected pathology.

With these electrodes a towel should be placed between the skin and the electrode, after clearing the area to be treated.

In the Curapuls 670, pulse duration, pulse frequency, power and treatment duration are adjustable. When two electrodes are used, the selected pulse frequency and treatment duration apply to both electrodes, while the other parameters can be set separately for each electrode.

Connect the necessary electrodes and switch the equipment on by means of the on/of switch. After switching on and performing the self test of the Curapuls 670, the information screen appears stating the software and language version.

After a few seconds the main menu appears with the following possibilities:

- Set parameters
- Treatment memory
- Treatment protocols

System settings

If you are in the main menu, you can set the standard settings of the Curapuls. To do so, keep the "arrow left" key pressed for 2 seconds. The sub screen appears with the following possibilities:

Language

By turning the central control you can choose the desired language.

By pressing the "arrow down" key, you go to the next option.

Contrast

By turning the central control you can choose the desired contrast.

Sound

You can choose between "on" and "off" by turning the central control. If you choose "on", you will be hearing a beep tone every time you press an "arrow" key, the Enter key or the Start key.

Mode

Here you can choose between a conventional screen and a screen only indicating the mean power.

Screen

You can choose between "normal", i.e. dark letters on a light background or "inverse" (light letters on a dark background).

After you have set your preferences, use the "arrow down" key to return to the main menu. Press the Enter key.

Setting the parameters

You can reach "Setting the parameters" in the main menu by placing the selection bar on it by means of the "arrow" keys and then pressing the "Enter" key.

The letters A and B in the left and right upper corner of the upper bar indicate that 2 induction electrodes have been connected. If only 1 electrode had been connected, only the letter A or B would be shown.

In the middle section of the screen you find the set parameters. You can adjust the desired parameters by means of the central control. With the 4 "arrow" keys, you can go to the various parameters and finally to the lower bar. If no Circuplode had been connected, the left and right column of this middle section will remain blanc.

If an electrode is insufficiently loaded during treatment, this is indicated by the fact that the cursor bar will light up and a beep will be heard. No power is applied until the loading is correct. A minimum "detector" power is used for measuring the load.

Below right is the number of the memory place. How you store a program, is described hereunder.

The presented parameters can be changed by first going to the parameter section with the "arrow up" key. There you can go to each parameter by means of the "arrow" keys and you can change them by turning the central control.

The frequency of 800 Hz can only be set if only one Circuplode-E is connected.

Should a 100 W Circuplode be connected and should the peak power be set to 200 W, the Curapuls 670 will automatically reduce the peak power to 100 W as soon as the treatment is started.

The intensity can be set in 5 Watt divisions.

The time can be set per 30 seconds, with a maximum of 30 minutes.

If you go to the "<<" sign in the lower left corner by means of the "arrow left" key and press the Enter key, you return to the main menu.

Storing a program

If you wish to store the program after you have changed the parameter settings, go to the "store" symbol in the lower right corner by means of the "arrow" keys. With the central control you can choose the desired number for storage. A "*" behind a number means that a program has already been stored under this number. If you push the Enter key, the already stored program will be overwritten, without further warning.

To clear a memory place number, you put all setting to their minimal values (65 μ s, 26 Hz, 0W and 0 min.) and you store these settings in the desired memory place.

With the "arrow down" key you go back to the main menu, from where you can proceed to the treatment memory or the treatment protocols. Push the Enter key.

Treatment memory

You choose a program number by turning the central control.

Only the numbers of the memory places where a program has been stored are shown. There is a differentiation between conventional mode and mean power. Should you, for example, have programmed 24 programs for conventional mode and 1 for mean power mode, you only get to see one program number when you wish to call a program for mean power.

Also the parameters of already stored programs can be changed and re-stored.

If you change a parameter, the text in the upper bar changes in "Treatment parameters" and next to the protocol number the "store" symbol appears, with behind it a number.

You can overwrite the present program (go with the "arrow" keys to the most right "store" symbol and push the Enter key) or assign a new memory place number (first choose a new number with the central control and then push the Enter key).

Treatment protocols

By means of the "arrow down" key or the central control you can choose from;

- Arthralgies
- Neuralgies
- Dermatology
- Vasomotoric insufficiency
- Insertion tenopathies
- Myalgies
- Soft tissue lesions

For an overview of the protocols, we refer to the appendix.

If you wish to go back to the protocol screen, you can press the "arrow left" key at any time.

You choose the desired protocol by going to it with the "arrow up" or "arrow down" key or with the central control and then pushing the Enter key. This also applies to choosing an underlying protocol.

The name of the protocol chosen by you appears in the top bar. It is also indicated which Circuplode connections are in use.

In the middle section you find the parameters.

In the middle of the lower bar, an advice is given which Circuplode to connect.

CPL90 stands for a Circuplode 90 mm.

CPL140 stands for a Circuplode 140 mm.

CPL-E stands for the ellipse-shaped Circuplode-E.

In the lower right corner you find the number of the protocol chosen by you. You can also look at other protocols by making sure that you stand on this symbol and then turning the central control.

If you go to the "<" sign in the lower left corner by means of the "arrow left" key and press the Enter key, you return to the main menu.

The treatment starts after you push the Start key.

The clock starts counting down as soon as one Circuplode makes contact.

Indications

Specific indications for pulsed shortwave therapy are:

Post-traumatic conditions, for example

- A sprain or dislocation
- Contusion
- Rupture
- Fracture
- Hematoma
- A graze

Post-operative conditions, for example

- jaw-, foot-, and hip operations

Inflammation, for example

- chronic osteitis
- bursitis (possibly with calcification)
- sinusitis

Peripheral blood circulation disorders

Woundhealing

Absolute contra-indications

- Malignant tumors
- Pacemakers
- Pregnancy
- Tuberculosis
- Fever
- Bacterial infections

Relative contra-indications

- Implanted metals
- Serious arterial and venous circulation disorders such as atherosclerosis, thrombosis etc.
- Heart disorders
- Acute infectious diseases, acute inflammation
- Thermic sensibility disorders

Before cleaning and other maintenance of the Curapuls 670, first switch off the unit and disconnect it from the mains supply.

The control panel and the housing can be cleaned with a damp cloth. The electrodes can be disinfected with a cloth moistened with 70% alcohol. The electrodes must not be sterilized by exposure to steam (autoclaving) or gas.

Check the electrodes and electrode cables regularly for damage and poor contacts. Damaged electrodes or electrode cables must not be used. The length of the cable is of critical importance. For this reason, only use cables supplied for the purpose by Enraf-Nonius.

End life span apparatus and accessories

Your Curapuls and accessories contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart these articles and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

Please ensure that you are well informed of the local rules and regulations regarding to the removal of apparatus and accessories.

Self-test

Immediately after switching on, the unit automatically carries out a self-test in which the built-in microprocessor checks a number of vital functions of the equipment.

After the self-test, the information screen appears which states the software and language version. In case of a fault you can not remedie yourself, please mention these data to the service department.

If a fault condition is detected during the self-test or during treatment, this will be signalled on the display, and the equipment will cease to function. A fault condition can occur as the result of external factors, such as an interruption of the mains supply. In such cases, the fault condition can be cleared by switching the unit off and on again. After switching on, the unit carries out another self-test after which, if the vital functions are found to be in order, the unit is again ready for use.

If the Curapuls 670 does not function correctly, or at all, the following can be checked to possibly trace the fault.

Display does not light up

- Check that the mains plug is correctly inserted in the wall socket.
- Check the mains supply to the wall socket.
- Check the fuses in the mains output and replace if necessary.

Treatment cannot be started

- Check that the treatment time has been set.
- Check that an electrode has been connected.
- Try another cable or electrode.

If these checks do not reveal the cause of the fault condition, the unit is probably defective. In this case, please contact your Enraf-Nonius authorized supplier.

Replacing fuses

Defective fuses can be replaced by the user. Switch off the unit. Remove the mains plug from the wall socket, and remove the mains connector from the input socket. Unscrew the fuse holder, e.g. using a coin. Replace the defective fuse with one of the same type (T6.3H250V). Screw in the fuse holder until it is finger-tight, and reconnect the mains cable.

Opening of the equipment and attempts at repair by unauthorized persons can adversely affect the safety of the equipment and are therefore not permitted. In such cases Enraf-Nonius reserve the right to cancel the manufacturer's guarantee.

Frequency	: 27.12 MHz ($\pm 0.6\%$)
Channels	: 2
Pulse duration	: 65, 82, 110, 150, 200, 300, 400 μs ($\pm 5\%$)
Pulse frequency	: 26, 35, 46, 62, 82, 110, 150, 200, 300, 400, 600 and 800 Hz ($\pm 5\%$)
Pulse power	: 0 - 200 W ($\pm 20\%$) into 50 ohm.
Mean power	: 0 - 64 W ($\pm 20\%$) into 50 ohm.
Treatment time	: 0 - 30 minutes ($\pm 5\%$)

TECHNICAL DATA

Mains voltage	: 110-120 V, 220-240 V ($\pm 10\%$)
Mains frequency	: 50-60 Hz
Current consumption	: max. 1.6 A (at 230V)
Electrode cable	: 50 ohm, 150 cm
Medical device classification	: IIb (according to Medical Device Directive (93/42/EEC))
Safety class	: I, type BF
Earth resistance	: typically 0.12 ohm (IEC requirement ≤ 0.2 ohm.)
Earth leakage current	: typically 200 μA (IEC requirement ≤ 500 μA)
Earth leakage current single fault condition	: typically 400 μA (IEC requirement ≤ 1000 μA)
Dimensions	: 39 x 40 x 93 cm (wxdxh)
Weight (unit)	: ca. 45 kg (excluding arm)
Electrode arm	: ca. 2.5 kg

Environment conditions for transport and storage

Environment temperature	: -10° till $+50^{\circ}$ C
Relative humidity	: 10 till 95 %
Atmospheric pressure	: 500 till 1060 hPa

Environments conditions for normal use

Environment temperature	: 10° till 40° C
Relative humidity	: 10 till 90 %
Atmospheric pressure	: 500 till 1060 hPa

Classification**Medical class IIb**

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

International safety standard IEC 60601-1

Non-ionising radiation

I

The unit must be connected to a grounded wall socket.

The leakage current and the ground connection of the equipment meets the requirements of IEC 60601-1 and IEC 60601-2-3.

For the ordering data of the Curapuls 670, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

English

Arthralgias

- arthritis, acute
- arthritis, chronic
- capsulitis, acute
- capsulitis, acute
- carpal tunnel syndrome
- cartilago inflam. patellae
- discopathie
- frozen shoulder
- meniscopathie
- sacro iliacale syndromes
- spondylitis
- whiplash syndroom

Neuralgias

- brachialgia
- causalgia
- entrapment neuropathy
- phantom pain
- herpes zoster
- intercostal neuralgia
- meralgy paraesthetica
- thoracic outlet syndrome
- pseudo radricula syndrome
- ischialgy

Dermatology

- decubitis, stage 1&2
- decubitis, stage 3&4
- keloid forming
- scarpain
- skin-flap
- ulcus crurisselectionar

Vasomotoric insufficiency

- arteritis
- intermittent claudication
- raynauds phenomenon
- perfusion problems
- veinary insufficiency
- secondary lymfoedema
- symp. reflex dystrophy

Insertion tendopathies

- achillodiny
- epicondylitis med/lat
- ligamentary, acute
- ligamentary, chronic
- tendinitis, acute
- tendinitis, chronic

Myalgias

- hypertonic
- myalgia, acute, deep
- myalgia, acute, superficial
- myalgia, chronic, deep
- myalgia, chronic, superf.
- myofascial triggerp., dull
- myofascial triggerp., sharp
- tension headache

Soft tissue laesions

- bursitis, acute
- bursitis, chronic
- muscle rupture, acute
- muscle rupture, chronic
- lumbago, acute
- lumbago, chronic
- surmenage syndrome, acute
- surmenage syndrome, chronic

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