# Eltrac 471

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo




Copyright:

Enraf-Nonius B.V.  
P.O. Box 810  
2600 AV DELFT  
The Netherlands  
Tel.: +31(0)15 - 26 98 400  
Fax: +31(0)15 - 25 61 686

Article number : 1471.750-44  
April 2001

# Eltrac 471

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo

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**Chapter 2 - Safety ..... 17**  
 2.1. Preliminary notes ..... 17  
 2.2. General ..... 17  
 2.3. Product liability ..... 17  
     2.3.1. *Limitations of liability* ..... 18

**Chapter 3 - Indications and Contra-indications ... 18**  
 3.1. Indications ..... 18  
 3.2. Contra-indications ..... 18

**Chapter 4 - Installation ..... 19**  
 4.1. Connection ..... 19  
 4.2. Connection mains cable ..... 19  
 4.3. Switching on and self test ..... 19  
 4.4. Installation ..... 19  
 4.5. Electromagnetic interference ..... 19

**Chapter 5 - Operation ..... 20**  
 5.1. Controls ..... 20  
 5.2. Introduction to operation ..... 23  
     5.2.1. *Parameters* ..... 23  
     5.2.2. *Displays* ..... 23  
     5.2.3. *Parameter setting* ..... 24  
 5.3. Sequence for setting parameters ..... 24  
 5.4. Operating sequence example ..... 25

**Chapter 6 - Maintenance (by the user) ..... 25**  
 6.1. Technical maintenance ..... 25  
 6.2. Cleaning of the apparatus ..... 25  
 6.3. Nylon cord ..... 25  
 6.4. Patient stop switch ..... 25  
 6.5. Environment information ..... 25

**Chapter 7 - Fault conditions ..... 26**  
 7.1. Self-test ..... 26  
 7.2. Displays do not light ..... 26  
 7.3. Dashes (-- --) appear on the display ..... 26  
 7.4. Treatment cannot be started ..... 26  
 7.5. Traction force [2] cannot be adjusted to a higher value/ LED flashes ..... 26  
 7.6. The base force [3] cannot be set ..... 26  
 7.7. Early termination of treatment ..... 26  
 7.8. Replacing fuses ..... 26

**Chapter 8 - Specifications ..... 27**  
 8.1. Product specifications ..... 27  
 8.2. Technical data ..... 27

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## 1.1. General

The ELTRAC 471 is a microcomputer controlled unit for continuous and intermittent traction. The equipment is only meant to be used by competent personnel in physiotherapy, rehabilitation or adjacent areas.

The newly designed control panel has integrated push buttons for selecting such parameters as traction force, hold time and treatment time. The selected parameters are set by rotating a single control knob. All parameters can be set before or after the start of treatment; the microcomputer maintains the correct treatment sequence. Parameters can also be changed during treatment making it possible, for example, to change from continuous to intermittent traction.

The cord used to apply the traction can be easily adjusted to the required length. A special tensioning system ensures that a minimum tension of 15 N (1.5 kg) is always applied to the cord. This enables the cord to be paid out easily, and ensures that it rolls up automatically when released.

Operation is, in all respects, versatile and simple.

## 1.2. Finally

You have made a wise choice in selecting the ELTRAC 471. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorised distributor.

## 2.1. Preliminary notes

It is important that you read these operating instructions carefully before using the ELTRAC 471. Please make sure that these instructions are available to all personnel who operates the equipment.

Pay attention to the following before using the ELTRAC 471:

1. Keep yourself informed of the contra-indications (see chapter 3).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave equipment.
3. Never leave the patient unattended during treatment.
4. The apparatus may not be used in so-called "wet rooms" (hydrotherapy rooms).

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

## 2.2. General

In the ELTRAC 471 the traction force is controlled and monitored by the built-in microprocessor. The tension in the cord is measured electronically, and compared with the set value. If these values differ, the motor equalizes the difference. In this way, the patient's movements are taken into account, and the tension in the cord remains constant.

A further advantage of the electronic traction control is that the speed of applying and releasing the traction force can be adjusted. In addition to the electronic safety system, the unit is also provided with a mechanical traction limiting system, and a patient stop switch which enables the patient to terminate the treatment.

All functions of the unit are controlled and monitored by the built-in microcomputer, ensuring a high degree of reliability and safety.

## 2.3. Product liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 7 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

**2.3.1. Limitations of liability**

To the maximum extent permitted by applicable law, in no event will Enraf-Nonius or its suppliers or resellers be liable for any indirect, special, incidental or consequential damages arising from the use of or inability to use the product, including, without limitation, damages for loss of goodwill, work and productivity, computer failure or malfunction, or any and all other commercial damages or losses, even if advised of the possibility thereof, and regardless of the legal or equitable theory (contract, tort or otherwise) upon which the claim is based. In any case, Enraf-Nonius's entire liability under any provision of this agreement shall not exceed in the aggregate the sum of the fees paid for this product and fees for support of the product received by Enraf-Nonius under a separate support agreement (if any), with the exception of death or personal injury caused by the negligence of Enraf-Nonius to the extent applicable law prohibits the limitation of damages in such cases.

Enraf-Nonius can not be held liable for any consequence resulting from incorrect information provided by its personnel, or errors incorporated in this manual and / or other accompanying documentation (including commercial documentation)

The opposing party (product's user or its representative) shall disclaim Enraf-Nonius from all claims arising from third parties, whatever nature or whatever relationship to the opposing party.

**3.1. Indications**

- Spondylosis
- Spondylarthrosis
- Spondylosis thesis
- Arthrosis deformans
- Disc prolapse
- Disc protrusion
- Disc degeneration
- Subluxation
- Cervical syndrome
- Shoulder-hand syndrome
- Costoclavicular syndrome
- Scalenus syndrome
- Scoliosis
- Transmandibular joint dysfunction
- Myalgia
- Osteochondrosis
- Osteophytosis

**3.2. Contra-indications**

- Malignant disorders
- Medullar compression
- Inflammation
- Osteoporosis
- Hypertension
- Rheumatoid arthritis
- Pregnancy

#### 4.1. Connection

- Mains supply connections must comply with the national requirements regarding medical rooms.
- Prior to connection of this apparatus to the mains supply, check that the voltage and frequency stated on the type plate correspond with the available mains supply.

#### 4.2. Connection mains cable

- Connect the supplied mains cable to the connection mains cable [21].
- Connect the mains cable to a grounded safety socket.

#### 4.3. Switching on and self test

- Switch on the unit using the On/Off button [1].
- Immediately after switching on, the unit carries out a self test.

##### Check that:

- All LED's and displays light up briefly.
- At the end of the test a beep can be heard.

Contact your supplier if this is not the case.

#### 4.4. Installation

- Do not install the unit in a location near to a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- Cordless telephones should not be used in close proximity of the unit.
- This unit should not be used in so-called "wet-rooms" (hydrotherapy rooms).
- The unit has to be installed in such way that liquid cannot enter.
- Connection of accessories other than ones specified by the manufacturer can adversely affect the safety of the patient and correct functioning of the equipment, and is therefore not permitted.

Should any liquid ingress into the housing of this equipment, unplug the unit from the wall socket (if connected) and have it checked by an authorised expert.

#### 4.5. Electromagnetic interference

- Operation in close proximity (i.e. less than 2 metres) to working shortwave or microwave therapy units may produce instability in the ELTRAC 471 output.
- To prevent electromagnetic interference, we strongly advise using separate mains groups (phases) for the ELTRAC 471 and the shortwave or microwave equipment, and a distance of at least 2 metres is kept between the ELTRAC 471 and the shortwave or microwave equipment.
- Ensure that the mains cable of the shortwave/microwave unit does not come near the Endomed or the patient.

If problems with electromagnetic interference persist, please contact your authorised distributor.

## 5.1. Controls

**[1] Mains (Power) switch**

Switch ON: press the switch to the right (I).  
Switch OFF: press the switch to the left (0).

**[2] Selector for traction force**

When this push-button is pressed, a green LED in the push-button lights to indicate that the traction force parameter has been selected.

The traction force can now be set by rotating the central control knob [13] and the set values indicated on display [8]. The traction force can be set between 15 N (Newton) and 900 N. If the traction force to be set exceeds 200 N, the > 20 acceptance push-button [10] must be pressed.


**Operation during treatment**

The LED in the selector [2] is lit: the traction force has been selected by the user and can be set as required.

The LED in the selector [2] flashes: the traction force has not been specifically selected, and can only be set to a lower value.

**Remark**

This latter situation occurs after another parameter has been set. 10 seconds hereafter, the unit automatically switches back to the traction force parameter but, for safety reasons, the (flashing) value can only be reduced. Any attempt to increase the traction force value in this mode will be blocked. This will be indicated on the display and by an acoustic warning signal.

**[3] Selector for traction force hold time**


When this push-button is pressed, a green LED in the push-button lights to indicate that the traction force hold time parameter has been selected. The

hold time can now be set by turning the central control knob [13] and the set value is indicated on the display [8]. The hold time can be set between 0 and 60 seconds in steps of 2 seconds between 2 and 10 seconds, and in steps of 5 seconds between 10 and 60 seconds.

**[4] Selector for base force (intermittent traction)**


When this push-button is pressed, a green LED in the push-button lights to indicate that the base force parameter has been selected. Selection of the base force serves to select the intermittent traction mode. The base force is set by turning the central control knob [13] and the set value is indicated on the display [9].

Since the base force can only be set to a value lower than the traction force [2], the traction force must be set first. Therefore, the display will initially show the traction force value.

**[5] Selector for base force hold time**

When this push-button is pressed, a green LED in the push-button lights to indicate that the base force hold time parameter has been selected.

The base force hold time can now be set by turning the central control knob [13] and the set value is indicated on the display [9]. The hold time can be set between 2 and 60 seconds in steps of 2 seconds between 2 and 10 seconds, and steps of 5 seconds between 10 and 60 seconds.

**[6] Selector for changeover speed**

When this push-button is pressed, a green LED in the push-button lights to indicate that the changeover speed parameter has been selected.

The changeover speed can now be set by turning the central control knob [13] and the set value is indicated on the bar display [7]. During intermittent traction the changeover speed is the speed at which the traction force changes to base force and vice versa. During static traction the changeover speed is the speed at which the traction force is built up at the beginning, and reduced at the end of treatment.

The changeover speed can be set to a choice of 10 levels. If no changeover speed is set after pressing the push-button [6], only the lowest LED of the bar display will light. This corresponds to a minimum changeover speed. If a higher level is set, more LED's will light. At maximum setting, all 10 LED's are lit.

**[7] Bar display****Setting changeover speed**

After selection of changeover speed [6] the bar display shows the changeover speed set. If only the lowest LED is lit, the changeover speed is set to the minimum; if all 10 LED's are lit, the changeover speed is set to the maximum.

**Display mode during treatment***Static traction*

In static traction, only the top LED [8] is lit during treatment. During the build up of traction, directly after the beginning of treatment, the top LED flashes until the traction in the cord has reached the set value.

### *Intermittent traction*

In intermittent traction, the bar display indicates the traction in the cord in relation to traction force and base force. If the traction in the cord is equal to the traction force set (indicated on display [8]), the top LED is lit.

If the traction in the cord is equal to the base force set (indicated on display [9]), the lowest LED is lit. All intermediate values are indicated by the LED's in between.

### **Changing the changeover speed during the treatment**

During treatment, the bar display normally indicates the traction phase. The changeover speed set can be indicated by pressing the changeover speed push-button [6], and can be adjusted by turning the central control knob [13]. The changeover speed value will remain on the display for 10 seconds, after which it automatically returns to indicating the traction phase.

#### **[8] Display for traction force/hold time**

This display indicates the traction force set, or the hold time. Display of the traction force set can be selected with push-button [2]; display of the hold time of the traction force can be selected with push button [3].

The indicated traction force must be multiplied by 10 to give the value in Newtons, i.e. '20' on the display is equivalent to  $20 \times 10 = 200$  N.

### **Display mode during treatment**

During treatment, the display normally shows the set traction force [2]. If the hold time push-button is pressed, the hold time will be displayed for 10 seconds. The hold time can now be adjusted by turning the central control knob [13]. The new hold time value will be shown for 10 seconds, after which the display automatically returns to indicating the traction force.

#### **[9] Display for base force/hold time**

In intermittent traction, this display indicates the base force set or the hold time. Display of the base force can be selected with push-button [4]; display of the hold time of the base force can be selected with push-button [5].


The indicated base force must be multiplied by 10 to give the value in Newtons, i.e. '20' on the display is equivalent to  $20 \times 10 = 200$  N.

The display does not light after switching on, or after termination of treatment. It only lights during intermittent traction, when the push-buttons base-force or base-force hold time are pressed (see [4], [5]).

### **Display mode during treatment**


During treatment, the display normally shows the set base force. If the hold time push-button [5] is pressed, the hold time will be displayed for 10 seconds. The hold time can now be adjusted by turning the central control knob [13]. The new hold time value will be shown for 10 seconds, after which the display automatically returns to indicating the base force.

#### **[10] Selector >20(selection of forces exceeding 200 N)**

This push-button is used for setting traction forces exceeding 200N. These can only be set if the > 20 push-button has been pressed. In other words, a high traction force has to be specifically selected. To indicate that the > 20 selector push-button has to be pressed, a yellow LED in the push button flashes and the buzzer sounds briefly. After selection, the LED is lit continuously.

#### **[11] Treatment time selector**

When this push-button is pressed, a green LED in the push-button lights to indicate that the treatment time parameter has been selected.

The treatment time can now be set by turning the central control knob [13] and the set value is indicated on display [12]. The treatment time can be set between 0 and 60 minutes in steps of 1 minute.

The treatment time must be set before treatment can start.

#### **[12] Treatment time display**

This display indicates the treatment time set (or remaining) in minutes. The maximum treatment time that can be set is 60 minutes. A flashing dot indicates that the timer is running. After the treatment time has elapsed, a buzzer sounds for 4 seconds and the traction force is reduced to the minimum value at a speed corresponding to the changeover speed [6], [7].

**[13] Central control knob**

All parameters are set with the rotary central control knob:

traction force, traction force hold time, base force, base force hold time, changeover speed, treatment time.

The value can be reduced by turning the knob to the left, and increased by turning the knob to the right. The knob has no start or stop positions and is coupled to the parameter indicated by a lit green LED. The set value of the parameter is indicated on the corresponding display. The range of settings is determined by the parameter selected (see also remark under [2]).

**[14] Start/stop selector**

This selector is used to start and stop treatment.

Treatment can only start when the following conditions are met:

- the patient stop switch is connected,
- the treatment time is set.

The start of treatment is indicated by the LED in the selector [15] changing from red to green. If condition a is not met, the buzzer sounds when the start/stop push-button is pressed.

Treatment can be interrupted before the end of the treatment time by pressing the start/stop push-button. The traction force is then reduced to the minimum at a speed corresponding to the changeover speed. The LED in the push-button then changes from green to red.

**[15] Start/stop LED**

Before treatment starts, or if treatment is interrupted, the LED in the start/stop selector push-button [14] lights red. When treatment starts the LED lights green.

**[16] Socket for patient stop switch**

The patient stop switch must be connected to this socket (see [20]).

Treatment cannot start unless the patient stop switch has been connected. The switch should be handed to the patient, and its use should be explained.

If it is pressed during treatment, the traction force will be reduced to the minimum value of 15 N at the maximum changeover speed, and the buzzer will sound. The unelapsed treatment time continues to be indicated on the display,

Connection of a patient stop switch other than that prescribed by ENRAF-NONIUS can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

**[17] Socket for remote control**

A remote control unit, which replaces the function of the central control knob [13], can be connected to this socket. The remote control unit is an optional accessory (please consult our catalogue).

Connection of a remote control other than that prescribed by ENRAF-NONIUS can adversely affect the safety of the patient and the functioning of the unit, and is therefore not permitted.

**[18] Pulley**


The pulley ensures optimum guidance of the cord for traction at any required angle.

**[19] Cord**

The cord can be wound into the ELTRAC 471, or paid out, thus increasing or decreasing the traction force when the cord is attached to the patient. The free end of the cord is provided with a carabin hook for easy attachment of accessories.

The cord used to apply the traction can be easily adjusted to the required length. A special tensioning system ensures that a minimum tension of 15 N (1.5 kg) is always applied to the cord. This enables the cord to be paid out easily, and ensures that it rolls up automatically when released.

The cord must be attached to the carabin hook by means of a knot capable of withstanding considerable traction, without coming undone. A suitable knot is shown in the figure below.

### [20] Patient stop switch

The patient stop switch must be connected to socket [16]. Treatment cannot start unless the patient stop switch is connected!

The switch should be handed to the patient and its use should be explained.

The patient stop switch is an emergency stop. Make sure that this emergency stop can be used by the patient during the whole treatment.

If it is pressed during treatment, the traction force will be reduced to the minimum of 15 N at the maximum change over speed, and the buzzer will sound. The remaining treatment time is still visible on the display

The patient stop switch must be tested every day before use.

The procedure is as follows:

- Connect the cord to the 'frame' on which the ELTRAC is mounted (couch, wall bar bracket or mobile traction frame).
- Adjust the traction force to 500 N
- Adjust a long enough treatment time (e.g. 5 min)
- Start the "treatment", so that the force can go to 500 N.
- Push the patient stop switch.
- The cord must now loosen itself.

### [21] Mains input and mains fuses

The mains cable provided should be connected to this socket. Two fuses (T0.315AH250V) are used (CSA 901 type: *only one fuse*). For replacement of faulty fuses, see paragraph 7.8. Replacing fuses


### [22] Potential equalisation point

For use in areas where local regulations require connection of a potential equalisation cable, an earth (ground) cable and plug can be supplied.

### [23] Type number plate

This plate provides information for identifying the equipment by type, version and serial number (for service, warranty, etc.), together with data such as mains voltage and current consumption.

### [24] Warning sticker cord guidance


Make sure the cord is guided over the pulley as shown in the left drawing.

## 5.2. Introduction to operation

In the following sections, which are devoted to the operation of the ELTRAC 471, particular attention is paid to the control panel (Paragraph 5.1.) and the parameter settings (Paragraph 5.3.). As the operating sequence is dependent on the therapy and accessories employed, only an example is given (Paragraph 5.4.).

### 5.2.1. Parameters

The control panel of the ELTRAC 471 has push buttons for selecting the following parameters:

Parameter indication in intermittent traction

- traction force [F1]
- traction force hold time [t1]
- base force [F2]
- base force hold time [t2]
- changeover speed [v]
- treatment time
- push button for forces exceeding 200 N push button for starting and stopping treatment

### 5.2.2. Displays

Displays are provided for:

- traction force or traction force hold time
- base force or base force hold time
- treatment time
- changeover speed or traction phase.

Traction force and base force are indicated in kilograms\*). Hold times for traction force and base force are indicated in seconds.

The treatment time is in minutes.

Before treatment, the changeover speed is indicated on the bar display, while during treatment this display indicates the current traction force.

### Table of conversion

kg	N	lbs
1,5	15	3,3
2	20	4,4
3	30	6,6
4	40	8,8
5	50	11
6	60	13
7	70	15
8	80	18
9	90	20
10	100	22
20	200	44
30	300	66
40	400	88
50	500	110
60	600	132
70	700	154
80	800	176
90	900	198


1 pound (lb) = 0.454 kgs.

\*) Strictly speaking, the kilogram is a unit of mass, while the unit of force is the Newton (N). However, as the force in traction therapy is still generally expressed in kilograms, the display of the ELTRAC 471 reads in kilograms. To obtain the value in Newtons, simply multiply the displayed figure by 10. For this reason, the relevant push-buttons are also marked 'x 10 N'.

### 5.2.3. Parameter setting

The sequence for parameter setting is as follows:

1. select the parameter to be set by pressing the appropriate push-button,
2. set the selected parameter by turning the central control knob.
3. the parameter setting is shown on the relevant display.

1 parameter setting (selection, 2 setting, 3 reading)

### 5.3. Sequence for setting parameters

The sequence for setting parameters is, in principle, determined by the user. However, there are two main possibilities:

- All parameters are preset before treatment begins. The set values can be adjusted during treatment. See figure a.
- Treatment is started first, and all parameters are set later. See figure b.


figure a


figure b

In *intermittent traction*, the traction force has to be set first, followed by the base force.

#### **Changeover from static to intermittent traction**

To change from static to intermittent traction, additionally the base force parameter must be set and, if required, the hold times (traction force and base force) and the changeover speed.

#### **Changeover from intermittent to static traction**

To change from intermittent to static traction during treatment, it is sufficient to set the base force to the same value as the traction force. The ELTRAC 471 will automatically switch to static traction after 10 seconds.

#### **Termination of treatment**

- a. automatically at the end of the treatment time,
- b. by pressing the start/stop push-button,
- c. by pressing the patient stop switch.

### **5.4. Operating sequence example**

As the operating sequence is dependent on the therapy and accessories employed, only one example is given.

#### **Requirements**

ELTRAC 471 and traction table

Pelvic belt

Thoracic belt

Flexistool

Adjustable traction guide

Fixation clamps (2x)

Fixation straps (2x)

- a. Slide the rolling traction section of the table up against the head section, and lock it in position.
- b. Lay the pelvic belt against the cranial edge of the traction section.
- c. Lay the thoracic belt against the caudal edge of the head section.
- d. Slide the two fixation clamps into the tube underneath the head section.
- e. Fit the fixation straps to the fixation clamps.
- f. Adjust the fixation straps to length so that when they are stretched out, the thoracic belt still lies against the caudal edge of the head section.
- g. Place the patient in position, and fasten the belts.
- h. Attach the cord of the ELTRAC 471 to the pelvic clamps.
- i. Give the patient stop switch to the patient and explain its operation.
- j. Set the treatment time and, if desired, also the other parameters.
- k. Release the rolling traction section.
- l. Start the treatment.
- m. Set the other parameters (see item j.).

## **6.1. Technical maintenance**

The ground connection of this equipment (by means of the mains cord) is important because of the safety provisions. Therefore it is essential that this connection is checked annually. A complete check up is recommended annually. This may be done by your supplier, or by another agency, authorised by the manufacturer. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

Maintenance and all repairs should only be carried out by an authorised agency. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorised persons.

The check up and/or technical maintenance must be carried out conform the procedure described in the service manual of the unit.

Opening of the equipment by unauthorised agencies is not allowed and will terminate any claim to warranty.

## **6.2. Cleaning of the apparatus**

Switch off the unit and disconnect it from the mains supply. The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).

## **6.3. Nylon cord**

The nylon cord should be regularly checked for damage and wear. For replacement please contact your local supplier.

## **6.4. Patient stop switch**

The patient stop switch must be disinfected using a cloth moistened with 70% alcohol.

## **6.5. Environment information**

Your ELTRAC 471 contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart the ELTRAC 471 and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

**Please ensure that you are well informed of the local rules and regulations regarding to the removal of equipment and accessories.**

**7.1 Self-test**

Immediately after switching on, the unit carries out a self-test in which the built-in microcomputer checks a number of vital functions of the equipment. If a fault condition is detected during the self-test or during treatment, this will be signalled by the appearance of dashes (-- --) on the display, and the equipment will cease to function.

**7.2. Displays do not light**

- Check that the mains plug is correctly inserted in the wall socket,
- Check that the connector is correctly inserted in the mains input socket,
- Check the mains cable (use a spare cable),
- Check the mains supply to the wall socket,
- Check the fuses.

If all the above are in order, the unit is probably defective, and should be checked by your local supplier.

**7.3. Dashes (-- --) appear on the display**

If dashes appear repeatedly after switching on the ELTRAC 471, the unit is probably defective, and should be checked by your local supplier.

**7.4. Treatment cannot be started**

Treatment can only be started when the following conditions are met:

- a) a patient stop switch is connected,
- b) the treatment time has been set.

**7.5. Traction force [2] cannot be adjusted to a higher value/ LED flashes**

This situation can be rectified quite simply by pressing the traction force selector [2].

This situation occurs when a parameter other than traction force has been adjusted after the start of treatment, i.e. ten seconds after the adjustment, the ELTRAC 471 automatically switches back to the traction force parameter. Then, for safety reasons, the traction force can only be set to a lower value and not to a higher value. See also remark paragraph 5.1. point [2].

**7.6. The base force [3] cannot be set**

Ensure that the traction force [2] has been set before attempting to set the base force. The base force can only be set to a lower value than the traction force.

**7.7. Early termination of treatment**

The ELTRAC 471 has a built-in, electronic safety device which is activated when large, unintended changes in the traction force of the cord occur. These can be caused by an abrupt movement of the patient, and/or by play in the accessories. We advise taking this into account when fixing the cord (application of the knot), as well as when using new straps or new (extra) cords.

When the safety device is activated, the traction force will be reduced to the minimum of 15 N at the maximum changeover speed, and the buzzer will sound. The unelapsed treatment time continues to be indicated on the display.

**7.8. Replacing fuses**

Defective fuses can be replaced by the user. Remove the mains plug from the wall socket, and remove the mains connector from the input socket. Unscrew the fuse holder. Replace the defective fuse with one of the same type (T0.315AH250V). Screw in the fuse holder until it is finger tight, and reconnect the mains cable.

### 8.1. Product specifications

Types of therapy	:	Static traction, Intermittent traction.
Traction force, minimum	:	15 N continuously adjustable.
Traction force, maximum	:	900 N continuously adjustable.
Treatment time	:	1-60 minutes, in steps of 1 minute.
Intermittent traction:		
Changeover speed	:	Choice of 10 levels.
Hold time, base force	:	0-60 s., 0-10 s in 2-s steps, 10-60 s in 5-s steps.
Safety:		<ul style="list-style-type: none"> <li>• Electronic monitoring of traction force.</li> <li>• Mechanical maximum force limitation.</li> <li>• Patient stop switch.</li> <li>• Selection of forces &gt;200 N.</li> </ul>

### 8.2. Technical data

Mains supply	:	230 V / 50 - 60 Hz.
Max. permissible variation	:	+/- 15%.
Max. current consumption	:	0.22 A (230 V).
Mains fuses	:	2 x T0.315AH250V
Medical device classification	:	IIb (according to Medical Device Directive (93/42/EEC))
Safety class	:	I type B, according to IEC 60601-1.
Earth resistance	:	Typically 0.12 Ohm. (IEC requirement $\leq$ 0.2 Ohm).
Earth leakage current*)	:	Typically 25 $\mu$ A (IEC requirement $\leq$ 500 $\mu$ A).
Earth leakage current under single fault condition	:	Typically 50 $\mu$ A (IEC requirement $\leq$ 1000 $\mu$ A).
Dimensions	:	31 x 27 x 16 cm (w x d x h).
Weight	:	12 kg.

#### Environment conditions for transport and storage

Environment temperature	:	-10° till +70° C
Relative humidity	:	10 till 90 % (not condensing)
Atmospheric pressure	:	500 till 1060 hPa

#### Environments conditions for normal use

Environment temperature	:	10° till 40° C
Relative humidity	:	10 till 90 % (not condensing)
Atmospheric pressure	:	500 till 1060 hPa

#### Classification

##### Medical class IIb

This equipment complies with all essential requirements of the Medical Device Directive (93/42/EEC).

I

##### Safety class I (IEC 60601-1)

The unit is equipped with a safety earth and must be connected to a mains plug with earth contact.

##### Type B

The leakage current and safety earth correspond to standard IEC 60601 -1.

\*) USA AND CANADA ONLY:

Earth leakage current:  $\leq$  100  $\mu$ A

For the ordering data of the ELTRAC 471, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

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