# Radarmed 650+

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo

Copyright:

**Enraf  
Nonius**

Enraf-Nonius B.V.

P.O. Box 810

2600 AV DELFT



The Netherlands

Tel.: + 31(0)15 - 26 98 400

Fax: + 31(0)15 - 25 61 686

Article number : 1435.750-41

March 2000

# Radarmed 650+

Bedieningshandleiding  
Operating Instructions  
Gebrauchsanweisung  
Mode d'emploi  
Instrucciones de manejo

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## TABLE OF CONTENTS

<b>1. Introduction</b> .....	<b>17</b>
1.1. General .....	17
1.2. Therapy possibilities .....	17
1.3. Radiators .....	17
1.4. Finally .....	17
<b>2. Preliminary notes</b> .....	<b>17</b>
<b>3. Safety</b> .....	<b>17</b>
3.1. Safety .....	17
3.2. Product Liability .....	17
<b>4. Indications and contra-indications</b> .....	<b>18</b>
4.1. Indications .....	18
4.2. Absolute contra-indications .....	18
4.3. Relative contra-indications .....	18
4.4. Special groups of patients .....	19
4.5. Note .....	19
<b>5. Rules for microwave treatment</b> .....	<b>19</b>
5.1. General .....	19
5.2. Preliminary measures of the unit and personnel .....	19
5.3. Preparation of the patient .....	20
5.4. Dosage .....	20
<b>6. Controls</b> .....	<b>21</b>
6.1. Control panel .....	21
6.2. Apparatus .....	21
<b>7. Installation</b> .....	<b>22</b>
7.1. Installation of the unit .....	22
7.2. Connection .....	22
7.3. Moving and transport .....	22
7.4. Electromagnetic interference .....	22
7.5. Connection mains cable .....	22
7.6. Connection radiator arm .....	22
7.7. Connection radiator .....	23
7.8. Connection HF-cable .....	23
7.9. Switching on and self test .....	23
<b>8. Operation</b> .....	<b>24</b>
8.1. Set up of the Radarmed 650+ .....	24
8.2. Arrangements regarding treatments .....	25
<b>9. Maintenance</b> .....	<b>25</b>
9.1. Technical maintenance .....	25
9.2. Cleaning of the apparatus .....	25
9.3. Cleaning of the radiators and HF-cable .....	25
9.4. End life span apparatus and accessories .....	25
<b>10. Fault conditions</b> .....	<b>26</b>
10.1. Display's fail to light up .....	26
10.2. Error code Er XXX* on the display .....	26
10.3. Power cannot be set .....	26
10.4. During treatment the HF-power switches off and a sound signal is audible ( 3 beeps of 1 sec.) .....	26
<b>11. Specifications</b> .....	<b>26</b>
11.1. Microwave therapy .....	26
<b>12. Technical specifications</b> .....	<b>27</b>
12.1. Unit .....	27
12.2. Environment conditions for transport and storage .....	27
12.3. Environments conditions for normal use .....	27
12.4. Classification .....	27
<b>13. Ordering data</b> .....	<b>27</b>

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---

### 1.1. General

Dit apparaat is bedoeld voor gebruik door vakkundig personeel in de fysiotherapie, revalidatie en aangrenzende gebieden. The Radarmed 650+ is a unit for continuous and pulsed microwave therapy with a frequency of 2450 MHz (12 cm wavelength). Safe dosage control based on the patient's sensitivity to heat is a big advantage of microwave therapy. Microwave therapy can easily be applied thanks to the unipolar application. Three different radiators are sufficient for all applications occurring in daily practice.

### 1.2. Therapy possibilities

The therapeutic effect of warmth has been put to use ever since the earliest beginning of medical science. The Radarmed 650+ is extremely suitable for localized application of warmth to deep located tissues. The Radarmed 650+ operates at a frequency of 2450 MHz, corresponding to a wavelength of 12 cm. At this wavelength, the high-frequency energy is absorbed particularly well in every kind of tissues containing a high percentage of water.

The possibilities of microwave therapy is significantly widened by the use of the pulsed mode. In this mode, microwave pulses with a fixed, high intensity are generated. In contrast with continuous microwave the higher intensity of the pulses briefly warms up the deeper lying tissues. At the same time, the same dosage rules can be used as for treatment with continuous microwave, since the sensation of warmth felt by the patient is the same with either form of treatment. Pulsed microwave can be used in cases where only a low dosage is permissible, but nevertheless deeper lying regions have to be reached.

### 1.3. Radiators

#### *Local field radiator*

The local field radiator enables greater depth effects to be obtained over smaller areas. It is the most universal and therefore most used radiator.

#### *Longitudinal radiator*

The longitudinal radiator is designed for the treatment of the extremities.

#### *Large field radiator*

The large field radiator is used for treating large areas of the body. The form of the radiator is adapted to the contours of the body. This radiator is very suitable for treating isolated areas such as the shoulder or knee joint, as well as larger muscle groups in the back and shoulder regions.

### 1.4. Finally

You have made a wise choice in selecting the Radarmed 650+. We are confident that your unit will continue to give satisfaction over many years of use. Nevertheless, if you have any queries or suggestions, please contact your authorised distributor.

---

## 2. PRELIMINARY NOTES

---

It is important that you read these operating instructions carefully before using the Radarmed 650+. Please make sure that these instructions are available to all personnel who operates the equipment.

Pay attention to the following before using the Radarmed 650+:

1. Keep yourself informed of the contra-indications (see chapter 4).
2. The apparatus may not be used in close proximity (i.e. less than 2 metres) to shortwave equipment.
3. The apparatus may not be used in so-called "wet rooms" (hydrotherapy rooms).

The manufacturer cannot be held responsible for the results of using this apparatus for any purposes other than described in these operating instructions.

---

## 3. SAFETY

---

### 3.1. Safety

The safety of the unit is guaranteed by a built in microprocessor which monitors all important functions of the unit, every 50 seconds. If during treatment the HF-cable is mistakenly disconnected (from apparatus or radiator), the unit automatically switches of the HF-power and 3 beeps of 1 second is audible. Continuation of the treatment is possible by setting a new intensity after the connection is restored.

### 3.2. Product Liability

A law on Product Liability has become effective in many countries. This Product Liability law implies, amongst other things, that once a period of 10 years has elapsed after a product has been brought into circulation, the manufacturer can no longer be held responsible for possible shortages of the product.

**4.1. Indications***Dermatology*

- Abscess
- Alopecia
- Carbuncle
- Frostbite, local
- Furuncle
- Gonorrhoeal arthritis
- Herpes zoster
- Paronychia (initial stage)
- Post-herpetic neuralgia

*Gynaecology*

- Amenorrhoea
- Dysmenorrhoea
- "Low back pain"
- Mamma abscess
- Mastitis
- Menstrual irregularity

*Internal Medicine*

- Bronchial asthma
- Bronchiectasis
- Bronchitis chronica
- Cholecystitis
- Circulation disorders of the acra
- Pleuritis chronica sicca
- Raynaud's disease

*Neurology*

- Intercostal neuralgia
- Ischialgia
- Neuralgia in general
- Neuritis
- Occipital neuralgia
- Root pain (radiculitis)

*Surgery and Orthopaedics*

- Arthritis
- Arthropathy
- Arthrosis deformans
- Articular rheumatism
- Bechterew's disease (symptomatic only)
- Brachialgia
- Bursitis
- Calf cramp
- Contusions
- Coxarthrosis
- Coxitis
- Distortions
- Epicondylitis
- Haematoma
- Ischialgia
- Lumbago
- Muscle cramps
- Meniscopathy
- Myalgia
- Periostitis
- Polyarthritis rheumatica
- Spondylosis
- Tendinitis
- Torticollis

**4.2. Absolute contra-indications**

- Patients with implanted heart pacemakers.
- Metallic endoprostheses in the treatment region.
- Malignant tumours. An increase in temperature may result in an increase in the rate of growth of the tumour or metastatic secondaries.
- Active tuberculosis or other active specific disorder.

**4.3. Relative contra-indications**

- Acute infectious disorders, especially when accompanied by high fever.
- Acute inflammatory processes of the bones and joints. Warmth may, in certain circumstances, induce or accelerate necrosis of the tissue. Furthermore, local heating may be more severe than anticipated in areas where the blood supply is reduced or inhibited as a result of lacking normal heat dissipation by blood circulation. Particular care is necessary when deciding the dosage if diathermy is used in the presence of ischemia or severe tissue damage.
- Traumata and vascular lesions for the first 24 to 36 hours after the incident, in order to prevent the danger of internal and external bleeding.
- Disorders in sensitivity to heat.
- Anaesthetized areas.
- In the region of the abdomen during menstruation.
- Patient with intra-uterine pessaries.
- Osteoporosis.
- Increased tendency towards bleeding and use of anticoagulants.
- Thrombophilia.
- The epiphyseal disks must also be exempted with children until the cessation of growth.
- The application of diathermy in the region of the eyes. When treating the head e.g. nose or frontal sinuses, it is recommended to take into account the reduced vascularization of the eyes and provide a pair of protective goggles in order to avoid heat stress. Microwave treatment of the eyes should be carried out by an eye specialist.
- Wearers of contact lenses. Contact lenses may inhibit heat dissipation and thus cause overheating of the eyes. They should be removed before treatment in the head region.
- Treatment in the region of the testicles. Special protective measures must be taken here since the testicles are extremely sensitive to heat, which may cause permanent damage or even sterility.
- Pregnancy. Special care must be taken to protect the embryo or foetus during microwave treatment in the region of the abdomen.

## 5. RULES FOR MICROWAVE TREATMENT

### 4.4. Special groups of patients

- Infants are best undressed completely before treatment. Their physical volume calls for special care in calculating the dosage and repeated checking of the skin temperature by hand.
- Never leave infants or sensitive adults, especially elderly patients, unattended during treatment.

### 4.5. Note

This list makes no claim to completeness. For the application of microwave therapy we refer to the booklet "Treatment tables" with reference number 1435.770. In this you will find information on the radiators, their applications, the recommended dosage and the number of treatments. The statements are related to the in paragraph 8.1 mentioned distances of the radiator to the skin and in respect of contra-indications and other treatment criteria from case to case.

### 5.1. General

Please note the following before a treatment is started.

- Patients with heart pacemakers or electrodes are contra-indicated and must be kept away from all rooms where microwave units are operated.
- Persons who are not being treated must not remain within a distance of 1,5 m of the radiator when the unit is being operated.

### 5.2. Preliminary measures of the unit and personnel

- During treatment a high-frequency field is present in the radiator's direction of radiation which may heat up metallic objects and cause malfunctions in electronic equipment. Persons within this zone should therefore remove metal objects (such as paging units, hearing aids, watches, jewellery and piercings).
- The radiator must only be aligned when there is no power output. In all cases, a dosage must be set only when there is a patient in front of the emitter. In this way, the patient properly absorbs the major part of the microwave energy, thus reducing the degree or interference to a minimum.
- Only original accessories may be used. Treatment with damaged radiators or cables (e.g. external deformation) is not permissible.
- Handle the radiators with care. Negligent treatment can influence the electromagnetic field of the radiator.
- Treatment chairs, tables or beds within the zone of radiation should not be made of metal or contain any metal parts or semi-conducting covering materials, as these may heat up. Chairs or couches made of wood with textile or plastic coverings are recommended.

### 5.3. Preparation of the patient

- Before the start of treatment, the patient should remove all articles of clothing containing metal (lurex) and metal items in pockets or any other such articles which are liable to come within the radiation zone. Such articles include rings, chains, clasps, bracelets, jewellery, watches, piercings, articles with metal fittings such as bras, suspender belts, corsets, hearing aids and metallic and electronic prostheses, etc.
- Moisture concentrations heat up more intensively in the treatment zone of the radiator. In leather or synthetic clothing occurs easily increased perspiration. It is therefore recommended to unclasp the body parts which are to be treated and, in the case of strongly perspiring, dried off (folds in the skin). Dry bandages can be irradiated quite safely. Diathermy should not be applied in regions covered by moist bandages or adhesive tape.
- Parts of the body containing metallic endoprotheses such as medullary pegs, metal splinters, metallic joints, pins, spirals, skull plates, dental fillings and crowns, wires and so on, should not be subjected to treatments or only at extremely low dosage levels.
- When treating narrow parts of the body (e.g. a wrist) it should be ensured that the radiator is set up so that sensitive parts of the body (e.g. eyes, testicles) are not in the irradiated area. In these cases it is recommended to use a smaller radiator.
- It should be made possible for the patient to adopt a relaxed and calm posture during the entire treatment.

### 5.4. Dosage

As long as the rules of treatment are adherent to, no harm will be incurred by patients as the result of microwave radiation. As the blood circulation increases strongly in the region the generated warmth is quickly dissipated and there is no built-up of heat. Injuries can be sustained only as the result of gross errors of treatment (heavy overdose), exactly as with any other form of physical or medical therapy.

If goes without saying that markedly aqueous tissue which is poorly vascularized, such as the eyes and testicles, should be treated only at very low dosage levels.

The decisive factor for all dosage calculations must be the patient's own sensitivity to heat. During the first treatment of a series, the patient should be asked immediately about his subjective feeling.

- Always set the dosage according to the individual patient and never in accordance with a fixed scheme.
- Every patient has a different susceptibility to warmth and this generally reaches a steady state only after about 5 minutes.
- The patient's sensitivity to warmth may vary during the course of a treatment (adaptation).
- The patient must be asked to report an increase of the warmth sensation. The dosage must then be reduced.
- The dosage must not be increased under any circumstances as a result of any lessening of the feeling of warmth during treatment.
- When sensibility disorders is suspected the subjective warmth perception of the patient is checked by means of the skin temperature. In order to avoid overheating in such cases, the dosage must be adjusted with extreme care. Apply only small doses for short periods to the treatment zone.

**Microwave therapy should never be applied to parts of the body under local anaesthetic because this might lead to local burns.**

For further information with regard to dosage we refer to the booklet "Treatment tables" with reference number 1435.770.

## 6.1. Control panel

(See jacket flap)

### [1] Indication LED therapy mode

The illuminated LED indicates the selected mode; continuous or pulsed.

### [2] Selector therapy mode: continuous or pulsed microwave

The unit automatically switches to continuous mode when switched on or at the end of treatment. It is possible to switch from continuous to pulsed mode and visa versa during treatment.

### [3] Display treatment time

The display shows the set c.q. remaining treatment time in minutes. When the time is expired a zero appears on the display, a sound signal is audible and the HF-power switches off. A flashing dot in the display means the clock is running and treatment started.

### [4] Treatment time setting selectors

These selectors are used to set the treatment time.

### [5] Power display

The display shows the set power. A flashing zero on the display indicates that no treatment time has been set. In this case it is not possible to set the power.

### [6] Power setting selectors

These selectors are used to set the power. It is only possible to set the power if also a treatment time has been set.

### [7] Stop button

Button to end a treatment premature. After operation of this button the treatment time runs down to zero, a sound signal is audible and the HF-power is shut off.

## 6.2. Apparatus

### [8] On/off switch

With this switch the Radarmed is turned on/off.

### [9] Ventilation slots

These openings must be kept clear for ventilation.

### [10] Connection points electrode arm

The electrode arm can be connected at choice at the left or right side of the unit (see paragraph 7.6 of this manual).

### [11] Hook for mains cable

Hook to hang up the mains cable when the unit is out of operation.

### [12] HF-connection

Connection for the HF-cable (see paragraph 7.8 of this manual).

### [13] Connection mains cable

Connect the Radarmed to a grounded safety socket of a voltage corresponding to the data mentioned on the type plate of the unit.

### [14] Connection for potential equalization cable

A special cable can be supplied for areas where local regulations require connection of a potential equalization cable.

### [15] Type plate

The type plate gives details of the model, series- and serial number (for service, warranty etc.) as well as ratings such as mains voltage and power consumption.

### 7.1. Installation of the unit

- Do not install the unit in a location near to a heat source such as a radiator.
- Avoid exposure to direct sunlight, rain, excessive dust, dampness, mechanical vibrations and shocks.
- Electronic equipment such as cordless telephones should not be used in close proximity of the unit.
- This unit may not be used in so-called "wet-rooms" (hydrotherapy rooms).
- The unit has to be installed in such way that liquid cannot enter.
- Always use the by Enraf-Nonius for this unit prescribed original accessories.
- Do not cover the ventilation slots.

Should any liquid ingress into the housing of this equipment, unplug the unit from the wall socket (if connected) and have it checked by an authorised expert.

### 7.2. Connection

- Mains supply connections must comply with the national requirements regarding medical rooms.
- Prior to connection of this apparatus to the mains supply, check that the voltage and frequency stated on the type plate correspond with the available mains supply.
- For safety reasons the brakes on the wheels of the unit should be locked before the beginning of a treatment.

### 7.3. Moving and transport

When moving or transporting of the Radarmed 650+, never pull at the radiator or radiator arm, to avoid the risk of tilting of the unit.

### 7.4. Electromagnetic interference

- Operation in close proximity (i.e. less than 2 metres) to working shortwave units may produce instability in the Radarmed output.
- In order to avoid interference, shortwave units should not be allowed to radiate in the direction of the microwave units, either in the same or in adjacent rooms.
- To prevent electromagnetic interference, we strongly advise using separate mains groups (phases) for the Radarmed and the shortwave equipment. Ensure that the mains cable of the shortwave unit does not come near the Radarmed or the patient.
- Electroacoustic equipment (e.g. radios, intercom systems, telephones systems, etc.) should be screened against microwave irradiation if interference is experienced. It may be sufficient to change the radiators direction of radiation.


If problems with electromagnetic interference persist, please contact your authorised supplier.

### 7.5. Connection mains cable

- Connect the supplied mains cable to the connection mains cable [13].
- Connect the mains cable to a grounded safety socket

### 7.6. Connection radiator arm

- Insert screw [1] and partly screw in.
- Place the arm against the side of the cabinet (connection point [10] left or right) so that the notch fits over screw [1].
- Tighten screw [2].
- Tighten both screws with the hexagon key supplied.

## 7.7. Connection radiator

To attach or remove the radiator, push the black locking knob (on top the radiator arm). When fitting the radiator, the locking knob must click into place.

## 7.8. Connection HF-cable

### 7.8.1. Connection to the unit


The HF-cable is connected through connection point [12] to the unit by means of a screw thread connection.

The length of the HF-cable is critical. Therefore only use the by Enraf-Nonius supplied cable.


### 7.8.2. Connection to the radiator

The HF-cable is connected to the radiator by means of a "fast" coupling.

#### To plug on the connector


Push on connector


Check to ensure that connector is secure

#### To unplug the connector

Pull back ring

Pull out connector

## 7.9. Switching on and self test

- Switch on the unit using the on/off switch [8].
- Immediately after switching on, the unit carries out a self test

Check that:

- all LED's light up briefly;
- at the end of the test a beep can be heard.

Contact your supplier if this is not the case.

## 8.1. Set up of the Radarmed 650+

**Switching on**

- Connect the radiator.
- Connect the HF-cable to the radiator and Radarmed.
- Switch on the unit using the On/off switch.

Immediately after switching on, the unit carries out a self test.

**Choice of therapy mode**

- Select continuous or pulsed mode

The illuminated LED indicates the selected mode.

**Alignment of the radiator**

- Align the radiator to the area which must be treated. Align the radiator carefully so that other parts of the body are minimal in the direction of radiation.
- Inform the patient.

For the alignment of the radiators take notice of the following distances of the radiator to the body surface.

- Local field radiator approx. 10 cm.
- Longitudinal radiator approx. 5 cm.
- Large field radiator approx. 1 cm.

**Treatment time**

- Select the treatment time.

The set time is indicated on the display.

**Power adjustment**

- Set the power.

The set power is indicated on the display. At the end of the treatment time a sound signal is audible and the HF-power switches off.

- The treatment time runs down as soon as a power is set. A flashing dot in the display indicates that the clock is running.

Treatment starts immediately after operation of the power setting selectors.

**Therapy adjustment**

- During treatment it is possible to switch from continuous to pulsed mode and visa versa.
- During treatment It is possible to increase and decrease the set power.

Every patient has a different susceptibility to warmth and this generally reaches a steady state only after about 5 minutes. In between adjustment of the power is therefore necessary.

**Preliminary ending of the treatment**

Treatment can be ended by operation of the Stop button. Treatment stops automatically when the treatment time has elapsed.

## 8.2. Arrangements regarding treatments

Also see chapter 5 “Rules for microwave treatment” of this manual.

### Before treatment

- Check the patient on the presence of possibly absolute and relative contra-indications.
- Articles of clothing containing metal or clothing with metal zippers should be taken off (empty pockets). Jewellery, watches and metallic prostheses should also be taken off.
- Test the sensitivity to warmth of the area to be treated.

### During treatment

- The patient must be in sight of the therapist at all times.
- The patient is regularly asked to report his/hers findings. As the patient to report an increase of the warmth. If necessary the treatment can be adjusted.

### After treatment

- The expected effects are checked.
- The patient is asked to comment subsequently on any reactions that may occur.

## 9.1. Technical maintenance

We recommend to have the unit checked annually. This may be done by your supplier, or by another agency, authorised by Enraf-Nonius. It is also recommended that a record of the service history is kept for all activities relating to service and maintenance. In some countries this is even obligatory.

Maintenance and all repairs should only be carried out by an authorised agency. Enraf-Nonius will not be held responsible for the results of maintenance or repairs by unauthorised persons.

The check up and/or technical maintenance must be carried out conform the procedure described in the service manual of the unit.

Opening of the equipment by unauthorised agencies is not allowed and will terminate any claim to warranty.

## 9.2. Cleaning of the apparatus

Switch off the unit and disconnect it from the mains supply.

The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution).

## 9.3. Cleaning of the radiators and HF-cable

The radiators can be disinfected using a cloth moistened with 70% alcohol. The radiators and HF-cable may not be sterilized by expose to steam or gas. Clean the radiators after each use.

The radiators and HF-cable should be regularly inspected for damage and bad contacts. Damaged radiators and a damaged HF-cable may not be used.

## 9.4. End life span apparatus and accessories

Your Radarmed and accessories contains materials which can be recycled and/or are noxious for the environment. At the end of the duration of life, specialised concerns can take apart these articles and sort out the noxious materials and materials for recycling. By doing so you contribute to a better environment.

**Please ensure that you are well informed of the local rules and regulations regarding to the removal of apparatus and accessories.**

---

## 10. FAULT CONDITIONS

---

### 10.1 Display's fail to light up

- Check if the mains cable is connected to mains supply.
- Check if the mains supply is available (power supply) for instance with a lamp.
  - If the lamp does not light up, check the mains fuses.
  - if the lamp lights up, check the fuses of the unit.

### 10.2. Error code Er XXX\* on the display

The unit has discovered a fault using its "self-test". Switch the unit back on. If the code re-appears, contact your supplier. The unit is probably defective.

\* Er XXX: a number is mentioned on position X

### 10.3 Power cannot be set

- HF-cable is not connected.
- No radiator is connected.
- No treatment time is set.

### 10.4 During treatment the HF-power switches off and a sound signal is audible ( 3 beeps of 1 sec.)

- The HF-cable is mistakenly disconnected from radiator or apparatus.
- The HF-cable is damaged and must be replaced.
- Bad contact occurred in the connection of the HF-cable to the apparatus and/or radiator.

---

## 11. SPECIFICATIONS

---

### 11.1. Microwave therapy

Frequency	: 2450 MHz $\pm$ 50 MHz
HF power, continuous	: 0 - 250 Watt at 50 Ohm
HF power, pulsed	: 0 - 250 Watt with a fixed peak power of 1500 Watt +/- 30%
Operating modes	: continuous and pulsed
Clock	: 0 - 30 minutes

**12.1. Unit**

Mains voltage	: 220-240 Volt
Frequency	: 50 Hz
Current consumption	: 1000 VA
Current rating	: 4,5 Ampere
Safety class	: 1 type B (according to IEC 601-1)
Medical device classification	: IIb (according to Medical Device Directive (93/42/EEC))
Weight	: 45 kg (including arm)
Dimensions	: 93x39x49 cm (hwxwd including castors)
Safety tests	: IEC 601-1, IEC 601-2-6, IEC 601-1-2, IEC 601-1-4 (TÜV Rheinland)
Marking	: 

**12.2. Environment conditions for transport and storage**

Environment temperature	: -10° till +50° C
Relative humidity	: 10 till 95 %
Atmospheric pressure	: 500 till 1060 hPa

**12.3. Environments conditions for normal use**

Environment temperature	: 10° till 40° C
Relative humidity	: 10 till 90 %
Atmospheric pressure	: 500 till 1060 hPa

**12.4. Classification****Medical class IIb**

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

**International safety standard IEC 601-1**

This equipment meets the requirements of the International Standard for the safety of medical electrical equipment IEC 601-1 +A1, +A2 (general standard), IEC 601-2-6 (standard for microwave equipment), IEC 601-1-2 (standard for EMC) and IEC 601-1-4 (standard for programmable systems).

**Non ionising radiation**

**Safety class 1**

The unit is equipped with a safety earth and must be connected to a mains plug with earth contact.

**Type B**

The leakage current and safety earth correspond to standard IEC 601-1 and IEC 601-2-6.

*Technical modifications reserved.*

**13. ORDERING DATA**

For the ordering data of the Radarmed 650+, standard accessories and additional accessories we refer to the Catalogue for Physiotherapy.

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