# Service Manual

High Flow CO<sub>2</sub> – Pneu  
2231

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## General Instructions

### ☞ **IMPORTANT!**

*Richard Wolf GmbH as the manufacturer and the seller of this product does not accept any liability for direct or consequential damage or injury, caused by improper use or operation of these instruments, in particular due to non-compliance with this instruction manual or inadequate re-processing or maintenance.*

*The fact that we supply technical documentation does not mean that we authorize the adjustment, repair or modification of the products.*

*This documentation may not be copied nor may its contents be passed on and commercially utilized, unless expressly agreed otherwise.*

### ☞ **NOTE!**

*The manufacturer is only responsible for the safety features of the products within the scope of the legal provisions provided that repairs and changes are carried out by himself or a third party explicitly authorized by him.*

*If non-authorized persons perform repairs or modifications of the products, the manufacturer does not accept any liability and the warranty becomes void. Components influencing safety may only be replaced with genuine parts.*

*After any repair or modification, the technical safety of the product must meet the applicable device standard.*

*Each repair, modification or test must be recorded in a report.*

*Proper care and maintenance are required to ensure safe operation of medical products.*

*We therefore recommend that the products be thoroughly checked for proper functioning and completeness upon receipt and prior to each use to prevent injury to the patient or user.*

*New and repaired products must be reprocessed in accordance with the instruction manual prior to their first use.*

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# ENGLISH

# Contents

<b>1</b>	<b>Maintenance</b> .....	<b>1</b>
1.1	Important Instructions .....	1
1.2	Maintenance of device and accessories .....	1
1.3	Preventive maintenance .....	1
1.4	Corrective maintenance .....	1
1.5	Visual Check .....	2
1.6	Electrical safety test .....	2
1.7	Function test .....	2
<b>2</b>	<b>Technical Description</b> .....	<b>3</b>
2.1	General .....	3
2.2	Program function .....	3
2.2.1	Main program routine .....	3
2.2.2	Read-in of device settings .....	3
2.2.3	Adjusting the volume .....	3
2.2.4	Checking the system tightness for smoke gas recirculation .....	3
2.2.5	Text output on video monitor .....	4
2.2.6	Resetting the CO <sub>2</sub> consumption indicator .....	4
2.2.7	Control of gas recirculation motor .....	4
2.2.8	Checking the remaining CO <sub>2</sub> supply .....	4
2.2.9	Starting / stopping insufflation .....	4
2.2.10	Preselecting the nominal pressure .....	5
2.2.11	Measuring the CO <sub>2</sub> flow .....	5
2.2.12	Calculation and display of flow value .....	5
2.2.13	Measuring the actual pressure .....	5
2.2.14	Calculation and display of actual and nominal pressures .....	5
2.2.15	Integration of the CO <sub>2</sub> consumption .....	5
2.2.16	Generating the insufflation cycle .....	6
2.2.17	Activation of valves, pushbuttons, LEDs, alarms and indicators .....	6
2.2.18	Valve 1 .....	6
2.2.19	Valve 2 .....	6
2.2.20	Valve 3 .....	6
2.2.21	Valve 4 .....	7
2.2.22	Valve 5 .....	7
2.2.23	Valve 6 .....	7
2.2.24	Valve 7 .....	7
2.3	Typical layout of a control using an I <sup>2</sup> C bus .....	8
<b>3</b>	<b>Setting Instructions (64 640.121–R01)</b> .....	<b>9</b>
3.1	Power supply PCB (64 151.716) .....	9
3.1.1	Test and measuring equipment .....	9
3.1.2	Fuses .....	9
3.1.3	Supply voltages on the power supply PCB .....	9
3.2	Mother board (64 151.715) .....	9
3.2.1	Test and measuring equipment .....	9
3.2.2	Function check .....	9
3.2.3	Adjusting the reference voltage .....	9
3.2.4	Adjusting the headline text, time and date .....	9
3.2.5	Function check of valves and adjustable pressure reducer .....	10
3.2.6	Alarm function .....	10

## Contents

3.2.7	Control of smoke gas pump .....	11
3.2.8	Volume (loudness) adjustments .....	11
3.3	PCB for video data insertion (64 151.746) .....	11
3.3.1	Test and measuring equipment .....	11
3.3.2	Preparation .....	11
3.3.3	Test setup .....	12
3.3.4	Adjustment of video signals on oscilloscope .....	12
3.3.5	Oscilloscope image .....	12
3.4	Adjustment of pneumatic section .....	13
3.4.1	Test and measuring equipment .....	13
3.4.2	Adjustment of pressure reducer (for location of pressure reducers see location chart 64 220.022) .....	13
3.4.3	Adjustment of servomotor (only after replacment of servomotor!) .....	13
3.4.4	Adjustment of flow .....	13
3.4.5	Adjustment of pressure switch (up to SNo 200) .....	14
3.4.6	Adjustment of pressure switch (as of SNo 201) .....	14
<b>4</b>	<b>Reparaturteile / Repair parts .....</b>	<b>15</b>
<b>5</b>	<b>Anhang / Annex .....</b>	<b>16</b>
5.1	Lageplan / Component drawing .....	16
5.2	Pneumatikpläne / Pneumatic diagrams .....	16
5.3	Verdrahtungspläne / Wiring diagrams .....	16
5.4	Stromlaufpläne / Circuit diagrams .....	16
5.5	Bestückungsdrucke / Component diagrams .....	16
<b>6</b>	<b>Protokolle / Reports .....</b>	<b>17</b>
6.1	Prüfprotokoll / Test Report .....	17
6.2	Wartungsprotokoll / Maintenance report .....	19

# 1 Maintenance

## 1.1 Important Instructions

This service manual describes the service measures specified for this product.

**IMPORTANT!**

*In order to carry out the service measures it is absolutely necessary to know the instruction manual of the product.*

## 1.2 Maintenance of device and accessories

Maintenance is subdivided into preventive and corrective maintenance:

- ◇ Preventive maintenance: Regular checking of device and accessories.
- ◇ Corrective maintenance: Modifications and repairs.

## 1.3 Preventive maintenance

To protect the person testing, the maintenance jobs and the checking of the device must be carried out in the described order.

- ◇ Visual check
- ◇ Electrical safety check
- ◇ Functional check

The required electrical test and measuring procedures are described in the EN 60601 / IEC 601 and DIN 57751 / VDE 751 standards.

The values indicated in the section entitled "safety check" relate to testing in accordance with EN 60601 / IEC 601.

## 1.4 Corrective maintenance

Corrective maintenance includes repairs and modifications to device and accessories.

After repairs or modifications, check all measurement values (basic values) of the device by means of the setting instructions and readjust the values if there are any deviations.

**NOTE!**

*Document all maintenance and checking work as well as modifications to the device or the accessories.*

## 1.5 Visual Check

Designation	Check for
Device and accessories	<ul style="list-style-type: none"> <li>◆ Contamination endangering safety and general cleanliness</li> <li>◆ Mechanical damage</li> <li>◆ Loose or missing parts</li> </ul>
Controls	◆ Mechanical function and easy operation
Marking / symbols	◆ Completeness and legibility
Safety-relevant labelling (such as warnings)	◆ Completeness and legibility
Fuses	◆ For value indicated on the identification plate by the manufacturer (nominal current and melting characteristics)
Lamp label (light sources only)	◆ intact and legible

## 1.6 Electrical safety test

Designation	Check for
Protective earth connection	<p><b>Measured values:</b></p> <ul style="list-style-type: none"> <li>◆ without supply line connected <math>\leq 0.1 \text{ Ohm}</math></li> <li>◆ with supply line connected <math>\leq 0.2 \text{ Ohm}</math></li> </ul> <p><b>Test conditions:</b>  <math>10 \text{ A} \leq I_{\text{meas}} \leq 25 \text{ A}</math> , <math>V_{0\text{meas}} \leq 6 \text{ V}</math> , <math>t_{\text{test}} \leq 5 \text{ s}</math></p> <p><b>To be checked:</b>                      The resistance between protective earth, protective earth connection and the conductive parts connected to them.</p>
Insulation resistance	<p><b>Measured values:</b></p> <ul style="list-style-type: none"> <li>◆ Insulation resistance <math>\geq 2 \text{ MOhm}</math></li> </ul> <p><b>Test conditions:</b>                      Test voltage 500V DC</p> <p><b>To be checked:</b>                      All mains / power lines relative to earth</p>
Leakage current	<p><b>Measured values:</b></p> <ul style="list-style-type: none"> <li>◆ Device leakage current <math>\leq 40 \text{ }\mu\text{A}</math></li> <li>◆ Patient leakage current <math>\leq 10 \text{ }\mu\text{A}</math></li> </ul> <p><b>Test conditions:</b>                      MD as per IEC 601</p> <p><b>To be checked:</b>                      The leakage current which can flow to earth from each pole of the power supply via the insulation through the protective earth or, via the patient, from the component in contact with the patient's body.</p>

## 1.7 Function test

See Instruction Manual section 3.5

## 2 Technical Description

### 2.1 General

The functions of the device are microprocessor-controlled and monitored. The High Flow CO<sub>2</sub> Pneu 2231 automatic insufflator software is written completely in assembler language.

The software problem definition is subdivided into tasks which are dealt with in subroutines. The relationships between the individual subroutines are as simple as possible. The various subroutines only have access to the respective memory cells which are enabled for the corresponding subroutine.

### 2.2 Program function

After reset the internal RAM receives the default value  $\phi H$ . All interrupts which are not used receive default values and are inhibited.

The serial RS 232 interface and the I<sup>2</sup>C bus interface in the device are initialized and enabled.

Effective parameters of the device such as nominal pressure, alarm signal volume and text mask for video data insertion are read from the non-volatile memory (EEPROM).

The video text controller is initialized and the text mask displayed on the monitor.

A self-test of the device completes the initialization phase: In this self-test the flow and pressure sensors are checked and a sound sequence checks the function of the acoustic alarm.

#### 2.2.1 Main program routine

The main program is structured in such a way that with every complete program run all subroutines contained in the repetition routine are called up.

Depending on the device status, the corresponding task is processed or the subroutine is exited immediately.

#### 2.2.2 Read-in of device settings

All LED status indicators and keys (pushbuttons) are read in via the I<sup>2</sup>C bus and the status of the LEDs (on/off) and keys (pushbuttons) are (actuated/not actuated) is stored in memory cells.

#### 2.2.3 Adjusting the volume

The volume (loudness) of the alarm sounds can be adjusted in 16 steps. For this, actuate the "alarm reset" key and in addition, press one of the "< / >" keys to increase or decrease the volume. With each pressing of the key the volume is changed by one step. The value adjusted in this way is stored in an EEPROM and via an I/O extension and different resistors it is connected to an amplifier.

#### 2.2.4 Checking the system tightness for smoke gas recirculation

To guarantee that no outside air can get into the sterile filter system for gas recirculation it is necessary to check the hose system for leaks.

The leakage test is started when insufflation is switched off and the user actuates the "pump off" key for more than 2 seconds. The hoses of the system must be connected in accordance with the illustration in the instruction manual.

Valve 5 and valve 6 are opened and a sound signal marks the start of the test. The 15L flow sensor measures the flow after 3 seconds; the flow is displayed on the digital display. During the leakage test the nominal pressure indicator is not updated. If the flow measured is  $< 0.2$  l/min the pump is enabled and a sound sequence signals the end of the test. If the flow measured is  $> 0.2$  l/min the pump cannot be started and the test can be repeated.

### 2.2.5 Text output on video monitor

In this subroutine the entire video data insertion is controlled and monitored via a video processor. When a key on the connected PC keyboard is actuated, the corresponding keyboard code is read in via interrupt and a validity check is carried out.

When operating function key F7 a screen mask appears allowing the input of the headline text of 27 characters. Pressing RETURN transfers the headline text to the video processor.

Function keys F8 and F9 generate screen masks for setting date and time. These masks serve to input date and time in a fixed format via the keyboard. Actuating the RETURN key loads the data via I<sup>2</sup>C bus into the software clock.

The device-specific data (pressure, flow, consumption) as well as date and time are updated online on the screen.

The text input in the text field is lost either after pressing F10 or when the device is switched off.

### 2.2.6 Resetting the CO<sub>2</sub> consumption indicator

Actuating the "consumption reset" key resets the 3-digit digital display for CO<sub>2</sub> consumption to zero. The gas consumption previously indicated is lost.

### 2.2.7 Control of gas recirculation motor

The pulse width output integrated in the microprocessor 80C552 is used to control the speed of the gas recirculation motor. Depending on the pushbutton depressed (I, II, III), the speed is adjusted in three stages. The speed stage selected is indicated by illumination of the corresponding pushbutton. The "pump off" pushbutton switches off the motor.

### 2.2.8 Checking the remaining CO<sub>2</sub> supply

The filling level of the CO<sub>2</sub> cylinder connected to the device is checked. This is effected via two pressure switches with different switching thresholds (30 bar and 5 bar); as a consequence the following messages are possible:

- Cylinder is full: green LED on, yellow LED off, no alarm
- Cylinder half empty: green LED off, yellow LED on, alarm low
- Cylinder is empty: green LED off, yellow LED blinks, alarm medium

### 2.2.9 Starting / stopping insufflation

When one of the three mode buttons (1L, 5L, 10L) is pressed, the corresponding valve combination is switched to generate the desired flow. The selected CO<sub>2</sub> flow is indicated by illumination of the corresponding selector pushbutton.

For safety reasons insufflation can only be started by actuating the 1L/min button. Only then can the CO<sub>2</sub> flow be increased to 5L/min or 10L/min.

Actuation of the "flow off" button switches all valves back to the initial status and the illumination of the selected flow button goes off.

#### 2.2.10 Preselecting the nominal pressure

The nominal pressure changes when one of two ("<" or ">") buttons are pressed. Valve 4 is energized and the servomotor on the pressure reducing valve turns clockwise or counterclockwise, depending on the button depressed.

As long as one of the control buttons remains depressed the pressure reducer continues to be adjusted while the nominal pressure is updated on the illuminated bar indicator. The nominal pressure can only be adjusted between the fixed limits of 5 and 25 mmHg.

When the control button is released the servomotor stops immediately and – after a brief delay – valve 4 switches back to the initial status. Then the displayed nominal value is deposited in the RAM and stored in the EEPROM. This saves the value when the device is switched off.

#### 2.2.11 Measuring the CO<sub>2</sub> flow

The voltage value at the flow sensor is read in via an AD converter input of microprocessor 80C552 and deposited in the RAM for further evaluation.

#### 2.2.12 Calculation and display of flow value

The flow value measured is calculated in such a way that the maximum possible measurement value (15L/min) corresponds with the maximum value of the respective indicator/display.

In the H-FLOW-ANZ subroutine the display driver protocol is generated and output via I<sup>2</sup>C bus.

#### 2.2.13 Measuring the actual pressure

The voltage value at the pressure sensor for actual pressure is read in via an AD converter input of microprocessor 80C552 and stored in the RAM for further evaluation.

#### 2.2.14 Calculation and display of actual and nominal pressures

The value measured actual pressure or nominal pressure is calculated in such a way that the maximum possible measurement value corresponds with the maximum value of the respective indicator/display.

The calculated measurement values are output via I<sup>2</sup>C bus in accordance with the display driver protocol.

#### 2.2.15 Integration of the CO<sub>2</sub> consumption

To obtain a precise indication of the CO<sub>2</sub> volume already consumed, the current flow measurement value is added up every 71.07 ms.

The display value is output via I<sup>2</sup>C bus.

### 2.2.16 Generating the insufflation cycle

This subroutine defines the duration of the insufflation time and the duration of the measuring period. For this, the actual pressure is compared with the nominal pressure.

If the nominal pressure exceeds the actual pressure, 1.991 s is loaded into the counter register for insufflation time and 0.995 s for the measuring time. When the difference between actual pressure and nominal pressure has decreased to 3 mmHg, 1.137 s is loaded into counter register for insufflation time and 0.996 for measuring time.

Identification bits are set as soon as the actual pressure has reached the nominal pressure, and also when the actual pressure exceeds the nominal pressure by 5 mmHg. Before the insufflation time has elapsed the CO<sub>2</sub> flow measurement is carried out and the insufflation pressure measured. Before the measuring time has elapsed the intraabdominal pressure (actual pressure) and the measuring flow are measured. If the measuring flow is below the defined limit value, CO<sub>2</sub> insufflation is interrupted for the time during which the measuring flow remains below this value, and the actual pressure is permanently updated.

A started insufflation cycle will be completed even if overpressure occurs during the insufflation time. This prevents possible oscillations of the closed-loop control circuit.

### 2.2.17 Activation of valves, pushbuttons, LEDs, alarms and indicators

These elements are controlled by I/O extensions via I<sub>2</sub>C bus.

#### 2.2.18 Valve 1

In the case of a fault which cannot be repaired, the valve can be operated manually during the intervention. (emergency insufflation). A flow of approx. 5 l/min is achieved.

#### 2.2.19 Valve 2

Valve 2 switches over between a flow rate of 1 L/min and a flow rate of 5 L/min, the switch position "energized (electrically on)" representing 1 L/min flow rate and "de-energized (electrically off)" representing 5 L/min flow rate.

With increased flow (>2 L/min or >7 L/min) valve V2 will not switch. Alarm E04 is displayed.

Valve 2 can only be energized when:

- the motor for nominal pressure preselection has not been energized
- the 1 L/min flow rate is or has been selected.

#### 2.2.20 Valve 3

Valve 3 switches over between a flow rate of 10 L/min and a flow rate of 5 L/min or 1 L/min, the switch position "energized" (electrically on) representing a flow rate of 10 L/min and "de-energized" (electrically off) representing a flow rate of 1 L/min or 5 L/min.

Valve 3 is switched when:

- the motor for for nominal pressure preselection has not been energized
- the flow rate of 10 L/min is or has been selected.
- there is no error message E04.

When a flow rate of 10 L/min is selected, the actual flow must at least be 6 L/min. If this flow rate is not reached, the unit automatically switches to a flow rate of 5 L/min; in this case the 10 L/min flow indicator will continue to light up. After 2 insufflation cycles the unit switches back to 10 L/min if a flow rate of more than 3.5 L/min is reached.

The switch-over to 10 L/min is effected in two steps.

#### 2.2.21 Valve 4

In the switch condition "electrically on", valve 4 bleeds the measuring line during nominal pressure adjustment.

Valve 4 is switched when:

- the motor for nominal pressure preselection is energized.
- writing into the EEPROM has been prepared
- The alarm button has not been pressed.

#### 2.2.22 Valve 5

Valve 5 generates the insufflation cycle and – depending on its switching state – changes between measuring time and insufflation time.

Valve 5 is energized (electrically on) when:

- it has been on before (locking effect)
- the nominal pressure is not yet reached
- the measuring flow exceeds the nominal value
- insufflation time has not yet elapsed
- one of the modes (1L, 5L or 10L of flow) is or has been selected
- the actual pressure exceeds the nominal pressure by no more than 5 mm Hg
- the motor for nominal pressure preselection is not energized
- no "actual pressure > 30 mmHg" alarm is triggered.

#### 2.2.23 Valve 6

Valve 6 acts as a safety valve. It closes when the actual pressure exceeds the nominal pressure by more than 5 mmHg and when there is an error message (E01, E04).

Valve 6 is energized (electrically on) when:

- one of the mode buttons (1L, 5L or 10L/min) has been selected
- the motor for nominal pressure preselection is not energized
- the actual pressure exceeds the nominal pressure by no more than 5 mmHg
- the "actual pressure > 30 mmHg" alarm is not triggered.
- when a smoke gas recirculation leakage test is being carried out.


#### 2.2.24 Valve 7

Valve 7 opens when the actual pressure exceeds the nominal pressure by 5 mm Hg. This allows for a release of the excessive pressure via valve 7.


Valve 7 is energized (electrically on) when:

- it has been on before (locking effect)
- or the actual pressure exceeds the nominal pressure by more than 3 mmHg
- or the actual pressure exceeds the nominal pressure by more than 5 mmHg
- or the actual pressure exceeds 30 mmHg
- and the insufflation cycle is "ON".

### 2.3 Typical layout of a control using an I<sup>2</sup>C bus

Data transmission

## 3 Setting Instructions (64 640.121–R01)

### 3.1 Power supply PCB (64 151.716)

#### 3.1.1 Test and measuring equipment

- Multimeter

#### 3.1.2 Fuses

Fuses: Si1 T 2.5 A slow  
Si2 T 315 mA slow  
Si3 T 1.25 mA slow  
Si4 T 2.5 A slow

#### 3.1.3 Supply voltages on the power supply PCB

- ◇ Switch on power switch, all LEDs on the PCB must light up.
- ◇ Check the following voltages with the multimeter:
  - ST2 Pin7 = + 5V  $\pm 0.3V$  (GND1 = Pin 8)
  - ST2 Pin6 = + 5V  $\pm 0.3V$  (GND3 = Pin 10)
  - ST2 Pin1 = + 12V  $\pm 0.3V$  (GND3 = Pin 9)
  - ST2 Pin4 = + 23V  $\pm 1.0V$  (GND3 = Pin 4)
  - ST2 Pin5 = + 25V  $\pm 1.0V$  (GND4 = Pin 4)

### 3.2 Mother board (64 151.715)

#### 3.2.1 Test and measuring equipment

- Multimeter
- Digital pressure gauge
- Monitor with S–video and FBAS connection (BNC)
- CCD camera with S–video and FBAS connection (BNC)
- PC keyboard, with switch–selectable XT mode (5050.801)

#### 3.2.2 Function check

- ◇ Establish all pneumatic and electric connections, either S–video or FBAS.  
(camera "OUT" to CO2 Pneu "IN", CO2 Pneu "OUT" to monitor "IN", FBAS–BNC connection)
- ◇ Connect PC keyboard.
- ◇ Switch on power switch.
- ◇ Open gas connection, gas supply indicator must light up green.

#### 3.2.3 Adjusting the reference voltage

- ◇ Connect multimeter to MP1 and MP2 (GND) on the motherboard
  - ◆ Use TP1 to adjust a voltage of 5 V  $\pm 0.05$  V.

#### 3.2.4 Adjusting the headline text, time and date

- ◇ Switch on monitor and camera.
  - ◆ The "text" button (TA9) must light up.
- ◇ Press the F7 function key on the PC keyboard.

- ◆ "Headline text input" appears on the monitor.
- ◇ Enter the required text with the PC keyboard and press the return key.
- ◇ Press the F8 function key on the keyboard and set the current date.
  - ◆ Sequence: Day XX ; month XX ; year XX.
- ◇ Press the return key.
  - ◆ The date must appear on the monitor on the upper right hand side in the first line.
- ◇ Press function key F9 on the keyboard and adjust the current time
  - ◆ Sequence: Hour XX ; minute XX ; second XX.
- ◇ Press the return key.
  - ◆ The time must appear on the monitor on the right hand side in the second line.
- ◇ Switch on the CO<sub>2</sub> Pneu automatic insufflator and switch it off again after 5 seconds.
  - ◆ The headline text, the date and the time must appear on the monitor as previously set.
- ◇ In the lower right hand corner of the monitor the measurement values for pressure (mmHg), flow rate (l/min) and consumption (l) must be displayed.
  - ◆ The values must correspond with the values displayed on the front PCB.

### 3.2.5 Function check of valves and adjustable pressure reducer

- ◇ Press the pressure preselection ">" button on the front PCB.
  - ◆ Adjust a value of 25mmHg on the preselected pressure display.
- ◇ Connect the pressure gauge to the "insufflation" output and press the "flow 1 l/min" button (TA1).
  - ◆ The pressure displayed on the pressure gauge must be 33 mbar  $\pm$  2mbar. The preselected pressure and the pressure indicated on the two bar indicators and the pressure display on the digital display must correspond.
  - ◆ Permissible deviation:  
Preselected pressure versus actual pressure  $\pm$  2mmHg.
- ◇ Remove pressure gauge.
- ◇ Select flow rate of 5 l/min (TA2).
  - ◆ A flow of approx. 6 l/min must be achieved.
- ◇ Select flow rate of 10 l/min (TA3).
  - ◆ A flow of approx. 11 l/min must be achieved.
- ◇ Select flow rate of 1 l/min (TA2).
  - ◆ A flow of approx. 1 l/min  $\mp$  0.2 l/min must be achieved.
- ◇ Switch off the flow (TA4)

### 3.2.6 Alarm function

- ◇ Actuate the "manual insufflation" button on the rear panel of the device with a pin and close the "insufflation" output.
  - ◆ An alarm sound must be heard and valve V7 must release the pressure.
  - ◆ On the "intraabdominal pressure" indicator a brief pressure drop is indicated.

### 3.2.7 Control of smoke gas pump

- ◇ Close the "insufflation" output and switch on a flow rate of 1 L/min.
- ◇ Press the "pump off" button (TA1) for approx. 1 s until an intermittent sound is heard.
  - ◆ after approx. 10 s a signal sound indicates that the pump is enabled to start.
- ◇ Press the stage I (TA II), II (TA4) and III (TA3) gas recirculation buttons one after the other.
  - ◆ The pump must work at three different speeds.
- ◇ Press the "pump off" button (TA1) briefly.
  - ◆ The LED of button TA3 will go out.

### 3.2.8 Volume (loudness) adjustments

- ◇ Press the "alarm reset" button (TA8) and keep it depressed until an intermittent sound is heard.
- ◇ Push the "plus" pressure selection (TA6) several times briefly.
  - ◆ The loudness is increased to the the maximum.
- ◇ Push the "minus" pressure selection (TA5) several times briefly.
  - ◆ The loudness is decreased to the the minimum.

## 3.3 PCB for video data insertion (64 151.746)

### 3.3.1 Test and measuring equipment

- Oscilloscope (e.g. Philips PM 3350)
- Test pattern generator (e.g. Philips PM 5515)
- Keyboard 5050.801
- Monitor (e.g. Sony PVM 1443MD)
- Camera with S–video output
- 1 BNC T–connector
- 3 BNC cables
- 2 S–video cables

### 3.3.2 Preparation


- ◇ Connect keyboard 5050.801
- ◇ Switch on test generator and select "chess board pattern"
- ◇ Switch on monitor (input A)
- ◇ Switch on oscilloscope (PM3350): 0.2V/Div ; T= 2ms ; "Auto" trigger mode ; couple "TVLP" ; source "A".

**NOTE!**

*The camera must not be switched on during the first part of the adjustment.*


### 3.3.3 Test setup

### 3.3.4 Adjustment of video signals on oscilloscope

- ◇ All potentiometers to center position.
- ◇ Press function key F1 on keyboard; the box must appear on the monitor (box mode "non-transparent")
  - ◆ Write "EEEEEEEEEEEE" on 2 lines.
- ◇ Use TP2 to adjust the black level for "text" to the same level as the black level for "chess board pattern".
- ◇ Use TP1 to adjust the white level for "text" to the same level as the white level for "chess board pattern".
- ◇ Switch box mode to "transparent" with the F3 key (keyboard)
- ◇ Press the F3 key.
  - ◆ The white characters must be set (if not, select with F4)
- ◇ Use TP5 to adjust the "box" black level to the same level as the black level for "chess board pattern".
- ◇ Use TP3 to adjust the "box" video amplitude to 350 mV from the black level.
  - ◆ With this adjustment the position of on the black level will generally be shifted; readjust until the correct setting is achieved.
- ◇ Switch to black characters with F4.
- ◇ Use TP4 to adjust the "box" white level to the same level as the white level for "chess board pattern".
- ◇ Switch off test pattern generator.
- ◇ Switch on camera.
- ◇ Switch monitor input to Y/C (S-video)
  - ◆ text insertion must appear on the monitor.

### 3.3.5 Oscilloscope image

### 3.4 Adjustment of pneumatic section

#### 3.4.1 Test and measuring equipment

- Multimeter
- Digital pressure gauge (0 – 150mmHg)
- Flow meter 0 – 10 l/min CO<sub>2</sub>
- Pressure gauge 0–6 bar, category 2

#### 3.4.2 Adjustment of pressure reducer

(for location of pressure reducers see location chart 64 220.022)

- ◇ Open valve on gas cylinder.
- ◇ Adjust pressure reducer DM 1 to 3 bar – 3.5 bar.
- ◇ Pull off hose at the needle valve inlet.
  - ◆ Connect digital pressure gauge and adjust a pressure of 115 mmHg at DM 2.
- ◇ Connect digital pressure gauge to the "insufflation" outlet and press the "manual insufflation" button on the rear panel.
  - ◆ Adjust a pressure of 45 mmHg at DM 3 and remove the pressure-gauge.

#### 3.4.3 Adjustment of servomotor (only after replacment of servomotor!)

- ◇ Connect cable of servomotor (with servomotor removed).
- ◇ Press pressure preselection button ">" to operate the servomotor until it has reached its outer stop.
- ◇ Disconnect plug and turn motor one revolution into the pressure reducer.
- ◇ Press pressure preselection button ">" and screw in motor until a pressure of approx. 20 mmHg is displayed on the "preselected pressure" bar indicator.
- ◇ Then press the pressure preselection button "<" and continue to screw in motor until a pressure of 30 mmHg is read.
- ◇ Lock motor with (lock) nut.
- ◇ Reconnect plug of servomotor.
- ◇ Press the pressure preselection button "<" until a pressure of approx. 20 mmHg is displayed
- ◇ Then press the pressure preselection button ">".
  - ◆ The pressure must be limited to 25 mmHg.

#### 3.4.4 Adjustment of flow

- ◇ Switch on a flow rate of 10 L/min
  - ◆ A flow of 12 l/min  $\pm$ 2 l/min must be achieved. (Digital display "gas flow")
- ◇ Switch on a flow rate of 5 l/min.
  - ◆ A flow of 6 l/min  $\pm$ 1 l/min must be achieved.
- ◇ Connect flowmeter and compare value.
- ◇ Switch on a flow rate of 1 l/min.
  - ◆ Use adjusting screw on DM 2 to adjust the throttle to a flow rate of 1 l/min  $\pm$ 0.2 l/min in accordance with the "gas flow" digital display. Switch off flow.

### 3.4.5 Adjustment of pressure switch (up to SNo 200)

- ◇ Remove plug from valve V7 and switch on the "CO2 Pneu" automatic insufflator.
- ◇ Set multimeter to voltage measurement and connect it to contacts 1 and 3 of the pressure switch.
  - ◆ No voltage must be measured (contact closed)
- ◇ Adjust the pressure switch to 40 mbar + 2 mbar.
- ◇ Press the "manual insufflation" switch on the rear panel of the unit with a pin.
  - ◆ After 10s the multimeter must read a voltage of 24 V (contact open)
- ◇ Release button and reconnect plug to valve V 7.
- ◇ Secure adjusting screw on pressure switch with glyptal lacquer.

### 3.4.6 Adjustment of pressure switch (as of SNo 201)

- ◇ Open CO<sub>2</sub> cylinder.
- ◇ Adjust the pressure switch to 40 mbar +2 mbar.
- ◇ Secure adjusting screw on pressure switch with glyptal lacquer.
- ◇ Switch on the "CO2 Pneu" automatic insufflator.  
 (Before this, the device must have been switched off for at least 1 min.!)
- ◇ When the device is switched on, the pressure switch is automatically checked for proper functioning.
  - ◆ The gas supply indicator lamps (green, yellow) flash alternately, and the check is acknowledged with a short acoustic signal.
  - ◆ A faulty pressure switch is indicated by the error message "E- 03".

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
Type		2231	
Name		Lageplan component drawing	
Date		30.09.01	
Author		V. Müller	
Checked		01.10.01	
Approved			
Scale		1:1	
Start Page		64 220.022	
End Page			
Name		Name	
Date		Date	
Checked		Checked	
Approved		Approved	

These drawings are valid until, in all cases, they are replaced by more recent ones, especially those of competing companies, which are accessible to our customers and remain our property.

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