Service Manual

Arthro Pump
2202

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...

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Important general instructions for use

Ensure that this product is only used as intended and described in the instruction manual, by adequately trained and qualified personnel, and that maintenance and repair is only carried out by authorized specialized technicians.

Operate this product only in the combinations and with the accessories and spare parts listed in the instruction manual. Use other combinations, accessories and wearing parts only if they are expressly intended for this use and if the performance and safety requirements are met.

Reprocess the products before every application and before returning them for repair as required by the instruction manual in order to protect the patient, user or third parties.


Subject to technical changes!

Due to continuous development of our products, illustrations and technical data may deviate slightly from the data in this manual.

CAUTION - USA only:

Federal law restricts this unit to be used or sold, except under the supervision of a medical doctor.

Safety instructions and levels of danger

Symbol	Level of danger
	WARNING! Failure to observe can result in death or severe injury.
	CAUTION! Failure to observe can result in slight injury or damage to the product.
	IMPORTANT! Failure to observe can result in damage to the product or surrounding.
	NOTE! Tips for optimum use and other useful information.

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ENGLISH

Contents

1	Maintenance	1
1.1	Important notes	1
1.2	Maintenance of device and accessories	1
1.3	Visual check	2
1.4	Electrical safety check	2
1.5	Functional Check	3
1.5.1	Introduction	3
1.5.2	Basic Functional Test	3
1.5.3	Pressure sensors	3
1.5.4	Seating of Pressure Chamber	4
1.5.5	Safety function in case of overpressure	4
1.5.6	Device shut-down in case of failure to detect rotary pulse	4
1.5.7	Pressure Differential Monitoring	4
1.5.8	Checking the Pressure Sensor Zero Point	5
1.5.9	Checking the Pedal Switch	5
1.5.10	Checking Automatic Instrument Recognition	5
1.5.11	Checking the cable remote control	5
2	Assemblies	6
2.1	Service Options	6
2.1.1	Calling the service menu	6
2.1.2	Pressure Measurement	7
2.1.3	Limit Value Calibration - Pressure Measurement	7
2.1.4	Offset Voltage of Motor	7
2.1.5	Speed Calibration	7
2.2	Error Messages	8
2.3	Replacing Assemblies	9
2.3.1	Replacing the Power Supply Components (panel-mounted 3-pole power plug)	9
2.3.2	Replacing the Pump Wheel	9
2.3.3	Replacing the Motor	9
2.3.4	Replacing the Control Module	10
2.3.5	Replacing the Switching Power Supply	10
2.3.6	Replacing the Line Filter	10
2.3.7	Replacing the UNI PCB	10
2.3.8	Replacing the KEY PCB	11
2.3.9	Replacing the microswitch	11
3	Anhang / Annex	12
3.1	Reparaturteile / Repair parts	12
3.2	Explosionszeichnungen / Exploded view	13
3.2.1	Legende / Legend	13
3.3	Funktionsschema / Operational Diagram	14
3.3.1	Legende / Legend	14
3.4	Lageplän E-Karten / Component drawing PCB	15
3.5	Blockschaltbild / Block Diagram	16
4	Protokolle / Reports	17
4.1	Prüfprotokoll / Test Report	17
4.2	Wartungsprotokoll / Maintenance report	19

1 Maintenance

1.1 Important notes

This service manual describes the external service measures defined for the product.

IMPORTANT!

Follow the instruction manual of the product when carrying out the service measures.

1.2 Maintenance of device and accessories

In order to protect the person performing the test, carry out maintenance and checking or testing must be carried out in the order specified below.

- ◇ Visual check
- ◇ Electrical safety check
- ◇ Functional check

The required electrical checking and measuring procedures are described in standards EN 60601 / IEC 601 and DIN 57751 / VDE 751.

The values indicated under the section entitled safety check relate to the test to EN 60601 / IEC 601.

After a repair job all measuring values (specified values) of the device must be checked in accordance with the instructions given in the Service Manual and readjusted if deviations are observed.

NOTE!

Any maintenance, checking or testing work on the device or the accessories must be documented.

1.3 Visual check

The disinfectants used in the operating theatre as well as the ozone freed when UV radiation sources are used, may change the surfaces of device parts.

Item to be checked	Check for
Device and accessories	<ul style="list-style-type: none"> ◆ Any soiling/contamination jeopardizing safety and general cleanliness ◆ Mechanical damage ◆ Loose or missing parts
Controls	◆ Mechanical function and easy operation
Front plate and housing	<ul style="list-style-type: none"> ◆ Any damage having an effect on function and safety ◆ Function of control panels and LEDs
Roller wheel	◆ Perfect condition, correct seating and easy operation
Membranes and covers	<ul style="list-style-type: none"> ◆ Undamaged membranes of pressure sensors ◆ Undamaged cover of pressure chamber switch
Labeling / Symbols	◆ Completeness and good legibility
Safety-relevant inscriptions (such as warnings)	◆ Completeness and good legibility
Fuse inserts	◆ Values indicated by the manufacturer on the identification plate (nominal current and melting characteristics)
Cabling	◆ Perfect condition (position, insulation and possible brittleness)
Printed circuit boards	◆ Corrosion or other damage
Lamp plate (only light sources)	◆ Existence and good legibility
	◆ Replace damaged parts immediately!

1.4 Electrical safety check

Item to be checked	Check for
Protective earth (PE) connection	<p>Measurement values:</p> <ul style="list-style-type: none"> ◆ With supply cable not connected $\leq 0.1 \text{ Ohm}$ ◆ With supply cable connected $\leq 0.2 \text{ Ohm}$ <p>Test conditions: $10 \text{ A} \leq I_{\text{meas}} \leq 25 \text{ A}$, $V_{0\text{meas}} \leq 6 \text{ V}$, $t_{\text{test}} \geq 5 \text{ s}$</p> <p>Check: The resistance between protective earth (PE) wire, protective earth (PE) connection and the conductive parts connected.</p>
Insulation resistance	<p>Measurement values:</p> <ul style="list-style-type: none"> ◆ Insulation resistance $\geq 2 \text{ MOhm}$ <p>Test conditions: Test voltage 500V DC (continuous voltage)</p> <p>Check: All mains (power) cables relative to earth.</p>
Leakage current	<p>Measurement values:</p> <ul style="list-style-type: none"> ◆ Device leakage current $\leq 100 \mu\text{A}$ <p>Test conditions: Measuring setup (MD) in accordance with IEC 601</p> <p>Check: The leakage current which flows from each pole of the power supply unit via the insulation through the protective earth (PE) wire or from the application part via the patient to earth.</p>

1.5 Functional Check

1.5.1 Introduction

This introduction should help you to better understand the functionality of this device in order to avoid operator's errors and measuring errors.

1.5.2 Basic Functional Test

In the basic functional test the indicators, pushbuttons and the delivery rate of the device are checked.

Preparation:

- ◇ Install the flexible tube set.
- ◇ The pressure chamber must click into place.
- ◇ Connect the irrigation tube to the fluid bag. The fluid bag must be suspended 1.0 m above the device.
- ◇ Place the instrument tube into a measuring vessel.
- ◇ Switch on the device with the mains/power switch. Select the following values:

Nominal Pressure:	150 mm Hg
Nominal Flow:	1.00 l/min
- ◇ Press the start/stop button. The green start/stop LED lights up. The pump wheel starts rotating.
- ◇ Fill the entire tube set with liquid and let the pump run for 1 min. to complete automatic instrument recognition.
- ◇ Then cross-clamp the instrument tube. Do not press the start/stop button.

Carry out the following:

- ◇ Empty the measuring vessel.
- ◇ Remove the cross clamp from the instrument tube and time 1 min.
- ◇ After 1 min press the start/stop button. The measuring vessel must contain 1.0 l ($\pm 10\%$) water.

If you obtain these values the basic functional test has been successfully completed.

1.5.3 Pressure sensors

- ◇ Switch on the device.
- ◇ Use a syringe to generate pressures of 50, 100 and 150 mm Hg (mercury column).
- ◇ The values of the pressure gauge and the actual pressure reading may deviate with a tolerance of $\pm 6\%$.

In the case of a wrong reading the device needs calibration (Chapter 2.1.2). If the device won't work properly after calibration, replace the AM3-UNI PCB.

1.5.4 Seating of Pressure Chamber

- ◇ Use the flexible test tube set to check the installation of the pressure chamber in the pump head.
- ◇ The pressure chamber must be easy to install. The pump head and the pressure chamber must be flush relative to each other.

In the case of malfunction replace the control module assembly.

1.5.5 Safety function in case of overpressure

- ◇ Select a nominal pressure of 150 mm Hg and a nominal flow of 1.5 l/min.
- ◇ Actuate the start/stop button.
- ◇ Use a syringe to generate a pressure of more than 150 mm Hg (<250 mm Hg).
- ◇ The pump wheel must turn backwards.
- ◇ After approx. 5 s a warning is sounded that is repeated every 5 s.
- ◇ Press the start/stop button again and release the pressure in order to complete this test.

In the case of malfunction the device needs calibration (Chapter 2.1.3). If the device still won't work after this, replace the AM3-UNI PCB.

1.5.6 Device shut-down in case of failure to detect rotary pulse

- ◇ Select a nominal pressure of 150 mm HG and a nominal flow of 0.20 l/min.
- ◇ Actuate the start/stop button.
- ◇ Carefully hold the pump wheel with your hand while it is rotating slowly.
- ◇ A warning signal is sounded and an error message 'E02' appears on the nominal pressure display.
- ◇ Switch off the device.

In the case of malfunction the device needs calibration (Chapter 2.1.5). If the device won't work properly after this, replace the AM3-UNI PCB. a

1.5.7 Pressure Differential Monitoring

- ◇ Remove the flexible tube set and switch on the device.
- ◇ Press the start/stop button to start the pump.
- ◇ Carefully tap a pressure sensor with your finger.
- ◇ Error message 'E01' is displayed on the actual pressure indicator and a warning signal is sounded
- ◇ Release the pressure sensor. The error message disappears.
- ◇ Repeat the test with a different pressure sensor.
- ◇ Switch off the device.

In the case of failure the device requires calibration (Chapter 2.1.2). If the device won't work after this, replace the control module assembly and/or the AM3-UNI PCB.

1.5.8 Checking the Pressure Sensor Zero Point

- ◇ Call menu item P2 in the user menu (Chapter 6). In P2 the values measured by the pressure sensors are displayed on the nominal pressure and actual pressure displays. These values must be approximately the same (± 2 mm Hg).
- ◇ Switch off the device. This completes the function test.

In the case of malfunction the device requires calibration (Chapter 2.1.3). If the device won't work after that, the control module assembly and/or the AM3 UNI PCB must be replaced.

1.5.9 Checking the Pedal Switch

- ◇ Switch on the device.
- ◇ Connect the fluid container and a flexible tube with LL.
- ◇ Connect the pedal switch to the corresponding connector on the rear panel of the device.
- ◇ Press the pedal switch. The nominal flow indicator will read '1.50' l/min and the nominal pressure indicator will display the value preselected in the User Manual. The 'Pedal switch/Wash' indicator lights up.
- ◇ Press the pedal switch again and the nominal flow and nominal pressure readings will return to the previous value.

1.5.10 Checking Automatic Instrument Recognition

**2202.001 as of serial number 98231,
2202.011 as of serial number 98100**

- ◇ Select a nominal pressure of 100 mm Hg (MC) and a nominal flow of 1.50 l/min. Touch the start/stop button.
- ◇ The device will automatically recognize the connected instrument. During this procedure the pump wheel may run briefly at various speeds and may stop several times.
- ◇ After this the pump wheel will run at high speed. The device detects the lack of pressure, i.e. that no instrument is connected to the end of the tube.

1.5.11 Checking the cable remote control

**2202.001 as of serial number 98231,
2202.011 as of serial number 98100**

- ◇ Connect the cable remote control to the connecting socket on the front panel of the device.
- ◇ The remote control can be used to control all functions of the device.
 - ◆ Pressure: Increasing and decreasing of the nominal pressure.
 - ◆ Flow: Increasing and decreasing of the nominal flow.
 - ◆ Wash: Quick change of pressure and flow ('Wash' function).
 - ◆ Start/stop: Starting /Stopping insufflation.

2 Assemblies

2.1 Service Options

Via the service menu you can change and display device parameters. There are 4 menu items (L).

L0 Pressure Measurement

L1 Limit value calibration for pressure measurement


L2 Offset voltage of motor

L3 Speed calibration

 **IMPORTANT!**

Always perform the menu items again after replacing an assembly.

2.1.1 Calling the service menu


◇ Switch on the device.

◇ While the actual pressure reading ' - - - ' blinks (1), touch the p_{minus} - and p_{plus} buttons (2. Observe actuation sequence!). L1 appears on the nominal flow display. The 'L' blinks (3).

◇ Push the start / stop button (4). The 'L' stops blinking. The nominal pressure display (5) (and if applicable the actual pressure display) read the current parameter of the menu level activated.


◇ With the p_{plus} or the p_{minus} button (6) you can change the parameter of the menu level.

◇ Press the start / stop button. The 'L' starts blinking. The selected parameter is stored.

◇ With the Q_{plus} or Q_{minus} button (7) a different menu level can be selected.


◇ Switch off the device. After switching on again you may work with the new parameters.

2.1.2 Pressure Measurement

- ◇ For this menu item the pressure chamber must not be in the pump head.
- ◇ Select menu level L0.
- ◇ The nominal pressure and actual pressure displays read the current pressure values detected by the pressure sensors. The two readings must be '115' (± 2).
- If this is not the case, adjust the two readings to '115' (± 2) by means of potentiometers P1 and P2 on the AM3-UNI PCB.
- ◇ After readjustment of the zero points calibrate menu items L1 - L3.
- ◇ Press the start / stop button to exit menu level L0.

2.1.3 Limit Value Calibration - Pressure Measurement

- ◇ Place the flexible test tube set into the pump head.
- ◇ Select menu level L1.
- ◇ Select a nominal pressure of 200 mm Hg (MC).
- ◇ Generate a pressure of 200 mm Hg (267 mbar) externally.
- ◇ Adjust the nominal pressure and actual pressure readings to the same value.
- ◇ The instrument display reads a digital value.
- ◇ If the values are the same, press the start / stop button in order to exit menu level L1.

2.1.4 Offset Voltage of Motor

- ◇ Remove the flexible test tube set.
- ◇ Select menu level L2.
- ◇ The device will automatically look for the motor offset voltage.
- ◇ When the 'L' in the nominal flow display blinks, the check is completed. You may now select a new menu level.

2.1.5 Speed Calibration

- ◇ Place a genuine tubing set into the pump head. This test will be performed without fluid delivery.
- ◇ Select menu level L3.
- ◇ The pump roller wheel turns. The nominal pressure display reads the rotations per minute (rpm).
- ◇ Potentiometer P3 on the AM3-UNI PCB must be adjusted in such a way that a speed of 355 rpm (± 5) is achieved.
- ◇ This menu item cannot be exited by means of the start / stop button. Switch off the device with the mains/power switch.

2.2 Error Messages

Error Messages and Warnings	Cause	Corrective Action
Service warning lights, 'E01' appears on the actual pressure display A warning signal is sounded	The pressure differential between the pressure sensors is greater than 50 mm Hg (Mc)	Calibrate the device (Chapter on Service menu). Then check the user menu item P2.
Service warning lights, 'E02' appears in the actual pressure display, A warning signal is sounded	The motor control circuit is defective.	Check all supply cables to the motor (loose cable etc.).
Service warning lights, 'E03' appears on the actual pressure display	The device is not calibrated.	Calibrate the device (Chapter on Service menu).
Service warning lights, 'E05' appears on the actual pressure display	Error on AM3-UNI PCB	Replace AM3-UNI PCB.
Service warning lights , 'E07' appears on the actual pressure display	Hardware error on the AM3-KEY PCB	Replace AM3-KEY PCB.
Two short warning signals are sounded	Timeout error	Replace AM3-UNI and AM3-KEY PCBs.
A long and a short warning signal are sounded consecutively	The AM3-KEY is defective.	Replace AM3-KEY PCB.
Display is dark	AM3-KEY PCB defective	Replace AM3-KEY PCB.
Pump wheel is stopped	Microswitch jammed or rollers on pump wheel do not turn easily	Replace microswitch or pump wheel.
Two short warning signals are sounded consecutively	Internal data error	Replace AM3-UNI PCB.

CAUTION!

For Model 2202.001 as of Serial Nr. 98231 and 2202.011 as of Serial Nr. 98100 a new key card is used.

2.3 Replacing Assemblies

2.3.1 Replacing the Power Supply Components (panel-mounted 3-pole power plug)

- ◇ The power supply components consist of the panel-mounted 3-pole power plug, plastic cover and cable.
- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Disconnect the screws and cables from the power supply unit.
- ◇ After having installed a new power supply unit perform the following tests:
 - ◆ Electrical safety tests

2.3.2 Replacing the Pump Wheel

- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Turn the pump wheel by hand until the grub screw faces upward. Loosen the grub screw.
- ◇ Pull the pump wheel off the shaft and install a new one.
- ◇ Tighten the grub screw insuring that it rests against the flat surface provided on the shaft. After tightening the grub screw check the pump wheel for proper alignment with the pump head. If it is not properly aligned, realign the pump wheel.
- ◇ After installing the new pump wheel perform the following tests:
 - ◆ Electrical safety tests

2.3.3 Replacing the Motor

- ◇ Make sure that the device is disconnected from the power supply (mains) Disconnect the power cable from the device.
- ◇ Remove the pump wheel as described under item 6.3.2.
- ◇ Unscrew the mounting screws and disconnect the cables from the motor.
- ◇ After having installed the new motor perform the following tests:
 - ◆ Electrical safety tests
 - ◆ Service menu
 - ◆ Functional check

2.3.4 Replacing the Control Module

- ◇ The control module consists of the design foil, front panel, upper part of pump head and lower part of pump head and is only available as a pre-assembled unit.
- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Remove the power switch, pump (roller) wheel with motor and micro switch.
- ◇ Unscrew the mounting screws and disconnect the cables of the control module.
- ◇ After installation of the new control module perform the following tests:
 - ▶ Electrical safety tests
 - ▶ Service menu
 - ▶ Functional check

2.3.5 Replacing the Switching Power Supply

- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Remove the cover of the switching power supply.
- ◇ Unscrew the mounting screws and disconnect the power supply cables.
- ◇ After installation of a new switching power supply perform the following tests:
 - ▶ Electrical safety tests

2.3.6 Replacing the Line Filter

- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Unscrew the mounting screws and disconnect cables from the line filter.
- ◇ After having installed the new line filter, perform the following tests:
 - ▶ Electrical safety tests

2.3.7 Replacing the UNI PCB

- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Unscrew the mounting screws and disconnect the cables from the PCB.
- ◇ After installation of the new PCB perform the following tests:
 - ▶ Electrical safety tests
 - ▶ Service menu
 - ▶ Functional check

2.3.8 Replacing the KEY PCB

- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the cable from the device.
- ◇ Unscrew the mounting screws and disconnect the cables from the PCB.
- ◇ After having installed the new PCB perform the following tests:
 - ▶ Electrical safety tests
 - ▶ Service menu
 - ▶ Functional check

2.3.9 Replacing the microswitch

- ◇ Make sure that the device is disconnected from the power supply (mains). Disconnect the power cable from the device.
- ◇ Unscrew the mounting screws and cables from the microswitch.
- ◇ After having installed the new microswitch perform the following tests:
 - ▶ Electrical safety tests
 - ▶ Functional check

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RICHARD WOLF Belgium S.A.
RICHARD WOLF France S.A.R.L.