# AESCULAP®




- D** **Gebrauchsanweisung/Technische Beschreibung**  
XENON Lichtquelle OP 930
- GB** **Instructions for use/Technical description**  
XENON Light Source OP 930
- F** **Mode d'emploi/Description technique**  
Source lumineuse au XÉNON OP 930
- E** **Instrucciones de manejo/Descripción técnica**  
Fuente de luz XENÓN OP 930
- I** **Istruzioni per l'uso/Descrizione tecnica**  
Fonte di luce XENON OP 930
- NL** **Gebruiksaanwijzing/Technische beschrijving**  
XENON lichtbron OP 930

Seiten  
Pages  
Pages  
Páginas  
Pagine  
Bladzijden

Deutsch	2 – 7
English	8 – 12
Français	13 – 18
Español	19 – 24
Italiano	25 – 30
Nederlands	31 – 35

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Read these instructions for use carefully before putting the instrument into operation. This will enable you to avoid damage caused by improper assembly or operation and therefore not covered by the guarantee and warranty.

## 1 Explanation of symbols

Potential equalization connection

Type BF unit

Optical cable connection

CE mark according to Directive 93/42/EEC


Read instructions for use

LED indicator for lamp change

XENON lamp indicator

Fuses

Dangerous electrical voltage

AC voltage

Marking of electric and electronic devices according to directive 2002/96/EG (WEEE), see chapter Disposal

## Contents

- 1 Explanation of symbols
- 2 Safety information
- 3 Indications for use
- 4 Mode of operation
- 5 Installation instructions and information
- 6 Installation
- 7 Preparation for operation
- 8 Function test
- 9 Changing the optical cable socket
  - 9.1 Replacing the optical cable socket
- 10 Commissioning
  - 10.1 Switch on light source
  - 10.2 Brightness control
- 11 Replacing the XENON lamp
  - 11.1 Removing the XENON lamp
  - 11.2 Installing the XENON lamp
  - 11.3 Resetting the operating hours counter
- 12 Disinfection and sterilization
- 13 Changing the fuses
- 14 Troubleshooting
- 15 Repair
- 16 Technical data
- 17 Accessories and spare parts
- 18 Disposal

## Legend for fold-out page

- ① ON switch
- ② OFF switch
- ③ Brightness control
- ④ Optical cable connection
- ⑤ Fast-action locks
- ⑥ Optical cable
- ⑦ Light source connection
- ⑧ Instrument connection
- ⑨ Unit feet
- ⑩ Potential equalization connection
- ⑪ Mains plug
- ⑫ Unit connector
- ⑬ Fuse holder
- ⑭ XENON lamp

## 2 Safety information

### Warnings and symbols

The safety information and instructions for use must be read and observed before commissioning the unit. In these instructions for use, all parts which affect your safety have been indicated by the following designations and symbols: Warning, Caution and Note.

#### Warning:

A warning indicates a danger. Failure to observe a warning may result in very severe injuries.

#### Caution:

Caution indicates possible danger. Failure to exercise caution may lead to injuries or damage to the unit.

#### Note:

Notes provide you with important information for optimal use of the unit.

#### Attention:

„Attention“ draws attention to the fact that certain maintenance or safety measures must be taken in order to avoid damage to the unit.

### Safety information

#### Warning:

Do not open the unit. There is the danger of an electric shock. Service work must be carried out only by authorized or qualified personnel.

Danger of burns! The lamp and light outlet of the optical cable or endoscope may be hot. Do not place connected optical cables on combustible objects, such as operating theatre drapes.

Danger of explosion! The unit must not be used in areas where there is risk of explosion. The mains plug may be connected to the power supply only outside areas with a risk of explosion.

#### Caution:

Danger of overheating! Do not cover lateral ventilation slots.

## 3 Indications for use

The light source is used for illuminating operation areas and body cavities in medical applications. Its use in other fields is not permissible.

## 4 Mode of operation

An electronic lamp ballast supplies the xenon lamp, which has a high light output. By means of a cold light reflector, this lamp produces a light beam which is passed through an optical cable to the field of view and results in homogeneous illumination there.

## 5 Installation instructions and information

#### Warning:

Before the unit is commissioned, the safety information and instructions for use must be carefully read and understood. Failure to observe the safety information and instructions for use may result in very severe injuries.

#### Do not open the unit!

There is a danger of an electric shock. Service work may be carried out only by authorized or qualified personnel.

#### Danger of explosion!

The unit must not be used in areas where there is risk of explosion. The mains plug may be connected to the power supply only outside areas with a risk of explosion.

#### Caution:

Electrical installations must comply with the national regulations. Follow the national regulations especially in connection with fire and explosion prevention.

#### Attention:

Ensure that the mains voltage agrees with the voltage specified on the rating plate (100 V - 240 V).

When several units are operated together, we recommend connecting them with the potential equalization rail.

For optimal heat removal, you must maintain a distance of 50 mm from the side and 60 mm from the back panel.

## 6 Installation

#### Caution:

Danger of overheating! Do not cover lateral ventilation slots.

- Ensure that the unit is placed on a level surface.
- Ensure that the support (table, ceiling bracket, etc.) is sufficiently stable.
- Ensure that the unit is not subjected to any impacts or vibrations during operation.

## 7 Preparation for operation

#### Warning:

Danger of burns! The lamp and light outlet of the optical cable or endoscope may be hot. Do not place connected optical cables on combustible objects, such as operating theatre drapes.

Danger of dazzling! Never look into the free end of connected optical cables or endoscopes.

#### Attention:

- Connect the optical cable to the light source only when the endoscope has been connected to the optical cable. Accordingly, first disconnect the optical cable from the light source after the operation.
- The front panel firmly when pulling out the optical cable.
- First connect the mains cable plug (11) to the unit connector (12). Then connect the mains cable plug to the socket.
- Optical cables (7) may be damaged by excessive mechanical stress. Ensure that the optical cables (7) are not crushed, buckled or wound too tightly (minimum diameter 15 cm).

## 8 Function test


#### Warning:

Each time before use, the operability of the unit must be checked. Faulty plugs or cables may not be used and must be replaced.

With proper operation, the green LED indicator above the ON switch and the green LED indicator below the symbol „XENON lamp“ light up after switching on.

#### Important:

If the red LED indicator (lamp change) lights up, you should replace the XENON lamp after the operation (see Replacing the XENON lamp).

## 9 Changing the optical cable socket

This unit is delivered with an AESCULAP/Storz connection as standard. If you require an optical cable socket from another manufacturer (Wolf, Olympus, ACMI), you can replace the optical cable socket.

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### 9.1 Replacing the optical cable socket

- Unscrew the optical cable socket (4) in the counterclockwise direction.
- Screw in the new optical cable socket in the clockwise direction as far as it will go.

## 10 Commissioning

### 10.1 Switch on light source

Press the ON switch (1). Both green LED indicators light up.

### 10.2 Brightness control

You can control the brightness continuously using the rotary knob (3).

#### Note:

Switching the cold light source on and off frequently leads to increased wear of the xenon lamp. For short stoppages, the brightness should therefore be turned down on the unit.

## 11 Replacing the XENON lamp

#### Warning:

- Owing to the high internal pressure of the lamp, an explosion may occur when the lamp is in the cold or hot state. The intense UV radiation developed by the lamp during the operation is dangerous for skin and eyes. Owing to the high luminous density, severe damage to the retina may occur on observation of the arc. The lamp may therefore be put into operation only after the housing has been closed.
- You must always wear protective gloves and safety goggles when installing and removing the lamp. The ignition device of the unit generates life-threatening high voltage (up to 25 kV). Service work may be carried out only by authorized or qualified personnel.
- Before opening the unit, disconnect mains plug, since plug, cable or other accessible parts in the interior of the light source may carry life-threatening mains voltages.

#### Caution:


- The unit has no spare lamp. To minimize the probability of failure of the lamp during an operation as low as possible, change the lamp after max. 500 operating hours. Since it is possible that frequent use may result in premature lamp failure in a light source, we recommend having a spare lamp available all the time.

#### Attention:

- Do not change the xenon lamp (14) until the unit has completely cooled. Otherwise, there is a danger of burns from the hot housing or from the hot lamp.

### 11.1 Removing the XENON lamp

- Open the fast-action locks (5) of the cover by a quarter turn in the counterclockwise direction with a screwdriver or a coin.
- Hold the cover at the lower end and pull it slightly outwards. Now lift the cover upwards. Ensure that there is no tension in the earthing cable of the cover.
- Press the lateral plug tongues and pull the plug out of the socket.
- Hold the lamp by the rear ceramic base and press the wire clips at the upper ends simultaneously together. Then move the clips downwards. You can now pull out the lamp backwards.


### 11.2 Installing the XENON lamp

- Position the new XENON lamp so that the centring pins rest exactly against the ceramic holder and the positioning pin engages the groove.
- Lock the lamp in this position by engaging the two wire clips again.
- Now connect the lamp plug to the socket. In doing this you must ensure that the plug tongues snap into the socket projections.
- Ensure that the connection between cover and earthing cable is satisfactory.
- Now carefully mount the cover again and fix the cover by pressing in the fast-action locks and by locking by a quarter turn in the clockwise direction.

#### Important:

- After each lamp change, the hours counter must be reset as described under 11.3.

### 11.3 Resetting the operating hours counter

The xenon light source has an integrated operating hours counter which measures the effective operating time of the lamp. If the lamp exceeds the maximum operating time of 500 hours, the red LED indicator (lamp change) on the foil keypad lights up. After the new xenon lamp has been inserted, the operating hours counter must be reset. To do this, switch on the unit. Press both LED indicators (red and green) simultaneously for approx. 3 sec. With correct use, the red and green LEDs light up four times alternately. The green LED then lights up constantly.

### 12 Disinfection and sterilization

#### Caution:

To avoid penetration of liquid into the unit, do not store liquids on or above the unit.

Disinfection of light source and optical cable:

- Use only cleaning agents and disinfectants which are approved nationally for surface disinfection.
- Wipe away the cleaning agent and disinfectant from the light source and optical cable with a moist cloth. Then dry the unit and optical cable with lint-free cloth.
- Keep the optical cable interfaces clean.
- Do not scratch surfaces. Rub with cotton wool pad.

#### Sterilization of the optical cable:

The optical cable is sterilized at 2 bar (134°C) in an autoclave.

#### Note:

Do not wind optical cables too tightly, in order to avoid breaking the thin glass fibres in the optical cable. The optical cable must not be wound to a diameter of less than 15 cm.

#### Note:


If the light losses are too high owing to broken glass fibres, you should replace the optical cable.

### 13 Changing the fuses

#### Caution:

You should disconnect the mains plug before changing the fuses.

- Release the catch with a small screwdriver and pull the fuse holder out of the fuse housing.
- Insert the new fuses into the fuse holder. For safety reasons, only original Aesculap fuses – Art. No. TA 020436 – should be used.
- Push the fuse holder back into the fuse housing until it snaps in audibly.

### 14 Troubleshooting

Fault	Cause	Fault detection	Rectification
Lamp and fan do not operate	- No voltage to unit	- Mains plug not connected	- Connect mains plug
	- Fuses blown	- Unit does not function in spite of proper operation	- Change fuses
	- Power unit faulty	- Unit does not function even after fuse change	- Repair by manufacturer
Lamp does not burn	- ON switch not pressed	- No light output at end of optical cable	- Press ON switch
	- Lamp faulty	- No light output at end of optical cable	- Lamp change (see Section 10)
	- Power unit faulty	- No light output at end of optical cable	- Repair by manufacturer
Ventilator does not run	- Ventilator faulty	- No air output at right ventilation grid	- Repair by manufacturer
Brightness markedly poorer	- Optical cable not fully inserted	- Optical cable can be easily pulled out	- Optical cable connection (see Section 6)
	- Heat filter or lens dirty	- Visual inspection of the filter and of the lens	- Cleaning of filter and/or lens
	- Heat filter or lens faulty	- Visual inspection of the filter and of the lens	- Repair by manufacturer

## 15 Repair

### Note:

In the event of inquiries, please state the serial number and type number of the unit as shown on the rating plate at the rear.

As the manufacturer of this unit, we consider ourselves responsible for the effects on safety, reliability and performance of the unit only when:

- assembly, extensions, readjustments, modifications or repairs are carried out by persons authorized by us.
- the electrical installation of the relevant rooms complies with the IEC requirements
- the unit is used in conformity with the instructions for use.

Unauthorized opening of the unit or repairs and modifications to the unit by unauthorized persons absolves us from any liability for the operational safety of the unit. All warranty rights during the guarantee time will be lost as a result.

Only original spare parts from AESCULAP AG & Co. KG must be used for repairs.

You are required to request the repairer to provide a list showing the type and extent of the work performed, together with company stamp, date and signature.

## 16 Technical data

Design:	Maintenance-free XENON light source in stainless steel/aluminium housing with pluggable mains cable.
Power supply:	100 – 240 V / 50 – 60 Hz
Power consumption:	240 W
Fuses:	2 x 4.0 A (250 V) slow-blow
Protection class:	I
Type:	BF
Lamp:	14.5 V / 180 W XENON short arc reflector lamp
Dimensions (W/H/D):	305 x 125 x 305 mm
Weight:	6.5 kg
Radio interference suppression	IEC 601-1-2
Applicable standards:	IEC 601-1

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### Mains cable

Art. No. TE 730	Germany	Denmark
	France	Austria
	Belgium	Spain
	The Netherlands	Portugal
	Sweden	Italy
	Norway	Finland
Art. No. TE 731	Switzerland	Liechtenstein
Art. No. TE 734	United Kingdom	Ireland
Art. No. TE 735	USA/Japan	Canada
Art. No. TE 736		Cold unit plug Isolating transformer

Technical modifications may be made without notification.

## 18 Disposal

### Note

The user institution is obliged to process the product before its disposal, see chapter Processing.

Always adhere to national regulations when disposing of or recycling the product or its components!

Products carrying this symbol are subject to separate collection of electric and electronic devices. Within the European Union, disposal is carried out without fees by the manufacturer.

If you have any questions concerning the disposal of the product, please contact your national B.Braun/Aesculap agency, see chapter Technical Service.

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AESCULAP AG & CO. KG