STORZ

KARL STORZ — ENDOSKOPE

205225 20 AUTOCON® 200

205225 20 AUTOCON® 200

205225 20 AUTOCON® 200

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Instruction Manual

Manual de instrucciones

 205225 20 AUTOCON® 200

 205225 20 AUTOCON® 200

 205225 20 AUTOCON® 200

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
Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the AUTOCON® 200. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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Controls, displays, connectors, and their uses

- ① Power switch
- ② Safety panel with warning lamps
- ③ Cutting mode display panel
- ④ Digital display of power limitation for cutting
- ⑤ Activation symbol for monopolar cutting
- ⑥ Coagulation effect display panel for incision
- ⑦ Coagulation mode display panel
- ⑧ Digital display of power limitation for coagulation
- ⑨ Activation symbols for monopolar and bipolar coagulation
- ⑩ Autostart mode display panel
- ⑪ Socket for neutral electrode
- ⑫ Socket for monopolar instruments
- ⑬ Socket for bipolar instruments
- ⑭ Autostart mode preselection button
- ⑮ ± buttons for preselecting power limitation in coagulation mode
- ⑯ M button for preselecting coagulation mode
- ⑰ Button for preselecting coagulation effect in cutting mode
- ⑱ ± buttons for preselecting power limitation in cutting mode
- ⑲ M button for preselecting cutting mode (Function selection only possible on 205225 20-010)

Rear panel:

- ① Volume control for acoustic signals
- ② Socket for footswitch
- ③ Potential equalization connector
- ④ Power supply cord receptacle
- ⑤ Power fuse holder

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Socket for bipolar electrodes, footswitch or Autostart operation

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Socket for monopolar electrodes, footswitch or handgrip operation

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The neutral electrode may remain applied to the patient during defibrillation
Applied part of type CF as per IEC 60601-1

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Potential equalization connector

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Environmental protection use period of 50 years (China RoHS)

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DANGER: Risk of explosion if used in the presence of flammable anesthetics.

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CAUTION: To avoid the risk of electrical shock, do not remove cover. Refer servicing to qualified service personnel.

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Keep out of reach of patients.

Do not store liquids on or above the unit.

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Photographs of the unit	A 1
Controls, displays, connectors, and their uses	A 2
Symbols employed	A 3

General information

The principle of high frequency surgery	2
The AUTOCON® 200 unit	8

Safety instructions

Warnings and cautions	10
Intended use	15
User qualification	15
Safety precautions at the site of installation ...	17
Safety precautions when using the unit	17
Important notes on the use of the monopolar technique	18
Safety devices	23

Installation and operating instructions

Unpacking the equipment	26
Basic equipment	26
Installing and connecting up the unit	26
Operating the unit	29
Test for proper operation	34
Program functions	36

Maintenance

Fuse replacement	45
Cleaning, disinfection and sterilization	46
Maintenance	51
Safety checks and inspections	51
Servicing and repair	52
Repair program	53
Limitation of liability	54
Manufacturer's warranty	54

Technical description

Troubleshooting	55
Technical data	60
Technical documentation	62
General circuit diagram	63

Spare parts, recomm. accessories

Spare parts list	64
Accessories	65

Appendix

Cleaning agents and disinfectants	67
---	----

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The principle of high frequency surgery

Biological tissue contains a more or less high concentration of electrolytes, making it sufficiently conductive to be treated electrosurgically. The thermal effect of HF current is used for separating (cutting) and coagulating tissue (desiccation of tissue and/or hemostasis).

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
HF currents must be used on the patient since low-frequency currents can stimulate nerve and muscle cells in the current flow due to electrochemical processes (electrolysis). These effects are small enough to be disregarded with frequencies above 100 kHz. In HF surgery, a distinction is made between the following application techniques:

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
Monopolar operation techniques

With the monopolar operation technique, the electric current flows from the small-face active electrode to a large-face passive or neutral electrode, the patient's body always forming part of a closed circuit. Tissue separation or coagulation takes place at the active electrode.

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
Bipolar operation techniques

With the bipolar operation technique, the current flows between two paired electrodes, the faces of which are either the same or different sizes, depending on the application (cutting or coagulation). Only a small part of the tissue is included in the circuit between the electrodes during application.

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Separation of tissue

When separating tissue (cutting), the cell sap of the tissue is vaporized. Due to the resulting vapor pressure, the tissue cells tear and the escaping vapor pressure causes the electrode to be insulated. If the voltage is high enough, the current flow is maintained by the commencement of arcing between the electrode and tissue. The thermal effect of the current causes the separation of tissue and the formation of surface slough, depending on the size of the arc.

Arcing commences at voltage values of around 200 Vp. If voltages above 500 Vp are used, the result is substantial formation of slough due to the energy released in the arc.

For cutting, automatic control loops allow an optimal, reproducible cutting quality to be achieved, which is to the greatest possible extent independent of cutting depth, incision speed, tissue resistance and suchlike. Power output is automatically adjusted to the level required for cutting.

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ENDOCUT operation

(only on 205225 20-010)

A special problem which occurs during endoscopic operations, particularly during polypectomies and papillotomies, is that the electrodes which are to be used for cutting must be fed through the narrow working channels of flexible endoscopes using long wire pulls. Therefore the operating surgeon has no direct control over the cutting procedure. However a controlled incision is a requirement for avoiding complications, especially during polypectomies and papillotomies. A too quickly made incision can lead to bleeding due to insufficient coagulation of the cut edges. An incision made too slowly can cause thermal damage, for example, to the intestinal wall.


Using the special ENDOCUT cutting control system the incision is automatically fractionated in such a way that intermittently brief, automatic arc-controlled cutting intervals result with defined pauses. In this way, a polypectomy snare cannot, for example, cut through a polyp at any desired speed. The speed of incision and degree of coagulation on the cut edges is more uniform.

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Tissue coagulation

To coagulate biological tissue, it must be heated to approx. 70 °C. To prevent the current from causing tissue separation, the power supply, and thereby the amount of heat, must be used specifically to produce denaturation of the tissue.

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The thermal effects which take place in biological tissue, depending on temperature, can be described as follows:

- Up to approx. 40 °C: No visible cell damage
- Approx. 40 °C and above: Reversible cell damage depending on the duration of current flow
- Approx. 49 °C and above: Irreversible cell damage
- Approx. 70 °C and above: Conversion of collagens into glucose – shrinking of tissue containing collagen – hemostasis
- Approx. 100 °C and above: Dehydration/desiccation: phase transition of intracellular and extracellular water from liquid into gas. Adhesive effect of glucose after dehydration. Shrinking of the coagulum after dehydration.
- Approx. 200 °C and above: Carbonization and fourth-degree pathological burns. Unpleasant odor of burnt tissue.

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Due to the inhomogeneity of the electrical and thermal properties of tissue and the resulting uneven current density distribution within the tissue, the temperature rises with an irregular distribution when power is supplied.

The current density, and thus the temperature, is highest in the vicinity of the coagulation electrode. The temperature decreases very rapidly with increasing distance from the electrode so that coagulation is at first limited to the area surrounding the coagulation electrode. The effective contact area between the tissue and coagulation electrode has a decisive influence on the maximum achievable area of the coagulation zone.

Due to the spatial extent of the coagulation zone being dependent on the current density distribution, the technique with which the HF current is applied is also of importance. Here, a distinction is made between monopolar and bipolar surface application, monopolar and bipolar puncture application as well as monopolar and bipolar grip application.

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Coagulation modes

The spatial extent of the coagulation zones also depends on the type of emergency supply and application technique.

It is useful to differentiate between 2 coagulation modes:


- Soft coagulation
- Forced or standard coagulation

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
Soft coagulation

With soft coagulation, the generator voltage is kept so low that no arcing is able to take place between the electrode and tissue during the coagulation process. Soft coagulation is recommended for all instances of coagulation in which monopolar or bipolar coagulation electrodes are used which are brought into direct contact with the tissue that is to be coagulated. Carbonization of the tissue is avoided.

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Forced or standard coagulation

With forced or standard coagulation, electric arcs are intentionally permitted between the coagulation electrode and tissue to enable deeper coagulations to be achieved than are possible with soft coagulation. This particularly applies if only small-face electrodes are available for coagulation. To avoid the risk of a cutting effect, a modulated HF voltage must be used. Here, carbonization of the tissue will occur. A typical example of applications is TUR, in which the cutting electrodes are also used for coagulation purposes.

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
The AUTOCON® 200 unit

The operator bases his assessment of the quality of the tissue cutting on the various parameters of the entire system. Such parameters include the geometry and size of the surgical probe used, the consistency of the tissue at the site of the operation, the penetration depth and cutting speed of the surgical probe, the conductivity of the irrigant used, and the application of the neutral electrode.

These 'external' parameters, which are extremely variable, influence the properties of conventional HF-equipment to such a degree that the operator must constantly readjust the power or generally set a higher power level. The excess power, however, produces severely necrotized cut edges and exposes the patient to unnecessarily high current and the associated risks.

Fig.: Power output (P) in relation to time (t) when cutting tissue.

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The AUTOCON® 200 is equipped with automatic control circuitry with which the peak HF voltage is kept constant. This control system ensures that the type of incision is highly elastic and the quality of the incision is largely consistent, i.e. without altering the settings the operating surgeon can guide the cutting electrode the tissue at virtually any depth, whereby the coagulation depth of the cut surfaces remains largely constant irrespectively.

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The AUTOCON® 200 HF surgical unit is a universally applicable unit which permits the precise setting of various incision and coagulation qualities as well as reliable reproducibility of these settings. Due to automatic open-loop and closed-loop control of the HF voltage, the HF voltage and intensity required for achieving the desired incision quality are provided in each phase of the incision.

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The AUTOCON® 200 permits monopolar incisions with voltage stability control and – as an optional extra – monopolar incisions in the ENDOCUT setting. 3 different modes are available for coagulation, 2 monopolar and 1 bipolar. The bipolar coagulation mode additionally permits operation with the Autostart function. In conjunction with automatic open-loop and closed-loop control of incision and coagulation quality, the instrument's various safety devices provide an exceptionally high level of user and patient safety. An error documentation facility as well as software-supported test programs ensure easy and rapid servicing.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, the accompanying text should be read carefully to ensure the safe and effective operation of this product. To make the words **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

Definitions

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

Note: A Note indicates special information about operating the product, or clarifies important information.

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Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

Warning: The electrical installation of the operating room must comply with the applicable IEC requirements.

Warning: Keep out of reach of patients.

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Warning: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

Warning: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this.

Warning: Grounding reliability can only be achieved when the equipment is connected to a "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly. The unit should not be used if any damage is evident.

Warning: Risk of explosion if used in the presence of flammable anesthetics.

Warning: Always unplug the unit before carrying out any maintenance work on it (e.g. cleaning).

Warning: To avoid the risk of electrical shock, do not open the unit. Refer servicing to the manufacturer or to qualified personnel. Removal of covers by unauthorized personnel will void the equipment's warranty.

Warning: Contact with noninsulated instruments by the active high frequency electrode can lead to unintentional burns. We therefore recommend using insulated instruments in the vicinity of the active electrode.

Warning: Avoid currents flowing via tissue bridges. The resulting current concentrations at the base adjoining healthy tissue may cause burns.

Warning: Use standard-compliant endoscopes with an insulated eyepiece to reduce the risk of capacitive leakage currents flowing to the user or to the attached camera system.

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Warning: Please always observe the maximum permitted high frequency peak voltage for operation of the instruments. Failure to observe the maximum permissible high frequency peak voltage for the operation of endoscopic accessories to high frequency generators can lead to destruction of the insulation and thus to possible burns at the contact points of the instrument to the patient or user.

Warning: The arc produced at the electrode during operation of the high frequency surgical device can impair the function of other electrical devices.

Warning: During simultaneous use of high frequency surgery and physiological monitoring equipment on a patient, every monitoring electrode should be positioned as far away as possible from the active and neutral electrodes of the high frequency surgical device. Needle electrodes for monitoring are not recommended. In all cases monitoring electrodes which contain devices for limiting the high frequency current are recommended.

Warning: A double-faced neutral electrode that has short-circuited precludes monitoring.

Warning: The unit may only be operated under the environmental conditions specified.

Warning: Please ensure that there is sufficient air circulation for your unit.

Warning: If liquid enters the unit, it must be given enough time to evaporate.

Warning: Patient contact with grounded metal parts must be avoided.

Warning: In the event of a sudden power loss, check that the neutral electrode is correctly positioned before increasing the power output.

Warning: Due to low-frequency electric currents, nerves and muscles of the patient can be stimulated during HF surgery.

Warning: Patients with heart pacemakers must be treated with special care. High frequency energy can interfere with heart pacemakers. The pacemaker manufacturer's instructions must be strictly followed.

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Use of several HF units on one patient

Warning: Only use units with 'floating' application parts. Always operate each unit with its own neutral electrode.

Avoid high frequency electric circuits, i.e. ensure adequate spacing (e.g. at the top, bottom, in the center, no switching of active to neutral electrodes).

All units should be operated according to the instructions provided in the respective instruction manuals.

Injury due to incorrect preselected high frequency parameters

Warning: Always select a high frequency output power setting suitable for the particular procedure in question in order to avoid:

- thermal invasion of the tissue, which can be caused by too high a power or effect setting;
- insufficient coagulation resulting in excessive bleeding, which can be caused by too low a power setting.

Warning: Avoid explosive gas concentrations that may be present in the area where a high frequency surgical procedure is being conducted. The following situations, among others, must be taken into consideration:

Gas which may support combustion is sometimes present in the gastrointestinal tract of an unprepared patient. Certain substances used for the preparation of the lower GI tract can enhance methane production. This is particularly relevant in colonoscopy, but has also been recorded in the upper GI tract. In addition, during transurethral resection of the prostate, it has been recorded that hydrogen can accumulate in the bladder above the irrigation solution, this is particularly true when a strong arc is used.

Working with explosive gas mixtures may lead to deflagration/explosion.

Caution: Only operate the unit with the voltage stated on the unit's identification plate.

Caution: Use only fuses of the correct rating.

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Caution: Avoid allowing liquids to enter the unit. Do not store liquids on or above the unit. It must also be ensured that no liquid bottles or tubing are placed near the unit.

If liquid enters the unit, despite the precautions, sufficient time should be planned for evaporation (this also applies for the formation of condensation).

Caution: The signal input and signal output parts of this equipment are designated by the manufacturer for exclusive connection to equipment which complies with IEC 60601-1.

Caution: It is essential to follow the instruction manual closely when cleaning the unit and preparing/sterilizing unit accessories.

Note: In the case of HF generator accessories, attention must be paid to the maximum HF peak voltage allowed. Please observe the instruction leaflet or instruction manuals provided to prevent damage to the material.

Note: Inspect high frequency cord, the connecting cable to the neutral electrode and endoscopic accessories regularly for damage to the insulation.

Note: Make certain to observe national regulations/laws regarding disposal.

Note: There have been various reports in the press about possible interference to medical instruments caused by nearby mobile telephones. We therefore recommend avoiding or prohibiting the use of mobile telephones near medical devices and pointing this out by appropriate signs.

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Please read these safety instructions carefully before using the unit. Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The AUTOCON® 200 is a high frequency surgical unit for monopolar cutting as well as monopolar and bipolar coagulation with suitable instruments. Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The AUTOCON® 200 may only be used with accessories which have been designated by KARL STORZ as suitable for the unit or the safe use of which is proven. Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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User qualification

The AUTOCON® 200 may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in the use of the instrument. Such instruction should only be provided by persons equipped to do so on the basis of their skills.

Your local representative or responsible KARL STORZ member of staff is available to provide training on your premises.

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Patient profile

The use of AUTOCON® 200 is not limited to a certain patient profile (sex, age, weight etc.). The attending physician must assess a patient's health for the relevant application and declare it suitable.

Affected body part or tissue type

The field of application is not confined to certain parts of the body or types of tissue.

User profile

- Adequate and recognized skills in the relevant application
- Adequate powers of comprehension to rationally assess the surgical situation in hand
- Adequate language skills in at least one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the unit
- Knowledge of the contents of the instruction manual
- Adequate knowledge of the purpose and whereabouts of the medical products log-book
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic)

Use

The unit may only be used in hospitals and doctors' offices in the ambient conditions given in the technical data.

Other conditions of use:

- Frequency of use up to several times a day
- Length of use from a few minutes to several hours a day
- Position it on a vibration-free, fixed base
- Can be moved only by installing it on an equipment trolley
- Can be swiveled only by installing it on ceiling supply units
- Intermittent operation:
25% duty cycle (10 sec. on, 30 sec. off)

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
Safety precautions at the site of installation

The unit may only be used in medical rooms whose electrical systems have been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the demarcated hazard zone. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using it. During treatment with the AUTOCON® 200 the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

The leads to the active electrodes should be placed so that they touch neither the patient nor other leads.

Any treatment may only be performed if there is visual observation of the action of the unit.

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Important notes on the use of the monopolar technique

Risk of burns

Certain risks are fundamentally attached to HF surgery for both user and patient. Above all, contact burns can occur due to unintentional thermal damage of tissue. With a knowledge of the possible causes and by observing the safety instructions below, burns can, however, be safely avoided.

Burns due to incorrect application of the neutral electrode

With incorrect or no application of the neutral electrode, there is a high risk of thermal tissue damage both at the point of application of the neutral electrode and at other points of the patient's body. When applying neutral electrodes, it is therefore imperative that care be taken to ensure that the effective contact area or conductance between the neutral electrode and the patient is adequate for the HF power used or intensity of the HF current. Here, effective contact area refers to the area of the neutral electrode in electroconductive contact with the patient's skin during HF surgery.

When selecting the neutral electrode take the following empirical values into account:

Neutral electrode made of electroconductive silicone:	1.5 watts/cm ² or 10 mA _{rms} /cm ²
Neutral electrode made of metal without conductive jelly:	3.0 watts/cm ² or 15 mA _{rms} /cm ²
Disposable neutral electrode:	5.0 watts/cm ² or 20 mA _{rms} /cm ²

Warning: Care must be taken to position the neutral electrode near the operating field on shaven tissue with good blood circulation.

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
The neutral electrode is secured to either the upper arm or thigh, depending on which is nearest the operating field. To ensure that the entire face of the neutral electrode is applied, it must be securely fastened with rubber straps and then wrapped with a bandage so that it cannot slip, even if the patient should move (note the risk of decubital necrosis). The neutral electrode must always be in perfect condition. Grease, oxidation, dirt etc. on the contact surface reduces the size of the electroconductive area and can lead to dangerously high current densities.

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Warning: A short circuit fault is detected when using a double-faced NE. Power is disconnected until optimal conditions return. Subsequently, a renewed HF activation is required.

Warning: The attachment of a neutral electrode made of conductive silicone (27805) prevents monitoring.

Note: The HF current flows to the neutral electrode through the tissue. It should be remembered that the distribution of current at the proximal corners and edges to which the current flows may be higher than at the distal corners and edges. For this reason, when arranging the neutral electrode care must be taken to ensure that the longer side of the neutral electrode is facing the operating site.

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
Patients must be correctly positioned before performing any operations using HF surgical instruments. An electrically insulating, thick and dry pad must be placed between the patient, operating table and clamps. If moisture, perspiration etc. are anticipated during the operation, penetration of moisture through this pad must be prevented by waterproof sheeting or other suitable measures. The patient's body tissue must not be allowed to come into contact with nearby electrically conductive objects (examples: operating table, metal clamps, wet towels, etc.).

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Since contact areas will generally be small in such cases, burn-necroses could arise from high frequency leakage currents due to high local current densities. In cases where contact between the neutral electrode and patient is poor this effect is increased.

Areas subject to more severe perspiration, extremities resting against the body or skin-to-skin contact should be dried by placing cover cloths in-between. Urine should be drained off by a catheter. If an injection cannula is used as an ECG electrode, the metal cone must not rest on the skin; the same applies to the leads of monitoring equipment.

- A Electrically insulated pad
B Grounded operating table

Warning: Skin-to-skin contact should be avoided.

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Burns due to unintentional activation of the HF generator

Patients can also suffer burns due to unintentional activation of the HF generator. Unintentional activation of the active electrode can occur, for example, if the footswitch is inadvertently depressed, if there is a short circuit in one of the cables, if electrically conductive liquids enter the handpiece or the footswitch, or if the HF generator is defective.

Warning: To prevent such occurrences:

- never leave electrodes lying on or near the patient in such a way that could result in direct or indirect (through the mediation of electrically conductive objects) contact between electrode and patient,
- ensure that the acoustic warning signal indicating the operational status of the HF unit is always set loud enough to be easily heard,
- ensure that the visual warning signal can always be seen by all attending personnel.

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Burns from contact with hot electrodes

Cutting and coagulation electrodes can become very hot during cutting or coagulation procedures. Burns can result from active or inactive electrodes coming into contact with tissue immediately following cutting or coagulation procedures.

Warning: Electrodes which are not currently in use must not come into contact with the patient.

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Burns due to improper use of HF power output

The optimal HF power level setting is that which yields the best results in terms of the desired effect (cutting/coagulation) with minimal incidence of undesirable side effects (i.e. charring of tissue). Since undesirable side effects occur almost exclusively at excessively high power levels, HF power levels should always be kept as low as possible. If the effects desired suddenly no longer be attainable at power levels experience has shown to be optimal, always check to be certain that everything is in proper working order and properly connected before increasing the HF power level:

- Is the neutral electrode properly connected?
- Is the active electrode properly connected?
- Are active and neutral electrodes clean?
- Is the HF generator power level correctly set?

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Burns arising from ignition of combustible media

Skin cleansing agents, disinfectants, and anesthetic gases, as well as endogenic gases, are generally readily combustible and can even explode. Since sparking inevitably occurs in HF surgery, vapors or gases evolved by such media could be ignited and result in burn injuries to patients. Thus, before switching on the HF generator, make certain that all combustible media and all objects that have been in contact with them (tissues, swabs, etc.) have been removed from the operating area.

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Warning: With TUR, H_2O molecules may be split into H_2 and O_2 by the arc between the resection loop and flushing liquid. These gases may collect at the top of the bladder in the form of a highly explosive gas mixture. Resection in this gas mixture can lead to dangerous explosions.

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Burns due to unsuitable application

If, during monopolar application, HF current flows through tissue constrictions, it must be remembered that a coagulation effect can occur at the constricted point. If tissue sections have the same cross section over an extended dis-tance, a uniform coagulation effect can occur over the entire distance.

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Stimulation of nerves and muscles

Due to low-frequency electric currents, nerves and muscles of the patient can be stimulated during HF surgery. The arcs between the elec-trode and tissue, which are unavoidable for cut-ting processes, spray coagulation and forced coagulation, have the effect of rectifying part of the HF alternating current, whereby more or less heavily modulated, low-frequency current com-ponents occur which are able to cause nerve stimulation. When using HF surgery on electri-cally stimuable structures, contractions of the muscles concerned must therefore be reckoned with.

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Heart pacemakers

Irreparable damage to heart pacemakers and interference with pacemaker operation, which could lead to ventricular fibrillation, are likely to result if radiofrequency surgery methods are used on patients with implanted heart pacemakers or pacemaker electrodes.

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Implants

In the case of patients with electrically conductive metal implants it must be ensured that with all applications involving HF surgery the current flow or current path of the HF electric current does not pass through or over the implant. At the current entry and exit points into and from the bio-tissue to the implant, current concentration may cause irreparable internal burns.

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Safety devices**Protective circuit in case of misdosage caused by the unit**

Deviations from the setpoint of each of the HF output parameters set are recognized by the unit and signalled to the user. If the deviation is so large that the required quality of incision or coagulation is no longer guaranteed, the HF generator is switched off.

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Monitoring the ON time

When used on the patient, the HF generator is only operated briefly when used as prescribed. The ON time of the HF generator is therefore constantly monitored. If a predetermined ON time (20 sec.) is exceeded, this is indicated by a visual signal on safety panel ②. If the HF generator is operated further, an acoustic warning signal is additionally sounded and the generator switched off.

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This protective circuit primarily serves the purpose of avoiding damage due to unintentional activation of the electrodes. The HF generator can be restarted at any time after being stopped by the protective circuit.

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Monitoring the neutral electrode

The AUTOCON® 200 has a protective circuit which monitors the electrical connection between the neutral electrode and unit and, when using double-faced electrodes, correct application on the patient.

If a single-faced neutral electrode is used, the indicator lamp ⑭ lights up green when the unit and the neutral electrode are correctly connected.

If the connection is interrupted, the LED goes out and all monopolar modes are blocked. If an attempt is made to activate the HF generator in a monopolar operating mode, the LED lights up red and an acoustic warning signal is given. At the same time, the corresponding warning signal lights up in safety panel ②, and an error message is shown on displays ④ and ⑧.

When using double-faced neutral electrodes, the electrical conductance between the two contact areas of the neutral electrode and the patient's skin is measured and compared with the intensity of the HF current flowing through the neutral electrode. If the intensity of the HF current is greater than it should be according to the conductance measured, the user is warned by visual and acoustic signals. If the conductance measured is too low, the monopolar operating modes are blocked.

With double-faced neutral electrodes, the application direction of the contact area is also monitored relative to the direction of the current. Here the intensity of the two partial currents which flow through the two partial faces of the neutral electrode is compared. Should the partial currents differ, the warning signal (exclamation mark) lights up on safety panel ②. If the partial currents differ too greatly, the warning signal flashes, a simultaneous acoustic warning sound is given, and the HF generator is switched off.

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Self-test program

Each time the instrument is switched on, a self-test program is run which recognizes and signals errors in the instrument's operator controls as well as errors in the accessories connected (hand grip, footswitch). This takes place by an audible signal and visual display with error code on displays ④ and ⑧.

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Unpacking the equipment

Carefully unpack the AUTOCON® 200 and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately. If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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
Basic equipment

- 1 AUTOCON® 200
- 1 Adaptor, 5/4 mm
- 1 Power cord
- 1 Instruction manual
- 1 Cover for instruction manuals
- 1 Operator's quick reference guide (self-adhesive)

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Installing and connecting up the unit

Note: The AUTOCON® 200 HF surgical unit including the accessories connected may be used only in medical facilities whose electrical systems have been installed in accordance with applicable national regulations. This unit is not intended for use in hazardous zones. Do not operate the unit within demarcated hazard zones while explosive anesthetic gases are in use.

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Place the unit on a level surface.

Note: HF surgery equipment generates high frequency and electromagnetic fields that can interfere with the operation of sensitive electronic equipment. HF surgical equipment, and particularly its cabling, should thus never be installed close to any equipment sensitive to electromagnetic interference.

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The unit is equipped with a connector ③ for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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Operate the unit only with the system voltage stated on the rating plate.

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
Connect power cord. Insert power cord into power cord receptacle ④ as far as it will go.

Note: The unit should only be operated with the power cord delivered by the manufacturer or a power cord of similar quality, which has been provided with a national inspection seal.

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! Warning: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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Insert the connection of the footswitch in socket ② and tighten the union nut.

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**Connecting the neutral electrode
(only with monopolar applications)**

Warning: It is imperative that the directions given on page 18 ff be observed before applying the neutral electrode.


Connect the cord of the neutral electrode to socket ⑪.

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Both single-faced and double-faced neutral electrodes can be connected up to the AUTOCON® 200. If connected correctly, the indicator next to socket ⑪ has a green light; if the indicator is not lit, the neutral electrode has not been connected properly.

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Note: When connecting double-faced neutral electrodes, the electrical connection between the neutral electrode and instrument as well as correct application to the patient is checked by the unit.

When connecting single-faced neutral electrodes, only the connection between the neutral electrode and unit can be checked.

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Connecting up the active electrode

Two different sockets are available for connecting active electrodes: Socket ⑫ is for connecting monopolar electrodes. Electrodes connected to this socket can be activated both via the hand grip and footswitch.

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Socket ⑬ is for connecting bipolar coagulation electrodes. Electrodes connected to this socket can only be activated via the footswitch.

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⚠️ Warning: Two electrodes can be simultaneously connected to sockets ⑫ and ⑬ and activated. In this case, care must be taken to ensure that the electrodes are put down in such a way as to prevent any confusion. If in doubt, test activation of the electrodes must precede introduction into the patient.

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Operating the unit

Turn on power switch ①. All symbols and display fields briefly light up on the front panel. The basic setting of the front panel is then displayed flashing.

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Softkoagulation:

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Forcierte Koagulation:

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The AUTOCON® 200 is delivered ex works with the following basic setting:

Monopolar cutting, voltage-controlled:

Power limitation 120 W,
coagulation effect stage 3.

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This basic setting is displayed each time the power switch is turned on, unless the user has programmed a different basic setting. If a basic setting programmed by the user becomes invalid, e. g. due to failure of the memory circuit or due to battery power loss, this basic setting is again adopted. To change the basic setting, see “Programming the basic setting of the front panel” on page 37.

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Cutting mode


Preselect monopolar operating mode*. To do this, keep the M button (M) pressed until the symbol for the required cutting mode lights up:

- Monopolar incision with voltage stability control
- ENDOCUT operation (*only the 205225 20-010 is equipped with this add-in).

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If necessary, set an HF power output different to the basic setting indicated on the front panel using +/- buttons (±). Limitation of the power output can be set between 1 and 200 watts. The power output selected is indicated on display (4).

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If necessary, set a coagulation effect for the incision edges different to the basic setting. To do so, preselect one of stages 1 (minimum coagulation effect) to 4 (maximum coagulation effect) by pressing button ⑦. The stage selected is illuminated on the display panel ④.

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The monopolar cutting electrodes can now be activated by pressing the yellow footswitch pedal or yellow release on the hand grip.

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
As soon as the electrodes are activated, the corresponding activation symbol ⑤ lights up. In addition, an acoustic signal is sounded during the activation phase. The volume of the acoustic signals can be adjusted with volume control ① on the rear of the unit.

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Coagulation mode

Preselect coagulation mode by keeping M button ⑩ pressed until the symbol for the coagulation mode required lights up on display panel ⑦:

- soft coagulation
- forced coagulation
- bipolar coagulation

Please refer to the section entitled “Coagulation modes” on page 6 ff for more detailed information regarding the individual coagulation modes.

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If necessary, set a power output different to the basic setting on the front panel using the ± buttons (15). Limitation of the power output can be set between 1 and 120 watts and is indicated on display (8).

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Caution: When using pointed bipolar coagulation forceps, the tips can suffer thermal damage due to electrical currents being too high. It is therefore recommended that the power limitation be set as low as possible and/or that care be taken to ensure that the tips of the coagulation electrodes do not touch each other.

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Autostart function

In “Bipolar coagulation” mode, the Autostart function enables automatic activation of the coagulation current as soon as both branches of the coagulation electrode touch the tissue.

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The Autostart function allows activation of coagulation with delay between 0.1 and 1 second. A value of 0.5 seconds is preset ex works. The ON delay times can be reset by the user according to his own judgment. Please refer to page 40.

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To activate the Autostart function, keep button ⑭ pressed until the symbol for the Autostart function lights up.

Warning: With the Autostart function switched on, it must be ensured that the coagulation electrodes are not put down on the patient or on an electroconductive surface.

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If the Autostart function is no longer required, it should be switched off for safety reasons.
If the foot or handswitch is operated with the Autostart function switched on, the autostart function is deactivated.

Note: The Autostart function can be activated only in bipolar coagulation mode. If button ⑭ is pressed in a monopolar coagulation mode, an intermittent acoustic error signal is sounded.

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- unipolare Koagulation
- bipolare Koagulation.

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The coagulation modes can now be activated with the blue footswitch or blue button on the handgrip.

As soon as the electrodes are activated, the corresponding activation symbol ⑨ lights up

- monopolar coagulation
- bipolar coagulation.

In addition, an acoustic signal is sounded during the activation phase.

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The volume of the acoustic signal can be adjusted with volume control ① on the rear of the unit.

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**Test for proper operation
Monopolar operating modes**

Connect the neutral electrode to socket ⑪, a monopolar active electrode to socket ⑫.

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Select a monopolar cutting mode on M-button ⑲ or a monopolar coagulation mode on M-button ⑳ and set the power limitation to maximum with ± buttons ⑱ or ⑲.

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Place a wet sponge, swab or suchlike onto the neutral electrode and bring into contact with the active electrode. Activate the HF generator by pressing the yellow foot pedal or the yellow button on the hand grip.

Water must evaporate on the active electrode. The corresponding “Monopolar incision” or “Monopolar coagulation” activation symbol must light up and an acoustic signal be sounded.

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Bipolar operating mode


Connect the bipolar instrument to socket ⑬ and bring into contact with a wet sponge, swab or suchlike.

Select bipolar coagulation mode with M-button ⑱ and set the power limitation to maximum with ± buttons ⑲.

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Activate the HF generator by pressing the blue footpedal. If the generator supplies power, water will evaporate at the point of contact. The corresponding “Bipolar coagulation” activation symbol must light up and an acoustic signal be sounded.

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After completing the functional test, reset the power limitation to the previous value.

Note: Should the results differ from those described above, the instrument may be put into operation only after inspection by an authorized technician.

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Program functions

The AUTOCON® 200 has several program functions which allow the user to carry out setting, testing or interrogation. The respective program can be called up by entering the corresponding program number after summoning the program mode.

The following programs can be called up:


- Program 1:
Programming the basic setting of the front panel
- Program 2:
Retrieving the stored error signals
- Program 3:
Setting the intensity of the seven-segment display
- Program 4:
Front panel display test
- Program 5:
Display of the software version number
- Program 6:
Programming the Autostart ON delay time
- Program 7:
Monitoring the neutral electrode
- Program 8:
Setting the no-load voltage for the "Forced coagulation" operating mode
- Program 9:
Setting the ON time limitation

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Calling up the program mode

While switched off, press M-button (19) in the "Cutting" function field and simultaneously turn on power switch (1).

"Pr." (program) appears on "Cutting" display (4) and program number "1" on "Coagulate" display (8).

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
Ending program mode

Program mode can be ended by turning off the power switch or by setting test program number "0".

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**Program 1:
Programming the basic setting of the front panel**

With this program the user can himself define and set the basic setting of the front panel. The basic setting is displayed flashing each time the power switch is turned on.

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After calling up program mode with ± buttons ⑮ in the coagulation field, select program number "1". Then press M-button ⑰. The preset value for the "Monopolar incision" operating mode lights up on the "Cutting mode" display panel ③.

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Enter the power limitation required with ± buttons ⑱. The value entered is indicated on display ④. Enter the required coagulation effect of the incision edges with button ⑰. The coagulation effect set is indicated by the corresponding symbol on display panel ⑥.

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Press M-button (19) again. The symbol for the ENDOCUT operating mode lights up (only if the unit is equipped with this additional function).
Enter the values required for power limitation and coagulation effect as described above.

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Now enter the basic settings required for the various coagulation modes. To do this, select individual coagulation modes with M-button (16).
Enter the respective power limitation with ± buttons (15). If necessary, activate the autostart function in the bipolar coagulation mode with button (14).

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Once you have entered the parameters required for each coagulation mode, turn the instrument off again at power switch (1). The basic setting is now stored and will be automatically set each time the instrument is put into operation.

Note: In the event of memory loss, e.g. due to an extended power failure, the basic setting entered ex works will again be adopted.

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**Program 2:
Retrieving the error signals stored**

The AUTOCON® 200 is equipped with an automatic error documentation facility. Should an error occur while switching on or during operation, the corresponding error number will be stored in the instrument. The last ten error numbers can be retrieved at any time via program 2.

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To do so, enter program number 2 with ± buttons 15 after calling up program mode.

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
Press M-button 19. "E" (error) and the number of the memory location appear on "Cutting" display 4. The error stored last will be found at memory location 1. The error number appears on "Coagulation" display 8. The error can be identified in the error list on page 55 ff of these instructions by means of the error number. The errors stored can be retrieved in sequence by pressing ± buttons 15. The errors stored can be erased by pressing M-button 19.

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
**Program 3:
Setting the intensity of the seven-segment
display**

With this program the light intensity of the seven-segment display can be changed. To do so, enter program number 3 with ± buttons 15 after calling up program mode.

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Press M-button 19.
"Int." (intensity) appears on "Cutting" display 4 and a number for the degree of intensity set (1 = minimum intensity, 10 = maximum intensity) on "Coagulation" display 8.
Set the degree of intensity required with ± buttons 15.

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**Program 4:
Front panel display test**

With this program the correct functioning of all visual signals on the front panel can be checked. The seven-segment displays each show 888.

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
To do so, enter program number 4 with ± buttons (15) after calling up program mode. Press M-button (19). All visual signals are switched on.

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**Program 5:
Display of the software version number**

With this program the number of the software version installed in the instrument can be retrieved.

To do so, enter program number 5 with ± buttons (15) after calling up program mode and press M-button (19).

“Snr” (software version number) appears on the “Cutting” display (4) and the number of the software installed on “Coagulation” display (8).

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**Program 6:
Programming the Autostart ON delay times**

With this program the delay time programmed by the manufacturer in Autostart mode can be redefined.

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To do so, enter program number 6 with ± buttons (15) after calling up program mode. Press M-button (19). "set" appears on "Cutting" display (4) and the delay time for Autostart on the "coagulation" display (8).

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Now enter the ON delay required for Autostart between 0 and 10 seconds, using the ± buttons (15). Press M-button (19) to store the value entered.

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
**Program 7:
Monitoring the neutral electrode**
With this program the instrument can be adjusted to use various types of neutral electrodes.
To do so, enter program number 7 with the ± buttons (15) after calling up the program mode.

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Press M-button (19). "Ne." and Setting 1 (standard setting) appear on the "Cutting" display, and a corresponding graphic representation appears on the "Coagulation" display.
Select the desired setting with the ± buttons (15).

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The following settings are possible:

- 1: Standard version, for single-faced and double-faced neutral electrodes (II I)
- 2: Only single-faced neutral electrodes (I)
- 3: Only double-faced neutral electrodes (II)
- 4: Only double-faced neutral electrodes, continuous visual and acoustic current density alarm, reduced alarm time constants (II)

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**Program 8:
Setting the no-load voltage in the “Forced Coagulation” operating mode**

With this program the no-load voltage in the “Forced Coagulation” operating mode can be changed in accordance with the adjusted power limitation.

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To do so, enter program number 8 with the ± buttons (15) after calling up program mode.

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Press M button (19).
“For” (Forced) appears on the “Cutting” display, and the adjusted voltage version appears on the “Coagulation” display.

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Select the desired voltage version (1– 3) with the ± buttons **15**.

Version 1:

From 1 W to 30 W power limitation, the no-load voltage continuously increases up to 1300 Vp. Above the 30 W power limitation, the no-load voltage is limited to 1300 Vp.

Characteristics: Minimal image distortion during video operation, minimal danger of burns in case of contact with the clamp, good coagulation without major generation of sparks even at a high power setting. In order to move the energy into the tissue, contact must be made with the tissue.

Version 2:

From 1 W to 30 W power limitation, the no-load voltage constantly increases up to 1300 Vp, and the no-load voltage continues to increase above the 30 W power limitation.

Characteristics: The characteristics are between versions 1 and 3.

Version 3:

From 1 W to 30 W power limitation, the no-load voltage continuously increases up to approx. 2300 Vp. Above the 30 W power limitation, the no-load voltage is limited to approx. 2300 Vp.

Characteristics: Possible image distortion during video operation, danger of burns in case of contact with the clamp, heavy generation of sparks during coagulation, even at the 30 W power setting. In order to move the energy into the tissue, there need not necessarily be contact with the tissue.

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**Program 9:
Setting the ON time limitation**

With this program the ON time limitation can be changed in the cutting and coagulating modes. To do so, enter program number 9 with the ± buttons (15) after calling up program mode.

Press M-button (19).
The set ON time limitation for the Cutting operating mode appears on the “Cutting” display. The value for the ON time limitation for Coagulation operating mode appears on the “Coagulation” display (each indicated in seconds).

The ON time limitation can be reset with the ± buttons (15) or (16) in a range between 3 and 900 seconds for any operating mode.

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Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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Remove the power fuse holder ⑤ with a screwdriver or other tool.

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! Caution: Only use fuses of the correct rating.

Insert new fuses of the appropriate rating.

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Refit the power fuse holder ⑤.
Reconnect the power cord.
Test the instrument for proper operation.

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Cleaning, disinfection and sterilization

! Warning: Always pull out power plug before cleaning!
Caution: Avoid allowing fluids to enter the unit.

The exterior of the unit and footswitch can be surface disinfected with a disposable cloth soaked in disinfectant (preferably alcohol-free).

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! Caution: After immersing the cable in solution carefully dry it out in order to avoid any malfunctions due to residual moisture at the mating contacts. After cleaning, steam sterilization should preferably be performed at 134 °C by the fractionated prevacuum method.

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Manual preparation

Note: Use a cleaning/disinfecting agent approved by KARL STORZ for cleaning/disinfection (see Appendix, page 67). You will also find the current list of approved products on our Web site (www.karlstorz.com).

! Caution: When preparing and using the solutions, follow the manufacturer's instructions paying close attention to proper dilution and working times. Prolonged immersion can lead to material changes.

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Machine preparation

Due to the partially strong motion of the water, it is necessary that all parts be fixated in holders which prevent damage. For sufficient rinsing out of the instruments, it is necessary that instrument channels are correctly connected.

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Sterilization

Note: Before sterilization, the components must be cleaned (see table, page 46):

- mechanically (manual) and/or
- washer/disinfector (thermal disinfection at 93°C)

Then pack the components into appropriate containers.

Warning: Sterilization is only possible on clean surfaces.

Caution: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and calibrated.

Caution: Any deviations from the recommended parameters for sterilization should be validated by the user.

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Steam Sterilization

Sterilization trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage.

Load the sterilizer according to its instruction manual.

At completion of the steam sterilization cycle, the instruments must cool down slowly.

Warning: Burn hazard! Parts are hot after steam sterilization. Allow them to cool down.

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Validated steam sterilization methods

KARL STORZ has validated the following steam sterilization methods:

- Prevacuum
- Fractionated prevacuum
- Gravity displacement

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Vacuum-/Prevacuum

The prevacuum method consists of four phases.

- **Conditioning phase**
The conditioning phase removes air from the chamber by pulling a vacuum and then warms the instruments by injecting steam.
- **Sterilization phase**
Sterilization takes place at a temperature of 134 °C and a pressure (absolute) of 3 bar over a period of 5-8 minutes.
- **Evacuation phase**
The steam is released.
- **Drying phase**
Drying takes place under renewed application of a vacuum over a period of approx. 5-20 minutes.

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Fractionated prevacuum procedure

The fractionated prevacuum procedure also consists of the four previously mentioned phases. In the conditioning phase, a vacuum is generated up to four times. Steam is then blown in and sterilization takes place at 121 °C or 134 °C in 20 or 5 minutes respectively (minimum times). Before the final drying phase, the steam is removed from the chamber during the evacuation phase. Drying takes place under vacuum.

Validated parameters for the fractionated prevacuum method:


Temperature:	134 °C (± 3 °C)
Pressure (p _{absolut}):	3 bar
Exposure time:	5 minutes

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Gravitation method

Warning: This method is only suitable for smooth and accessible surfaces, not for instruments with channels, gaps or stopcocks.

The gravitation method is also made up of 4 phases:

- **Conditioning phase**
Steam is pressed into the chamber and heats the instruments.
- **Sterilization phase**
Sterilization takes place at a temperature of 121°C and a pressure (absolute) of 2 bar over a period of 30 - 45 minutes. Sterilization may also take place at 134 °C and 3 bar, however the time required is 5 - 10 minutes.
- **Evacuation phase**
Steam is released through the outlet.
- **Drying phase**
Drying takes place under ambient air pressure by heating through the heating jacket of the autoclave.

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
Ethylene oxide (EtO) gas sterilization

All sterilization cycles must include a pre-conditioning cycle:

Temperature: 54 ± 2°C
Time: 30 minutes

KARL STORZ has validated ethylene oxide (EtO) sterilization using the following parameters:

Gas Mixture	EtO: HCFC*
Temperature	54 ± 2°C
Relative Humidity	60 ± 20%
Pressure	0.56–0.70 bar
Exposure Time	120 minutes
EtO concentration	600 ± 30mg/l

* Gas mixture composed of 10% by weight ethylene oxide and 90% by weight chlorotetrafluoroethane (HCFC-124)

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The airing time is 12 hours at a temperature of 50-55°C. This increases to 24 hours at room temperature (20°C).

Maximum permissible levels (in accordance with 21 CFR 812.100; US Federal Register, Vol. 43, No. 122, June 23, 1978) for residues following EtO sterilization are as follows:

Ethylene oxide: 25 ppm* / 250 ppm**

Ethylene chlorohydrine: 25 ppm* / 250 ppm**

Ethylene glycol: 250 ppm* / 5000 ppm**

** for instruments in contact with skin and mucous membrane,
* for instruments in contact with blood.

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Gas sterilization with formaldehyde (FO)

Caution: During formaldehyde sterilization the temperature must not exceed 65°C.

When using formaldehyde the instruments can be used immediately after the sterilization cycle has been completed. An additional airing time is not necessary.

Note: Observe national laws and regulations for the gas sterilization methods.

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Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.

Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.

Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal, Vol. 62, No. 1, pp. 23-29, 1995.

References for Cleaning and Sterilization

Gruendemann, B.J. and Meeker, M.H. Alexander's Care of the Patient in Surgery, 7th edition. The C.V. Mosby Company, St. Louis, Mo. 1983.

Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR No. 12-1994.


Descoteaux, J-G, Poulin, E.C., Julein, M. and Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal Vol. 62, No. 1, pp. 23 - 29, 1995.

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Association for the Advancement of Medical Instrumentation. Designing, Testing and Labeling Reusable Medical Devices for Reprocessing in Health Care Facilities: A Guide for Device Manufacturers. AAMI TIR Nº 12-1994.

Descoteaux, J-G, Poulin, E.C., Julein, M. y Guidoin, R. Residual Organic Debris on Processed Surgical Instruments. AORN Journal Vol. 62, Nº 1, págs. 23-29, 1995.

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Maintenance

Regular maintenance can contribute to identifying potential problems before they become serious, enhancing the reliability of the unit and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer.

Maintenance servicing (functional and safety test) is recommended at least once a year, regardless of the accident prevention regulations or testing intervals for medical devices applicable in the various countries.

For detailed instructions please refer to the latest version of the Service Manual.

Safety checks and inspections

Regular safety checks are recommended to prevent accidents due to aging, wear or failure of medical devices.

The high-frequency surgical unit must be inspected at least once a year. Please observe the accident prevention regulations and testing intervals applicable in the various countries.

Warning: Servicing may only be performed by specialists authorized by KARL STORZ.

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If any defects or shortcomings which could endanger patients, operating personnel or others are ascertained during these safety checks, the unit must not be operated until such time as these defects or shortcomings have been eliminated through qualified technical servicing.

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In accordance with the Medical Device Directive (MDD) this instruction manual contains information on the type, extent and frequency of safety checks.

When carrying out safety checks, the following items must be tested:

- Visual inspection of the unit and accessories for proper condition
- Checks of protective grounding and insulation, and for leakage currents
- Testing for proper operation of all visual and acoustic signals
- Testing of the monitoring devices
- Testing of the function monitor
- Testing of the automatic start mode
- Checking the activation of the HF power using the dual-pedal footswitch
- Testing the HF power in various operating modes
- Measurement of power consumption when operating at no-load and under load

Note: For detailed information on the range and implementation of the safety checks please refer to the Service Manual SV 3349 (English version).

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen.

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and used meet the applicable laws and standards, and
- the unit has been used in accordance with its operating instructions at all times.

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Manufacturer's warranty

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen, Germany

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 Tuttlingen, Germany

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Warning: Always unplug the unit beforehand.

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Troubleshooting

- Error description:*
When activating CUT or COAG, another activation signal is recognized by the unit (only for Spanish version). The national version is set using jumpers on the front panel.
- Remedy:*
Do not double-activate.

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- Error description:*
Resistance load on unit was too low.
- Remedy:*
Either the contact area of the active electrode was too large or the power limitation set too low or the active electrode was in contact with a metallic instrument, e.g. trocar sheath.

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- Error description:*
The HF generator of the unit was overloaded for too long. The maximum power output has been automatically reduced.
- Remedy:*
The unit can supply more than 200 watts of power for a short period of time. More than 300 watts, averaged over 1 second, are not permitted for safety reasons.

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Note: This troubleshooting list only includes error signals which are relevant for the user and which can be rectified by him.

Error numbers not included in this list require the instrument to be inspected by a qualified service technician.

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
Diagrams
AUTO CUT mode

Fig.: Peak value of the HF output voltage U_{HF} depending on adjustable power limit P_{max} and Effect as a parameter.

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Soft coagulation
Bipolar coagulation

Fig.: Peak value of the HF output voltage U_{HF} depending on adjustable power limit P_{HF}

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Standard compliance (for 205225 20)

According to IEC 60601-1 and IEC 60601-2-2:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type CF 

Information on electromagnetic compatibility

HF units are assigned to CISPR 11, Group 1 since, as per IEC 60601-2-2, they can only be measured when the output switch is not activated. The AUTOCON® 200 complies with the EMC guidelines of IEC 60601-1-2 : 1993, as stated in IEC 60601-2-2 : 1998.

Directive compliance (for 205225 20)

According to Medical Device Directive (MDD):

This medical product belongs to Class II b
This medical product bears the CE mark in accordance with MDD 93/42/EEC.
A code number after the CE mark indicates the responsible notified body.

Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the instrument that have been designated as repairable by their respective manufacturers.

Supply of such technical documentation relating to the instrument shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the instrument.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)
- (A)= Only for the flexible endoscopes of the new generation (year of manufacture > 2005)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets
Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Plus Sekusept Pulver Classic Ultrasemit AF	1, 2, 3 1, 2, 3 1, 2, 3
Ester Pharma GmbH	Ultrasseptin Aktiv Ultrasseptin Classic	1, 2, 3 1, 2, 3
	Afid	1, 2
	Afid plus	2
Fresenius AG	Sporcid FF	1, 2
Holifa Polska Sp. zo.o Indeba Industria E Commercio LTDA	Poisept Holifa Letahdeido	1, 2 2
	Darodor 4000 Liquido	1
Jose Collado S.A.	Darodor 3000 Darodor Sinaldehyd 2000 Aldasan 2000	2 2 1, 2
	Almyrol	1, 2
	Desoform	1, 2, 3
Lysoform Dr. Rosemann GmbH	Lysoformin 3000	1, 2
	MediDis	2
	MediZyme	1
	Edisonite Super	1, 3
Merz Hygiene GmbH	Mucadont-IS Mucadont-Zymaktiv Mucocit-T neu	2, 3 1, 3 1, 2, 3
METASYS Medizintechnik GmbH (formerly LONZA)	ID 50	2
Orochemie	A 20	1, 3, A
Promagent AB	Wavacide Descoton Extra Descoton forte	2 1, 2 2
Schuhmacher, Dr. GmbH	Perfektan Endo Perfektan Neu Perfektan TB	2, 3 1, 2 1, 2, 3, A
	Gigasept Gigasept AF forte	2, 3 1, 2
	Gigasept FF Gigasept Instru AF Gigasept Med forte	1, 2 1, 2, 3, A 1, 2
	Gigazyme Lyselvol	1, 3 2
Steris	EnzyCare 2 neodisher LM 2 neodisher mediclean	1, 3 1, 3 1
	neodisher mediclean forte	1, 3
Weigert, Dr. GmbH & Co.	neodisher medizym neodisher SeptoClean neodisher SeptoPreClean	1, 3 1, A 1, 2
	Aidal Plus	1, 2
Whiteley Industries PTY.Ltd	Matrix	1

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For flexible endoscopes, only immerse the examination sheath, not the housing/handle part.

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schülke & Mayr GmbH	Thermosept alka clean	1
	Thermosept DK	2
	Thermosept ED	2
	Thermosept ER	1
	Thermosept RKN-zym	1
Steris	Instru-Klenz	1, 3
	neodisher DuoClean	1, A
Weigert, Dr. GmbH & Co.	neodisher FA	1
	neodisher FA forte	1, 3
	neodisher Oxivario	1
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher SeptoClean	1, 2, A
	neodisher Septo DA	2
	neodisher Septo DN	2

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GARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

To be completed by the supplier/importer:
Company stamp/signature:

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Tipo del aparato: Número de serie:

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Firma/Fecha:

To be filled out by instrument owner:
Return address/company stamp:

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KARL STORZ Endoscopy South Africa, (Pty) Ltd.
8th Floor Convention Tower
Cnr Heerengracht & Coen Steytler, Foreshore
Cape Town 8001, South Africa
P.O.Box 6061, Roggebaai Cape Town,
8012, South Africa
Phone: +27 21 417 2600, Fax: +27 21 421 5103
E-Mail: info@karlstorz.co.za

TOO KARL STORZ ENDOSCOPY Kasachstan
Khodjanova 17
050060 Almaty, Kazakhstan
Phone/Fax: +7 72 72 49 43 63, +7 72 72 49 41 00
e-mail: KarlStorzKazakhstan@itte.kz

KARL STORZ Endoscopy India Private Ltd.
C-126, Okhla Industrial Area, Phase-1
New Delhi 110020, India
Phone: +91 11 26815445-51,
Fax: +91 11 268129 86
E-Mail: karlstorz@vsnl.com

KARL STORZ GmbH & Co. KG
Resident Representative Office
80/33 (44/19) Dang Van Ngu,
F:10-Q. Phu Nhuan
Ho Chi Minh City, Vietnam
Phone: +848 991 8442, Fax: +848 844 0320

KARL STORZ Endoscopy (S) Pte. Ltd.
#05 – 08 San Centre, 171 Chin Swee Road
Singapore 0316, Singapore
Phone: +65 6376-1066, Fax: +65 6376-1068
E-Mail: karlstorz@pacific.net.sg

KARL STORZ Endoscopy Australia Pty. Ltd.
174 Parramatta Rd.
Camperdown, N.S.W. 2050, Australia
Phone: +61 02 8594 9150, Fax: +61 02 8594 9199
E-Mail: info@karlstorz.au

KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay,
Hong Kong, People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020, People's Republic of China
Phone: +86 10 8525 3725, Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030, People's Republic of China
Phone: +86 21 6113-1188, Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016, People's Rep. of China
Phone: +86 28 8620-0175, Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281, Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Asia Marketing Pte Ltd
3791 Jalan Bukit Merah
06-11 e-Centre @ Redhill
Singapore 159471, Singapore
Tel. No. +65 63761066, Fax. No. +65 63761068
Email : infoasia@karlstorz.com.sg
Email : serviceasia@karlstorz.com.sg

KARL STORZ Endoscopy Singapore Sales Pte Ltd
3791 Jalan Bukit Merah
10-18 e-Centre @ Redhill
Singapore 159471
Tel. No. +65 65325548, Fax No. +65 65323832
Email : infosgpp@karlstorz.com.sg
Email : servicesgpp@karlstorz.com.sg

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6F1, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966, Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

STORZ
KARL STORZ—ENDOSKOPE

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com