STORZ

KARL STORZ — ENDOSKOPE

263310 20-1

ENDOMAT® n. HAMOU®

263310 20-1

HAMOU® ENDOMAT®

263310 20-1

ENDOMAT® según HAMOU®

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Instruction manual

Manual de instrucciones

	263310 20-1	ENDOMAT® n. HAMOU®
	263310 20-1	HAMOU® ENDOMAT®
	263310 20-1	ENDOMAT® según HAMOU®

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Thank you for your expression of confidence in the KARL STORZ brand name. Like all of our other products, this product is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality item of equipment from KARL STORZ.

This instruction manual is intended to serve as an aid in the proper setup, installation, and operation of the HAMOU® ENDOMAT®. All essential details of the equipment and all actions required on your part are clearly presented and explained. We thus ask that you read this manual carefully before proceeding to work with the equipment. Insert this manual in its protective wallet and keep it available for ready reference in a convenient and conspicuous location near the equipment.

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Controls, displays, connectors, and their uses

- ① Main power switch
- ② Mode of operation indicator lamps
- ③ Flow-rate indicator
- ④ Flow-rate setpoint indicator
- ⑤ Irrigation pressure level indicator
- ⑥ Irrigation pressure limit indicator
- ⑦ Suction pressure limit indicator
- ⑧ Suction pressure setpoint indicator
- ⑨ Pushbutton, releases pump rollers
- ⑩ Label strip (replaceable)
- ⑪ Pressure transducer
- ⑫ Pilot lamp
- ⑬ Irrigation START/STOP button
- ⑭ Suction line connector
- ⑮ Suction START/STOP button
- ⑯ ± buttons for adjusting suction pressure setpoint
- ⑰ ± buttons for adjusting irrigation pressure setpoint
- ⑱ ± buttons for adjusting flow-rate setpoint
- ⑲ Footswitch connector
- ⑳ **SCB**[®] * connectors
- ㉑ Potential equalization connector
- ㉒ Power fuses
- ㉓ Power supply cord receptacle
- ㉔ Fluid collection bottle
- ㉕ Fluid collection bottle vacuum line fitting (with bacterial filter)
- ㉖ Fluid collection bottle instrument line fitting

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

* The KARL STORZ-SCB[®] interface (KARL STORZ Communication Bus), based on the CAN field bus, permits remote control of equipment functions, as well as remote display of equipment parameters.

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
Note: The reusable HYS 263300 40 and LAP 263300 41 tubing sets are available as an optional extra (see 'Recommended accessories').

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Device description

The HAMOU® ENDOMAT® is a combination irrigation/suction pump for use in gynecological and laparoscopic surgery. Use of the most advanced technologies allows its application in both hysteroscopic and laparoscopic surgery. Adaptation to the type of surgery intended via provision of the relevant optimal operating parameters occurs automatically when the correct type of tubing set for the hysteroscopic or laparoscopic application is inserted. Inserting pressure lines has been greatly simplified through use of a redesigned reel pump. The unit electronically monitors the correct positioning of pressure lines.

Operator adjustments and unit operation are made using precision pushbutton switches. Displays for setpoints and current actual performance data, which are positioned next to each other, allow optimal monitoring of unit operational status. Electronic safety circuitry interrupts irrigation/suction operation in the event that unit performance persistently deviates from the setpoint, and issues an audible alarm to warn the operator. An electronic self-check routine that checks all unit components each time the unit is switched on and notifies the operator of any failures detected further enhances the HAMOU® ENDOMAT®'s high standards of safe operation.

Property rights

This product is protected in the USA by (at least one of the following) US Patent No(s) 5,788,688; 6,397,286; 6,484,221; 6,824,539.

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Warnings and cautions

Please read this manual and follow its instructions carefully. The words **Warning**, **Caution**, and **Note** convey special meanings. Wherever they are used in this manual, they should be carefully reviewed to ensure the safe and effective operation of this product. To make the words, **Warning** and **Caution** stand out more clearly, they are accompanied by a pictogram.

Warning: A Warning indicates that the personal safety of the patient or physician may be involved. Disregarding a Warning could result in injury to the patient or physician.

Caution: A Caution indicates that particular service procedures or precautions must be followed to avoid possible damage to the product.

Note: A Note indicates special information about operating the product, or clarifies important information.

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Warning: Read this instruction manual thoroughly and be familiar with its contents prior to using this equipment. Before using the unit, read the following safety instructions carefully to avoid putting your patients, personnel or yourself at risk.

Note: Any damage to the unit resulting from incorrect operation is not covered by the manufacturer's warranty.

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Warning: The electrical installations in the operating room in which the unit is connected and operated must comply with the applicable IEC standards.

Warning: Keep out of reach of patients.

Warning: The instructions and interface specifications for medical devices and/or system components used in combination must be observed precisely.

Warning: Combinations of medical devices are only assured to be safe if

- they are identified as such in the respective instruction manuals or
- the intended purpose and interface specifications of the devices used in combination permit this (cp. IEC 60601-1-1).

Warning: Grounding reliability can only be achieved when the equipment is connected to "Hospital Only" or "Hospital Grade" outlet (i.e., approved for use in an operating room environment). Routinely inspect electrical plug and cord. Do not use if inspection reveals damage.

Warning: Test this equipment prior to each surgical procedure to ensure that it functions correctly.

Warning: DO NOT use the HAMOU® ENDOMAT® in the presence of flammable anesthetics.

Warning: Always unplug the unit before all maintenance work (e. g. cleaning).

Warning: To avoid the risk of electrical shock, do not open the unit. Refer servicing to qualified personnel.

Warning: Pressure-controlled pump motor. When the operating switch is on, never reach into the pump rollers even if the pump is stationary. The pump may start up at any time when the operating switch is on.

Warning: The tubing set must be vented before use on a patient.

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Warning: It is important to carefully monitor the fluid level in the irrigation bag to ensure that the irrigant does not completely run out and result in air being pumped into the patient. Another bag of irrigant should be started prior to emptying the bag being used.

Warning: The HAMOU® ENDOMAT® does not measure intracorporeal pressure. The pressure setting is the pressure at which the irrigant is being pumped into the system. Irrigation line pressure and flow can be read on displays ⑤ and ③, respectively. An intermittent "beep" tone will be heard if the pressure limit setting is briefly exceeded. For additional safety there is an overpressure valve integrated into the pressure measuring chamber to prevent maximum delivery pressure from being exceeded.

Warning: Use of electromechanical medical devices can cause electromagnetic interference in other devices, particularly cardiac pacemakers. Precautions must be taken to ensure that the patient's well being is maintained in the event of such interference.

Warning: DO NOT reuse the disposable tubing set. DISCARD after use.

Warning: Reusable tubing must be cleaned, disinfected and sterilized before use and subsequent reuse using validated preparation procedures.

Warning: For all cleaning and disinfection procedures on contaminated instruments, the guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed.

Warning: The tubing set must be vented before every examination to remove any air bubbles.

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Caution: Only operate the unit at the line voltage stated on the manufacturer's identification plate.

Caution: When replacing fuses, use only fuses of the correct rating.

Caution: The HAMOU® ENDOMAT® may only be operated with tubing sets and accessories which have been designated as suitable for the unit by KARL STORZ.

Caution: Avoid allowing fluids to enter the unit. Do not store liquids on or above the unit.

Caution: Always position the HAMOU® ENDOMAT® at or above the position of the patient to assure proper function of the unit.

Caution: If a power failure occurs after the tubing set has been installed on the unit, the self-test procedure can be by-passed by depressing the footswitch when turning on the power switch.

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Note: This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE). At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point. Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

Note: Make certain to observe national regulations/laws regarding disposal (e.g. used irrigation fluid, tubing sets, irrigation bags, etc.).

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Before using the unit on the patient it is imperative that you be acquainted with how the unit operates and is controlled.

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Intended use

The HAMOU® ENDOMAT® is a combined suction/irrigation pump for use in hysteroscopic and laparoscopic interventions. Application is performed using special suction/irrigation instruments.

Use of the unit in fields other than those indicated above is not allowed for safety reasons.

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The HAMOU® ENDOMAT® may only be operated with tubing sets and accessories which have been designated as suitable for the unit by KARL STORZ.

No responsibility for the safe operation of the instrument can be accepted if any tube systems other than those specified are used. The following serious defects may arise in particular:

- Inadequate grip on the tube connectors means a loss of sterility and may lead to the unit being damaged by ingress of liquid.
- Premature breakage of the pump tube segment and consequent loss of sterility and damage to the unit caused by ingress of liquid.
- Display of an incorrect actual irrigation rate due to deviations in tube diameter.

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Unauthorized conversions or modifications to the unit are not allowed for safety reasons.

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
User qualification

The HAMOU® ENDOMAT® may only be used by physicians and medical assistants who have a corresponding specialized qualification and who have been instructed in use of the unit.

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Safety precautions at the site of installation

The unit may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.

It is not intended for use in hazardous zones. This means, for example, that when using easily combustible and explosive inhalation anesthetics or mixtures thereof, the unit must not be operated inside the hazard zone shown in the diagram. This also applies to highly combustible and potentially explosive chemicals, e.g. skin disinfectants and fast-acting surface disinfectants.

The unit is equipped with a connector for attaching a ground line. It should be connected up in accordance with national regulations.

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Safety precautions when operating the unit

It is the user's responsibility to make sure the unit is safe and operates properly before using it.

During treatment with the ENDOMAT® the patient must be treated and kept under observation with the usual medical care. This includes keeping a check on the progress of treatment, as well as monitoring the vital levels and the anesthetic.

Any treatment may only be performed if there is visual observation of the action of the unit.

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Safety features

The HAMOU® ENDOMAT® is equipped with the following safety features:

Self-Check

A self-check procedure is performed each time the unit is switched on. The components checked are as follows:

- internal self-check of the electronic circuitry
- tube detection test
- pressure measuring chamber detection test
- pressure sensor test

If the results of any of these tests are negative, it will not be possible to start the unit. In this case the unit will issue an intermittent warning sound for 10 seconds. The error display will continue to flash until the error is remedied.

Note: In order to insert the tube set before putting the unit into operation, the self-check can be bypassed by pressing the footswitch during power up.

Make sure the unit performs a self-check procedure at least once every day of operation.

Display element test

On switching on the instrument all the display elements will light up. This is to enable the user to check the display for any faulty segments.

Monitoring the irrigation pressure

Irrigation will be stopped as soon as the actual pressure value exceeds the specified setpoint. If the deviation lasts for more than several seconds, there will be an audible warning signal.

As an additional safety measure the pressure measuring chamber was fitted with an integrated overpressure valve which prevents the maximum irrigation pressure from being exceeded.

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Unpacking the equipment

Carefully unpack the HAMOU® ENDOMAT® and remove the unit and its accessories from their packaging. Check for missing items and evidence of shipping damage. File any complaints with the manufacturer or supplier immediately.

If possible, retain the original packing materials for later use; these can come in handy if the unit has to be transported.

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Installation

Note: The HAMOU® ENDOMAT® including the accessories connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations.


The unit is not intended for use in hazardous zones. Do not operate the unit within the hazard zone shown in the diagram while explosive anesthetic gases are in use.

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Fitting a label strip to the front panel

Undo the left-hand front side cover with a screwdriver.

Note: The appropriate controls/indicators on the unit are designated on the label strip.

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Remove any unmarked label strip. Select a label strip ⑩ (currently 6 languages available) and slide it between the plastic and the metal front. Fold the strip along the broken line to an angle of 90° and approx. 2 mm further on (toward the end of the strip) fold it back to an angle of 90° so that the side cover can be refitted without any obstructions. Replace the side cover.

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The unit is equipped with a connector ② for attaching a ground line. The unit's ground line should be installed by a qualified electrician.

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Installation and operating instructions

Connect power cord. Insert power cord into power cord receptacle 23 as far as it will go.

Warning: Operate the unit only with the system voltage stated on the rating plate.

Warning: Only insert the power plug into and remove it from electrical outlets located outside areas subject to explosion hazards.

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KARL STORZ-SCB®

Note: To avoid mistakenly pulling out the SCB® connecting cable, the SCB® connector possesses a protection device.

Pull back the protection device of the SCB® connector and insert the connector into one of the SCB® sockets 20. Connect the other end of the cable to the KARL STORZ-SCB® (KARL STORZ Communication Bus) control unit or other SCB® units (see KARL STORZ-SCB® System Instruction Manual).

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If necessary, connect footswitch (see Accessories) to the socket 19.

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
Mode of operation selection

The HAMOU® ENDOMAT® unit is suitable for both hysteroscopic and laparoscopic applications. The appropriate mode of operation will be automatically selected when the correct type of tubing for the application is inserted. The unit determines the intended application from the differing 'encodings' of tubing sets and then automatically selects the appropriate mode of operation.

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Operating the unit

Switch on the unit at its main power switch ①. The red pilot lamp ⑫ will illuminate. The unit will then conduct a self-check procedure, checking all system components for proper operation.

The indicators on the setpoint bargraphs ④, ⑥ and ⑧ for flow, irrigation pressure and suction pressure are flashing.

Note: Once switched on, the unit awaits confirmation or a change of the last set values used (see note, p. 17). This is signalled by flashing of the corresponding LEDs in the bargraph displays. The LED in the irrigation START/STOP button ⑬ also flashes. As long as the displays are flashing, the START function is blocked.

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Set value confirmation

To confirm all set values, press one of any ± buttons or the START/STOP button ⑬.

The LED in the START/STOP button ⑬ goes out and the LEDs in the bargraph displays are continuously illuminated.

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Power-up with irrigation START/STOP button pressed down

At device power-up with the footswitch pressed down or with the irrigation START/STOP button ⑬ pressed down, tests for the following routines are not performed: pressure measuring chamber detection test, irrigation sensor test and irrigation sensor calibration.

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Inserting the tubes

Note: The tubing set must not be inserted until the unit has been switched on with the power switch ①. If the tubing set is inserted before switching on, the self-check procedure blocks the equipment. If the footswitch or the irrigation START/STOP button ⑬ is pressed simultaneously while switching on the power switch, the self-check procedure may be avoided. Make sure the unit performs a self-check procedure at least once every day of operation.

When inserting the tubing set, push button ⑨ to release the pump rollers. The pump rollers spring back.

Caution: Use only tubing sets designed for the HAMOU® ENDOMAT®.

Note: The pump tube connectors have different widths so that the pump tube can be inserted in one direction only.

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Disposable tubing set

Place a length of tubing on the tube holder. Position the pressure measuring chamber on the pressure transducer ⑪ and turn clockwise until the pressure measuring chamber locks into place (bayonet-lock).

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Install the reusable tubing set for LAP 263300 41 or HYS 263300 40, available as an optional extra, in the same way as you would the disposable tubing set.

The instruction manual delivered with the tubing sets must be observed.

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Connect the correct end of the length of tubing to the irrigation-liquid bag (penetration fitting), and to the 'LUER-lock' irrigation-line fitting on the instrument to be employed.

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Once the tubing set and pressure measuring chamber have been correctly installed, the pilot lamp ⑫ will switch from red to green.

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The mode of operation selected will be indicated by lamp ⑫.

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Caution: If the pilot lamp ⑫ remains red, either there is a failure somewhere in the unit, or the tubing set has been incorrectly inserted. If the pilot lamp ⑫ fails to switch from red to green after the tubing set has been removed and reinserted, then the unit will have to be checked by an authorized technician before it can be used.

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Removing the tubing set

Press pushbutton ⑨ to release pump rollers.
Remove the tubing from the tube holder.
Loosen and detach the pressure measuring chamber from the pressure transducer by turning it counterclockwise.

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Note: For further information on the maintenance/preparation of pressure chambers/tubing sets, please consult the instruction manual delivered with the reusable tubing set.

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**Connecting the tubing set
for suction (optional)**

Screw the bacterial filter into the vacuum line fitting 25 on the suction bottle and attach the suction line to the filter housing.

Caution: To avoid contaminating the unit a bacterial filter must be fitted between the fluid collection bottle and the unit.

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Connect the other end of the suction line to the fitting 14.

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Use a length of tubing to interconnect the suction line fitting of the fluid collection bottle 26 and the suction line fitting on the instrument to be employed.

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
Note: Pressure, flow rate, and suction pressure settings are retained in the unit's memory while the unit is switched off, and will be redisplayed flashing when the unit is switched on. To confirm the set values, **one-time activation** of any ± buttons or the START/STOP button ⑬ is enough, or reset them as described following.

Then the LED in the START/STOP button ⑬ goes out and the LEDs in the bargraph displays are continuously illuminated.

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
Presetting irrigation parameters

Enter the flow rate setpoint using the associated ± buttons ⑱.
The value set will be indicated on the associated display ④.

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
Enter the limiting irrigation pressure using the associated ± buttons ⑰.
The value set will be indicated on the associated bargraph display ⑥.

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Presetting suction pressure

Enter the suction pressure setpoint using the associated ± buttons ⑲. The value set will be indicated on the associated bargraph display ⑧.

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Test for proper operation and air-relief of the tubing system

 **Warning:** Test this equipment prior to each surgical procedure to ensure that it functions correctly.

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Note: Have a collecting container ready to catch irrigation liquid as it runs out.


Press the irrigation START/STOP button ⑬ or the footswitch until all the air has been let out of the tubing system (i.e. no bubbles) (see Warning, page 5).

The flow rate indicated on bargraph display ③ must agree with the flow-rate setpoint indicated on display ④.

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Gradually close the metering valve on the instrument. As soon as the pressure indicated on bargraph display ⑤ rises above the limiting value indicated on display ⑥, the reel pump should be halted. If the actual value exceeds the set value for a number of seconds there is an audible warning signal.

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
Pinch off the suction line with your fingers and check whether suction line underpressure is formed. The suction line pressure indicated on display ⑦ must agree with the suction line pressure setpoint indicated on display ⑧.

 **Warning:** If there are any substantial deviations from these results ($> \pm 10\%$), the unit must be inspected by an authorized technician before further use.

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Initiating irrigation

Irrigation can now be initiated by actuating the irrigation START/STOP button ⑬, or by depressing the footswitch. Irrigation line pressure and flow rate can be read on displays ⑤ and ③, respectively.

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If the set pressure should be briefly exceeded the pump will be halted and an audible alarm will be heard after a few seconds. For passive safety there is an overpressure valve integrated into the pressure measuring chamber to prevent maximum irrigation pressure from being exceeded.

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Initiating suction

Suction can be initiated by actuating the suction START/STOP button ⑮. Suction pressure can be read on display ⑦.

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
Fuse replacement

Switch off the unit and remove the power plug from the electrical outlet.

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
Remove the power fuse holder ② with a screwdriver or other tool.

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! Caution: Only use fuses of the correct rating.


Insert new fuses of the appropriate rating.

	100...240 V~
Power fuse	2 x T2AL250V

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Refit the power fuse holder ②.
Reconnect the power cord.
Test the unit for proper operation.

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Cleaning, disinfection and sterilization

Note: The instructions for “Cleaning, Sterilization and Care of KARL STORZ Instruments” must be followed. There the procedures for Cleaning, Disinfection and Sterilization are explained in detail.

Warning: Always pull out power plug before cleaning!

Caution: Avoid allowing fluids to enter the unit.

The exterior surfaces of the unit must be cleaned with a disposable cloth moistened with disinfectant (preferably free of alcohol).

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Cleaning, disinfection and sterilization of reusable tubing

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Warning: When cleaning and disinfecting contaminated instruments (tubing sets), the local guidelines of the Employer’s Liability Insurance Association and equivalent organizations must be observed.

Warning: Prior to the sterilization of reusable tubing, particularly in cases of manual preparation, it must be ensured that all chemical residue has been rinsed off or removed (see also preparation steps).

Warning: Disposable tubing sets, locks and hydrophobic bacteria filters must not be reused.

Note: For the cleaning/disinfectant solution, use chemicals that have been specially approved by KARL STORZ for the preparation of medical devices (see Appendix, pages 37, 38).

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Preliminary cleaning of reusable tubing

Place the reusable tubing in a container filled with a cleaning solution (as per manufacturer's instructions) immediately after use to prevent contaminants from drying on the surface of the tubing set.

Caution: When preparing and using the solutions, follow the manufacturer's instructions, paying close attention to proper concentration and exposure times. Prolonged immersion can lead to material changes.

Water quality requirements

Tap water can be used to prepare the cleaning and disinfectant solutions. When rinsing for the last time, for example to remove any chemical residue, microbiologically pure/sterile water must be used.

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Cleaning/Preparation

Warning: The local guidelines of the Employer's Liability Insurance Association and equivalent organizations must be observed when performing all cleaning and disinfection procedures on contaminated instruments.

Manual cleaning of the suction bottle

The suction bottle must be cleaned thoroughly, as must the overflow protection unit. Clean the channels with a brush, taking care to avoid spraying. Then rinse thoroughly, if necessary repeat this procedure several times. Ensure that the floats in the overflow protection unit can move freely and do not become stuck due to residue.

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Manual preparation of reusable tubing

1. Disconnect the reusable tubing from the connectors as necessary and then immerse it in a cleaning solution (made up according to the instructions provided by the manufacturer).
2. Using a syringe, rinse out the tubing underneath the (water) surface. Repeat the procedure several times, if necessary. Remove any contaminants using brushes, sponges, soft cloths or a cotton applicator.
Cleaning accessories are available from KARL STORZ (see catalog).
3. After cleaning, rinse the tubing and place it in a suitable container with disinfectant solution, making sure that all air bubbles have been expelled. Ensure that the tubing is completely covered in the solution for the duration of exposure.
4. Cleaning brushes must be cleaned and disinfected immediately.
5. Remove the tubing after the exposure time and rinse with microbiologically pure/sterile water. Ensure that all chemical residue in and on the tubing is removed.
6. Dry the tubing with sterile compressed air.
7. After preparation, inspect the tubes for cleanliness and damage. If soiling is still visible, the preparation procedure must be repeated.

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Machine preparation of reusable tubing

To prepare using a machine, connect the reusable tubing to the retaining mechanisms in such a way that the instrument channels can be completely rinsed through. Visible, heavy soiling must be removed in preliminary cleaning prior to machine preparation.

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! Warning: Check the entire tubing set for leaks after preparation. As a result, it is possible to detect a leak early, through which liquid could escape and enter the unit.

To check for leaks use leakage tester 13242 XL for example. For this purpose connect one end to the leakage tester and block off the other end, with your finger for example. If the tubing set develops a leak, a pressure drop will be observed on the pressure gauge of the leakage tester.

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Sterilization

The suction fluid collection bottle and lid, including the overflow protection unit, can be either steam sterilized or gas sterilized.

We recommend steam sterilizing the reusable tubing before use at 134°C (+3°C) using the fractionated prevacuum procedure.

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! Warning: Successful sterilization is only possible on clean, disinfected surfaces.

Warning: The recommended sterilization parameters are only valid with sterilization equipment that is properly maintained and validated.

Warning: Any deviations from the recommended parameters for sterilization should be validated by the user.

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Steam Sterilization

KARL STORZ has validated the following sterilization method:

Fractionated prevacuum (steam sterilization)

Warning: Tubing sets may only be sterilized using the fractionated prevacuum procedure.

Trays should be positioned in the sterilizer so that there is adequate circulation and penetration of steam, air removal and condensate drainage. Load the sterilizer according to its instruction manual.

At the completion of the steam sterilization cycle, all instruments should remain untouched until adequately cooled.

Warning: Burn hazard! Parts are hot after steam sterilization. Allow them to cool down.

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Fractionated prevacuum procedure

The fractionated prevacuum procedure consists of four phases.

- **Conditioning phase**
The conditioning phase removes air from the chamber by pulling a vacuum up to four times. Steam is then blown in.
- **Sterilization phase**
Sterilization takes place at a temperature of 134°C (+3°C) and an pressure (absolute) of 300 kpa (3 bar) over a period of 5 minutes (minimum time).
- **Evacuation phase**
The steam is released.
- **Drying phase**
Drying takes place under renewed application of a vacuum over a period of approx. 5–20 minutes.

Validated parameters for the fractionated prevacuum method:

Temperature: 134°C (+3°C)

Pressure (p_{absolute}): 300 kpa (3 bar)

Exposure time: 5 minutes for all instruments.

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
Maintenance

Preventive maintenance is not essential. Regular maintenance can, however, contribute to identifying potential problems before they become serious, thus enhancing the instrument's reliability and extending its useful operating life. Maintenance services can be obtained from your local representative or from the manufacturer. Regardless of the accident prevention regulations or testing intervals for medical instruments prescribed in different countries, we recommend a functional or safety test of the unit at least once a year. For detailed instructions please refer to the version of the Service Manual currently valid, catalog number SV 3373 (English version).

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Servicing and repair

Defective items of equipment must be serviced and repaired exclusively by persons authorized by us; all repair work must employ original KARL STORZ parts only.

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Disposal

This unit has been marked in accordance with the European Directive 2002/96/EC on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the unit as electronic scrap.

Please ask either KARL STORZ GmbH & Co. KG, a KARL STORZ subsidiary or your specialist dealer for information on your local collection point.

Within the scope of application of this Directive, KARL STORZ GmbH & Co. KG is responsible for the proper disposal of this unit.

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Repair program

For fiberscopes and equipment, individual repair is necessary. Usually to bridge the repair period, you will receive a unit on loan which you then return to KARL STORZ as soon as you receive the repaired unit.

In Germany you can refer repairs direct to

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 TUTTLINGEN

In other countries please contact your local KARL STORZ branch or authorized dealer.

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Important information

In order to prevent infection, it is strictly forbidden to ship contaminated medical devices. All medical devices must be decontaminated on site, to avoid contact and aerogenous infections (among staff).

We reserve the right to return contaminated instruments/equipment to the sender.

Repairs, changes or expansions which are not carried out by KARL STORZ or by experts authorized by KARL STORZ will invalidate all guarantee rights. KARL STORZ gives no guarantee on the correct functioning of equipment or instruments which have been repaired by unauthorized third parties.

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Limitation of liability

KARL STORZ shall be liable for failure or deterioration in the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the unit is connected and operated meet the applicable laws and standards; and
- the unit has been used in accordance with its operating instructions at all times.

Manufacturer's warranty

For a period of two years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship. In so doing we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Please fill out the attached warranty card on the last page and return it as soon as possible to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 TUTTLINGEN/GERMANY

Note: The unit must always be sent to the following address (also during the warranty period; where applicable, with warranty card):

KARL STORZ GmbH & Co. KG
Abt. Reparaturservice
Dr. Karl-Storz-Str. 34
78532 TUTTLINGEN/GERMANY

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such opening, repair, or modification performed during the warranty period shall void all warranty.

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Troubleshooting

Warning: Always unplug the unit before all maintenance work.

Symptom:

- Complete failure of the unit.

Possible causes:

- No power from the power line.
- Defective power fuse.
- Power cord is not properly connected to power cord connector.

Remedy:

- Check that there is electricity to the wall outlet.
- Change fuses as described in the instruction manual. Use only fuses of the correct rating.
- Push power cord firmly into power cord connector.

Symptom:

- Poor suction.

Possible causes:

- Leakage occurring in suction lines.

Remedy:

- Check all suction lines, and replace any exhibiting evidence of leakage.
- Check the seating of the sealing cap.

Symptom:

- No suction.

Possible causes:

- Floats of overflow protection unit blocking suction inlet.

Remedy:

- Check the fluid level in the vessel, and empty, if required.
- Clean the floats and check them for freedom of movement.

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Error codes

If an error occurs when the unit is powered up, the unit will issue an intermittent warning sound for 10 seconds. The error display will continue to flash until the error is remedied. The number of flashing pulses indicates the error code.

Only one error code will be displayed: the one which was detected first. In addition, the error code is shown on display ③ 'Irrigation flow act. value'. The error code corresponds to the number of LEDs lit up. This code is for the user's information or provides Service with preliminary information.

Error	Description of the error codes
1	Pressure measuring chamber attached while switching on
2	Checksum error of the EPROM (checksum ≠ 0)
3	Error during RAM test
4	Error during initial flow measurement, flow value > 0 (≠ 0)
5	Offset error of pressure sensor, p measured < -20 mmHg
6	Offset error of pressure sensor, p measured > 40 mmHg
7	Motor control module too hot > 135°C
8	Excessive motor current, irrigation pump > 200 mA ± 50%
9	Excessive motor current, suction pump > 200 mA ± 50%
10	Offset error, suction pressure too low > 40 mbar
11	Offset error, suction pressure too high < -40 mbar
12	Tubing detection operated
13	Motor control too hot (during operation)
14	Excessive irrigation-pump-motor current (during operation)
15	Excessive suction-pump-motor current (in operation)
16	Maximum number of SIO errors reached

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Standard compliance (for 263310 20-1)

According to IEC 60601-1, UL 2601.1, CAN/CSA C22.2 No. 601.1-M90:

- Type of protection against electric shocks: Protection Class I
- Degree of protection against electric shocks: Applied part of type BF

According to IEC 60601-1-2:

Please read the Electromagnetic Compatibility Information in the appendix (p. 39-51).

Directive compliance (for 263310 20-1)

According to Medical Device Directive (MDD):

This medical product belongs to Class II b

This medical product bears the CE mark in accordance with MDD 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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Technical documentation

On request, the manufacturer will provide those circuit diagrams, itemized parts listings, descriptions, sets of adjustment instructions, and other items of available documentation to suitably qualified user personnel duly authorized by the manufacturer for their use in repairing those components of the unit that have been designated by their respective manufacturers as repairable.

Supply of such technical documentation relating to the unit shall not be construed as constituting manufacturer's authorization of user's personnel, regardless of their level of technical training, to open or repair the unit.

Explicitly exempted herefrom are those maintenance and repair operations described in this manual.

We reserve the right to make engineering modifications in the interest of promoting technological progress and generating performance improvements without obligation on our part to submit prior notice thereof.

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
Extract from the list of approved chemicals for the preparation of telescopes and endoscopic instrument sets
Manual preparation

The following instruments and telescopes produced by KARL STORZ must not be fully immersed in a liquid: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, prismatic light deflectors. For motors and their handles the appropriate manual provides detailed information.

The following instruments produced by KARL STORZ must not be placed in an ultrasound cleaning bath: rigid telescopes, flexible telescopes, other instruments with optical components made of glass.

When preparing and using the solutions, follow the solution manufacturer's instructions, paying close attention to concentration and exposure times. Regarding the microbiological effectiveness please ask the manufacturer. The following cleaning and disinfectant solutions have been approved regarding material compatibility:

Ecolab GmbH & Co. OHG	Sekusept Pulver Classic	1, 2, 3
	Sekusept Pulver Classic + Aktivator	1, 2, 3
Estee Pharma GmbH	Ultrasmit AF	1, 2, 3
	Ultrasseptin Aktiv	1, 2, 3
	Ultrasseptin Classic	1, 2, 3
Fresenius AG	Afid	1, 2
	Afid plus (neu)	1, 2
	Sporcid FF	1, 2
Indeba Industria E Comercio LTDA	Letahdeido	2
Jose Collado S.A.	Darodor 4000 Liquido	1
	Darodor 9000	2
	Darodor Sinaldehyd 2000	2
Lonza AG	ID 50	2
	Aldasan 2000	1, 2
Lysoform Dr. Rosemann GmbH	Almyrol	1, 2
	Desoform	1, 2
	Lysoformin 3000	1, 2
Medichem International	MedDis	2
	MediZyme	1
	Edisonite Super	3
Merz Consumer Care GmbH	Mucadont-IS	2, 3
	Mucadont-Zymaktiv	1, 3
	Mucocit-T	1, 2, 3
Orochemie	A 10	2
	ID 210	2
Promagent AB	Wavacide	2
	Gigasept	2, 3
	Gigasept FF	2, 3
Schülke & Mayr GmbH	Gigasept Instru AF	1, 2, 3, A
	Lysefol V	2
Schuhmacher, Dr. GmbH	Descoton forte	2
	Perfektan Endo	2
Steris	EnzyCare 2	1, 3
	neodisher LM 2	1
	neodisher mediclean	1
Weigert, Dr. GmbH & Co.	neodisher mediclean forte	1, 3
	neodisher medizym	1, 3
Whiteley Industries PTY.Ltd	Aidal Plus	2

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*** Caution:** The use of CAUSTIC SODA can cause changes to the surface of aluminium parts (including coated parts), plastics and soldered connections, and may reduce the service life of the instrument set. For **flexible endoscopes**, only immerse the examination sheath, **not** the housing/handle part.

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Typo

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- A = Not suitable for flexible endoscopes (according to manufacturers' instructions)

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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Machine preparation

The following instruments and telescopes produced by KARL STORZ must not be prepared in a cleaning and disinfecting machine: HAMOU® telescopes with ocular focusing, magnifiers with focusing, IMPERATOR drill handles, standard, tungsten carbide and diamond-tipped drills, prismatic light deflectors, and magnetic instrument rack. For motors and their handles the appropriate manual provides detailed information.

Selection of a method for cleaning and disinfection must be discussed with the machine manufacturer and the manufacturer of the chemical reagents. Only special methods must be used which have been verified for this purpose. Regarding the microbiological effectiveness please ask the manufacturer. Thermal disinfection is preferable. The following agents for machine cleaning and disinfection have been approved regarding material compatibility:

Schuhmacher, Dr. GmbH	Thermoton Endo	2
Steris	Instru-Klenz	1, 3
	neodisher FA	1
Weigert, Dr. GmbH & Co.	neodisher Oxivario (only in combination with neodisher FA)	1
	neodisher FA forte	1, 3
	neodisher mediclean	1
	neodisher mediclean forte	1, 3
	neodisher mediklar	5
	neodisher medizym	1
	neodisher Septo DA	2
	neodisher Septo DN	2

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Type

- 1 = Cleaner
- 2 = Disinfectant
- 3 = Suitable for ultrasound baths
- 4 = Neutralizer
- 5 = Rinse aid

Note: You will also find the current list of approved products on our Web site (www.karlstorz.com).

Note: When agents other than those named are used, KARL STORZ assumes no liability for any damage which may occur.

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WARNING: Medical electrical equipment needs special precautions regarding Electromagnetic Compatibility (EMC). Observe the EMC instructions in this appendix during installation and commissioning.

The HAMOU® ENDOMAT® model 263310 20-1 corresponds to EN/IEC 60601-1-2 :2001 [CISPR 11 Class B] and therefore meets the EMC requirements of the Medical Device Directive (MDD) 93/42/EEC.

These limits are designed to provide reasonable protection against the typical electromagnetic interference to be expected in a medical environment. The HAMOU® ENDOMAT® model 263310 20-1 is a Group 1 unit (as per CISPR 11). Group 1 contains all the 'equipment and systems which generate or use RF energy only for their internal functioning'.

NOTE: The tables and guidelines that are included in this appendix provide information to the customer or user that is essential in determining the suitability of the equipment or system for the electromagnetic environment of use, and in managing the electromagnetic environment of use to permit the equipment or system to perform its intended use without disturbing other equipment and systems or non-medical electrical equipment. If this equipment does cause harmful interference with other devices, the user is encouraged to try to correct the interference by one or more of the following measures:

- reorient or relocate the receiving device
- increase the separation between the equipment
- connect the equipment into an outlet on a circuit different from that to which the other device(s) is connected.

If you have any further questions, please contact your local representative or our service department.

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WARNING: The HAMOU® ENDOMAT® model 263310 20-1 should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

WARNING: The use of portable and mobile RF equipment may have an impact on this and other pieces of medical equipment.

WARNING: The use of accessories and cables other than those specified in the KARL STORZ instruction manual may result in increased emissions or decreased immunity of the HAMOU® ENDOMAT® model 263310 20-1. The accessories and cables listed below have been shown to comply with the requirements of EN/IEC 60601-1-2. When using accessories and cables other than those specified here, it is the responsibility of the user to ensure that they comply with EN/IEC 60601-1-2.

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Table 200				
Accessories and cables which have been shown to comply with EN/IEC 60601-1-2:				
Type	Shielded	Length [m]	Ferrite	Used for
200103 30	Yes	3	No	Footswitch model 200103 30
200900 30	Yes	4	No	SCBcom model 200900 30
PA	No	> 3	No	Potential equalization
Power cord	No	3	No	Power supply

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Table 201		
Guidance and manufacturer's declaration – electromagnetic emissions		
HAMOU® ENDOMAT® model 263310 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the HAMOU® ENDOMAT® model 263310 20-1 should ensure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	HAMOU® ENDOMAT® model 263310 20-1 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The HAMOU® ENDOMAT® model 263310 20-1 is suitable for use in all establishments (e.g. hospitals and doctors' offices) including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

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Table 202			
Guidance and manufacturer's declaration – electromagnetic immunity			
HAMOU® ENDOMAT® model 263310 20-1 is intended for use in the electromagnetic environment specified below. The customer or user of the HAMOU® ENDOMAT® model 263310 20-1 should ensure that it is used in such an environment.			
Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	Complies ± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	Complies ± 2 kV for power supply lines ± 1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	Complies ± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T * (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 sec	Complies <5% U_T * (>95% dip in U_T) for 0.5 cycle Complies 40% U_T (60% dip in U_T) for 5 cycles Complies 70% U_T (30% dip in U_T) for 25 cycles Complies <5% U_T (>95% dip in U_T) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the equipment or system requires continued operation during power mains interruptions, it is recommended that the equipment or system be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	Complies 3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
* Note: U_T is the a.c. mains voltage prior to application of the test level.			

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Table 204

**Guidance and manufacturer's declaration – electromagnetic immunity –
for medical electrical equipment and systems that are not life-supporting**

HAMOU® ENDOMAT® model 263310 20-1 is intended for use in the electromagnetic environment specified below.
The customer or user of the HAMOU® ENDOMAT® model 263310 20-1 should ensure that it is used in such an environment.

Immunity test	EN/IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
<p>Conducted RF IEC 61000-4-6</p> <p>Radiated RF IEC 61000-4-3</p>	<p>$3 V_{rms}$ 150 kHz to < 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>$3 V_{rms}$</p> <p>$3 V/m$</p>	<p>Portable and mobile RF communications equipment should be used no closer to any part of the HAMOU® ENDOMAT® model 263310 20-1, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = [3.5/3]\sqrt{P} \quad 150 \text{ kHz to } < 80 \text{ MHz}$ $d = [3.5/3]\sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = [7/3]\sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts [W] according to the transmitter manufacturer and d is the recommended separation in meters [m].</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a, should be less than the compliance level in each frequency range^b.</p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection by structures, objects and people.

a Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than **3 V/m**.

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Table 206

**Recommended separation distances between portable and mobile RF communications equipment and the
HAMOU® ENDOMAT® model 263310 20-1**

The HAMOU® ENDOMAT® model 263310 20-1 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] according to frequency of transmitter		
	150 kHz to < 80 MHz	80 MHz to < 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance **d** in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where **P** is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

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WARRANTY

For a period of **two** years after delivery to the end-user, we shall replace free of charge equipment that can be proven to have defects in material or workmanship.

In doing so we cannot accept to bear the cost of transportation or the risk of shipment. The warranty referred to in our Standard Conditions of Business shall apply.

Complete warranty card upon purchase/delivery and return immediately to:

KARL STORZ GmbH & Co. KG
Postfach 230
78503 Tuttlingen/Germany

To be completed by the supplier/importer:
Company stamp/signature:

To be filled out by instrument owner:
Return address/company stamp:
Field of application:
Type of instrument: Serial no.:
Purchase date:
Signature/Date:

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KARL STORZ Endoscopy China Ltd.
Hong Kong Representative Office
Unit 1601, Chinachem Exchange Square
1 Hoi Wan Street, Quarry Bay, Hong Kong,
People's Republic of China
Phone: +8 52 28 65 2411, Fax: +8 52 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Beijing Representative Office
Room 610, China Life Tower
No. 6, Chaowai Street
Beijing, 100020,
People's Republic of China
Phone: +86 10 8525 3725,
Fax: +86 10 8525 3728
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Unit 3901-3904, Tower 1 Grand Gateway
No. 1 Hong Qiao Road
Shanghai, 200030,
People's Republic of China
Phone: +86 21 6113-1188,
Fax: +86 21 6113-1199
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Chengdu Representative Office
F-5, 24/F., Chuanxing Mansion,
No. 18 Renming Road South
Chengdu, Sichuan, 610016,
People's Rep. of China
Phone: +86 28 8620-0175,
Fax: +86 28 8620-0177
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy China Ltd.
Guangzhou Representative Office
Room 1119-20, Dongshan Plaza
69 Xianlie Road Middle, Dongshan District,
Guangzhou, Guangdong, 510095,
People's Rep. of China
Phone: +86 20 8732-1281,
Fax: +86 20 8732-1286
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy Japan K. K.
Bois Hongo Building 6Fl, 3-42-5 Hongo
Bunkyo-ku, Tokyo 113-0033, Japan
Phone: +81 3 5802-3966,
Fax: +81 3 5802-3988
E-Mail: info@karlstorz.co.jp

www.karlstorz.com

STORZ
KARL STORZ — ENDOSKOPE

KARL STORZ GmbH & Co. KG

Mittelstraße 8
78532 Tuttlingen

Postfach 230
78503 Tuttlingen
Germany

Telefon: +49 (0)7461 708-0
Telefax: +49 (0)7461 708-105
E-Mail: info@karlstorz.de
Web: www.karlstorz.com