


**BW 485/L**

**Blut- und Infusionswärmer  
Blood and infusion warmer**



- anwenderfreundlich
  - mikroprozessorgesteuert
  - hohe Durchflußraten
  - geringer Energieverbrauch
  - Untertemperaturalarm
- 
- user friendly
  - microprocessor controlled
  - high transfusion rates
  - minimized energy consumption

**Service Manual**

**EG  
(DR) 5**



# 1. CONTINUOUS TEST

## 1.1 General

- Start-up
- After approx. 15 min. measure the temperature at the heat exchanger
- After 24 h measure the temperature again

During this test the initial and final values must be within the tolerance range and the measurements must be taken under the same ambient conditions.

The final temperature may deviate from the initial value by  $\pm 0.5^{\circ}\text{C}$ .




Fig. 1.1.1 Measurement set-up

- 1 IR thermometer
- 2 Heat exchanger (BW ring)

Measuring instrument used: HORIBA infrared thermometer IT 330  
 Accuracy:  $\pm 1\%$ ,  $\pm 1$  digit at full scale  
 Emission: black body

| Date           | Created by | Release | File                      | Revision no. | Page |
|----------------|------------|---------|---------------------------|--------------|------|
| February, 1997 | Riedl      |         | BW 485 Anweisung Englisch | 01           |      |

## Temperature Adjustment

- Connect adapter box to computer and start terminal program (e.g. Procomm)
- Connect adapter box to BW 485 CON2 (on PCB)
- Set switch to „TEST“ position. Device then emits an audio signal at regular intervals
- Preset spindle trimmer:

|    |                  |                         |
|----|------------------|-------------------------|
| P3 | Gain 1           | right-hand end position |
| P6 | Gain 2           | left-hand end position  |
| P5 | Alarm limit 41°C | left-hand end position  |
| P4 | Set temperature  | right-hand end position |

- Switch on BW 485
- Set temperature controller (P4) to required temperature (display on PC: „Set“). For fast presetting, press „SET“ key on adapter box in order to display set temperature on bar chart.
- Set current externally measured heater temperature with spindle trimmer P6 (display on PC: „SENS2“)
- Turn spindle trimmer P3 slowly to left until device starts to heat up. When heating is on L13 is illuminated. Slowly increase heater temperature to set value.

The two temperatures displayed on the PC as „Sens1“ and „Sens2“ must be as close to each other as possible ( $\pm 0.3$ ) otherwise an alarm is triggered.

The difference between the two temperatures is displayed on the PC under „Diff1“. Because measurements are taken at different points, the internally and externally measured temperatures may deviate slightly from one another.

### Setting the high temperature alarm limit:

1. Set switch in adapter box to „Test“ position
2. Set temperature to 41,5°C with P4 and wait until temperature is constant
3. Turn P5 s l o w l y to right until relays drop out and alarm is triggered
4. Let device cool down. Heat it up again and check cut-off point again
5. Increase set temperature (39°C) s l o w l y to 41°C and control alarm shut off
6. Reset at previous temperature

## 4. Temperature Adjustment

### Service interface:

The BW 485 is equipped with a service interface. This serial interface outputs a range of data which can be very useful for trimming and trouble shooting.

Connect the BW485 to your PC or terminal with the adapter box. The data can be displayed using any common terminal program e.g. Procomm.

Data format: 4800,n,8,1

### Key to serial output data:

|        |  |
|--------|--|
| Set    | Set value  |
| Sens2  | Temperature measured by Sensor2 (temperature signal for display and 41°C safety cut-off) |
| Sens1  | Temperature measured by Sensor 1 (control temperature)                                   |
| Diff   | Difference between set and actual (only when actual < set, otherwise 00.0)               |
| P%     | Current heating output in % of maximum power-rating                                      |
| Diff 1 | Difference between Sens1 and Sens2   |
| Err    | Error counter  |
| t      | Time until low temperature alarm is activated (0=on)                                     |


### Error codes:


The BW485 indicates various errors and displays them on the bar chart as follows:

|          |  |
|----------|--|
| ■□□□□□□□ | Low temperature alarm (1 <sup>st</sup> LED flashes)    |
| □■□□□□□□ | Temperature difference between Sens1 and Sens2 > 0.5°C |
| □□■□□□□□ | Sens1 > 41°C   |
| □□□■□□□□ | Sens2 > 41°C   |

BW 485

**BIEGLER**  
MEDIZINELEKTRONIK






BW 485

**BIEGLER MEDIZINELEKTRONIK**

|              |              |
|--------------|--------------|
| Printv. VI 2 |              |
| No. 100:     | No. 101:     |
| 10.12.1996   | F.NETAUSCHKE |



## PARTS LIST

|                     |   |   |
|---------------------|---|---|
| BW 485 / L (385 V7) | 1 | 3 |
|---------------------|---|---|

| PART            | NO. | DESCRIPTION                   | PART NO.   | ITEM       |
|-----------------|-----|-------------------------------|------------|------------|
| diode           | D1  | 1N4148                        | BD 1004148 | B3         |
| diode           | D2  | 1N4148                        | BD 1004148 | C3         |
| diode           | D3  | 1N4148                        | BD 1004148 | B3         |
| diode           | D4  | 1N4148                        | BD 1004148 | B3         |
| diode           | D5  | 1N4148                        | BD 1004148 | B3         |
| RFI capacitor   | C6  | 1.40.00MKTX2 0,47 $\mu$ F     | BB 3011447 | E1-2       |
| rectifier       | GL1 | DF04                          | BE 1000004 | C4         |
| LED             | L12 | $\varnothing$ 3mm red         | BD 4020004 | B4         |
| LED             | L13 | $\varnothing$ 3mm red         | BD 4020004 | C1         |
| LED             | L1  | $\square$ 5mm green           | BD 4205300 | A1-2       |
| LED             | L2  | $\square$ 5mm green           | BD 4205300 | A2         |
| LED             | L3  | $\square$ 5mm green           | BD 4205300 | A2         |
| LED             | L4  | $\square$ 5mm green           | BD 4205300 | A2         |
| LED             | L5  | $\square$ 5mm green           | BD 4205300 | A2         |
| LED             | L6  | $\square$ 5mm green           | BD 4205300 | A3         |
| LED             | L7  | $\square$ 5mm green           | BD 4205300 | A3         |
| LED             | L8  | $\square$ 5mm green           | BD 4205300 | A3         |
| LED             | L9  | $\square$ 5mm green           | BD 4205300 | A3         |
| LED             | L10 | $\square$ 5mm orange          | BD 4305300 | A3-4       |
| microcontroller | U1  | microcontroller ST62T15       | AG 1000625 | A1-2, B1-2 |
| optocoupler     | U2  | MOC3041                       | AI 1003041 | C1-D1      |
| oper. amplifier | U10 | ICL7621DCPA                   | AT 1007621 | A1         |
| oper. amplifier | U7  | ICL7611DCPA                   | AT 1007611 | B1         |
| PCB             |     | SKM                           | JP 1204857 |            |
| quartz          | Q1  | 4 MHz                         | AJ 1000046 | A3-4       |
| relay           | RE2 | JW2SN                         | CB 1212000 | B2,C2      |
| relay           | RE1 | JW2SN                         | CB 1212000 | B2-3,C2-3  |
| capacitor       | C8  | 15 pF RM2,5                   | BB 1010151 | A3         |
| capacitor       | C9  | 15 pF RM2,5                   | BB 1010151 | A3         |
| capacitor       | C7  | electrolytic 470 $\mu$ F/40 V | BC 1154477 | B4,C4      |
| capacitor       | C10 | foil 100 nF                   | BB 1020015 | B2         |
| capacitor       | C12 | foil 100 nF                   | BB 1020015 | B3         |
| capacitor       | C13 | foil 100 nF                   | BB 1020015 | C4         |
| capacitor       | C14 | foil 100 nF                   | BB 1020015 | C1         |
| capacitor       | C15 | foil 100 nF                   | BB 1020015 | C1         |
| capacitor       | C3  | foil 100 nF                   | BB 1020015 | B1         |
| capacitor       | C4  | foil 100 nF                   | BB 1020015 | A1         |
| capacitor       | C5  | foil 100 nF                   | BB 1020015 | B1,C1      |
| capacitor       | C11 | tantalum 4 $\mu$ 7/35 V       | BC 3302475 | C3         |

| Datum     | Erstellt | Freigabe | Datei                   | Änderungsindex | Seite |
|-----------|----------|----------|-------------------------|----------------|-------|
| Dec, 1996 | Riedl    |          | BW485 teillist englisch | 01             | 1     |



## PARTS LIST

|                     |   |   |
|---------------------|---|---|
| BW 485 / L (385 V7) | 2 | 3 |
|---------------------|---|---|

| PART     | NO. | DESCRIPTION   | PART NO.   | ITEM   |
|----------|-----|---------------|------------|--------|
| resistor | R25 | 100 $\Omega$  | BA 1020101 | C1     |
| resistor | R14 | 10k $\Omega$  | BA 1020103 | A4     |
| resistor | R29 | 10k $\Omega$  | BA 1020103 | B4, C4 |
| resistor | R30 | 10k $\Omega$  | BA 1020103 | B3, C3 |
| resistor | R34 | 10k $\Omega$  | BA 1020103 | A1     |
| resistor | R22 | 1k $\Omega$   | BA 1020102 | A4, B4 |
| resistor | R24 | 1k $\Omega$   | BA 1020102 | A2     |
| resistor | R26 | 1k $\Omega$   | BA 1020102 | C3     |
| resistor | R28 | 22k $\Omega$  | BA 1020223 | B1     |
| resistor | R32 | 22k $\Omega$  | BA 1020223 | A1, B1 |
| resistor | R11 | 27 $\Omega$   | BA 1020270 | D1     |
| resistor | R12 | 270k $\Omega$ | BA 1020274 | B1     |
| resistor | R20 | 270k $\Omega$ | BA 1020274 | A1     |
| resistor | R15 | 27k $\Omega$  | BA 1020273 | B1     |
| resistor | R16 | 27k $\Omega$  | BA 1020273 | A1     |
| resistor | R17 | 27k $\Omega$  | BA 1020273 | A1, B1 |
| resistor | R18 | 27k $\Omega$  | BA 1020273 | A1     |
| resistor | R23 | 2k2 $\Omega$  | BA 1020222 | B1     |
| resistor | R31 | 2k2 $\Omega$  | BA 1020222 | A1     |
| resistor | R13 | 330 $\Omega$  | BA 1020331 | D1     |
| resistor | R1  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R10 | 470 $\Omega$  | BA 1020471 | A3     |
| resistor | R2  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R3  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R4  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R5  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R7  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R8  | 470 $\Omega$  | BA 1020471 | A2-3   |
| resistor | R9  | 470 $\Omega$  | BA 1020471 | A3     |
| resistor | R6  | 470 $\Omega$  | BA 1020471 | A2     |
| resistor | R21 | 680 $\Omega$  | BA 1020681 | C4     |
| resistor | R33 | 8k2 $\Omega$  | BA 1020822 | A1     |

| Datum     | Erstellt | Freigabe | Datei                   | Änderungsindex | Seite |
|-----------|----------|----------|-------------------------|----------------|-------|
| Dec, 1996 | Riedl    |          | BW485 teillist englisch | 01             | 2     |

## PARTS LIST

BW 485 / L (385 V7)

3

3

| PART              | NO.  | DESCRIPTION      | PART NO.   | ITEM                |
|-------------------|------|------------------|------------|---------------------|
| pin-connector     | CON1 | 8p SPOX 5281 08A | CC 1852818 | C1, D1              |
| pin-connector     | CON2 | 6p/2,54          | CC 1000001 | B4, C4              |
| pin-connector     | CON4 | 3p SPOX 5281 03A | CC 1352813 | D1, E1              |
| pin-connector     | CON5 | 4p SPOX 5281 08A | CC 1452818 | D2-3                |
| pin connector     | JP7  | 3p/2,54          | CC 1000002 | B1                  |
| fuse              | S11  | 1,6 A T          | CE 1201600 | E2-3                |
| fuse holder       | S11  |                  | CE 9510062 | E2-3                |
| fuse holder       | S11  |                  | CE 9510062 | E2-3                |
| fuse              | S12  | 1,6 A T          | CE 1201600 | D2-3                |
| fuse holder       | S12  |                  | CE 9510062 | D2-3                |
| fuse holder       | S12  |                  | CE 9510062 | D2-3                |
| fuse              | S13  | 250 mA T         | CE 1200250 | C3                  |
| fuse holder       | S13  |                  | CE 9510062 | C3                  |
| fuse holder       | S13  |                  | CE 9510062 | C3                  |
| voltage regulator | U3   | LM78L05          | AB 1034006 | B2                  |
| voltage regulator | U4   | LM7805           | AB 1034105 | C3                  |
| spindle trimmer   | P3   | 3296 1k          | BA 7120102 | B1-2                |
| spindle trimmer   | P4   | 3296 20k         | BA 7120203 | B3                  |
| spindle trimmer   | P5   | 3296 20k         | BA 7120203 | B1-2                |
| spindle trimmer   | P6   | 3296 1k          | BA 7120102 | B1-2                |
| buzzer            | SU1  | KPE-214A         | CH 1214001 | A2,B2-3             |
| supressor-diode   | D6   | P6KE 15A C549    | BD 6615549 | C4                  |
| triac             | TR1  | BTA 08-700       | AA 3008700 | D1, E1              |
| driver            | U5   | ULN2003A         | AT 1002003 | A2,B2               |
| driver            | U6   | ULN2003A         | AT 1002003 | A2-3,B2-3           |
| transformer       | TR2  | Myra 9V 4VA      | BG 1002094 | C3-4, D3-4,<br>E3-4 |

| Datum     | Erstellt | Freigabe | Datei                   | Änderungsindex | Seite |
|-----------|----------|----------|-------------------------|----------------|-------|
| Dec, 1996 | Riedl    |          | BW485 teillist englisch | 01             | 3     |

EC Certificate of Conformity

No. 0001 V1 / 08.96

Authorized agent: E. BIEGLER G.m.b.H.  
Medical Electronics

Address: Allhangstrasse 18 a  
3001 Mauerbach  
Austria

Manufacturer: E. BIEGLER G.m.b.H.  
Medical Electronics  
3001 Mauerbach  
Austria

Product description: BW 485  
from serial no. V7 2001

BW 485L  
from serial no. V7 7001

The product described above is a medical product (as per 93/42/EEC Section 1, Para. 2a) and conforms to the provisions of the following European Directives:

93/42/EEC Council Directive on the harmonization of statutory provisions for medical products.

E. Biegler G.m.b.H.

Mauerbach, Austria,


28.05.97

(Date)

F. Netauschek/Development Manager  
(Name, Position)

  
(Signature)




I. Biegler/Managing Director  
(Name, Position)

  
(Signature)

This document certifies conformity with the above directives but does not constitute warranty of a quality as defined by the Austrian product liability law.


The safety regulations in the product documentation supplied must be observed.

**STAATLICH AUTORISIERTE  
VERSUCHSANSTALT FÜR RADIOTECHNIK  
PRÜFSTELLE FUNKENTSTÖRUNG PFE und FUNKGERÄTE-SICHERHEIT PFS**

| <b>Mitteilung von Prüfergebnissen</b>   |   |  |         |
|---|---|--|---------|
| geprüft im Auftrag von  | Biegler Medizin Elektronik<br>Allhangstr. 18a<br>3001 Mauerbach                     |  |         |
| Prüfgut   | <b>Blut- und Infusionswärmer<br/>BW 485 (L)</b>                                     |  |         |
| Eingelangt am   | 22. August 96   | Auftrags-Nr.   | 24013/R |
| Beendet am  | 22. August 96   |  |         |
| Geprüft nach  | EN 60601-1-2  |  |         |
| <br>_____<br>Ing. Thun<br><br>Sachbearbeiter |  | <br>_____<br>Dipl. Ing. Dr. techn. K. Tögl<br><br>Gutachter |         |
| <b>Test Report</b>  |   |  |         |



# QUALITY ASSURANCE



CERTIFICATE OF COMPLIANCE OF THE  
QUALITY SYSTEM WITH  
STANDARDS MENTIONED BELOW

Certificate No.: Q1 95 03 14553 005

The Certification Body for Quality Systems  
of TÜV PRODUCT SERVICE GMBH certifies that

**Biegler GmbH**  
**Allhangstr. 18a**  
**A-3001 Mauerbach**

in the facility:  
**Biegler GmbH**  
**Allhangstr. 18a**  
**A-3001 Mauerbach**

for the following area:

Ventilators and respiratory training apparatus,  
Blood and infusion warmers,  
Fibrin sealant application systems,  
Cooling systems for medical solutions,  
Olfactometers, Electrical stimulation therapy devices  
Laboratory and diagnosis equipments

has established and is operating a quality system  
which meets the requirement(s) of

ISO 9001 : 1994 / EN 46 001 : 10/93

as documented in the audit report no. 01/338-5-003/02

This certificate is valid until March 1998,  
provided that a periodical surveillance is conducted.

Munich, March 24, 1995

**TÜV PRODUCT SERVICE GMBH**  
**ACCREDITED CERTIFICATION BODY**  
for Quality Systems



*Z. B. Leue*

TÜV PRODUCT SERVICE GMBH, Ridlerstrasse 31, 80339 Munich, Germany

