Service Instructions

DrägerService
Medizintechnik
Medical Systems

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11.87 Lenke No.1

Incubator 8000 service strategy

Service at customers by Inspection Service engineer: - Checking of unit in accordance with Test Card and fault list
- Replacement of assemblies/p.c. board replacement if necessary

Service at Branch/Agency: - P.c. board replacement
- Repair at component level following corresponding training.

DrägerService


Detlev Lenke

11.87 Lenke No.2

Shifting of motor capacitor

As a result of unfavourable dimensional tolerances, the 68 μ F motor capacitor of the unit may make contact with the valve block and cause a short-circuit in the 24 V circuit. The 24 V circuit supplies the valves, the load relays for the heaters, the safety relay and the motor with voltage. In the event of such a fault, the incubator would assume the INOP status and give an audio and visual alarm; there is thus no danger for the patient.

When you next have access to the incubator, please move the 68 μ F motor capacitor towards the fan motor (see sketch).

DrägerService
Medizintechnik
Medical Systems

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Following parts will be sent free of charge and referring to device to each subsidiary from VM-DS-PB:

- 1 Conversion instruction
- 1 EPROM 83 00963
- 1 Operating instruction, German, 90 27300, 5th edition/January'88
- 1 Diode 83 00251 (1N5625) for reduction of the vaporizer performance
- 1 Shrink hose
- 2 Hinge assembly 2M 19542
- 3 Sealing ring 2M 8777 for water bottle holder
- 1 Sealing CH 9039 for water pipe
- 1 Ring M 16206

- 2 Diodes 83 00021 (1N4148) Horn suppression
- 1 Capacitor 18 01376 (Tantalum 10µF)
- 1 Spring V 4647 for brief operating instruction
- 1 Metal screw B 4.8 x 9.5 for water bottle holder
- 1 Steel-water pipe 2M 19619 (does not concern all devices)
- 1 Hose 2M 19620 (- - -)
- 1 Screw 2M 17288 (Hose holder, does not concern all devices)
- 1 Report form
- 1 Description of software modification

Additionally with costs if exported

- 1 Folio keyboard for skin temperature
- 1 P.c. board Display
- 1 Enclosure, English, to the operating instruction (only for English-speaking countries)

The possible difference between delivered Incubators and Incubators converted with skin temperature regulation is separately charged by the sales department.

We would like to ask you to perform the conversion until the end of march.


DrägerService

Dräger *Lenke*

Faulty O₂ calibration

02.88 Lenke No. 4

Cause: When effecting calibration with high oxygen concentrations inside the Incubator, oxygen-enriched air from the interior may flow past the O₂ sensors if there is a gap between the climate sensor and the metal block below the climate sensor. This will bias the calibration value for 21 vol. % O₂.

Precise working instructions:

1. Loosen the two socket-head cap screws and fold down the flap beneath the electronics unit.
2. Detach the cable and pneumatics connections and then remove the electronics unit.
3. Strip incubator down to unit module.
4. Unscrew water pipe to unit module.
5. Detach power cord at terminals (bottom left) on p.c. board "unit".
6. Loosen the two screws (from above) and then remove the unit.
7. Move the capacitor holders towards the fan motor and secure with 2 tie bands.
8. Installation is to be performed in the reverse order.

This measure affects the first 25 incubators with the serial nos. 2043-1 to 2043-25, delivered as of the 42nd calendar week.

DrägerService

**Conversion**

01.88 Lenke No.3

All 170 Incubators 8000 delivered up to now are converted to latest state. The conversion contains following items:

- New software (version no. 0.3) with improved regulation for temperature and humidity, possibility of skin temperature regulation and changed warning philosophy. There are 3 additional Service-modes for the Service:
 - Mode 6 = Measuring value of the air temperature sensor.
 - Mode 7 = Measuring value O₂-sensor A (outer sensor)
 - Mode 8 = Measuring value O₂-sensor B (inner sensor)The valid range for 21 Vol.% of the sensors is 306 - 812
- New operating instruction (German)
- New hinge assembly for front flap
- Reduction of vaporizer performance by an additional diode in the heating circuit (concerns 75 devices)
- Set regulating thermostat of the vaporizer from 135°C to approx. 120°C, in order to prevent the non-desired switching off of the vaporizer-heating by the safety thermostat.
- New steel water pipes with hose for devices which were delivered still with alloy water pipe (conversion concerns 75 devices)
- Horn suppression in 10-minutes-check by soldering in 2 additional diodes and one capacitor on the p.c. board "power supply" (concerns 75 devices)
- Slimming down a housing edge using a file in the interior of the Incubator to prevent motor noise
- Additional sealing rings and conversion of the bottle holder to prevent water leakages.
- Stronger spring for the flap of the brief operating instruction
- Retrofitting of the skin temperature regulation
- Retrofitting of the screwing on of the hose holder (concerns 75 devices)

DrägerService
Medizintechnik
Medical Systems

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cont. of IDM No. 4

Remedy: Set the gap between the climate sensor and the metal block to a minimum. This may involve enlarging the two holes in the Plexiglas hood using a needle file.

Please check the climate sensor when you next visit hospitals. Only a few units are affected by this measure.

DrägerService

Dieter Lenke

INOP-LED

02.88 Lenke No. 5

Reason: Lightening up of INOP-LED without audio alarm and during normal operation of the unit.


Cause: Unfavourable tolerances of several elements on the p.c. board "motherboard".

Remedy: In the event of a fault, the transistor V 27 on the p.c. board "motherboard" is to be replaced with a different type of transistor.



old: BCY59/IX

new: 2N3904 (Order No. 1801805)

BCY59:

2N3904:

The resistance of R52 and R53 are likewise to be changed:

R52: old 14.7 k Ω

new 7.5 k Ω (order No. 1806882)

R53: old 38.3 k Ω

new 10.0 k Ω (order No. 1805339)

This modification has been incorporated into series production. General retrofitting is not necessary.

DrägerService

Dieter Lenke

INOP due to CPU reset

02.88 Lenke No. 6

Reason: The Incubator sporadically assumes the INOP status. The Incubator usually functions normally again once it has cooled down and been switched on again.

Cause: The 5 V voltage monitoring on the p.c. board "CPU" is incorrectly calibrated and, although there is no undervoltage, triggers a reset if the electronics heat up. A particular characteristic feature of this fault is that the LED on the p.c. board "CPU" lights up.

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
cont. of IDM No. 6

Remedy: Re-calibrate 5 V voltage monitoring on p.c. board "CPU".

Tools and aids required: Voltmeter, small screwdriver for potentiometers and locking compound.

Working instructions:

- Fold down flap beneath electronics module after loosening the two hexagon socket-head screws.
- The p.c. boards required for calibration are now accessible from below.

- Using potentiometer R19 on the p.c. board "power pack" adjust the supply voltage - measured between TP1 (+) and TP2 (GND) on the p.c. board "CPU" - from approx. 5 V to $4.75 \text{ V} \pm 0.02 \text{ V}$. Caution: Never adjust to a setting in excess of 5.5 V!!!
- Two initial conditions are possible for calibration of the voltage monitoring on the p.c. board "CPU":
 - a) If the LED V11 lights up, turn potentiometer R1 in a clockwise direction until the LED V11 goes out or only lights up faintly. Then turn in the opposite direction until the LED V11 just lights up again.
 - b) If the LED V11 does not light up or only lights up faintly, turn potentiometer R1 in an anticlockwise direction until the LED V11 just lights up.
- Using potentiometer R19 (p.c. board "power pack"), now slowly increase the operating voltage (caution: never to a setting in excess of 5.5 V!) until the LED V11 goes out or lights up faintly.
Test value: $U = 4.89 \text{ V} \pm 0.03 \text{ V}$.
- Set operating voltage with R19 (p.c. board "power pack") to $5.0 \pm 0.05 \text{ V}$.
- Secure both potentiometers with locking compound.


This measure affects just a few units from the serial number range 2030-XX to 2080-XX. There is no need to re-calibrate units where this fault has not yet occurred. Repair units from the a.m. range of serial numbers must always be re-calibrated.

Castors for Incubator 8000

Reason: Fracture of castor pins

Cause: Bending stress caused by slipping of block 2M 19266

Remedy: To reduce the bending stress, the block 2M 19266, which previously was only placed in position, should be bonded to the trolley and thus prevented from slipping. Bonding is to be effected by using double-sided adhesive tape "Tesafix 4968". This requires a piece 50 x 35 mm with a hole of approx. 12 mm dia. in the centre. The threaded pin of the castor is to be secured with Loctite 221.

This work applies to all Incubators up to serial no. 2140-0000 and also to 2140-0008. The necessary work has already been performed on 2130-0011, -0022, -0023, -0024, -0028 and -0030 in Lübeck. All Branches/Subsidiaries will automatically be receiving sufficient "Tesafix 4968" from VM-DS-PB. The work outlined above should be performed at the next opportunity or on the occasion of the conversion campaign II (new water jar mount and noise attenuation for swivel windows).

DrägerService

E. W. Seeger *Lenke*

Retrofitting of new swivel-window seal

Reason: The old swivel-window seal does not make allround contact with the swivel window.

Cause: The swivel-window seal is round and becomes corrugated when fitted in the non-round opening in the hood.

Remedy: Pre-shaped seals:

Seal, left 2M 19469

Seal, right 2M 19470

The new seals are white (previously: transparent)

This retrofitting work affects all Incubators up to serial no. 2140-0000 and also 2140-0008. New seals have already been fitted on 2130-0011, -0022, -0023, -0024, -0028 and -0030.

DrägerService
Medizintechnik
Medical Systems

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Continuation IDM no. 8

All Branches/Subsidiaries will automatically be receiving sufficient quantities of new seals from VM-DS-PB (left 2M 19469, right 2M 19470).

This retrofitting work should be performed at the next opportunity or on the occasion of conversion campaign II (new water jar mount and noise attenuation for swivel windows).

The urgency of this work is to be agreed upon with Sales.

DrägerService

B. W. Sauer *Lenke*

3.88 Lenke No. 9

Temperature safeguard for air heater

Repair of the temperature safeguard 18 23 809 is not envisaged. In the event of a fault a new temperature safeguard is always to be fitted.


DrägerService

B. W. Sauer *Lenke*

6.88 Lenke No. 10


Regulating thermostat for water deficiency signal

The adjustable capillary thermostat for water deficiency signal is replaced by the thermostat $\pm 15^{\circ}\text{C}$ (18 25 992). A conversion is always then necessary if in case of repair a new heater (subassembly) 2M 19 709 (aluminium block and heating cartridge) is fitted.

All units are prepared for the installation of the new thermostat. During installation attention is absolutely to be paid to the following:

- Screw in thermostat only lightly by hand
- The connecting cables have to be coated with a PVC-hose (per m 11 97 797), see diagram) in the range of the air heating and the mains filter because of the 4 KV-potential voltage.

DrägerService
Medizintechnik
Medical Systems

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Continuation IDM no. 10

The overlap is bigger with the new thermostat than with the old one. After releasing of the water deficiency signal and subsequent filling of the water supply several minutes may pass, until the water deficiency signal switches off. This will be incorporated into the new Test Card.

DrägerService

B. Reyer *J. Lenke*

6.88 Lenke Nr. 11

Conversion campaign II

All Incubators 8000 which were delivered up to now, are retrofitted to latest state during this last retrofitting campaign.

The retrofitting includes the following items:

- new water jar mount
- noise attenuation for swivel windows

Each subsidiary/agency will receive free of charge and referring to device from VM-DS-PB following parts:

- 1 x Conversion instruction
- 1 x Water jar mount 2M 19 806
- 1 x Lock washer for water distributor
- 4 or 6 x Retainer with buffer (according to number of swivel windows)

We like to ask you to agree upon with sales the urgency of this work and to perform this retrofitting campaign, including the measures in accordance with IDM no. 7 and 8, until the end of October at the latest.

DrägerService

B. Reyer *J. Lenke*

Retrofitting Package III on Incubator 8000

Retrofitting work consists of the following:

1. Improved cooling behaviour by means of additional fan. Whenever the air heating is not switched on, an additional fan is switched in - irrespective of the processor - which dissipates excess heat from the electronics. The installation of such an additional fan does not lead to an increase in the noise level inside the incubator.
2. New software (EPROM) and larger data memory (8K-RAM). The most important improvements brought about by the new software are as follows:
 - Improved control loops for all parameters without permanent, steady-state system deviations.
 - No calibration prompt for O₂ measurement function following power failure.
 - Switching of valves in 10 minute test dispensed with.
 - No evaluation of difference voltage between the two O₂ capsules as criterion for used-up sensor capsule. Different O₂ capsules are also accepted, as long as the voltages of the individual sensors are within the specified voltage window when calibrating with 21 vol. %.
 - Monitoring of O₂ concentration following switch-off of O₂ module. If the O₂ concentration then increases by 3 vol. % within 3 minutes, this results in an INOP alarm (fault 67) with additional flashing of the instantaneous O₂ value.

- O₂ setpoint adjustment expanded to 75 vol. %.
- The humidity control is active even under 27° C and with warning + 1.5° C air temperature.
The humidity display range is now 0 % to 99 %. The instantaneous humidity level is indicated in the additional service mode 9.
- The initial air-temperature setpoint after switch-on of the incubator is now 33.0° C.
- ~~RAM test expanded for 8K-RAMs. 8K-RAMs (previously 2K-RAM)~~
are mandatory with the new software.
- Alarms for setpoint deviations (± 1.5° C for air temperature, ± 0.5° C for skin temperature, ± 5 vol. % O₂) are only given if the permissible deviation is exceeded. Previously, the alarm was given if the permissible deviation was reached.
- Following switch-on of the incubator or one of the modules, the lettering "SEt" flashes alternately with the instantaneous actual-value display (with O₂ module following calibration) as a prompt for acknowledging the setpoints. Acknowledgement is effected by pressing the button "+ setpoint" or "- setpoint" in the respective module. If no acknowledgement is given, the module continues to function as before; the flashing of "Set" is merely informative.
- The following initial conditions apply when switching from air to skin temperature:
 - a) Cold start (initial warm-up, audio alarm tone suppression -1.5° C is still active)
The instantaneous skin-temperature actual value is taken over as setpoint. An air-temperature setpoint of 33.0° C is set as zero crossing for the skin-temperature characteristic curve.

b) Warm start

The instantaneous skin-temperature actual value is taken as setpoint. The instantaneous air-temperature set value is taken as zero crossing for the skin-temperature characteristic curve.

It is therefore advisable to first set the correct temperature via the air temperature and then to switch to skin-temperature regulation. No or only minor setpoint changes for the skin temperature are then necessary and the skin-temperature characteristic curve is automatically in the right temperature range.


As regards the skin-temperature sensor alarm, the following are set irrespective of starting mode (cold or warm):

A skin-temperature setpoint of 36.3° C when switching, a setpoint of 35.0° C in the case of skin-temperature actual values < 35.0° C and a setpoint of 37.0° C for skin-temperature actual values > 37.0° C.

- When switching from skin to air-temperature regulation, the air temperature specified by the skin-temperature controller is taken as setpoint for the air-temperature regulation.
- After swivelling in the climate sensor, the heating is switched off for 90s and the audio alarm $\pm 1.5^{\circ}$ C is suppressed for 10 minutes.

3. Renewal of evaporator connection (see IDM No. 14)

4. Wind insulating tape around connecting cables of keypads to provide protection against moisture and mechanical damage. This is no longer necessary with newer incubators where the connecting cables are laid internally between keypad and p.c. board "display" and thus do not project.

5. Assorted new miscellaneous items such as 2 new stickers, replacement of a screw at the climate sensor hinge and a new hose for the evaporator.


The following tools and auxiliary items are required for retro-fitting purposes:

- various screwdrivers and hexagon-socket wrenches
- 8 mm pipe wrench
- caliper gauge
- drill with 3.5 mm bit
- soldering iron and tin-lead solder
- Wacker Silicon E 41 transparent or heat-sealing adhesive
- various tie bands
- waterproof fibre-tipped pen

Retrofitting instructions

These retrofitting instructions describe the retrofitting sequence.

1. Note down data from service modes 1, 2 and 3.
(For description refer to Test Card, item 3.). The subsequent replacement of the RAM would otherwise cause the data to be lost.
2. Fold down flap beneath slide-in electronics module after loosening the 2 screws or the 4 screws in the case of newer models.
3. Remove slide-in electronics module.
4. Remove p.c. board "CPU" from slide-in electronics module.
 - Replace EPROM
 - Replace RAM
 - Reconnect jumper next to RAM

5. Unscrew left side panel (next to p.c. board "CPU") from electronics slide-in module.
 6. Drill side panel with the aid of the drilling template 82 90 387 in accordance with drawing "conversion kit, p.c. board fan" 82 90 390. Alter type number of electronics slide-in module on side panel:

82 00 901	=	82 90 401	or
82 00 902	=	82 90 402	
 7. Screw on p.c. board "fan" on inside of side panel.
-
8. Remove humidity keypad (only if connecting cable of keypad protudes) and wind insulating tape around connecting cable; in doing so, do not wind more than twice around connector. Then install humidity keypad again.
 9. Solder on wiring harness 82 90 389 on motherboard and on p.c. board "fan" in accordance with drawing "conversion kit, p.c. board fan" 82 90 390. To do so:
 - Scrape off protective lacquer (transparent) and solder resist (green) on conductor for connection of red and black cable and tin conductor.
 - Scrape off protective lacquer on IC connection for connection of green cable and tin.
 - Solder on wiring harness.
 - Fit side panel and solder on wiring harness on p.c. board "fan"; in doing so, secure soldering points with shrink-down tubing.
 - Secure soldering points on p.c. board "motherboard" with Wacker Silicon E 41 or hot-sealing adhesive.

DrägerService
Medizintechnik
Medical Systems

10. Fit CPU using module location in centre of p.c. board "motherboard".
11. Wind adhesive tape around connecting cables of other keypads (as Item 8).
Exception: Do not remove calibrated skin-temperature keypad (West Germany only) as otherwise the calibration seal would be broken.
Then install keypads again.

12. Remove unit.

13. Renew cap on evaporator of unit (old cap white, new cap black).

14. Unit fuses

Move F6 = 200 mA to position F5 and renew F6 with 500 mA.
Please also refer to "fuses" sticker in retrofitting set.


To date:

New:

F5 = 125 mA
F6 = 200 mA

F5 = 200 mA
F6 = 500 mA

Stick "fuses" sticker from inside onto panel beneath electronics slide-in module (see diagram)

15. Installation of additional fan

- Unscrew plate with connector for climate sensor.
- Unscrew sensor plug with wiring harness from plate.
- Screw sensor plug with wiring harness to new plastic section; in doing so, make use of the longer countersunk screws provided.
- Screw on plastic section with fan in place of plate; the black protective cap can be dispensed with.
- Re-attach any loosened tie bands.

16. Install unit. New seals and a new hose are provided in the retrofitting set for the water pipe.

17. Install slide-in electronics module and connect all cables. The plug connection of the fan is polarized.

On account of the checks to be performed in accordance with the new Test Card do not close the flap beneath the slide-in electronics module.

18. Renew hexagon socket-head screw at hinge of climate sensor. The new screw is made of stainless steel (does not apply to all units).

19. Assemble unit (with exception of flap beneath slide-in electronics module) ready for operation.


20. Check engaging function of swivel window

The swivel windows must close securely even in the event of slight closing pressure and subsequently pulling hard on the outer edge of the swivel window. If this function is not guaranteed:

DrägerService
Medizintechnik
Medical Systems

- Remove washer beneath rubber cone of swivel-window latch if there is one and check function again.
- If there is no washer beneath rubber cone of swivel-window latch, or if removing the washer does not lead to any improvement, gradually shorten rubber cone until secure engaging is guaranteed.

21. Stick on sticker "distilled water" next to jar holder.

22. Enter data from service modes 1, 2 and 3 into unit again (refer also to R in microfiche).
23. Check incubator in accordance with new Test Card.
24. Enter serial no. of incubator in new operating instructions (last page) and destroy old operating instructions.
25. Instruct customer on the basis of the description of software modifications given in these retrofitting instructions and the new operating instructions.

26. Document retrofitting of incubator and provision of instruction in customer's equipment manual (form D/1 or A/1). (Applies only to West Germany and West Belin)
 - Existing EPROM replaced with new EPROM 82 90 370
 - New operating instructions 90 27 341
 - Additional fan fitted
 - New sticker 2M 19 841
 - New sticker 2M 20 056
27. Fill in acknowledgement sheet for retrofitting work and send it to VM-DS-PB (absolutely necessary on account of need for registration).

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Abtlg. VM-DS-PB

Postfach 13 39

2400 Lübeck 1

West Germany

Acknowledgement sheet for retrofitting of Package III in
Incubator 8000

1. Customer's address:

2. Date on which retrofitting work performed:

3. Equipment data, rating plate:

Item No.:

Serial No.:

4. Equipment data following retrofitting:

Serial No. EPROM 82 90 370: _____

Serial No.,

operating instructions 90 27 341: _____

5. The old EPROM is to be destroyed.

6. Branch/Agency, Subsidiary, Sales Partner:

Date

Signature


5. Technical safety check
 To do so, fold down flap beneath slide-in electronics unit after loosening the two 2 or with newer units the 4 screws.

5.1. Protective conductor test
 Test points: Protective conductor terminal on left of unit, trolley, screws on height-adjustable column, O₂ connection socket for flowmeter, flap beneath slide-in electronics unit.


Test value: $R \leq 0.2 \Omega$ 0

5.2. Equivalent leakage current measurement
 For measurement purposes the switch-on relay and the safety relay must be jumpered.
 Caution: Detach mains connection beforehand!
 Without humidity module -> connection a
 With humidity module -&; connection a + b

p.c. board "unit" 8200600 - 3:

P.c board "unit" 8200600 - 4:

USE FIRST MEASURED VALUE AS A REFERENCE

SUBSEQUENT MEASUREMENTS MAY EXCEED 1ST MEASUREMENT BY 100%

5.2.1. Equivalent unit leakage current measurement
The subsequent measurements may exceed the first-measured value in the test card by max. 100%, however the stated value must not be exceeded.


DOES NOT EXCEED 0.75 mA

Test value: $I_A < 0.75 \text{ mA}$
Advisory: The first measured value is to be taken over, when a new test card is started..
First measured value: $I_A = \dots \text{mA}$

V

5.2.2. Equivalent leakage current measurement to the housing of the skin temperatur sensor.
(only with skin temperature regulation)

Connect shortened circuit mains plugs (L and N, not protective conductor) to socket "Ab1" and test cable to the socket "Com" of the VDE tester.
(with converted VDE testers socket "SL" instead of "Com").
Equivalent leakage current measurement to the contact surface of the skin temperature sensor.

Test value: $I_A = 0.0+0.1 \text{ mA}$

V

If the first measured value exceeds this value, the insulation of the sensor is faulty.

Important: After having performed measurement the connections of the soldering terminals are to be removed again.

5.3. Visual check of the mains voltage conducting cables.
The power cable and all visible connecting cables to the unit must be fixed safely.

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


Exchange ICs
install PCB with Loctite 221

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Back Casings of 8000
loosen plate with sensor plug
remove black cap
Mount plug on new fan plate (Loctite)
Put cable to right of transformer
install fan

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1. General scope

These testing and adjustment specifications apply to the following equipment:

- Slide-in module, air temperature, Incubator 8000 and air unit or air/water unit.
- Slide-in module, air and skin temperature, Incubator 8000 and air unit or air/water unit.

2. General description

2.1 Intended use, usage conditions

The air-temperature slide-in module is an assembly designed to regulate the air temperature in the Incubator 8000 as well as to optionally regulate the O₂ concentration and humidity level.

The air and skin temperature slide-in module also makes it possible to regulate the air temperature as a function of the child's skin temperature.

2.2 Schematic design, component block diagram

All electrical assemblies can be broken down in 5 basic components:

1. Unit

With mains cut-in relay with auxiliary transformer, mains transformer, fuses, fan motor, air heating with semiconductor relay, evaporator with semiconductor relay and thermostat for lack-of-water alarm (only with humidity regulation option), 2 safety relays and the connection for 2 O₂ valves.

2. Electronics slide-in unit

With p.c. board "CPU Standard 2", p.c. board "analog" (2 versions possible - > with and without skin-temperature measurement amplifier), p.c. board "power pack", p.c. board "motherboard", displays and keyboards for air temperature, skin temperature (skin-temperature regulation option), oxygen (oxygen regulation option), humidity (humidity regulation option) and the pneumatic pressure switches for electrical height adjustment (option).

Some As EUITA.

3. Climate sensor

With 2 air-temperature sensors, 2 microswitches for detection of sensor setting, 1 connection socket for skin-temperature sensor, 2 O₂ sensors (oxygen regulation option) and a humidity sensor (humidity regulation option).

4. Valve block

With 2 solenoid valves (oxygen regulation option) with fixed metering for 14 L/min. in each case and a check valve for external oxygen supply via a flowmeter.

5. Electrically adjustable height setting

With pneumatic actuation via pressure switch or fixed column.


Both versions feature the mains switch and an overload circuit-breaker in the column.

2.3 Description of Incubator 8000 block diagram

p.c. board designations:

- 1 Motherboard
- 2 Analog
- 3 Unit
- 4 Powerpack
- 5 Switch
- 6 Air display
- 7 Skin-temperature display
- 8 O₂ display
- 9 humidity display
- 10 Fan

- identical p.c. boards


Inkubator 81


VM-DS-PB-le-3525

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Mappen-Nr.
6141.20

Ausgabe
7.87


Block diagram incubator	
PIA1	PIA2
RAM/ROM	Timer
P.W.P. Failure	Relay
LED PWR-fail	Relay
Horn	Relay
Heater (Air)	Relay
Heater (H ₂ O)	Relay
Valve 1	Relay
Valve 2	Relay
O ₂ -valve	Relay
O ₂ -valve	Relay
Additional Fan	Relay
Blower	Relay
Feedback fan	Relay
I/O	Relay
Emergency Stop	Relay
Cut Out Water-level	Relay
Cut Out	Relay

All contact settings are shown in idle position

Dräger

Block diagram incubator

Components of Incubator 8000 with all options

Component designations:

D IC, digital
N IC, analog
V Individual semiconductor
C Capacitor
R Resistor
S Switch
K Relay
F Fuse
T Transformer

Definition of terms:

Block : Numbering as per block diagram
Action in normal operation : Function of block when there is no malfunction
Action in the event of malfunction : Function of (also not defective) block in the event of a fault
Indication : Indication of status

Logic-level definitions:

L = low (= 0 volts)
H = high (= +5 V or +15 V)
Y = high impedance or input
X = don't care

Establishment of "safe condition" in the event of malfunction:

Air heating off
H₂O heating off
O₂ valve 1 off
O₂ valve 2 off

SEE PAA 2.3
- main "UNIT."

Block No. 1 air heating 400 W (3)

- Action in normal operation : Air heating on/off
- Actuated by : Mains relay (3/K1)
Excess-temperature relay (3/K2), Semiconductor relay 1 (3/HLR1), Air-heating thermoswitch (3), Fuse (3/F1, F2)
- Action in the event of malfunction : Air heating off
- Actuated by : a) Switch-off semiconductor relay 1
- by means of PIA, see Item 13
- by means of elapsed monoflop, see Item 21
- by means of switch-off 24 V, see Item 17
b) switch-off by excess temperature, see Item 12
c) Switch-off by air-heating thermoswitch, see Item 9
- Indication : By way of air-heating current monitoring, see Item 8

Block No. 2 H₂O heating 200 W (3)

- Action in normal operation : H₂O heating on/off
- Actuated by : Mains relay (3 K1),
Excess-temperature relay (3/K2), Semiconductor relay 2 (3/HLR2), Thermoswitch H₂O (3), Fuse (3/F1, F2)
- Action in the event of malfunction : H₂O heating off
- Actuated by : a) Switch-off semiconductor relay 2
- by means of PIA, see Item 14
- by means of elapsed monoflop, see Item 21
- by means of switch-off 24 V, see Item 17

b) Switch-off by excess-temperature relay, see Item 12

c) Switch-off by H₂O thermo-switch, see Item 10

Indication : By way of indication from semiconductor relay 2, see Item 14

Block No. 3 O₂ valve 1

Action in normal operation : Valve on/off

Actuated by : PIA1-PB2 = H (on)
= L (off)
(1/D1, D7C, V5)

Action in the event of malfunction : Valve off

Actuated by : 1) PIA1-PB2 = L
(1/D1, D7C, V5)
2) Monoflop 1 (1/3A) or Monoflop 2 (1/3B) elapsed, see Item 21
3) Switch-off 24 V, see Item 17

Indication : PIA2-PB2 = H (on)
= L (off)
By means of electr. potential at coil (1/V6, D2)

Block No. 4 O₂ valve 2

Action in normal operation : Valve on/off

Actuated by : PIA1-PB3 = H (on)
= L (off)
(1/D1, D7D, V7)

Action in the event of malfunction : Valve off

Actuated by : 1) PIA1-PB3 = L
(1/D1, D7D, V7)
2) Monoflop 1 (1/3A) or
Monoflop 2 (1/3B) elapsed,
see Item 21
3) Switch-off 24 V, see Item 17

Indication: : PIA2-PB3 = H (on)
= L (off)
by means of electr. potential
at coil (1/V8, D2)

Block No. 5 fan (3)

Action in normal operation : Air circulation

Actuated by : 24 V = always runs with "power
on"

Indication : See Item 6, fan failure alarm

Block No. 6 fan failure alarm

Action in normal operation : PIA1-PA6 = L
(fan running)

Actuated by : Solenoid in impeller induces
pulses in coil of fan failure
alarm (3/coil) (1/V18, N1,
V13, D6A, V14, D1)

Action in the event of malfunction : PIA1-PA6 = H

Actuated by : No retrigger pulses for
monoflop (1/D6A), (1, R25,
C11, D6A, V14, D1)

Indication : 1) Fan running
PIA1-PA6 = L
2) Fan running and test:
PIA1-PA6 = H
3) Fan not running:
PIA1-PA6 = H

Test by : PIA2-CB2 = H
Suppresses fan pulses
(1/D2, V19)

Block No. 7 thermostwitch for lack of water 115 °C (3)

Action in normal operation : Water present:
PIA2-PA1 = L
Lack of water:
PIA2-PA1 = H

Actuated by : H₂O heating and water supply

Block No. 8 current monitoring heating air (3)

Action in normal operation : Heating on:
PIA2-PA7 = L
Heating off:
PIA2-PA7 = H

Actuated by : Current transformer
(3/RWP 209G),
Monoflop (1, N3)

Block No. 9 air thermostwitch (3)

Action in normal operation : "On"

Actuated by : T < 315 °C

Action in the event of malfunction : Air heating "off"

Actuated by : T > 315 °C

Indication : Current monitoring
Heating air, see Item 8

Block No. 10 H₂O thermostwitch (3)

Action in normal operation : "On"

Actuated by : T < 135 °C

Action in the event of malfunction : H₂O heating "off"

Actuated by : T > 135 °C

Block No. 11 mains switch (5/S1), auxiliary transformer (3/T54)

Mains relay (3/K1)

Action in normal operation : 1) Mains switch "on"

- a) $\overline{\text{NMI}} = \text{H}$
- b) Mains relay on (3/K1)
- c) Emitter (4/V21) on GND

2) On switch-off
NMI = L

3) Mains switch "off"

- a) Mains relay off (3/K1)
- b) Emitter (4/V21) open

Actuated by : Re 1)
a) (5/S 1.3) closed (4/V8, V9, V22)
b) (5/S 1.1) closed (3/T54, V2)
c) (5/S 1.2) closed
Re 2) (5/S 1.3) open, (4/V8, V9, V22)
Re 3)
a) (5/S 1.1) open
b) (5/S 1.2) open

Indication : 1) $\overline{\text{NMI}} = \text{H}$
PIA2-PB6 = H
and
PIA2-PB7 = H (on)

2) $\overline{\text{NMI}} = \text{L}$

PIA2-PB6 = L

and

PIA2-PB7 = L (off),

(5/S 1.3) open,

+ 15 V (4/V7, R4, V9, V22)

Block No. 12 excess-temperature relay (3/K2)

Action in normal operation : Enabling of mains voltage for
air heating and H₂O heating

Actuated by : Excess-temperature relay
bistable (1/K1.2), see Item 16

Action in the event of malfunction : Switch-off of mains voltage for
air heating and H₂O heating

Actuated by : Excess-temperature relay
bistable (1/K1.2) see Item 16

Indication : Current monitoring heating
air, see Item 8

Block No. 13 semiconductor relay 1 (3)

Action in normal operation : Air heating on/off

Actuated by : PIA1-PB0 = H (on)
PIA1-PB0 = L (off)

Action in the event of malfunction : Air heating off *

Actuated by : 1) PIA1-PB0 = L
2) Monoflop 1 (1/D3A) or
Monoflop 2 (1/D3B) elapsed,
see Item 21
3) Switch-off 24 V, see Item 17

Indication : Current monitoring heating
air, see Item 8

Block No. 14 semiconductor relay 2 (3)

- Action in normal operation : H₂O heating on/off
- Actuated by : PIA1-PB1 = H (on)
PIA1-PB1 = L (off)
- Action in the event of malfunction : H₂O heating off
- Actuated by : 1) PIA1-PB1 = L
2) Monoflop 1 (1/D3A) or
Monoflop 2 (1/D3B) elapsed,
see Item 21
3) Switch-off 24 V, see Item 17
- Indication : Potential at control input of
semiconductor relay 2
PIA1-PB1 = H (on)
PIA2-PB1 = L (off)

Block No. 15 power pack

- Action in normal operation : Supplies:
+ 5 V
+ 15 V
- 15 V
+ 24 V
+ 9 V (battery buffered)
- Action in the event of malfunction : 1) Switch-off of transformer
primary
2) Fusing of secondary
- Actuated by : Re 1)
a) 2-pole fuse (3/F1, F2)
b) Transformer fuse, primary
(3/F4)
c) Efen sensor in trans-
former (3)
- Re 2)
a) Fusing +5 V (3/F7)
b) Fusing +15 V (3/F6)
c) Fusing -15 V (3/F5)

Indication : Via voltage monitoring:
+24 V, see Item 17
+ 5 V, see Item 19
+/- 15 V, see Item 35
Mains voltage, see Item 23

Block No. 16 excess-temperature relay, bistable (1/K1A, B)

Action in normal operation : Activation of excess-temperature relay by closed contact (1/K 1.2), see Item 12

Actuated by : Leading edge of
PIA1-PB4 = L > H
via (1/C16), (1/V9) to
coil (1/K1A)

Action in the event of malfunction : Deactivation of excess-temperature relay via open contact of (1/K 1.2), see Item 12

Actuated by : 1) Excess-temperature comparator = H (2N5)
2) PIA1-PB6 = H
on coil (1/K1B)

Indication: Via 2nd contact of (1/K 1.1)
PIA1-PA0 = H (normal)
PIA1-PA1 = L (excess temperature)
Current monitoring heating air, see Item 8
Indication, semiconductor relay 2, see Item 14

Block No. 17 relay, monitoring 24 V (3/K4)

Action in normal operation : 24 V actuation voltage for:

- Semiconductor relay 1 (heating air)
- Semiconductor relay 2 (heating H₂O)
- O₂ valve 1
- O₂ valve 2

Actuated by : +5 V voltage monitoring, see Item 19

- Action in the event of malfunction : Disconnection of 24 V actuation voltage for:
- Semiconductor relay 1
 - Semiconductor relay 2
 - O₂ valve 1
 - O₂ valve 2
- Actuated by : a) +5 V voltage monitoring, see Item 19
- b) Failure of 24 V voltage supply
- Indication : a) Current monitoring, heating air, see Item 8
- b) Indication
Semiconductor relay 2,
see Item 14
- Test : By way of 5 V voltage monitoring, see Item 19

Block No. 18 relay O₂ (3/K3)

- Action in normal operation : Enabling of 24 V supply for O₂ valves
- Actuated by : Relay coil (3/K3) de-energized
O₂ sensor comparator = L,
see Item 30 and PIAL-PB6 = L
- Action in the event of malfunction : Switch-off of 24 V supply for O₂ valves
- Actuated by : Relay coil (3/K3) in O₂ sensor comparator = H or PIAL-PB6 = H
- Indication : Via potential of valve coils, see Items 3 and 4

Block No. 19 5 V voltage monitoring (1/N3)

Action in normal operation	:	Relay 24 V on (3/K4), see Item 17
Actuated by	:	PIA1-CA2 = L and PIA1-CB2 = H (1/V23, N3, V21)
Action in the event of malfunction	:	Relay 24 V off (3/K4), see Item 17
Actuated by	:	1) 5 V supply: Overvoltage \rightarrow -5.5 V $+ 0.2/- 0.1$ V 2) 5 V supply: Undervoltage < 4.4 V $+/- 0.2$ V 3) PIA1-CA2 = H 4) PIA1-CB2 = L 5) Reset
Indication	:	Via current monitoring heating air, see Item 8 Via indication, semiconductor relay 2, see Item 14
Test by	:	1) Test, comparator, undervoltage: PIA1-CB2 = L PIA1-CA2 = L 2) Test, comparator, overvoltage PIA1-CB2 = H PIA1-CA2 = H

Block No. 20 power failure monitoring +15, LED power failure (6/V26)

Action in normal operation	:	Power failure = Y (open collector)
Actuated by	:	+ 15 V OK (4/R11, V13, V25, V26)

- Action in the event of malfunction : Power failure = L,
LED (6/V26) lights up
- Actuated by : a) Failure +15 V supply
(4/rechargeable battery,
V25, V26)
- b) Mains switch "on" and
failure of +15 supply
(4/V5, R1)
- Test by : e.g. pulling mains connector

Block No. 21 watchdog 1 and 2 (1/D3A, B)

- Action in normal operation : Both mono-outputs = H
= Enabling of actuation for
- Semiconductor relay 1
(heating, air)
- semiconductor relay 2
(heating, H₂O)
- O₂ valve 1
- O₂ valve 2
(1/D9B, D9C, D7A-D)
- Actuated by : Punctual retriggering of mono-
flops by addressing
corresponding address:
Mono 1 = 8X3XH
Mono 2 = 8X4XH
- Action in the event of malfunction : One or both mono outputs = L
= blocking of actuation for:
- semiconductor relay 1
- semiconductor relay 2
- O₂ valve 1
- O₂ valve 2
- Actuated by : No punctual retriggering of
monos
- Indication : Current monitoring heating
air, see Item 8
Indication semiconductor
relay, see Item 14
Indication of potentials of
coils of O₂ valves, see Item
3, 4

Block No. 22 horn (4)

Action in normal operation : No action

Action in the event of malfunction : Sounds

Actuated by : 1) Power failure and mains switch "on" (4/storage battery, V5, V2), 15 V power supply = 0 V (S 1.3)

2) Elapse of watchdog 1 and/or 2 (1/D3A, D3B, D9B, R3), (4/V3, R2, V4, V21) see Item 17

3) PIA1-PB7 = H (1/D1, R2), (4/V3, R2, V4, V21)

Indication : Audible

Block No. 23 mains voltage monitoring, monoflop (4/D1)

Action in normal operation : $\overline{\text{NMI}} = \text{H}$

Actuated by : Retrigger pulses via (4/V15-V18, R12, V27)

Action in the event of malfunction : $\overline{\text{NMI}} = \text{L}$

Actuated byq : Power failure for more than one half-wave of the rectified voltage

Indication : $\overline{\text{NMI}} = \text{L}$
PIA2-PB6 = L with power failure and PIA2-PB7 = H (mains on), (5/S 1.3) closed (4/V8, V9, V22) see Item 11

Block No. 24 watchdog 3, INOP-LED (1/D6B)

- Action in normal operation : INOP-LED off
- Actuated by : Punctual retriggering of monoflop by addressing address 8X5XH
- Action in the event of malfunction : INOP-LED lights up
- Actuated by : No punctual retriggering of monoflop, (1/D6B), (6/R16, V17, V34, R14), +5 V supply
- Indication : Visual

Block 25 CPU Standard 2

Action in normal operation:

The p.c. board CPU Standard 2 is a microprocessor system featuring the microprocessor 6802, which contains the entire data and program storage area and various peripheral modules. Voltage monitoring and watchdog logic is also provided.

The address, data and control lines of the CPU are connected to the operating voltage by way of pull-up resistance networks, so as to obtain level matching with respect to the subsequent HCMOS modules. The address and control lines are buffered directly at the CPU with D2, D3 and D5, in order to enable the assembly to be operated in DMA mode. These driver modules are controlled by the signal BA.

The address decoder consists of the programmable module D6 and the subsequent decoder D7. Chip select signals for the memory modules D8 ... D12 and one select signal each for the internal I/O modules, the external I/O modules and the data bus driver are generated in D6. The programmability of D6 makes it possible - together with the jumper fields X2 ... X5, X8 - to have extremely flexible assembly with various memory modules (see 3.2.1 and 3.2.2). The signal BA at D7 is used to prevent an external CPU (DMA mode) having memory access to the I/O modules.

A serial interface can be realized with the ACIA D14. The interrupt output of D14 is connected via the diode V2 to the interrupt lines of the other peripheral modules. Moreover, the interrupt lines of D14 is routed via the inverter D28 to the printed circuit connector X1. Such special treatment of the ACIA interrupt makes it possible to have a hardware handshake with rapid data transfer via the serial interface.

The timer module D13 can be operated with an internal and an external block. One of the three clock pulses 100 Hz, 1 kHz, 10 kHz can be supplied (see Table 1) to each clock input C1, C2, C3 via the jumper field X6. The three gate inputs G1, G2, G3 and the timer outputs Q1, Q2, Q3 are routed to the printed circuit board connector X1.

The REAL TIME CLOCK D15 has its own clock generator and is battery-buffered. As this module has a bus interface for a multiplexed address/data bus, it must be addressed in two cycles in each case via the inputs AS and DS. The CS signal must therefore be permanently active and is only deactivated via the transistor V12 in the event of too low an operating voltage. The chip select signal CSRTC is contained in the circuitry of the inputs AS and DS.

The module D22 monitors the operating voltage and generates a RESET pulse for the CPU D1. When the operating voltage U is switched on, the RC-element R5C3 is initially still discharged. As long as U is less than 4.9 V, the output OUT2 of D22 remains on lo and discharges the RC-element R5C3 via V1. Once U is in excess of 4.9 V, OUT2 switches to hi and R5C3 is charged. The buffered RAM is now activated via R7, the operating voltage VCCBAT is switched to RAM via V6 and V10 and the REAL TIME CLOCK is activated via V12. As soon as the voltage at C3 increases to a level in excess of the internal reference voltage of D22 (approx. 1.3 V), the RESET pulse at OUT1 is terminated and the CPU can start.

If the operating voltage drops below 4.75 V, RAM and REAL TIME CLOCK are inhibited and a CPU-RESET is triggered. A RESET pulse can also be triggered during operation via the diodes V3 and V4 due to the discharge of R5C3.

The other timer D19 is reset in exactly the same manner as in mode A.

In the N mode the counter D20 can be set to a defined initial status by way of write access to a fixed address (CLRCNT = low). The sequence of addresses specified by the wiring of D20 and D21 is CSCNT + n, with n = 0, 4, 1, 5, 8, 12, 9, 13, 2, 6, 3, 7, 10, 14, 11, 15.

The free input U7 of D17 has no function and can be used for subsequent expansion.

If the processor is halted with the line HALT, the watchdog oscillator is also stopped.

Action in the event of malfunction: Reset

Action in the event of malfunction : Reset

Actuated by : No punctual retriggering of the two monoflops D18 and D19 with the addresses 9000H and 9001H

Block No. 25a 5 V voltage monitoring CPU

Action in the event of malfunction : Generates reset

Actuated by : Undervoltage of 5 V supply of < 4.75 V

Block No. 26a display interface, air (6/D1)

Action in normal operation : Read-in of keyboard, Actuation of 7-segment displays and LEDs: Control/operation (6/V18), Alarm (6/V19), Silence (6/V20), Fan failure (6/V21), Deviation > 1.5 °C (6/V22), Air heating (6/V23), > 39 °C (6/V24), > 37 °C (6/V25)

Actuated by : Data : 8X00H
Control : 8X01H

Indication : Data read-back potential from display RAM

Block No. 26b display interface, skin temperature (7/D1)

Action in normal operation : Actuation of 7-segment displays and LED: Sensor (7/V17) Control (7/V18) +/-0.5 °C (7/V21)

Actuated by : Data : 8X02H
Control : 8X03H

Indication : Data read-back potential from display RAM

Block No. 26c display interface O₂ (8/D1)

Action in normal operation : Actuation of 7-segment displays and LED:
Sensor (8/V17)
Control (8/V19)
> 40 vol. % (8/V19)
Cal. (8/V20)
+/-0.5 °C (8/V21)

Actuated by : Data : 8X04H
Control : 8X05H

Indication : Data read-back potential from display RAM

Block No. 26d display interface, humidity (9/D1)

Action in normal operation : Actuation of 7-segment displays and LED:
Sensor (9/V17)
Control (9/V18)
H₂O (9/V21)

Actuated by : Data : 8X06H
Control : 8X07H

Indication : Data read-back potential from display RAM

Block No. 27 A/D converter (2/D2)

Action in normal operation : 12 bit A/D conversion

Actuated by :- SOC: (CS1) 8X0CH.
B0-B7: (CS2) 8X0AH.
B8-B11, Busy, Data Valid:
(CS3) 8X08H.


Test : By way of +/-15 V voltage supply, see Item 35

Block No. 28 Multiplexer (2/D1)

Action in normal operation : Switching of voltage to be digitalized to A/D converter.


Actuated by : Address PIA1-PA: 8X1XH,
(1/D1), (2/V24 to V31)

Channel 0: Air temperature	- LXLL	LXXX
1: Excess temperature	- LXLL	HXXX
2: Skin temperature	- LXLH	LXXX
3: Skin temperature Sensor alarm	- LXLH	HXXX
4: O ₂ -Sensor A	- LXHL	LXXX
5: O ₂ -Sensor B	- LXHL	HXXX
6: Humidity	- LXHH	LXXX
7: A/D test	- LXHH	HXXX

Indication : Via PIA2-PA: 8X2XH, (1/D2)

Channel 0: Air temperature	- XXLL	LXXX
1: Excess temperature	- XXLL	HXXX
2: Skin temperature	- XXLH	LXXX
3: Skin temperature Sensor alarm	- XXLH	HXXX
4: O ₂ -Sensor A	- XXHL	LXXX
5: O ₂ -Sensor B	- XXHL	HXXX
6: Humidity	- XXHH	LXXX
7: A/D test	- XXHH	HXXX

Test : By way of +/-15 V voltage supply, see Item 35

Block No. 29 excess-temperature comparator (2/N5)

Action in normal operation	:	No action, output = L
Actuated by	:	Air temperature < 39.8 °C
Action in the event of malfunction	:	Comparator output = H (+15 V) switches off the excess-temperature relay (see Item 12) and thus the valves and heating via the relay UET (see Item 16)

Actuated by : Air temperature > 40.0 °C
Indication : Via relay excess-temperature,
see item 16
Test by : Test resistance (2/R107)

Block No. 30a/30b comparator O₂ (2/N10)

Action in normal operation : No action
Actuated by : O₂ sensor voltages > 9 mV
Action in the event of malfunction : Switches off valves via O₂
relay, see Item 18
Actuated by : O₂ sensor voltage of one
capsule < 9 mV (sensor
short-circuit or no O₂
capsule connected)
Indication : O₂ valve coil voltage
indication, see Item 3, 4
Test by : e.g. removing O₂ capsules

Block No. 31a Excess-temperature amplifier (2/N3, N4)

Action in normal operation : Bridge amplifier for excess
temperature sensor,
35.0 °C = $R_{\text{sensor}} = 3266\Omega$
= 625 mV at TP UT1 =
3.500 V at TP UT2

Block No. 31b air-temperature amplifier (2/N1, N2)

Action in normal operation : Bridge amplifier for air-
temperature sensor,
35.0 °C = $R_{\text{sensor}} = 3266\Omega$
625 mV at TP LT1 =
2.917 V at TP LT2

Block No. 32 skin-temperature amplifier (2/N6, N7, N8)

Action in normal operation : Bridge amplifier for skin-
temperature sensor,
36.0 °C = $R_{\text{sensor}} = 1412\Omega$
= 308.6 mV at TP HT1 =
2.857 V at TP HT2

Block No. 33a sensor amplifier A (2/N10)

Action in normal operation : Amplifies O₂ sensor voltage

Block No. 33b O₂ sensor amplifier B (2/N10)

Action in normal operation : Amplifies O₂ sensor voltage

Block No. 34 humidity amplifier

Action in normal operation : Supplies via (2/V16, R54) the +12 V supply voltage for the humidity module and steps down its output voltage (2/R52, R53)

Block No. 35 +/-15 V voltage monitoring (2/R55, R83)

Action in normal operation : No action

Action in the event of malfunction : Test A/D converter
Test +/-15 V

Actuated by : Voltage divider (2/R55/R83) between + and -15 V,
Check for 4.85 +/-0.75 V

Block No. 36 relay, test comparator, excess temperature (2/K2)

Action in normal operation : Switches a test resistor (2/R107) in accordance with > 40 °C to excess-temperature sensor input

Actuated by : PIA2-PB5 = H (1/D2)

Indication : Current monitoring heating air, see Item 8
Indication semiconductor relay 2, see Item 41
Indication by way of electr. potential at valve coils, see Item 3, 4

Block No. 41a/b O₂ sensor A/B

Action in normal operation : Output voltage proportional to oxygen content

Block No. 42 humidity sensor

Action in normal operation : Output voltage of humidity module proportional to relative humidity in limit range between 3.5 V = 30 % rel. humidity and 7.75 V = 90 % rel. humidity

Block No. 43 sensor switch

Action in normal operation : Indicates position of sensor head to PIA2-PBO and PIA2-PB4

Block No. 44 additional fan

Action in normal operation : Fan on: PIA2-PA7 = H
Fan off: PIA2-PA7 = L

Actuated by : Current transformer (3/RWP 209G)
Monoflop (1, N3),
Monoflop (10, N1)

Action in the event of malfunction : Fan on (except +15 V is missing)

2.4 Alarms

The incubator features, a stepped alarm behaviour.

Faults occurring are indicated as a function of their significance.

The alarms for set-point deviations, overtemperature and lack-of-water feature an intermittend alarm tone which can be deactivated for 10 minutes. A horn with a non-deactivatable continuous tone is assigned to all other faults. The corresponding alarm LED also flash in each case. Each individual alarm likewise leads to an indication on the central alarm LED:

- Three lines (air and skin temperature) or two lines (O₂, humidity) flash in the respective actual-value display.
- In the event of a fan fault, the air-temperature actual value flashes.
- The corresponding alarm LED flashes.
- Central alarm LED flashes.
- Continuous horn tone (cannot be deactivated).

As regards the sensor alarms for oxygen concentration and humidity, the alarm is reset following switch-off of the module concerned and further operation can be continued with the other functions. The same applies to the skin-temperature module, however, the defective sensor must first be removed in this case.

2.4.3 Module fault

Hardware faults affecting a single module result in the following alarm messages for the skin-temperature, oxygen-concentration and humidity modules:

- The letters "Err" flash up on the corresponding actual-value display if the module is switched on or appear in steady form if the module is switched off.
- INOP-LED lights up if module is switched on.
- Continuous horn tone (non-deactivatable) if module is switched on.

The functions of the other modules which are working properly are retained.

2.4.4 INOP

Serious hardware or software faults result in complete failure of the incubator. The following alarm messages are given:

- The INOP-LED lights up.

- Continuous horn tone (non-deactivatable).

The equipment is not serviceable.

3. Service-Software

3.1 Service-mode

Two different service modes are implemented in the Incubator. The first mode can only be reached by actuating simultaneously the buttons "Audial alarm off" and "-Air" and by switching on the device.

For further information see "List of errors".

While the device is operating the second mode can be reached as follows:

Actuate buttons "Audial alarm off" and "-Air" simultaneously for approx. 8 s permanently. ^{DFC ↓}

The device then starts service-mode 0:

°C Air/Air		
X	Y	/
/	/	0


Service-mode 0 indicates the version number of the software

Version number is X, Y
e.g. 0.4

The corresponding number of the service mode is flashing in the right digit of the "Air". Display for the desired value. By actuating button "+Air" the next mode can be activated. Mode 0 re-appears after the last service mode. Escape from service mode can be achieved by pressing a control-button.

-S- appears in the set-value-displays of the modules "Skin" "O₂", "Humidity" where applicable.

Mode 1:

Service mode 1
indicates the
operational hours dar.

Elapsed time meter:
ABCDE in hours

Setting elapsed time meter:

By pressing button "Audible alarm off" the corresponding variable digit is selected (flashing indicator).

Actuating button "> 37 °C", the value of the corresponding digit will increase by 1.

In case of a replacement of the complete electronics or the p.c. board "CPU" the last value of the elapsed time meter has to be transferred on the new unit by calling the set mode as described above. If this figure is not available for some reason, the respective value is to be taken from the last entry in the test certificate.

DrägerService
Medizintechnik
Medical Systems

Mode 2:

°C Air		
A	B	C
D	E	2

Service mode 2 indicates the number of the actuations of button "check".

ABCDE number of actuations

Again recommended to check once / Day.

The counter for actuations is variable as described in mode 1.

The value of the elapsed time meter must be entered in the test certificate.

*change by Horn cancell and > 3.9
↑
10/100, 1000's
when NO*

Mode 3:

°C Air		
A	B	C
D	E	3

Service mode 3 indicates the number of the actuations of button "check 36 °C".

ABCDE number of actuations

To check to see if this has been carried out regularly.

The counter for actuations is variable as described in mode 1.

The value of the elapsed time meter must be entered in the test certificates.

on skin Temp.

Mode 4:

°C Air		
/	X	/
/	/	4

Service mode 4 indicates the position of the sensor head.

- X: 1: Swivelled in
2: Intermediate position
3: Swivelled out
4: Sensor disconnected

Alarm sound in 3

DrägerService
Medizintechnik
Medical Systems

Mode 5:

°C Air		
X	Y	Z
/	/	5

Service mode 5 indicates set-value of overtemperature sensor

Overtemp.:
XY, Z in °C

30.4
-5

DrägerService
Medizintechnik
Medical Systems

Mode 6:

°C Air		
X	Y	Z
/	/	6

Service mode 6 Mode means the
air temperature sensor

Air temperature
XY, Z in °C

Mode 7:

°C Air		
X	Y	Z
/	/	7

Service mode 7 indicates the measured value of the O₂-sensor A (outer sensor)

O₂ A-Value:
XYZ measured value of the A/D converter in decimal numbers

Valid range for 21 Vol.% O₂:

306 to 812

$$= 9.5 \text{ mV}$$

Mode 8:

°C Air		
X	Y	Z
/	/	8

Service mode 8 indicates the measured value of the O₂-sensor A (inner sensor)

O₂ A-Value:
XYZ measured value of the A/D converter in decimal numbers

Valid range for 2l Vol.% O₂:

306 to 812

3.2 List of errors Incubator 8000

This list of errors contains all possible errors occurring to Incubator 8000, only in case of INOP error it is to be additionally checked, if the Reset-LED on the p.c. board CPU lights. The numbers of errors can be indicated via operating system of Incubator 8000 as follows:

1. Actuate simultaneously buttons "Audible alarm off" and "Decrement set value air".
2. Switch on Incubator.
3. Keep buttons pressed for approximately 5 sec. perpetually. Premature disengagement will lead to a system cold start.
4. The device responds with Character "F" or the "Air" display for the actual value.

By actuating button "Inc. air" the search routine continues for further possible errors and their indication.

If an indicated error is to be cleared, the button "Reset" is to be pressed followed by button "Inc. air" in addition. The respective error is then cleared, and the searching routine continues. If no further error occurs the "Air"-display for the desired value indicates "- - -".

5. The indication of errors is terminated by switching off the Incubator.

The number of the error is identical with the physical address, e.g.: F1 = \$ 000,1, F10 = \$ 000A.

Explanations:

Column "May occur during operation":
Indicated errors are concerning internal and external hardware errors. e.g. if an O₂ sensor is worn out, error "F 59" is indicated.

"Self test": All errors being recognized during self test (5 sec. after power up and then restoring every 10 minutes).

"Permanent feedback": Errors being discovered in various routines while programme is running.

Column "Result":

Indicates the effects of the error. In the event of INOP-Errors the Incubator cuts off the main line. Errors concerning the optional modules (skin temp., humidity and O₂) can be disabled on the respective unit itself which responds with the display code "ERR". The name of the variable characters written in brackets indicating the programme part where the error (=FF) is stored.

DrägerService
Medizintechnik
Medical Systems

Number of error	Set in routine	May occur during operation	Reason	Result
1	AD-Manage		A/D converter has not finished conversion after 33 m sec. (permanent feedback)	Inop
2	AD-Manage		A/D converter is active before being started. (permanent feedback)	Inop
3	Test 10 min	x	Water heating unable to be activated. If humidity option does not exist this error code may appear. (self test)	Error Humidity-module
4	Test 10 min	x	Measured value of skin temperature test resistor exceeding specified range from 35.8 °C to 36.2 °C. If skin temperature option does not exist, error code may appear! (self test)	Error Skin-module
5	Test 10 min	x	Valve 1 unable to be activated in spite of existing O ₂ -sensor. If O ₂ option error code will appear. Error message 5 is not applicable as of EPROM 8290371 (software 0.4) (permanent feedback)	Error O ₂ -module
6	Skin sens	x	Skin sensor short-circuit or measuring amplifier faulty. If skin temperature does not exist, error code may appear. (permanent feedback)	Error Sensor Skin-module

Number of error	Set in routine	May occur during operation	Reason	Result
7	Dis-unit		RAM-Reading of display controller(s) malfunctioning. (permanent feedback)	Inop
8	Test 10 min		Heating of water nebulizer unable to be switched off via PIA.	Inop
9	Test 10 min		Fan failure test unable to be activated. Pulse signals of fanwheel cannot be disabled or Mono faulty. (self test)	Inop
10	Test 10 min		Air heater cannot be activated or feedback-line malfunctioning. (self test)	Inop
11	Test 10 min		Air heater cannot be switched off via PIA. (self test)	Inop
12	Test 10 min		Over-temperature signal via PIA does not cut air heater. (self test)	Inop
13	Test 10 min		Watchdog cycle (Mono 1 or 2) too fast. (30 - 60 ms) (self test)	Inop
14	Test 10 min		Watchdog Background (Mono 1) does not run or runs too slow or does not cut off slow air-heater. (self test)	Inop

Number of error	Set in routine	May occur during operation	Reason	Result
15	Test 10 min		Watchdog IRQ (Mono 2) does not run or runs too slow or does not cut off slow air-heater. (self test)	Inop
16	Test 10 min		Watchdog IRQ (Mono 2) does not run or runs too slow or does not cut off heater for water nebulizer (self test)	Inop
17	Test 10 min		Watchdog IRQ (Mono 2) does not run or runs too slow or does not cut off valve 1. (self test)	Inop
18	Test 10 min		Watchdog IRQ (Mono 2) does not run or runs too slow or does not cut off valve 2. (self test)	Inop
19	Test 10 min		A/D reference voltage >4.79 V (self test)	Inop
20	Test 10 min		Test voltage check upper limit: Air-heater does not cut off. (self test)	Inop
21	Test 10 min		Test voltage check lower limit: Water heating does not cut off. (self test)	Inop
22	Test 10 min		Test voltage check lower limit does not cut off valve 1. (self test)	Inop

DrägerService
Medizintechnik
Medical Systems

Number of error	Set in routine	May occur during operation	Reason	Result
23	Test 10 min		Test voltage check lower limit does not cut off valve 2. (self test)	Inop
24	Test 10 min		Test voltage check lower limit does not cut off air-heater. (self test)	Inop
25	Test 10 min		Inspite of non-resettable overtemperature signal from HW-Comparator (>40 °C) air the heater does not cut off. (self test)	Inop
26	Test 10 min		When using test resistor (40 °) HW-Comparator does not switch off bistable relay. (self test)	Inop
27	Test 10 min		Reset of overtemperature (releasing the actuators) does not work. (self test)	Inop
28	Test 10 min		O ₂ -valve 1 does not switch off via PIA or feedback wrong (self test)	Inop
29	Test 10 min		Over-temperature signal via PIA does not cut valve 1. (self test)	Inop

Number of error	Set in routine	May occur during operation	Reason	Result
30	Test 10 min		O ₂ -valve 2 does not switch off via PIA or feedback wrong. (self test)	Inop
31	Test 10 min		Over-temperature signal via PIA does not switch off valve 2. (self test)	Inop
32	IRQ Timer		Air heater does not switch on. (permanent feedback)	Inop
33	IRQ Timer		Air heater is switched on permanently. (permanent feedback)	Inop
34	IRQ Timer		Valve 1 or valve 2 are permanently switched on, although not desired. (permanent feedback)	Inop
35	IRQ Timer		The humidity boiler is switched on permanently. (permanent feedback)	Inop
36	Fan	x	The fan does not turn due to faulty motor or fuse; fan wheel may not be mounted at all.	Fan Failure
37	AD-Manage		The feedback lines of the multiplexer addressing are not in accordance with the commands. (self test)	Inop

Number of error	Set in routine	May occur during operation	Reason	Result
38	IRQ Timer		The humidity boiler does not switch on. (permanent feedback)	Error Humidity-module
39	IRQ Timer		Valve 1 or Valve 2 does not switch on. (permanent feedback)	Error O ₂ -module
40	O ₂	x	Voltage for sensor 1 too high or corresponding amplifier faulty. (permanent feedback)	Malffunction of sensor O ₂
41	O ₂	x	Sensor differential voltage O ₂ exceeding permitted range; possibly sensor worn out or measuring amplifier faulty. (permanent feedback)	Malffunction of sensor O ₂
42	O ₂ Cal	x	No valid calibration value was found. Sensor(s) may be worn out or calibration trial when O ₂ content was varying. (permanent feedback)	Malffunction of sensor O ₂
43	O ₂	x	Average value of the two sensors is not in permitted range; possibly sensor(s) worn out or measuring amplifier faulty. (permanent feedback)	Malffunction of sensor O ₂
44	Air sensor	x	Air temperature sensor high-resistive. Error code is indicated when sensor is not connected. (permanent feedback)	Malffunction of sensor O ₂

Number of error	Set in routine	May occur during operation	Reason	Result
45	Humidity	x	Humidity sensor exceeding the permitted voltage range: Error code is indicated when sensor is not connected. (permanent feedback)	Malfunction of sensor Humidity
46	Air sensor		Difference between over-temperature settings exceed the permitted value. (permanent feedback)	Malfunction of sensor Air
47	Test 10 min	x	The connected skin temperature sensor is changing its values by less than 0.007 °C within 20 min. Possible cause: Contact of reed relay on p.c. board analog is sticking. If skin option does not exist error code may appear. (self test)	Option Error Skin
48	H ₂ O	x	Humidity boiler indicates lack of water. (permanent feedback)	Lack of water
49	Boot		Logic circuit "boot programme" is faulty.	Dead-end loop back-ground
50	Back-		Logic circuit "background programme controlling" is faulty.	Dead-end loop back-ground
51	CRC 16 Bit		ROM contents faulty.	Dead-end loop ROM-CRC

Number of error	Set in routine	May occur during operation	Reason	Result
52	Walk-Pat		Result of RAM test negative	Dead-end loop RAM-Walkpat
53	PIA Test		PIA's cannot be initialized.	Dead-end loop PIA-Unit
54	Timer Test		Timer cycles are incorrect.	Dead-end loop Timer
55	IRQ Timer		Logic circuit "IRQ controlling" is faulty.	Dead-end loop IRQ
56	RAM Test		Result of RAM test negative	Dead-end loop RAM-CRC
57	Air sensor		Air temperature sensor short circuit. (permanent feedback)	Malfunction of sensor Air
58	Skin sensor	x	Skin temperature sensor high resistive or measuring amplifier faulty. If skin option does not exist, error code may appear. (permanent feedback)	Malfunction of sensor Skin
59	O ₂	x	Cell voltage for Sensor 1 not high enough or corresponding amplifier faulty. (permanent feedback)	Malfunction of sensor O ₂

Number of error	Set in routine	May occur during operation	Reason	Result
60	O ₂	x	Voltage for Sensor 2 too high or corresponding amplifier faulty. (permanent feedback)	Malfunction of sensor O ₂
61	O ₂	x	Voltage for Sensor 2 not high enough or corresponding amplifier faulty. (permanent feedback)	Malfunction of sensor O ₂
62	Test 10 min	x	Valve 2 does not switch on inspite of existing O ₂ -sensors. If O ₂ -option does not exist, this error code appears. Error message 5 is not applicable as of EPROM 8290371 (software 0.4) (self test)	Error Option O ₂
63	ROM-Exor		The ROM is faulty. (permanent feedback)	Dead-end loop ROM
64	Reference		Data in ROM is faulty. (self test)	Constants
65	Test 10 min		A/D-reference voltage < 3.99 Volt. (self test)	Inop
66	INOP		Contact of the power switch is open or PIA D2. Port B7 is faulty. (permanent feedback)	Inop
67			After switching off the O ₂ module the O ₂ concentration increased within 3 minutes by more than 3 vol. %. Possible malfunction causes: O ₂ valve leaky or external O ₂ supply	Inop and additional O ₂ display

3.3 Description of 10-minute test

10-minute test sequence

A software-controlled check of the hardware of relevance to safety is carried out following switch-on of the incubator, following a mains failure and every 10 minutes.

Sequence in program

- Salvaging of input statuses
- Salvaging of setting of relay, overtemperature, bistable
- Is there an overtemperature situation currently present ($> 40\text{ °C}$)?

Yes: In the event of overtemperature ($> 40\text{ °C}$) the overtemperature comparator switches off the overtemperature relay by way of the relay, overtemperature, bistable. In this case, if an overtemperature reset is not possible on account of the fact that an overtemperature situation is still present ($> 40\text{ °C}$), a 10-minute test cannot be performed and a test is merely carried out to check whether the air heating is switched off.

Check-back signal: current through air heater?

Yes: INOP + fault 25

No: No fault + end of 10-minute test

No: Start of 10-minute test

"Switch-on of air heating" (fault 10)

- Switch on air heating only via PIA
- Check-back signal: current through air heater?
Yes: No fault + continue
No: INOP + fault 10 + continue

"H₂O heating on and O₂ valves off" (faults 3, 28, 30)

- Switch on H₂O heating only via PIA
- Check-back signal: current through H₂O semiconductor relay?
Yes: No fault + continue
No: Humidity module fault + fault 3 + continue

- Check-back signal: current through valve 1?
Yes: INOP + fault 28 + continue
No: No fault + continue
- Check-back signal: current + through valve 2?
Yes: INOP + fault 30 + continue
No: No fault + continue

"O₂ valves on and air/H₂O heating off"
(faults 5, 8, 62)

- Switch on O₂ valves 1 and 2 only via PIA (does not apply of EPROM 8290371)
- Check-back signal: current through valve 1? (does not apply as of EPROM 8290371)
Yes: No fault + continue
No: O₂ module fault + fault 5 (only if sufficiently large O₂ cell voltages are measured) + continue
- Check-back signal: current through valve 2? (does not apply as of EPROM 8290371)
Yes: No fault + continue
No: O₂ module fault + fault 62 (only if sufficiently large O₂-cell voltages are measured) + continue
- Check-back signal: current through H₂O semiconductor relay?
Yes: INOP + fault 8 + continue
No: No fault + continue

"Blocking fan pulses (fault 9)

(The fault message "fan not running" is not generated here, but rather during continuous monitoring of the fan in normal operation. This test only checks the monitoring circuit of the fan failure alarm.)

- Block fan pulses via PIA (mono is not retriggered)
- Check-back signal: fan pulses blocked (fan failure)?
Yes: No fault + continue
No: INOP + fault 9
- Enable fan pulses

"Air heating off" (fault 11)

- Air heating is still switched off via PIA
- Check-back signal: current through air heater?
Yes: INOP + fault 11 + continue
No: No fault + continue

"Testing monoflops" (faults 13, 14)

- Retriggering of monos 1, 2 and 3
- Switch on air heating via PIA
- Waiting loops, monos 1 and 2 must not elapse and block air heating
- Check-back signal: current through air heater?
Yes: No fault + continue
No: INOP + fault 13 + continue
- Retriggering of monos 2 and 3
Mono 1 must elapse 140 ms after retriggering and block air heating.
- Check-back signal: current through air heater?
Yes: INOP + fault 14 + continue
No: No fault + continue

"Comparator 40 °C test (overtemperature)" (fault 26)

- Switch off heaters and valves via PIA
- Switch test resistor 40 °C to overtemperature amplifier
- Check-back signal: has relay, overtemperature, bistable switched to overtemperature?
Yes: No fault + continue
No: INOP + fault 26 + continue

"Overtemperature by way of PIA" (faults 12, 27, 29, 31)

- Trigger overtemperature and O₂ relay by way of PIA
- Switch off test resistor 40 °C (normal setting), overtemperature relay does however remain switched off.
- Trigger overtemperature and O₂ relay by way of PIA and simultaneously attempt to trigger an overtemperature reset (normally not possible, overtemperature must remain).
- Switch on air heating and valves via PIA (must not switch on)
- Check-back signal: current through valve 1?
Yes: INOP + fault 29 + continue
No: No fault + continue
- Check-back signal: current through valve 2?
Yes: INOP + fault 31 + continue
No: No fault + continue
- Check-back signal: current through air heater?
Yes: INOP + fault 12 + continue
No: No fault + continue
- Switch off air heating and valves via PIA
- Overtemperature reset (enabling of relay, overtemperature, bistable)

- Switch on air heating only via PIA
- Check-back signal: current through air heater?
Yes: No fault + continue
No: INOP + fault 27 + continue

"Testing mono 2" (faults 15, 16, 17, 18)

- Retriggering of monos 1, 2 and 3
- Switch on heaters and valves via PIA
- Waiting loop for testing mono 2
- Retriggering of monos 1 and 3
- Waiting loop; mono 2 must elapse 140 ms following last retriggering and block both heaters and valves
- Check-back signal: current through air heater?
Yes: INOP + fault 15 + continue
No: No fault + continue
- Check-back signal: current through H₂O semiconductor relay?
Yes: INOP + fault 16 + continue
No: No fault + continue
- Check-back signal: current through valve 1?
Yes: INOP + fault 17 + continue
No: No fault + continue
- Check-back signal: current through valve 2?
Yes: INOP + fault 18 + continue
No: No fault + continue

"Testing A/D converter reference" (faults 19, 65)

- Read-out of reference-channel voltage
- Reference-channel voltage less than 3.99 V?
Yes: INOP + fault 65 + continue
No: No fault + continue
- Reference-channel voltage greater than 4.79 V?
Yes: INOP + fault 19 + continue
No: No fault + continue

"Testing 5 V voltage monitoring" (faults 20, 21, 22, 23, 24)

- Switch off heaters and valves via PIA
- Test upper limit of voltage monitoring via PIA (24 V relay must switch off heaters and valves)
- Switch on air heater via PIA
- Check-back signal: current through air heater?
Yes: INOP + fault 20 + continue
No: No fault + continue
- Switch off air heater via PIA
- Cancellation of testing of upper limit of voltage monitoring

- Test lower limit of voltage monitoring via PIA (24 V relay must switch off heaters and valves)
- Switch on heaters and valves via PIA
- Check-back signal: current through H₂O semiconductor relay?
 - Yes: INOP + fault 21 + continue
 - No: No fault + continue
- Check-back signal: current through valve 1?
 - Yes: INOP + fault 22 + continue
 - No: No fault + continue
- Check-back signal: current through valve 2?
 - Yes: INOP + fault 23 + continue
 - No: No fault + continue
- Check-back signal: current through air heater?
 - Yes: INOP + fault 24 + continue
 - No: No fault + continue
- Cancellation of testing of lower limit of voltage monitoring (normal operation)

"Testing 36.0 °C skin-temperature simulation and evaluation of relay contact sticking" (faults 4, 47)

- Switch to 36.0 °C skin-temperature equivalent resistance
- Read-out of skin-temperature-amplifier voltage
- Temperature less than 35.8 °C or greater than 36.2 °C?
 - Yes: Skin module fault + fault 4 + continue
 - No: No fault + continue
- Has bit 2 switched in the last 20 minutes in normal operation (= change in temperature)?
 - Yes: No fault + continue
 - No: Skin module fault + fault 47 (only if skin temperature sensor connected and measured value between 33 and 38 °C) + continue

"End of 10-minute test"

- Restore input statuses
- Restore position of relay, overtemperature, bistable
- In the event of INOP fault situation: INOP-LED lights up, permanent horn tone, heaters and valves switched off.

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

PIA/port	Bit; < Input, > Output; CS: check-back signal	" 0 "	" 1 "
PIA D1 Port A (8010H)	0: < CS, relay, OT, bistable 1: < Skin-temperature identifier 2: < Identifier 02 3: > Multiplexer address A0 4: > Multiplexer address A1 5: > Multiplexer address A2 6: < Fan malfunction 7: > Multiplexer enable E-transverse	Overtemperature Skin module present 02 module present Running- enabled	o.k. Version without skin module without 02 module Not running/test disabled
PIA D1 Port B (8012H)	0: > Air heater 1: > H2O heater (evaporator) 2: > O2 valve 1 3: > O2 valve 2 4: > Relay, OT bistable - RESET - 5: > Relay, test 36° 6: > Relay, OT bistable - SET - 7: > Horn	"Off" "Off" "aus" "aus" 0 => 1: Air and H2O heater enable Normal operation Normal operation No horn tone	"On" "On" "On" "On" Check 36° Overtemperatur Horn sounds
PIA D2 Port A (8020H)	0: < Humidity identifier 1: < Thermoswitch, H2O deficiency 2: > DW-ext. alarm 3: < CS Multiplexer A0 4: < RM Multiplexer A1 5: < RM Multiplexer A2 6: < CS, relay 36° 7: < CS, air heater	Humidity module present o.k. o.k. Check "Off"	Version without humidity module Lack of water ALARM Normal operation "Off"
PIA D2 Port B (8022H)	0: < Sensor calibration switch 1: < CS, H2O heater 2: < CS, O2 valve 1 3: < CS, O2 valve 2 4: < Sensor operating switch 5: > Test overtemperature, 40° 6: < NMI-transverse 7: < Mains	"On" "On" "On" Normal operation NMI Mains "off"	"Off" "Off" "Off" Ot test No NMI Mains "on"
PIA-CONTR		Data word	Data word
PIA D1 CA2 (8011H)	> 5 V - Voltage monitoring	Normal operation: 34H Output: 0	Upper limit test: 3CH Output: 1
PIA D1 CB2 (8013H)	> 5 V - Voltage monitoring	Normal operation: 3CH Output: 1	Lower limit test: 34H Output: 0
PIA D2 CB2 (8023H)	< > Fan-pulse suppression	Normal operation: 04H Input (int. pull-up)	Test: 3CH Output: 1

5. Calibration specifications

Contains all calibration specifications for the basic components described under 2.2.

5.1 Unit

5.1.1 Thermostat, lack-of-water alarm

The thermostat of the lack-of-water alarm is to be set to 120 °C in line with the scale. This thermostat must switch off the water heating prior to the safety thermostat (135 °C) via the processor system and it should not switch in normal operation with water to the maximum specified set value. In future, this adjustable thermostat will be replaced by a permanently set thermostat (design similar to that of safety thermostat). Installation potential is provided with all units.

5.1.2 Fan motor


5.1.2.1 Clearance of the motor shaft

Permissible clearance of the motor shaft = 0.1 mm

5.1.2.2 Fan motor removal and installation

- Remove unit.
- Unscrew coil of fan-failure alarm from aluminium block of motor mount.
- Detach cable connection between motor and p.c. board "unit".
- Unscrew retaining bracket of motor.
- Remove motor after loosening the 4 retaining nuts.
- Fit new motor springs to set screws and screw in retaining nuts approx. 3.5 mm measured from upper edge of threaded rod to nuts.

- Find adjustment: The shaft of the motor must protrude 43.0 ± 0.5 mm over the upper edge of the aluminium block and be at right angles to it. The motor must hang freely in the springs. The nuts are then to be secured with Locite 221.

- Installation is to be performed in reverse order.

When connecting the earthing studs, care is to be taken to ensure that these do not protrude over the edge of the spring-mounted motor retaining plate. The countersunk screws of the motor retaining bracket are to be secured with locking compound.

- Following installation, a check is to be made as to whether the impeller catches on the through and as to whether it moves freely in the spring mount.

5.2 p.c. board "CPU Standard 2"

5.2.1 Component configuration


D 9	- Not configured
D 10	- Not configured
D 11	- Not configured
D 12	- EPROM
D 13	- Timer 6840
D 14	- Not configured
D 15	- Not configured
D 6	- Programmable address decoder 8304920
D 8	- RAM (2K x 8) up to software 0.3 or (8K x 8) as of software 0.4

See 5.2.2 for layout plan

5.2.2 Assignment of jumper fields

Fitting of jumper fields X2 ... X8 with shorting links.
The pin numbering on the p.c. board is counterclockwise.

X 2	- Not fitted
X 3	- Not fitted
X 4	- Not fitted
X 5	- 1-2, 3-4 (32 K x 8) EPROM
X 6	- 1-18, 4-15, 8-11 (ext. clock for timer)
X 7	- Not fitted
X 8	- 3-4 (2 K x 8 RAM) up to software 0.3 or 1-2 (8 K x 8 RAM) as of software 0.4

5.2.3 Adjustment of voltage monitoring

If the operating voltage drops below 4.75 V, the RAM is blocked and a CPU reset is triggered. The voltage monitoring circuit can be adjusted with potentiometer R1 on the p.c. board "CPU".

- Using potentiometer R19 on the p.c. board "power pack" (see 5.3) adjust supply voltage, measured between TP1 (+) and TP2 (GND) on p.c. board "CPU", from approx. 5 V to $4.75 \text{ V} \pm 0.02 \text{ V}$.

Caution: Never set to in excess of 5.5 V!

- 2 initial statuses are possible for adjustment of the voltage monitoring on the p.c. board "CPU":
 - a) If the LED V11 lights up, turn potentiometer R1 in a clockwise direction until the LED V11 goes out or glows weakly. Then turn in opposite direction until LED V11 just lights up again.
 - b) If the LED V11 does not light up or glows weakly, turn potentiometer R1 in an anti-clockwise direction until the LED V11 just lights up.
- Then slowly increase operating voltage with potentiometer R 19 (p.c. board "power pack") (caution: never in excess of 5.5 V!) until the LED V11 goes out or glows weakly.

Test value: $U = 4.89 \text{ V} \pm 0.03 \text{ V}$

- Set operating voltage with potentiometer R 19 (p.c. board "power pack") to $5.00 \text{ V} \pm 0.05 \text{ V}$.
- Secure both potentiometers with locking compound.

5.3 p.c. board "power pack"

The 5 V supply voltage, measured between TP1 (+) and TP2 (GND) on the p.c. board "CPU", is to be set with the potentiometer R19 on the p.c. board "power pack" to $U = 5.00 \text{ V} \pm 0.05 \text{ V}$. The potentiometer is then to be secured with locking compound.

Refer also to Item 5.5.3 "Checking of operating voltages".

5.4 p.c. board "motherboard"

No adjustment necessary.

5.5 p.c. board "analog"

5.5.1 Measuring instruments and auxiliaries required


- Multimeter
- 2 resistance decades (skin-temperature adjustment only)
- Incubator 8000 sensor simulator
Order No. 7901240
- Skin-temperature sensor simulator
Order No. 7901236
- Incubator 8000 measurement adapter (not absolutely necessary)
Order No. 7901473

5.5.2 Important Information

All measuring points of importance for adjustment are combined in a 26-pole plug connector.

X4/1	TP O ₂ A	X4/14	TP HT 1B
2	TP Humidity	15	TP HT Analog
3		16	TP D/A Ref.
4	TP LT 2	17	TP HT 1
5	TP O ₂ B	18	TP +15 V
6	TP OT 1	19	TP GNDA
7	TP OT	20	TP +5 V Ref.
8	TP OT 2	21	TP GNDA
9	TP LT Analog	22	TP -15 V
10	TP A/D-Converter	23	TP GNDD
11	TP HT 2	24	TP -5V
12	TP O ₂ -Alarm	25	TP +5 V
13	TP O ₂ -Analog	26	TP GNDD

DrägerService
Medizintechnik
Medical Systems

5.5.3 Checking of operating voltages

Switch on unit.

The supply voltages are checked using a digital voltmeter at the plug connector and the "analog" p.c. board.


TP 25 -> TP 26 (GNDD): $\pm 5.0 \text{ V} \pm 0.05 \text{ V}$

TP 18 -> TP 26 (GNDD): $+ 15.0 \text{ V} \pm 0.6 \text{ V}$

TP 22 -> TP 26 (GNDD): $- 15.0 \text{ V} \pm 0.6 \text{ V}$

TP 20 -> TP 21 (GNDA): $+ 5.0 \text{ V} \pm 0.4 \text{ V}$

TP 24 -> TP 21 (GNDA): $- 5.0 \text{ V} \pm 0.4 \text{ V}$

5.5.4 Adjustment of air-temperature amplifier

- Connect Incubator 8000 sensor simulator (7901240) to X2.
- Switch on incubator.
- Wait approx. 35 s for self-test
- Set sensor simulator to 35.0 °C (R = 3266 Ohm).
- Connect digital voltmeter to X4/3 (+) and X4/21 (GNDA).
- Use potentiometer R7 to set voltage to 625 mV \pm 0.5 mV.
- Connect digital voltmeter to X4/4 (+) and X4/21 (GNDA).
- Use potentiometer R4 to set voltage to 2.917 V \pm 5 mV.
- Display on air actual-value display => 35.0 °C \pm 0.1 °C.
- Set sensor simulator to 25.0 °C (R = 5000 Ohm).
- Check voltage at X4/4
U_{des.} = 2.083 V \pm 0.05 V.
- Display on air actual-value display => 25.0 °C \pm 0.2 °C.

TEMP. SENSOR Alarm

- Set sensor simulator to 45.0 °C (R = 2184 Ohm).
- Check voltage at X4/4
U_{des.} = 3,750 V ± 0.05 V.
- Display on air actual-value display => 45.0 °C
± 0.2 °C.

5.5.5 Adjusting overtemperature amplifier

- Connect Incubator 8000 sensor simulator (7901240) to X2.
- Switch on incubator.
- Wait approx. 35 s for self-test.
- Call up service mode 5 = measured value, overtemperature amplifier (see Item 3.1, service mode).
- Set sensor simulator to 35.0 °C (R = 3266 Ohm).
- Connect digital voltmeter to X4/6 (+) and X4/21 (GNDA).
- Use potentiometer R 15 to set voltage to 625 mV ± 0.5 mV.
- Connect digital voltmeter to X4/8 (+) and X4/21 (GNDA).
- Use potentiometer R12 to set voltage to 3.500 V ± 5 mV.
- Display on air actual-value display => 35.0 °C ± 0.1 °C.
- Set sensor simulator to 25.0 °C (R = 5000 Ohm).
- Check voltage at X4/8
U_{des.} = 2.500 V ± 0.05 V.
- Display on air actual-value display => 25.0 °C ± 0.2 °C.
- Set sensor simulator to 45.0 °C (R = 2184 Ohm).

*Audible Alarm when A
No Led on over Temp*

*only alarms
because
in seller
of Temp set
NOT over
40°*

- Check voltage at X4/8
 $U_{des.} = 4.500 \text{ V} \pm 0.05 \text{ V}.$
- Display on air actual-value display $\Rightarrow 45.0 \text{ }^\circ\text{C} \pm 0.2 \text{ }^\circ\text{C}.$

Note: If the measured values between the air-temperature amplifier and the overtemperature amplifier differ by more than $0.5 \text{ }^\circ\text{C}$, a sensor alarm is given.

The two measured values may deviate by a maximum of $0.3 \text{ }^\circ\text{C}$ from one another in the range of adjustment of the sensor simulator and in normal operation.

Et OK. on simulator but more on normal operation the sensors may be faulty

5.5.6 Adjustment of overtemperature comparator

- Connect Incubator 8000 sensor simulator (7901240) to X2.
- Switch on incubator. *TAKE out*
- Wait approx. 35 s for self-test. *of Service*
- Set sensor simulator to $40.0 \text{ }^\circ\text{C}$ ($R = 2663 \text{ Ohm}$). *MODE*
- Connect digital voltmeter to X4/7 (+) and X4/21 (GNDA).
- Set potentiometer R19 such that voltage switches from negative to positive. *A.C.W. $\rightarrow +$*
- Set sensor simulator to $39.5 \text{ }^\circ\text{C}$ ($39.7 \text{ }^\circ\text{C}$)
 $U \leq -12 \text{ V}$ *more -ve -14V*
- Set sensor simulator to $40.5 \text{ }^\circ\text{C}$ ($40.2 \text{ }^\circ\text{C}$)
 $U \geq +12 \text{ V}$ *+14V*

5.5.7 Adjustment of skin-temperature amplifier

- Connect skin-temperature sensor simulator (7901236) to climate sensor or to Incubator 8000 sensor simulator. *Plug onto ada simulator*
- Set skin-temperature sensor simulator to $36.0 \text{ }^\circ\text{C}$ ($R = 1412 \text{ Ohm}$).
- Connect digital voltmeter to X4/17 (+) and X4/21 (GNDA).

- normally around 54K
50-60K Range*
- Connect a resistance decade to the soldering terminals of R24 and set such that the voltage at the digital voltmeter is $308.6 \text{ mV} \pm 0.3 \text{ mV}$.
 - Measure out the set resistance value to the nearest 10 Ohm and solder in the resistance value determined (if applicable, series or parallel connection of 2 resistors).
 - Check adjustment.
 - Connect digital voltmeter to X4/11 (+) and X4/21 (GNDA).
 - Connect a resistance decade to the soldering terminals of R31 and set such that the voltage at the digital voltmeter is $2.857 \text{ V} \pm 10 \text{ mV}$.
 - Measure out the set resistance value to the nearest 10 Ohm and solder in the resistance value determined (if applicable, series or parallel connection of 2 resistors).
 - Check adjustment.
 - Display on ~~air~~ *skin Temp.* actual-value display $\Rightarrow 36.0 \text{ }^\circ\text{C} \pm 0.1 \text{ }^\circ\text{C}$.
 - Press button "check 36.0 °C"
 \Rightarrow Display on skin-temperature actual value display
 $\Rightarrow 36.0 \text{ }^\circ\text{C} \pm 0.1 \text{ }^\circ\text{C}$.
skin Temp. wont change because skin connects to Rummy
 - Set skin-temperature sensor simulator to $33.2 \text{ }^\circ\text{C}$
 \Rightarrow Display on skin-temperature actual value display
 $\Rightarrow 33.2 \text{ }^\circ\text{C} \pm 0.2 \text{ }^\circ\text{C}$.
 - Set skin-temperature sensor simulator to $37.8 \text{ }^\circ\text{C}$
 \Rightarrow Display on skin-temperature actual value display
 $\Rightarrow 37.8 \text{ }^\circ\text{C} \pm 0.2 \text{ }^\circ\text{C}$.

5.5.8 Checking of O₂ input amplifiers and O₂connectors

- Connect Incubator 8000 sensor simulator (7901240) to X2.
- Set sensor simulator to $35 \text{ }^\circ\text{C}$.

- Switch on incubator.
- Wait approx. 35 s for self-test
- Deactivate any overtemperature alarm by pressing the "reset!" button.
- Set voltage for the O₂ amplifiers at the Incubator 8000 sensor simulator to 50 mV with the two O₂ potentiometers. *Remove leads and set vol least on Test Box.*
- Connect digital voltmeter to X4/1 (+) and X4/21 (GNDA) and check output voltage of O₂ amplifier
A: $U = 1.97 \text{ V} \pm 0.05 \text{ V}$.
- Connect digital voltmeter to X4/5 (+) and X/21 (GNDA) and check output voltage of O₂ amplifier
B: $U = 0.72 \text{ V} \pm 0.02 \text{ V}$.
- Set voltage for the O₂ amplifiers at the Incubator 8000 sensor simulator to 9.5 mV with the two O₂ potentiometers.
- Connect digital voltmeter to X4/12 (+) and X4/21 (GNDA): $U = < 0.1 \text{ V}$.
- Slowly reduce voltage for O₂ amplifier A DVM reading: $U = > 11 \text{ V}$.
- Set voltage for O₂ amplifier A to 9.5 mV again
DVM reading: $U = < 0.1 \text{ V}$. *100 mV* *1st Reading*
- Slowly reduce voltage for O₂ amplifier B DVM reading: $U = > 11 \text{ V}$. *12V* *2-8V*
- Set voltage for O₂ amplifier B to 9.5 mV again
DVM reading: $U = < 0.1 \text{ V}$. *100 mV* *Reset over*

5.5.9 Adjustment of air-temperature analog output (not used at present)

- Set sensor simulator to 25.0 °C.
- Connect digital voltmeter to X4/9 (+) and X4/21 (GNDA).
- Use potentiometer R67 to set voltage to $2.5 \pm 10 \text{ mV}$.

*Not
Registered*

- Not
used*
- Set sensor simulator to 40.0 °C.
 - Use potentiometer R69 to set voltage to 10 V ± 10 mV.
 - Set sensor simulator to 35.0 °C
DVM reading: 7.5 V ± 10 mV.

5.5.10 Adjustment of skin-temperature analog output (not used at present)

- not
used*
- Additionally connect skin-temperature sensor simulator and set to 37.8 °C.
 - Connect digital voltmeter to X4/15 (+) and X4/21 (GNDA).
 - Connect a resistance decade to the soldering terminals of R60 and set such that the voltage on the DVM is 9.60 V ± 10 mV.
 - Measure out set resistance value to the nearest 10 Ohm and solder in the resistance value determined.
 - Check adjustment.
 - Set skin-temperature sensor simulator to 33.2 °C
DVM reading: 0.40 V ± 10 mV.
 - Set skin-temperature sensor simulator to 36.0 °C
DVM reading: 6.0 V ± 10 mV.

5.5.11 Adjustment of reference voltage of D/A converter (not used at present)

- not
used*
- Connect digital voltmeter to X4/16 (+) and X4/21 (GNDA).
 - Use potentiometer R81 to set voltage to -10 V ± 10 mV.

5.5.12 Adjustment of A/D converter

This is only necessary if the voltage adjustment does not tally with the actual value displays or if the A/D converter or N19 is replaced.

- Enter service mode 5 (measured value, overtemperature sensor).
- Connect a resistance decade to the sensor simulator for overtemperature and switch to R_{EXT.}
- Connect digital voltmeter to X4/8 and X4/21 (GNDA).
- Set voltage via resistance decade to 0.000 V + 1 mV. *27.72 K.*
- Set zero point of A/D converter with potentiometer R 75 such that the actual value display just switches from 0.0 °C to 0.1 °C.
- Set voltage to 5.000 V ± 1 mV. *1.819 K*
- Connect a further resistance decade to the soldering terminals of R71 and adjust such that the actual value display just switches from 49.9 °C to 50.0 °C. *43.2 K* *MAKE SURE IN SERVICE mode*
- Measure out the set resistance and solder in the resistance value determined (if necessary, series or parallel connection of 2 resistors). *5*
- Check adjustment at
 - 5.00 V = 50.0 ± 0.1 °C
 - 4.00 V = 40.0 ± 0.1 °C
 - 3.00 V = 30.0 ± 0.1 °C
 - 2.00 V = 20.0 ± 0.1 °C
 - 1.00 V = 10.0 ± 0.1 °C
 - 0.00 V = 0.0 ± 0.1 °C

Remove Dummy sensor.

DrägerService
Medizintechnik
Medical Systems


5.6. Adjustment and functional control of humidity sensor

5.6.1. Necessary test tools

- Temperature measuring unit, compl. order no. 79 01 148
- Humidity and temperature sensor order no. 79 01 476

The accuracy of the humidity and temperature sensor must be guaranteed by regular check and adjustment with the adjustment set (order no. 79 01 477) *this checks*

5.6.2. Humidity sensor, complete

5.6.3. Adjustment and functional control when assembled

- Unscrew cover from climate sensor
- Connect climate sensor without cover to the Incubator.
- Fasten humidity sensor 79 01 476 max. 60 mm away from the humidity sensor of the Incubator.

For this, lay the connecting cable of the humidity sensor over the holder of the double side wall and fix humidity sensor at the temperature sensor (black disk) of the climate sensor using adhesive tape or wire.

- Seal the air gap which exists, because the cover of the climate sensor is unscrewed, with adhesive tape.
- Set Incubator to following desired values:

Air temperature: 33.0 to 36.0°C

Humidity: 65 to 75% rel. humidity

The indication of the Testoterm humidity sensor is decisive for the humidity. In case of deviations the desired value at the Incubator is to be changed.

- The desired value at the Incubator must be constant for at least 30 minutes before performing adjustment and functional check.
- Set humidity display at the Incubator to the value of the Testoterm humidity sensor with poti R 2. When performing a functional check the permissible deviation is 6% rel. humidity.
- Secure adjustment poti R 2 with safety lacquer.

5.7. Fresh air suction

A functional check of the fresh air suction is always necessary, if:

- a) with external oxygen supply the O_2 -concentration is not attained in the Incubator.
- b) the basic housing or filter housing 2M 19 592 with metal 2M 19 573 is renewed.

By twisting of the metal 2M 19 573 at the filter housing 2M 19 592 the diameter in the suction port behind the fresh air filters can be changed, this makes an adjustment possible.

By using a sharp thing the metal can almost be changed without disassembling the filter housing. If not, remove filter housing, reduce suction opening and when re-assembled enlarge it so long until the desired value of the fresh air suction is attained.

5.7.1. Functional check and pre-adjustment

Measurement with a Volumeter at the bore of the metal 2M 19 538 before the fresh air filters.


Assumption: New fresh air filters

Desired value of fresh air suction = 30 ± 3
-5 L/min.

5.7.2. Final adjustment

Assumption: new fresh air filters

Test set-up:

This test set-up makes a very accurate adjustment of the fresh air suction possible, without inaccuracy in measurement by power loss of a measuring instrument. For this, measure by changing the flow at the second open bore of the buffer volume with a smoke tube or with the smoke of a glowing cigarette, when the actual value of the fresh air suction falls below or is exceeded. The test value of the fresh air suction is 30 ± 3 L/min.

7. Trouble-shooting with the Penbug system

Basic setting of Penbug system: S1 = 2
S2 = 1
RAM
Penbug

7.1. Read-out of fault list

In addition to keyboard interrogation, the fault list can also be read out via the Penbug.

Work sequence: Basic Penbug setting
Connect and switch on Gercom and
Penbug
Switch on incubator

Press "return" button:

=> > - - - - enter start address 0001

Address \$0001 corresponds to fault 1,

Address \$000A corresponds to fault 10 etc

The address range can be run through with the button " ". A memory content of \$FF corresponds to a fault message.

7.2. Actuation of air heating, water heating, valve 1 and valve 2

Prerequisite:

The climate sensor with O₂ sensors must have been connected with adequate voltage. There must be no overtemperature currently present. The Incubator 8000 sensor simulator (79 01 240) can also have been connected in place of the climate sensor.

Actuation is effected via the PIA D1 (address \$8X10 - \$7X13):

PORT B	=	Outputs
CA 2	=	Output = low (voltage monitoring)
CB 2	=	Output = high (voltage monitoring)
PB 0	=	Air heating (High = on)
PB 1	=	Water heating (High = on)
PB 2	=	Valve 1 (High = on)
PB 3	=	Valve 2 (High = on)
PB 4	=	Overtemperature reset (high = overtemperature reset)

The following watchdogs (monos) must be retriggered:

WD 1 CPU	=	\$9000	
WD 2 CPU	=	\$9001	
Mono 1	=	\$8X3X	(INOP-LED)
Mono 2	=	\$8X4X	
Mono 3	=	\$8X5X	

The program is stored in the battery-buffered RAM of the CPU.

DrägerService
Medizintechnik
Medical Systems

Address (hex)	Memory unit	Remarks
1. CPU Watchdogs retriggeren		
00	F7	STA B
01	90	Address WD1 CPU
02	00	
03	F7	STA B
04	90	Address WD2 CPU
05	01	
2. PIA D1 initialisieren (complete)		
06	C6	LDA B
07	40	01000000 -> without Bit 6 all inputs
08	F7	STA B
09	88	Address Port A
0A	10	
0B	C6	LDA B
0G	34	XX1101XX Low Output (CA 2)
0D	F7	STA B
0E	88	Address Control register A
0F	11	
10	C6	LDA B
11	FF	11111111 -> all outputs
12	F7	STA B
13	88	Address Port B
14	12	
15	C6	LDA B
16	3C	XX1111XX High Output (CB 2)
17	F7	STA B
18	88	Address Control register B
19	13	

3. Actation of Peripherie (heaters and valves)

1A	C6	LDA B
1B	1F	00011111
		<pre> Air heater Water heater Valve 1 Valve 2 Overtemperature- Reset </pre>
1C	F7	STA B
1D	88	Address Port B
1E	12	

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Over Test-Programm

The heaters and valves can be addressed individually by changing the memory content of address \$001B.

1B	11	Air heating only
1B	12	Water heating only
1B	14	Valve 1 only
1B	18	Valve 2 only

The check-back signals from the periphery (see also 6.1., description of PIA ports), which can be measured directly with an oscillograph or a voltmeter at the PIA ports, are very often of importance for trouble-shooting.

7.3. Actuation of display elements

This test program can be used to switch all the display elements of the individual displays on or off.

Exception: The mains failure LED is not actuated by the processor system.

Addresses: Air display = \$8X00, \$8X01
 Skin display = \$8X02, \$8X03
 O₂ display = \$8X04, \$8X05
 Humidity display = \$8X06, \$8X07
 INOP-LED = \$8X3X

The test program is stored in the battery-buffered RAM of the CPU.

Example: All display elements of the air display light up with the exception of the INOP and the mains failure LED:

All displays can be blanked by altering the memory content of address \$21 from \$00 to \$FF. The INOP-LED lights up if the address \$8833 is not addressed (e.g. by changing to \$8844 or by inserting NOP commands). The other modules can be addressed with this test program by changing the module addresses.
For example:

Air module	-	O ₂ module
Address \$8X00	-	\$8X04
Address \$8X01	-	\$8X05

Address (hex)	Memory unit (hex)	Remarks
---------------	-------------------	---------

1. All Watchdogs retrigger

00	F7	STA B
01	90	Address WD1 CPU
02	00	
03	F7	STA B
04	90	Address WD2 CPU
05	01	
06	F7	STA B
07	88	Address Mono 3
08	55	(INOP-LED)
09	F7	STA B
0A	88	Address Mono 1
0B	33	(Horn)
0C	F7	STA B
0D	88	Address Mono 2
0E	44	(Horn)

2. Initialization of display module and actuation of display elements

	Air	Skin	O ₂	Humidity	
0F	C6				LDA B
10	0C				= 16x8bit character display-left entry + encoded scan-sensor-matrix
11	F7				STA B
12	88				Address Air-Display
13	01	03	05	07	
14	C6				LDA B
15	20				00100000 = clock-time
16	F7				STA B
17	88				Address Air-Display
18	01	03	05	07	
19	C6				LDA B
1A	90				10010000 = write display
1B	F7				STA B
1C	88				Address Air-Display
1D	01	03	05	07	
1E	86				LDA A
1F	08				= Counter
20	C6				LDA B
21	00				= all Display on
22	F7				STA B
23	88				Address Air-Display
24	00	02	04	06	
25	4A				DEC A (counter-1)
26	26				BNE
27	F8				
28	7E				JMP
29	00				Return to address \$0000
2A	00				

Circuit-Diagramms

Circuit and Assemble-Diagramms

- **Mainunit Incubator 8000**
- **Valveblock**
- **Aggregat**
- **Electronicinsertion Air + Air/Skin**
- **PC Motherboard**
- **PC Powersupply**
- **PC Analog Air + Air/Skin**

- **PC CPU**
- **Indicator + Keyboard**
- **Climatesensor**

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

Drägerwerk AG Lubeck

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Original of the document/contract is kept in the archive of the company.

Use of the document / contents is limited without express authority. All rights reserved.
DRÄGERWERK AG, Lübeck

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

Hugo Kkmann, Kiel 03 6663 83

Use of this document / contents is forbidden without express authority. All rights reserved.
DRÄGERWERK AG, Lübeck

Use of the document / contents is forbidden without express authority. All rights reserved.
DRÄGERWERK AG, Lübeck

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Use of the document/contents is forbidden without express authority. All rights reserved.
without special instruction. All rights reserved
without express consent. All rights reserved.
Without express authority. All rights reserved. DRÄGERWERK AG, Lübeck