

1. Contents

1.	Contents	1
2.	Preface	2
3.	Important remarks	3
3.1	Safety aspects	3
3.2	Data registration.....	3
4.	Technical data	5
4.1	General technical data.....	5
4.2	Environmental conditions.....	5
4.3	Transport and storage	5
5.	Circuit description	6
5.1	Power supply	6
5.2	Patient Current Routing	6
5.3	Vacuum	6
5.3.1	Control Loop	6
5.3.2	Pulsed Vacuum.....	6
5.3.3	End of Treatment.....	7
5.3.4	Water management.....	7
6.	Safety inspection list (yearly check-up).....	8
7.	Hardware labels.....	9
7.1	Upgrading the hardware code.	9
8.	Replacement procedures	10
8.1	Opening the unit.....	10
8.2	Replacement of power supply PCB	11
8.3	Replacement of the Main PCB	11
9.	Upgrades	12
9.1	Vacuum pump switches of after 2 minutes even when current is flowing	12
10.	Spare parts list.....	13
11.	Schematics and PCB layout.....	13

2. Preface

The Phyaction V is a Uniphy “low noise” vacuum unit with following characteristics :

- Vacuum unit with 2 independent channels.
- 2 or 4 poles
- Electronic vacuum control
- Continuous and pulse mode

The Phyaction V unit can be connected to the following Gymna electrotherapy units :

- Phyaction C
- Phyaction E
- Phyaction I

As a member of the Uniphy family it is built according to the requirements of the MDD security standards. All functions of the unit are controlled by a microprocessor, ensuring a high degree of reliability and safety.

This service manual gives a complete and accurate technical picture of the Phyaction V. In doing so, it will hopefully help you to reach your goal: to form a correct diagnosis and to solve the client's problem as thoroughly as possible.

If you have any questions or if you need additional information about this manual or about the use of the unit, please do not hesitate to contact us.

3. Important remarks

3.1 Safety aspects

To understand and practice all procedures described in this manual a good technical background is a must.

GymnaUniphy cannot be held responsible for any actions executed on the unit by unauthorised persons, or for executing any procedures not prescribed in this manual.

All information in this manual has a GymnaUniphy nv. copyright.

3.2 Data registration

The distributor must be able to provide the following data for each unit:

Instrument data: Part and serial number.

Gymna has the original configuration of each unit.

The **Phyaction V** has following critical parts:

- Ref 324597 : Main PCB Phyaction V
- Ref 138103 : Pump Vaco 200
- Ref 324826 : Power supply SMPS Vaco 200

All these parts have a serial-or batchnumber that has to be filed in case of a configuration change :

- Customer data: Name, full address and date of delivery.
- Service activities: All service activities must be filed.
If any critical parts are changed, we also expect the distributor to file the new serial numbers of these part(s).

To file this data, a data registration document is available on next page.

Data registration document

Information

Distributor name	Distributor address

Customer name	Customer address

Instrument name:	Serial number:
Date of delivery:	

Service activities

Date:		
Error complaint:		
Service action		
Critical spare parts exchange list		
Critical part name	Old serial number	New serial number

Date:		
Error complaint:		
Service action		
Critical spare parts exchange list		
Critical part name	Old serial number	New serial number

4. Technical data

This is a summary of the technical data as described in the user manual.

4.1 General technical data

● Dimensions Phyaciton V (w x h x d)	267 x 95 x 270 mm
● Weight Phyaciton V	3,5 kg
● Weight including accessories	4,6 kg
● Mains voltage	100 - 240 VAC, 50 - 60 Hz
● Maximum power consumption	30 VA
● Safety class	Class I (earthed socket required)
● Insulation	Type BF (floating patient circuit)
● Fuses	2 x T2AL250V
● Volume water reservoir	± 180 ml
● Working pressure continuous vacuum	38 - 320 hPa
● Working pressure pulsation vacuum	58 - 480 hPa
● Vacuum rhythm	1,5/1,5 - 1,5/4,5 s (on/off time)

4.2 Environmental conditions

● Temperature:	+10 °C to +40 °C
● Relative humidity	30% to 75%
● Atmospheric pressure	700 hPa to 1060 hPa

4.3 Transport and storage

● Transport weight	6,5 kg
● Storage temperature	-20 °C to +60 °C
● Relative humidity	10% to 100%, including condensation
● Atmospheric pressure	200 hPa to 1060 hPa
● Transport classification	Single piece by mail

5. Circuit description

5.1 Power supply

The device is mains powered. It sports an all-range off-line switched mode power supply , meaning it can be used all over the world without having to select a working voltage or frequency. But of course the the plug of the power cord has to fit into the local outlet type.

5.2 Patient Current Routing

The electrotherapy outputs of the main unit are connected to the vacuum unit, at the rear, and the electrodes are connected to the front of the vacuum unit. As long as the vacuum is not used the traditional electrodes function as before, even when the vacuum unit is not switched on. As soon as the vacuum is turned on the patient current is switched from the traditional electrodes to the connectors for the hoses at the other end of which the suction cups are. In the hoses is a wire that feeds the patient current to a metal plate inside the cup. Contact with the skin of the patient is established by a moisturized sponge. The switchover always happens for channel A but the operator can add channel B by pressing a button that toggles the choice between traditional and vacuum electrodes.

5.3 Vacuum

When the vacuum control is turned up, air is pumped out to create the vacuum until the set level is reached. The hose connectors limit the air flow into each hose so one can put the vacuum electrodes one at a time in their place. When channel B is not used, its connectors should be plugged by connecting them to each other by a hose, otherwise the pump keeps running and deeper vacuum levels cannot be reached.

5.3.1 Control Loop

A pressure sensor measures the vacuum level and a control loop powers the pump if the set level is deeper or opens a valve to let air in when lower . In continuous mode this results in only the pump being run very lightly, making up up for the small air leakage between the cups and the skin and at the connectors.

5.3.2 Pulsed Vacuum

When the electrodes are in place the operator can pulsate the vacuum to mitigate the effect vacuum has on patients with a vulnerable skin. When pulsating, the vacuum is raised by 50% during 1.5 second and then returns for 1.5 to 4.5 seconds to the level set for continuous operation. As a result the pump will run during the 1.5 second period of deepened vacuum after which the valve will be opened until the base level is reached again.

5.3.3 End of Treatment

Prolonged exposure to the vacuum under the suction cups can lead to injuries, e.g. when the patient is forgotten by the operator (actually has happened!). To prevent this, the unit senses the patient current in both channels when the vacuum is turned up. To keep the patient circuits floating this is done by optocouplers blocking a timer from running while the vacuum is switched on. As soon as no current is flowing in both channels, this timer will start running and, on timing out after a couple of minutes, switch off the vacuum. The yellow output indicator light(s) next to the vacuum connectors will now flash until the vacuum control knob is turned back to zero.

5.3.4 Water management

Because the sponges in the vacuum cups must be moistened for conductivity, unavoidably droplets are sucked into the system. These droplets do not only contain water but also fibers from the sponges, flakes of skin and other substances present under the cups. To prevent these to reach the pump and pollute it, the air crosses a water container first where the not gaseous particles are precipitating. The now clean air is processed by the pump and released into the ambient through a silencer. After a while this container will get filled with water and if not emptied, water would reach the pump and be sprayed out through the silencer. To prevent this happening there is a level detector inside the water container. When too high a level is detected, this is shown to the operator by a blinking light next to the 'water container full' symbol on the front. A running session will not be ended by this situation but when trying to turn up the vacuum from zero for a new treatment, nothing will happen. The operator first has to empty the water container by detaching the hose at the rear side of the unit and letting out the water.

6. Safety inspection list (yearly check-up)

Customer identity

Unit:

Serial number:

Date:

Test 1: General inspection

Pass Fail N.A

- Results of previous safety inspections are available
- The user manual is present
- The device label and suppliers label are clearly readable
- The enclosure, the intensity controls, the keys and the display are undamaged
- Mains entry and power cord are undamaged
- The output connectors are undamaged
- The electrode connectors and cables are undamaged.
- All the operation functions work
- The measured output signals are correct
- All the alarm functions work.

Remarks:

Test 2: electrical safety test according to VDE 0751

Parameter	Measured	Limit	Comments	Pass	Fail
Protective earth resistance	Ω	$< 0.2 \Omega$			
Enclosure leakage current	μA	$< 1000 \mu\text{A}$			
Patient leakage current	μA	$< 5000 \mu\text{A}$			

The unit With serial number has passed/failed the test

Name service engineer

Signature

7. Hardware labels.

Each PCB (except the display PCB and the power supply) has a label with the hardware code which GymnaUniphy sticks on the PCB during production.

The main PCB consists of an eeprom in which the actual hardware code of the used PCB is saved.

The hardware code consists of 4 characters and is composed as follows (from left to right):

0...	From 1 to 9.	Increase in case of a principle technology change.	Hardware is not compatible anymore. All previous changes made in hardware shall be inherited in the new layout. The last 3 alphanumerical numbers will begin with A
.A..	A to Z	Increase in case of a principle technology change.	Components are not anymore compatible with the previous layout. All previous changes made in hardware shall be inherited in the new layout. The last 2 alphanumerical numbers will begin with A
..A.	A to Z	Increase in case of component change.	
...A	A to Z	Increase in case of a PCB track change.	

Example:

Suppose that resistor R68 changes from 22 K to 100 K.
The actual hardware version 1AAA will be upgraded to 1ABA.

7.1 Upgrading the hardware code.

Every hard-or software change will be presented in a technical bulletin. In case of a hardware change GymnaUniphy will provide the necessary components as well as a new hardware label. This label must be put on the upgraded PCB.

8. Replacement procedures

The most frequent service operations will be exchanging defective boards rather than repairing them.

When handling static sensitive devices such as PC boards of the Phyaction V, the following precautions should be observed:

- Persons should be earthed by means of a wrist strap.
- Ground all electrical equipment, work bench, soldering iron...

This chapter will explain how to open the unit and remove the PCB's for exchange.

8.1 Opening the unit


Step 1 : Remove the 2 bolt at the back of the unit



Step 2 : Take of the cover.




8.2 Replacement of power supply PCB



Unscrew the 4 screws of the power supply PCB and replace with a new PCB

8.3 Replacement of the Main PCB



The main PCB is one board in production that is divided over the front of the unit.
To remove one of these boards it is necessary to loosen the bolts of the metal vacuum outlets.

9. Upgrades

9.1 Vacuum pump switches off after 2 minutes even when current is flowing


The Phyaction V has an automatic current sense circuit. If the unit doesn't detect any current for more than 120 sec, the vacuum is switched off and the vacuum cups will be released from the patient.

We have noticed that in some cases, the vacuum is switched off, even when there is current detected. In order to solve this, the unit needs to be modified.


Step 1: open the unit

Step 2: remove the main PCB from the unit

Step 3: remove R247 from the main PCB (photo below). Be careful about the 2 adjacent capacitors.



Step 4: install a 1Mohm resistor between C233p+ and CON216p7 (photo below).



Step 5: change the hardware code into 1AAB


Step 6: close the unit and perform a functional and safety test

10. Spare parts list

138103	PUMP, PM16170-86, G/PHY.V/VACO.500
324597	MAIN PCB., PHYACTION V / VACO 200
324650	CABLE ASSY ET-IN, VACO500/PHYACTION V
324727	CABLE ASSY MAINS INLET, SMPS, 200-SERIES
324793	BOTTOM BOX PHYACTION V
324804	METAL TOP COVER PHYACTION V
324837	CABLE ASSY MAINS INLET, SMPS, 200-SERIES
352578	VACO OUTLET CONNECTOR, PHYACTION V


11. Schematics and PCB layout

See next pages





COMPONENTS MARKED WITH [1] TO BE INSERTED FOR VAC0500
 COMPONENTS MARKED WITH [2] TO BE INSERTED FOR PHYACTION V

Property of GymnaUniphy N.V. Pasweg 6A, B-3740 Bilzen		Rev.
TITLE : VACUUM I/O & UI		2.1
DOCUMENT NUMBER :	ELVA003-SD201	1
DATE 1.0 :	23/03/2004	Sheet
DATE 2.0 :	02/06/2004	1/3
DATE 2.1 :	08/12/2005	Size
LAST SAVE :	12-12-2005-18:49	A4



COMPONENTS MARKED WITH [1] TO BE INSERTED FOR VACO 500
 COMPONENTS MARKED WITH [2] TO BE INSERTED FOR PHYACTION V

Property of GymnaUniphy N.V.		Pasweg 6A, B-3740 Bilzen	
TITLE :	VACUUM CONTROL		
DOCUMENT NUMBER :	ELVA003-SD202	1	Rev. 2.1
DATE 1.0 :	23/03/2004		Sheet 2/3
DATE 2.0 :	02/06/2004	DESIGN : L.MEIJER/N.H.BEUN	Size A4
DATE 2.1 :	08/12/2005	DRAWING : L.MEIJER/N.H.BEUN	
LAST SAVE :	3-1-2006_16:51		



Property of **GymnaUniphy N.V.** Pasweg 6A, B-3740 Bilzen

TITLE : VACO 500 & PHYACTION V

DOCUMENT NUMBER : ELVA003-SD200 1

DATE 1.0 = 23/03/2004

DATE 2.0 = 02/06/2004

DATE 2.1 = 08/12/2005

DESIGN : L.MEIJER/N.H.BEUN


DRAWING : L.MEIJER/N.H.BEUN

LAST SAVE : 12-12-2005_18:26


Rev.	2.1
Sheet	3/3
Size	A4

Signature for agreement

	Name	Signature	Date
<i>Review Moderator</i>			
<i>R&D Manager</i>			



0 1000 2000 3000 4000 5000 6000 7000 8000



8000 7000 6000 5000 4000 3000 2000 1000 0

