

NEUROTHERM
RADIO FREQUENCY
LESION GENERATOR
MODEL JK3

Service Manual

Document No. R28.503

Using the NeuroTherm Stimulation Test Kit

Introduction

There are times when the cannula and thermocouples look to be properly positioned in the patient but the patient does not feel any stimulation.

The Stimulation Test Kit provides a positive test that the electrode and RF lesion generator are operating correctly. The test can be performed within the sterile field.

Preparation

Ensure that the Stimulator Test Kit is kept sterile and always available.

Contents	Red cable
	Circular test block
	Blue sterilising tray
	Instructions for use

Sterelisation Instructions

After use the Stimulation Test Kit should be cleaned as per institutional policy using any normally available cleansing agent.

Sterilize by autoclave as for porous materials. Maximum permitted temperature 140°C.

Check Thermocouple

When the thermocouple is plugged into the JK3 the red LED labelled THERMOCOUPLE will turn green. When inserted into the cannula the temperature reading on the JK3 in TEST or in LESION mode will read between 35° and 38°C.

Check the Reference Plate

Ensure that the reference plate is properly connected. With the thermocouple and cannula in place in the patient check the IMPEDANCE.

Check the Stimulation

Plug the red test lead into the JK3 and the other end into the test round block. Place the test block on the sterile trolley.

Switch the JK3 into STIMULATOIN mode, set to 100 Hz or 3 Hz, turn up the amplitude and note the meter increases up to 2 Volts approximately.

Remove the sterile thermocouple from the patient, the cannula can stay in place, and touch the end of the thermocouple onto the sterile block.

With the amplitude turned up a buzz (100 HZ) or a tick (3 Hz) will be heard.

This is proof positive that the stimulation voltage is actually being delivered to the tip of the cannula.

Repositioning the Electrode

If all is correct with the machine and electrodes then the position of the cannula is suspect. Continue to reposition until the patient feels the stimulation.

If a satisfactory threshold cannot be found turn up the stimulation voltage to 0.5 Volt (or 1 Volt) and keep the stimulation on while the needle is slowly moved around. Ensure an assistant is ready to turn the amplitude down or off as soon as stimulation is felt as it can be painful for the patient. When stimulation is felt properly test for the sensory threshold turning up from zero

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

1. CONTENTS
 - 1.1 Contents
 - 1.2 List of Figures

2. WARNINGS AND CAUTIONS

3. INTRODUCTION AND APPLICABILITY OF THIS MANUAL
 - 3.1 Introduction and applicability of this manual
 - 3.2 Summary of Equipment Revision Changes
 - 3.3 Summary of Board Revision Changes
 - 3.3.1 Revision A - All Boards
 - 3.3.2 Revision B - August 1993
 - 3.3.3 Revision C – December 1998 Model JK3
 - 3.4 Manual Updates
 - 3.4.1 Neurotherm Radio Frequency Lesion Generator Unit Manual Changes
 - 3.4.2 Record of Manual Updates carried out

4. GENERAL DESCRIPTION
 - 4.1 Specification
 - 4.2 Principle of Operation
 - 4.2.1 Theory of Neuro-Surgical Radio Frequency Lesioning
 - 4.2.2 The Active and Passive Electrodes
 - 4.2.3 The Surgical effects of Radio Frequency Lesioning
 - 4.2.4 Lesioning
 - 4.2.5 Monopolar and Bipolar Lesioning
 - 4.3 General Block Diagram
 - 4.4 Wiring Diagram
 - 4.5 Rear Connector Configuration
 - 4.6 Front Panel Connections

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

- 5.0 DETAILED DESCRIPTION OF MODULES
 - 5.1 Power Supply Unit Board No. 2800808001C
 - 5.2 Fuse Board No. 2800808001D
 - 5.3 Impedance Board No. 2800808002
 - 5.4 Stimulate Board No. 2800808003
 - 5.5 RF Generator Board No. 2800808004
 - 5.6 RF Voltage and Current Metering Board No. 2800808005
 - 5.7 Temperature Board No. 2800808006
 - 5.8 Timer and Interlock Board No. 2800808007
 - 5.9 Counter Interface Board No. 2800808008
 - 5.10 System and Layout and Wiring (2800809001)

- 6.0 CIRCUIT DIAGRAMS AND PARTS LISTS
 - 6.1 Power Supply Unit Board – 2800808001C
 - 6.2 Fuse Board – 2800808001D
 - 6.2 Impedance Board - 2800808002
 - 6.3 Stimulate Board - 2800808003
 - 6.4 RF Generator Board - 2800808004
 - 6.5 RF Voltage and Current Metering Board - 2800808005
 - 6.6 Temperature Board - 2800808006
 - 6.7 Timer and Interlock Board - 2800808007
 - 6.8 Counter Interface Board - 2800808008
 - 6.9 General Systems Schematic - 2800809001

- 7.0 SERVICE AND TROUBLESHOOTING
 - 7.1 General Service Information
 - 7.2 Disassembly and Reassembly
 - 7.3 Troubleshooting
 - 7.3.1 Basic Troubleshooting
 - 7.3.2 Further Troubleshooting

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.0 ADJUSTMENTS

- 8.1 Adjustments after component replacements
- 8.2 Impedance Board (2800808002) Adjustments
- 8.3 Stimulate Board (2800808003)
- 8.4 Voltage and Current Metering (2800808005) Adjustments
- 8.5 Temperature Board (2800808006) Adjustments
- 8.6 Timer and Interlock Board (2800808007) Adjustments

9.0 FUNCTIONAL FIELD CHECK PROCEDURES

10.0 SPARE PARTS

- 10.1 Spare Parts

11.0 EARLIER REVISIONS

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

LIST OF FIGURES

- | | | |
|------|--|---------------------------------|
| 4.1 | Output Power/Load Resistance Curves | |
| 4.2 | Maximum available Stimulate Output Voltage/Load Resistance Curve | |
| 4.3 | Simplified Block Diagram | |
| | | |
| 5.1 | Power Supply Unit - Circuit Diagram | |
| 5.2 | Block Diagram - Impedance Circuit | |
| 5.3 | Impedance Board - Circuit Diagram | |
| 5.4 | Stimulate Board - Circuit Diagram | |
| 5.5 | Typical Stimulate Waveform | |
| 5.6 | RF Generator Board - Circuit Diagram | |
| 5.7 | Typical Lesioning Waveform | |
| 5.8 | RF Voltage And Current Metering Board - Circuit Diagram | |
| 5.9 | Temperature Board - Circuit Diagram | |
| 5.10 | Relationship between Thermocouple Input and Amplifier Output Signals | |
| 5.11 | Timer and Interlock Board - Circuit Diagram | |
| 5.12 | Counter Interface Board - Circuit Diagram | |
| 5.13 | System Wiring - Circuit Diagram | |
| 5.14 | Voltages on Circuit Boards related to position of Main Function Switch | |
| | | |
| 6.1 | Power Supply Board | - 2800808001 - Circuit Diagram |
| 6.2 | Power Supply Board | - 2800808001 - Component Layout |
| 6.3 | Impedance Board | - 2800808002 - Circuit Diagram |
| 6.4 | Impedance Board | - 2800808002 - Component Layout |
| 6.5 | Stimulate Board | - 2800808003 - Circuit Diagram |
| 6.6 | Stimulate Board | - 2800808003 - Component Layout |
| 6.7 | RF Generator Board | - 2800808004 - Circuit Diagram |
| 6.8 | RF Generator Board | - 2800808004 - Component Layout |
| 6.9 | RF Voltage and Current Metering Board | - 2800808005 - Circuit Diagram |
| 6.10 | RF Voltage and Current Metering Board | - 2800808005 - Component Layout |
| 6.11 | Temperature Board | - 2800808006 - Circuit Diagram |
| 6.12 | Temperature Board | - 2800808006 - Component Layout |
| 6.13 | Timer and Interlock Board | - 2800808007 - Circuit Diagram |
| 6.14 | Timer and Interlock Board | - 2800808007 - Component Layout |
| 6.15 | Counter Interface Board | - 2800808008 - Circuit Diagram |
| 6.16 | Counter Interface Board | - 2800808008 - Component Layout |

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.17	Fuse Board	- 2800808001D – Circuit Diagram
6.18	Fuse Board	- 2800808001D – Component Layout
6.19	General System Schematic	- 2800809001 - Circuit Diagram

2.0 WARNINGS AND CAUTIONS

Warning indicates a potentially harmful situation to yourself or others.

HAZARDOUS ELECTRICAL OUTPUT: The equipment is for use **ONLY** by qualified medical personnel.

Do **NOT** under any circumstance, perform any testing or maintenance on the equipment while it is being used on a patient.

Do **NOT** use extension cords or adapters of any type. The power cord and plug must be intact and undamaged.

Should the power cord or plug become cracked, frayed, broken or otherwise damaged, it must be replaced immediately.

If the equipment has in any way suffered mechanical damage it should be returned to the Supplier for Inspection and Test before further use.

Unplug the power cord before cleaning or service.

The operator should not perform any servicing of the equipment. Any servicing should only be carried out by qualified personnel.

EXPLOSION HAZARD: Never use this equipment in the presence of flammable anesthetics.

ELECTRIC SHOCK HAZARD: Always turn the equipment off before cleaning and **DO NOT** allow **ANY** fluid to enter the ventilation holes or sockets.

ELECTRIC SHOCK HAZARD: Do not touch any exposed wiring or conductive surface while cover is off and the equipment is energised. The voltage present when the electric power is connected to the equipment can cause injury or death. **Never** wear a grounding wrist strap when working on energised equipment.

FUSE REPLACEMENT: For continued protection against fire hazard, replace only with same type and rating of fuse as displayed on the rear Serial Number Plate.

2.0 WARNINGS AND CAUTIONS (continued)

A CAUTION indicates a condition that may lead to equipment damage or malfunction

Servicing of the equipment in accordance with this service manual should never be undertaken in the absence of proper tools, test equipment and the most recent revision of this service manual which is clearly and thoroughly understood.

When repairing circuit boards, great care should be taken in handling boards as all boards contain static sensitive devices. Before servicing a board, ground yourself and the relevant tool to discharge any accumulated static charge by wearing a wrist strap and placing the board on a static mat. If a board has to be returned, use anti-static bags or containers.

The tests and repairs outlined in this manual should only be attempted by trained personnel. Unauthorised service may void the warranty of the unit.

Check the voltage rating on the rear Serial Number Plate before connecting the equipment to AC Mains Power. The equipment must never be operated at the wrong mains voltage.

Use insulated tools when adjusting the internal controls on the equipment.

When cleaning the outer casing or display panel of the equipment do not use abrasive agents or solvents.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

3.0 INTRODUCTION AND APPLICABILITY OF THIS MANUAL

3.1 Introduction and applicability of this manual

This service Manual (Document No. R28.503) gives the information required to maintain and repair the Neurotherm Radio Frequency Lesion Generator Unit Model JK3. The main body of this manual deals with the present production revision of the equipment. Differences between equipment revision are summarised in Section 3.2. Section 3.3 lists the technical changes made to the equipment.

The revision of the whole equipment is changed if such technical changes are made which make some spare parts incompatible with earlier units. The initial equipment numbering, as shown on the rear Serial Number Plate contains no revision letter (e.g. Serial No. R18-101-98) however later revisions will contain a revision letter (e.g. Serial No. R18-101-98 Rev A). If the whole machine is upgraded such that early machines cannot be easily amended, then the upgraded machines will start from a nominated serial number as indicated in Section 3.2.

Within the equipment, function units such as Printed Circuit Boards will be changed or updated from time to time, these may or may not introduce a revision of the whole equipment. Each printed circuit board contains an identity number and a serial number of the board together with the issue number - designated by a letter. The initial issue letter was A for all Boards. In all cases the spare parts order code is also the Board Identity Number and Issue Number (e.g. 2800808002).

In order to understand this manual it is necessary to have a complete understanding of the function and operations of the Lesion Generator Machine. This information can be obtained from the Operators Manual which contains full operating instructions.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

3.2 Summary of Equipment Revision Changes

Initial production revision of this equipment – Model JK2

Revision - August 1993

Serial No's from R18-009 onwards

Incorporation of Digital displays and digital timing

Upgraded Boards. Major internal rewiring

Power Supply Board	2800808001 B
Impedance Board	2800808002 B
Stimulate Board	2800808003 B
Timer and Interlock Board	2800808007 B
General System Schematic	2800809001 B

New Board - Counter Interface Board 2800808008 A

In December 1998, the equipment was changed to Model JK3.

The Serial Nos started from R18-101-98.

The major differences from Model JK2 were

- a) New case
- b) New Fuse Board 2800808001D
- c) New Stimulate Board 28008003C
- d) New Temperature Board 28008006D
- e) New Timer and Interlock Board 28008007D
- f) Modified Counter Interface Board 28008008C

The board changes added a number of features to the functions of the JK2 machine

- 1 Two more temperature limits, 70°C and 80°C were added to the existir one, 90°C on the Lesion Function.
- 2 A Pulsed RF function was added
- 3 The Stimulation Rate LED is now controlled by the Interlock Circuit
- 4 A test socket and a test block is provided to give audio confirmation th stimulation pulses are present.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

3.3 Summary of Board Revision Changes

3.3.1 Revision A - All Boards Model JK2

Initial production revision of this equipment

3.3.2 Revision B - August 1993 Model JK2

Power Supply Board 2800808001 B

5 volt supply added for Lamps and Meters
40 volt rectifier moved onto Board

Impedance Board 2800808002 B

221 and C220 added to smooth impedance display

Stimulate Board 2800808003 B

Digital Display Interface added
Interlock added to prevent stimulate output when initially switched to stimulate function

Timer and Interlock Board 2800808007 B

New Board Design to operate digital timer

Counter Interface Board 2800808008 A

3.3.3 Revision C – December 1998 Model JK3

Power Supply Board 2800808001C

5 Volt supply changed to Switch Mode
Fuses removed from Board

Fuse Board 2800808001D

The fuses for the 4 supplies are on this Board

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Stimulate Board 2800808003C

Detector Circuit added to give audio confirmation of presence of pulse

Temperature Board 2800808006D

Two extra temperature limits added (70°C and 80°C)
Pulsed RF added

Timer and Interlock Board 2800808007D

Several relays removed from Board and replaced by integrated circuits
Functions unchanged

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

3.4 Manual Updates

3.4.1 Neurotherm Radio Frequency Lesion Generator Unit Manual Changes

This is the JK3 manual applicable to the machines from serial number shown below:-

Issue No.	Page	Change	Date
1		As Issued	May 1999

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

3.4.2 Record of Manual Updates carried out

Update Number	Carried out by Name	Date
---------------	---------------------	------

4.0 GENERAL DESCRIPTION

4.1 Specifications

Size:

Width	400mm (15.7")
Height	172mm (6.8")
Depth	430mm (16.9")

Weight:

8.9kg	(19.5lbs)
-------	-----------

Electrical:

UK	240 volts 50Hz Fused 1 ampere on live and neutral
USA	110 volts 60Hz Fused 2 ampere on live and neutral

The power supply is built to Class 2 standard. The mains transformer and all mains related parts are doubly insulated from the Main Enclosure. The mains transformer has separate isolated bobbins for mains and low voltage windings. It also has a thermal fuse built into the primary winding that will fail at 125°C should overheating occur. The instrument is not connected to mains earth. (Class 2).

<u>Standards:</u>	This instrument complies with BS 5724, EC 601 – 1 and MDD 93/42/EEC.
-------------------	--

Impedance:

Measuring frequency	53KHz (\pm 5KHz)
Measuring source voltage	Less than 500 mV AC
Digital Display	0-2000 ohms (in 1ohm steps)
Meter Reading	Meter reads biological impedance
Lamp Indicator	Amber LED lights when this facility is in circuit
Self Test	500 ohms resistor built into machine

4.1 Specifications (continued)

Stimulator:

Voltage Amplitude	Continuously Variable between 0 and 2 volts. The voltage supplied is displayed on the Digital Meter in 0.01 volt steps.
Pulse Rates	3 and 100 pulses per second Accuracy $\pm 5\%$ of reading
Pulse Width	1 mSec at 100 Hz
Wave Form	Biphasic pulses. Negative Pulse leading
Lamp Indicator	Green LED flashes at pulse rate of 3Hz and is continuous at a pulse rate of 100Hz.
Safety Interlock	Stimulation voltage cannot be delivered unless the voltage control is first set to the off position; this prevents any accidental application of stimulation voltage. If the stimulation section is selected and voltage control is on, the display and output will stay at zero, and the LED will not flash.

RF Lesion Power:

RF Frequency	300KHz ($\pm 10\%$) sine wave.
Power Output	Continuously variable. Maximum power output 8 watts into 200 ohm
Voltmeter	0 - 40 RF volts (RMS)
Milliamperere	Meter 0 - 200 RF milliamps (RMS)
Self Test	Into 200 ohm dummy load resistance built into machine
Lamp Indicators	Red LED lights when this facility is in circuit. Second red LED flashes when Lesion Power is switched on.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.1 Specifications (continued)

Lesion Timer:

Range Timer counts down from preset times of 60,90, 12 and 180 seconds. Timer stops when lesion power is turned off. Timer can be restarted when lesion power is established. Turning the Function switch to any position other than "Lesion" resets the timer to its set point.

Lamp Indicator Clock indicates amount of time remaining for lesioning. At the end of the preset period an alarm tone sounds and power is removed from the needle. Tone can be stopped by turning the Lesion power to the OFF position. Timer resets to set point for the next lesion.

Audible Indicator A one second "tick" indicates the clock is running.

Temperature Monitor:

Meter Range 30°C to 100°C

Probes:

Use only Rocket-Ker Thermocouple Probes No individual adjustment necessary.

Safety Features:

Safety Cut Out Lesion Power is automatically cut off if Temperature reaches the level set by the maximum Temperature Control. Lesion Power will automatically cut back in once the temperature at the tip of the probe drops 2°C less than the maximum temperature selected.


Safety Interlock Lesion Power cannot be delivered unless the RF Power control is first set to the OFF position, this prevents any accidental application of the RF power. If the lesion power is selected and power control is on, a high pitched warning tone is given out.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.1 Specifications (continued)

Earth Leakage Data	Typical	Maximum Allowable
1	Enclosure leakage current	
	Normal	3 microamps
	Reverse	4 microamps
	Single fault condition	
	Normal	0 microamps
	Reverse	10 microamps
2	Patient leakage current	
	Normal	4 microamps
	Reverse	5 microamps
	Single fault condition	
	Normal	0 microamps
	Reverse	10 microamps
3	Patient Auxilliary Current	
	Normal	1 microamp
	Reverse	2 microamps
	Single fault condition	
	Normal	0 microamps
	Reverse	5 microamps
4	Patient leakage with applied mains	
	Normal	0 microamps
	Reverse	0 microamps
5	Patient Leakage with mains applied antiphase	
	Normal	0 microamps
	Reverse	0 microamps


Figure 4.1 Output Power/Load Resistance Curve



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 4.2 Maximum available Stimulate Output Voltage/Load Resistance Curve

Load Resistance	Stimulate Output Voltage
0	
50	1.7
100	1.9
200	2
500	2.08
800	2.10
1000	2.10



4.2 Principle of Operation

4.2.1 Theory of Neuro Surgical Radio Frequency Lesioning

Definition

RADIO FREQUENCY LESIONING IS DEFINED AS THE HEATING EFFECT PRODUCED BY THE PASSAGE OF HIGH FREQUENCY ELECTRIC CURRENT THROUGH A BODY

Surgical RF Lesioning uses the properties of high frequency (RF) current to produce controlled heating capable of producing limited lesions in or around nervous tissue. Further, this radio frequency lesion generator has the capacity of accurate temperature measurement of the actuator tip making precise temperature and time controlled lesions.

The body will conduct any electrical current, however, the amount of stimulation the tissue receives - particularly muscle and nerve fibres is generally inversely proportional to the frequency. It follows therefore that the passage of mains electricity: 240V at a frequency of 50 cycles per second (50Hz) through the body will cause intense muscle reaction, indeed if passed directly through the heart, it may cause Cardiac Arrest.

However the passage of a current of 2000V at a frequency of 2.5 million cycles per second (2.5MHz) would only cause local heating at the point of contact. High frequency currents (over 1.0MHz) give good lesion effect but are difficult to confine to wires and are prone to interfere with surrounding electrical equipment such as ECG monitors etc. For these reasons the system operates at 300KHz, as a compromise between these two extremes.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.2.1 Theory of Neuro Surgical Radio Frequency Lesioning (continued)

History

The use of radio frequency (rf) heat in lesion making has proved to be the most popular method of local destruction of bodily nervous tissue although a number of other methods have been extensively researched.

The reason that RF lesioning is preferred to other methods such as cryogenics, focused ultrasound and chemical injection is because:-

- 1) RF heat lesion volumes tend to have smooth and well-circumscribed bounding surfaces;
- 2) RF lesion sizes can be reproducibly quantified using temperature monitoring of the tip of the lesion electrodes;
- 3) unwanted side effects such as boiling, charring and sticking can be avoided by temperature monitoring;
- 4) because of the easy application of stimulation, impedance monitoring and recording, target localisation and stability of RF electrode positioning are easier to accomplish;
- 5) RF lesion electrodes are robust and easily made in different configurations to suit specific anatomical usages.

4.2.1 Theory of Neuro Surgical Radio Frequency Lesioning (continued)

Operation

RF voltage is connected between an 'active' electrode and a 'dispersive' electrode with the bodily tissue completing the electrical circuit, the RF current then flows through the tissue. The electric field created by the RF voltage creates an electric force on the ions in the tissue electrolyte causing them to move back and forth at the RF. Frictional dissipation of the ionic current within the fluid medium causes tissue heating. This is the origin of the RF heat lesion.

The temperature varies as a function of distance from the electrode. It is at between 45 - 50°C that boundary neural tissue is killed. Assuming a uniform medium, the two most important parameters in determining the temperature distribution are the temperature next to the electrodes surface and the radius of the electrode. The electrode's tip length determines the length of the lesion.

The RF heat is generated in the tissue, not in the electrode tip, and because the tip lies in the tissue, it absorbs heat from it. The electrode does not absorb much heat from the tissue, therefore at equilibrium the tip temperature nearly equals that of the hottest tissue adjacent to it.

Temperature is the fundamental lesioning parameter. By monitoring the electrode temperature and choosing the proper electrode size, more consistent and safer lesions can be made, and adverse effects which are caused by the overheating of the bodily tissue, such as boiling, charring and sticking can be avoided.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.2.2 The Active and Dispersive Electrode

As is common with most electrical circuits Radio Frequency units have two electrodes - the output or ACTIVE electrode and the PASSIVE or DISPERSIVE electrode for returning the current to the machine.

The level of heat generated is dependent upon two factors:-

- i) The amount of electrical power applied
- ii) The surface area of the electrodes in contact with the body.

The ratio between these two factors is termed: CURRENT DENSITY

It follows that the smaller the surface area of the electrodes, the greater the current density will be for a given amount of power applied. The greater the current density - the greater the localised heating effect.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.2.3 The Surgical effects of Radio Frequency

As the high frequency current is passed through the body tissue, heat is generated in the tissue around the point of contact. As the TISSUE becomes heated, the intracellular water is slowly driven out of the tissue and localised destruction of the nervous tissue occurs. Within the temperature range 45 - 50°C neural tissue is killed.

The Lesion size is critically dependent on the temperature isotherm around the electrode tip and the electrode tip size.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.2.4 Lesioning

A high frequency current is passed through the body, as all animal tissues present a natural resistance to the flow of current, heat is generated in the tissue around the point of contact. When the tissue becomes heated, the intracellular water is slowly driven out of the tissue, and the tissue dies.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.2.5 Monopolar and Bipolar Lesioning

There are two types of surgical radio frequency lesioning in use, they are termed: MONOPOLAR and BIPOLAR.

Although in general for pain work monopolar lesioning is used, there is no practical difficulty in using Bipolar if the need should arise.

An example of Bipolar operation is when lesioning the lumbar communicating branch of the posterior primary ramus. Here a second 2mm bare tip insulated probe is inserted and this probe is connected to the Dispersive Plate with a separate sterile lead. If following careful testing satisfactory differential is apparent RF power can be delivered.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.3 General Block Diagram

The Neurotherm Radio Frequency Lesion Generator consists of the following modular parts (see Figure 4.3):-

Mains Input Transformer

Fuse Board (2800808001D)

Power Supply Board (2800808001C)

Impedance Board (2800808002)

Stimulate Board (2800808003)

RF Generator Board (2800808004)

RF Voltage and Current Metering Board (2800808005)

Temperature Board (2800808006)


Timer and Interlock Board (2800808007)

Counter Interface Board (2800808008)

Front Panel Assembly and Circuit Board Rack Wiring (2800809001).


Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 4.3 Simplified Block Diagram



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3


4.4 Wiring Diagram



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.5 Rear Connector Configuration

Rear Panel Mains Input Connection



1 MAINS ON/OFF SWITCH

This is a rocker type switch, combined with an I.E.C connector socket with two 'in-line' anti-surge fuses in a single unit to BS 4265.

2 MAINS IEC CONNECTOR


The three pin plug of the mains lead must be pushed into this socket. This cannot be done incorrectly i.e with the live and neutral reversed because of the orientation of the unused earth pin.

3 FUSES

The mains transformer is protected by two in-line fuses, one on the live line and one on the neutral line. These fuses are located to the right hand side of the main on/off switch. The fuses are 20mm Anti Surge 2 amp to BS4265 for a machine supplied to the U.S. for the U.K. or Europe the fuses are 20 mm Anti Surge 1 amp to BS 4265.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

4.6 Front Panel Connection



1 LEFT HAND 4MM SOCKET


This socket is for the lead of the Dispersive Patient Pate, which should be of at least 100 square centimetres, (10 square inches).

2 RIGHT HAND 4 WAY LEMO SOCKET

This socket is for the connection of the probe to the lesion generator and for carrying the thermocouple signal.

3 CENTRAL 4 MM SOCKET

This socket is to connect the Test Block with which an audio indication of the presence of stimulation pulses can be obtained.



CONNECTOR
PIN

1	Screen
2	R.F.
3	Thermocouple +
4	Thermocouple -

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.0 DETAILED DESCRIPTION OF MODULES

5.1 Power Supply Unit (2800808001C)

Mains power is connected into the rear of the equipment and goes via two fuses to double pole ON/OFF switch. See Figure 5.1. The output of the switch is connected to the main power Isolating Transformer TR101.

The Isolating Transformer is of a double bobbin type and has four secondary windings namely:-

Winding 1	30V	1 amp
Winding 2	15V	1 amp
Winding 3	15V	1 amp
Winding 4	12V	1 amp


All Windings go directly to the Power supply Card (2800808001C).

The Power Supply Card provides the following supplies:-

- a) Approximately 40V DC for the RF Lesioning Power. The supply is used by the RF Generator Board (2800808004);
- b) +12V DC Regulated Supply for all boards;
- c) -12V DC Regulated Supply for the following Boards:-
 - Stimulate Board (2800808003)
 - Voltage and Current Metering Board (2800808005)
 - Temperature Board (2800808006)
- d) 5V DC Regulated Supply for Indicator lamps on the front panel and Digital Displays.

All four supplies are fused on the Fuse Board with 1 amp T Type Antisurge fuses.

Figure 5.1 Power Supply Unit - Circuit Diagram



REPLACE ROBBINS FOR PRIMARY AND SECONDARY WINDINGS
 W1 111 TUBUS 1-000 30V
 W2 88 TUBUS 0-21000 15V
 W3 56 TUBUS 0-21000 15V
 W4 45 TUBUS 0-21000 15V

DRAWING NO.		DATE	
2800808001		APRIL 88	
DESIGNED BY	CHKD BY	DATE	
JEM	WMC	APRIL 88	
1	Completed	2	Updated
2	Updated	3	Updated
3	Updated	4	Updated
4	Updated	5	Updated
5	Updated	6	Updated

LESION GENERATOR
POWER SUPPLY UNIT

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.2 Fuse Board (2800808001D)

Each of the supply lines 12V, 5V, -12V, 40V is independently fused. Supplies come from the Power supply Unit (2800808001C) and pass through the Fuse Board before going to the rest of the machine.

5.3 Impedance Board (2800808002)

Biological Impedance is measured using an AC Signal, by driving a constant low value current through the unknown patient impedance and measuring the voltage developed.

A block diagram of the measuring system is shown in Figure 5.2 and the Impedance Board Circuit is shown in Figure 5.3.

Referring to the circuit diagram, Timer TM201 generates a square wave output whose frequency is set by P201. This frequency is set to approximately 53KHz and the output from Pin 3 of TM201 is a 5V peak to peak waveform. This waveform is modified by Inductor L201 and capacitors C206 and C207 to give a constant voltage sine wave at resistor R206. Since R206 at 10K ohms is considerably larger than any patient impedance the current through R206 is essentially constant. This current flows through the isolation transformer (T201) and the patient. A voltage is therefore developed across the patient impedance and is amplified by Amplifier OPA 201, OPA 202 is an active filter which filters out noise and other spurious signals, the cleaned up signal is then amplified by OPA 203 and then rectified and filtered by D202, D203 and C219 before being passed to the Digital Impedance Meter via series potentiometer P204.

The Digital Impedance Meter reads the voltage across resistor R221 (750 ohms) and is set to a range of 2 volts. A reading of 1999 ohms (1.999 volts) reads 1999 on the display. Capacitor C220 is connected across the meter to smooth the impedance reading.

Figure 5.2 Block Diagram - Impedance Circuit





Figure 5.3 Impedance Board - Circuit Diagram



METER SET TO 2VOLT
FULL SCALE
DISPLAYS 2000 AT 2VOLT

ALL CAPACITOR VALUES IN MICRO-FARADS
UNLESS OTHERWISE STATED
UNLESS OTHERWISE STATED

DRAWING NO.		2800808002			
DRAWN BY		JEM		CHKD BY	WJC
DATE		APRIL 88			
1	R222	10K	3-9-91		
2	R221	EC201	Updated		
3	Bomb	To B	9/1/93		
4	Updated		10/2/94		
5					
6					

LESION GENERATOR
IMPEDANCE CIRCUIT

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.4 Stimulate Board (2800808003)

When the Stimulate function is selected on the Lesion Generator the Stimulate Board is activated. The Circuit of the Stimulate Board is shown in Figure 5.4.

Referring to Figure 5.4, Timer TM301A produces a square wave output whose frequency is set by either resistors R301A and R301B for 100Hz pulse rate or resistors R302A and R302B for 3Hz pulse rate, which set of resistors are active is set by the Rate Switch (mounted on the front panel).

The output (Pin 5) from TM301A is fed to integrator Z301A and timer TM301B.

Z301A and Z301B convert the pulse train from a single sided pulse to a biphasic pulse as shown in Figure 5.5. P301 is used to set the balance between the upward and downward sections of the waveform. C309 is used to produce a flat top to the biphasic waveforms and Resistance P302 is used to set the maximum biphasic voltage.

Timer TM301B converts the 3mS pulse width on its input to 60mS on its output to give reasonable length flashes on the Stimulate Rate LED mounted on the front panel.

A typical stimulate waveform is shown in Figure 5.5. The actual peak to peak waveform is approximately twice the indication on the front panel Stimulate level Digital Display Meter. The actual maximum available stimulate pulse height will vary slightly with patient impedance. The exact relationship is shown on Figure 4.2.

The Stimulate Board also contains interlock relays and controls for setting up the Digital Display Meter.

The potentiometer (RV2) mounted on the front panel has two windings and an integral ON/OFF switch. One winding sets the stimulate output, and the other winding provides a voltage to the Digital Display Meter. The meter is set for the 20 volt range. Full scale on the winding is set to 2 volts (adjusted by P303) and this reads 2.00 on the display.


For stimulate output to be present, for the display to read , and for Stimulate Rate LED to flash Relays RLB1 and RLB2 need to be energised. This condition only occurs if RLA at some point is de-energised. RLA can only be de-energised if potentiometer RV2 is turned fully anti-clockwise to open circuit, the ON/OFF switch.

Once RLB1 and RLB2 are energised a contact of RLB1 holds them energised irrespective of whether RLA is energised or not.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

If the active lead is connected to the Test Socket, and there is a stimulate signal greater than 1 mA, this signal is detected by transformer TF302 and amplified by Z302A and Z302B. Timer Z303 then provides a pulse to the sounder S301 which gives an audio indication of the presence of a stimulate signal.

Figure 5.4 Stimulate Board - Circuit Diagram




1	8/1/60	EC 501 28A	6X4
2	8/1/60	EC 501 27A	6X5
3	7/1/63	1A8Z6	
4	7/1/63	B 504	
5	13/1/62	C 56	
6			

DRAWING NO.		DATE	
2800808003		APRIL 68	
DRAWN BY	CHG'D BY	DATE	
JEM	MMC		

LESION GENERATOR STIMULATE CIRCUIT	
------------------------------------	--

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.5 Typical Stimulate Waveform



Waveform shown for 100 pulses/second

$V = 0.2$ Volts nominal

$T = 10\text{mS}$ for 100 pulses/second

$T = 333\text{mS}$ for 3 pulses/second

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.5 RF Generator Board (2800808004)

The RF Generator Board whose circuit diagram is shown in Figure 5.6 is essentially a simple oscillator comprising FET401, Transformer TF401 and Capacitor C402. This oscillator has an oscillation frequency of approximately 300KHz. which is set by C402.


The oscillator circuit is powered from a DC voltage which can be varied from 1.2v to 40v by altering the offset adjustment of the variable voltage regulator VREG401. The offset adjustment is the Lesion Power control mounted on the front panel.

The output signal from the RF Board is an undistorted sine wave as shown in the Figure 5.7, whose shape is unaltered by changes in patient impedance.

In order to enable the RF power to be switched off if the maximum temperature (measured at the lesion point) is exceeded, a simple opto isolated shut off circuit is incorporated on the board. This operates such that when the "OSC OFF" signal goes low, the opto isolator OPTO 401 is energised and Zener Diode ZD402 is shorted out effectively tying the gate of FET401 to a fixed voltage which inhibits the oscillations.

Neurotherm Radio Frequency Lesion Generator Service Manual Model JK3


Figure 5.6 RF Generator - Circuit Diagram



NOTES:
 1. RESISTORS: Ω, KΩ, MΩ, UNLESS OTHERWISE STATED
 2. CAPACITORS: pF, nF, μF, UNLESS OTHERWISE STATED
 3. ALL CAPACITOR VALUES IN MICRO-FARADS
 4. UNLESS OTHERWISE STATED

DRAWING NO. 2800808004	1	2	3	4	5	6
	DRAWN BY JEM	CHKD BY JFC	DATE APRIL 88			
	LESION GENERATOR RF GENERATOR CIRCUIT					
	MORGAN AUTOMATION LTD					
	MAY 1999					
	5-12					

Figure 5.7 **Typical Lesioning Waveform**



Frequency 300KHZ

Max RMS Voltage 40V

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.6 RF Voltage and Current Metering Board (2800808005)

The Circuit Diagram of the RF Voltage and Current Metering Board is shown in Figure 5.8.

The Board has two basic functions, these are:-

a) to convert RF Current and Voltage to DC levels so they can be monitored by the DC Meters,

and


b) to provide a 200 ohm load under test conditions.

Relays RL501 and RL502 are normally de-energised, and the Active and Dispersive Signals from the RF Generator Board are internally connected across the 200 ohm load R508 and R509. When the Lesion Generator is in the "Lesion" condition, 12 volt is applied to the board on Pins 24A and C, and the two relays are activated allowing the RF Generator Board Signal to be connected to the patient.

Regarding the metering side of the circuit, the current transformer TF501 monitors the in-line current on the active line. 1/28 of it is rectified by D501 and D502, the DC output is amplified by a variable gain amplifier OPA 501B. The gain of the amplifier is set so that for 200mA current, the 1mA full scale meter reads 200mA. In a similar way TF502 measures a proportion of the current flowing between the active and dispersive lines. The actual impedance between the two lines is in the order of 1600 ohms, made up of C501, R501 and TF502, so about 25mA flow through this impedance chain. 7/28 of this current is rectified by D503 and D504 and amplified by variable gain amplifier OPA 501A. The gain of the amplifier is set so that 40v between Active and Dispersive gives 1mA (40v) full scale on the RMS Voltmeter.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.8 RF Voltage and Current Metering - Circuit Diagram



NOTES:
ALL CAPACITOR VALUES IN MICRO-FARADS
ALL RESISTOR VALUES UNLESS OTHERWISE STATED
TF 501 PRIMARY 1.00U HAN
TF 502 PRIMARY 1.00U HAN
SEC 25 TURNS
BOTH ON B7C 5E75346 (PK 5008) FORMER

DRAWING NO.		1	2	3	4	5	6
DRAWN BY		JCP					
CHECKED BY		JPC					
DATE		APRIL 88					
DRAWING NO.		2800308005					
TITLE		LESION GENERATOR RF VOLTAGE & CURRENT METERING					

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.7 Temperature Board (2800808006)

The temperature measuring section of the Lesion Generator is accommodated on the Temperature Board, the circuit of which is shown in Figure 5.9.

The external thermocouple (Type K - chromel-alumel) is connected to thermocouple amplifier with integral cold junction (Z601). The connection is via multipole low pass filter made up of R604, R605, R606, R607 and C604, C606, C609, C610, C611, which has an attenuation of over 100dB at 300KH enabling temperature to be taken during the application of RF Voltages. The amplifier provides both cold junction compensation and linearisation of the incoming signal as shown in Figure 5.10. The output from the amplifier is approximately 10mV/°C.

Output from the thermocouple amplifier goes via amplifier Z 603 when an offset voltage of approximately 300mV is added (via Potentiometer P601) to cater for the temperature meter being scaled from 30°C. A gain control (Potentiometer) P602 is also provided to set the 100°C point on the meter.

The output voltage from Z 603 also provides one of the inputs to comparator Z 604. The other input to the comparator is provided by a voltage from the resistor chain comprising R612, R613 and P603. With the temperature limit switch set to 90°C P603 is set so that when the meter is reading approximately 90°C, the comparator output goes high, TR601 switches on and the EXCESS TEMPERATURE signal goes low. This signal is used to disable the RF Lesioning signal generated by the RF Generator Board.


Similarly P605 and P606 control the excess of 70°C and excess of 80°C settings.

Relay RL601 is energised in the RF Pulsed mode and P604 sets the 42°C limit. Z604 provides the Pulse Mode envelope which is used to drive the Excess Temperature signal and hence allow the RF signal from the RF Generator Board to be provided in controlled bursts.

The Thermocouple Amplifier Z601 provides an output which indicates whether thermocouple is in circuit or whether a thermocouple is open circuit. This output is used to drive relay RL602 which indicates via a Green/Red LED on the Front Panel.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.9 Temperature Board - Circuit Diagram



DRAWING NO		2800B08005	
TITLE		LESION GENERATOR TEMPERATURE CIRCUIT	
1	TRAP POINTS added	DATE	APRIL 88
2	Rechecked	CHKD BY	MTC
3	Recorder Output No	DRAWN BY	JCH
4	3 Wire	DATE	APRIL 88
5	2 Wire	DATE	
6		DATE	

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.10 Relationship between Thermocouple Input and Amplifier Output Signals

Thermocouple Temperature °C	Thermocouple Voltage mV	Output Voltage from AD 595 in mV
0	0	2.7
10	0.397	101
20	0.798	200
30	1.203	300
40	1.611	401
50	2.022	503
60	2.436	605
70	2.851	707
80	3.266	810
90	3.681	912
100	4.095	1015

5.8 Timer and Interlock Board (2800808007)

The Timer and Interlock Board is shown in Figure 5.11.

The Board contains two main sections which relate to:

- a) One second timing
- b) Machine interlocks

The one second timing pulses are generated by the Programmable Timer Z702, the one second period being set by potentiometer RV701. Output from Z702 is buffered by transistor T701 and drives the Tick relay on the Counter Interface Board.

The Timer is disabled if the Lesion Power Control is off or if the counter registers zero energising RLE. Both these conditions put an 'H' on Z703 Pin 13.


The Interlock part of the circuit consists of Relays A,D,G and Integrated Circuits Z704 and Z705. RF Active and Dispersive lines from the RF Generator Board and Metering Board come into the Interlock Board on Pins 30C (Active) and 32C (Dispersive). When Relay A is activated, these signals are connected to Pins 30A and 32A which connects through to the Active and Dispersive Sockets on the Lesion Generator. The Interlock circuit is controlled by the switch on the Lesion Power Potentiometer which is open when the potentiometer is in the off position.

When the Lesioning Function is selected on the main Function Switch + 12 volts is provided to the Board, and when Relay G is in an energised state - the power potentiometer is not in its off position - this condition disables Relay A and sounds the alarm.

Once the Lesion Power Potentiometer is switched off Z704 causes the Alarm to switch off Relay A then become energised allowing signal through to the external sockets of the Lesion Generator. Relay D is also energised enabling the signal from the LED Flasher Z705 to flash the Lesion Power Lamp. Relay G switches off the impedance display if the Lesion Power Potentiometer is away from its OFF position. This situation occurs when the main Function Switch is in the Test position.

Relay RLD is also used to disable the RF output of the Generator. RLD needs to be energised for the machine to give power to the patient.

Figure 5.11 Timer and Interlock Board - Circuit Diagram



DRAWING NO.		2800308007	
ORIGIN BY		MJC	
DATE		APRIL 08	
1	Number of PCB's	2	Revision
2	3	3	4
3	5	6	7
4	8	9	10
5	11	12	13
6	14	15	16
7	17	18	19
8	20	21	22
9	23	24	25
10	26	27	28
11	29	30	31
12	32	33	34
13	35	36	37
14	38	39	40
15	41	42	43
16	44	45	46
17	47	48	49
18	50	51	52
19	53	54	55
20	56	57	58
21	59	60	61
22	62	63	64
23	65	66	67
24	68	69	70
25	71	72	73
26	74	75	76
27	77	78	79
28	80	81	82
29	83	84	85
30	86	87	88
31	89	90	91
32	92	93	94
33	95	96	97
34	98	99	100

LESION GENERATOR
TIMER & INTERLOCK CIRCUIT

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.9 Counter Interface Board (2800808008)

The Counter Interface Board is located in the Digital Timer and provides the interfacing of signals from the Timer and Interlock Board, Timer Selector Switch and Digital Counter.


The circuit diagram of the Board is shown in Figure 5.12.

The Timer Selector Switch has three Banks, two banks are used to select the preset times for Lesioning, selection is made via diodes. The third bank de-energises RL4 during switching from one selection position to another, this causes the Load Count input of the counter to go to 5v for a short period and loads the preselected counts into the Digital Counter. When the clock on the Timer and Interlock Board is enabled, one second pulses are sent to the Counter Interface Board. These pulses cause RL1 to switch ON and OFF (to give an audible 'Tick') and RL1 contact to provide a count down pulse onto the Digital Counter.

The Digital Counter counts down to zero, and on zero energises Relay RL2. This provides pulse to the Timer and Interlock Board briefly energising relay RLE. The contacts of RL2 are used to energise RL3 on the Counter Interface Board and provide a Load pulse to the counter to reset the count to the preset value selected by the Timer Selector Switch.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.12 Counter Interface Board - Circuit Diagram



COUNTER INTERFACE BOARD 2800808008

7-1-94
UPDATED 10-5-99

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

5.10 System Layout and Wiring (2800809001)

The Wiring of the Lesion Generator is shown in Figure 5.13.

The seven boards which make up the individual functions of the machine are mounted in a card rack which is accessible by removing the rear panel.

The main function switch SW1, provides power to the various boards when particular functions are selected as detailed in Figure 5.14.

The various Meters and Displays on the machine are illuminated when particular functions that relate to them are selected.

The Counter Interface Board is located in the Digital Counter which displays Lesion Time.

The various Digital Meters are powered by a 5 volt supply. The 0 volts for this supply is floating in relation to the 0 volts for the +12,-12,+40 volt supplies.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.13 System Wiring - Circuit Diagram

Only available on request

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 5.14 Voltages on Circuit Boards related to position of Main Function Switch

	FUNCTION SWITCH				
	Off	Test	Impedance	Stimulate	Lesion
Impedance		+12V	+12V		
Stimulate				+12V -12V	
RF Generator.		40V +12V			40V +12V
V/I Metering		+12V -12V			+12V -12V
Temp.		+12V -12V			+12V -12V
Timer/Interlock		+12V			+12V
<hr/>					
<u>Meters/Lamps/ Displays</u>					
Temperature		X			X
RF Volts		X			X
RF Milliamps		X			X
Impedance		X	X		
Timer		X			X
Stimulate				X	

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.0 CIRCUIT DIAGRAMS AND PARTS LISTS

This section contains details of all circuit boards, and general wiring:-

- a) Parts Lists
- b) Circuits
- c) Board Layouts

Details are included for the following Boards:-

- | | |
|--|-----------------|
| a) Power Supply Unit Board | No. 2800808001C |
| b) Fuse Board | No. 2800808001D |
| c) Impedance Board | No. 2800808002 |
| d) Stimulate Board | No. 2800808003 |
| e) RF Generator Board | No. 2800808004 |
| f) RF Voltage and Current Metering Board | No. 2800808005 |
| g) Temperature Board | No. 2800808006 |
| h) Timer and Interlock Board | No. 2800808007 |
| i) Counter Interface Board | No. 2800808008 |
| j) System Schematic - Interwiring | No. 2800809001 |


Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.1 Power Supply Unit Board – 2800808001C

Components:-

R101	2R7	2.5W	5% Wirewound
R102	4K7	¼W	5% Carbon Film
R103	4K7	¼W	5% Carbon Film
R104	4K7	¼W	5% Carbon Film
C101	4700MF	63V	Electrolytic
C102	4700MF	25V	Electrolytic
C103	470nF	63V	Polyester
C104	470nF	63V	Polyester
C105	100nF	63V	Polyester
C106	4700MF	25V	Electrolytic
C107	470nF	63V	Polyester
C108	470nF	63V	Polyester
C109	100nF	63V	Polyester
C110	4700MF	25V	Electrolytic
C111	470nF	63V	Polyester
C112	330 MF	16V	Electrolytic
C113	100nF	63V	Polyester
Z101	BZY 93C 39R		39V Zener 25W (on Heat Sink)
REG 101	7812	12V Regulator on Heat Sink	
REG 102	7912	-12V Regulator on Heat Sink	
REG 103	MIC 4576	5V Regulator on Heat Sink	
REC 101	KBPC 602	6A Bridge Rectifier	
REC 102	KBPC 602	6A Bridge Rectifier	
REC 103	KBPC 602	6A Bridge Rectifier	
REC 104	KBPC 602	6A Bridge Rectifier	
D101	Schottky Diode	IN5821RL	
L101	Inductor 40µH		
L102	Inductor 40µH		
EC101	32 way connector DIN 41612 Type B		
	Handle Black		

Figure 6.1 Power Supply Board - 2800808001B - Circuit Diagram




TR 101
SEPARATE BOBBINS FOR PRIMARY AND SECONDARY WINDINGS
W1 111 TUBUS 1.0MVA 50V
W2 88 TUBUS 0.810MVA 15V
W3 58 TUBUS 0.810MVA 15V
W4 45 TUBUS 0.4MVA 15V

DRAWING NO.		2800808001	
1	Connector Pin-out	1	Updated
2	Updated to 3 Vias	2	Updated
3	C Vias	3	Updated
4	Updated	4	Updated

LESION GENERATOR
POWER SUPPLY UNIT

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.2 Power Supply Board - 2800808001 - Component Layout



6.2 Impedance Board - 2800808002B

Components


R201	270	¼W 5% Carbon Film
R202	12K	¼W 5% Carbon Film
R203	4K7	¼W 5% Carbon Film
R204	1K	¼W 5% Carbon Film
R205	Link	
R206	10K	¼W 5% Carbon Film
R207	1M	¼W 5% Carbon Film
R208	100	¼W 5% Carbon Film
R209	51K	¼W 5% Carbon Film
R210	750K	¼W 5% Carbon Film
R211	75K	¼W 5% Carbon Film
R212	1K	¼W 5% Carbon Film
R213	10K	¼W 5% Carbon Film
R214	1K1	¼W 5% Carbon Film
R215	2K2	¼W 5% Carbon Film
R216	1M	¼W 5% Carbon Film
R217	100	¼W 5% Carbon Film
R218	750K	¼W 5% Carbon Film
R219	1K	¼W 5% Carbon Film
R220	10K	¼W 5% Carbon Film
R221	750	¼W 5% Carbon Film
P201	2K	25 turn Trimmer Type 94P
P202	1K	25 turn Trimmer Type 94P
P203	10K	25 turn Trimmer Type 94P
P204	2K	25 turn Trimmer Type 94P
C201	220MF	16V Electrolytic
C202	100nF	63V Polyester
C203	22MF	16V Electrolytic
C204	1200pF	160V Polystyrene
C205	10nF	100V Polyester
C206	3300pF	160V Polystyrene
C207	3300pF	160V Polystyrene
C208 A + B	390pF	Silver Mica
C209	10nF	100V Polyester
C210	100nF	63V Polyester
C211	47pF	Silver Mica
C212	100nF	63V Polyester
C213	1200pF	160V Polystyrene

6.2 Impedance Board - 2800808002B (continued)

Components
(continued)

C214	1200pF	160V Polystyrene
C215	100nF	63V Polyester
C216	47pF	Silver Mica
C217	100nF	63V Polyester
C218	100nF	63V Polyester
C219	100nF	63V Polyester
C220	1000MF	10V Electrolytic
D201	6V2	Zener Diode
D202	OA90	Diode
D203	OA90	Diode
D204	IN4148	Diode
OPA 201	CA3130E	
OPA 202	CA3130E	
OPA 203	LM380N	
TM201	NE555N	Timer
T201	Type PT4	Pulse Transformer
L201	10 μ H	Inductor
EC 201	64 way connector DIN 41612 Type B	
White Handle		

Figure 6.3 Impedance Board - 2800808002B - Circuit Diagram



METER SET TO 200V
FULL SCALE
DISPLAYS 2000 AT 200V

ALL CAPACITANCE VALUES IN MICRO-FARADS
UNLESS OTHERWISE SPECIFIED
ALL RESISTANCE VALUES IN OHMS
UNLESS OTHERWISE SPECIFIED

1	Q223	3-9-91
2	R221, C222	10/2/91
3	Board to B. 2800808002	10/2/91
4	Updated	10/2/91
5		


DRAWING NO. 2800808002

DRAWN BY: CHND BY: DATE: ...

LESION GENERATOR
IMPEDANCE CIRCUIT

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.4 Impedance Board - 2800808002 - Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.3 Stimulate Board – 2800808003C

Components:-


R301A	16K	¼W 5% Carbon Film
R301B	Not Fitted	
R302A	33K	¼W 5% Carbon Film
R302B	620K	¼W 5% Carbon Film
R303	3K	¼W 5% Carbon Film
R304	6K2	¼W 5% Carbon Film
R305	470K	¼W 5% Carbon Film
R306	560	¼W 5% Carbon Film
R307	4K7	¼W 5% Carbon Film
R308	10K	¼W 5% Carbon Film
R309	15K	¼W 5% Carbon Film
R310	2K2	¼W 5% Carbon Film
R311	150K	¼W 5% Carbon Film
R312	3K9	¼W 5% Carbon Film
R313	33	¼W 5% Carbon Film
R314	10K	¼W 5% Carbon Film
R315	10K	¼W 5% Carbon Film
R316	300K	¼W 5% Carbon Film
R317	10K	¼W 5% Carbon Film
R318	30K	¼W 5% Carbon Film
R319	27K	¼W 5% Carbon Film
R320	10K	¼W 5% Carbon Film
R321	330	¼W 5% Carbon Film
P301	10K	25 turn Trimmer Type 94P
P302	100K	25 turn Trimmer Type 94P
P303	2K	25 turn Trimmer Type 94P
C301	47MF	16V Tantalum
C302	680nF	100V Polyester
C303	10nF	100V Polyester
C304	470MF	16V Electrolytic
C305	100nF	63V Polyester
C306	10nF	100V Polyester
C307	100nF	63V Polyester
C308	10nF	100V Polyester
C309	470pF	100V Ceramic
C310	1000MF	16V Electrolytic
C311	470nF	63V Polyester
D301	IN4148	Diode
TR301	2N3704	Transistor
TM301	NE556	Dual Timer

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Z301	LM747	Dual Operational Amplifier
Z302	LM358	Dual Operational Amplifier
Z303	HCF4047	Multivibrator
RLA	12V	Double Pole Subminiature Relay
RLB1	10V	Double Pole Subminiature Relay
RLB2	10V	Double Pole Subminiature Relay
S301	Sounder	IMO
TF301	1+1:2+2	Transformer
TF302	1+1:2+2	Transformer
EC 301	64 way connector DIN 41612 Type B	
Green Handle		

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.5 Stimulate Board - 2800808003BC- Circuit Diagram



1	01510 Z301
2	01510 Z301
3	7413 1N4748
4	7413 B 5816
5	13114 C 54


DRAWING NO. 2800808003

LESION GENERATOR STIMULATE CIRCUIT

DRAWN BY CHMO BY DATE

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.6 Stimulate Board - 2800808003 - Component Layout




Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.4 RF Generator Board - 2800808004A

Components

R401	24	¼W 5% Carbon Film	
R402	2R7	¼W 5% Carbon Film	
R403	10	¼W 5% Carbon Film	
R404	5K1	¼W 5% Carbon Film	
R405	200	¼W 5% Carbon Film	
R407	1K	¼W 5% Carbon Film	
C401	470pF	Silver Mica	} CAPACITORS } IN } PARALLEL
C402A	33nF	63V Polystyrene	
C402B	33nF	63V Polystyrene	
C402C	33nF	63V Polystyrene	
C402D	Not fitted		
C403	100nF	63V Polyester	
C404	1MF	250V Polyester	
C405	10000pF	Polypropolene	
C406	100nF	63V Polyester	
D401	IN4148	Signal Diode	
FET401	IRF 530	MOSFET Transformer	
ZD401	BZY 88	15V Zener Diode	
ZD402	BZY 88	15V Zener Diode	
OPTO401	MCA 2231X	Opto Isolator	
TF401	Transformer with 3 windings		
VREG401	LM317T	Variable Voltage Regulator	
Heatsink	Type 07WN00750A200		
Red Handle			
EC401	32 way connector DIN 41612 Type D		

Figure 6.7 RF Generator - 2800808004A - Circuit Diagram




NOTES
ALL RESISTORS 0.25W UNLESS OTHERWISE STATED
ALL CAPACITORS 50V UNLESS OTHERWISE STATED
TF401/5 TURNS 1000
TF401/2 TURNS 1000
TF401/3 TURNS 1000
TF401/4 TURNS 1000
TF401/5 TURNS 1000
UNLESS OTHERWISE STATED

DRAWING NO.		2800808004			
DRAWN BY		JEM		DATE	APRIL 88
CHECKED BY		MJC		DATE	
1	2	3	4	5	6

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.8 RF Generator Board - 2800808004 - Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.5 RF Voltage and Current Metering Board - 2800808005A


Components

R501	100	¼W 5% Carbon Film
R502	1K2	¼W 5% Carbon Film
R503	4K3	¼W 5% Carbon Film
R505	1K2	¼W 5% Carbon Film
R506	3K3	¼W 5% Carbon Film
R508	100	¼W 5% Carbon Film
R509	100	¼W 5% Carbon Film
P501	500	25 turn Trimmer Type 94P
P502	500	25 turn Trimmer Type 94P
C501	330pF	Silver Mica 630V
C502	100nF	63V Polyester
C503	100nF	63V Polyester
C504	100nF	63V Polyester
C505	100nF	63V Polyester
D501	IN4148	
D502	IN4148	
D503	IN4148	
D504	IN4148	
D505	IN4148	
OPA501	TL072N	
TF501	1:28 Transformer	on FX3008 Former
TF502	7:28 Transformer	on FX3008 Former
RL501	FBR211 24V	Single Pole Relay
RL502	FBR211 24V	Single Pole Relay
EC501	32 way connector	DIN 41612 Type D

Black Handle

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.9 RF Voltage and Current Metering Board 2800808005A - Circuit Diagram




NOTES
ALL CAPACITOR VALUES IN MICRO-FARADS
UNLESS OTHERWISE STATED
RESISTOR VALUES IN OHMS UNLESS
OTHERWISE STATED
TF 501 RELAY 11 1/2" X 1 1/2" X 1 1/2"
TF 502 RELAY 7 1/2" X 1 1/2" X 1 1/2"
BOTH ON S7C 527244 (FX 3006) FORMER

DRAWING NO.		1	2	3	4	5	6
2800808005		1	2	3	4	5	6
DRAWN BY		JEM		DATE		APRIL 88	
CHECKED BY		JMC		DATE		APRIL 88	
LESION GENERATOR RF VOLTAGE & CURRENT METERING							

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.10 RF Voltage and Current Metering Board - 2800808005 - Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.6 Temperature Board – 2800808006D

Components


R601	1K	¼W 5% Carbon Film
R602	6K8	¼W 5% Carbon Film
R603	82K	¼W 5% Carbon Film
R604	2K2	¼W 5% Carbon Film
R605	2K2	¼W 5% Carbon Film
R606	2K2	¼W 5% Carbon Film
R607	2K2	¼W 5% Carbon Film
R608	12K	¼W 5% Carbon Film
R609	3K	¼W 5% Carbon Film
R610	1M	¼W 5% Carbon Film
R611	22K	¼W 5% Carbon Film
R612	6K8	¼W 5% Carbon Film
R613	3K9	¼W 5% Carbon Film
R614	5K1	¼W 5% Carbon Film
R615	Not fitted	
R616	10K	¼W 5% Carbon Film
R617	5K1	¼W 5% Carbon Film
R618	330K	¼W 5% Carbon Film
R619	13K	¼W 5% Carbon Film
R620	4K7	¼W 5% Carbon Film
R621	10K	¼W 5% Carbon Film
R622	5K1	¼W 5% Carbon Film
R623	10K	¼W 5% Carbon Film
R624	1K	¼W 5% Carbon Film
P601	20K	25 turn Trimmer Type 94P
P602	10K	25 turn Trimmer Type 94P
P603	1K	25 turn Trimmer Type 94P
P604	2K	20 turn Trimmer Type 94P
P605	2K	20 turn Trimmer Type 94P
P606	10K	20 turn Trimmer Type 94P
C601	100nF	63V Polyester
C602	100nF	63V Polyester
C604	100nF	63V Polyester
C605	10nF	100V Polyester
C606	10nF	100V Polyester
C607	100nF	63V Polyester
C608	100nF	63V Polyester
C609	100nF	63V Polyester
C610	10nF	100V Polyester
C611	10nF	100V Polyester

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

C612	2.2MF	63V Polyester
C613	10nF	100V Polyester
C614	100nF	63V Polyester
D601	IN4148	Diode
D602	IN4148	Diode
D603	IN4148	Diode
D604	IN4148	Diode
D605	TC04BCZM	1.26V Precision Bandgap Voltage Reference
TR601	2N3705	Transistor
TR602	BC337	Transistor
Z601	AD595CQ	Thermocouple Amplifier
Z602	LM311N	Voltage Comparator
Z603	TL071	Operational Amplifier
Z604	NE555N	Timer
RL601	Double Pole Relay 12V	
RL602	Double Pole Relay 12V	
EC601	64 way Connector DIN 41612 Type B	
Blue Handle		

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.11 Temperature Board – 2800808006D - Circuit Diagram




DRAWING NO		2800808006			
DRAWN BY		JEM		CHKD BY	WMC
DATE		APRIL 88			
1	TRF Points added				
2	Rechecked	3/9			
3	Recorder Output - N1				
4	Printer Output - N1				
5	Printer Output - N1				
6	Printer Output - N1				

LESION GENERATOR
TEMPERATURE CIRCUIT

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.12 Temperature Board - 2800808006 - Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.7 Timer and Interlock Board – 2800808007D

Components

R701	1K	¼W 5% Carbon Film
R702	1K5	¼W 5% Carbon Film
R703	Not fitted	
R704	10K	¼W 5% Carbon Film
R705	10K	¼W 5% Carbon Film
R706	4K7	¼W 5% Carbon Film
R707	10K	¼W 5% Carbon Film
R708	33K	¼W 5% Carbon Film
R709	10K	¼W 5% Carbon Film
R710	10K	¼W 5% Carbon Film
R711	1K	¼W 5% Carbon Film
R712	3K3	¼W 5% Carbon Film
R713	3K3	¼W 5% Carbon Film
R714	10K	¼W 5% Carbon Film
R715	22K	¼W 5% Carbon Film
R716	22K	¼W 5% Carbon Film
R717	47K	¼W 5% Carbon Film
R718	22K	¼W 5% Carbon Film
R719	47K	¼W 5% Carbon Film
R720	22K	¼W 5% Carbon Film
R721	47K	¼W 5% Carbon Film
R722	1K	¼W 5% Carbon Film
R723	1K	¼W 5% Carbon Film
C701	220MF	16V Electrolytic
C702	Not fitted	
C703	Not fitted	
C704	100nF	63V Polyester
C705	100nF	63V Polyester
C706	100nF	63V Polyester
C707	100nF	63V Polyester
C708	1MF	50V Electrolytic
C709	1MF	50V Electrolytic
D701	IN4001	Signal Diode
D702	IN4001	Signal Diode
D703	IN4001	Signal Diode
D704	IN4001	Signal Diode
D705	In4148	Diode
RV 701	10K	25 turn Trimmer Type 94P

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Z701	LM 3909N	LED Flasher
Z702	4541	Programmable Timer
Z703	4001B	Quad 2 Input CMOS NOR Gates
Z704	4013B	
Z705	4001B	

T701	2N 3705	Transistor
T702	BC337	Transistor
T703	BC337	Transistor
T704	BC337	Transistor


RLA	BT Double Pole Relay	10V Coil
RLD	BT Double Pole Relay	12V Coil
RLE	BT Double Pole Relay	12V Coil
RLG	BT Double Pole Relay	12V Coil

Connector DIN 41612 Type C Connectors A+C Even Contacts

Handle - Yellow

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.13 Timer and Interlock Board – 2800808007D - Circuit Diagram



1	Number of ICS	3
2	Revision	1/12
3	Version	1/12
4		
5		


DRAWING NO. 2800808007
 DRAWN BY: CMM/7 DATE: ...

LESION GENERATOR
 TIMER & INTERLOCK CIRCUIT

⑤ WIRE NO

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.14 Timer and Interlock Board - 2800808007 - Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3


6.8 Counter Interface Board – 2800808008C

Components

R1	10K	1/8W 5%Metal Film
R2	10K	1/8W 5% Metal Film
C1	47MF	6.3V Electrolytic
D1	IN4148	Signal Diode
D2	IN4148	Signal Diode
D3	IN4148	Signal Diode
D4	IN4148	Signal Diode
D5	IN4148	Signal Diode
D6	Not fitted	Signal Diode
D7	IN4148	Signal Diode
D8	IN4148	Signal Diode
D9	IN4148	Signal Diode
TR1	2N 3704	Transistor
RL1	12V	Miniature Single Pole Relay
RL2	5V	Miniature Single Pole Relay
RL3	12V	Miniature Single Pole Relay
RL4	5V	Miniature Single Pole Relay

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.15 Counter Interface Board – 2800808008C - Circuit Diagram




COUNTER INTERFACE BOARD 2800808008C

7-1-94
UPDATED 10-5-99

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.16 Counter Interface Board - 2800808008 - Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.9 Fuse Board – 2800808001D

R101 1K 6 watt Vitreous

F101 Fuseholder with
 cover

F102 Fuseholder with
 cover

F103 Fuseholder with
 cover

F104 Fuseholder with
 cover


Fuses 1 AMP T 101 to 104
 (Bussmann)

Black Handle

Edge connector DIN41612 Type B


Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Fig 6.17 Fuse Board 2800808001D - Circuit Diagram



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Fig 6.18 Fuse Board 28008001D – Component Layout



Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.10 General System Schematic - 2800809001

Components

RV1	4K7 1W Potentiometer with integral 2 pole switch
RV2	2 gang potentiometer, each gang 1K,1W with integral ON/OFF switch
C1	1000MF 25V Electrolytic
SW1	Function Switch – 6 positions - 9 pole
SW2	Start Time Switch - 2 positions - 1 pole with spring return
SW3	Stimulator Rate Switch - 2 positions - 1 pole
SW4	Lesion Time - 4 positions - 3 pole
SW5	Max Temperature – 3 positions – 1 pole
M1	0-1mA DC Movement Meter - Temperature
M2	0-1MA DC Movement Meter - 0 to 200mA RMS
M3	0-1mA DC Movement Meter - 0 to 40V RMS
DIG 1	Digital Panel Meter 3½ digit - to display impedance
DIG 2	Digital Panel Meter 3½ digit - to display stimulate voltage
DIG 3	Digital Counter 4 digit - to display time
IND 1	Meter Illumination 6mm x 24mm 6V Festoon Bulb
IND 2	Meter Illumination 6mm x 24mm 6V Festoon Bulb
IND 3	Meter Illumination 6mm x 24mm 6V Festoon Bulb
IND 4	Green LED with Bezel - Mains
IND 5	Green LED with Bezel - Stimulate Rate
IND 6	Amber LED with Bezel - Lesion Power On
IND 7	Green/Red LED with Bezel – Thermocouple Open Circuit
D 1	IN 4002 Diode
AL 1	2 Tone Alarm Sounder

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

6.10 General System Schematic - 2800809001 (continued)

SKT 1 LEMO ECG OB-304C
SKT 2 4mm Shrouded Socket - Black
SKT3 Test Socket 4mm Black

Mains Input Connection

Corcom Power Input Unit Type 6V MIS complete with Mains On/Off Switch and 2 fuses.

Transformer

Primary 0-120V, 0-120V 50/60Hz

Double Bobbin Section

Secondary 1 30V RMS at 1A RMS
Secondary 2 15V RMS at 1A RMS
Secondary 3 15V RMS at 1A RMS
Secondary 4 12V RMS at 1A RMS

Thermal Fuse fitted in primary winding

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

Figure 6.19 General System Schematic - 2800809001 - Circuit Diagram

Diagram available on request

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.0 SERVICE AND TROUBLESHOOTING

7.1 General Service Information

There are several levels of service available these are:-

24 hour Loan Service.

A Loan machine will be despatched to any UK customer within 24 hours of receipt of order. The Loan Unit will be sent in a strong carrying case so that the unit for repair can be returned to RDG Medical via prepaid carriage. The hospital unit will be returned after repair in the same type of carrying case. The Loan Unit can then be returned to RDG Medical.

Loan Units will be charged on a daily basis with a minimum charge for carriage and return.

24 hour Board Replacement Service.

For those customers who have established EBME repair facilities, RDG Medical will provide a full PCB replacement service. RDG Medical will despatch within 24 hours of receipt or order, a fully tested replacement board for delivery the next day.

Usually field service is limited to swapping the suspect circuit board etc., the faulty board can then be repaired at RDG Medical.

If you need to seek service advice, please quote the unit serial number, product revision letter and a detailed fault description. After any replacement see Section 8.1 (adjustments) and after any service perform the functional field check procedure as outlined in Chapter 9.

Depending on the capabilities of the EBME repair facility some repairs may not be possible. In such cases the faulty parts should be sent to RDG Medical for repair.

The tests and repair outlined in this section should only be attempted by trained personnel with the appropriate equipment. Unauthorised service may void the warranty and its associated guarantees.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.2 Disassembly and Reassembly

- 1 Access to the Boards is achieved by removing the name plate panel on the back of the machine, secured by 4 M4 screws.
- 2 The machine is disassembled by:-
 - a) Removing the 2 outer M4 screws on the rear panel
 - b) Removing the 6 M4 screws along the two sides of the machine

The top housing of the machine can now be separated from the bottom housing however the cable from the card rack to the front connectors will restrict how far the two can be separated.

The cable harness length to the front panel is so designed that the top housing can be placed on its right hand side next to the right hand side of the machine, and the front panel can be connected in this position.

The front panel is connected to the main rack by a 50 way connector.

All boards within the machine are then totally accessible and can individually be removed from their rack or operated on extender boards.

It is possible to fully operate and monitor the internal workings of the machine with the cover off. There are no voltage in excess of 40V DC on the rack, cards or front panel.

Assembly is essentially reversing what was described above, however some caution is necessary. Before replacing the top housing check that all internal connectors are in place, the 50 way front panel connector has its retainers fully secured, and the cables to the front connectors are not snagged or trapped in any way.

The top housing is then placed over the base and secured with the 10 screws.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3 Troubleshooting

7.3.1 Basic Troubleshooting

This section deals with very basic Troubleshooting which should always be carried out prior to attempting major repair.

SYMPTOM	POSSIBLE CAUSE AND REMEDY
MAINS LAMP DOES NOT LIGHT	<p>Unit is not receiving power Check Mains Cable is correctly inserted into a wall socket and into the rear socket of the machine.</p> <p>Check the mains switch on wall socket and rear of machine is in the 'ON' position. Check fuses in rear of machine or fuse in mains plug (if fitted).</p>
ON IMPEDANCE FUNCTION IMPEDANCE SHOWS OPEN CIRCUIT	<p>Problems with Dispersive Plate Lead. Check the Dispersive Plate Lead is correctly connected and plugged into the Front Panel Socket.</p> <p>To test that the Impedance Function is operating correctly turn FUNCTION SWITCH to TEST position and check that displayed reading is approximately 500 ohms</p>
NO OUTPUT AT ACTIVE PROBE	<p>Problems with Ancillary Cables. Check LESION selected on Function Switch and Lesion Power control turned off and then on. Check that the Temperature Probe is correctly plugged into front panel connector.</p> <p>If in doubt, replace Temperature Probe.</p> <p>To test whether the machine is internally generating RF Power, turn FUNCTION SWITCH to TEST Position and turn Lesion Power Control clockwise. Observe meter which both should move to full scale deflection at the fully clockwise position of the control.</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.1 Basic Troubleshooting

SYMPTOM	POSSIBLE CAUSE AND REMEDY
LOW RF POWER AT ACTIVE PROBE	Problems with Dispersive Electrode or lead. Check Dispersive Plate is securely applied. If in doubt change for a new one. (If patient impedance is above 200 ohms as power is applied RF voltage will increase faster than RF Current - this can be observed on the two meters). If impedance is too high the lesion power available will not be sufficient to achieve a high temperature at the lesion point.

If this Basic Troubleshooting fails to remedy the fault either:-

- a) Report the fault to RDG Medical for a machine exchange;
- or
- b) Contact an authorised Service Engineer;
- or
- c) Carry out the further troubleshooting as detailed below

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting

Further troubleshooting or field service of this machine is intended to be done as board swapping where possible. Thus the most important part of troubleshooting is to pin-point the faulty module.

In order to carry out this further troubleshooting it is necessary to remove the name plate cover at the back of the machine, or disassemble the machine as detailed in Section 7.2.2.

Leave the Front Panel connected to the Card Rack and put the top housing on its right hand side next to the right hand side of the machine. Ensure all Circuit Boards are pushed fully home in the Card Rack.

Connect Mains Lead.

Follow through the procedure outlined below, omitting sections as necessary:-

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
<p>MAINS LAMP ON FRONT PANEL DOES NOT LIGHT</p>	<p>Switch on. Check rear fuses and mains plug fuse, check that mains voltage is available at machine. With meter and prods check that AC voltages are present at the pins of the secondary side of the mains transformer</p> <p>If these voltages are not present, either the mains input socket is faulty or transformer has blown its thermal fuse.</p>	<p>REAR FUSES MAIN INPUT UNIT</p> <p>MAINS POWER TRANSFORMER</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
	Put Fuse Board on an extender board or lead and check DC Voltages at the four 1Amp Fuses.	
	'0 volts' is available on the rear of the Board on Terminals 11,18,24,28. The '0 volts' for the 5V supply is on Terminal 11. The output side of the four fuses is near the front edge of the board.	INTERNAL FUSES
	From the top of the board, fuses are for +12V, 5V -12V, 40V.	
	If a fuse has blown, it is likely that either the power supply is faulty or another board is shorting out one of the supply lines.	POWER SUPPLY BOARD FAULTY
	Turn off machine and leave for some minutes for capacitors to discharge. Check the impedance of each power line (at the fuses) to 0 volts. For the 5V supply check to the '0 volts' of this supply (i.e. Terminal 11). If any power line is low impedance, disconnect boards one at a time to find offending board.	SHORT CIRCUIT ON CIRCUIT BOARD
	If power supply is at fault and it is obvious where the fault lies, replace component or request a replacement board.	

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
<p>DISPLAYS AND METERS FAIL TO LIGHT WHEN FUNCTION IS SELECTED</p>	<p>If all meter lights are inoperative fault is either at the transformer, power supply board or Function Switch SW1. Check transformer as above, and 5v supply on Power Supply Board. If these are intact then check wiring to Function Switch (Wire 48) or meter commons (Wire 47).</p> <p>If a single meter lamp fails to come on, check internal Festoon bulb or associated wiring.</p>	<p>POWER SUPPLY BOARD/FUSE</p> <p>FUNCTION SWITCH</p> <p>LAMP WIRING</p> <p>METER BULBS</p>
<p>RF POWER CANNOT BE GENERATED WITH FUNCTION SWITCH IN TEST POSITION</p>	<p>The fault either lies on the RF Generator Board or on the RF Voltage and Current Metering Board.</p> <p>Put RF Generator Board on an Extender Lead and check DC output voltage of VREG 401 goes from 3 to 37V when Lesion Power Control is turned clockwise, if not fault lies with either external potentiometer or Regulator.</p> <p>If DC voltage control is correct, connect an oscilloscope between Pins 10 and 12 of edge connector of RF Generator Board. Signal should be 110V peak to peak at 300KHz sine wave with Power Control Potentiometer fully clockwise.</p>	<p>RF GENERATOR OR RF VOLTAGE AND CURRENT METERING BOARD</p> <p>VREG 401 ON RF GENERATOR BOARD OR LESION POWER CONTROL</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE AREA	FAULT
	<p>If this sine wave is present check the RF Voltage and Current Metering Board. Put Board on an Extender Lead and check that RL501 and RL502 are energised when Test Position is selected on the Function Switch, i.e. check 24V is present across diode D505.</p> <p>Check 110V 300KHz sine wave is present across R508 and R509.</p> <p>If voltage is present then check metering circuit and associated wires to meters.</p>	<p>RF VOLTAGE AND CURRENT METERING BOARD</p>	
<p>RF POWER NOT AVAILABLE IN LESION POSITION</p>	<p>Connect a 200 ohm load between Dispersive socket and Temperature Probe, and check meters move up together to full scale with increasing clockwise rotation of Lesion Power Potentiometer. If they do not, check RF Generator Board and Voltage and Current Metering Boards as above. If these Boards operate correctly check 300KHz sine wave across Pins 32C and 30C on Timer and Interlock Board and between 30A and 32A on the Board. If the signal is available on the 32C and 30C Terminals but not on 30A and 32A then interlock system on the Board should be investigated. If signal is available check wiring to front connectors.</p>	<p>TIMER AND INTERLOCK BOARD</p> <p>WIRING TO FRONT CONNECTOR</p>	

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
<p>IMPEDANCE MEASUREMENT FAILS TO OPERATE</p>	<p>With Function Switch in the Impedance Position connect a 500 ohm resistor across Pins 10A and 14A of IMPEDANCE BOARD. Check whether Impedance Meter reads. If it does, then check wiring from Card to Function Switch. If Display fails to read check wiring to Display. Check voltages (0 and +12v) are available on board, if not check Function Switch. If voltages are present request a replacement Impedance Board.</p>	<p>WIRING TO FUNCTION SWITCH</p> <p>IMPEDANCE BOARD</p>
<p>STIMULATE FUNCTION FAILS TO OPERATE</p>	<p>With Function Switch in the Stimulate Position check voltage across Pins 28 and 29 of Stimulate Board Edge connector using an oscilloscope. With Stimulate rate set to 100Hz.</p> <p>Positive and Negative going pulses should be observed (1mS wide) whose height varies with the Stimulate Voltage Control. If no pulses are observed, check -12V and +12V supplies - if these are not present check wiring to Function Switch. If pulses are present but not getting to front panel connectors check wiring to Function Switch and to front panel connectors. See System Schematic Wiring Diagram Figure 6-19.</p>	<p>WIRING TO FUNCTION SWITCH</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
	If voltages are present but there is no output from board, check timer output Pins 5 and 9 of TM301A for pulses. Replace timer if necessary or request a replacement board.	TIMER TM301A STIMULATE BOARD
TEMPERATURE FUNCTION FAILS TO OPERATE	<p>With Function Switch in Test Position, connect either a probe or signal injector into the front panel connector and vary effective input temperature. If no movement is observed on meter, check whether signal is arriving at the Temperature Board by observing voltage across Pins 1 and 14. Note voltages are very small at these points (1mV per 25°C) - see Figure 5-10 for expected signal levels. If voltage is not present check wiring to front panel.</p> <p>If voltage is present, check output of Z601A at Test Point A and 0V.</p> <p>Signal should be 10mV/°C as detailed in Figure 5.10. If voltage is not present Z601 is suspect.</p> <p>If voltage is present check output of Z603. If this voltage is present, check wiring to meter.</p> <p>If the Board appears faulty first check that the +12V and -12V supplies are present before requesting a replacement Board.</p>	<p>FRONT CONNECTOR WIRING</p> <p>Z601 ON TEMPERATURE BOARD</p> <p>TEMPERATURE METERING WIRING</p> <p>TEMPERATURE BOARD</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
<p>SYSTEM FAILS TO REDUCE POWER WHEN TEMPERATURE EXCEEDS 90°C</p>	<p>This fault is either caused by the Temperature Board not sending an Excess Temperature Signal to the RF Generator Board, or the Opto Isolator (OPTO401) not operating.</p> <p>With Function Switch in Test Position, apply some Lesion Power so that RF Voltage and Current Meters read. Connect a signal injector into the front connectors and turn it up so that the temperature meter moves above 90°C.</p> <p>Meters should show decreasing output as temperature moves above 90°C.</p> <p>If meters fail to show decreasing power, check function of Z602. Pin 7 should change state when voltage on Pin 2 exceeds that on Pin 3.</p> <p>Check voltage on Pin 3 is set at correct level for 90°C Cut out.</p> <p>If Z602 operates correctly check that Pin 25 on edge connector operates. If not TR601 is suspect. If it is correct check wiring between Temperature and RF Generator Boards or replace OPTO401 on RF Generator Board.</p>	<p>Z602 ON TEMPERATURE BOARD</p> <p>TR601 ON TEMPERATURE BOARD OPTO401 ON GENERATOR BOARD</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

7.3.2 Further Troubleshooting (continued)

SYMPTOM	PROCEDURE	PROBABLE FAULT AREA
CLOCK FAILS TO START OR OPERATES INCORRECTLY	<p>If clock fails to start when the Function Switch is in the LESION position and Lesion Power Control is on, then check state of Z702 Programmable Timer and whether a pulse train of one second pulse are coming from Pin 8. Check that Pin 1 of Z703 goes high when start timer is held clockwise. Check Pin 13 is low, if not check appropriate switches.</p> <p>If clock appears to misoperate check Counter and Counter Interface Board. Replace if necessary.</p> <p>Check status of all interlocks by putting the Timer and Interlock Board on an extender card.</p> <p>If interlock conditions are correct and Board still fails to operate correctly - replace the Board.</p>	<p>LESION POWER CONTROL SWITCH</p> <p>TIMER START SWITCH Z702 OR Z703 ON TIMER AND INTERLOCK BOARD</p> <p>COUNTER AND COUNTER INTERFACE BOARD</p> <p>TIMER AND INTERLOCK BOARD</p>

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.0 ADJUSTMENTS

Adjustments are limited to 5 boards namely.

Impedance Board	(2800808002)
Stimulate Board	(2800808003)
Voltage and Current Metering Board	(2800808005)
Temperature Board	(2800808006)
Timer and Interlock Board	(2800808007)

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.1 Adjustments after component replacements

After any component replacement check from the following table which adjustments have to be performed:-

Component Replaced	Adjustments
Transformer or Power Supply Board	All controls may need adjustment on: Impedance Board Stimulate Board Voltage and Current Metering Board Temperature Board
RF Voltage or Current Meters	Voltage and Current Metering Board
Temperature Meter	Temperature Board
Impedance Display	Impedance Board
Impedance Board or any component on it	Impedance Board
Stimulate Board or any component on it	Stimulate Board
Voltage and Current Metering Board or any component on it	Voltage and Current Metering Board
Temperature Board or any Component on it	Temperature Board
Timer and Interlock Board or any Component on it	Timer and Interlock Board

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.2 Impedance Board (2800808002) Adjustments

1. Connect Mains Power to the machine.
2. Connect the Test Box to the Active and Dispersive sockets on the front of the machine and set the load to 500 ohms.
3. Turn Main Function Switch to Impedance Position.
4. Connect Impedance Board on an extender lead so that it is clear of the card rack.
5. Connect an oscilloscope between 0V (Pin 1) and Pin 3 of TM201 (one end of R204).
6. Switch on Machine.
7. Adjust Frequency Potentiometer (P201) to obtain 53KHz (19 μ S) 5V peak to peak square wave signal with ON time (5V) = 3 x OFF time (0V).
8. Connect scope between 0V and wiper of P202 (one end of R206). Adjust Measuring Level Potentiometer P202 to get 1V peak to peak 53KHz signal (rounded saw tooth) with small high frequency component. This signal sits 200mV above 0V line.
9. Connect scope between 0V and junction of D202 and D203. Adjust Amplifier Gain (P203) to get 2.5V peak to peak 53 KHz signal with no ringing for 500 ohm load.

Equivalent voltage for other loads:-

200 ohm	1 volt
500 ohm	2.5 volts
1000 ohm	4 volts

10. Adjust P204 (Adjust Meter Potentiometer) to get correct reading at 500 ohms. Observe reading at 1K ohms and slightly adjust P203 and P204 both clockwise or anti-clockwise to get best fit over impedance range. Ensure display goes overscale on open circuit output..
11. Check that amplifier does not ring or fold over by connecting 1K5 ohm load, monitoring signal on junction of D202 and D203, and disconnecting and reconnecting dispersive lead. If there is any sign of signal oscillating or if amplifier folds over, reduce the output by turning P203 clockwise and restore meter scale reading by turning P204 clockwise.
12. Lock the four trimmers with varnish.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.3 Stimulate Board (2800808003) Adjustments

1. Connect Mains Power to the Machine.
2. Connect the Test Box to the Active and Dispersive Sockets on the front of the machine and set the load to 500 ohms.
3. Switch on and turn Main Function Switch to Stimulate Position.
4. Turn Stimulate Rate Switch to 100Hz and Stimulate Amplitude Control fully clockwise.
5. Connect an oscilloscope across the output terminals of the machine.
6. Adjust set shape potentiometer on Stimulate Board (P301) to get positive and negative pulses of equal amplitude.
7. Adjust Set Gain potentiometer on Stimulate Board (P302) to get 2 volts positive and negative pulses.
8. Adjust Set Display Potentiometer (P303) to get 2 volt display on Front Panel Display.
9. Lock the trimmers with varnish

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.4 Voltage and Current Metering Board (2800808005) Adjustments

1. This board is operational in the "TEST" and "LESION" positions of the Function Switch.
2. With the Function Switch in the LESION Position connect a Test Load of 200 ohms between the Active and Dispersive sockets on the front of the machine.
3. Connect an oscilloscope across the output terminals of the machine and turn Lesion Power Control fully clockwise. Signal obtained should be a 300KHz sine wave of about 110V peak to peak.
4. Measure maximum peak to peak value and calculate RMS voltage value as below:-

$$\text{RMS Voltage} = \frac{\text{Peak to Peak Voltage}}{2.828}$$

5. Adjust P502 (Adjust Voltage) to get RMS Voltage meter to read the correct RMS Voltage on the meter (should be near full scale).
6. Calculate RMS Current as below:-

$$\text{RMS Current} = \frac{\text{RMS Voltage}}{200}$$

7. Adjust P501 (Adjust Current) to get RMS Current Meter to read correct value (should be near full scale).
8. Lock the two trimmers with varnish.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.5 Temperature Board (2800808006) Adjustments

1. This board is operational in the "TEST" and "LESION" positions of the Function Switch.
2. With the Function Switch in the LESION Position connect a Test Load of 200 ohms between the Active and Dispersive sockets on the front of the machine.
3. Connect the Temperature Board on an extender so it is clear of the Card Rack.
4. With the Function Switch in the Lesion Position inject a signal of approximately 1.2mV into the thermocouple input pins of the Active socket to give an output of 310mV when measured on Test Point A on the Board. Adjust P601 (set 30°C) to get meter to read 30°C.
5. Inject a signal of approximately 4mV into the thermocouple input pins of the Active Socket to give an output of 1025mV when measured on Test Point A on the Board. Adjust P602 (set 100°C) to get meter to read 100°C.
6. With Temperature Limit Switch in 90°C position adjust voltage injector signal to get a meter reading of 90°C. Measure voltage on Pin 2 Z602. Adjust P603 (set Cut off) to get the same voltage on Test Point B.
7. Test Cut Out control by turning down temperature by reducing voltage injector signal. Turn Lesion Power On to get full scale readings on the RF Voltage and Current Meters. Increase voltage injector signal and check that when temperature reading increases above 90°C, the RF Current and Voltage Signals decrease.
8. With Temperature Limit Switch in 80°C position set meter reading to 80°C and repeat #6 by adjusting P606.
9. Repeat above at 70°C adjusting P605.
10. Switch Mode to Pulsed RF and repeat above at 42°C with P604.
11. Lock all trimmers with varnish.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

8.6 Timer and Interlock Board (2800808007) Adjustments

1. This board is operational in the 'TEST' and 'LESION' positions of the Function Switch.
2. With the Function Switch in the 'TEST' position, set time on the Digital Counter on the Front Panel to 60 seconds.
3. Turn the Lesion Power Switch ON.
4. Start the timer and observe the countdown from 60 seconds to zero.
5. Check the time taken by the counter against a stopwatch.
6. If the two times are not in agreement, adjust RV701 as appropriate and repeat stages 3 to 5 until they are the same.
7. Repeat procedure with 120 seconds on the counter.
8. When satisfied that the two times are within one second over a 120 second period, lock the trimmers with varnish.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

9.0 FUNCTIONAL FIELD CHECK PROCEDURES

The following test should be carried out after any service to ensure proper operation of the equipment

1. Switch the equipment ON.
2. Check that the green MAINS Indicator is illuminated.
3. Connect Test Box to Active and Dispersive connectors.
4. Switch Function Switch to TEST Position, check Impedance Display reads approximately 500 ohms. Check two RF Meters and Temperature Meter light up. Turn power up and check RF Voltage and RF Current go to Full Scale on their meters.
5. Turn up voltage on Test Box and observe Temperature moves up scale on Temperature Meter. Set Temperature Limit Switch to 90°C and check as temperature increases above 90°C RF Voltage and RF Current decrease.
6. Turn Function Switch to IMPEDANCE and note that Impedance Display lights up. Select Impedance on Test Box and check that Impedance Meter displays the correct impedance. $\pm 20\%$
7. Turn Function Switch to STIMULATE and check that Green Rate LED flashes at both 3Hz and 100Hz. Check with scope on Test Box that stimulate voltage is the correct shape and frequency. Range should be 0 to 2 volts.
8. With Test Box set to 200 ohm LOAD and Lesion Power Potentiometer turned on, switch Function Switch to LESION Position. Lesion power lamp should light together with RF Voltage, RF Current and Temperature Meters. A High Pitch sound should be heard.
9. Turn Lesion Power Potentiometer off, High pitch sound should cease. Turn Lesion Power Potentiometer On, Lesion Power On, LED should flash, and RF Voltage and RF Current Meters should move upscale with clockwise rotation of the Lesion Power Potentiometer.
10. Start Timer by rotating Start Timer Switch and releasing. Turn off Lesion Power using the potentiometer, note that clock stops. Turn on Lesion Power and when clock is restarted it should start from where it left off. Let Timer run to zero, Lesion Power should switch off and a high Pitch Tone should be heard. Tone ceases when Lesion Power Switch is turned off.

9.0 FUNCTIONAL FIELD CHECK PROCEDURES (continued)

11. Monitor shape of RF Lesion signal on Test Points on Test Box. At Maximum Power signal should be 110V peak to peak sine wave at approximately 300KHz.
12. Check that 60,90,120 and 180 seconds can be set on the Clock Timer by rotating the Set Time Switch.
13. Check the leakage current of the unit. Check the condition of all cables and connectors.

Neurotherm Radio Frequency Lesion Generator Service Manual
Model JK3

10.0 SPARE PARTS

10.01 Spare Parts

NOTE: Accessories are Listed in Operators Manual.

Item	Item Description	Order No.
1A	Primary Fuse EUR 1 AMP T TYPE	2800-020*
1B	Primary Fuse USA 2 AMP T TYPE	2800-021*
2	Power Supply Fuses 1 AMP T TYPE	2800-020*
3	Mains Entry Module with Switch	2800-023
4A	Mains Cable UK	2800-024*
4B	Mains Cable EUR	2800-025*
4C	Mains Cable USA	2800-026*
5	Transformer	2800-027
6	Alarm Buzzer	2800-028
7	Power Supply Board 2800808001	2800-001*
8	Impedance Board 2800808002	2800-002*
9	Stimulate Board 2800808003	2800-003*
10	RF Generator Board 2800808004	2800-004*
11	RF Voltage and Current Metering Board 2800808005	2800-005*
12	Temperature Board 2800808006	2800-006*
13	Timer and Interlock Board 28008008007	2800-007*
14	Front Panel complete with switches and meters and displays	2800-008
15	Enclosure	2800-029
16	Replacement Lamps for Meters	2800-030*
17	Kit of Fixings	2800-031*
18	Set of 4 feet with fixings	2800-032
19	Test Box	2800-033
20	Fuse Board	2800-009

* This part is recommended for stock

11.0 EARLIER REVISIONS

This section lists any technical details relevant to earlier versions of the Lesion Generator.

{The main body of this technical manual deals with the present production revision of the equipment}.

Serial Nos.

Relevant Technical Details

NEUROTHERM RADIO FREQUENCY LESION GENERATOR
MODEL JK3
SERVICE MANUAL

12-1



*The Intelligent
Application of
Technology*

Rake Heath House
London Road
Hillbrow
Liss
Hants
GU33 7NT

Telephone
+44 (0) 1730 895900
Facsimile
+44 (0) 1730 895922

EC DECLARATION OF CONFORMITY


Morgan Automation Ltd. declares that the apparatus known as :-

Neurotherm – Model JK3

is in conformity with the following Standards and Requirements :-

BS EN 60601-1-2 : 1993
ISO 13485 : 1996
BS EN 46001 : 1994
BS 5724 Section 2:3 : 1983
BS 5724 Section 1:1 :1992 IEC 601-1-1:1992
Directive 93/42/EEC

Is subject to the procedure set out in Annex II of Directive 93/42/EEC under the supervision of Notified Body Number 0120, SGS Yarsley International Certification Services Ltd.

Signed :-  Date : 3-12-98

Name ...H.M. CLARKE..... (Technical Director)

Ref: TECHFILE/NTDEC.doc

Morgan Automation Ltd. Reg. Office: Rake Heath House, London Road, Hillbrow, Liss, Hants GU33 7NT Reg. No. 2174066



NEUROTHERM RADIO FREQUENCY LESION GENERATOR
MODEL JK3
SERVICE MANUAL

12-2



429 Brighton Road, South Croydon, Surrey, CR2 6EU, UK Tel: 0181-660 4374 Fax: 0181-660 9417

YEAR 2000 CONFORMITY

16th December 1997

General Statement

This is to confirm that the NeuroTherm RF Lesion Generator does not contain a date sensitive microprocessor or date sensitive programme.

ROD GEMMELL
Managing Director



International Tel: +(44) 181 660 4374

International Fax: +(44) 181 660 9417

E-mail: admin@rdgmedical.com

www.rdgmedical.com

A division of Croydon Industries Limited

Registered No: 1390305 (England)