



Model TED 200-T7

Oxygen Monitor

Service Manual



-  **TYPE B EQUIPMENT:** Equipment providing a particular degree of protection against electric shock, particularly regarding—
- allowable LEAKAGE CURRENT
 - Reliability of the protective earth connection (if present).

Teledyne Analytical Instruments
16830 Chestnut Street
City of Industry, CA 91749

P/N SM70027
04/09/99
ECO#: 99-0159

Table of Contents

| | |
|-------------------------------------|----|
| Description of the TED 200T7 | 3 |
| Sensor | 3 |
| Alarms | 3 |
| Display and Touch Panel | 3 |
| | |
| Circuit Description | 4 |
| Power Supply | 4 |
| | |
| Repair/Replacement | 4 |
| Battery Replacement | 4 |
| T-7 Sensor Replacement | 4 |
| Opening the Case | 4 |
| Sensor Simulator | 5 |
| | |
| Troubleshooting | 6 |
| | |
| Manufacturing Test Procedures | 8 |
| | |
| Specifications | 11 |
| | |
| Drawing List | 12 |
| | |
| Parts List | 13 |

Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

This manual is intended to provide information to help qualified maintenance personnel service and repair the TED Model 200T7 oxygen monitor. The ability to read and follow schematics is assumed, as is basic electronic knowledge of the functions and characteristics of operational amplifiers.

The only equipment needed to service the monitor is a digital voltmeter; however, you may find that a means of simulating the sensor is useful in troubleshooting. A schematic for a simple cell simulator circuit is provided on page 6.

Description of the TED 200T7

The TED 200T7 uses a galvanic fuel cell sensor to detect the oxygen level in the gas immediately surrounding the tip of the sensor. The sensor signal is processed and displayed as a percentage.

Sensor

The T-7 Micro-Fuel Cell puts out a minute current in proportion to the partial pressure of the oxygen in the gas surrounding the sensor. The current is zero at zero oxygen pressure and is a maximum at 100% oxygen.

Alarms

The two alarms are set through the touch panel. The HI alarm can be set anywhere from 21 to 100, or OFF. The LO alarm can be set from 17 to 99. There is an alarm test option as well as an alarm silence key, which defeats the audible/visible alarm in increments of 30, 90, and 180 seconds.

Display and Touch Panel

The display is a digital liquid crystal display (LCD) showing oxygen level in percent.

All functions are microprocessor-controlled. Commands are input through the touch panel and processed in a ROM.

Circuit Description

The TED 200T7 is a high precision digital and analog circuit implemented in fine geometry, surface mount technology. Due to the highly integrated nature of this circuitry, Teledyne does not recommend that the circuitry or the board be troubleshot in the field.

Power Supply

The TED 200T7 uses 4 "AA" alkaline batteries which provide 6 volts and a ground. An integrated circuit chip regulates this voltage to 4.61 volts.

Repair/Replacement

Battery Replacement

1. Turn the instrument off.
2. Slide battery holder out.
3. Remove the battery holder and take out the old batteries (if any).
4. Install 4 "AA" alkaline batteries (other types of batteries will give erroneous battery test readings) into the holder, observing proper polarity.
5. Snap the battery holder back onto the leads and slide the holder back into the compartment.

T-7 Sensor Replacement

1. Turn off the instrument.
2. Unplug the coiled cable from the old sensor.

3. Remove the new sensor from its protective plastic bag.
4. Plug the coiled cable into the new sensor and calibrate.

Opening the Case

1. Turn the instrument off
2. On the bottom of the unit are four rubber bumper feet. Remove the screw in the center of each foot.
3. Carefully lift off the cover. If desired, turn the unit over, holding the two case halves together, and remove the top cover.

To close the instrument, reverse this procedure, being careful not to pinch any wires between the halves of the case.

Sensor Simulator

Occasionally, it is necessary to simulate a sensor for testing. The sensor connector and simulator are wired as shown below.

The sensor simulator shown will only give a general indication that the TED 200T7 electronic circuitry is operating normally. Actual calibrations should be in air with the sensor attached.





Figure 1a & 1b: the sensor simulator circuit.

1.0 INTRODUCTION

This functional test is performed as the final check of the TED 200 T7 Portable Oxygen Monitor prior to shipment.

2.0 APPLICABLE DOCUMENTS

- 2.1 Schematic D-68003 (1 of 2)
- 2.2 Schematic D-68003 (2 of 2)
- 2.4 Instruction Manual.

3.0 TEST OBJECTIVE

The objective of this test is to exercise all the functions of the monitor to ensure it performs to the specifications as outlined in the above referenced documents.

4.0 TEST EQUIPMENT

4.1 Power Requirements - As stated per the running sheet, otherwise, power shall be 6.3 Vdc.

4.2 Equipment Requirements - As follows:

- 5.2.1 AA Battery, 1.5 Vdc (4 ea.) or adjustable power supply.
- 5.2.2 Keithley Current Source
- 5.2.3 "Dummy" (sensor) cell.

5.0 TEST SETUP

Circuit oscillations may be generated when using the Keithley. If this happens, place a 10k or 20k ohm resistor in series with the input connector.




Figure 1. Test Setup
TED200-T7

TED 200T7 Service Manual

| STEP # | OPERATOR ACTION | EXP. RESULT | TOLERANCE |
|------------------|--|---|--------------|
| 1 | Install four AA 1.5 Vdc batteries or adjustable DC power supply. | N/A | |
| 2 | Press the ON/OFF key once. | Unit LCD displays "High Alarm 100" "Low Alarm 17" and flashing message "AIR CAL" | Exact |
| 3 | Install the dummy sensor cell, then connect the Keithley instrument to the dummy cell. | N/A | |
| 4 | Adjust Keithley to read zero μ a. | N/A | |
| ZERO ADJUSTMENT: | | | |
| 5 | Connect DMM BLACK lead to TP2 and RED lead to TP3. | N/A | |
| 6 | Move a Shorting block JP4-1 and JP4-2 to JP4-2 and JP4-3 | N/A | |
| 7 | Adjust R7 located on the analog input board so that the DMM read -5.5mv. | -5.5mv | ± 0.5 mv |
| 8 | Return shorting block JP4-2 and JP4-3 to JP4-1 and JP4-2. | N/A | |
| 9 | Adjust the Keithley to 112 μ a. | N/A | |
| 10 | Remove DMM. | N/A | |
| 11 | Press the CAL key once | LCD display counts down from 10 to 0, then displays 21. | ± 1 |
| | | LCD displays message "CAL IN 100%". | Exact |
| 12 | Adjust the Keithley for 535.885 μ a. | LCD displays 100 | $\pm 5\%$ |
| 13 | Press the CAL key twice. | LCD displays message "CAL IN 100%" and counts down from 15 to 0, then displays 100. | $\pm 1\%$ |
| | | Alarm buzzer sounds and red LED flashes. | Exact |



| STEP # | OPERATOR ACTION | EXP. RESULT | TOLERANCE |
|--------|--|---|-----------|
| 14 | Adjust the Keithley output to 267.9 μ a. | LCD displays 50 | Exact |
| 15 | Press the SET HI ALARM key, then press the UP arrow key. | LCD displays "OFF for HI ALARM setting and flashing message "ALARM DEFEATED". | Exact |
| 16 | Using the DOWN arrow key, adjust the HI ALARM setting to a value less than 100, e.g., "88". | LCD displays value set. | Exact |
| 17 | Press the SET LO ALARM key, then press the UP arrow key until the LO ALARM setting reads "21". | LCD displays value set. | Exact |
| 18 | Adjust the Keithley to 112 micro amps. | LCD displays 21. Alarm buzzer sounds. | Exact |
| 19 | Press the ALARM SILENCE key once. | Red LED flashes. After 5 seconds, the LCD display counts down from "30" in the lower right hand corner. Alarm buzzer stops. Red LED continues flashing. | Exact |
| 20 | Press the SET LO ALARM key | N/A | |
| 21 | Using the DOWN arrow key, set the LO ALARM to "19". | Red LED goes off. | Exact |
| 22 | Press the ALARM SILENCE key | "30" is displayed in lower right hand corner of LCD display. | Exact |
| 23 | Press the ALARM SILENCE key a second time. | "90" is displayed in lower right hand corner of LCD display. | Exact |
| 24 | Press the ALARM SILENCE key a third time. | "180" is displayed in lower right hand corner of LCD display. | Exact |

TED 200T7 Service Manual

| STEP # | OPERATOR ACTION | EXP. RESULT | TOLERANCE |
|--------|--|--|-----------|
| 25 | Press the ALARM SILENCE key a fourth time. | "00" is displayed in lower right hand corner of LCD display. | Exact |
| 26 | Press the ALARM TEST key twice. | HI ALARM setting counts up from 19 to the value set in Step 9, e.g. "88". | Exact |
| | | Alarm buzzer sounds twice. | Exact |
| | | Red LED flashes twice. Then the setting counts back down to the LO ALARM setting of 19. | Exact |
| | | The buzzer sounds twice again and the red LED flashes twice. The LCD then displays 21. | ±1 |
| 27 | Press the BATT TEST key. | LCD displays "BATTERY" HOURS LEFT 999". | Exact |
| 28 | Press the ON/OFF key twice. | The unit shuts off. | Exact |
| 29 | Press the ON/OFF key once. | LCD displays 21. | ±1 |
| 30 | Adjust Keithley output to 0.000 µa. Allow meter to stabilize. | LCD displays 0 | ±1% |
| 31 | Adjust Keithley output to 535.885 µa. | LCD displays 100. | Exact |
| 32 | Verify the unit is in normal operating mode (e.g., displaying an O ₂ concentration). Then move the shorting block from E1-1/E1-2 to E11-22/E11-3. | After 5 to 15 seconds, the unit's LCD goes blank the alarm buzzer sounds continuously, & the red LED stays on. | Exact |
| 33 | Return the shorting block from E1-2/E1-3 to E1-1/E1-2. | N/A | |
| 34 | Turn the unit under test POWER switch OFF. Disconnect all test leads and test equipment. | N/A | |

Specifications

| | |
|--|--|
| Range: | 0-100% oxygen |
| Accuracy: | (after 2-point calibration) $\pm 2\%$ of full scale for 8 hours at constant temperature |
| Response Time: | 90% in less than 8 seconds (typically 6-8 sec.) at 25°C. |
| Display Resolution: | Digital LCD with battery life indication. |
| Calibration Time: | @ 21%: 10 sec. @ 100%: 15 sec. |
| Battery Life: | Approximately 1,000 hrs. continuous use in non-alarm conditions using 4 "AA" alkaline batteries. |
| Humidity Range: | 0-95% RH |
| Expected Sensor Life: | 1 year in most applications |
| Sensor Type: | Class T-7 (Galvanic) |
| Dimensions: | 7"W x 5"D x 2-1/2"H (178 mm x 127 mm x 63½ mm) |
| Weight: | Less than 2 lbs. |
| Cable Length: | Retracted: 2 ft. Extended: 10 ft. |
| Storage Temp.: | 0-50°C (Recommended Temp. 10-30°C) |
| Operating Temp.: | 0-40°C |
| Alarm Indicators: | Audible/Visible Audio: Pulsating .5-1 KHZ tone @ 70 to 90 DBA, 50% duty cycle. Visual: Red high-brightness LED. |
| Stability: | Less than 1% O ₂ over an 8 hour period at constant pressure and temperature |
| O₂ level alarm ranges: | Hi alarm: 18 to 100 Percent Low alarm: 17 to 100 percent |
| O₂ level alarm accuracy: | Alarm thresholds digitally stored and alarm effectivity is digitally compared against the value calculated by the micro processor. Hence, the accuracy of the alarms is the same as the accuracy of the unit as specified above. |
| Warm up time: | Typically 20 minutes or less. |

RFI Protected

Drawing List

D-68003 Schematic Diagram, Surface Mount PCB
D-68004 PCB Assembly, Board Set
D-68987A TED 200T7 Final Assembly

Optional Accessories

1 B34102 Mounting Clamp
1 A51589 T-7 adapter cap, female (22 mm)
1 A51588 T-7 adapter cap, male (22 mm)
1 C53790 Calibration assembly
1 A284 Universal adapter set for pediatric circuits (15 mm)
1 A274 Tee adapter, autoclavable
1 A283 Tee adapter, metal

A minimum charge is applicable to spare parts orders.

IMPORTANT: Orders for replacement parts should include the model number, serial number, and range of the analyzer for which the parts are intended.

Orders should be sent to:

TELEDYNE ELECTRONIC TECHNOLOGIES

Analytical Instruments

16830 Chestnut Street

City of Industry, CA 91749-1580

Phone (626) 934-1500

FAX (626) 961-2538

TWX (910) 584-1887 TDYANLY COID

or your local representative


Parts List for TED 200T7

| Qty. | P/N | Item Description |
|------|--------|---|
| 1 | A51329 | Micro-Fuel Cell T-7 w/Flow Diverter P/N A50057 |
| 4 | B99 | "AA" Size Alkaline Battery |
| 1 | B551 | Battery Holder |
| 1 | B69934 | Cable Assembly |
| 1 | A268 | Tee Adapter (22mm) |

NOTE: UNLESS OTHERWISE SPECIFIED

1. PREPARED FOR DOCUMENTS LISTED
2. FOR ASSEMBLY SEE DRAWING
3. RESISTANCE VALUES ARE IN OHMS, 1/100, 1/10, 1/2, 1, 10, 100, 1000
4. CAPACITANCE VALUES ARE IN MICROFARADS, 100, 10, 1, 0.1, 0.01, 0.001

| | | | |
|-----|----------------------------|----------|----|
| REV | DESCRIPTION | DATE | BY |
| 1 | INITIAL DESIGN FOR DRAWING | 10/10/68 | AM |
| 2 | FOR APPROVAL AND RELEASE | 11/24/71 | AA |



| | | | |
|--|--|--|----------|
| PREPARED FOR DOCUMENTS LISTED 1. PREPARED FOR DOCUMENTS LISTED 2. FOR ASSEMBLY SEE DRAWING 3. RESISTANCE VALUES ARE IN OHMS, 1/100, 1/10, 1/2, 1, 10, 100, 1000 4. CAPACITANCE VALUES ARE IN MICROFARADS, 100, 10, 1, 0.1, 0.01, 0.001 | | TELETYPE ANALYTICAL INSTRUMENTS A DIVISION OF THE ELECTRONIC TECHNOLOGY | |
| SCHEMATIC DIAGRAM TEU 20017 SMT PCB | | DATE 11/24/71 | BY AA |
| D-68003 | | 1 | 1 |

0 7 6 5 4 3 2 1

NOTE: UNLESS OTHERWISE SPECIFIED.


1. INTERPRET DIMENSIONS AND TOLERANCES IN ACCORDANCE WITH MIL-STD-208.
2. FOR SCHEMATIC SYMBOLS SEE DRAWING.
3. SOLDER JET CONTACT MEASUREMENTS.

▲ MARK POSITION AND SERIAL NUMBER PERMANENT AND LEGIBLE IN AREA INDICATED.

▲ MARK CONFIGURATION "A", "B" OR "C" PERMANENT AND LEGIBLE.
 A - ENGLISH VERSION, B - FRENCH VERSION OR C - GERMAN VERSION IN AREA INDICATED.

| REV. | DESCRIPTION | DATE | BY | CHK. |
|------|------------------------|------|----|------|
| 1 | ISSUED FOR DESIGN | | | |
| 2 | ISSUED FOR MANUFACTURE | | | |
| 3 | ISSUED FOR ASSEMBLY | | | |

- STEPS FOR INSTALLING JS ASSEMBLY**
1. SOLDER JS P/W JUMPS FIRST ONTO PCB WITH NO GAPS.
 2. SOLDER WIRES FROM INSIDE OF BOARD BRIDGING SW P/W AIRWAYS TO PCB WITHIN SILICONE-RESIN OUTLINE. USE IDENTIFIED COLOR SILICONE FOR P/W. THE WIRE IS P/W MARKED "YEL" AND IS ON.
 3. SOLDER BRUSH BRIDGING ONE PIN ONLY TO PCB WITH NO GAPS BETWEEN SW AND PCB.
 4. SOLDER COLORED WIRES ON OUTSIDE OF BOARD BRIDGING SW TO CORRESPONDING P/W WITHIN PCB WITH WIRE SOLDERED TO P/W MARKED "RED" AND IS ON.



6. IF PRESENT, REMOVE LED PROTECTIVE FILM BEFORE POSITIONING LED ON BOARD. SAVE THE FILM TO REUSE AFTER THE LED HAS BEEN INSTALLED ON THE PCB.
7. INSTALL BRIDGE FOR LED ON SOLDER SIDE UNDER PCB.
8. INITIAL LED ON BOARD. WITH ORIENTATION AS SHOWN ON SOLDER SIDE OF PCB. AFTER ALL SOLDERING IS COMPLETE.
9. REPLACE THE PROTECTIVE FILM OVER LED FACE WHEN THE LED IS INSTALLED ON THE PCB.
10. TRIM ALL COMPONENT LEADS TO .00 INCH AFTER SOLDERING.
11. INSTALL BRUSH JUMPER ON S1 PINS 1 AND 2.
12. INITIAL BRUSH JUMPER ACROSS PINS 1 AND 2 OF JS FOR NORMAL OPERATION. JUMPER INSTALLATION ON JS PINS 3 AND 4 IS RESERVED FOR TEST USE ONLY.

▲ SOLDERING INSTRUCTIONS FOR SW P/W LEADS. LEAD TEMPERATURE:
 TYPICAL: 300 DEGREES C
 MAXIMUM: 325 DEGREES C



WARNING: ANTI-STATIC HANDLING PROCEDURES ARE REQUIRED IN ASSEMBLY OF THIS PCB. REF. A29973

| PCB PART NO. | | REV. | | DATE | |
|--------------|--------|------|---|------|---|
| 1 | 100004 | 1 | 1 | 1 | 1 |
| 2 | | 2 | 1 | 1 | 1 |
| 3 | | 3 | 1 | 1 | 1 |
| 4 | | 4 | 1 | 1 | 1 |
| 5 | | 5 | 1 | 1 | 1 |
| 6 | | 6 | 1 | 1 | 1 |
| 7 | | 7 | 1 | 1 | 1 |
| 8 | | 8 | 1 | 1 | 1 |
| 9 | | 9 | 1 | 1 | 1 |
| 10 | | 10 | 1 | 1 | 1 |
| 11 | | 11 | 1 | 1 | 1 |
| 12 | | 12 | 1 | 1 | 1 |
| 13 | | 13 | 1 | 1 | 1 |
| 14 | | 14 | 1 | 1 | 1 |
| 15 | | 15 | 1 | 1 | 1 |
| 16 | | 16 | 1 | 1 | 1 |
| 17 | | 17 | 1 | 1 | 1 |
| 18 | | 18 | 1 | 1 | 1 |
| 19 | | 19 | 1 | 1 | 1 |
| 20 | | 20 | 1 | 1 | 1 |
| 21 | | 21 | 1 | 1 | 1 |
| 22 | | 22 | 1 | 1 | 1 |
| 23 | | 23 | 1 | 1 | 1 |
| 24 | | 24 | 1 | 1 | 1 |
| 25 | | 25 | 1 | 1 | 1 |
| 26 | | 26 | 1 | 1 | 1 |
| 27 | | 27 | 1 | 1 | 1 |
| 28 | | 28 | 1 | 1 | 1 |
| 29 | | 29 | 1 | 1 | 1 |
| 30 | | 30 | 1 | 1 | 1 |
| 31 | | 31 | 1 | 1 | 1 |
| 32 | | 32 | 1 | 1 | 1 |
| 33 | | 33 | 1 | 1 | 1 |
| 34 | | 34 | 1 | 1 | 1 |
| 35 | | 35 | 1 | 1 | 1 |
| 36 | | 36 | 1 | 1 | 1 |
| 37 | | 37 | 1 | 1 | 1 |
| 38 | | 38 | 1 | 1 | 1 |
| 39 | | 39 | 1 | 1 | 1 |
| 40 | | 40 | 1 | 1 | 1 |
| 41 | | 41 | 1 | 1 | 1 |
| 42 | | 42 | 1 | 1 | 1 |
| 43 | | 43 | 1 | 1 | 1 |
| 44 | | 44 | 1 | 1 | 1 |
| 45 | | 45 | 1 | 1 | 1 |
| 46 | | 46 | 1 | 1 | 1 |
| 47 | | 47 | 1 | 1 | 1 |
| 48 | | 48 | 1 | 1 | 1 |
| 49 | | 49 | 1 | 1 | 1 |
| 50 | | 50 | 1 | 1 | 1 |
| 51 | | 51 | 1 | 1 | 1 |
| 52 | | 52 | 1 | 1 | 1 |
| 53 | | 53 | 1 | 1 | 1 |
| 54 | | 54 | 1 | 1 | 1 |
| 55 | | 55 | 1 | 1 | 1 |
| 56 | | 56 | 1 | 1 | 1 |
| 57 | | 57 | 1 | 1 | 1 |
| 58 | | 58 | 1 | 1 | 1 |
| 59 | | 59 | 1 | 1 | 1 |
| 60 | | 60 | 1 | 1 | 1 |
| 61 | | 61 | 1 | 1 | 1 |
| 62 | | 62 | 1 | 1 | 1 |
| 63 | | 63 | 1 | 1 | 1 |
| 64 | | 64 | 1 | 1 | 1 |
| 65 | | 65 | 1 | 1 | 1 |
| 66 | | 66 | 1 | 1 | 1 |
| 67 | | 67 | 1 | 1 | 1 |
| 68 | | 68 | 1 | 1 | 1 |
| 69 | | 69 | 1 | 1 | 1 |
| 70 | | 70 | 1 | 1 | 1 |
| 71 | | 71 | 1 | 1 | 1 |
| 72 | | 72 | 1 | 1 | 1 |
| 73 | | 73 | 1 | 1 | 1 |
| 74 | | 74 | 1 | 1 | 1 |
| 75 | | 75 | 1 | 1 | 1 |
| 76 | | 76 | 1 | 1 | 1 |
| 77 | | 77 | 1 | 1 | 1 |
| 78 | | 78 | 1 | 1 | 1 |
| 79 | | 79 | 1 | 1 | 1 |
| 80 | | 80 | 1 | 1 | 1 |
| 81 | | 81 | 1 | 1 | 1 |
| 82 | | 82 | 1 | 1 | 1 |
| 83 | | 83 | 1 | 1 | 1 |
| 84 | | 84 | 1 | 1 | 1 |
| 85 | | 85 | 1 | 1 | 1 |
| 86 | | 86 | 1 | 1 | 1 |
| 87 | | 87 | 1 | 1 | 1 |
| 88 | | 88 | 1 | 1 | 1 |
| 89 | | 89 | 1 | 1 | 1 |
| 90 | | 90 | 1 | 1 | 1 |
| 91 | | 91 | 1 | 1 | 1 |
| 92 | | 92 | 1 | 1 | 1 |
| 93 | | 93 | 1 | 1 | 1 |
| 94 | | 94 | 1 | 1 | 1 |
| 95 | | 95 | 1 | 1 | 1 |
| 96 | | 96 | 1 | 1 | 1 |
| 97 | | 97 | 1 | 1 | 1 |
| 98 | | 98 | 1 | 1 | 1 |
| 99 | | 99 | 1 | 1 | 1 |
| 100 | | 100 | 1 | 1 | 1 |



MODEL / SERIES: TED 200T7

A- ENGLISH VERSION

DRAWING TITLE: TED 200T7 BOARD SET


DATABASE
FILENAME
D68004.XLS

DATE: 04/20/99
A-A- 4/21/99
ECO NO.: 99-0167

| ITEM NO | QUANTITY | | | | | | | | | YET PART NUMBER | REF./PCB LOCATION | ITEM DESCRIPTION | REMARKS |
|---------|----------|---|---|---|---|---|---|---|---|-----------------|----------------------------|------------------------------|-------------------------------------|
| | A | B | C | D | E | F | G | H | I | | | | |
| 1 | 1 | | | | | | | | | R2191 | R17 | RES., 1.5M Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 2 | 1 | | | | | | | | | R2195 | R27 | RES., 100 Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 3 | 9 | | | | | | | | | R2129 | R1,15,24,26,36,42,43,46,49 | RES., 100K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 4 | 1 | | | | | | | | | R2179 | R32 | RES., 10K Ohms, 1/8W, 0.1% | (SMT 1206 PKG) |
| 5 | 4 | | | | | | | | | R2272 | R5, 6, 12, 48 | RES., 10M Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 6 | 2 | | | | | | | | | R2273 | R19, 30 | RES., 15K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 7 | 1 | | | | | | | | | R2274 | R33 | RES., 195K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 8 | 6 | | | | | | | | | R2130 | R4, 14, 21, 47, 50 | RES., 1M Ohms, 1/8W, 5% | (SMT 1206 PKG) |
| 9 | 1 | | | | | | | | | R2275 | R8 | RES., 2.26M Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 10 | 2 | | | | | | | | | R2235 | R22, 51 | RES., 2.74K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 11 | 1 | | | | | | | | | R2276 | R35 | RES., 200K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 12 | 1 | | | | | | | | | R2277 | R45 | RES., 210K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 13 | 1 | | | | | | | | | R2278 | R34 | RES., 3.32M Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 14 | 1 | | | | | | | | | R2279 | R38 | RES., 301K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 15 | 1 | | | | | | | | | R2280 | R25 | RES., 3.01M Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 16 | 1 | | | | | | | | | R2281 | R13 | RES., 402K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 17 | 1 | | | | | | | | | R2282 | R2 | RES., 437K Ohms, 1/8W, 0.25% | (SMT 1206 PKG) |
| 18 | 1 | | | | | | | | | R2283 | R31 | RES., 47.5K Ohms, 1/8W, 1% | (SMT 1206 PKG) PFC-W1206R-02-4373-B |
| 19 | 2 | | | | | | | | | R2284 | R16, 29 | RES., 493K Ohms, 1/8W, 0.1% | (SMT 1206 PKG) |
| 20 | 2 | | | | | | | | | R2285 | R16, 28 | RES., 523K Ohms, 1/8W, 1% | (SMT 1206 PKG) PFC-W1206R-02-4933-B |
| 21 | 1 | | | | | | | | | R2286 | R44 | RES., 542K Ohms, 1/8W, 0.1% | (SMT 1206 PKG) |
| 22 | 1 | | | | | | | | | R2287 | R11 | RES., 576K Ohms, 1/8W, 1% | (SMT 1206 PKG) PFC-W1206R-02-5423-B |
| 23 | 1 | | | | | | | | | R2288 | R39 | RES., 715 Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 24 | 4 | | | | | | | | | R2182 | R20, 23, 40, 41 | RES., 8.25K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 25 | 1 | | | | | | | | | R2126 | R3 | RES., 82.5K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 26 | 1 | | | | | | | | | R2289 | R37 | RES., 9.53K Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 27 | 0 | | | | | | | | | R2290 | R9, R10 | RES., XXXX Ohms, 1/8W, 1% | (SMT 1206 PKG) |
| 28 | 1 | | | | | | | | | P366 | R7 | POT., 50K, WW | (THRU HOLE) |
| 29 | 2 | | | | | | | | | T1232 | Q1, Q4 | XSTR., MMB13640LT1 | (SMT SOT-23) |
| 30 | 4 | | | | | | | | | T1191 | Q2, 3, 5, 6 | XSTR., MMBT2222ALT1 | (SMT SOT-23) |
| 31 | 7 | | | | | | | | | D462 | D1-D7 | DIODE, FULL4148.TR | (SMT MLL-34) |
| 32 | 1 | | | | | | | | | D475 | D9 | DIODE, MLL4001 | (SMT MLL-34) |
| 33 | 1 | | | | | | | | | Z25 | D8 | DIODE, ZENER, 114968, 130V | (THRU HOLE D07) |
| 34 | 1 | | | | | | | | | CP1916 | C19 | CAP., 047 uF, 50V, CERAMIC | (SMT 0605) |
| 35 | 10 | | | | | | | | | CP1753 | C6-10, 14, 16, 17, 25, 26 | CAP., 0.1 uF, 50V, CERAMIC | (SMT 0605) |




NOTES: UNLESS OTHERWISE SPECIFIED.

- 1 STICK A 4.0 X .75 X .13 THICK FOAM RUBBER PAD ON TOP OF COVER UNDER SURFACE AS SHOWN.
- 2 BREAK ITEM 3 INTO TWO PIECES, THEN TRIM BOTH EDGES, ON PERFORATED TABS. INSTALL ITEM 3 IN PLACE AS SHOWN.
- 3 LOCATE APPROX. AS SHOWN
- 4 OVERALL WEIGHT: WITH BATTERY 1 Pound 6 Ounces
WITH-OUT BATTERY 1 Pound 4 Ounces
- 5 KIT WITH UNIT
- 6 APPLY E-100 TO ITEM 41 AS SHOWN.



| LEGEND | |
|--------|--------------------|
| SYMBOL | DESCRIPTION |
| A | 1TED 200 T7 |
| B | MDL 850 (SECURIST) |
| C | OM-100 (NEWPORT) |

| REV | | REVISIONS | | | |
|-----|-----------------|-----------|------|------|--|
| REV | DESCRIPTION | DATE | APP | REV | |
| 1 | INC ECO 88-0207 | 11/17/88 | R.H. | R.H. | |
| 2 | INC ECO 88-0180 | 4/13/88 | R.J. | V.F. | |



| ITEM | QTY | PART NO. | DESCRIPTION |
|------|-----|-----------|--|
| 42 | REF | D68987-RS | ROUTE SHEET |
| 41 | 4 | ---- | WASHER, FIBER #4 (.031 THK) |
| 40 | 1 | A51327 | T7 SENSOR |
| 38 | 1 | B-68984 | WARNING INFORMATION SHEET |
| 38 | 1 | A-70130 | |
| 37 | 1 | A-57342 | BOX LABEL |
| 36 | 1 | A-57483 | |
| 35 | REF | D-88003 | SCHEMATIC DIAGRAM |
| 34 | 1 | 1/68987C | OPERATING MANUAL |
| 33 | 1 | 1/68987A | OPERATING MANUAL |
| 33 | REF | 1/68987A | OPERATING MANUAL |
| 33 | REF | CP1922 | ROBOL CABLE 8-WIRES |
| 32 | 4 | CP172 | COMP. TERMINAL |
| 31 | 1 | CP1725 | MSX-IT CONNECTOR |
| 30 | 1 | A 768 | T-ADAPTER, (NOT SHOWN, INCLUDED IN SHIPPING BOX) |
| 28 | 1 | B 99 | 7AA BATTERY, (NOT SHOWN, INCLUDED IN SHIPPING BOX) |
| 28 | 1 | A-68403 | LABEL, SERIAL NUMBER (TED 200 T7) |
| 27 | 1 | A-57484 | LABEL, SERIAL NUMBER (MDL 850) |
| 27 | 1 | A-89013 | LABEL, SERIAL NUMBER (OM-100) |
| 26 | 1 | A-41151 | LABEL, FCC |
| 25 | 1 | BB192 | BOX, WITH FOAM |
| 24 | 1 | A-70089 | LABEL, CSA |
| 23 | 1 | 065899-CV | LABEL, MODEL AND CLASS |
| 22 | 1 | A-45582 | LABEL, WARNING |
| 21 | 1 | ---- | FOAM RUBBER 4" X 3/4" X 1/8" |
| 20 | 4 | SS 272 | SCREW, 4-40 X 1-3/8" LG BINDER HEAD |
| 19 | 4 | B 330 | BUMPER |
| 18 | 1 | D-89934 | COILED CABLE, ASSEMBLY |
| 17 | 2 | ---- | KEY'S NUT, #8 |
| 16 | 2 | ---- | SCREW, 8-32 X 3/8" |
| 15 | 1 | CP 1253 | PROBE CLIP |
| 14 | 2 | ---- | SCREW, 7-56 X 1/4" LG BINDER HEAD |
| 13 | 2 | ---- | LOCKWASHER #2 EXTERNAL TOOTH |
| 12 | 1 | A 178 | ALARM |
| 11 | 2 | S 823 | SPACER, 3/16" O.D. 2-56 X 8/16" LG |
| 10 | 2 | ---- | SCREW, 8-32 X 3/8" FH UNDERCUT |
| 9 | 1 | B 551 | 4 CELL "AA" BATTERY HOLDER WITH ORATOR |
| 8 | 2 | LL 48 | LOCKWASHER #4 EXTERNAL TOOTH |
| 7 | 2 | SS 78 | SCREW, BINDER HEAD, 4-40 X 1/4" |
| 6 | 1 | H 240 | HANDLE, VEMALINE PRODUCTS /SABA-B |
| 5 | 1 | A-87452 | |
| 5 | 1 | A-88822 | |
| 5 | 1 | A-88823 | BACK PANEL SILKSCREEN DETAIL |
| 4 | 1 | A-88885 | LABEL, L100 |
| 3 | 1 | D-88004 | P.C. ASSEMBLY, BOARD SET, TED 200 T7 |
| 2 | 1 | B-88980A | FRONT PANEL, ASSEMBLY |
| 1 | 1 | B-55014 | CASE (PER C-35015 SHEET 1 AND 2) |
| 1 | 1 | B-57568 | CASE (PER C-35015 SHEET 1 AND 2) |
| 1 | 1 | B-55014 | CASE (PER C-35015 SHEET 1 AND 2) |

| SIGNATURES | | DATE | TITLE | SCALE |
|------------|------|----------|-------|-------|
| BY | DATE | 11-11-88 | | 1/16" |
| CHKD | | | | |
| APP'D | | | | |
| ENGR | | | | |
| CB | | | | |

| | | | | |
|---|--|--|---|----------------|
| DR SCALE ONE UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES DECIMALS ARE TO TWO PLACES FRACTIONS ARE TO SIXTEENTHS UNLESS OTHERWISE SPECIFIED | ALL DIMENSIONS ARE TO UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE TO UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE TO UNLESS OTHERWISE SPECIFIED | ELECTRONIC TECHNOLOGIES ANALOGIC INSTRUMENTS 10000 WOODBRIDGE BLVD NEWPORT NEWS, VA 23606 | TITLE TED 200 T7 FINAL ASSEMBLY | SCALE 1/16" |
|---|--|--|---|----------------|

AS SHOWN D-68987 1-2