


**M**ALLINCKRODT


**PURITAN-BENNETT**

# Infant Star<sup>®</sup> 500/950

*Ventilator  
Systems*



*Service and  
Repair  
Instructions*



---

Copyright information

Copyright 1999 Mallinckrodt™ Inc. All rights reserved. The Infant Star® 500/950 ventilator systems are manufactured in accordance with Mallinckrodt proprietary information, covered by one or more of the following U.S. Patents and foreign equivalents: 4,727,871; 4,747,403; 4,838,259; 5,237,987; 5,452,714; 5,385,142; 5,542,415. Infant Star is a registered trademark of Mallinckrodt Inc.

The information contained in this manual is the sole property of Mallinckrodt Inc. and may not be duplicated without permission. This manual may be revised or replaced by Mallinckrodt Inc. at any time and without notice. You should ensure that you have the most current applicable version of this manual; if in doubt, contact the Technical Publications Department of Mallinckrodt Inc. at the Carlsbad, CA location. While the information set forth herein is believed to be accurate, it is not a substitute for the exercise of professional judgment.

The ventilator should be operated and serviced only by trained professionals. Mallinckrodt's sole responsibility with respect to the ventilator, and its use, is as stated in the limited warranty provided.

Nothing in this manual shall limit or restrict in any way Mallinckrodt's right to revise or otherwise change or modify the equipment (including its software) described herein, without notice. In the absence of an express, written agreement to the contrary, Mallinckrodt has no obligation to furnish any such revisions, changes, or modifications to the owner or user of the equipment (including its software) described herein.

*Manufactured by:*

**Nellcor Puritan Bennett Inc.**  
A subsidiary of Mallinckrodt Inc.  
2200 Faraday Avenue  
Carlsbad, CA 92008-7208 USA

Phone: 760.603.5300  
Fax: 760.603.5902

*Corporate Headquarters:*

**Mallinckrodt Inc.**  
675 McDonnell Boulevard  
P.O. Box 5840  
St. Louis, MO 63134 USA

Phone: 314.652.2000

*Authorized Representative:*

**Mallinckrodt Europe BV**  
Hambakenwetering 1  
5231 DD 's-Hertogenbosch  
The Netherlands

Phone: +31.73.6485200  
Fax: +31.73.6410915

---

## Authorized Service

This manual is intended for use by Mallinckrodt trained and authorized Service Personnel. Any adjustments or procedures that exceed the scope of this manual should be referred to Mallinckrodt.

Mallinckrodt does not condone or approve of service activity on its products by anyone other than Mallinckrodt trained and authorized personnel. Mallinckrodt is not responsible for any repairs made by unauthorized procedures or personnel.

Use of the incorrect part or failure to exercise due care in the installation, removal, servicing, or calibration of parts and equipment may result in damage to the equipment and possible malfunction of the equipment, which may in turn result in damage to property and injury (including death) to persons. The purchaser and installer of these parts shall bear full responsibility and liability for the above.

For specific operating instructions and clinical theory of operation, refer to the *Infant Star 500/950 Ventilator's Operating Instructions Manual*. Service personnel should become familiar with the Operating and Repair techniques before servicing this equipment.

## Electromagnetic Susceptibility

The *Infant Star 500/950 Ventilator Systems* comply with the requirements of IEC 601-1-2 (EMC Collateral Standard), including the E-field susceptibility requirements at a level of 10 volts per meter, at frequencies from 26 MHz to 1 GHz, and the ESD requirements of this standard up to 6 kV direct and 15 kV air discharge. However, even at this level of device immunity, certain transmitting devices (cellular phones, walkie-talkies, cordless phones, paging transmitters, etc.) emit radio frequencies that could interrupt ventilator operation if their proximity is located too close to the ventilator. It is difficult to determine when the field strength of these devices becomes excessive. Practitioners should be aware that radio frequency emissions are additive, and that the ventilator must be located a sufficient distance from transmitting devices to avoid interruption. Do not operate the ventilator in a magnetic resonance imaging (MRI) environment. Section 9 describes possible ventilator alarms and what to do if they occur. Consult with your institution's biomedical engineering department in case of interrupted ventilator operation, and before relocating any life-support equipment.

## Definitions

This manual uses three special indicators to convey information of a specific nature. The text under each heading states the reason the heading is used. They are listed below:

---

### Warning

Means there is a possibility for causing a potential safety hazard, which can result in an adverse reaction, or injury to the patient, you, or others.

---

### Caution

Means that special care is to be exercised to protect against the possibility of damage to the instrument or other property.

---

### NOTE:

Indicates points of particular interest for more efficient and convenient operation.

---

---

## Warnings, Cautions, and Notes

Please take the time to familiarize yourself with the following caveats as they cover safety considerations, special handling requirements, and regulations that govern the use of the Infant *Star* 500/950 ventilator systems.

- To ensure proper servicing and avoid the possibility of physical injury, only qualified personnel should attempt to service or make authorized modifications to the ventilator.

The user of this product shall have sole responsibility for any ventilator malfunction due to operation or maintenance performed by anyone not trained by Mallinckrodt staff.

- To avoid an electrical shock hazard while servicing the ventilator, be sure to remove all power to the ventilator by disconnecting the power source and turning off all ventilator power switches.

An electrical shock hazard is present inside the ventilator cabinet when the Vac power is connected, even though the power switch is in the OFF position. Always take appropriate safety precautions to avoid accidental shock.

The ventilator has a "hospital grade" Vac Power cord and plug. Grounding reliability can only be assured when connected to a tested receptacle identified as "hospital grade." Tampering with or reconfiguring the cord or plug supplied with the unit may void the factory warranty, cause damage to the unit, or be injurious to those associated with its use.

- To avoid a fire hazard, keep matches, lighted cigarettes, and all other sources of ignition (for example, flammable anesthetics and/or heaters) away from the Infant *Star* 500/950 ventilator systems and oxygen hoses.

Do not use oxygen hoses that are worn, frayed, or contaminated by combustible materials such as grease or oils. Textiles, oils, and other combustibles are easily ignited and burn with great intensity in air enriched with oxygen.

In case of fire or a burning smell, immediately disconnect the ventilator from the oxygen supply, facility power, and backup power source.

- When handling any part of the Infant *Star* 500 or Infant *Star* 950 ventilator system, always follow your hospital infection control guidelines for handling infectious material.

Mallinckrodt recommends that users of its products that require cleaning and sterilization/disinfection, consider the *National Standards and Recommended Practices for Sterilization* published by the Association for the Advancement of Medical Instrumentation (AAMI), as well as the following Center for Disease Control (CDC) publications: *Guidelines for Maintenance of In-use Respiratory Therapy Equipment and Guidelines for Prevention of Nosocomial Pneumonia*.

- Patients on life-support equipment should be appropriately monitored by competent medical personnel and suitable monitoring devices.

The Infant *Star* 500/950 ventilator systems are not intended to be comprehensive monitoring devices and do not activate alarms for all types of dangerous conditions for patients on life-support equipment

- For a thorough understanding of ventilator operations, be sure to thoroughly read this manual before attempting ventilator service.
- Before activating the ventilator, be sure to check the equipment for proper operation and, if appropriate, run a Quick Checkout as described in Section 3 of this manual.
- Federal law (U.S.) restricts the sale of this device to, or by the order of, a physician.
- Check the ventilator periodically as outlined in this manual; do not use if defective. Immediately replace parts that are broken, missing, obviously worn, distorted, or contaminated.

- 
- An alternate source of ventilation should always be available when using the *Infant Star* 500/950 ventilator systems.
  - The *Infant Star* 500/950 ventilator systems are a member of the *Infant Star* family of products.

### Intended Use

The *Infant Star* is designed to provide continuous ventilatory support for neonatal and pediatric patients less than 40 lbs. (18 kg), with an adjustable FLOW RATE range of 4 to 40 liters per minute (LPM), and BACKGROUND FLOW range of 2 to 30 LPM (2 to 32 LPM software versions 107 and above).

### Functional Considerations

Verify that the ventilator is Y2K compliant. If the ventilator is noncompliant, it may display incorrect dates.

### Shipping instructions

A shipping container for the ventilator is available. To ensure proper preparation and shipping, contact Mallinckrodt for instructions.

### Applicability

This service manual is applicable to *Infant Star* 500/950 ventilators only.

Mallinckrodt reserves the right to change specifications and procedures without notice.

# CONTENTS

<b>1</b>	<b>Repair Statement/Program Implementation</b>	
1.1	Infant <i>Star</i> 500/950 Ventilators . . . . .	1-1
1.2	Preventive Maintenance . . . . .	1-1
<b>2</b>	<b>Specifications</b>	
2.1	Infant <i>Star</i> 500/950 Ventilator Features . . . . .	2-1
2.2	System Specifications . . . . .	2-2
<b>3</b>	<b>Quick Checkout</b>	
3.1	Quick Checkout . . . . .	3-1
3.1.1	Infant <i>Star</i> 500 Quick Checkout Procedure . . . . .	3-1
3.1.2	Infant <i>Star</i> 950 Quick Checkout Procedure . . . . .	3-4
<b>4</b>	<b>Theory of Operation</b>	
4.1	Introduction . . . . .	4-1
4.1.1	Microprocessor Controlled Gas Delivery System . . . . .	4-1
4.1.2	Adjustable Background Flow . . . . .	4-3
4.1.3	Inadvertent PEEP Automatic Correction . . . . .	4-3
4.1.4	Built-In Patient Monitoring . . . . .	4-4
4.1.5	Tubing Obstruction Alarms . . . . .	4-4
4.1.6	Patient Data Output . . . . .	4-5
4.1.7	Battery Back Up . . . . .	4-5
4.1.8	Updatable Software . . . . .	4-5
4.1.9	Ventilator Functions . . . . .	4-5
<b>5</b>	<b>Preventive Maintenance</b>	
5.1	Preventive Maintenance/Extended Test . . . . .	5-1
5.1.1	When to Run . . . . .	5-1
5.2	Test Equipment and Service Material Required . . . . .	5-2
5.3	Preliminary Ventilator Cleaning and Inspection . . . . .	5-3
5.4	Preliminary Ventilator Set Up . . . . .	5-3
5.5	Preventive Maintenance/Extended Test Procedure . . . . .	5-4
5.5.1	Electrical Safety Test . . . . .	5-4
5.5.2	Power-On Self Test . . . . .	5-5
5.5.3	A05 Test (High Peak Inspiratory Pressure) . . . . .	5-5
5.5.4	Low Air and Low Oxygen Pressure Switch Test . . . . .	5-6
5.5.5	Check Valve Test . . . . .	5-7
5.5.6	Proximal Pressure Calibration . . . . .	5-7
5.5.7	Background Flow Control Test . . . . .	5-8
5.5.8	Ventilator Rate Control Test . . . . .	5-9
5.5.9	Peak Inspiratory Pressure . . . . .	5-10
5.5.10	Inspiratory Time Control Test . . . . .	5-11

# Contents

5.5.11	Duration of Positive Pressure Test	5-12
5.5.12	Expiratory Time and I:E Ratio Display Test	5-12
5.5.13	PEEP/CPAP Control Test	5-13
5.5.14	Oxygen Blender Test	5-15
5.5.15	Oxygen Blender Potentiometer Test	5-16
5.5.16	Crossover Test	5-17
5.5.17	Air/Oxygen Crossover Network Test	5-17
5.5.18	Blended Gas Outlet Test	5-18
5.5.19	Low Inspiratory Pressure Alarm, Audible and Alarm Silence Tests	5-18
5.5.20	Low PEEP/CPAP Alarm Test	5-19
5.5.21	Obstructed Tube Alarm Tests	5-20
5.5.22	High CPAP Alarm Test	5-21
5.5.23	Overpressure Pressure Relief Valve Test	5-21
5.5.24	Vent Inop Test	5-22
5.5.25	Airway Leak Alarm Test	5-23
5.5.26	HFV Rate Control Calibration and Test (Infant Star 950)	5-23
5.5.27	HFV Flow Test (Infant Star 950)	5-24
5.5.28	HFV Amplitude Test (Infant Star 950)	5-25
5.5.29	HFV PEEP/CPAP Control Test (Infant Star 950)	5-26
5.5.30	IMV + HFV Inadvertent PEEP Test	5-27
5.5.31	Internal Battery Operation Test	5-27
5.5.32	Hour Meter	5-28
5.6	Troubleshooting	5-29

## 6 Repair and Adjustment

6.1	Oxygen Regulator (R2) Adjustment	6-1
6.1.1	Tools, Equipment, and Supplies	6-1
6.1.2	Balance Air and Oxygen Pressures	6-1
6.1.2.1	Oxygen Regulator (R2) Adjustment	6-4
6.2	Check Valve (CV5) Leak Test and Safety Regulator (R8) Calibration	6-5
6.3	Proportional Valve Flow Solenoids/Repair and Maintenance	6-8
6.3.1	Tools, Equipment and Supplies	6-8
6.3.2	Background Flow Control Tests	6-8
6.3.3	Flow Solenoid Malfunctions	6-9
6.3.4	Maintenance/Procedure	6-10
6.4	Oxygen Blender Calibration	6-11
6.4.1	Tools, Equipment and Supplies	6-11
6.4.2	Ventilator Set Up	6-11
6.4.3	Oxygen Blender Quick Check	6-13
6.4.3.1	Check the 0.015-inch Restrictor for an Occlusion	6-14
6.4.4	Blender Flow Control Operation	6-15
6.4.5	O2 Control Knob Stop Alignment	6-18
6.4.6	Setting Oxygen Thresholds	6-18
6.4.6.1	Setting 21% Oxygen	6-18
6.4.6.2	Setting 100% Oxygen	6-19
6.4.7	Verify Oxygen Readings	6-20
6.5	Oxygen Blender Potentiometer Calibration	6-20
6.6	PEEP/CPAP Needle Valve Adjustment	6-21

6.7	Proximal Pressure Transducer Calibration .....	6-21
6.7.1	Tools, Equipment and Supplies .....	6-21
6.8	Venturi Needle Valve Adjustment .....	6-23
6.9	Analog Meter Potentiometer Adjustment .....	6-24
6.10	Replace Alarm Battery .....	6-25

## 7 Schematics

7.1	Printed Circuit Board Functionality, Schematics, and Circuit Diagrams .....	7-1
7.1.1	Introduction .....	7-1
7.1.2	Power Supply PCA .....	7-2
7.1.3	Analog PCA .....	7-5
7.1.4	MP2 PCA .....	7-7
7.1.5	MP101 PCA .....	7-9
7.1.6	Pressure Transducer PCA .....	7-11
7.1.7	Front Panel PCA .....	7-13
7.1.8	I/O (Input/Output) PCA .....	7-16
7.1.9	Alarm Driver PCA (Original) .....	7-18
7.1.10	Differential Pressure PCA .....	7-20
7.1.11	CE Mark Ventilators .....	7-22
7.1.11.1	Power Monitoring .....	7-22
7.1.11.2	Alarm Driver PCA - Current Limiting .....	7-22
7.1.11.3	Proximal/Differential Transducer PCA .....	7-25
7.1.12	Pneumatic Assembly Functional Descriptions .....	7-26
7.1.12.1	IMV Regulator (R5) .....	7-26
7.1.12.2	PEEP Regulator (R6) .....	7-26
7.1.12.3	Air Regulator (R1) .....	7-27
7.1.12.4	Oxygen Regulator (R2) .....	7-27
7.1.12.5	O2 Blender .....	7-27
7.1.12.6	Mini Regulator (R7) .....	7-27
7.1.12.7	Proportional Flow Manifold .....	7-27
7.1.12.8	Jet Venturi .....	7-27
7.1.13	Pneumatic Diagrams .....	7-28
7.1.14	Electrical Interconnect Diagrams .....	7-32
7.1.15	CE Mark Electrical Interconnect Diagrams .....	7-36

## 8 Hardware Identification

8.1	Hardware Identification .....	8-1
-----	-------------------------------	-----

## 9 Alarms and Troubleshooting

9.1	Software Diagnosed Alarms .....	9-1
9.2	Infant Star Troubleshooting Guide .....	9-4
9.2.1	IMV Standard Test .....	9-4
9.2.2	Obstructed Tube Alarms A01 Test, A01, A02, A05 .....	9-5
9.2.3	PEEP 0 to 24 cmH2O .....	9-6
9.2.4	PEEP Circuit .....	9-7
9.2.5	Venturi Solenoid (SV9) .....	9-8

# FIGURES

Figure 3-1.	Infant <i>Star</i> 500 Ventilator .....	3-3
Figure 3-2.	Infant <i>Star</i> 950 Ventilator .....	3-6
Figure 6-1.	Infant <i>Star</i> 500 Pneumatics Base Assembly .....	6-2
Figure 6-2.	Air and Oxygen Pressures Applied to Blender .....	6-3
Figure 6-3.	Removing Pipe Plug from Oxygen Regulator .....	6-4
Figure 6-4.	Air (R1), Oxygen (R2), and Safety (R8) Regulator Assemblies .....	6-7
Figure 6-5.	Pneumatics Base (rear view) .....	6-11
Figure 6-6.	Pneumatics Base with O <sub>2</sub> Analyzer .....	6-12
Figure 6-7.	Oxygen Blender 0.015-inch Restrictor .....	6-14
Figure 6-8.	Oxygen Blender and Flow Control Assembly (SV13 and Pilot Valve) .....	6-16
Figure 6-9.	Front Panel Display/Pneumatic Controls .....	6-17
Figure 6-10.	Proximal Pressure Transducer Adjustment .....	6-22
Figure 6-11.	Analog PCA and Analog Meter Gain Potentiometer .....	6-25
Figure 7-1.	Power Supply PCA .....	7-3
Figure 7-2.	Power Supply PCA Block Diagram .....	7-4
Figure 7-3.	Analog PCA .....	7-5
Figure 7-4.	Analog PCA Block Diagram .....	7-6
Figure 7-5.	MP2 PCA .....	7-7
Figure 7-6.	MP2 PCA Block Diagram .....	7-8
Figure 7-7.	MP101 PCA .....	7-9
Figure 7-8.	MP101 Block Diagram .....	7-10
Figure 7-9.	Pressure Transducer PCA .....	7-11
Figure 7-10.	Pressure Transducer PCA Block Diagram .....	7-12
Figure 7-11.	Front Panel PCA .....	7-14
Figure 7-12.	Front Panel Block Diagram .....	7-15
Figure 7-13.	I/O (Input/Output) PCA .....	7-17
Figure 7-14.	I/O (Input/Output) PCA Block Diagram .....	7-17
Figure 7-15.	Alarm Driver PCA .....	7-18
Figure 7-16.	Alarm Driver PCA Block Diagram .....	7-19
Figure 7-17.	Differential Pressure PCA .....	7-20
Figure 7-18.	Differential Pressure PCA Block Diagram .....	7-21
Figure 7-19.	Alarm Driver PCA (current limiting) .....	7-23
Figure 7-20.	Alarm Driver PCA (current limiting) Block Diagram .....	7-24
Figure 7-21.	Proximal/Differential PCA Block Diagram .....	7-26
Figure 7-22.	Infant <i>Star</i> 500 Inspiration Pneumatic Diagram .....	7-28
Figure 7-23.	Infant <i>Star</i> 500 Exhalation Pneumatic Diagram .....	7-29
Figure 7-24.	Infant <i>Star</i> 950 Inspiration Pneumatic Diagram .....	7-30
Figure 7-25.	Infant <i>Star</i> 950 Exhalation Pneumatic Diagram .....	7-31
Figure 7-26.	Original Infant <i>Star</i> Models 500/950 (1 of 4) .....	7-32
Figure 7-27.	Original Infant <i>Star</i> Models 500/950 (2 of 4) .....	7-33

# Figures

---

Figure 7-28.	Original Infant <i>Star</i> Models 500/950 (3 of 4)	7-34
Figure 7-29.	Original Infant <i>Star</i> Models 500/950 (4 of 4)	7-35
Figure 7-30.	Infant <i>Star</i> Models 500/950 230 V (1 of 4)	7-36
Figure 7-31.	Infant <i>Star</i> Models 500/950 230 V (2 of 4)	7-37
Figure 7-32.	Infant <i>Star</i> Models 500/950 230 V (3 of 4)	7-38
Figure 7-33.	Infant <i>Star</i> Models 500/950 230 V (4 of 4)	7-39
Figure 8-1.	Yoke Assembly	8-2
Figure 8-2.	Front Cover Assembly	8-4
Figure 8-3.	Rear Cover Assembly	8-5
Figure 8-4.	Pneumatics Base Assembly	8-7
Figure 8-5.	Pneumatic Console Assembly	8-10
Figure 8-6.	Relief Valve Assembly	8-12
Figure 8-7.	O <sub>2</sub> Blender Assembly	8-14

# TABLES

Table 2-1.	Infant <i>Star</i> 500/950 Comparison Chart .....	2-2
Table 2-2.	System Specifications .....	2-2
Table 4-1.	Solenoid Combinations to Set Flow .....	4-2
Table 5-1.	Performance Verification Intervals .....	5-1
Table 5-2.	Test Equipment Tools .....	5-2
Table 6-1.	Oxygen Regulator Calibration Chart .....	6-4
Table 6-2.	Safety Regulator Outlet Pressure .....	6-6
Table 6-3.	Background Flow Allowable Tolerances .....	6-9
Table 8-1.	Yoke Assembly Parts List .....	8-3
Table 8-2.	Front Cover Assembly Parts List .....	8-4
Table 8-3.	Rear Cover Assembly Parts List .....	8-5
Table 8-4.	Pneumatics Base Assembly Parts List .....	8-8
Table 8-5.	Pneumatic Console Assembly Parts List .....	8-10
Table 8-6.	Relief Valve Assembly Parts List .....	8-12
Table 8-7.	O <sub>2</sub> Blender Assembly Parts List .....	8-14
Table 9-1.	"Ventilator Inoperative" Alarm Troubleshooting Chart .....	9-1

# Repair Statement/Program Implementation

## 1.1 Infant Star 500/950 Ventilators

### Caution

Always adhere to electrostatic discharge control (ESD) procedures when servicing a ventilator.

Infant Star<sup>®</sup> 500/950 ventilators are designed, and technical support is provided, to allow service and repair by qualified personnel. This means personnel who routinely repair sophisticated electronic, mechanical, and pneumatic equipment; have been specifically trained to repair Mallinckrodt products.

If a medical facility and/or distributor does not provide service and repair, Mallinckrodt offers a complete service training program on Infant Star ventilators. For more information contact:

Mallinckrodt Inc.

1-800-635-5267

press 2 for Technical Support.

To support all levels of service and repair the following projects are on going:

- Comprehensive service manuals specialized tooling and testing instrumentation.
- Most Infant Star major subassemblies are available on an exchange basis.
- Mallinckrodt offers service and repair training classes to medical facility and distributor personnel. For more information contact Mallinckrodt Technical Support at the number listed above.

## 1.2 Preventive Maintenance

To reduce mechanical failures, it is important that a routine preventive maintenance program be followed to maintain a functional ventilator.

- Perform the Quick Checkout Procedure (described in Section 3) each time the ventilator is used for a new patient or, at a minimum, once a month.
- Perform the Extended Preventive Maintenance Test (described in Section 5) twice a year.

Regular service extends the life of the ventilator. Overhaul service is recommended after 10,000 hours or 7 years; whichever comes first.

Replace the pneumatic chassis at the 10,000 hour overhaul. During the 10,000-hour overhaul, other components are replaced based on age and usage, and install-free update. Electrical components are replaced or rebuilt to current specifications. For more information on overhaul programs, contact Mallinckrodt Technical Support at the number listed above.

# Specifications

## 2.1 Infant *Star* 500/950 Ventilator Features

The Infant *Star* 500/950 ventilators are microprocessor controlled. The Infant *Star* 500 is a conventional neonatal ventilator whereas the Infant *Star* 950 offers both conventional neonatal and high-frequency ventilation (HFV) modes. The Infant *Star* 950 ventilator has an expanded pneumatic system and added control functions. These additional features provide conventional ventilation, high-frequency oscillatory type of ventilation, or a mixing of both.

The design and construction of the Infant *Star* 500/950 ventilators are identical, with the exception of the following features (see Table 2-1):

- The Infant *Star* 950 control panel MODE selector switch has two additional selections:
  - HFV ONLY (high-frequency ventilation)
  - HFV + IMV (high-frequency ventilation plus intermittent mandatory ventilation breaths)
- The Infant *Star* 950 control panel has two additional ventilator control potentiometers and an LED display:
  - HFV RATE (operator control for setting HFV breathing rate)
  - HFV AMPLITUDE (operator control for setting flow amplitude in HFV mode)

The Infant *Star* 500 Control Panel does not have the additional components and graphic overlay covers that support HFV.

- The Infant *Star* 950 pneumatic system has four additional solenoid valves on the flow manifold in the manifold block: SV11 valves A and B, SV12 valves A and B. The solenoid low-voltage wiring is connected to the existing terminal block.
- To prevent unauthorized clinical use of the HFV mode, the Infant *Star* 950 has a keyed ON/OFF switch on the rear of the pneumatics compartment.
- The jet orifice, located in the exhalation block, is larger on the Infant *Star* 950 ventilator (non-CE mark) to accommodate HFV, as compared to the *Star* 500 ventilator (non-CE mark) which does not have HFV. Conversely, Infant *Star* 500/950 ventilators that are CE marked have the same size jet orifice. (CE mark ventilators have a needle valve on the venturi solenoid that adjust for HFV.)

**Table 2-1: Infant Star 500/950 Comparison Chart**

Feature	Infant Star 950	Infant Star 500	Comments
HFV Mode ON/OFF key switch	Available -- Harness connection to Front Panel PCB	Not available	Switch and harness deleted, chassis hole plugged on Infant Star 500
Mode Selector: HFV Only HFV + IMV	Available	Not available	
Controls: HFV Rate HFV Amplitude	Available	Not available	Potentiometers are deleted from existing front panel display PCB on Infant Star 500
LED Display: HFV Amplitude	Available	Not available	LED deleted from existing front panel display PCB on Infant Star 500
Solenoids: SOL11/SOL12 (2 ea.)	Installed on flow manifold	Not available	Valve ports plugged on flow manifold of Infant Star 500
Exhalation Block Venturi Jet Orifice	CE mark and Non-CE mark ventilators: 0.032-in. ID orifice jet (blue)	Non-CE mark ventilators: 0.020-in. ID orifice jet (silver) CE-mark ventilator: 0.032-in. ID orifice jet (blue)	Infant Star 500 has <i>standard</i> Venturi jet

## 2.2 System Specifications

Table 2-2 below lists the Infant Star 500/950 ventilator system specifications.

**Table 2-2: System Specifications\***

Control/Parameter	Range
<b>FRONT CONTROL PANEL Alarms</b>	
High Insp. Pressure	5 to 105 cmH <sub>2</sub> O (Infant Star 950 only)
Low Insp. Pressure	3 to 60 cmH <sub>2</sub> O
Low PEEP/CPAP	
PEEP/CPAP Setting (cmH <sub>2</sub> O)	Pressure difference for 25 seconds or longer between setting and proximal required to activate alarm (cmH <sub>2</sub> O)
0 to 5	2
6 to 8	3
9 to 12	4
13 to 24	5

Table 2-2: System Specifications\* (continued)

Control/Parameter	Range
<b>Alarms (continued)</b> Airway Leak Obstructed Tube  Insufficient Exp. Time  Ventilator Inop Low Oxygen Pressure Low Air Pressure Low Battery Ext Power Loss Alarm Silence Visual Reset	Leak makeup flow exceeds Background Flow by 13 Lpm for 4 seconds Detects any problems that cause: <ul style="list-style-type: none"> <li>• Proximal Pressure to rise to the HIP</li> <li>• Prolonged or impeded exhalation</li> <li>• PEEP/CPAP pressure to rise above the set point by 6 cmH<sub>2</sub>O for 5 seconds</li> <li>• Internal pressure to rise above set HIP by 10 cmH<sub>2</sub>O</li> <li>• Internal pressure to rise above set HIP by 5 cmH<sub>2</sub>O</li> </ul> 0.3 seconds at rates 100 bpm 0.2 seconds at rates 100 bpm  System malfunction or fully depleted internal battery < 45 psig (310 kPa) (35 psi for Infant <i>Star</i> 500 with CE mark) < 45 psig (310 kPa) (35 psi for Infant <i>Star</i> 500 with CE mark) Activates within 5 minutes prior to full-battery discharge Loss of ac power (battery life 30 minutes minimum - if fully charged) 60 seconds Push to clear alarm indicators (except Power Loss and Ventilator Inop)
<b>VENTILATOR SETTINGS</b>	
Modes	<ul style="list-style-type: none"> <li>• IMV</li> <li>• CPAP</li> <li>• HFV ONLY (Infant <i>Star</i> 950 only)</li> <li>• HFV+IMV (Infant <i>Star</i> 950 only)</li> </ul>
Ventilator Rate	1 to 150 bpm
Peak Insp. Pressure	5 to 90 cmH <sub>2</sub> O
Inspiratory Time	0.1 to 3 seconds
PEEP/CPAP	0 to 24 cmH <sub>2</sub> O
Flow Rate	4 to 40 Lpm
Background Flow	2 to 30 Lpm 2 to 32 Lpm (software version 107 and above)
High Insp. Pressure Alarm	5 to 105 cmH <sub>2</sub> O
Low Insp. Pressure Alarm	3 to 60 cmH <sub>2</sub> O
Manual Breath	Delivers single mandatory breath
HFV Rate (950 only)	2 to 22 Hz
HFV Amplitude (950 only)	Varies pulse flow from 12 to 120 Lpm

Table 2-2: System Specifications\* (continued)

Control/Parameter	Range
<b>PATIENT MONITOR</b>	
Proximal Airway Pressure	-10 to +120 cmH <sub>2</sub> O
Proximal PEEP/CPAP	-9 to +99 cmH <sub>2</sub> O
Peak Insp. Pressure (PIP)	0 to 110 cmH <sub>2</sub> O
Proximal Mean	-9.9 to +99.9 cmH <sub>2</sub> O
Selected Data	Push button control allows a single display to be used for multiple measurements
Dur Pos Press	0.02 to 5.0 seconds
Expiration Time	0.2 to 59.9 seconds
I:E Ratio	1:0.1 to 1:99.9
<b>LOWER SECTION</b>	
Pressure Relief Valve	5 to 120 cmH <sub>2</sub> O
Exhalation Block	Houses exhalation valve diaphragm and jet venturi
Oxygen %	21 to 100
<b>REAR PANEL CONTROL</b>	
Alarm	72 to 88 decibels
Hours (Elapsed Hours Meter)	0.0 to 99999.9 hours
To <i>Star Sync</i>	9-pin DIN connector
ac Circuit Breaker (100/115 V)	100 V – 2.0 amps 115 V – 1.2 amps
Mains/Battery Charger (220 V)	220 V – 0.6 amps
dc Circuit Breaker	12 V – 8.0 amps
12 Vdc input	11.8 to 18 V – 5 amp battery power source
ON/OFF Switch	Turns ventilator ON and OFF
RS-232 Serial Output Connection	1200 baud
Analog Pressure Output	10 mV / cmH <sub>2</sub> O – 0.1 to 1.2 V
Remote Alarm	30 V - 0.25 amps
<b>INPUTS</b>	
Electrical – Nominal	
100 V	90 - 105 Vac, 50/60 Hz, 100 W
115 V	105 - 123 Vac, 50/60 Hz, 100 W
220 - 240 V	198 - 264 Vac, 50/60 Hz, 100 W
Pneumatics	Air and oxygen, 45 to 90 psig (309 to 619 kPa)

Table 2-2: System Specifications\* (continued)

Control/Parameter	Range
<b>PHYSICAL CHARACTERISTICS</b>	
Ventilator	16.25" high x 11.25" wide x 19.25" deep (41 cm x 29 cm x 49 cm)
Weight	50 lbs (23 kg) – does not include pedestal stand
Pedestal Stand	44" high x 25" wide (112 cm x 46 cm)
Weight	14 lbs (6.4 kg)
* Specifications subject to change without notice.	

## Quick Checkout

### 3.1 Quick Checkout

#### 3.1.1 Infant *Star* 500 Quick Checkout Procedure

---

**NOTE:**

Reference Figure 3-1 when performing this procedure.

---

1. Connect both source gases and electrical power; attach the breathing circuit, humidifier and test lung.
  2. Turn circuit breaker (MAINS/Battery Charger) to the ON (I) position.
  3. Turn ventilator ON/OFF switch to ON and note the software revision number in the SELECTED DATA display.
- 

**NOTE:**

If during start up, any number under 103 appears in the SELECTED DATA display, the software level should be updated. Contact your local distributor or Mallinckrodt at 1-800-635-5267, press 2 for Technical Support.

---

4. All LED indicators are illuminated except Low Battery, Ventilator Inop., and Ext Power Loss. (The Ventilator Inop. may blink during power-up.)
5. The audible alarm is momentarily activated, verifying its function.
6. Set standard test conditions:

<b>Mode</b>	IMV
<b>Ventilator Rate</b>	30 bpm
<b>Peak Inspiratory Pressure</b>	40 cmH <sub>2</sub> O
<b>Inspiratory Time</b>	1.0 seconds
<b>PEEP/CPAP</b>	0 cmH <sub>2</sub> O
<b>Flow Rate</b>	20 Lpm
<b>Background Flow</b>	4 Lpm
<b>High Inspiratory Pressure</b>	45 cmH <sub>2</sub> O
<b>Low Inspiratory Pressure</b>	30 cmH <sub>2</sub> O
<b>Selected Data</b>	Expiratory Time (seconds)
<b>Oxygen Percent</b>	60
<b>Pressure Relief Valve</b>	Fully clockwise

**NOTE:**

After 30 seconds the MEAN display, located in the Patient Monitor portion of the control panel, is  $19 \pm 2$  cmH<sub>2</sub>O and the PIP display is  $40 \pm 1$  cmH<sub>2</sub>O.

7. During a mechanical breath, kink the From Patient tube to prevent gas flow. Verify that the proximal pressure drops abruptly when the A03 OBSTRUCTED TUBE alarm activates.
8. Disconnect the To Patient tube at the ventilator and block the connector located under the pressure-relief valve. Verify that the HI PIP, A05 OBSTRUCTED TUBE, HIGH INSP. PRESSURE and LOW INSP. PRESSURE alarms activate and the displays illuminate. Reconnect the tube, cancel the audible alarm, then push VISUAL RESET button. Verify PEEP/CPAP display is  $0 \pm 2$  cmH<sub>2</sub>O.
9. Enter the following settings:

<b>Mode</b>	CPAP
<b>PEEP/CPAP</b>	10 cmH <sub>2</sub> O
<b>Flow Rate</b>	8 Lpm
<b>Background Flow</b>	8 Lpm


**NOTE:**

After 30 seconds the Proximal Airway Pressure meter, PEEP/CPAP, MEAN, and PEEP/CPAP setting displays should read  $10 \pm 1$ .

10. After checking at 10 cmH<sub>2</sub>O, reset PEEP/CPAP to 20 cmH<sub>2</sub>O and again verify displays. After 30 seconds, the displays (mentioned in the note above) should read  $20 \pm 2$  cmH<sub>2</sub>O.
11. If the unit passes the Quick Checkout procedure, it is ready for patient use. Set the ventilator according to the physician's orders before connecting it to a patient.
12. To verify functionality of the audible/visual alarms for low oxygen and low-air pressure, momentarily close off gas sources. The alarms should activate when gas pressures are removed and reactivate when gas pressures are reapplied.


**Warning**

If a mechanical or electrical problem is discovered while running the Quick Checkout procedure or while operating the ventilator, the ventilator must be removed from use and serviced by a certified biomedical technician or person with equivalent experience and training in this type of equipment.



**NELLCOR PURITAN BENNETT**  
**Infant Star 500**

*Lower Best Console 1155409*



S-00110

Figure 3-1. Infant Star 500 Ventilator

### 3.1.2 Infant Star 950 Quick Checkout Procedure

#### Warning

If a mechanical or electrical problem is discovered while running the Quick Checkout procedure or while operating the ventilator, the ventilator must be removed from use and serviced by a certified biomedical technician or person with equivalent experience and training in this type of equipment.

#### NOTE:

The quick checkout procedure is based on using the following Mallinckrodt products:

- high-frequency breathing circuit (P/N – N-501023)
- standard test lung (P/N – N-1101262)

#### NOTE:

Reference Figure 3-2 when performing this procedure.

1. Connect the two source gases and electrical power, attach the patient circuit, humidifier, and test lung. Turn the ventilator power switch to the ON position.

Verify:

- All segments on all digital displays are illuminated and bright.
- All LED indicators are illuminated except EXT POWER LOSS, VENTILATOR INOP, and LOW BATTERY. VENTILATOR INOP may blink during power-up.
- An audible alarm momentarily sounds.
- A software revision number appears in the SELECTED DATA window.

2. Set the following standard test conditions:

Mode	IMV
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
Inspiratory Time	0.5 seconds
PEEP/CPAP	0 cmH <sub>2</sub> O
Flow Rate	20 Lpm
Background Flow	4 Lpm
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
HFV Rate	12 Hz
HFV Amplitude	Fully counterclockwise
Selected Data	Expiratory Time (seconds)
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Key	ON


Verify:

- No HFV oscillations can be heard
  - PEEP/CPAP is 0 cmH<sub>2</sub>O
  - Peak Inspiratory Pressure is 40 ±1 cmH<sub>2</sub>O
  - MEAN Airway Pressure is 9 ±1 cmH<sub>2</sub>O
  - Expiratory Time is 1.5 seconds
1. During a mechanical breath, kink the From Patient tube to restrict gas flow. Verify that the proximal pressure drops abruptly when the AO3 OBSTRUCTED TUBE alarm activates.
  2. Disconnect the To Patient tube at the ventilator connection and block the outlet. Verify the A05 OBSTRUCTED TUBE and LOW INSP. PRESSURE alarms and displays. Reconnect the tube; the alarm should stop. Push the VISUAL RESET button to clear the visual alarms.
  3. Turn the Mode switch to HFV + IMV. Rotate the HFV Amplitude control fully clockwise. High-frequency pulses should be visible on the analog gauge and audible between each mechanical IMV breath. High-frequency pulses should not occur during an IMV breath. HFV Amplitude display should register at least 34 cmH<sub>2</sub>O. HFV Amplitude levels will vary based on compressible volumes in the humidifier, breathing circuit, and test lung. For example, if the water level in the humidifier is low, the compressible volume increases and the amplitude decreases.
  4. Rotate the HFV Rate control knob to 20 Hz. HFV pulses should be rapid (too fast to count). Turn the HFV Rate control knob back to 2 Hz; HFV should be relatively slow (approximately 3 pulses between each IMV breath). Set the HFV Rate back to 12 Hz.
  5. Turn the HFV Amplitude control knob fully counterclockwise. The amplitude should incrementally decrease to a minimum of 6 to 13 cmH<sub>2</sub>O.
  6. Set the following:

Mode	CPAP
PEEP/CPAP	10 cmH <sub>2</sub> O
Background	8 Lpm
HFV Key	OFF


Verify: *LOW INS PRESS* Fully Counter Clockwise

- 10 cmH<sub>2</sub>O of PEEP appears in a total of five displays:
    - Proximal Airway Pressure gauge
    - PEEP/CPAP display in the Ventilator Settings portion of the control panel
    - MEAN Airway Pressure monitor display
    - PEEP/CPAP display in the Patient Monitor portion of the control panel
    - Low Inspiratory Pressure setting display
  - There should be no high-frequency oscillations.
7. After checking at 10 cmH<sub>2</sub>O, set the PEEP/CPAP control to 20 cmH<sub>2</sub>O; again verify that the five displays (listed above) read 20 cmH<sub>2</sub>O.



NELLCOR PURITAN BENNETT

**Infant Star 950**



5-00109

Figure 3-2. Infant Star 950 Ventilator

# Theory of Operation

## 4.1 Introduction

---

**NOTE:**

The following information provides an overview of the Infant *Star* 500/950 ventilators. For an operational description, reference Section 7 for pneumatic and electronic schematics, circuit diagrams, and their functional descriptions.

---

The Infant *Star* 500/950 ventilators have eight major functions that ensure safety and provide flexibility:

- full microprocessor control of gas delivery
- adjustable background-flow ventilation to reduce work of breathing
- automatic correction of inadvertent PEEP
- built-in patient monitoring
- five automatic alarms for tubing obstruction
- data output to a personal computer
- 30 minutes of internal battery power
- updateable software

### 4.1.1 Microprocessor Controlled Gas Delivery System

The Infant *Star* 500/950 ventilator gas-delivery system is microprocessor controlled. The microprocessor and software program precisely shape the flow and amount of pressure delivered to the patient.

The Infant *Star* 500/950 ventilators have solenoid valves. These solenoid valves are pneumatic switches that are electronically controlled. Gas from the blender enters the proportional manifold, where 6 solenoid valves (SV) on the Infant *Star* 500, produce the following liters per minute (Lpm) flow rates:

- SV1 = 2 Lpm
- SV2 = 4 Lpm
- SV3 = 8 Lpm
- SV4 = 16 Lpm
- SV5, SV6 = 32 Lpm

To produce the additional flow required by the Infant *Star* 950 for HFV, 4 additional solenoid valves in the proportional manifold operate with the other 6 solenoid valves:

- two SV11 valves A / B and two SV12 valves A / B = 64 Lpm

See Table 4-1 for solenoid combinations to control flow.

Table 4-1: Solenoid Combinations to Set Flow

2 Lpm	4 Lpm	6 Lpm	8 Lpm	10 Lpm	12 Lpm	14 Lpm	16 Lpm	30 Lpm	32 Lpm	34 Lpm	62 Lpm	64 Lpm
SV1	SV2	SV1 SV2	SV3	SV1 SV3	SV2 SV3	SV1 SV2 SV3	SV4	SV1 SV2 SV3 SV4	SV5/6	SV1 SV5/6	SV1 SV2 SV3 SV4 SV5/6	SV11 (2 ea) SV12 (2 ea)

When flow is required, the microprocessor activates one or more solenoid valves.

A mandatory breath is often examined by looking at the proximal airway pressure waveform.

A typical pressure curve is produced by a combination of flow rate, peak inspiratory pressure (PIP), and inspiratory time. The flow rate establishes the angle, or the slope, of the inspiratory pressure waveform for the mandatory breath. A high-flow rate produces a square-shaped waveform, while a lower-flow rate produces a sloped waveform.

In the Infant Star 500/950, PIP is an operator-selected value. When flow is established, the gas flow continues at the set rate until the pressure reaches 75% of PIP. At the 75% point, the microprocessor gradually reduces the gas flow, preventing pressure overshoot.

At 10 Lpm (operator-selected value), the flow waveform shows 8, 6, 4 and 2 Lpm and then goes to 0 (zero) flow. Since the endotracheal tube is not cuffed, the 2 Lpm solenoid turns ON and OFF to hold the pressure at the PIP point during the inspiratory plateau. Following inspiration, the exhalation valve opens and expiration occurs immediately. In the *Infant Star 500/950*, the pressure waveform is software modeled to optimize ventilation, regardless of the flow or breaths per minute (bpm) settings.

#### 4.1.2 Adjustable Background Flow

The adjustable background flow feature produces a low-level flow of 2 Lpm which is applied as a background flow to eliminate the possibility of rebreathing. A sensitive pressure transducer monitors pressures at the patient's airway. When the patient produces a spontaneous inspiratory effort, the signal is transmitted back to the microprocessor, which then activates one or more solenoids to provide additional gas to the airway (within 30 milliseconds).

---

**NOTE:**

For software version 105 and lower, the adjustable background flow can be set up to 30 Lpm  
For software version 107 and higher, flow can be set up to 32 Lpm. Flow is set in 2 Lpm increments for all software versions.

---

With the *Infant Star 500/950*, software versions 105 and lower, additional gas flow is in direct proportion to the patient's effort. For example, if the patient effort produces  $-1$  cmH<sub>2</sub>O at the airway, the patient will receive approximately 6 Lpm, a combination of the 2 and 4 Lpm solenoids. An effort of 5 cmH<sub>2</sub>O produces up to 34 Lpm. *Infant Star 500/950* software version 107 or higher has increased the responsiveness of the demand flow system. The demand flow system supplies extra flow to the circuit when the patient's spontaneous breathing exceeds the background flow setting. When the patient effort exceeds 1 cmH<sub>2</sub>O for 30 ms, flow is incremental. For every 30 ms that the pressure remains  $-1$  cmH<sub>2</sub>O, flow is incremented until the sub-baseline pressure is gone or 35 - 40 Lpm flow is reached.

The operator can set flow rates that satisfy the requirements of the mandatory breath, while at the same time, reducing the potential for increased expiratory work of breathing.

#### 4.1.3 Inadvertent PEEP Automatic Correction

The combination of the background flow, PEEP/CPAP, and ventilator rate settings (stored in microprocessor) predict the amount of inadvertent PEEP produced.

Based on this prediction, it increases or decreases the ON cycle of a solenoid directing gas flow through the Jet Venturi system which is built into the exhalation valve block. This Jet Venturi system produces a subambient pressure at the expiratory connection of the exhalation leg of the breathing circuit. When the jet flow is increased, more subambient pressure is produced; when it is decreased, there is less subambient pressure. Within 5 seconds of the pressure correction, the microprocessor checks the pressure at the patient's airway and compares that value to the PEEP/CPAP operator-set value. Based on this comparison, automatic compensation occurs until the inadvertent PEEP/CPAP is eliminated.

As the ventilator settings are adjusted, the microprocessor works with the pressure transducer to monitor the patient's airway, make corrections and hold the PEEP/CPAP level constant.

#### 4.1.4 Built-In Patient Monitoring

While other systems may require add-on equipment for patient monitoring, the Infant *Star* 500/950 ventilators have a built-in monitoring system.

Digital displays simultaneously show pressure measurements of PEEP/CPAP, mean airway pressure (which is a time weighted average), and PIP. The pressure signal is taken from the infant's airway through the proximal pressure tube to the pressure transducer and then to the microprocessor that does the calculations, integrates them, and controls the display.

The integration of ventilator pressure monitoring reduces the duplication of common components when monitors are purchased separately. This approach also allows us to include 2 pressure transducers: one at the airway and one inside the ventilator. Other systems, using pressure gauges or stand-alone monitors, are vulnerable to undetected high airway pressure caused by disconnection or blockage of the proximal pressure tube. Using two transducers allows a comparison of pressures, alerting the clinician, and preventing undetected high pressure. This system protects against proximal pressure tube blockage, disconnection, or even a single transducer failure. The electronics also work to protect the patient's airway in case of a circuit blockage.

#### 4.1.5 Tubing Obstruction Alarms

The operator-set peak inspiratory pressure (PIP) value, programs the microprocessor to limit pressure and compare the set value to what is actually occurring at the airway. Five different detection comparisons are checked, including 3 levels of Pressure Overshoot, the shape of the Expiratory Curve, actual CPAP versus the setting, and high resistance in the inspiratory leg of the breathing circuit. Under all of these conditions, the Infant *Star* 500/950 ventilator systems senses an obstructed tube and reports a specific error message. At the same time, the patient airway is protected from sustained pressure by shutting off gas flow.

- The A01 message indicates high inspiratory pressure (HIP), indicating proximal pressure has reached the HIP setting pressure. If the HIP pressure setting is momentarily reached, gas flow stops, the exhalation valve opens, and the patient airway is vented to PEEP/CPAP. The obstructed tube indicator activates both an audible alarm and visual signal. This type of alarm is most frequently caused by a patient fighting the ventilator, hiccuping, or crying.
- If the exhalation valve opens, but the pressure continues to rise due to a blockage in the expiratory leg of the breathing circuit, a second alarm occurs at HIP +5 cmH<sub>2</sub>O. An A02 alarm displays and an audible and a visual signal is activated. When this occurs, the circuit vents through the inspiratory leg of the circuit, thus protecting the patient's airway from sustained pressure.
- A third alarm occurs if something partially blocks the expiratory leg of the breathing circuit and prevents exhalation from occurring quickly. (Under normal operation, PIP drops more than 50% in the first 300 milliseconds.) If something is impeding the patient's exhalation, and pressure is not dropping quickly, the A03 message and an audible alarm activates, indicating elevated expiratory resistance. A complete pneumatic shutdown occurs for 1 second to vent the patient to ambient.
- The A04 message indicates CPAP 6 cmH<sub>2</sub>O greater than what has been set, for up to 5 seconds. This can occur when water, tubing fatigue or some other factor is blocking gas flow in the expiratory leg of the breathing circuit. If there is a complete pneumatic shutdown for 3 seconds, the patient is vented to ambient.
- The fifth detection point, A05, indicates a blockage in the inspiratory leg of the circuit or in the proximal pressure line. The ventilator alarms when the internal pressure transducer records a pressure 10 cmH<sub>2</sub>O above the set HIP.

Under the A01 alarm condition, the gas delivery stops and both the visual and audible obstructed tube alarms occur. The exhalation valve opens to vent the airway to PEEP. Under the A02, A03 and A04 alarm conditions, gas delivery stops and both the visual and audible

obstructed tube alarms activate. The exhalation valve and the internal vent system open to vent the airway pressure to ambient. This automatic response protects the patient's airway, while the error messages identify the most likely problem.

#### 4.1.6 Patient Data Output

Patient data can be sent directly to an IBM compatible personal computer and then to a printer. A Standard RS-232 cable, compatible with all PCs, connects to the output on the back panel of the ventilator. Control settings, airway pressures, and respiratory trending can be automatically recorded.

#### 4.1.7 Battery Back Up

The Infant *Star* 500/950 ventilators have a built-in battery backup, and operate on the battery at all times. When connected to an electrical outlet, a self-contained charger keeps the battery at full power. In the event of a power failure, the Infant *Star* 500/950 ventilators continue operating for a minimum of 30 minutes on battery power.

#### 4.1.8 Updatable Software

Microprocessor control of display calculations and pneumatics, allows major changes to be made to the software by replacing EPROMs installed on printed circuit boards.

#### 4.1.9 Ventilator Functions

The upper module contains the microprocessor electronics and the controls to establish the ventilator parameters.

The lower compartment contains the pneumatics.

Two modes of ventilation are available on the Infant *Star* 500:

- CPAP
- IMV for Infant *Star* 500

Four modes of ventilation are available on Infant *Star* 950:

- CPAP
- IMV
- HFV ONLY
- HFV IMV

Three additional patient-triggered modes are available for both the Infant *Star* 500 and 950 by attaching the *Star Sync* Interface:

- Assist/Control
- SIMV
- CPAP/backup

Dedicated displays identify the control settings:

Flow Rate	4 to 40 Lpm
Background Flow	2 to 30 Lpm (2-32 Lpm 107 software and higher)
Ventilator Rate	1 to 150 bpm
Peak Inspiratory Pressure	5 to 90 cmH <sub>2</sub> O
Inspiratory Time	1 to 3 seconds
Low Inspiratory Pressure	3 to 60 cmH <sub>2</sub> O
HIP Alarm	5 to 105 cmH <sub>2</sub> O
PEEP	0 to 24 cmH <sub>2</sub> O
Infant <i>Star</i> 950 ventilator	additional HFV Rate of 2 to 22 Hz and Amplitude control settings.

A push-button control allows a single display to be used for observing the I:E ratio, expiration time, or duration of positive pressure.

The three proximal pressure displays provide information about the actual pressure being applied at the patient's airway. The displays include PEEP/CPAP, MEAN Airway Pressure, and PIP.

The intensity of the audible alarm is adjustable. Visual displays show specific alarm violations. The alarm silence button provides a 60-second silence for all violations with the exception of continuous power loss and ventilator inoperative. When an alarm self corrects, the audible part of the alarm automatically stops while the visual indicator continues to provide information on what caused the problem. Pressing the Visual Reset button clears the visual alarm indicator.

Alarms are activated by:

- low inspiratory pressure
- low PEEP/CPAP
- airway leak
- obstructed tube
- insufficient expiratory time
- low oxygen pressure
- high inspiratory pressure
- low air pressure
- external power loss
- low battery
- ventilator inoperative

The lower chamber of the ventilator contains the pneumatic systems, including a mechanical pop-off valve, block assembly (exhalation valve, heater, jet venturi), and oxygen blender control.

The pneumatics base has both air and oxygen inlets on the back panel. Each is protected with water traps and filters. A single outlet provides the blended gas from the operator-set concentration (set via the oxygen control knob). The blended gas may be connected to an accessory flowmeter mounted on a rail system for hand resuscitation or aerosol delivery.

The Infant *Star* 950 has an HFV ON/OFF switch. The rear of the electronics module includes:

- ventilator ON/OFF switch
- hour meter
- RS-232 serial output port for a PC connection
- *Star Sync*<sup>™</sup> and remote alarm connector
- analog pressure outputs for connecting to a strip chart recorder
- external battery connection
- circuit breakers
- adjustable baffle for changing the audible alarm intensity

# Preventive Maintenance

## 5.1 Preventive Maintenance/Extended Test

### 5.1.1 When to Run

Table 5-1 lists the service performed and the Preventive Maintenance/Extended Test and calibration procedure that are performed after service. Run the electrical safety tests, which are a part of the extended test, according to your organization's requirements.

**Table 5-1: Performance Verification Intervals**

Service Action Performed	Performance Verification Test Requirements
For any service	Test 5.5.1
<b>Electronic Section</b>	
Rear Panel Assembly	Tests 5.5.1,20,32
Power Cord	Tests 5.5.1
Hour Meter	Test 5.5.32
Alarm	Tests 5.5.20
ac Circuit Breaker	Tests 5.5.1
<b>Yoke Assembly</b>	
MP-2 PCA	Test 5.5.32
Analog PCA	Tests 5.5.6,8,9,10,12,20
Alarm Driver (A/D) PCA	Tests 5.5.7,24
Power Supply (P/S) PCA	Tests 5.5.1
Transformer	Tests 5.5.1
10 V Battery Pack (Ext. Battery Operation)	Test 5.5.31
Front Panel (Infant Star 500)	Tests 5.5.6,8,9,10,12,20
Front Panel (Infant Star 950)	Tests 5.5.6,8,9,10,12,20,26
Analog Meter	Adjust R76 - Gain Adjust, Test 5.5.6

Table 5-1: Performance Verification Intervals (continued)

Service Action Performed	Performance Verification Test Requirements
<b>Pneumatic Section</b>	
High Pressure Gas Filter Assembly	R8 Regulator adjustment, Test 5.5.7
Air Regulator Assembly	Tests 5.5.4,5,7,15
Oxygen Regulator Assembly	Oxygen Regulator Calibration, Tests 5.5.4,5,7,15
Oxygen Blender Assembly	Oxygen Blender Calibration, Tests 5.5.7,15
Proportional Flow Manifold (Infant <i>Star</i> 500)	Tests 5.5.4,7,19,22.1
Proportional Flow Manifold (Infant <i>Star</i> 950)	Tests 5.5.4,7,19,22.1,26,27,28
IMV Regulator Assembly	Verify output, Test 5.5.9
PEEP Regulator Assembly	Test 5.5.13
Mini Regulator	Verify output is 13 psi, Tests 5.5.4,13
Transducer PCA	Tests 5.5.4,7,19,22.1,26,27,28
SV8 PEEP, SV9 Venturi Solenoids (Infant <i>Star</i> 500)	SV8 - PEEP Calibration, Tests 5.5.13,14
SV8 PEEP, SV9 Venturi Solenoids (Infant <i>Star</i> 950)	SV8 - PEEP Calibration, Tests 5.5.13,14,26
SV7 Safety Vent Valve Solenoid	Test 5.5.22.1

## 5.2 Test Equipment and Service Material Required

The tools listed in Table 5-2 are required for total performance verification and calibration. Refer to specific tests for individual tool requirements.

Table 5-2: Test Equipment Tools

Gauge, Capsuhelic -5 - +5 cmH <sub>2</sub> O (P/N - N-4300010)
<i>PTS 2000 (or equivalent) OR RT200</i>
Adjuster Tool, Oxygen Regulator (P/N - N-97101016)
Tee Assembly, (2 required) (P/N - N-9710161)
<b>Manometers: (or equivalent)</b>
0 to 100 cmH <sub>2</sub> O
Accuracy: ±1% of full scale
0 to 10 cmH <sub>2</sub> O (Positive or Negative Pressure)
Accuracy: ±1% of full scale
Flowmeter: 0 to 200 Lpm, Float or Ball type. Dwyer RMC 103 (or equivalent)
Flowmeter: 0 to 50 Lpm, Float or Ball type. Dwyer RMC 101 (or equivalent)
Oxygen Analyzer: MiniOx 3000 (or equivalent)
Tubing: 1/8 in., and 10 mm, as required
Syringe: 50 ml
Tee Adapter: 1/8 in. pneumatic

For 950  
HF

Table 5-2: Test Equipment Tools (continued)

Tee Adapter:	10 mm pneumatic
Adapter:	1/8 in. to 10 mm
Frequency Counter (use Fluke Volt Meter) for Infant <i>Star</i> 950 HFV Rate	
Digital Volt Meter:	3-1/2 Digit Fluke 87 or similar
<b>Patient Circuit:</b>	
Standard circuit:	P/N – N-501011 (Infant <i>Star</i> 500/950)
HFV circuit:	P/N – N-501023 (Infant <i>Star</i> 950)
Infant <i>Star</i> Test Lung:	P/N – N-1101262
Rubber Test Lung:	P/N – 4-000612-00 (or similar)
Analyzer:	Electrical Safety
<b>Regulators:</b>	
With gauge 0 - 50 psi air supply	
With gauge 0 - 50 psi oxygen supply	
Static Control System	
<b>Miscellaneous Supplies (or equivalent)</b>	
Isopropyl Alcohol for cleaning	
Cotton Swabs	
Loctite Adhesive Sealant, #271, Red	
GE RTV 118 Silicone	
Torque Seal	
Devcon 14210 5-Minute Epoxy	
Dow Corning 340 Silicone Heat Sink Compound	
Fisher & Paykel low compressible volume humidifier chamber (P/N – N-501209)	

### 5.3 Preliminary Ventilator Cleaning and Inspection

To prevent disease transmission, use personal protective equipment when handling contaminated bacterial filters or other patient accessories.

**NOTE:**

If any problems are found during the preliminary inspection, correct them before proceeding with the performance verification. The failure to correct such problems now may affect the remainder of the performance verification.

Clean and inspect the ventilator for obvious problems such as missing broken parts, loose assemblies or disconnected wires, connectors or tubing. Repair as needed.

### 5.4 Preliminary Ventilator Set Up

Set up the ventilator for the Extended Test:

1. Install a complete infant circuit and Infant *Star* test lung:
  - P/N – N-501011 for the Infant *Star* 500/950 patient circuit
  - P/N – N-501023 for the Infant *Star* 950 HFV patient circuit
  - P/N – N-501209 for the Fisher & Paykel low compliance humidifier chamber
  - P/N – N-1101262 for the Infant *Star* test lung
2. Connect the ventilator to both an air and oxygen source. (Ventilator hour meter reading should be observed at this time.)

## 5.5 Preventive Maintenance/Extended Test Procedure

To ensure systematic performance verification and logical fault diagnosis, perform these tests in the order given. If a test must be repeated, the current control settings are completely defined at the beginning of each test procedure.

---

**NOTE:**

- If running the performance verification tests in the order written, only the **boldface** ventilator settings must be changed.
  - If a test fails, identify the cause of a malfunction by noting the item number at the end of the step; then refer to that item number in Section 5.6. Perform tasks in the referenced item to identify and correct the causes of malfunctions.
- 

Follow these general guidelines when running the performance verification:

- If a problem occurs during the performance verification, check to verify that the procedures were correctly performed before attempting to repair the ventilator.
- Verify all data entered.
- Use alarm silence and alarm reset keys to restore the ventilator to the test condition. Using these keys may be necessary because of the nonstandard set of conditions that may cause an alarm condition to occur.
- When making ventilator settings, be aware that because of the inter-relationship between some ventilator settings, you may not always be able to make all settings in the indicated sequence.
- Except for the alarm silence and alarm-reset keys, do not change the control settings during these procedures, unless specifically instructed.
- When repairs are completed, repeat the test. When the test is successful, proceed with the next test, as applicable.

---

**Warning**

Follow accepted safety procedures for electrical equipment when making connections, adjustments, or repairs.

---

### 5.5.1 Electrical Safety Test

---

**Warning**

If the ventilator fails an electrical safety test, do not proceed to the next electrical safety test until the problem is corrected and the test is repeated.

---

**NOTE:**

There are no specific ventilator setting requirements for the electrical safety test.

---

1. Make sure the unit is turned OFF.
2. Verify ground resistance is 0.1  $\Omega$  or less (See Section 5.6, item 1.)
3. Turn ON the ventilator.
4. Verify that forward and reverse current leakage to ground is 25 mA on 115 Vac units and 40 mA on 220 Vac units. (See Section 5.6, item 2.)

### 5.5.2 Power-On Self Test

1. Turn ON the power switch.
2. Set the unit to the following Standard Test Conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	30 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

3. Verify that all displays illuminate (except Low Battery and Ext Power Loss). Vent Inop may blink at initial power-up. The software level momentarily appears in the Selected Data window. (See Section 5.6, item 3.)

### 5.5.3 A05 Test (High Peak Inspiratory Pressure)

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	30 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Disconnect To Patient tube at manifold and block the connector on the bottom of the Pressure Relief valve.

3. Verify HI-PIP, A05, Obstructed Tube and Low Inspiratory Pressure Alarms and Displays. (See Section 5.6, item 4.)
4. Reconnect tube, allow audible alarm to self-cancel, then push Visual Reset.
5. Set the following test conditions:
 

<b>Mode</b>	<b>CPAP</b>
<b>PEEP/CPAP</b>	<b>10 cmH<sub>2</sub>O</b>
<b>Flow Rate</b>	<b>8 Lpm</b>
<b>Background Flow</b>	<b>8 Lpm</b>
6. Verify that PEEP/CPAP reads 10 cmH<sub>2</sub>O in the following four displays:
  - Proximal Airway Pressure gauge
  - PEEP/CPAP window
  - MEAN Pressure window
  - PEEP/CPAP setting display
7. After checking at 10 cmH<sub>2</sub>O, reset PEEP/CPAP to 20 cmH<sub>2</sub>O, and again verify displays. (See Section 5.6, item 4.)

#### 5.5.4 Low Air and Low Oxygen Pressure Switch Test

1. Set the following test conditions:

<b>Mode</b>	<b>CPAP</b>
<b>PEEP/CPAP</b>	<b>Fully counterclockwise</b>
<b>Flow Rate</b>	<b>20 Lpm</b>
<b>Background Flow</b>	<b>4 Lpm</b>
<b>Ventilator Rate</b>	<b>30 bpm</b>
<b>Peak Inspiratory Pressure</b>	<b>40 cmH<sub>2</sub>O</b>
<b>High Inspiratory Pressure</b>	<b>45 cmH<sub>2</sub>O</b>
<b>Inspiratory Time</b>	<b>1.0 second</b>
<b>Low Inspiratory Pressure</b>	<b>30 cmH<sub>2</sub>O</b>
<b>Display Select</b>	<b>Expiration Time</b>
<b>Oxygen Percent</b>	<b>60</b>
<b>Pressure Relief Valve</b>	<b>Fully clockwise</b>
<b>HFV Rate</b>	<b>12 Hz</b>
<b>HFV Amplitude</b>	<b>Fully clockwise 3 - 4 turns (no hard stop)</b>

---

#### NOTE:

If using a *PTS 2000*, the HIGH PRESSURE function is required for the measurements in this section; all readings are the "current" reading.

---

2. Reduce air pressure to 40 psi.
3. Verify that the Low Air Pressure alarm is ON.
4. Raise the air pressure to 45 psi.
5. Verify that the Low Air Pressure alarm is OFF. (See Section 5.6, item 5.)

6. Repeat steps 2 - 5 for oxygen pressure. (See Section 5.6, item 6.)  
Ensure that the Safety Vent Valve does not open.

### 5.5.5 Check Valve Test

1. Set up the ventilator as follows:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	30 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Disconnect air and oxygen high-pressure hoses alternately at the ventilator, and ensure that there is no leakage from the connection while the other gas system is applying pressure.
3. To test the check valves, insert the disconnected end of the high-pressure hose into a cup of water.
4. If there are more than fifty bubbles in 10 seconds, the leak is significant and the check valve should be replaced. (See Section 5.6, item 7.)

### 5.5.6 Proximal Pressure Calibration

1. Set the following test conditions:

Mode	CPAP
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	30 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude x	Fully clockwise 3-4 turns (no hard stop)

2. Auto-Zero the ventilator by disconnecting the proximal airway pressure line and power the ventilator OFF, then ON, 3 times.
3. Verify that the proximal airway meter and the PEEP/CPAP display read 0 cmH<sub>2</sub>O ±0.50. (See Section 5.6, item 8.)
4. Using the *PTS 2000* (or equivalent), tee a syringe to the + low-pressure port and the proximal airway pressure line. ~~Pressure Transducer~~
5. Apply pressure, using a syringe, until the desired pressure is displayed.

Apply (cmH <sub>2</sub> O)	Measure Limits
5	4 to 6
25	24 to 26
50	48 to 52
75	73 to 77
95	92 to 98

6. Verify that the PEEP/CPAP display and Analog meter readings are within specifications. (See Section 5.6, item 8.)  
Ignore the Obstructed Tube Alarm during this test.
7. When the ventilator digital readings are verified as correct, it will not be necessary to use an external pressure gauge in future tests.

### 5.5.7 Background Flow Control Test

1. Set the following test conditions:

Mode	CPAP
PEEP/CPAP	Fully counterclockwise
Flow Rate	40 Lpm
Background Flow	See below
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	30 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	21
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

#### NOTE:

The background control increments at 2 Lpm (2 - 30 Lpm for software version 105 or lower) and (2 - 32 Lpm for software version 107 or higher).

2. Using the *PTS 2000* (or equivalent), measure the flow from the To Patient port against the values listed in the chart below. Use the Air and High Flow function to take the current reading. (See Section 5.6, item 9.)
3. If running software 107, disconnect the proximal airway pressure line.

**NOTE:**

- The flow rate must be set higher than, or equal to, the background flow.
- The proximal PEEP must read 0 (zero) for the flow measurement to be accurate.

Background Flow Rate Control (Lpm)	Measure Limits (Lpm)
2	1.8 to 2.2
4	3.6 to 4.4
8	7.2 to 8.8
16	14.4 to 17.6
30	27.0 to 33.0
32	(Software 107 or higher) 28.8 to 35.2

Disconnect  
Prox airway  
Tubing Press

**5.5.8 Ventilator Rate Control Test**

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	See below
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	30 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Verify that Ventilator Rate display increments:
  - Between 1 and 60, the unit increments in 1 bpm steps.
  - Between 60 and 130, the unit increments in 2 bpm steps.
  - Between 130 and 150, the unit increments in 5 bpm steps.
3. Set VENT RATE to 30 bpm, and measure the time for 15 breaths.  
Test Limits: 29 to 31 seconds. (See Section 5.6, item 10.)

### 5.5.9 Peak Inspiratory Pressure

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	See below
Background Flow	See below
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	See below
High Inspiratory Pressure	105 cmH <sub>2</sub> O (blinking)
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	21
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Check that PIP increments in 1 cmH<sub>2</sub>O steps, from 5 to 90.
3. Using the chart below set the PIP and Flow Rate. Set Background Flow to 4 Lpm. Repeat the test using a Background Flow of 20 Lpm. (See Section 5.6, items 11, 12.)

PIP (cmH <sub>2</sub> O)	Flow Rate (Lpm)	Background Flow (Lpm)	PIP Limits (cmH <sub>2</sub> O)
10	20	4 then 20	9 to 11
20	20	4 then 20	18 to 22
40	20	4 then 20	38 to 42
60	20	4 then 20	58 to 62
80	20	4 then 20	78 to 82
20	10*	4 only	18 to 22

**NOTE:**

If the difference between the PEEP/CPAP and PIP settings is < 5 cmH<sub>2</sub>O, the ventilator will not accept the PIP setting and the PIP value display will blink. Perform the steps below to verify this range.

4. Set PEEP/CPAP to 10 cmH<sub>2</sub>O.
5. Set PIP to 14 cmH<sub>2</sub>O (PIP value should blink).
6. Adjust PIP to 18 cmH<sub>2</sub>O (PIP value should stop blinking).

### 5.5.10 Inspiratory Time Control Test

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	15 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	See below
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Verify the Inspiratory Time display increments:

- 1 step between 0.10 to 0.60 seconds
- 2 steps between 0.60 to 1.0 seconds
- 10 steps between 1.0 to 3.0 seconds

3. Set mode:

Ventilator Rate	30 bpm
Inspiratory Time	1.0 seconds

Observe the Expiratory Time in the Selected Data window while measuring the Inspiratory Time at the Analog Output connector on the rear cover of unit (pins 1 and 3).

#### NOTE:

If using a DVM for the next two measurements, first verify the frequency is 0.5 Hz, then verify the duty cycle is between 40% and 60%. *USE J18 on Alarm Driver Board*

4. Verify the Expiratory Time display reads 1.0 seconds. (See Section 5.6, item 13.)
5. Verify the Inspiratory Time reads between 0.95 to 1.1 second on the measurement device. (See Section 5.6, item 13.)

### 5.5.11 Duration of Positive Pressure Test

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
<b>Display Select</b>	<b>Duration of Positive Pressure</b>
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Verify the duration of positive pressure reads between 1.13 and 1.19 in the Selected Data window. (See Section 5.6, item 14.) *Most units 1.12 - 1.13*

### 5.5.12 Expiratory Time and I:E Ratio Display Test

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
<b>Flow Rate</b>	<b>10 Lpm</b>
Background Flow	4 Lpm
<b>Ventilator Rate</b>	<b>See below</b>
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
<b>Inspiratory Time</b>	<b>0.4 second</b>
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
<b>Display Select</b>	<b>Expiration Time</b>
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Using the following chart, adjust the Ventilator Rate and verify the Expiration Time display and Insufficient Expiration Time alarm. (See Section 5.6, item 15.)

Ventilator Rate (bpm)	Expiration Time (seconds)	Insufficient Expiration Time Alarm
1	59.6	OFF
62	0.57	OFF
88-98*	0.30	ON
100	0.20	OFF
102-150*	0.20	ON

\*Ventilator rate settings will blink during test.

3. Change Data Select to I:E Ratio.
4. Adjust Ventilator Rate and verify that the I:E Ratio display corresponds to:  
(See Section 5.6, item 15.)

Ventilator Rate (bpm)	I:E Ratio
1	1:99.9 blinks
6	1:24.0
33	1:3.5
102-150	1:---

### 5.5.13 PEEP/CPAP Control Test

1. Set the following test conditions:

Mode	CPAP
PEEP/CPAP	See below
Flow Rate	40 Lpm
Background Flow	See below
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Verify the Proximal Airway meter and, while turning the PEEP/CPAP knob, verify that the PEEP/CPAP Range includes 0 to 20 cmH<sub>2</sub>O ±1. (In most cases PEEP will go to 24 cmH<sub>2</sub>O maximum, but 20 cmH<sub>2</sub>O is acceptable.)

3. Adjust Background Flow and PEEP/CPAP setting according to the table below.
4. Verify the PEEP/CPAP Proximal Display for 5 seconds after stabilization, and ensure that PEEP/CPAP setting, and the PEEP/CPAP values are within their specified ranges.  
(See Section 5.6, item 17.)

Background Flow (Lpm)	PEEP/CPAP Level (cmH <sub>2</sub> O)	Range
4	0	±1
4	2	±1
4	6	±1
4	15	±1
4	20	±1
20	2	±1
20	6	±1
20	15	±1
20	20	±1
30	2	±2
30	6	±2
30	15	±2
30	20	±2

*Calibrate  
Venturi  
if High*

5. Set the following test conditions:

<b>Mode</b>	<b>IMV</b>
<b>PEEP/CPAP</b>	<b>Fully counterclockwise</b>
<b>Flow Rate</b>	<b>40 Lpm</b>
<b>Background Flow</b>	<b>See below</b>
<b>Ventilator Rate</b>	<b>See below</b>
<b>Peak Inspiratory Pressure</b>	<b>40 cmH<sub>2</sub>O</b>
<b>High Inspiratory Pressure</b>	<b>45 cmH<sub>2</sub>O</b>
<b>Inspiratory Time</b>	<b>See below</b>
<b>Low Inspiratory Pressure</b>	<b>3 cmH<sub>2</sub>O</b>
<b>Display Select</b>	<b>Expiration Time</b>
<b>Oxygen Percent</b>	<b>60</b>
<b>Pressure Relief Valve</b>	<b>Fully clockwise</b>
<b>HFV Rate</b>	<b>12 Hz</b>
<b>HFV Amplitude</b>	<b>Fully clockwise 3 - 4 turns (no hard stop)</b>

6. Set Background Flow, Inspiratory Time, and Ventilator Rate values shown in the table below; then verify that PEEP/CPAP display reads less than the maximum allowed. (See Section 5.6, item 17.)

Flow Rate (Lpm)	Background Flow Lpm)	Inspiratory Time (seconds)	Ventilator Rate (bpm)	Maximum Allowed PEEP/CPAP (cmH <sub>2</sub> O)
40	30	0.50	60	5
30	30	0.33	90	2
20	20	0.15	120	2
10	10	0.10	150	2

### 5.5.14 Oxygen Blender Test

1. Blender Dial Alignment:
  - a. Turn the O<sub>2</sub> control knob fully counterclockwise.
  - b. Check that the mark on the O<sub>2</sub> control knob aligns with the 21 percent mark on the console base. (Mark is considered aligned if any part of one mark makes contact with the other.)
  - c. Turn O<sub>2</sub> control knob fully clockwise.
  - d. Check that the mark on the O<sub>2</sub> control knob aligns with the 100 percent mark on the console base. (Mark is considered aligned if any part of one mark makes contact with the other.)
2. With air and oxygen connected to the ventilator, connect the ventilator To Patient port to the High Flow port on the *PTS 2000* (or equivalent). Take readings in the % Oxygen and current function.
3. Set the following test conditions:

<b>Mode</b>	<b>CPAP</b>
<b>PEEP/CPAP</b>	Fully counterclockwise
<b>Flow Rate</b>	<b>20 Lpm</b>
<b>Background Flow</b>	<b>See below</b>
<b>Ventilator Rate</b>	<b>30 bpm</b>
<b>Peak Inspiratory Pressure</b>	40 cmH <sub>2</sub> O
<b>High Inspiratory Pressure</b>	45 cmH <sub>2</sub> O
<b>Inspiratory Time</b>	<b>1.0 second</b>
<b>Low Inspiratory Pressure</b>	3 cmH <sub>2</sub> O
<b>Display Select</b>	Expiration time
<b>Oxygen Percent</b>	<b>See below</b>
<b>Pressure Relief Valve</b>	Fully clockwise
<b>HFV Rate</b>	12 Hz
<b>HFV Amplitude</b>	Fully clockwise 3 - 4 turns (no hard stop)

4. With the ventilator connected to air and oxygen supply, attach a calibrated oxygen analyzer to the To Patient port, leaving one side of the analyzer's sensor with a 12-inch tube open to atmosphere. (See Section 6, Figure 6-6.)
5. If the ventilator has software 107 installed, disconnect the proximal airway pressure line and zero the Pressure Transducer PCA. *→ Nit Alarm Silence*
6. Set Background Flow to 16 Lpm and verify readings. (See Section 5.6, item 18.)
7. Repeat the test at 4 Lpm.

Oxygen Percent	Measured Value
21	21 to 22
30	27 to 33
60	57 to 63
90	87 to 93
100	99 to 100

### 5.5.15 Oxygen Blender Potentiometer Test

1. Set the following test conditions:

<b>Mode</b>	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
<b>Display Select</b>	<b>I:E Ratio</b>
<b>Oxygen Percent</b>	<b>60</b>
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

#### Caution

Be sure to perform this test while unit is plugged into ac power. Performing this test using battery power may give erroneous readings.

2. Be sure to perform this test while unit is plugged into ac power. Performing this test under battery power may give erroneous readings.
3. Verify: Diagnostics window 2 displays the Oxygen Blender Potentiometer setting. It should read  $\pm 3\%$  of what the O<sub>2</sub> control knob indicates. (See Section 5.6, item 19.)
4. To adjust the displayed value, loosen the Allen screw on the white gear on the potentiometer.

5. By pulling the white gear forward, the shaft of the Potentiometer can be rotated to display the appropriate number in the diagnostic window.
6. Push the white gear to its original position on the shaft and tighten Allen screw.

### 5.5.16 Crossover Test

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	40 Lpm
Background Flow	30 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	21
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Using a *PTS 2000* (or equivalent), measure high-pressure readings using the current reading for the High Pressure function.
  - a. Set air line pressure to 75 psig (or maximum possible considering environmental conditions).
  - b. Set oxygen line pressure to 40 psig.
3. Using the same test setup as in 5.5.15, verify that oxygen concentration does not rise above 22%. (Ignore the LOW O<sub>2</sub> Alarm.) (See Section 5.6, item 20.)

### 5.5.17 Air/Oxygen Crossover Network Test

1. Set the following test conditions:

Mode	CPAP
PEEP/CPAP	Fully counterclockwise
Flow Rate	40 Lpm
Background Flow	30 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Disconnect the Proximal Airway Pressure line.
3. Using a *PTS 2000* (or equivalent), connect a hose between the High Flow port and the To Patient port on the ventilator.
4. Set the Flow Rate to 40 Lpm and Background Flow to 30 Lpm.
5. Disconnect air pressure source and allow ventilator to stabilize for approximately 10 seconds.
6. Verify the High Flow reading is > 26 Lpm. (See Section 5.6, item 21.)
7. Reconnect air source, allow ventilator to stabilize.
8. Disconnect the oxygen pressure source and allow the ventilator to stabilize.
9. Verify the High Flow reading is > 26 Lpm. (See Section 5.6, item 21.)
10. Reconnect the oxygen pressure source.

### 5.5.18 Blended Gas Outlet Test

1. Turn Ventilator Power OFF.
2. Using the *PTS 2000* (or equivalent), connect a hose between the High Pressure port and the blended-gas port on the ventilator. - *Back of Unit*
3. Verify the current pressure is 18 psi  $\pm$  1 psi. (See Section 5.6, item 22.)
4. Using the *PTS 2000* (or equivalent), connect a hose between the High Flow port and the ventilator blended-gas port.
5. Verify the Blended Gas Outlet reads between 18 and 22 Lpm (or  $\pm$ 10% of the value stated on label). (See Section 5.6, item 22.)

### 5.5.19 Low Inspiratory Pressure Alarm, Audible and Alarm Silence Tests

1. Set the following test conditions:

<b>Mode</b>	<b>IMV</b>
<b>PEEP/CPAP</b>	Fully counterclockwise
<b>Flow Rate</b>	<b>20 Lpm</b>
<b>Background Flow</b>	<b>4 Lpm</b>
<b>Ventilator Rate</b>	30 bpm
<b>Peak Inspiratory Pressure</b>	40 cmH <sub>2</sub> O
<b>High Inspiratory Pressure</b>	45 cmH <sub>2</sub> O
<b>Inspiratory Time</b>	1.0 second
<b>Low Inspiratory Pressure</b>	See below
<b>Display Select</b>	Expiration Time
<b>Oxygen Percent</b>	60
<b>Pressure Relief Valve</b>	Fully clockwise
<b>HFV Rate</b>	12 Hz
<b>HFV Amplitude</b>	Fully clockwise 3 - 4 turns (no hard stop)

2. Check that Low Inspiratory Pressure changes are in increments of 1 cmH<sub>2</sub>O between 3 and 60 cmH<sub>2</sub>O, and the readings are stable. (Ignore alarm.)
3. Set the Low Inspiratory Pressure to 42 cmH<sub>2</sub>O.
4. The Low Inspiratory Pressure alarm should activate. (See Section 5.6, item 23.)
5. Press ALARM SILENCE.
6. Verify that ALARM SILENCE period equals 1 minute ( $\pm 3$  sec.)
7. Press ALARM SILENCE twice.
8. Verify that Audible turns OFF and back ON. (See Section 5.6, item 23.)
9. Open Volume Control baffle.
10. Verify loud "beeps" (back panel) and low steady tone (alarm driver PCA). (See Section 5.6, item 23.)
11. Close the baffle and verify that the alarm signal is audible from a distance of six feet, while facing the front of the ventilator.
12. Reduce the Low Inspiratory Pressure Alarm setting to 28 cmH<sub>2</sub>O.
13. Verify that audible alarm turns OFF, and Low Inspiratory Pressure LED stays ON until Visual Reset button is pressed. (See Section 5.6, item 23.)

#### 5.5.20 Low PEEP/CPAP Alarm Test

1. Set the following test conditions:

<b>Mode</b>	<b>CPAP</b>
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
<b>Background Flow</b>	<b>10 Lpm</b>
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
<b>Low Inspiratory Pressure</b>	<b>3 cmH<sub>2</sub>O</b>
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Disconnect the hose from the proximal airway pressure inlet on the front of the pneumatics unit.
3. Connect a syringe to the proximal airway pressure inlet fitting; adjust the syringe pressure to zero PEEP/CPAP.

4. Perform the following steps for each PEEP/CPAP setting listed below:

PEEP/CPAP Setting (cmH <sub>2</sub> O)	Alarm ON (cmH <sub>2</sub> O)	Alarm OFF (cmH <sub>2</sub> O)
0	-2	-1
2	0	1
6	3	4
9	5	6
13	8	9
20	15	16

5. Set syringe pressure (as shown on the PEEP/CPAP display setting) to the Alarm Off pressure setting.
6. Verify that the Low PEEP/CPAP alarm is OFF. (See Section 5.6, item 24.) If the alarm is ON, wait (up to 25 seconds) for the alarm to turn OFF, then press the Visual Reset button to clear the LED.
7. Decrease the syringe pressure to the alarm on setting (ignore the airway-leak alarm).
8. Verify that the low PEEP/CPAP Pressure Alarm activates within 30 seconds, and takes 20 seconds to deactivate. (See Section 5.6, item 24.)

### 5.5.21 Obstructed Tube Alarm Tests

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Crimp the tube at the test lung and verify HIGH PIP alarm. PIP should be > 40 cmH<sub>2</sub>O, but not higher than 50 cmH<sub>2</sub>O in the proximal airway digital display. *Usually over 60*
3. Remove crimp.
4. Squeeze the lung when the breath reaches peak pressure.
5. Verify that the unit can display HI PP A01, A02 (unit may not display the A02 alarm) and A05 alarms. (See Section 5.6, item 25.)

6. Kink the From Patient line hard enough to diminish flow and activate Obstructed Tube Alarm (HI PIP A03).
7. Verify that the Proximal Pressure drops abruptly when the alarm activates. Continue to kink the From Patient tube to activate a HI PP A04. (See Section 5.6, item 25.)

### 5.5.22 High CPAP Alarm Test

1. Set the following test conditions:

Mode	CPAP
PEEP/CPAP	10 cmH <sub>2</sub> O
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Turn Ventilator OFF.
3. Connect a syringe to the Proximal Airway Pressure inlet.
4. Turn Ventilator ON.
5. Increase syringe pressure until OBSTRUCTED TUBE (HI CP A04) Alarm activates.
6. Log the highest PEEP/CPAP display reading before the alarm activates (Test Limits: 15 to 25 cmH<sub>2</sub>O). (See Section 5.6, item 26.)

*Normally Alarms at 16*

### 5.5.23 Overpressure Pressure Relief Valve Test

1. Set the following test conditions:

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	12 bpm
Peak Inspiratory Pressure	90 cmH <sub>2</sub> O
High Inspiratory Pressure	95 cmH <sub>2</sub> O
Inspiratory Time	3.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz

*If pressure will not hold Reg 5 (IMV) may not be set to 62 cm H<sub>2</sub>O.*

*Check by removing outlet tubing from bottom of reg 1, put in extension tube & check low pressure on meter. To adjust remove top of reg 5 adjust to screwdriver*

- HFV Amplitude Fully clockwise 3 - 4 turns (no hard stop)
2. Turn the Relief Valve adjustment counterclockwise and note the Relief Valve setting on the PIP display. (Ignore alarms.)
  3. Increase the adjustment, and verify that it can be adjusted from 12 to 90 cmH<sub>2</sub>O. (See Section 5.6, item 27.)

#### 5.5.24 Vent Inop Test

1. Set the following test conditions:
 

Mode	IMV
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)
2. Remove the rear cover of ventilator head.
3. Disconnect J9 connector at the Analog PCA.
4. Verify that Vent Inop condition occurs. (See Section 5.6, item 28.)
5. Reconnect J9 connector.
6. Verify that the Vent Inop does not deactivate. If the Vent Inop does deactivate, the Alarm Driver PCA may have malfunctioned. (See Section 5.6, item 28.)
7. Switch the ventilator main power switch OFF and ON.
8. Verify normal conditions.
9. Disconnect J14 at the I/O PCA and repeat the procedure.

**5.5.25 Airway Leak Alarm Test**

1. Set the following test conditions:

Mode	CPAP
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
<b>Background Flow</b>	<b>8 Lpm</b>
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Disconnect the To Patient tube.
3. Turn PEEP/CPAP control clockwise to 3 cmH<sub>2</sub>O.
4. Verify that the leak alarm activates within five seconds. (See Section 5.6, item 29.)

**5.5.26 HFV Rate Control Calibration and Test (Infant Star 950)**

1. Set the following test conditions:

Mode	HFV ONLY
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
<b>Background Flow</b>	<b>4 Lpm</b>
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)

2. Using a storage oscilloscope or a frequency counter, observe the Proximal Airway pressure waveform. The measurement should be taken from the rear panel, on the Analog Output 3-pin Din connector pins (1 positive and 3 negative).
3. Remove the HFV Rate knob. While viewing the measurement device, adjust the HFV RATE shaft until the time from HFV pulse to HFV pulse equals 0.250 seconds.
4. Install the knob so that the indicator points to 4 Hz.
5. Adjust the following HFV RATE settings and record the time from HFV pulse to HFV pulse. Verify it to be within the specified limits listed below. (See Section 5.6, item 30.)

HFV Rate Setting (Hz)	HFV Rate Limits (seconds) ( $\pm 10\%$ )	Infant Star 950 HFV Rate (Hz)
10	0.090 to 0.110	9.09 to 11.11
14	0.064 to 0.079	12.66 to 15.62
20	0.045 to 0.055	18.18 to 22.22

### 5.5.27 HFV Flow Test (Infant Star 950)

**NOTE:**

This procedure is based on using the Mallinckrodt:

- high frequency breathing circuit (P/N – N-501011)
- standard test lung (P/N – N-1101262)

1. Connect a ball-type flowmeter (or equivalent) in series with the To Patient tube.
2. Set the following test conditions:

Mode	HFV ONLY
PEEP/CPAP	24 cmH <sub>2</sub> O
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully counterclockwise 3 - 4 turns (no hard stop)

3. Verify that the flow is between 6 and 13 Lpm. (See Section 5.6, item 31.)
4. Slowly turn the HFV Amplitude clockwise. Verify the flowmeter bouncing ball is normal and that the flow increases as the knob is rotated clockwise.
5. Set the HFV Amplitude fully clockwise.
6. Set the HFV Rate to 14 Hz.
7. Verify that the HFV Flow is between 34 – 43 Lpm. (See Section 5.6, item 31.)
8. Remove the flowmeter.

### 5.5.28 HFV Amplitude Test (Infant Star 950)

1. Install the HFV Breathing Circuit for this test.
2. Set the following test conditions:

Mode	HFV ONLY
PEEP/CPAP	Fully counterclockwise
Flow Rate	20 Lpm
Background Flow	4 Lpm
Ventilator Rate	30 bpm
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	1.0 second
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully counterclockwise 3 - 4 turns (no hard stop)

3. Verify that the HFV Amplitude display is between 6 - 12 cmH<sub>2</sub>O. (See Section 5.6, item 32.)
4. Set the HFV Amplitude fully clockwise and the Rate to 14 Hz.
5. Verify that the HFV Amplitude display is between 42 and 64 cmH<sub>2</sub>O. (See Section 5.6, item 32.)

### 5.5.29 HFV PEEP/CPAP Control Test (Infant Star 950)

**NOTE:**

This procedure is based on using the Mallinckrodt:

- High frequency breathing circuit (P/N – N-501011)
- Fisher & Paykel low compressible volume humidifier chamber (P/N – N-501209)
- Standard test lung (P/N – N-1101262)

**1. Set the following test conditions:**

<b>Mode</b>	<b>HFV ONLY</b>
<b>PEEP/CPAP</b>	<b>See below</b>
<b>Flow Rate</b>	<b>20 Lpm</b>
<b>Background Flow</b>	<b>4 Lpm</b>
<b>Ventilator Rate</b>	<b>30 bpm</b>
<b>Peak Inspiratory Pressure</b>	<b>40 cmH<sub>2</sub>O</b>
<b>High Inspiratory Pressure</b>	<b>45 cmH<sub>2</sub>O</b>
<b>Inspiratory Time</b>	<b>1.0 second</b>
<b>Low Inspiratory Pressure</b>	<b>3 cmH<sub>2</sub>O</b>
<b>Display Select</b>	<b>I:E Ratio</b>
<b>Oxygen Percent</b>	<b>60</b>
<b>Pressure Relief Valve</b>	<b>Fully clockwise</b>
<b>HFV Rate</b>	<b>12 Hz</b>
<b>HFV Amplitude</b>	<b>Fully counterclockwise 3 - 4 turns (no hard stop)</b>

2. Adjust PEEP/CPAP settings as listed below.
3. Verify the PEEP/CPAP setting and the PEEP/CPAP display agree for at least 5 seconds before recording the Servo Correction. (See Section 5.6, item 33.)
4. To view the Servo Correction, press and hold both the Alarm Silence and Visual Reset buttons until “:-1” appears in the Selected Data window.

PEEP/CPAP (cmH <sub>2</sub> O)	Limits (cmH <sub>2</sub> O)
3	-16 to +16
4	-16 to +16
6	-16 to +16
8	-16 to +16
10	-16 to +16
14	-16 to +16
18	-16 to +16
22	-16 to +16
24	-16 to +16

### 5.5.30 IMV + HFV Inadvertent PEEP Test

1. Set the following test conditions:

Mode	IMV + HFV
PEEP/CPAP	Fully counterclockwise
Flow Rate	See below
Background Flow	See below
Ventilator Rate	See below
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O
High Inspiratory Pressure	45 cmH <sub>2</sub> O
Inspiratory Time	See below
Low Inspiratory Pressure	3 cmH <sub>2</sub> O
Display Select	Expiration Time
Oxygen Percent	60
Pressure Relief Valve	Fully clockwise
HFV Rate	12 Hz
HFV Amplitude	Fully counterclockwise 3 - 4 turns (no hard stop)

2. Verify that the PEEP/CPAP display is within limits. (See Section 5.6, item 34.)

Flow Rate (Lpm)	Background Flow (Lpm)	Inspiratory Time (seconds)	Ventilator Rate (bpm)	PEEP/CPAP Limits (cmH <sub>2</sub> O)
40	30	0.50	60	0 to 5
30	30	0.33	90	0 to 5
20	20	0.15	120	0 to 5
10	10	0.10	150	0 to 5

### 5.5.31 Internal Battery Operation Test

The ventilator internal battery should be charged every 2 months to promote longer battery life (allow 2 hours to fully charge).

1. Charge the internal battery:
  - a. Turn the ventilator OFF.
  - b. Set the green circuit breaker switch ON.
  - c. Connect the internal battery to an outside ac source and charge for 2 hours.

## 2. Set the following test conditions:

Mode	IMV	
PEEP/CPAP	Fully counterclockwise	
Flow Rate	20 Lpm	
Background Flow	4 Lpm	
Ventilator Rate	30 bpm	<i>Windows in diagno</i>
Peak Inspiratory Pressure	40 cmH <sub>2</sub> O	<i>window 6 - 3V Ref</i>
High Inspiratory Pressure	45 cmH <sub>2</sub> O	<i>E31 - Vent INOP</i>
Inspiratory Time	1.0 second	<i>3UREF Bad +/- 10% of 3V</i>
Low Inspiratory Pressure	3 cmH <sub>2</sub> O	
Display Select	Expiration Time	<i>3UREF = 5V ÷ 256 = .01953</i>
Oxygen Percent	60	<i>154 usually seen in window in diagnostic.</i>
Pressure Relief Valve	Fully clockwise	
HFV Rate	12 Hz	<i>154 x .01953 = 3 Volts</i>
HFV Amplitude	Fully clockwise 3 - 4 turns (no hard stop)	

**NOTE:**

Once the battery is fully charged, it can be tested in the IMV mode.

3. Turn ventilator power ON, unplug the ac power cord and verify:
  - Displayed values do not change more than ±1.
  - EXT POWER LOSS LED is turned ON.
  - Ventilator delivers positive pressure breaths.
  - No active alarms.
4. The ventilator should run for a minimum of 30 minutes before it produces a low battery LED indication. Replace battery if failed. (See Section 5.6, item 35.)

**5.5.32 Hour Meter**

Note that the hour meter increased in time (tenths of an hour) during the unit checkout. (See Section 5.6, item 36.)

## 5.6 Troubleshooting

Use this troubleshooting information in conjunction with the preventive maintenance / extended tests.

1. **Ground line resistance 0.1  $\Omega$  or less failure**
  - a. Check for short circuits across or visible damage to power cord.
  - b. Verify secure connection of ground wires at terminal on weldment.
  - c. Try another ac outlet.
2. **Leakage current out of specification**
  - a. Check for damage to power cord.
  - b. Replace power supply.
3. **Unit power-up test fails**
  - a. Replace power switch. (Check all connections.)
  - b. If unit displays error codes refer to Section 9.
  - c. If unit fails upon power-up, replace MP2 PCA.
  - d. If unit completes a self test then fails, replace MP101 PCA.
  - e. Replace EPROMs.
4. **A05 test fails (high peak inspiratory pressure)**
  - a. Adjust pressure relief fully clockwise.
  - b. Replace Pressure Transducer PCA.
  - c. Replace Analog PCA.
5. **Low air pressure switch test fails**
  - a. Verify air source are within specs.
  - b. Verify wiring of pressure switch #1.
  - c. Recalibrate/replace pressure switch #1.
6. **Low oxygen pressure switch test fails**
  - a. Verify oxygen source are within specs.
  - b. Verify wiring of pressure switch #2.
  - c. Recalibrate/replace pressure switch #2.
7. **Check valve test fails**
  - Replace air/oxygen check valve.
8. **Proximal pressure test fails**
  - a. Verify functionality of manometer. Replace if necessary.
  - b. Check for static field on Front Panel. (Wipe overlay with moist towel.)
  - c. Re-zero manometer.
  - d. If analog gauge is off, recalibrate R76 on Analog PCA.
  - e. Check Proximal line for leaks.




Figure 7-6. MP2 PCA Block Diagram

### 7.1.5 MP101 PCA

The MP101 PCA (P/N – N-1101345) uses a 80631 microprocessor and its associated buss structure to perform program control of the operation of the Infant Star 500/950 ventilators in basically the same manner as the MP2 PCA.

A function of the MP101 programming is to control the operation of the Analog PCA and the decoding of the digital information obtained from it.

The programming of this PCA also provides for input and output port timing signals that are used on the I/O PCA for Front Panel Mode Selection/Display Information and output direction control of information and control signals.




Figure 7-7. MP101 PCA




Figure 7-8. MP101 Block Diagram

### 7.1.6 Pressure Transducer PCA

The Pressure Transducer PCA (P/N – N-1101163) is located in the pneumatic base assembly and provides the following outputs:

- PT2
- Prox + and Prox -
- PT1

The PT2 signal takes the proximal pressure reading of 0 - 100 cmH<sub>2</sub>O and converts it to output voltage 0 - 2.344 Vdc, which is directly proportional to the input pressure. This is accomplished by the use of a constant I source and a metering circuit (U4 and U2) that feeds the PT2 output amp to accomplish the signal conversion. The output of the PT2 gain amp is paralleled to the prox amp which converts the 0 - 100 cmH<sub>2</sub>O proximal pressure to a proportional 0 - 1 Vdc output used for oscilloscope viewing of the breathing pulses.

The PT1 signal takes the 0 - 100 cmH<sub>2</sub>O To Patient pressure and converts this signal to -0.298 Vdc to +2.126 Vdc in the same scenario as PT2.

#### Caution

To prevent water damage to the Pressure Transducer PCA, an in-line filter (P/N - N-4403017) can be added to the proximal pressure line.





Figure 7-9. Pressure Transducer PCA



S-00131

Figure 7-10. Pressure Transducer PCA Block Diagram

### 7.1.7 Front Panel PCA

The Front Panel PCA (P/N – N-1150012) assembly is divided into three sections:

- front panel controls
- mode and display select
- strobe decode and driver (includes display LEDs)

The front panel controls consist of eight separate potentiometers connected to the +5 VA supply. The output of each are:

- |                       |         |
|-----------------------|---------|
| • Ventilator Rate     | 0 - 5 V |
| • Peak Insp. Time     | 0 - 5 V |
| • Inspiratory Time    | 0 - 5 V |
| • PEEP/CPAP           | 0 - 5 V |
| • Flow Rate           | 0 - 5 V |
| • Background Flow     | 0 - 5 V |
| • High Insp. Pressure | 0 - 5 V |
| • Low Insp. Pressure  | 0 - 5 V |


These analog signals are routed to the analog PCA for conversion to digital numbers.

The Mode Select and Display consists of five individual switches that route the +5 Vdc to the MP2 and I/O PCAs. These settings are decoded and determine how the Infant *Star* 500/950 ventilators are being used.

The strobe decode takes the DB0 - DB7 buss from the MP2 PCA and decodes them into strobes 0 - 7 at the high-going edge of the DS signal to U6 and U7. The DB0 - DB7 buss is then multiplexed to block data latches 0 - 4 to provide seven segment drivers to the LEDs associated with each of the four respective blocks. When a strobe is active to a specific LED within a block that has segments decoded, these segments of the LED light cause the LED to provide the appropriate front panel indications, as required by the internal program of the Infant *Star* 500/950 ventilators. The strobe decode also routes the following signals across the PCA:

- Meter + to E1
- Meter - to E2
- Intlk 1 to Intlk 2

It also provides display information for Power Loss, Ext. Power Loss, and Vent Inop.



S-00132

Figure 7-11. Front Panel PCA




Figure 7-12. Front Panel Block Diagram

### 7.1.8 I/O (Input/Output) PCA

The I/O PCA (P/N – N-1101157) serves three functions:

- input/output buss routing
- watchdog timer circuits for MP101 and MP2 PCAs
- interlock routing of seven separate signals

The functions of the input/output are to gate information into and from AD0 - AD7 buss as follows:

- In 0 enables Vent Mode, Vent Type, Manual Breath and MP2 Hand Shake information from Front Panel to AD0 - AD7 buss.
- In 1 enables Display Select and Re-Zero information to AD0 - AD7 buss.
- In 2 enables *Star Sync* Trigger to the AD0 - AD7 buss from the front panel.
- Out 0 transfers information from the buss to the MP2 PCA through L0 - L7 (information from the A/D buss).
- Out 1 passes information from the MP101 PCA via the AD0 - AD7 buss to the Alarm Driver PCA controlling valve enabling lines (out of U7).
- Out 2 controls AD0 - AD7 information out of U12.
- Out 3 controls information out of U6 from the AD0 - AD7 buss which is used for valve control information and for gating of the MP101 watchdog timer. The watchdog timer circuits serve to provide resets to the appropriate memory PCAs when these circuits sense the loss of activity from either of these PCAs.

Interlock routing is provided for the following signals:

- Vent Inop LED
- Power Loss LED +
- Power Loss LED -
- Int. Batt LED
- Interlock 1
- PIP Interlock
- PEEP Value Interlock




Figure 7-13. I/O (Input/Output) PCA




Figure 7-14. I/O (Input/Output) PCA Block Diagram

### 7.1.9 Alarm Driver PCA (Original)


The Alarm Driver PCA (P/N – N-1101162) is divided into three operational sections:

- Alarm Monitor
- Valve Control
- Power Monitoring

The Alarm Monitor checks the status of five control signals and detects whether certain operational conditions exist. When any parameter is in the wrong state, the alarm monitor removes power from the solenoids and provides for Front Panel indication of trouble(s) by causing the Vent Inop LED to light. Once this condition has occurred, the power must be recycled to clear the Vent Inop thus restoring the 10 V SV to the solenoids.

The Valve Enable provides drive to the appropriate solenoids as required by the I/O information (U7 on I/O PCA).

Power Monitoring provides for both visual and audible indications whenever input voltage control changes occur. Visual indications appear on the Front Panel and the audible indication is provided by an internal system of buzzers.



S-00135

Figure 7-15. Alarm Driver PCA




Figure 7-16. Alarm Driver PCA Block Diagram

### 7.1.10 Differential Pressure PCA

The Differential Pressure PCA (P/N – N-1101186) controls the flow of blended gases from the oxygen blender and keeps the accumulator pressure at 38 psi, including:

- Upstream Pressure
- Downstream Pressure
- Low Flow Oxygen or Air Pressure Alarms
- Low Manifold Pressure Alarm (PS3) 8 psi

When the pressure in the Accumulator drops to approximately 28 psi, SV13 turns off, opening the pilot valve allowing flow and building up the pressure in the Accumulator. When pressure at the Accumulator increases to within 2 psi of R1 (air regulator pressure setting), SV13 energizes, closing the pilot valve.

The Differential Pressure PCA is connected to SV13, PS-3, 10 VS, 10 VA, and the Upstream and Downstream Pressure Tubing.

To measure the Upstream and Downstream pressure controls flow, the pressure is fed to the Differential Pressure PCA where a transducer (XDCR1) converts the pressure to a negative and positive voltage output. The resultant voltages are sent through two stages of amplification. This amplified output is then compared with a known Voltage Reference and is used to gate SV13 on and off. This gating can be observed by watching DS1 on the PCA which illuminates when SV13 energizes.




Figure 7-17. Differential Pressure PCA




Figure 7-18. Differential Pressure PCA Block Diagram

### 7.1.11 CE Mark Ventilators

**NOTE:**

CE marked ventilators can be identified by the CE labels affixed to the rear of the ventilator, the ferrite cores on the RS-232 cable battery, the analog output harnesses, and the capacitor between terminals 1 and 4 of Terminal Block 3 (TB3). Contact your local Mallinckrodt representative if these components are not installed.

#### 7.1.11.1 Power Monitoring

There is a green rocker switch, on the ventilator rear panel, whose function is to charge the ventilator battery. When the ventilator has power and is turned ON, and the green rocker switch is set to the ON (I) position, a light behind the switch illuminates. This indicates the internal battery is charging.

The rocker switch is equipped with an ac circuit breaker which protects the ventilator against possible shorts and line surges. When the circuit breaker trips, the rocker switch switches OFF (O position) thus interrupting ac power to the ventilator. To reset the circuit breaker, set the rocker switch back to the ON position.

If the ventilator has power and is turned ON, but the green rocker switch is turned OFF, the power supplied to the ventilator does not charge the internal battery. When this is the case, the External Power Loss LED, on the front panel display, illuminates. This indicates the ventilator is operating on the internal battery. After approximately 30 minutes of battery operation, the Low Battery light illuminates and an alarm sounds. If the ventilator continues operation without the rocker switch being reset (turned ON), the ventilator will declare a VENTILATOR INOPERATIVE condition after approximately 5 minutes.

#### 7.1.11.2 Alarm Driver PCA - Current Limiting

The Alarm Driver PCA (P/N – N-1175017) is divided into five operational sections (Figure 7-19):

- Alarm Monitor
- 12.5 V Power Supply Generation.
- Valve Control
- Power Monitoring
- Alarm Control

The Alarm Driver PCA checks the following status signals:

- VALVE ENA2
- MP2 OK
- VALVE ENA1
- MP101 OK
- PS OK
- over-current status from the OXYGEN/AIR gases blending SV13

When any of these parameters are in the wrong state, the alarm monitor will negate the SYS\_RDY and SYS\_GOOD signals. Setting the SYS\_RDY signal to a logic low state, stops the 10 V solenoid from supplying power to all solenoids. Setting the SYS\_GOOD signal to a logic low state, causes the safety valve (SV7) to open and the VENTURI INHIBIT valve (SV0) to deactivate, thus causing no venturi flow to reach the exhalation block. The Front Panel Vent Inop LED and the External and Internal audible alarms are activated simultaneously. Once


this condition occurs, power must be recycled to clear the Vent Inop alarm and thus restore power to the solenoids.

The 12.5 V power supply generates 12.5 Vdc whenever 10 V solenoid is present. The 12.5 V switching power supply is current limiting at approximately 3A and is precisely adjusting potentiometer R44 (in the voltage monitoring section) to set output voltage.

Valve Driving provides power to the appropriate solenoids as required by the I/O information (U7 on I/O PCA). All valve drivers are current limiting at 500 mA. If any valve exceeds this limit, the driver for that valve will shut off until power is recycled. All flow valves are driven by two stage drivers. At approximately 8 mS after activation, the voltage across the valve is 11 Vdc and then changed to 5.5 Vdc until the valve is deactivated. The PEEP, EXHALATION, VENTURI, and SPARE valves are driven with a constant 11 Vdc.

Power Monitoring monitors the supply of current to the heater mounted on the exhalation block and the O<sub>2</sub> blender. When the heater current exceeds approximately 0.62 A, the heater relay K1 is latched OFF until power is recycled. When the O<sub>2</sub> blender current exceeds approximately 0.6A, the SYS\_RDY signal is forced to a logic low state thus causing the 10 V solenoid to disconnect from the valves and activate the alarms as described above. It also monitors Power status: \*INT BATT and LOW BATT signals from the Power supply and I/O PCA. The corresponding LEDs are activated as a result of these signals.

Alarm Control drives the internal, external, and remote alarms with the same signal (\*BUZZ). This signal can be activated by the AUD ALM, VENT INOP2, VALVE ENA2, MP2 OK, VALVE ENA1, MP101 OK, and PS OK signals. The alarms are powered by 10 VP if the ac line voltage is present, or by the internal battery if the ac line voltage is disconnected. When there is a system power loss, the 4.8 V back-up battery provides power to the internal audible alarm. The oscillator with 0.5 seconds on-time and 2 seconds off-time intermittently drives the internal alarm.



S-00136

Figure 7-19. Alarm Driver PCA (current limiting)




Figure 7-20. Alarm Driver PCA (current limiting) Block Diagram

### 7.1.11.3 Proximal/Differential Transducer PCA

The Proximal/Differential Transducer PCA (P/N – N-1175018) is located in the pneumatic compartment base assembly and its function is to convert the proximal and internal pressure signals to electrical signals for pressure control and monitoring purposes; its other function is to control the O<sub>2</sub> blender gases. Following is a detailed description of each function of the PCA.

The Proximal and Internal Pressure Conversion part of the PCA provides the following outputs:

- PT2
- Prox + and Prox -
- PT1

The PT2 signal is converted from the proximal pressure reading of 0 - 100 cmH<sub>2</sub>O. The proportional voltage of 0 - 2.344 Vdc is accomplished by using a constant current source (CR1, CR2, U1A, and R10) and an instrument amplifier circuit (U3, U4, and U1D). The output of the differential amplifier (U3 and U4 which are part of the instrument amplifier) is connected to PT2 and Proximal op-amps which convert the differential signals to single-ended signals. The PT2 op-amp output is reference to 2 V. The output of the Proximal op-amp is reference to an analog ground. The signal is directly connected to an outside monitoring device. A high-impedance resistor network (R29, R30, and R31) divides the Proximal pressure signal down to provide internal circuit protection. The voltage produced by this circuit is 0 - 1 Vdc, which is directly proportional to 0 - 100 cmH<sub>2</sub>O.

The PT1 signal is converted from a 0 - 100 cmH<sub>2</sub>O pressure reading to a -0.298 Vdc - +2.126 Vdc electrical signal in the same manner as PT2.

The Controlling Air/Oxygen Gases Blender part of this PCA controls the flow of blended gases from the O<sub>2</sub> blender and keeps the Accumulator pressure at 38 psi including:

- Upstream Pressure
- Downstream Pressure
- Low Flow Oxygen or Air Pressure Alarms
- Low Manifold Pressure Alarm (PS3) 8 psi

When the Accumulator pressure drops to approximately 28 psi, SV13 is turned off, opening the pilot valve allowing flow and building up pressure in the Accumulator. When accumulator pressure is within 2 psi of the line Pressure, SV13 energizes and the pilot valve turns off.

The Proximal Differential Transducer PCA provides 10 V (10VSX) for SV13. The Upstream and Downstream Pressure Tubing is connected to the differential pressure transducer PT3.

Measuring the Upstream and Downstream pressure controls flow. This pressure is fed to the pressure transducer (PT3) which converts the pressure to a differential electrical signal. The resultant voltage is amplified by an instrument amplifier (U9A, U9B, U9C and associated discrete components) which amplifies and converts the signal to a single-ended signal. The output is then compared with a 3.75 V reference voltage (set up by R63 and R74) by comparator U10A. The output signal of this comparator is inverted by U10B whose output is connected directly to the SV13 driver (Q4). The driver status can be observed by watching DS1 that is lit when the solenoid energizes.




Figure 7-21. Proximal/Differential PCA Block Diagram

## 7.1.12 Pneumatic Assembly Functional Descriptions

### 7.1.12.1 IMV Regulator (R5)

R5 is gray and is located at the rear of the pneumatic base, on the bottom right side, adjacent the transducer PCA. It is calibrated to 65.2 cmH<sub>2</sub>O. During a delivered breath, the IMV regulator provides pressure to the back of the exhalation diaphragm through SV10 (IMV/PEEP select solenoid valve).

### 7.1.12.2 PEEP Regulator (R6)

R6 is located at the rear, bottom left side, of the pneumatic base. It is gray and is visibly different from the IMV regulator due to its dome port. The one-to-one regulator (a 1 cmH<sub>2</sub>O of pressure on the dome delivers 1 cmH<sub>2</sub>O through the output) is controlled by SV8 (PEEP solenoid valve). When the ventilator is set to PEEP/CPAP, the PEEP regulator applies pressure

to the back of the exhalation diaphragm when SV10 is deenergized. When a breath is *not* being delivered and there is no PEEP to apply pressure to the back of the exhalation diaphragm, the patient exhales to atmosphere.

#### 7.1.12.3 Air Regulator (R1)

R1 is located at the rear of the pneumatics base, on the top right side, near the transducer PCA. It is calibrated to 38 psi (35 psi for CE marked 500 models) and supplies 38 psi pressure to the oxygen blender, oxygen regulator dome, IMV, PEEP, mini regulators, and SV9. The R1 check valve allows flow into the ventilator, does not allow the air that has entered to return, and is connected by an attaching inlet manifold block.

#### 7.1.12.4 Oxygen Regulator (R2)

R2 is positioned at the rear of the pneumatics base on the top left side. It is calibrated to the output pressure of R1 using a differential capsuhelic gauge -5 to +5 cmH<sub>2</sub>O to approximately -0.9 cmH<sub>2</sub>O. It has a bleed on the bottom rear, left side. It provides the equivalent pressure to that delivered by R1 to the blender. If R1 drops in pressure, R2 drops equally.

#### 7.1.12.5 O<sub>2</sub> Blender

The O<sub>2</sub> Blender is located just behind the O<sub>2</sub> blender control knob located on the front of the pneumatic base. It is a long aluminum block with brass sleeves on at each end. It is calibrated to ±3 percent of the control knob setting, with the exception of 21% and 100% which are calibrated ±2%. When the O<sub>2</sub> control knob is positioned at 21, the oxygen plunger shaft (to which the knob attaches) is turned back, blocking any flow coming from the oxygen portion of the blender. When the control knob is positioned at 100, the oxygen plunger shaft is turned all the way in, blocking flow from the air side of the blender. All other settings allow a mixture of the gases producing different oxygen percentages. The blended gases are shared by the blended gas outlet, and the To Patient port.

#### 7.1.12.6 Mini Regulator (R7)

R7 is located at the front of the pneumatic base to the left of the blender assembly. It is calibrated to 13 psi to supply a 13 psi PEEP pressure to SV8 and the safety valve regulator (R8). There is a lock-down ring under the adjuster, which is pushed down and tie wrapped after calibration.

#### 7.1.12.7 Proportional Flow Manifold

The Proportional Flow Manifold is located in the front left portion of the pneumatic base and behind the To Patient port. It supplies the To Patient port and the blended gas outlet with gas from the blender. The regulator (R4) is calibrated to 18 psi and the flows are controlled by the use of a restrictor in the block under each solenoid. If the flow is 16 Lpm, a restrictor is not normally required.

#### 7.1.12.8 Jet Venturi

A lower pressure is produced inside the exhalation block assembly to reduce the inadvertent PEEP/CPAP in the patient circuit. (Inadvertent PEEP/CPAP is unwanted pressure that is higher in the patient's circuit than the PEEP/CPAP set on the ventilator.) The solenoid (SV9) duty cycle varies to control the flow to the venturi inhibit solenoid (SV0). SV0 is energized during exhalation to allow the flow to go to the jet venturi. The jet flow through the venturi creates an area of low pressure. This low pressure aids the exhalation flow and removes the inadvertent PEEP/CPAP.

7.1.13 Pneumatic Diagrams




Figure 7-22. Infant Star 500 Inspiration Pneumatic Diagram




Figure 7-23. Infant Star 500 Exhalation Pneumatic Diagram





### 7.1.14 Electrical Interconnect Diagrams




Figure 7-26. Original Infant Star Models 500/950 (1 of 4)




Figure 7-27. Original Infant Star Models 500/950 (2 of 4)







Figure 7-29. Original Infant Star Models 500/950 (4 of 4)

### 7.1.15 CE Mark Electrical Interconnect Diagrams





Figure 7-30. Infant Star Models 500/950 230 V (1 of 4)






Figure 7-32. Infant Star Models 500/950 230 V (3 of 4)

S-00154




Figure 7-33. Infant Star Models 500/950 230 V (4 of 4)

# Hardware Identification

---

## 8.1 Hardware Identification

The following assemblies and identifying components are included in this section:

- Yoke
- Front Cover
- Rear Cover
- Pneumatic Base
- Pneumatic console
- Relief Valve
- O<sub>2</sub> Blender/Flow control

---

**NOTE:**

Part numbers are provided for convenience of identifying components. They are not intended to be used to order parts, as they may not be accurate with the latest part number for a given component. When ordering parts, numbers should be verified in the latest Price List or by contacting Mallinckrodt at 1-800-635-5267, press 2 for Technical Support.

---




Figure 8-1. Yoke Assembly

**NOTE:**

The front panel ground has been moved to the bridge rectifier on latest Infant Star 500/950 ventilators

Table 8-1: Yoke Assembly Parts List

Reference	Description	Part Number
A	Alarm Driver PCA	N-1101162-SP
	Current Limiting Alarm Driver PCA	N-1175017-SP
B	Cable Assembly J-4	N-1101241
C	MP-2 PCA <i>responsible for display, shows Soft Rev.</i>	N-1101156-SP
D	I/O PCA <i>light low batt, ext power loss</i>	N-1101157-SP
E	MP101 PCA - <i>Brain of unit</i>	N-1101345-SP
F	Analog PCA	N-1175040-SP
G	Power Supply PCA	N-1101160-SP
H	Cable Assembly, J-7	N-1101239
J	Cable Assembly, J-28	N-1101235
K	100 Vac Transformer	N-1101350-SP
	Transformer Universal Assembly 115, 220 Vac	N-1150042-SP
L	Bracket, Battery Upper	N-2101307
M	Battery Assembly	N-1101230-SP
N	Bracket, Battery, Lower	N-2101306
P	Cover, Pneumatic Housing, ISV	N-2175025
R	Cable Assembly, J-10	N-1101240
T	Foam, Cover	N-2101389

*If K1 Relay opens it will not allow  
a signal to the solenoids*

*B & R connection 1A Vdc*

*E23-E30 pressure transducer to alarm board*

*E06 OR E08 - PIP window*




Figure 8-2. Front Cover Assembly

Table 8-2: Front Cover Assembly Parts List

Reference	Description	Part Number
A	Knob, .43 dia x .62L, Selection	N-2120020-2
B	Knob, .43 dia x .62L,	N-2120020-1
C	Plug, 3/16 Hold x 11/32 Head, Nylon	N-4800201
D	Overlay, Front Panel, Model 500 <i>109.28</i>	N-9910371
	Overlay, Front Panel, Model 950	N-9910764
E	Sub-Panel, ISV Front	N-2175018
F	Front Panel PCA, Model 500	N-1150012-SP
	Front Panel PCA, Model 950	N-1195012-SP
G	Gauge, 10 to 120 cmH <sub>2</sub> O	N-4300025
H	Meter Bracket	N-2150035
J	Front Cover/Bezel Assembly	N-1175015




Figure 8-3. Rear Cover Assembly

Table 8-3: Rear Cover Assembly Parts List

Reference	Description	Part Number
A	Clamp, Capacitor New Style	N-2150029
B	Capacitor, 2300 uf, 20V, 20%	N-6700006
C	CKT BKR, SPST, 8A/250V	N-7600006
D	Cable Assembly, Remote Alarm Connector J-37	N-1800010-SP
E	Hour Counter Assembly	N-1101255
F	Relay SPDT, 12Vdc, 10A, .19QCD	N-5150011
G	Cover, Rear, ISV	N-1175011
H	Circuit Breaker, DPST, 2A, 115Vac Lamp, GRN	N-7600043
	Circuit Breaker, 1.2A, 115V, Universal Transformer	N-7600060
	Circuit Breaker, .6A, 230V, Universal Transformer	N-7600059
K	Cable Assembly, J-17, RS-232	N-1101264
	Cable Assembly, J-17, RS-232 Shielded (CE model)	N-1101282
L	Switch, Rocker, DPDT, 3A @ 250Vac, Black	N-7600004
M	Connector Assembly, Wired Proximal	N-1101292

Table 8-3: Rear Cover Assembly Parts List (continued)

Reference	Description	Part Number
N	Cord, Power W/Hospital Grade Plug	N-8109001
	Cord, Universal Power	N-8109010
P	Alarm, Audio, Modified	N-2101301
R	Cable Assembly, Star Sync, Internal	N-1800077
S	Switch Seal, Push Button 3/8-32 Cover, Alarm Remote	N-7800169
T	Strain, Relief	N-7800026
U	Anchor Mount Cable Tie, 14.7 Large	N-3440026 N-3440027
V	Terminal Block, 4 POS	N-2120051
W	Terminal Block, 3 POS	N-2120052
X	Insulation Terminal Block	N-2120050
Y	Bracket, Relay, ISV (B/C)	N-2175008
Z	Torrid, Ferrite, CYLD, 196 OHM, 100MH <sub>Z</sub>	N-7800228




Figure 8-4. Pneumatics Base Assembly

Table 8-4: Pneumatics Base Assembly Parts List

Reference	Description	Part Number
A	PEEP Regulator	N-1700002-SP
B	Proportional Manifold Assembly, Model 500	N-1150022-SP
	Manifold Assembly, Proportional HFV, Model 950	N-1121275-SP
C	Double Spade Disconnect Adapter, Male	N-7700137
D	SV7 Solenoid	N-1101229-SP
E	PEEP Restrictor and Tubing Assembly	N-1101228
F	Foam Accumulator	N-2101288
	O <sub>2</sub> Regulator Assembly, 1/8-inch NPT, 35 psi (CE model and overhauled since January 1999)	N1195006-SP
G	O <sub>2</sub> Regulator Assembly	N-1101331-SP
	O <sub>2</sub> Regulator Assembly, 1/4-inch NPT, 38 psi (35 psi CE model 500)	N-1150068-SP
H	1/4" to 1/4" NPT Bushing Chrome Nipple Fitting	N-3600151
	1/8-inch NPT Chrome Nipple (CE model and overhauled since January 1999)	N-3600174
J	Chrome 1/4 x 1/4" NPT Elbow Fitting	N-3600150
	1/8-inch NPT 90° Elbow Fitting (CE model and overhauled since January 1999)	N-3600175
K	O <sub>2</sub> Body CGA 1240 x 1/4" Chrome	N-3600153
	O <sub>2</sub> Body Adapter dCGA 1240 x 1/8 inch (CE model and overhauled since January 1999)	N-3600176
L	Cable Tie 14.5 LG .19 Wide NYL	N-3440009
M	Accumulator Bundle Assembly	N-1101215
	Bundle Assembly	N-1121215
N	Air DISS FTG Body Adapter, 1/4" NPT Chrome	N-3600152
P	DIS Demand Check Valve 1/4" NPT O <sub>2</sub> Fitting Chrome	N-3600154
R	Cable Assembly, Pres. Sev (1,2,3)	N-1101303
S	XDCR PCA (non-CE mark ventilator)	N-1101163-SP
	Proximal/Differential Transducer PCA (CE model ventilator)	N-1175018-SP
T	IMV Regulator Assembly	N-1700056-SP
U	Sleeve, Pole Mount	N-2100112
V	Prt/restrictor Assembly Bleed	N-11012667
W	CBI tie 3.9 LG 0.10 Wide NYL	N-3440008
X	Muffler/Bleed part Assembly	N-1101268-SP

Table 8-4: Pneumatics Base Assembly Parts List (continued)

Reference	Description	Part Number
Y	Bracket Assembly SV-8, SV-9 solenoids, Model 500	N-1101227
	Bracket Assembly SV-8, SV-9 Solenoids, Model 950	N-1121227
AA	Mini REgulator Assembly	N-1101253-SP
AB	Inhibit Assembly, Venturi SV-0 Solenoid	N-1101258-SP
AC	Base, Pneumatic, Silk-screen	N-2190003-21
AD	Foot, Rubber, 1-18 D x 7/8 Tall Black	N-4800204
AE	Bracket, Mounting Terminal Block	N-2101338
AG	Fitting, barbed, 1/8 Tubing - 1/8 NPT	N-3610071
AH	Fitting, ELL, 1/8 NPT x 10-32	N-3600182
AJ	Muffler, 10-32, Sintered	N-4800118
AK	Knob, Thumb Screw, Knurled, N08	N-5400004
AM	Filter Assembly 1/4" NPT 5 Micron	N-4402005
	Filter Assembly 1/8" NPT 5 Micron (CE model and overhauled since January 1999)	N-4402068
	Filter Element, 5 Micron	N-4403006
AN	Air Regulator Assembly	N-1101202-SP
	Air Regulator Crossover Assembly (CE models only)	N-1101529-SP
AP	Blender Gas Outlet Assembly	N-1101221-SP
AR	O <sub>2</sub> Blender/Snap Assembly	N-1220068-SP
	O <sub>2</sub> Blender Assembly (CE models only)	N-1195007-SP
AS	Manifold Assembly, Proportional HFV, Model 950	N-1121275-SP
AT	Key Lock Assembly, Model 950	N-1195019-SP
	Key Lock Assembly (950 CE model only)	N-1800059-SP
AU	Replacement Key HFV	N-7600068




Figure 8-5. Pneumatic Console Assembly

Table 8-5: Pneumatic Console Assembly Parts List

Reference	Description	Part Number
A	Venturi Nozzle Assembly Model 500	N-1150016-SP
	Venturi Nozzle Assembly Model 950 or CE model 500	N-1195016-SP
B	Clamp Assembly, Venturi Block	N-2175013
C	Knob, Clamp, ISV	N-2175038
E	FTG, Short Coupling, #10-32, BRS	N-3600141
F	Dowel Pin, .125 DIA x .750L	N-3400062
G	Pivot Block, ISV	N-2175014
H	Screw, SKTHD, Set 6-32 x 3/16	N-3100511
J	Support, Heater Block	N-2175010
K	Block Assembly, Heater	N-1150019
L	Valve Assembly, PEEP/IMV, SV-10	N-1150021-SP
M	Dowel Pin, .125 dia x 1.50 L, SS	N-3400063
N	FTG, Barb, 90 degree 1/8 NPTF x 3/1	N-3610054

Table 8-5: Pneumatic Console Assembly Parts List (continued)

Reference	Description	Part Number
P	Bracket, Base Support, Right	N-2175023-2
S	Screw, Panhead, Phil, 8-32 x 3/8	N-3100036
T	Bracket, Base Support, Left	N-2175023-1
U	Plate, Heating, ISV	N-2175007
V	Thermistor, open @ 140/close @ 120	N-5500008
W	Resistor, Power, 30 ohms, 10W, 1%	N-6800053
X	Venturi Block Assembly, Threaded	N-1150024-SP
Y	Valve Assembly, Relief	N-1150009-SP
AA	Screw, Phil, Panhead, #6-32 x 3/8	N-3100026
AB	Retainer, Pivot Block ISV	N-2175015





Figure 8-6. Relief Valve Assembly

Table 8-6: Relief Valve Assembly Parts List

Reference	Description	Part Number
B	Collar, Shaft	N-2175052
C	Washer, Nylon, .505 ID x 1.25 OD	N-3220154
D	Washer, Wave, .578 ID x .995 OD	N-3220155
F	Bracket, Relief Valve MTC	N-2175027
H	Cam, Relief Valve, ISV	N-2175031
K	Spring, Relief Valve	N-2175053
L	Ball, 3/4 dia. Delrin	N-4800182
M	Retaining Ring, .25 OD, .219 DX .029 GRV	N-3700045
N	Shaft, Relief Valve	N-2175033
P	Valve Block, Relief	N-2175032
R	Fitting, Barb NYL, 1/8 NPT - 1/8 Hose	N-3610020

Table 8-6: Relief Valve Assembly Parts List (continued)

Reference	Description	Part Number
S	Fitting, Street-T, 1/8 NPT x 1/8 NPTF (2)	N-3600122
U	Washer Flat, #6 .156 ID x .375 OD	N-3220006
V	Washer, Split Lock, #6	N-3220032
W	Screws, Phil Panhead #6-32 x 3/8	N-3100026
X	Plunger, Relief Valve Cam	N-2175030
Y	Bushing, Flanged, .250 x .295 D x .25, NYL	N-4800195
AA	Fitting, 10 and 15 mm Barb	N-2175029

Figure 8-7. O<sub>2</sub> Blender AssemblyTable 8-7: O<sub>2</sub> Blender Assembly Parts List

Reference	Description	Part Number
A	Blender Assembly, O <sub>2</sub> , Model 500	N-1220069-SP
B	Housing, Restrictor	N-2101278
C	Fitting, Hoser 10-32 to 1/8 Barb Assembly	N-3600129
D	Restrictor, Jeweled Typer .015"	N-3540015
E	Fitting, Short Coupling, #10-32	N-3600141
F	Fitting O <sub>2</sub> Bleed	N-2101313
G	Pot, Single Turn	N-6300009
H	Bracket, Potentiometer	N-2101221
J	Flow Control Assembly, Blender	N-1101301-SP
K	Differential Pressure PCA	N-1101186-SP
L	Bumper, .44 DIA x .20 High	N-4800042
M	Insulation, Differential	N-2101365
N	Button, Mounting	N-7800024
P	Gear, Sput, Molded	N-4240001

# Alarms and Troubleshooting

## 9.1 Software Diagnosed Alarms

The ventilator microprocessor can diagnose problems with the functioning of the ventilator. When one of these conditions is identified, the phrase "Ventilator Inoperative" and an error code are seen in the "Selected Data" display window. A list of these error codes, the problem they identify and the probable cause for each appears in Table 9-1, Ventilator Inoperative Alarm Troubleshooting Chart.

**Table 9-1: "Ventilator Inoperative" Alarm Troubleshooting Chart**

Error Code	Problem	Probable Cause
None	10 V battery fully discharged fuse blown or battery damaged. Accompanied by a POWER LOSS LED. Power Supply or MP2 PCA failures.	10 V battery (if during internal Battery operation) MP2 PCA. Power Supply Assembly.
E01	Unable to perform automatic rezeroing of Proximal Pressure Transducer on power up or after 24 hours of continuous operation.	Pressure Transducer PCA. (May require rezero of Proximal Pressure Transducer.)
E02	Data error (checksum) error found in data from MP101. Two errors detected within 10 msec.	MP101 PCA I/O PCA MP2 PCA
E03	Internal checks	Position 4 of SW1 on MP2 PCA is not in open (off) position or MP2 PCA failures.
E04	Internal checks	MP101 PCA failure
E05	Error found in data sent to MP2 (PIP display shows index number of erroneous data).	MP101 PCA I/O PCA
E06, E07	Control panel potentiometer (Pot) out of range E06 = low, E07 = high PIP shows POT NUMBER Pot No.    Pot Name  0            Flow Rate 1            Ventilator Rate 2            PEEP/CPAP 3            Inspiratory Time 4            Low Inspiratory Pressure 5            Peak Inspiratory Pressure 6            Not used 7            Not used 8            HFV Rate 9            HFV Amplitude	Control Panel PCA Analog PCA

*Check J11 cable  
also look at pins of  
pot or solder splashes*

Table 9-1: "Ventilator Inoperative" Alarm Troubleshooting Chart (continued)

Error Code	Problem	Probable Cause
E08	MP2 waited too long for data from MP101.	MP101 PCA MP2 PCA
E09	Two consecutive A/D errors.	Analog PCA
E10	Exhalation valve (SV 10) closed for longer than 3.5 seconds as indicated by EXHAL (J5X-8) on Alarm/Driver PCA staying low for longer than 3.5 seconds. E10 is also displayed 3.5 seconds after a VENT INOP occurs when no other error code has been previously displayed.	Alarm Driver PCA Power Supply Assembly Various cable assemblies MP101 PCA MP2 PCA I/O PCA
E11	No flow greater than 4 Lpm for 66 seconds indicated by the FLOW signal on Alarm/Driver PCA staying low for more than 66 seconds.	Alarm Driver PCA The PIP was held higher than the set PIP by more than 15 cmH <sub>2</sub> O for 66 seconds.
E12	No HFV breath in 4 seconds as indicated by FLOW on Alarm/Driver PCA staying low for longer than 4 seconds in HFV mode.	Alarm Driver PCA The PIP was held higher than the set PIP by more than 15 cmH <sub>2</sub> O for 4 seconds during an obstructed tube alarm.
E13	MP2 failed to acknowledge MP101 data transmission.	MP2 PCA
E14	MP101 Startup RAM check error.	MP101 PCA (RAM IC)
E15	MP101 Background RAM check error.	MP101 PCA (RAM IC)
E16	MP101 ROM check error.	MP101 PCA (ROM IC)
E17	Not used	
E18	MP2 Startup RAM check error.	MP2 PCA (RAM IC)
E19	MP2 Background RAM check error.	MP2 PCA (RAM IC)
E20	MP2 Background ROM check error.	MP2 PCA (ROM IC)
E21	MP2 Background ROM check error.	MP2 PCA (ROM IC)
E22	MP2 detects error in its own software interrupt routine.	MP101 PCA
E23	Zero offset of 0 to 100 cmH <sub>2</sub> O A/D converter output >5 cmH <sub>2</sub> O.	Pressure transducer PCA. (May require rezero of Proximal Pressure transducer.)
E24	Zero offset of 0 to 100 cmH <sub>2</sub> O A/D converter output <-5 cmH <sub>2</sub> O.	
E25	Zero offset of 0 to -50 cmH <sub>2</sub> O A/D converter output >5 cmH <sub>2</sub> O	
E26	Zero offset of 0 to -50 cmH <sub>2</sub> O A/D converter output <-5 cmH <sub>2</sub> O.	
E27	Zero offset of -5 to +50 cmH <sub>2</sub> O A/D converter output >5 cmH <sub>2</sub> O	
E28	Zero offset of -5 to +50 cmH <sub>2</sub> O A/D converter output <-5 cmH <sub>2</sub> O	Pressure transducer PCA. (May require rezero of Proximal Pressure transducer.)
E29	Zero offset of METER VOLTS A/D converter output >6 cmH <sub>2</sub> O.	
E30	Zero offset of METER VOLTS.	

E23  
 E24  
 E25  
 E26  
 E27  
 E28  
 E29  
 E30  
 All  
 these  
 errors  
 are  
 caused  
 by  
 tubing

Table 9-1: "Ventilator Inoperative" Alarm Troubleshooting Chart (continued)

Error Code	Problem	Probable Cause
E31	Internal 3 V battery reference varies greater than 10%.	Power Supply PCA Analog PCA
E32	Front Panel Identification is NOT valid.	Front Panel is incompatible with installed software or identification is intermittent
E33	Delivery mode is NOT valid.	Ventilator delivery mode is not a valid selection
E34	Dead Battery.	10.8 V dc Battery Pack is bad.

## 9.2 Infant Star Troubleshooting Guide

This troubleshooting guide should be used with a Mallinckrodt Circuit P/N – N-501011, and Test Lung P/N – N-1101262.

The *Infant Star Service Manual* has the references for standard test conditions and regular testing equipment. Always check the external parts of the *Infant Star* for a problem before opening it for repair.

The following guides address conditions that may or may not occur with your ventilator. These guides state problems that may be encountered, they list a probable cause(s) for the identified problems, and they give suggestions for troubleshooting and correcting these conditions.

### 9.2.1 IMV Standard Test

Problem	PIP set to 40 cmH <sub>2</sub> O, but only reaches 38.
Possible Causes	<ol style="list-style-type: none"> <li>1. SV1 not energized or flows from 2 - 4 Lpm are low.</li> <li>2. I time was set to 0.1 but should be 1.0.</li> <li>3. Small hole in the diaphragm.</li> <li>4. Pop-off valve not seating, not screwed fully clockwise.</li> <li>5. Pop-off valve damaged or debris on ball or seat.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Check flows; if low, clean restrictor in the proportional manifold.</li> <li>2. Set I time to 1.0.</li> <li>3. Replace diaphragm.</li> <li>4. Turn pop-off valve fully clockwise.</li> <li>5. Clean or replace ball or seat.</li> </ol>
Problem	Display settings drift more than 2 increments.
Suggested Actions	Replace F/P PCA or Analog PCA.
Problem	PIP set to 10 cmH <sub>2</sub> O, but reaches 12 or higher.
Probable Causes	<ol style="list-style-type: none"> <li>1. Diaphragm installed backward.</li> <li>2. Leak in Flow Manifold - caused by very dirty or broken flow solenoid spider.</li> <li>3. Film on diaphragm making diaphragm stick or bad diaphragm.</li> <li>4. Diaphragm seat on venturi block not completely seated.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Reinstall diaphragm.</li> <li>2. Clean or replace solenoid.</li> <li>3. Clean or replace diaphragm.</li> <li>4. Replace venturi block.</li> </ol>
Problem	PIP set to 90 cmH <sub>2</sub> O, but only reaches 87.
Possible Causes	<ol style="list-style-type: none"> <li>1. IMV Regulator out of calibration.</li> <li>2. Pop-off valve leaking.</li> <li>3. Circuit on test lung leaking.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Calibrate IMV Regulator to 60 cmH<sub>2</sub>O.</li> <li>2. Clean, check the adjuster and spring and ball; replace as necessary.</li> </ol>

9.2.2 Obstructed Tube Alarms A01 Test, A01, A02, A05

Problem	A01 Alarm 5 cmH <sub>2</sub> O over PIP (set by control panel -- HIP)
Possible Causes	<ol style="list-style-type: none"> <li>1. Too short a circuit, and/or Proximal line.</li> <li>2. Possible leaky flow solenoid.</li> <li>3. SVO not closing.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Locate a proper circuit and test lung.</li> <li>2. Clean flow solenoids; check for broken spiders. Replace if necessary. If the problem still exists test for leaky solenoid:                      Turn off the ventilator.                      Connect air and oxygen.                      Plug SV7 safety vent valve hole (black, threaded cylinder under flow manifold facing toward outside of ventilator).                      Connect tube from To Patient port into cup of water. Bubbles indicate a leak coming from flow manifold. Clean or replace solenoids as necessary.</li> <li>3. Check electrical connection.</li> </ol>
Problem	A05 Alarm 10 cmH <sub>2</sub> O over PIP
Possible Causes	<ol style="list-style-type: none"> <li>1. Blocked inspiratory tube on patient circuit or To Patient port is plugged.</li> <li>2. Water in Proximal line or in proximal accumulator inside ventilator.</li> <li>3. Needle valve out of calibration.</li> <li>4. Large leak in flow solenoid valve.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Clear blockage.</li> <li>2. Disconnect tubing on both sides of accumulator and blow out with an air nozzle. Install proximal line filter P/N – N-4403017 water from entering ventilator.</li> <li>3. Calibrate needle valve.</li> <li>4. Check flow valves for broken spider.</li> </ol>
Problem	No A05 Alarm
Probable Causes	Faulty Transducer PCA
Suggested Actions	To test A05 Alarm: plug patient port. Replace Transducer PCA.

### 9.2.3 PEEP 0 to 24 cmH<sub>2</sub>O

Check PEEP calibration before continuing. If necessary, clean SV8.

Problem	PEEP set to 5 cmH <sub>2</sub> O but reaches 6 or higher
Possible Causes	<ol style="list-style-type: none"> <li>1. PEEP calibration may be incorrect.</li> <li>2. PEEP regulator cracking pressure too high.</li> <li>3. Diaphragm installed backwards.</li> <li>4. Low jet venturi operation.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Check or calibrate PEEP/CPAP needle valve adjustment.</li> <li>2a. Inspect 0.024 restrictor (downstream from SV8), clean if necessary.</li> <li>2b. Check R7 setting.</li> <li>3. Reinstall diaphragm or replace diaphragm.</li> <li>4a. Check or calibrate venturi needle valve adjustment, if installed.</li> <li>4b. Check R7 setting.</li> </ol>
Problem	PEEP set to 5 cmH <sub>2</sub> O but only reaches 4 or less
Possible Causes	<ol style="list-style-type: none"> <li>1. Diaphragm not seated properly.</li> <li>2. Damaged jet venturi O-ring.</li> <li>3. Cap on PEEP regulator loose, cracked, or broken.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1a. Check or calibrate PEEP/CPAP needle valve adjustment.</li> <li>1b. Reseat diaphragm.</li> <li>2a. Check O-ring between jet and venturi body.</li> <li>2b. Inspect venturi body for damage.</li> <li>3. Tighten regulator cap, if required or replace PEEP regulator (R6).</li> </ol>

### 9.2.4 PEEP Circuit

Problem	Cannot calibrate a PEEP of 20 cmH <sub>2</sub> O (Without any pressure from the PEEP circuit, a setting of 20 to 24 cmH <sub>2</sub> O will result in a 13 to 15 cmH <sub>2</sub> O PEEP.)
Possible Causes	<ol style="list-style-type: none"><li>1. Leak in internal PEEP circuit</li><li>2. SV8 not energizing.</li><li>3. No air flow through SV8.</li><li>4. Pop-off valve not fully closed.</li><li>5. Faulty internal PEEP circuit.</li><li>6. Faulty needle valve.</li><li>7. Faulty SV8.</li><li>8. Flow does not reach 8 Lpm.</li></ol>
Suggested Actions	<ol style="list-style-type: none"><li>1. Check all tube connections.</li><li>2. Check electrical connections.</li></ol>

**Caution**

Do not break the needle valve. Disconnect SV8 tubing only. Do not pull off needle valve.

3. Pull needle valve off SV8 and listen for air pulses.
4. Turn fully clockwise.
5. Check circuit for broken white elbow.
6. Clean or replace needle valve.
7. Clean or replace SV8.
8. Correct flow setting.

### 9.2.5 Venturi Solenoid (SV9)

Controls the baseline pressure when the PEEP circuit is not in use. Venturi and PEEP never operate at the same time.

---

**NOTE:**

See venturi calibration for newer Infant *Star* models 500 and 950.

---

Problem	A04 Alarm
Possible Causes	<ol style="list-style-type: none"> <li>1. Jet Venturi nozzle plugged.</li> <li>2. Pressure going to back of diaphragm.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Clean nozzle. Using an eyeloop, check nozzle's orifice for a crack; replace if necessary.</li> <li>2a. First, disconnect tube behind diaphragm, if problem corrects reconnect tube.</li> <li>2b. Second, disconnect the PEEP regulator output, if problem corrects check PEEP circuit.</li> <li>2c. Third, disconnect the IMV regulator output, if problem corrects check SV10 IMV select valve.</li> </ol>
Problem	Unstable baseline
Possible Causes	<ol style="list-style-type: none"> <li>1. Dirty SV9</li> <li>2. Faulty SV9</li> <li>3. Inlet pressures are 50 psi or higher.</li> </ol>
Suggested Actions	<ol style="list-style-type: none"> <li>1. Clean SV9.</li> <li>2. Replace SV9.</li> <li>3. Check static or 45 psi or above dynamic.</li> </ol>