1. Einleitung

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BiPAP 30 System Adjustment Data Sheet (Pre Run-In)

7.4.1 Blower Motor Speed

Hz	4780-4820 Hz
----	--------------

7.4.2 Flow Transducer Calibration

Flow Meter L/min	0	60
\dot{V}_{tot}		
Range (volts)	-0.015 to 0.014	0.985 to 1.014

7.4.3 Valve Calibration

Control Knob Setting	28
IPAP	
EPAP	
Pressure (cm H ₂ O)	27.7-28.3

7.4.5 Pressure Monitor Calibration (S/T-D 30 Only)

Manometer Pressure Settings	Measurement (volts)	Range (volts)
28.0		2.65 to 2.94

7.4.7 Trigger Level Verification

Time Meter Reading (End)

Pass	Fail
(Circle One)	

Serial Number _____ Signature _____ Date _____ PASS FAIL

BiPAP 30 System Final Data Sheet (Post Run-In)

7.4.1 Blower Motor Speed

Hz	4780-4820 Hz
----	--------------

7.4.2 Flow Transducer Calibration

Flow Meter L/min	0	30	60	90	120
\dot{V}_{tot}					
Range (volts)	-0.035 to 0.034	0.350 to 0.640	0.850 to 1.140	1.300 to 1.690	1.750 to 2.240

7.4.3 Valve Calibration

Control Knob Setting	Min	6	12	18	24	30
IPAP						
EPAP						
Pressure (cm H ₂ O)	<2.5	3.5 to 8.4	9.5 to 14.4	15.5 to 20.4	21.5 to 26.4	27.5 to 32.4

7.4.4 Pressure vs. Flow Calibration

Flow Valve Settings	50	100	
Manometer Pressure Settings			Pressure (cm H ₂ O)
6.0			4.5 to 7.5
12.0			10.5 to 13.5
18.0			16.5 to 19.5
24.0			22.5 to 25.5
29.0			27.5 to 30.5

7.4.5 Pressure Monitor Calibration (S/T-D 30 Only)

Manometer Pressure Settings	Measurement (volts)	Range (volts)
0.0		-0.15 to 0.14
5.0		0.35 to 0.64
10.0		0.85 to 1.14
15.0		1.35 to 1.64
20.0		1.85 to 2.14
25.0		2.35 to 2.64
28.0		2.65 to 2.94

7.4.6 Blocked and Unblocked Watts

WATTS			
Blocked		Unblocked	
	<110		<110

7.4.7 Trigger Level Verification

Pass	Fail
(Circle One)	

7.6.3 & 7.6.4 Spontaneous & Spontaneous/Timed

	Spont.	Spont/Timed
S/T		Pass/Fail
S/T-D		Pass/Fail
	4.40 to 5.95 sec	Circle One

7.6.5 Timed Mode Verification

	BPM min	%IPAP min	BPM max	%IPAP max
S/T				
S/T-D				
Range	10 to 30 sec	5 to 15%	1.8 to 2.2 sec	85 to 95%

7.6.6 DCP Verification (S/T-D 30 Only)

Pass	Fail
(Circle One)	

Time Meter Reading (End)

Serial Number _____ Signature _____ Date _____ PASS FAIL

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**Clinical
Manual**

Table of Contents

Chapter 1:	General Description	1
Chapter 2:	Warnings, Cautions, and Notes	3
Chapter 3:	Description of Controls	9
Chapter 4:	Modes of Operation	15
Chapter 5:	Setting Up the BiPAP 30 System	23
Chapter 6:	Administering Oxygen with the BiPAP 30 System	27
Chapter 7:	Administering Aerosolized Medication with the BiPAP 30 System	29
Chapter 8:	Performance Verification.....	31
Chapter 9:	Detachable Control Panel 30	33
Chapter 10:	Cleaning Instructions	37
Chapter 11:	Routine Maintenance.....	39
Chapter 12:	Troubleshooting.....	41
Chapter 13:	Abbreviations/Definitions	43
Chapter 14:	Specifications	47
Chapter 15:	BiPAP S/T-D 30 Outputs	49
	Limited Warranty.....	51

Chapter 1: General Description

The BiPAP 30 System Models S/T 30 and S/T-D 30

The BiPAP 30 System is a positive pressure, electrically driven ventilation system with electronic pressure control. The unit's pressure controls are adjusted to deliver pressure support for patient ventilatory assistance.

The BiPAP 30 System is primarily intended to augment patient ventilation by supplying pressurized air through a patient circuit. It senses the patient's breathing effort by monitoring airflow in the patient circuit and adjusts its output to assist in inhalation or exhalation. This assistance is provided by the administration of two levels of positive pressure. During exhalation, pressure is variably positive or near ambient. The inspiratory level is variably positive and is always higher than the expiratory level.

The Respironics BiPAP 30 System responds to changes in patient flow rates that indicate movement to inhalation and exhalation. It performs reliably in sensing the patient's breathing efforts, even with the presence of most leaks in the patient circuit. Automatic adjustment of the trigger threshold in the presence of leaks is what makes the BiPAP 30 System ideal for mask-applied ventilation assistance.

The BiPAP 30 System senses air flow in the patient circuit by using a flow transducer in conjunction with the patient air outlet. The flow signal is composed of flow into the patient's lungs and flow due to leaks in the patient circuit. This flow data is continuously processed, and adjustments are automatically made to trigger thresholds. The adjustments provide high sensitivity to inspiratory effort and rapid adjustment for changing leaks while minimizing auto-triggering.

The BiPAP 30 System can operate in the following four modes. In each mode, the prescribed setting for each functional control must be determined by the physician/clinician based on appropriate patient testing and monitoring of the necessary physiologic parameters.

SPONTANEOUS (S) MODE

The unit cycles between the Inspiratory Positive Airway Pressure (IPAP) and Expiratory Positive Airway Pressure (EPAP) levels in response to patient triggering. The patient determines the respiratory rate.

SPONTANEOUS/TIMED (S/T) MODE

The unit cycles between the IPAP and EPAP levels in response to patient triggering. If the patient fails to initiate an inspiration, the unit will cycle to IPAP based on a preset interval determined by the Breaths Per Minute (BPM) control.

TIMED (T) MODE

The unit cycles between the IPAP and EPAP levels based solely on the timing intervals as determined by the rate (BPM) and inspiratory time (% IPAP) controls. The patient may superimpose spontaneous breathing over the IPAP and EPAP levels.

CONTINUOUS POSITIVE AIRWAY PRESSURE (CPAP)

With the Function Selector Knob set in either the IPAP or EPAP position, the pressure set on the corresponding dial will be delivered continuously.

Chapter 2: Warnings, Cautions, and Notes

WARNING:	Indicates the possibility of injury to the patient or the operator
CAUTION:	Indicates the possibility of damage to the device
NOTE:	Places emphasis on an operating characteristic

2.1 WARNINGS

General

- Federal law restricts these devices to sale by or on the order of a physician.
- This manual serves as a reference. It should be used in conjunction with the instructions and protocol set by the physician at the institution where the device is being used. The instructions in this manual are not intended to supersede established medical protocols regarding the use of assist ventilators.

Safety

- The BiPAP 30 System is an assist ventilator and is intended to augment the ventilation of a spontaneously breathing patient. **IT MUST NOT BE USED AS A LIFE SUPPORT VENTILATOR.** It is not intended to provide the total ventilatory requirements of the patient.
- The BiPAP 30 System should not be applied to patients who are incapable of maintaining life-sustaining ventilation in the event of malposition of the patient interface.
- The BiPAP 30 System is intended for use on pediatric or adult patients with:
 - chronic respiratory failure;
 - acute respiratory insufficiency;
 - acute respiratory failure; and,
 - obstructive sleep apnea syndrome.
- The BiPAP 30 System is intended for use with a Respironics, Inc. approved patient circuit only. A Respironics approved circuit consists of :
 - smooth inner lumen tubing;
 - an exhalation device such as the Whisper Swivel[®] or Respironics Disposable Exhalation Port Assembly; and,
 - a noninvasive patient interface (nasal mask, mouthpiece, or full face mask).

Additional accessories may be added to the circuit to meet specific needs. Every time changes are made to the circuit configuration, the delivered pressures must be monitored at the patient connection with the unit cycling. If using other patient circuits, they must have a continuous leak exhalation port with flow comparable to Respironics exhalation ports.

- The continuous flow exhalation port (i.e., Whisper Swivel) is designed to exhaust patient exhalation from the patient circuit. Continuous flow is required for safe operation. **DO NOT BLOCK OR OTHERWISE TRY TO SEAL THE GAS VENTS ON THE EXHALATION PORT.**

WARNINGS continued

Service

- Wiring connections and adjustments should be performed by authorized service personnel ONLY.
- Repairs and adjustments MUST be performed by experienced personnel, fully acquainted with this equipment. Service done by inexperienced, unqualified personnel or installation of unauthorized parts could cause personal injury, invalidate the warranty, or result in costly damage. Contact Respironics, Inc. or your health care dealer if the unit fails to meet performance specifications.
- To avoid electrical shock, disconnect the electrical supply before changing the fuses.
- Replace fuses with those of the same type and rating only.

Cleaning

- To avoid electrical shock, unplug the BiPAP 30 unit before cleaning it.

2.2 Cautions

- The DCP 30 cable connector port is NOT a standard computer communications port. Connecting the output port to a computer could damage the computer and the BiPAP 30 unit.
 - On BiPAP S/T-D 30 models, output signal pin connections are active. Care must be exercised to prevent electrical shorting between the pin connections.
 - Only Respironics cables should be used. Use of any other cable types may cause damage to the BiPAP S/T-D 30 System.
 - For patient pressure monitoring, from an external monitor, use only tubing with a smooth inner lumen between the patient circuit and monitor.
 - When transporting or storing the BiPAP 30 System, precautions should be taken to avoid exposure to temperature extremes. If exposure to such temperatures has occurred, the unit should be allowed to come to room temperature before being turned on.
 - The unit should be positioned on its base or rear panel for proper operation.
 - The voltage selector switch must be set according to local line voltage, and the proper fuses used.
 - When using Respironics smooth inner lumen tubing without a heated humidifier, the overall length of the tubing should not exceed twelve feet.
 - Regularly inspect the DCP 30 cable.
 - Remove the mask before turning the BiPAP 30 unit off.
-

2.4 Contraindications

The following conditions may contraindicate the use of a BiPAP 30 System:

- a. Patients with or susceptible to pneumothorax or pneumomediastinum should be monitored closely when applying positive pressure. Pre-existing bullous lung disease may represent a relative contraindication.
- b. Hypotension induced by positive pressure ventilation.
- c. A history of allergy or hypersensitivity to the mask material where the risk from allergic reaction outweighs the benefit of ventilatory assistance.
- d. Aspiration of gastric contents. Exercise caution when applying a full face (nasal/oral) mask to patients because of the increased possibility of aspirating gastric contents. This is particularly true in patients who are at risk for or who have been vomiting. In these patients, placement of a nasogastric suction tube may be advisable prior to the application of a full face mask.
- e. Acute sinusitis or otitis media. Although it is not necessary in many patients, nasal positive pressure therapy may need to be temporarily suspended in some patients with acute sinusitis or otitis media.
- f. Patient's inability to maintain a patent airway or to adequately clear secretions.
- g. Life-threatening or potentially life-threatening epistaxis.
- h. Patient's without an intact respiratory drive.



2.5 Patient Warnings



- Advise the patient to immediately report any unusual chest discomfort, shortness of breath, or severe headache.
- If using a full face mask, the patient should be advised not to eat or drink two to three hours prior to bedtime.
- If skin irritation or breakdown develops due to the mask, refer to the User Circuit Accessory Guide for appropriate action.
- The following are potential side effects of noninvasive positive pressure ventilation:

- Ear discomfort
- Conjunctivitis
- Skin abrasions due to noninvasive interfaces
- Aerophagia (Gastric Distention)

Chapter 3: Description of Controls

BiPAP 30 System Models S/T 30 and S/T-D 30

FRONT PANEL	
 <p>BiPAP[®] VENTILATORY SUPPORT SYSTEM S/T-D 30 RESPIRONICS INC.[®]</p>	
	<p>POWER SWITCH</p> <p>Turns the electrical power to the BiPAP 30 unit on and off.</p>
	<p>TUBING CONNECTOR</p> <p>Accepts 22 mm (approximately 7/8") I.D. smooth inner lumen tubing and also accepts a mainflow bacteria filter.</p>
	<p>PRESSURE PORT</p> <p>S/T-D 30 units only: Accepts manometer tubing with a smooth inner lumen. Interfaces with an internal pressure transducer to provide circuit pressure monitoring with output to a recorder. Source of the pressure display for the Detachable Control Panel 30 (DCP 30). The range is 0 to 30 cm H₂O.</p> <p>CAUTION: Never connect oxygen or any positive pressure source to this port.</p>

CONTROLS	
	<p data-bbox="458 374 752 408">IPAP CONTROL</p> <p data-bbox="458 461 900 495">Inspiratory Positive Airway Pressure</p> <p data-bbox="458 525 1328 587">Active when the Function Selector Knob is in any of these four positions: IPAP, S, S/T, or T.</p> <p data-bbox="458 622 1370 656">When the BiPAP 30 unit is delivering IPAP, the adjacent LED is illuminated.</p> <p data-bbox="458 686 749 720">Range = 4 to 30 cm H₂O.</p>
	<p data-bbox="458 759 757 794">EPAP CONTROL</p> <p data-bbox="458 847 893 881">Expiratory Positive Airway Pressure</p> <p data-bbox="458 911 1328 973">Active when the Function Selector Knob is in any of these four positions: EPAP, S, S/T or T.</p> <p data-bbox="458 1007 1379 1042">When the BiPAP 30 unit is delivering EPAP, the adjacent LED is illuminated.</p> <p data-bbox="458 1072 749 1106">Range = 4 to 30 cm H₂O.</p> <p data-bbox="458 1136 1384 1317">NOTE: The IPAP and EPAP controls are electrically coupled. The unit will not deliver a higher EPAP than IPAP level. If the EPAP control is set higher than the IPAP control, the BiPAP 30 unit will be locked to the IPAP setting, and the IPAP LED will remain lit. This is active throughout the S, S/T, and T settings on the Function Selector Knob.</p>

REAR PANEL

Warning / Avertissements

WARNING: This is a NON-CONVENTIONAL ventilator and is intended to support patient breathing. IT MUST NOT BE USED AS A LIFE SUPPORT VENTILATOR. It is not intended to provide the total ventilatory requirements of the patient.

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Danger: Shocking Hazard. Do not use in the presence of flammable anesthetics.

Danger: Flaming Anesthetics. Do not utilize an anesthetic without proper instructions.

Caution: To avoid electrical shock, disconnect supply before changing line.

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Caution: Voltage Selector Switches, before, should be set for proper line voltage (North America - 115 V).

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Warning: For continued protection against risk of fire, replace with same type and rating of fuse.

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

Serial No. 000000

CONTROL PANEL DOOR

Contains the warnings, cautions, and electrical information necessary for safe operation of the BiPAP 30 System.

VOLTAGE SELECTOR SWITCH

Set this switch to the appropriate line voltage.
 U.S.A. and Canada—115 Volts
 European and other countries—Check the country's line voltage requirements for the correct setting.

ELAPSED TIME INDICATOR

Provides resolution to 1/10 hour. Runs only when the Power switch is turned ON (I).

Chapter 4: Modes of Operation

4.1 Continuous Positive Airway Pressure (CPAP)

CPAP can be delivered in either of two ways:


1. **Function Selector Knob: IPAP**
Active Control: IPAP
2. **Function Selector Knob: EPAP**
Active Control: EPAP

NOTE: For this combination of controls, the active pressure control being used actually indicates the level of CPAP being delivered.

In the CPAP mode, the patient breathes spontaneously at the indicated pressure level. The patient is in total control of both rate and tidal volume. The BiPAP 30 System responds as necessary to satisfy the patient's flow demand and circuit leak while maintaining the preset pressure at a stable level.


Example:

CPAP = 10 cm H₂O


Example:

IPAP = 15 cm H₂O

Vest


P

EPAP = 10 cm H₂O


Example:

IPAP = 15 cm H₂O

EPAP = 4 cm H₂O

Vest


P

Refer to the examples and diagrams on the following page.

NOTE: For S/T-D 30 Units only— When using a recorder, an event marker on the V_{est} channel identifies timer-triggered breaths in the S/T mode. The event marker is a deflection below baseline at the initiation of the inspiratory flow.


Example:

IPAP = 15 cm H₂O
 EPAP = 10 cm H₂O
 BPM = 10

Example:

IPAP = 15 cm H₂O
 EPAP = 4 cm H₂O
 BPM = 12


Example:

IPAP = 15 cm H₂O
 EPAP = 10 cm H₂O
 BPM = 10
 % IPAP = 30

Based on these settings:

Respiratory cycle = 60 seconds/10 BPM = 6 seconds
 % IPAP = 6 seconds x 30% = 1.8 seconds
 EPAP = 6 seconds - 1.8 seconds = 4.2 seconds


Example:

IPAP = 15 cm H₂O
 EPAP = 4 cm H₂O
 BPM = 12
 % IPAP = 40

Based on these settings:

Respiratory cycle = 60 seconds/12 BPM = 5 seconds
 % IPAP = 5 seconds x 40% = 2 seconds
 EPAP = 5 seconds - 2 seconds = 3 seconds

Chapter 5: Setting Up the BiPAP 30 System

- Step 1** See Figure 1. Connect the inlet of a bacteria filter to the airflow outlet port on the front of the BiPAP 30 unit. Connect one end of the patient tubing to the outlet of the bacteria filter.


Figure 1

To protect the unit from contamination, place a low resistance main flow bacteria filter on the outlet port of the unit. The filter must be positioned in the circuit so that it is filtering air directly from the unit. If a humidifier is used, the bacteria filter should not be placed in the circuit between the outlet of the humidifier and the patient, since this can cause an increase in filter resistance.

The addition of a low resistance main flow bacteria filter does not effect the trigger sensitivity of the system. Depending upon the degree of resistance imposed by the filter, the pressure stability (the ability of the unit to maintain the delivered pressure within specified tolerances) may be affected. The filter must remain clean and dry to minimize circuit resistance. The main flow bacteria filter should be replaced in accordance with the manufacturer's specifications and instructions.

The bacteria filter should be in-line when the pressure settings are determined. If a filter is added to a system already in use, observe the patient to determine whether the pressure settings should be adjusted to compensate for any pressure drop resulting from the filter.

NOTE: Heat moisture exchangers should not be used with the BiPAP 30 System. Testing has shown that the ability of the unit to maintain the prescribed pressure is significantly affected by the addition of a heat moisture exchanger to the patient circuit.

Step 3 Plug the electrical cord into the power connector on the back of the unit. See Figure 3.


Figure 3

Step 4 Make sure the voltage selector switch on the back of the unit is set correctly. See Figure 4.

United States and Canada setting: 115 V
European and other countries: Check the country's line voltage requirements for the correct setting.

Step 5 Make sure the ON/OFF switch is OFF (0), then plug the electrical cord into a three-prong outlet. See Figure 5.

NOTE: In the United States and Canada, the use of a three-prong to two-prong adapter or extension cord is not recommended.


Figure 4

Step 6 S/T-D 30 units only—Connecting a Recorder:

If you are using a Detachable Control Panel 30 (DCP 30), refer to the DCP 30 instruction manual for details on connecting a recorder.

Step 7 Turn the unit on, occlude the patient port of the circuit, and confirm the presence of flow from the exhalation port.


Figure 5

Chapter 6: Administering Oxygen with the BiPAP 30 System

WARNING: Oxygen should be administered only on the order of a physician.

WARNING: Oxygen supports combustion. Oxygen should not be used while smoking or in the presence of an open flame.

WARNING: Do not add oxygen to the front of the BiPAP 30 unit (at the inlet filter opening).



NOTE: When using oxygen, turn the BiPAP 30 System ON (I) before turning on the oxygen source. The oxygen should then be turned off before the unit is turned OFF (O).

To add oxygen to the circuit, a controlled, low - flow source is needed (less than 15 L/min).

- Attach the oxygen tubing to one of the sample ports on or near the patient's interface, as indicated in Figure 21a.
- Do not use a flow rate greater than 15 L/min.
- A bacterial filter must be inserted between the outlet on the unit and the circuit tubing (See Figure 6b).
- Turn the BiPAP 30 System on before turning on the oxygen supply.
- Turn off the oxygen supply before turning the BiPAP 30 System off.
- If power to the BiPAP 30 System has been turned off for any reason before the oxygen supply is turned off, allow the oxygen to dissipate from the unit before turning the BiPAP 30 System on.

When adding oxygen to the patient circuit, the liter flow is gradually adjusted until the patient's oxygenation needs are adequately met. If adequate oxygenation cannot be achieved, an alternate method of ventilatory support should be considered.

If the method of oxygen administration is changed at any time, or as patient status changes, the patient should be monitored to make sure that his or her oxygenation level is sufficient.

Chapter 7: Administering Aerosolized Medication with the BiPAP 30 System

Small volume nebulizers or adaptors, which allow the use of metered dose inhalers, may be added to the patient circuit. If a small volume nebulizer is used, testing has shown that the liter flow used to drive the nebulizer does not impact the functioning of the BiPAP 30 System.

The aerosol delivery system used should be added to the patient circuit on the patient side of the continuous flow exhalation port. This will minimize the impact of the continuous circuit flow on the aerosol treatment. If a main flow bacteria filter is not being used, one should be added during the treatment period to prevent any of the aerosolized medication from entering the BiPAP 30 System via the patient circuit.

Inspect the exhalation device after each aerosol administration to ensure that the ports have not been occluded by the medication. If necessary, clean or replace the exhalation device.

WARNING: Aerosolized medications may collect in the Exhalation Port. To avoid an obstruction to flow, inspection of the Exhalation Port is essential during and after medication delivery.

Chapter 8: Performance Verification

WARNING: Performance verification should be performed prior to each patient setup of the BiPAP 30 System.

NOTE: For performance verification, Respiration recommends use of a DCP 30, output to a recorder, or a manometer capable of reading pressures in centimeters of water.

Step 1 Open the control panel door on the rear of the unit. Set the Function Selector Knob to EPAP. Set EPAP to minimum.

Step 2 Connect the patient circuit following the Set-up Instructions in this manual. Connect the proximal pressure tubing to one of the pressure ports on the mask or exhalation valve.

Step 3 Adjust the settings on the control panel as follows:

Function Selector Knob: IPAP
 Controls: IPAP 4 cm H₂O
 EPAP 4 cm H₂O
 BPM 10
 %IPAP 30

NOTE: If you are using a BiPAP 30 Detachable Control Panel 30 (DCP 30), make all setting adjustments on the DCP 30. When the DCP 30 is connected, the remote control settings override the settings on the BiPAP 30 unit's control panel.

Step 4 For S/T-D 30 units with Pressure (P) output to a DCP 30: Connect the remaining end of the proximal pressure tubing to the BiPAP 30 pressure port on the front panel of the unit.

For S/T 30 units or S/T-D 30 units not using a DCP 30: Connect the remaining end of the proximal pressure tubing to the manometer.

Step 5 Turn the BiPAP 30 unit ON (I) and allow it to warm up for at least 15 minutes prior to performance verification.

Step 6 Completely occlude the patient circuit outlet. Turn the BiPAP 30 unit ON (I). Verify flow from the exhalation port.

If using a recorder, verify recorder calibration:

- a. Set the Function Selector Knob to Zero.
- b. Calibrate the recorder to "0" for V_T , V_{est} , and V_{tot} .
- c. Set the Function Selector Knob to Cal.
- d. Calibrate the recorder to +0.5 V, -0.5 V for V_T and V_{est} and +0.5 V, 0.0 V for V_{tot} .
- e. Disconnect the proximal pressure tubing from the BiPAP S/T-D 30 unit, calibrate P to "0."
- f. Reattach the proximal pressure tubing to the BiPAP 30 unit's pressure port.


Chapter 9: Detachable Control Panel 30

9.1 General Description

The Detachable Control Panel 30 (DCP 30) is an accessory available for use with the BiPAP S/T-D 30 Ventilatory Support System. The DCP 30 combines all of the functional controls found on the rear of the BiPAP 30 unit with a selectable display that provides pressure, estimated exhaled tidal volume, and estimated circuit leak data. The DCP 30 allows the operator to make adjustments to settings from a remote location. Output connectors that allow the patient output signals to be displayed on a recorder can be found on the side of the DCP 30.

CAUTION: Use only Respirationics cables designed for use with the BiPAP 30. Use of any other cable types may cause damage to the BiPAP 30 System.

9.2 Description of Controls

Refer to the Description of Controls section for complete descriptions of the panel controls. The function controls on the DCP 30 duplicate those of the BiPAP S/T-D 30 unit.

NOTE: Refer to the Detachable Control Panel 30 Operating Instructions for complete information.

9.4 Assembly Instructions

- Step 1** With the power to the BiPAP S/T-D 30 System turned off, use the BiPAP S/T-D 30/Remote Cable to connect the DCP 30 to the BiPAP S/T-D 30 System. When connecting the cable to the unit, make sure the cable is securely screwed into place. (NOTE: To disconnect the cable from the unit, unscrew the pins on the sides of the cable and pull.) See Figure 7.

WARNING: Damage to the BiPAP S/T-D 30 Remote Cable may affect the operation of the BiPAP S/T-D 30 System. If the cable is damaged, replace it immediately. Never disconnect by pulling on the cable. Always use the end connectors to release the cable.


Figure 7

- Step 2** Attach the patient circuit to the BiPAP S/T-D 30 unit and connect the proximal pressure line. See Figure 8.

NOTE: For pressure monitoring, use only tubing with a smooth inner lumen.

- Step 3** If patient output recordings are desired, use the Remote/Recorder cables to connect the DCP 30 to the recorder. See Figure 9.

CAUTION: Before connecting the remote output cables, first connect the BiPAP S/T-D 30 System to a grounded outlet, then connect the DCP 30. This sequence eliminates the potential for electro-static discharge damage.

NOTE: Use only Respiration recorder cables. Do not substitute any other type of cables.

- Step 4** Verify recorder calibration. Turn the BiPAP S/T-D 30 unit on.

- Set the DCP 30 Function Selector Knob to Zero.
- Set the recorder to "0" for V_T , V_{est} and V_{tot} .
- Set the Function Selector Knob to Span.
- Calibrate the recorder to +0.5 V, -0.5 V for V_{est} and V_{tot} and to 0.0 V, +0.5 V for V_T .
- Disconnect the proximal pressure line from the BiPAP S/T-D 30 unit and set the recorder pressure (P signal) to "0."
- Reattach the proximal pressure line to the BiPAP S/T-D 30 System pressure port.


Figure 8


Figure 9

Chapter 10: Cleaning Instructions

10.1 Cleaning the BiPAP 30 System

CAUTION: Unplug the BiPAP 30 System before cleaning it. Do not immerse it in water or allow any liquid to enter the cabinet. Do not get the filter wet.

CAUTION: DO NOT clean any parts of the system with alcohol or cleaning solutions containing alcohol. DO NOT clean the system by steam autoclave or gas sterilization methods. These cleaning processes may harden or deform the flexible plastic parts of the system and adversely affect their function.

Step 1 Using a cloth slightly dampened with water and a mild cleanser, wipe off the outside of the BiPAP 30 System enclosure.

Step 2 Let the BiPAP 30 unit dry before plugging it in again.

10.2 Cleaning the Detachable Control Panel 30

CAUTION: Unplug the BiPAP S/T-D 30 System before cleaning it.

CAUTION: DO NOT immerse the DCP 30 in water or allow any liquid to enter the unit.

CAUTION: DO NOT clean any parts of the system with alcohol or cleaning solutions containing alcohol. DO NOT clean the system by steam autoclave or gas sterilization methods. These cleaning processes may harden or deform the flexible plastic parts of the system and adversely affect their function.

Step 1 Using a cloth slightly dampened with water and a mild cleanser, wipe off the outside of the DCP 30 enclosure.

Step 2 Let the DCP 30 unit dry before use.

10.3 Cleaning the Patient Circuit

Refer to the instructions provided with the circuit components you are using.

Chapter 11 Routine Maintenance

Refer to the Service Manual for more complete maintenance information and requirements.

11.1 Changing the Filter

CAUTION: Failure to replace a dirty filter may cause high operating temperatures of the BiPAP 30 unit, may reduce the flow, and may reduce the output pressure.

The white filter on the front of the unit is disposable and needs to be replaced after thirty days of use or sooner if it seems especially dirty.

- Step 1 Turn OFF (O) the BiPAP 30 unit and unplug the electrical cord from the wall outlet and from the back of the unit.
- Step 2 Locate the button in the center of the filter holder beneath the surface of the filter. Push in on the button and pull the filter out from all four sides of the filter holder. See Figure 10.
- Step 3 Throw away the dirty filter. The filter is not washable.
- Step 4 Center a new filter over the filter holder. Carefully push in on the center button and tuck the filter in on all four sides. See Figure 11.
- Step 5 Release the button. The filter should be intact and fit securely, covering the entire holder. Remove and readjust the filter if necessary.

The filter must be in place at all times when the BiPAP 30 System is being used.


Figure 10


Figure 11

Chapter 12: Troubleshooting

The table below lists common problems you may have with your BiPAP 30 System and possible solutions to those problems. If none of the corrective actions listed seems to work, contact your medical equipment distributor.

PROBLEM	POSSIBLE CAUSES	CORRECTIVE ACTION
1. No air flow from the unit.	Voltage selector switch is set incorrectly.	Verify that the voltage selector switch is not set too high (i.e., your electrical line is 115 V but your switch is set on 230 V).
	Possible internal problem.	Call your medical equipment distributor. Obtain a replacement unit.
2. Unit does not operate when turned on.	Power cord not firmly connected to the unit or the wall outlet.	Check that the power cord is plugged into the unit and into the wall outlet.
	Unit not connected to a "live" outlet.	Check to ensure that the unit is plugged into a "live" outlet (e.g., plug a working lamp into the outlet).
	Unit fuse(s) has blown.	Replace the fuses. See Section 12.2. NOTE: Always replace both fuses at the same time. Verify that the voltage selector switch setting is correct. If the new fuses blow out, contact your medical equipment distributor.
3. Unit stops and starts.	Power cord not securely connected to the unit or the outlet.	Check that the power cord is completely connected at the rear of the unit and at the outlet.

Chapter 13: Abbreviations/Definitions

Assist Ventilator	A device designed to assist the patient's spontaneous breathing and augment ventilation.
Auto-trigger	Movement to the IPAP level in the S or S/T modes that was not caused by patient effort or the BPM control.
BPM	Breaths Per Minute. The rate at which the unit will trigger to the IPAP level in the S/T mode if the patient does not initiate a spontaneous trigger. The rate at which the unit will cycle to IPAP in the T mode.
CPAP	Continuous Positive Airway Pressure. The mode that is delivered when the Function Selector Knob is in either the IPAP or the EPAP position. A constant level of pressure that is maintained throughout the patient's respiratory cycle during spontaneous breathing.
EPAP	Expiratory Positive Airway Pressure. The pressure level that is maintained during the expiratory phase of the S, S/T, or T modes. The pressure level that is maintained continuously with the Function Selector Knob in the EPAP position.
Estimated Leak Flow Rate	The Total Flow Rate is analyzed to determine what portion of the flow through the circuit is due to intentional (Whisper Swivel or Respironics Disposable Exhalation Port) and unintentional (mask fit) leaks. This continuous unintentional leak level is referred to as the Estimated Leak Flow Rate (\dot{V}_{leak}).
Estimated Patient Flow Rate	The Total Flow Rate is measured, and the Estimated Leak Flow Rate is determined. The balance of the flow is due to the instantaneous changes that occur during patient inspiration/expiration. This is referred to as the Estimated Patient Flow Rate (\dot{V}_{est}).
Instantaneous Flow	Changes in the flow rate through the patient circuit due to inspiration/expiration.
Intentional Leak	Component of the Estimated Leak Flow Rate (\dot{V}_{leak}) caused by flow through a continuous bleed exhalation port such as the Whisper Swivel or the Disposable Exhalation Port.
IPAP	Inspiratory Positive Airway Pressure. The pressure level that is maintained during the inspiratory phase of the S, S/T, or T modes. The pressure level that is maintained continuously with the Function Selector Knob in the IPAP position.
%IPAP	The percent of the total respiratory cycle (the duration of which is set by the BPM control) that the unit remains at the IPAP level in the Timed mode.
P Output	S/T-D 30 units with DCP 30 only. Pressure. Data from an internal pressure transducer is transmitted to a recorder, providing a tracing of the delivered IPAP and EPAP levels.

V_{leak}	The Total Flow Rate is analyzed to determine what portion of the flow through the circuit is due to intentional (Whisper Swivel or Respironics Disposable Exhalation Port) and unintentional (mask fit) leaks. This continuous unintentional leak level is referred to as the Estimated Leak Flow Rate (V_{leak}). The V_{leak} is a calculated component of V_{tot} . Composed of intentional (Whisper Swivel or Disposable Exhalation Port) and unintentional (mask) leaks.
V_{tot} Output	S/T-D 30 units with DCP 30 only. Total Flow Rate. Data from the internal flow transducer is transmitted to a recorder, providing a tracing of the total flow from the unit.
V_T Output	S/T-D 30 units with DCP 30 only. Estimated Tidal Volume in liters. V_{est} is converted to a volume reading and transmitted to a recorder, providing a tracing of the patient's tidal volume.
VT_E	Expiratory Tidal Volume.
VT_I	Inspiratory Tidal Volume.

Chapter 14: Specifications

FUNCTIONS:

IPAP: Delivers set IPAP pressure continuously

EPAP: Delivers set EPAP pressure continuously

Modes: Spontaneous
Spontaneous/Timed
Timed
CPAP

PRESSURE AND FLOW CHARACTERISTICS:

Pressure Range: IPAP: 4 cm H₂O to 30 cm H₂O* at ±60 L/min flow
EPAP: 4 cm H₂O to 30 cm H₂O* at ±60 L/min flow
Pressure controls are coupled so that the unit will not deliver a higher EPAP than IPAP level when the unit is operated in the S, S/T, or T modes.

Pressure Calibration Accuracy: ± 2 cm H₂O of the indicated control setting

Pressure Setting Stability: ± 1 cm H₂O over 8 hours

NOTE: For specifications, all pressure measurements were made at the blower outlet unless stated otherwise.

*Measured at the end of a 6-foot section of smooth inner lumen tubing.

CYCLING:

Breaths per Minute (BPM): 4 to 30

Breathing Rate Accuracy: Timed mode ± 10% or 2 BPM of the indicated reading, whichever is greater

% IPAP: 10 to 90%

% IPAP Setting Accuracy: ± 10% of the indicated setting


Minimum IPAP Time: Approximately 300 msec
S, S/T mode only for the spontaneous breath
Approximately 200 msec for S/T mode only for timer triggered breath

Maximum IPAP Time: 3 seconds ± 0.5 second at a constant flow rate condition of 100 L/min
S, S/T modes only

Chapter 15: BiPAP S/T-D 30 Outputs

S/T-D 30 Units Only

Output Signal Pin Connections

CAUTION: Care must be exercised to prevent electrical shorting between the pin connections. Only Respironics cables should be used. Use of any other cable types may cause damage to the BiPAP S/T-D 30 unit.


Output Signals Of DCP 30: (All Signals DC)

Estimated Patient Flow Rate - V_{est}
 1.0 V = 60 L/min
 Range = ± 150 L/min

Total Flow Rate - V_{tot}
 1.0 V = 60 L/min
 Range = ± 150 L/min

Estimated Tidal Volume - V_T
 1.0 V = 2 liters
 Range = 0 to 4 liters

Pressure - P
 0.1 V = 1 cm H₂O
 Range = 0 to 35 cm H₂O

CAUTION: The output signals are not isolated from the BiPAP 30 power system. The signals must be used only when connecting to a Respironics approved system with a Respironics approved cable.

Limited Warranty

Respironics, Inc. warrants that the BiPAP 30 System, shall be free from defects of workmanship and materials and will perform in accordance with the product specifications for a period of one year from the date of sale by Respironics, Inc. to the dealer. If the product fails to perform in accordance with the product specifications, Respironics, Inc. will repair or replace, at its option, the defective material or part. Respironics, Inc. will pay customary freight charges from Respironics, Inc. to the dealer location only. This warranty does not cover damage caused by accident, misuse, abuse, alteration, and other defects not related to material or workmanship.

RESPIRONICS, INC. DISCLAIMS ALL LIABILITY FOR ECONOMIC LOSS, LOSS OF PROFITS, OVER-HEAD OR CONSEQUENTIAL DAMAGES WHICH MAY BE CLAIMED TO ARISE FROM ANY SALE OR USE OF THIS PRODUCT. SOME STATES DO NOT ALLOW THE EXCLUSION OR LIMITATION OF INCIDENTAL OR CONSEQUENTIAL DAMAGES, SO THE ABOVE LIMITATION OR EXCLUSION MAY NOT APPLY TO YOU.

THIS WARRANTY IS GIVEN IN LIEU OF ALL OTHER EXPRESS WARRANTIES. IN ADDITION, ANY IMPLIED WARRANTIES, INCLUDING ANY WARRANTY OF MERCHANTABILITY OR FITNESS FOR THE PARTICULAR PURPOSE, ARE LIMITED TO ONE YEAR. SOME STATES DO NOT ALLOW LIMITATIONS ON HOW LONG AN IMPLIED WARRANTY LASTS, SO THE ABOVE LIMITATION MAY NOT APPLY TO YOU. THIS WARRANTY GIVES YOU SPECIFIC LEGAL RIGHTS, AND YOU MAY ALSO HAVE OTHER RIGHTS WHICH VARY FROM STATE TO STATE.

To exercise your rights under this warranty, contact your local, authorized Respironics, Inc. dealer or Respironics, Inc. at:

1001 Murry Ridge Drive
Murrysville, Pennsylvania 15668-8550
(412/733-0200)

Index

E

Elapsed Time Indicator 13
EPAP Control 11

F

Front Panel
 ON/OFF Control 9
 Pressure Port 9
 Tubing Connector 9
Function Selector Knob 10
Fuse Holder 14

I

IPAP Control 11

L

Limited Warranty 51

M

Maintenance 39-40
 Changing the Filter 39-40
 Changing the Fuses 40
Model/Serial Number Label 14
Modes of Operation 15
 Continuous Positive Airway Pressure 15-16
 Spontaneous 16-17
 Spontaneous/Timed 18-19
 Timed 20-33

N

Notes 6

O

Power Control 9
Output Jack, DCP 34
Oxygen Administration 25-27

P

Patient Cautions 7
Performance Verification 31-32, 36-37
Power Cord Connector 14
Pressure Port 9

THE UNIVERSITY OF CHICAGO PRESS