
Neu ab:
New since: 04.2009
Nouveau depuis:
Nuevo desde:

HELIODENT^{PLUS}

Schaltunterlagen
Wiring References
Schémas de circuits
Colección de esquemas

! **ACHTUNG**

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! **ATTENTION**

*Interference with electromedical devices by radio telephones:
To guarantee the operational safety of electromedical devices, the operation of mobile radio telephones in the medical practice or hospital area is prohibited.*

! **ATTENTION**

*When opening the unit:
Please observe the usual precautionary measures for handling PCBs (ESD).
Touch a ground point to discharge static electricity before handling any components*

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
Legend

DX1 Generator board

S1	Main switch
F200	Power Fuse
F201	Power Fuse
F300	Fuse for switched mode power supply
V316	LED, Supply voltage +15V
V409	LED, Reference values apply
V413	LED, Base heater = on, Radiation = out
V414	LED, kV controller defective
V503	LED, kVactual cable not connected
V611	LED, Release activated
H400	Buzzer, radiation
D500	Glow lamp, output stage running

DX4 Display board

S1	Changeover switch, release internal/external
S2	Changeover switch, safety circuit ON/OFF
V108	LED, 3.3V Supply for display
V201	LED, Controller running
V600	LED, 8V Input voltage
V602	LED, 3.3V µController
H200	Buzzer, radiation

Membrane keyboard

S4	X-ray release key on the control membrane
LED	Signal LED for operational readiness

DX42 Display module

H1 X-ray tube assembly

Heiz-	Tube filament heating
Heiz+	Tube filament heating
Minus	Ground
kVactual	Measured high voltage
V	High-voltage transformer connection
W	High-voltage transformer connection

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** Property of Sirona. All rights reserved. **
Sirona Dental Systems GmbH