



SERVICE BULLETIN

KODAK Multiloader

Eastman Kodak Company, Health Imaging, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 242

May 1996

3053 *Kodak* Multiloader 700

Possible Missing Part from Modification Kit 26

Some of the Modification Kits shipped by Service Parts Management, Rochester NY, may be missing the BRAKE STOP, Part Number 928 0981.

If you are preparing to install this Modification, do not disassemble the Multiloader 700 until you have inventoried the kit and ensured that it is complete.

The kits currently being shipped from Service Parts Management, Rochester NY, have been updated with the missing part. Anyone whose kit did not contain the BRAKE STOP, Part Number 928 0981 can order the part from Service Parts Management.

Correction to Site Specification, Publication Number 8B7005

The correct dimensions for the location of the SEAL PLATE for the Through-the-Wall installation to allow proper draining of the HUMIDIFIER are given below:

On Page 34, referring to Figure 18:

- Change Dimension G to 220mm (8.7 in.).
- Change Dimension H to 430mm (16.9 in.).

On Page 36, referring to Figure 19:

- Change Dimension G to 220mm (8.7 in.).
- Change Dimension H to 430mm (16.9 in.).

On Page 38, in Figure 20, referring to the vertical location of the Opening for Plumbing Seal Plate in Plywood Panel A:

- Change 389mm (15.3 in.) to 220mm (8.7 in.).
- Change 599mm (23.6 in.) to 430mm (16.9 in.).

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

CES Personnel Contact:

Film Handling TAC
U.S.: 1-800-822-1414
Canada: 1-800-433-1414
International: 716-724-4675

Kodak and X-Omat are trademarks.

9B8330

