



SERVICE BULLETIN

Kodak X-Omat Processors

Kodak Industrex Processors

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 220

Revised: July, 1995

Kodak M35A-M X-Omat Processor
Kodak M35-M X-Omat Processor
Kodak Industrex Processor, Model M35-A
Kodak Industrex Processor, Model M35

Announcing the Availability of
Modification No. 5 for the M35A-M X-Omat Processor
Modification No. 2 for the M35-M X-Omat Processor
Modification No. 1 for the Industrex Processors, Models M35-A and M35

Deterioration may occur on the surface of the D ROLLER 1C7638 located in the DEVELOPER RACK. This condition may cause pickoff artifacts on *Kodak Min-R* Film.

Modification 5 for the M35A-M X-Omat Processor, Modification 2 for the M35-M X-Omat Processor, and Modification 1 for the *Industrex* Processors, Models M35-A and M35, help eliminate the potential for generating pickoff artifacts.

The MODIFICATION KIT, Part Number 8B6583, will be available at no charge until December 31, 1995.

All modifications are Type 1, Required. This modification **must** be installed on all DEVELOPER RACKS 914264 and 1C7315 purchased between September, 1994, and April, 1995, and on the DEVELOPER RACKS of all of the following processors:

Processor	Starting Serial Number	Ending Serial Number
M35A-M X-Omat	51981	52453
M35-M X-Omat	76000	76193
Industrex, Model M35-A	85030	85042
Industrex, Model M35	95030	95052

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg 20
Rochester, NY 14650-1131

1-(800)-336-4722

Please be prepared to give your Customer Number and your equipment serial number.

Kodak, Industrex, Min-R, and X-Omat are trademarks.

N-931

