



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging, Rochester, New York 14650

SERVICE BULLETIN NO. 259

February 1997

3027 KODAK RPX-Omat Processor, Model M7B

Possible Incorrect Connector

There is a possibility that some of the M7B Processors with the Serial Number range of 28378 - 28432 may have the incorrect CONNECTOR on the BUZZER ASSEMBLY, Part Number 594759 (Subs to 1C0033). See the Parts List, Publication Number, 981160, Figure 8-1.

When you do the next preventive maintenance, check that the BUZZER ASSEMBLY has the correct CONNECTOR, and that it operates properly. The BUZZER should sound 3 seconds after the film exits the DETECTOR ASSEMBLY. The correct CONNECTOR is red with clips to latch the CONNECTOR in position. The incorrect CONNECTOR is black, without clips to latch the CONNECTOR in position. Replace the incorrect BUZZER ASSEMBLY. There is no charge for the new BUZZER ASSEMBLY if the Serial Number of the Processor is provided when ordering the new BUZZER ASSEMBLY.

Ralph Peterson
Customer Equipment Services

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

CES Personnel Contact:

Film Handling TAC
U.S.: 1-800-822-1414
Canada: 1-800-433-1414
International: 716-724-4675

Kodak and X-Omat are trademarks.

9B8932

