



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 228

December 1995

Kodak AUTOMATIC FILM FEEDER, MODEL 11

Modification No. 2 for the *Kodak* AUTOMATIC FILM FEEDER, MODEL 11

The purpose of this modification is to prevent interruptions in film feeding caused by the effects of static and switch bounce on the FEEDER. It should only be installed on FEEDERS exhibiting problems of intermittently stopping feeding. The kit contains redesigned SUPPLY (200) and DISPLAY (300) PRINTED CIRCUIT BOARDS and a new SENSOR INPUT PRINTED CIRCUIT BOARD as well as TOOTHED WASHERS. A GROUND WIRE, other accessories and instructions necessary to complete the installation are also included.

This is a Type 1 Selective Modification. This modification is now available from Service Parts Management (1-800-431-7278), as Modification Kit 8B6580. This kit will be available at no charge for a period of 1 year.

FEEDERS with serial numbers 1597 and below were shipped without this modification installed. The modification has been incorporated into production so that units with serial number 1598 and above are shipped with the modification installed.

For more information, please contact:

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg 20
Rochester, NY 14650-1131

(800) 336-4722

Please be prepared to give your Customer Number and your equipment serial number.

Kodak is a trademark.

N-940

sb121495_3128.doc

Printed in USA

HEALTH SCIENCES DIVISION
EASTMAN KODAK COMPANY • ROCHESTER, N.Y. 14650

