



# SERVICE BULLETIN

## KODAK X-OMAT Processor

Eastman Kodak Company/Customer & Technical Services/Health Sciences Markets Division/Rochester, New York 14650

SERVICE BULLETIN NO. 208

November 1994

### *Kodak X-Omat 480 RA Processor* *Kodak X-Omat M7B Processor*

#### **Solid State Relay Change**

Starting with Serial Number 2377, the 480 RA Processors will be using new Solid State Relays that are better matched to the load they control.

#### Kodak X-Omat 480 RA Processor

SSR	LOAD	OLD	NEW
U1	Dryer Heater	20 A P/N 655290	45 A P/N 1C7329
U2	Developer Heater	20 A P/N 655290	45 A P/N 1C7329
U5	Fixer Heater	20 A P/N 655290	45 A P/N 1C7329
U3	Developer Replenishment Pump Motor	20 A P/N 655290	10 A P/N 1C7330
U4	Fixer Replenishment Pump Motor	20 A P/N 655290	10 A P/N 1C7330

#### **NOTE**

Make sure to use the appropriate amperage relay.

**See back page for Kodak X-Omat M7B Processor information.**

Starting with Serial Number 27,956 the M7B processors will be using new Solid State Relays that are better matched to the load they control.

**Kodak X-Omat M7B Processors**

SSR	LOAD	OLD	NEW
U1	Dryer Heater	20 A P/N 655290	45 A P/N 1C7329
U2	Developer Heater	20 A P/N 655290	45 A P/N 1C7329



Possible damage from electrostatic discharge. SOLID STATE RELAYS should be regarded as ESD-sensitive devices. Take appropriate precautions when handling them and during installation.

When ordering either of the new SOLID STATE RELAYS (10 A or 45 A), a METAL-OXIDE VARISTOR (MOV), Part No. 943378, is included with the relay.

Make sure to install the new MOV on the AC side of the SOLID STATE RELAY.

For more information, please contact:

Eastman Kodak Company  
Imaging Service Center  
Health Sciences Division  
100 Carlson Road  
Rochester, New York 14653-9013  
(716) 724-1789, Option 1

Please be prepared to give your Service Customer Number and your equipment serial number.

svcbul\_208.txt  
N-920

*Kodak and X-Omat are trademarks.*

