



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 223

SEPTEMBER 1995

KODAK X-OMAT 180 LP Processor KODAK X-OMAT 270 RA Processor KODAK X-OMAT Multiloader 300

Required Mod 6 Available for the KODAK X-OMAT 180 LP Processor

- Purpose:**
- Install new style HEAT EXCHANGER FITTINGS to minimize chemical leakage from the HEAT EXCHANGER connections.
 - Install new style RACK SPRINGS that are more resistant to corrosion than the existing RACK SPRINGS.
 - Install a new WASH TANK CLIP to increase the volume of water around the HEAT EXCHANGER in the WASH TANK to decrease false "Loss of Developer Cooling Ability" errors.

Mod 6 is a Type 1, required modification.

Parts are available at no charge from Service Parts Management, September, 1995.

Check the Serial No. on the Processor and in the table below. It may not be necessary to do the whole Mod on your Processor.

| 180 LP PROCESSOR Serial No. Range: | Parts Required: | Action Required: |
|---|--|---|
| 1367 and below | Order Mod Kit 5B6770 | Do Sections 2 through 6 of the Mod instructions |
| 1368 to 1429 | Order Mod Kit 5B6770 | Do Sections 2, 5 and 6 of the Mod instructions |
| 1430 to 1476 | None - Mod 6 was factory installed. The 6 was <u>not</u> circled on the Mod Label | Circle 6 on the MODIFICATION LABEL |
| 1477 and above | None - Mod 6 is factory installed | None |

Required Mod 18 Available for the KODAK X-OMAT 270 RA Processor

- Purpose:**
- Install new style RACK SPRINGS that are more resistant to corrosion than the existing RACK SPRINGS.
 - Install new style Developer and Fixer Replenishment FITTINGS to minimize chemical leakage from the current FITTINGS.
 - Install new style HEAT EXCHANGER FITTINGS to minimize chemical leakage from the HEAT EXCHANGER connections.
 - Install a new WASH TANK CLIP to increase the volume of water around the HEAT EXCHANGER in the WASH TANK to decrease false “Loss of Developer Cooling Ability” errors.

Mod 18 is a Type 1, required modification.

Parts are available at no charge from Service Parts Management, September, 1995.

Check the Serial No. on the Processor and in the table below. It may not be necessary to do the whole Mod on your Processor.

| 270 RA PROCESSOR Serial No. Range: | Parts Required: | Action Required: |
|---|---|--|
| 1930 and below | Order Mod Kit 5B6766 | Do Sections 2 through 7 of the Mod instructions |
| 1931 to 1980 | Order Mod Kit 5B6766 | Do Sections 2, 4, 5, 6 and 7 of the Mod instructions |
| 1981 to 2089 | Order Mod Kit 5B6766 | Do Sections 2, 5, 6 and 7 of the Mod instructions |
| 2090 to 2134 | Order Mod Kit 5B6766 | Do Sections 2, 6 and 7 of the Mod instructions |
| 2135 to 2220 | None - Mod 27 was factory installed. The 27 was <u>not</u> circled on the Mod Label | Circle 18 on the MODIFICATION LABEL |
| 2221 and above | None - Mod 18 is factory installed | None |

Required Mod 27 Available for the KODAK X-OMAT Multiloader 300

- Purpose:**
- Install new style RACK SPRINGS that are more resistant to corrosion than the existing RACK SPRINGS.
 - Install new style Developer and Fixer Replenishment FITTINGS to minimize chemical leakage from the current FITTINGS.
 - Install new style HEAT EXCHANGER FITTINGS to minimize chemical leakage from the HEAT EXCHANGER connections.
 - Install a new WASH TANK CLIP to increase the volume of water around the HEAT EXCHANGER in the WASH TANK to decrease false “Loss of Developer Cooling Ability” errors.

Mod 27 is a Type 1, required modification.

Parts are available at no charge from Service Parts Management, September, 1995.

Check the Serial No. on the Processor and in the table below. It may not be necessary to do the whole Mod on your Processor.

| MULTILOADER 300 Serial No. Range: | Parts Required: | Action Required: |
|--|---|---|
| 919 and below | Order Mod Kit 5B6766 | Do Sections 2 through 7 of the Mod instructions |
| 920 to 949 | Order Mod Kit 5B6766 | Do Sections 2, 5, 6 and 7 of the Mod instructions |
| 950 to 964 | Order Mod Kit 5B6766 | Do Sections 2, 6 and 7 of the Mod instructions |
| 965 to 984 | None - Mod 27 was factory installed. The 27 was <u>not</u> circled on the Mod Label | Circle 27 on the MODIFICATION LABEL |
| 985 and above | None - Mod 27 is factory installed | None |

For more information, please contact:

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg 20
Rochester, NY 14650-1131

(716) 724-1789, Option 1

Please be prepared to give your Customer Service Number and your equipment serial number.

Kodak and X-Omat are trademarks.

N-935

