



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 224

September 1995

KODAK X-OMAT Multiloader 300
KODAK X-OMAT 270 RA Processor
KODAK X-OMAT 180 LP Processors

Modification 29 for the KODAK X-OMAT Multiloader 300

Modification 19 for the KODAK X-OMAT 270 RA Processor

Modification 7 for the KODAK X-OMAT 180 LP Processor

All Kodak X-Omat 180 LPS Processors have the newly designed WASH RACK factory installed.

The purpose of this Mod is to improve the overall performance of the WASH RACK. If you are experiencing problems with film washing uniformity, corrosion of metal E-RINGS, or film transport failures, install this version of the WASH RACK for the following improvements:

- **More uniform washing.** A redesigned water MANIFOLD improves the distribution of water across the ROLLERS.
- **More corrosion resistance.** SCREWS now hold the SHAFTS of the ROLLERS in the WASH RACK, which eliminates the E-RINGS. The SPRINGS on the WASH RACK are still made of a special, high-quality stainless steel alloy.
- **Improved film transport.** A new GUIDE SHOE and an additional IDLER ROLLER improve the transport of film from the WASH RACK into the DRYER RACK. Pinned ROLLERS keep the drive and driven ROLLERS correctly engaged.

This is a Type 1 Selective Modification, which means that it may be needed on some processors to improve performance.

The modification is now available from Service Parts Management, as Modification Kit 8B6612. This kit will be available at no charge for a period of 6 months.

Following are the serial number ranges that do not have this Mod factory installed.

| | |
|------------------------|------------|
| Multiloader 300..... | 101 - 1040 |
| 270 RA Processor | 101 - 2397 |
| 180 LP Processor..... | 101 - 1644 |

This modification has been incorporated into production units now being shipped from the factory.

For more information, please contact:

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg 20
Rochester, NY 14650-1131

(800) 336-4722, Press 5

Please be prepared to give your Customer Number and your equipment serial number.

Kodak and X-Omat are trademarks.

N-936

