



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 249

July 1996

3434 KODAK X-OMAT 3000 RA Processor

Fuses in the Kodak X-Omat 3000 RA Processor

There is a possibility that some *Kodak X-Omat* 3000 RA Processors left the factory with fuses that have the incorrect current rating. The fuses in question are F1 and F5. The correct current rating for these fuses is 1.25 amps. The incorrect fuse that may have been installed has a current rating of 5 amps. This is not a safety or reliability issue. However, in order to meet agency approvals and for version consistency, the units within the serial number range of 3322 - 3465 should be checked on your next service call. If necessary, install the proper fuse, Part Number 454771.

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

CES Personnel Contact:

Film Handling TAC
U.S.: 1-800-822-1414
Canada: 1-800-433-1414
International: 716-724-4675

Kodak and X-Omat are trademarks.

9B8344

