



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 229

January 1996

Kodak X-Omat M43 Processor
Kodak X-Omat M43A Processor
Kodak X-Omat Clinic 1 Processor

Availability of the Anti-Biological Growth Upgrade Kit 8B6930

To reduce the potential of biological growth in the wash system, a Continuous Water Usage Mode feature has been developed. This feature includes a software and hardware modification to the M43, M43A, and Clinic 1 PROCESSORS. The Anti-Biological Growth Upgrade Kit 8B6930 will be available January 1, 1996.

Products with the serial numbers listed below will be manufactured with this feature.

| Serial No. | Processor |
|---------------|-----------|
| 467 and above | M43A |
| 526 and above | M43 |
| 534 and above | Clinic 1 |

Note

Please note that for the Anti-Biological Growth Upgrade Kit to operate, Modification 1, 1C7782, **must first be installed.**

Replacement of EXIT BAR 5B6291 with EXIT BAR 1C7748

The EXIT BAR 5B6291, located on the TOP COVER of the PROCESSOR, has a MYLAR STRIP that can cause scratches on film leaving the DRYER. The new EXIT BAR 1C7748, which prevents these scratches, is now available. No special tool is required to do this replacement. Products with the serial numbers listed below will be manufactured with the new EXIT BAR.

| Serial No. | Processor |
|---------------|-----------|
| 467 and above | M43A |
| 526 and above | M43 |
| 534 and above | Clinic 1 |

Replacement of PHENOLIC ROLLER 458423 with PHENOLIC ROLLER 8B6568

The PHENOLIC ROLLER 458423, located in the WASH RACK, may bow and cause film transport problems. The new PHENOLIC ROLLER 8B6568, which prevents film transport problems, is now available. See the SERVICE MANUAL for the *Kodak X-Omat* M43 and M43A PROCESSORS and the *Kodak X-Omat* Clinic 1 PROCESSOR, Publication No. 981090 for installation instructions. Products with the serial numbers listed below will be manufactured with the new PHENOLIC ROLLER.

| Serial No. | Processor |
|-------------------|------------------|
| 463 and above | M43A |
| 518 and above | M43 |
| 534 and above | Clinic 1 |

For more information, please contact:

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg 20
Rochester, NY 14650-1131

(800) 336-4722, Press 2

Please be prepared to give your Customer Number and your equipment serial number.

Kodak and X-Omat are trademarks.

N-941

