



# SERVICE BULLETIN

---

## KODAK X-OMAT Multiloader

---

Eastman Kodak Company, Health Imaging, Rochester, New York 14650

SERVICE BULLETIN NO. 257

December 1996

**3053 KODAK X-OMAT Multiloader 700**  
**3443 KODAK X-OMAT Multiloader 700 Plus**  
**3066 KODAK X-OMAT Multiloader 700 Side-By-Side Kit**

### Amendment to Site Specifications

#### Serviceability

When installing a Multiloader 700 or a Multiloader 700 Plus, it is important that consideration is given to the condition of the floor. It will be necessary from time to time to un-dock and separate the Multiloader from the Processor and/or Side-By-Side, to gain access to the Processor or the Multiloader 700/700 Plus for servicing.

The floor should be:

- Free from obstructions such as grills, drains, cables, etc.
- Level and flat
- Resistant to denting from the weight of the equipment
- Smooth enough to allow the free passage of the equipment castors

#### Personal Safety

When carrying out an un-docking procedure, it is necessary that the correct procedure is used to separate the units.

- Lift the TOP COVER and pull from the front of the Multiloader 700 in a straight line.
- Avoid twisting your back.
- Use your legs to apply the force necessary.
- Do not pull from an awkward angle.
- If you have a known weakness you should get assistance from a colleague.

These points are even more important where the equipment is installed close to, or at the minimum distance from the wall, as laid down in the Site Specifications, Part No. 8B7005.

Ralph Lipka  
Customer Equipment Services

**For More Information:**

**U.S. Distributors Contact:**

Health Imaging  
Customer Support Operations  
1-800-336-4722

**CES Personnel Contact:**

Film Handling TAC  
U.S.: 1-800-822-1414  
Canada: 1-800-433-1414  
International: 716-724-4675

*Kodak and X-Omat are trademarks.*

9B8926

