



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging, Rochester, New York 14650

SERVICE BULLETIN NO. 250

October 1996

3058 KODAK X-OMAT Multiloader 300
3059 KODAK X-OMAT 270 RA Processor
3198 KODAK X-OMAT 180 LPS Processor
3227 KODAK X-OMAT 180 LP Processor
3133 KODAK M35-AM X-OMAT Processor
3408 KODAK M35-M X-OMAT Processor

New Type 1 Required Modification Kit Available:

Modification No. 35 for the Kodak X-Omat Multiloader 300

Modification No. 21 for the Kodak X-Omat 270 RA Processor

Modification No. 8 for the Kodak X-Omat 180 LP and 180 LPS Processors

The purpose of the Modification is to improve the reliability of the Processor by:

- Minimizing the possibility of corrosion of the electrical connectors by moving the connections for the DRYER THERMISTOR and the DRYER OVER-TEMPERATURE THERMOSTAT to the outside of the TANK.
- Reducing the possibility of dilution of the fixer by installing a new YOKE ASSEMBLY which more accurately directs the water into the FIXER/WASH CROSSOVER TROUGH ASSEMBLY.
- Eliminating corrosion of the SPROCKETS in the DRYER RACK by installing plastic SPROCKETS.

Processor	Serial Number Range
<i>Kodak X-Omat Multiloader 300</i>	101 - 1499 (Serial Number is for the Processor)
<i>Kodak X-Omat 270 RA Processor</i>	ALL
<i>Kodak X-Omat 180 LP Processor</i>	101 - 3399
<i>Kodak X-Omat 180 LPS Processor</i>	10,001 - 13,199

Modification Kit, Part No. 9B8368, is now available from Service Parts Management. This kit will be available at no charge for a period of 6 months.

- In order to install this Modification, you will need the following tools:

Wire Crimping Tool TL-2020 (or the equivalent)

Wire Stripping Tool TL-2275 (or the equivalent)



Caution

Some of the existing SPROCKETS may be so corroded that you could cause damage to the SHAFT of the SQUEEGEE DRIVE ROLLER when you remove the SPROCKET.

To avoid another service call and lost time, have a new SQUEEGEE DRIVE ROLLER with you.

The part numbers are:

Processor	Part Number
<i>Kodak X-Omat Multiloader 300</i>	1C0985
<i>Kodak X-Omat 270 RA Processor</i>	
<i>Kodak X-Omat 180 LP Processor</i>	8B7110
<i>Kodak X-Omat 180 LPS Processor</i>	



Note

Both Part Numbers sub to Kit No. 9B8087. The Kit includes the new improved SQUEEGEE ROLLER, Part No. 8B7110, and the parts and instructions to install it.

New TANK and New Publications for the *Kodak X-Omat 270 RA Processor* in the *Kodak X-Omat Multiloader 300*

Starting with Processor, serial number 1301, all *Kodak X-Omat 270 RA Processors* in the *Kodak X-Omat Multiloader 300*, have a new style TANK and an improved wash system. This is the same TANK and wash systems that is used in the *Kodak X-Omat 3000 RA Processor*, and was recently introduced into the *Kodak X-Omat 180 LP* and *180 LPS Processors* (see Service Bulletin No. 239). As with the *180 LP* and *180 LPS Processors*, this TANK and wash system are not meant to be retrofitted into existing Processors in the field.

The updated Parts List, IP3058-2, dated 5/96 which supersedes IP3058-2, dated 3/94 covers the new TANK and associated parts as well as general updates for the *Kodak X-Omat 270 RA Processor* in the *Kodak X-Omat Multiloader 300*. Please scrap any old stock and replace it with IP3058-2, dated 5/96. This is especially important for customers who order a Parts List for a Processor with a serial number greater than 1300.

New Smooth GUIDE SHOES Available for the *Kodak X-Omat 270 RA Processor* and the *Kodak M35-M, M35A-M X-Omat Processors*

In applications where film is run emulsion side down, the developer turnaround GUIDE SHOES can put + density marks ('tic' marks) on the edges of the film. In Mammography applications, these artifacts are objectionable because the edge of the film is part of the image area. SMOOTH GUIDE SHOES, which are interchangeable with the existing GUIDE SHOES, have been designed to significantly reduce these artifacts.

As announced in Service Bulletin No. 244, *Kodak Min-R 2000 Film* may exhibit non-uniform development in *Kodak M35A-M X-Omat* or *Kodak M35-M X-Omat Processors*. This non-uniform pattern can be reduced by processing the film emulsion side down. If marks from the GUIDE SHOES appear on *Kodak Min-R 2000 Film*, obtain a SMOOTH GUIDE SHOE Kit through your Kodak Sales Representative.

In the *Kodak X-Omat 270 RA Processor* in a *Kodak X-Omat Multiloader 300* and in the *Kodak M35-M X-Omat* or *Kodak M35A-M X-Omat Processor* with a Mini-Loader (Europe only), mammography film is always run emulsion side down. Therefore, these GUIDE SHOES are now available separately as a service part for these Processors.

It is strongly recommended that the Installation Instructions be used in order to ensure proper orientation and adjustment of these GUIDE SHOES.

The following Part Numbers, for the individual SMOOTH GUIDE SHOE and Instructions are available through Service Parts Management:

Processor	Part Number	Installation Instructions
<i>Kodak X-Omat 270 RA Processor in a Kodak X-Omat Multiloader 300</i>	8B6559	9B8336
<i>Kodak M35-M X-Omat or Kodak M35A-M X-Omat</i>	5B6445	9B8335

Modification Kits Updated:

Modification No. 29 for the Kodak X-Omat Multiloader 300

Modification No. 19 for the Kodak X-Omat 270 RA Processor

Modification No. 7 for the Kodak X-Omat 180 LP Processor

Effective June 20, 1996, the FLOW CONTROL (Dole) VALVE in this kit was changed from a $\frac{3}{4}$ Gallon Per Minute (gpm) to a $\frac{1}{2}$ gpm VALVE to reduce the amount of splashing at the water inlet on the WASH RACK. If you are experiencing splashing problems or fixer dilution after installing the above Modifications with the $\frac{3}{4}$ gpm VALVE, order the $\frac{1}{2}$ gpm FLOW CONTROL VALVE, Part No. 5B6630.

Note

The following Type 1 Required Modification must already be installed on the Processor when $\frac{1}{2}$ gpm FLOW CONTROL VALVE is used, or developer cooling problems may occur:

Processor	Modification No.	Part No.
<i>Kodak X-Omat Multiloader 300</i>	27	5B6766
<i>Kodak X-Omat 270 RA Processor</i>	18	5B6766
<i>Kodak X-Omat 180 LP Processor</i>	6	5B6770

Larry Weishaar
Service Engineer
Customer Equipment Service

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

CES Personnel Contact:

Film Handling TAC
U.S.: 1-800-822-1414
Canada: 1-800-433-1414
International: 716-724-4675

Kodak, Min-R and X-Omat are trademarks.

9B8346

