



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 251

August 1996

3041 *Kodak M35 X-Omat Processor*

Height Adjustment Kit

When docking a *Kodak M35 X-Omat Processor* to a *Kodak Ektascan 1120 Laser Printer*, remove each of the 4 CASTERS from the Laser Printer and place a SHIM between the Laser Printer and each CASTER. Then bolt each CASTER and SHIM to the Laser Printer. This procedure will adjust the height of the Laser Printer and improve the clearance of the film coming from the Laser Printer to the Processor.

A Height Adjustment Kit No. 4B4840 is available from Service Parts Management. The Kit contains 4 SHIMS and 16 SCREWS.

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

CES Personnel Contact:

Film Handling TAC
U.S.: 1-800-822-1414
Canada: 1-800-433-1414
International: 716-724-4675

Kodak, Ektascan and X-Omat are trademarks.

9B8347

