



# SERVICE BULLETIN

---

## KODAK X-OMAT Processor

---

Eastman Kodak Company, Health Imaging Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 236

Revised May 1996  
Supersedes April 1996

### 3255 KODAK X-OMAT M43 Processor

#### Kodak X-Omat M43 Processor Circuit Breaker

An incorrect 20 AMP Main Circuit Breaker CB1 may have been installed in some *Kodak X-Omat* M43 Processors. There is no additional risk of personal injury, fire, or safety created by this error. The correct Circuit Breaker is 10 AMP, Part No. 918552.

An unscheduled site visit is not required. At the next scheduled site visit, inspect all *Kodak X-Omat* M43 Processors and do the following:

- Verify that the correct 10 AMP Circuit Breaker, Part No. 918552, is installed. The AMP rating is located on the Drive Side of CB1.
- If the Circuit Breaker is the incorrect 20 AMP Circuit Breaker, replace it with the 10 AMP Circuit Breaker, Part No. 918552.
- For *Kodak X-Omat* M43 Processors that have incorrect Circuit Breakers installed, send the serial number and the machine date code to your Regional Specialist for feedback to Rochester.
- *Kodak X-Omat* M43 Processors are sold as export only.

#### For More Information:

##### U.S. Distributors Contact:

Health Imaging  
Customer Support Organization  
1-800-336-4722

##### CES Personnel Contact:

Film Handling TAC  
U.S.: 1-800-822-1414  
Canada: 1-800-433-1414  
International: 716-724-4675

*Kodak and X-Omat are trademarks.*

8B7497

sb040496\_236\_r.doc

Printed in USA

**HEALTH IMAGING DIVISION**  
EASTMAN KODAK COMPANY • ROCHESTER, N.Y. 14650

