**Balance** | **S7600**  
**Scales**

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**GB** "FAZZINI" wishes to thank you for purchasing one of its range of products and you are invited to read this booklet thoroughly. This booklet contains all of the information necessary to use the device properly. Please adhere closely to the notices contained within and conserve the booklet, without making changes, in a place where it is easy to find and is accessible to the operator. The contents of this manual can be changed, without warning or additional obligations, so that the manufacturer can include variations and improvements. It is prohibited to copy or translate any of this booklet without the written permission of "FAZZINI".

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**THANKS FOR BUYING ANOTHER FAZZINI PRODUCT**

The scale have a large weighing range, it is widely used in hospitals, clinics and the domestic use. It is also the basic examination device for medical and health sectors.

**ASSEMBLING INSTRUCTIONS**

*Cautions*

This precision device is very easy to be fitted together because all the major parts are pre-assembled in the factory. Please unpack the carton carefully and remove all the packages and fillers before starting assemble.

*Description*

Before using it, please take out the wadding paper between the two sides of the scale body, tighten up the pressing screws on the pan using a screwdriver and place the scale on the floor or table horizontally.

If you would like to use it on the mattress or flannelette blanket, please put wood board under it.

Before weighting, rotate the knob in order to make the index aim at "zero", and the bring 1-2 Kg pressure to bear on the scale body for several times: if the index can aim at "zero" accurately, it can be used.

If you would like to weight the infant who is standing up, just uninstall the pan and make the index at "zero" again, and then it can be used.

*Maintenance*

Don't rotate the knob overwhelmingly in order to avoid drop.

Don't dismantle it freely. Don't place it near fire, heat and steam, and avoid direct sunshine.

If there is stain on the scale, please clean it with soft cloth and neural cleaning agent.

Don't use hot water and organic solvent.

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