Plethysmograph

### Plethysmograph

# Plethysmograph Diagnostics MeSH D010991 [1]

A **plethysmograph** is an instrument for measuring changes in volume within an organ or whole body (usually resulting from fluctuations in the amount of blood or air it contains).

### **Organs studied**

#### Lungs

Pulmonary plethysmographs are commonly used to measure the functional residual capacity (FRC) of the lungs—the volume in the lungs when the muscles of respiration are relaxed—and total lung capacity.

In a traditional plethysmograph, the test subject is placed inside a sealed chamber the size of a small telephone booth with a single mouthpiece. At the end of normal expiration, the mouthpiece is closed. The patient is then asked to make an inspiratory effort. As the patient tries to inhale (a maneuver which looks and feels like panting), the lungs expand, decreasing pressure within the lungs and increasing lung volume. This, in turn, increases the pressure within the box since it is a closed system and the volume of the box compartment has decreased to accommodate the new volume of the subject.

Boyle's Law is used to calculate the unknown volume within the lungs. First, the change in volume of the chest is computed. The initial pressure and volume of the box are set equal to the known pressure after expansion times the unknown new volume. Once the new volume is found, the new volume minus the original volume is the change in volume in the box and also the change in volume in the chest. With this information, Boyle's Law is used again to determine the original volume of gas: the initial volume (unknown) times the initial pressure is equal to the final volume times the final pressure.

The difference between full and empty lungs can be used to assess diseases and airway passage restrictions. An obstructive disease will show increased FRC because some airways do not empty normally, while a restrictive disease will show decreased FRC. Body



**Plethysmograph** or "body box" used in lung measurements.



A man undergoing whole body plethysmography.

plethysmography is particularly appropriate for patients who have air spaces which do not communicate with the bronchial tree; in such patients gas dilution would give an incorrectly low reading.

Newer lung plethysmograph devices have an option which does not require enclosure in a chamber.

Plethysmograph

#### Limbs

Some plethysmograph devices are attached to arms, legs or other extremities and used to determine circulatory capacity. In water plethysmography an extremity, e.g. an arm, is enclosed in a water-filled chamber where volume changes can be detected. Air plethysmography uses a similar principle but based on an air-filled long cuff, which is more convenient but less accurate. Another practical device is mercury-filled strain gauges used to continuously measure circumference of the extremity, e.g. at mid calf. Impedance plethysmography is a non-invasive method used to detect venous thrombosis in these areas of the body.

#### **Genitals**

Another common type of plethysmograph is the penile plethysmograph. This device is used to measure changes in blood flow in the penis. Although some researchers use this device to assess sexual arousal and sexual orientation, the data are usually not admissible in court cases in the United States. An approximate female equivalent to penile plethysmography is vaginal photoplethysmography, which optically measures blood flow in the vagina. [2]

### **Use in preclinical research**

Plethysmography is a widely used method in basic and preclinical research to study respiration. Several techniques are used:

# Respiratory parameters from conscious freely moving animals: whole-body plethysmography

Whole-body plethysmography is used to measure respiratory parameters in conscious unrestrained subjects, including quantification of bronchoconstriction.

The standard plethysmograph sizes are for the study of mice, rats and guinea pigs. On request, larger plethysmographs can also be manufactured for other animals, such as rabbits, dogs, pigs, or primates.

The plethysmograph has two chambers, each fitted with a pneumotachograph. The subject is placed in one of them (subject chamber) and the other remains empty (reference chamber).

The pressure change is measured by a differential pressure transducer with one port exposed to the subject chamber and the other to the reference chamber.

# Respiratory parameters from conscious restrained animals: double-chamber / head-out plethysmography

The double-chamber plethysmograph (dcp) measures respiratory parameters in a conscious restrained subject, including airway resistance and conductance. Different sizes of plethysmograph exist to study mice, rats or guinea pigs.

The head-out configuration is identical to the standard configuration described above except that there is no head chamber.

Of course the collar seal is still applied, so that the body chamber remains airtight. With only a thoracic signal, all parameters can be obtained except for specific airway resistance (SRaw) and specific airway conductance (Sgaw).

Plethysmograph 3

### Resistance/compliance from sedated animals

In anesthetized plethysmography, lung resistance and dynamic compliance are measured directly because the subject is anesthetized.

Depending on the level of sedation, the subject may be spontaneously breathing (SB configuration) or under mechanical ventilation (MV configuration). A flow signal and a pressure signal are required to calculate compliance and resistance.

### References

- [1] http://www.nlm.nih.gov/cgi/mesh/2011/MB\_cgi?field=uid&term=D010991
- [2] Daniel Bergner (2009-01-22). "What Do Women Want?" (http://www.nytimes.com/2009/01/25/magazine/25desire-t.html?hp). New York Times.

### **Further reading**

 Glaab T, Taube C, Braun A, Mitzner W (2007) Invasive and noninvasive methods for studying pulmonary function in mice. *Respiratory Research* 8:63 (free full text article (http://respiratory-research.com/content/8/ 1/63))

### **External links**

- Information on pulmonary function testing and physician reimbursement (http://www.cardiopulmonaryservices.com/)
- Whole Body Plethysmography for animal research (http://mednt.jp/whole+body+plethysmography)

### **Article Sources and Contributors**

Plethysmograph Source: http://en.wikipedia.org/w/index.php?oldid=432966136 Contributors: Andthu, Arcadian, Arthena, Bender235, Billymac00, Bryan Derksen, DSITelemetry, Datasci, Evands, Finchbook01, Gigacephalus, Gogo Dodo, Goplat, Hooperbloob, JP75, Jmabel, Johnpseudo, Jokestress, Lamronby, Michael Kaine, NCurse, Nraesmith, Pinethicket, Rrburke, Shureg, Stevenfruitsmaak, Vice regent, Vivonoetics, WhatamIdoing, Wikidemon, Xiaphias, Yamamoto Ichiro, 14 anonymous edits

## **Image Sources, Licenses and Contributors**

Image:Plethysmograph or BodyBox Empty.jpg Source: http://en.wikipedia.org/w/index.php?title=File:Plethysmograph\_or\_BodyBox\_Empty.jpg License: Creative Commons Attribution-Sharealike 3.0 Contributors: Original uploader was Finchbook01 at en.wikipedia (Original text: Finchbook01 (talk))

Image:Body plethysmography male subject.jpg Source: http://en.wikipedia.org/w/index.php?title=File:Body\_plethysmography\_male\_subject.jpg License: Creative Commons Attribution-Sharealike 3.0 Contributors: Finchbook01 (talk) Original uploader was Finchbook01 at en.wikipedia

### License

Creative Commons Attribution-Share Alike 3.0 Unported http://creativecommons.org/licenses/by-sa/3.0/