



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 227

October 1995

KODAK X-OMAT 480 RA Processor

Modification No. 1 for the *Kodak X-Omat 480 RA PROCESSOR*

The purpose of this Mod is to improve the reliability of the LEVEL PROBES. The new encapsulated PROBE addresses the problems of poor electrical connections caused mainly by corrosion of the CONNECTORS at the old PROBE assembly. The kit contains the new PROBE and cabling as well as the accessories and instructions necessary to install them.

This is a Type 1 Selective Modification, it should only be installed on processors that are having problems with the level sensing.

The modification is now available from Service Parts Management, as Modification Kit 5B6532.

Processors with serial numbers 2812 - 2815 and 2810 and below were shipped without this modification installed.

This modification has been incorporated into production units now being shipped from the factory.

For more information, please contact:

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg 20
Rochester, NY 14650-1131

(800) 336-4722, Press 5

Please be prepared to give your Customer Number and your equipment serial number.

Kodak and X-Omat are trademarks.

N-939

