



# SERVICE BULLETIN

## KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 226

October 1995

### *Kodak AUTOMATIC FILM FEEDER, MODEL RA*



#### **DIP SWITCH Location Changed on the ADJUSTMENT CIRCUIT BOARD**

To reduce the effect of electrostatic discharge on the FILM FEEDER, the DIP SWITCHES have been relocated on the ADJUSTMENT CIRCUIT BOARD. On current units, the SWITCHES are located on the ADJUSTMENT CIRCUIT BOARD so that they are inside the ELECTRICAL ENCLOSURE. See the illustration. The DIP SWITCHES are not accessible through the FRONT PANEL.

This information is covered in Publication Change Notice No. 2 (Pub No. 5B6743) for the Service Guide (Pub No. 240977) and can be ordered from Advertising and Distribution at 1-800-677-9933. The binder complete Publication Number for the *Kodak AUTOMATIC FILM FEEDER, MODEL RA* is Pub No. 240993.

#### **Note**

On FEEDERS with Serial No. 50,899 and above, the DIP SWITCHES are inside the ELECTRICAL ENCLOSURE. See the figures below.



**For more information, please contact:**

Eastman Kodak Company  
Health Sciences Division  
Imaging Service Center  
343 State Street, 3rd fl, Bldg 20  
Rochester, NY 14650-1131

(800) 336-4722, Press 5

Please be prepared to give your Customer Number and your equipment serial number.

*Kodak and X-Omat are trademarks.*

N-938

