



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company/Customer & Technical Services/Health Sciences Markets Division/Rochester, New York 14650

SERVICE BULLETIN NO. 209

October 1994

Kodak X-Omat PROCESSOR, Model M8

Ordering a New CAPACITOR for the RECIRCULATION PUMP MOTOR

In 1982 the supplier for the RECIRCULATION PUMP MOTOR used in the M8 PROCESSOR changed from MAC to Baldor. The assembly part number for the RECIRCULATION PUMP MOTOR, however, did not change. The part number for the PUMP MOTOR is still 594391.

The PUMP MOTORS supplied by Baldor after 1982 require a different CAPACITOR than the PUMP MOTORS supplied by MAC before 1982. So when ordering a new CAPACITOR, identify which PUMP MOTOR is used in your PROCESSOR and use the table below to determine which CAPACITOR to order.

Motor Supplier	Capacitance	Capacitor Part Number
MAC	3 μ F	480190
Baldor	4 μ F	477871

For more information, please contact:

Eastman Kodak Company
Imaging Service Center
Health Sciences Division
100 Carlson Road
Rochester, New York 14653-9013
(716) 724-1789, Option 1

Please be prepared to give your Service Customer Number and your equipment serial number.

svcbul_209.txt
N-921

Kodak and *X-Omat* are trademarks.

Printed in USA

Health Sciences Division
EASTMAN KODAK COMPANY • ROCHESTER, N.Y. 14650

