



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 233

April 1996

3229 KODAK X-OMAT M43A Processor
3255 KODAK X-OMAT M43 Processor
3406 KODAK X-OMAT Clinic 1 Processor

Announcement of Kit No. 9B8328 for the *Kodak X-Omat* M43, M43A, and Clinic 1 Processors

The insulation on the electrical box WIRES in Processors with the serial numbers listed below could be damaged by rubbing against the TANK SUPPORT. Installation of the Kit, Part No. 9B8328 available from Service Parts Management, will remove the potential for damaging the WIRES.

Serial No.	Processor
595 and below	M43
506 and below	M43A
703 and below	Clinic 1

If you have ideas for Service Bulletin Articles, contact either the Health Imaging Customer Support Organization at 1-800-336-4722 or TAC at 1-800-822-1414.

Kodak and X-Omat are trademarks.

8B7487

