



# SERVICE BULLETIN

## KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 219

JULY 1995

***Kodak X-Omat M43 PROCESSOR***  
***Kodak X-Omat M43A PROCESSOR***  
***Kodak X-Omat Clinic 1 PROCESSOR***

### **Availability of the Film Transport Enhancement Kit 5B6960**

To improve film transport in the WASH RACK and keep the SQUEEGEE ROLLER in place, a Film Transport Enhancement Kit 5B6960 is available as of June 26, 1995, for the products and serial numbers listed below. This kit is free for one year for the products and serial numbers listed.

<b>Below Serial No.</b>	<b>Processor</b>
457	M43A
480	M43
496	Clinic 1

### **Replacement for the Drain Valve Stem 918240 and Drain Cap Assembly 1C4841**

To improve the performance of the DRAIN VALVES, the existing STEM and CAP will be replaced with STEM ASSEMBLY 5B6460, available July 21, 1995. The STEM ASSEMBLY includes both the STEM and CAP. The new DRAIN VALVES use a  $\frac{1}{4}$ -turn design.

Thank you for your continued service support of Kodak products.

#### **For more information, please contact:**

Eastman Kodak Company  
Health Sciences Division  
Imaging Service Center  
343 State Street  
Rochester, NY 14650-1131  
(800) 336-4722, Press 5

Please be prepared to give your Customer Number and your equipment serial number.

*Kodak and X-Omat are trademarks.*

N-934

