



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Sciences Division, Imaging Service Center, Rochester, New York 14650

SERVICE BULLETIN NO. 216

May 1995

KODAK X-OMAT 480 RA Processor

Dryer Transport Enhancement Kit Available

The purpose of the Dryer Transport Enhancement Kit is to provide increased drive to the DRYER EXIT TURNAROUND ASSEMBLY on the *Kodak X-Omat* 480 RA PROCESSOR. This will be helpful for processors that have experienced film jams in the DRYER EXIT TURNAROUND ASSEMBLY, or have difficulty transporting roll film.

Installation time is approximately 1 hour.

The Kit will be available on April 17, 1995, and will be free of charge until September 1, 1995. Order Kit 5B6862 from Eastman Kodak Service Parts Management, Rochester, New York.

KIT 5B6862 includes the following:

Quantity	Part No.	Description
1	1C7648	BEARING BLOCK ASSEMBLY
1	5B6865	SCREW - Idler Pulley
1	852623	NUT
1	535009	LOCK WASHER
1	107485	WASHER
1	5B6028	DRIVE ROLLER ASSEMBLY
1	1C7649	CAP - Bearing Block
1	5B6350	DRIVE BELT
1	5B6349	Installation Instructions

For more information, please contact:

Eastman Kodak Company
Health Sciences Division
Imaging Service Center
343 State Street, 3rd fl, Bldg. 20
Rochester, NY 14650-1131

(716) 724-1789, Option 1

Please be prepared to give your Customer Service Number and your equipment serial number.

Kodak and X-Omat are trademarks.

N-216

