



SERVICE BULLETIN

KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging, Rochester, New York 14650

SERVICE BULLETIN NO. 274

October 1997

3058 KODAK X-OMAT Multiloader 300
3059 KODAK X-OMAT 270 RA Processor
3198 KODAK X-OMAT 180 LPS Processor
3227 KODAK X-OMAT 180 LP Processor
3255 KODAK X-OMAT M43 Processor
3229 KODAK X-OMAT M43A Processor
3406 KODAK X-OMAT Clinic 1 Processor
3434 KODAK X-OMAT 3000 RA Processor

New Tan BEARING, Part No. 9B8719 to Replace the White BEARING, Part No. 779119

We are beginning to use the new Tan BEARING in production line Processors and a replacement part is available in Service Parts Management for existing Processors. A new Part No., 9B8719 has been assigned for the new Tan BEARING and the old White BEARING, Part No. 779119 will sub over to it. The new material has shown a significant increase in SHAFT and BEARING life during the trade trials. This will be a direct replacement for the old White BEARING.

Although the new and old style BEARINGS can be intermixed in a RACK without any problem, we recommend the Service Provider replace all the old style BEARINGS at the next scheduled service call on a given Processor. You will see a significant increase in the service life over the original BEARINGS and reduced wear on the SHAFTS. When replacing the BEARINGS, the SHAFT should be checked for any signs of scoring which could cause premature wear of the BEARING. Replace the SHAFT if there is any significant wear. Please note this scoring would affect either type of bearing material.

Processors will be manufactured with the new BEARINGS starting with the following Serial Numbers:

<i>Kodak X-Omat M43 Processor</i>	Serial No. 840
<i>Kodak X-Omat M43A Processor</i>	Serial No. 575
<i>Kodak X-Omat Clinic I Processor</i>	Serial No. 65
<i>Kodak X-Omat 180 LP Processor</i>	Serial No. 3875
<i>Kodak X-Omat 180 LPS Processor</i>	Serial No. 13,405
<i>Kodak X-Omat 3000 RA Processor</i>	Serial No. 4080
<i>Kodak X-Omat Multiloader 300</i>	Serial No. 1720*

*Serial number is for the *Kodak X-Omat 270 RA Processor*

We are also in the process of upgrading the White BEARING, Part No. 911684 for the TOP EXIT RACK of the *Kodak X-Omat* 180 LP Processor, *Kodak X-Omat* 180 LPS Processor, *Kodak X-Omat* 270 RA Processor and the *Kodak X-Omat* 3000 RA Processor. We will implement this after we implement the above change. The original White BEARING, Part No. 911684 will then automatically sub to the new Tan BEARING in Parts Service.

This BEARING will be phased into production by December. We will not publish a serial number list for this change.

The Processors affected by this change are:

Kodak X-Omat 180 LP Processor

Kodak X-Omat 180 LPS Processor

Kodak X-Omat 270 RA Processor

Kodak X-Omat 3000 RA Processor

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

CES Personnel Contact:

Film Handling TAC
U.S.: 1-800-822-1414
Canada: 1-800-433-1414
International: 716-724-4675

Kodak and *X-Omat* are trademarks.

7C8284

