



# SERVICE BULLETIN

## KODAK X-OMAT Processor

Eastman Kodak Company, Health Imaging, Rochester, New York 14650

SERVICE BULLETIN NO. 267

August 1997

**3041 KODAK M35 X-OMAT Processor**  
**3047 KODAK M35A X-OMAT Processor**  
**3408 KODAK M35-M X-OMAT Processor**  
**3133 KODAK M35A-M X-OMAT Processor**

### SMOOTH GUIDE SHOE Modification Kit Now Available

The main purpose of the SMOOTH GUIDE SHOES in this Modification Kit is to reduce GUIDE SHOE artifacts when processing *Kodak Min-R* 2000 Film emulsion down.

#### **Modification Kit Part No. 7C7865**

Processor	Mod. No.	Type of Mod.	Serial No.
M35-M	3	Type 1 Required	76,995 and below
M35A-M	6	Type 1 Required	53,244 and below
M35	2	Type 1 Selective	All
M35A	2	Type 1 Selective	All

The Type 1 Required Modification is for the **M35-M** and **M35A-M** Processors. All service providers with these processors are expected to obtain the Modification Kit and have it installed by September 1, 1998. The Kit will be free from Service Parts Management and upon return of the Modification Reimbursement Card, their account will be credited \$45.00 for 1 hour of labor for each Modification installed.

The Type 1 Selective Modification is for the **M35** and **M35A** Processors. All service providers with these processors are expected to obtain the Modification kit if they are using or plan to use *Kodak Min-R* 2000 Film. It will remain a free Modification until September 1, 1998. The Kit will be free from Service Parts Management and upon return of the Modification Reimbursement Card, their account will be credited \$45.00 for 1 hour of labor for each Modification installed.

The **M35** and **M35A** Processors are **not** being modified in current production and the Modification Labels will not have a number circled. However, if a unit is modified in the field, be sure to circle the number of the Modification.

Please order only as many Kits as you can install at this time to prevent a backorder situation so that all service providers will have access to the current availability.

Christopher Sanzone  
Health Imaging Services

**For More Information:**

**U.S. Distributors Contact:**

Health Imaging  
Customer Support Operations  
1-800-336-4722

**CES Personnel Contact:**

Film Handling TAC  
U.S.: 1-800-822-1414  
Canada: 1-800-433-1414  
International: 716-724-4675

*Kodak, Min-R, and X-Omat are trademarks.*

7C7889

