



SERVICE BULLETIN

KODAK Products

Eastman Kodak Company, Health Imaging, Rochester, New York 14650

SERVICE BULLETIN NO. 264

SEPTEMBER 1997

KODAK X-OMAT Processors KODAK MULTILOADERS KODAK Side-By-Side

List of Current Software

- This publication is not intended as a notification to upgrade your equipment. It is meant to be a handy reference of the latest software.
- All software is available at Service Parts Management.

PART NO.	DESCRIPTION	VERSION
Common RA Software		
990029	5600 BOARD - MAIN PROGRAM (U12)	2.50
5B6278	PROCESSOR DIAGNOSTICS DISK - except 3000RA & 5000RA	2.55
699551	500 BOARD - BOOTSTRAP (U20) - except 180 LP and 180 LPS	2.50
480 RA		
2B6500	Software Download DISK	2.50
990034	PROM SETS (U17, U18, U19) - except Japan	2.50
1C7613	Japanese PROM SETS (U17, U718, U19)	2.51
8B7098	MIN R 2000 Film Download DISK	3.00
8B7094	MIN R 2000 Film PROM SETS (U17, U18, U19) - except Japan	3.00
8B7095	MIN R 2000 Film PROM SETS (U17, U18, U19) - Japan	3.00
460 RA		
2B6501	Software Download DISK	2.50
990033	PROM SETS (U17, U18, U19) - except Japan	2.50
1C7608	Japanese PROM SETS (U17, U18, U19)	2.51
9B8530	MIN R 2000 Film Download DISK	3.00
9B8526	MIN R 2000 Film PROM SETS (U17, U18, U19) - except Japan	3.00
9B8528	MIN R 2000 Film PROM SETS (U17, U18, U19) - Japan	3.00

PART NO.	DESCRIPTION	VERSION
270 RA		
1C7618	Software Download DISK	2.57
5B6283	PROM SETS (U17, U18, U19) - except Japan	2.57
1C7622	Japanese PROM SETS (U17, U18, U19)	2.57
9B8537	MIN R 2000 Film Download DISK	3.00
9B8533	MIN R 2000 Film PROM SETS (U17, U18, U19) - except Japan	3.00
9B8535	MIN R 2000 Film PROM SETS (U17, U18, U19) - Japan	3.00
180 LP and 180 LPS		
5B6623	Software Download DISK - units with SORTER must have SORTER Mod 1 installed	2.24
5B6618	PROM MAIN PROGRAM (U19) - units with SORTER must have SORTER Mod 1 installed	2.24
5B6607	PROM BOOTSTRAP (U20) - units with SORTER must have SORTER Mod 1 installed	2.52
M43, M43A, Clinic 1		
1C7790	M43, M43A PROM MAIN PROGRAM (U5) - for S/N 349 and below, without Mod 1 installed	1.72
1C7791	Clinic 1 PROM MAIN PROGRAM (U5) - for S/N 349 and below, without Mod 1 installed	1.72
8B6944	M43, M43A PROM MAIN PROGRAM (U5) - for S/N 350 and above or with Mod 1 installed	1.90
8B6943	Clinic 1 PROM MAIN PROGRAM (U5) - for S/N 350 and above or with Mod 1 installed	1.90
3000 RA		
9B8449	Software Download DISK (no PROM SETS available for 5000 or 6000 BOARD)	3.31
8B6711	5000 BOARD - BOOTSTRAP (U6)	2.00
8B6712	6000 BOARD - BOOTSTRAP (U18)	2.00
5000 RA		
9B8754	Software Download DISK (no PROM SETS available for 5000 or 6000 BOARD)	4.15
8B6711	5000 BOARD - BOOTSTRAP (U6)	2.00
8B6712	6000 BOARD - BOOTSTRAP (U18)	2.00
Side-By-Side		
9264176	Side-by-Side Software	1.5
MULTILOADER 300		
9219216	Operating Software - Mod 46	3.31
MULTILOADER 700 Plus		
9285016	AØ PCB - U33, U34, U35, A2PCB-U3	5.1
MULTILOADER 700		
9193716	AØ PCB - U33	4.2
9193726	AØ PCB - U34	4.2
9193786	AØ PCB - U35	4.2
9196766	A2 PCB - U3	3.2

Larry Weishaar, Service Engineer

For More Information:

U.S. Distributors Contact:

Health Imaging
Customer Support Operations
1-800-336-4722

GCS&S Contact:

Film Handling TAC

U.S.: 1-800-822-1414

Canada: 1-800-433-1414

International: 716-724-4675

Kodak and X-Omat are trademarks.

7C7876

