




Technical Service Information 2003

Models 8500, 8500M, 8500A, and 8500MA



English



REV	EN	DATE	ORIG	APPD
J	1576	12/22/97	MJI	KEW
K	1608	04/09/98	RSR	KEW
L	1797	12/04/98	TDA	KEW

NONIN MEDICAL, INC. CONFIDENTIAL
 THE DESIGN AND CONSTRUCTION OF THIS DEVICE IS THE PROPERTY OF NONIN MEDICAL, INC. AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE EXPRESS WRITTEN PERMISSION OF NONIN MEDICAL, INC.
 TITLE 8500 HANDHELD WITH MEMORY
 REF 85500
 PART NUMBER 0825-000
 REV L
 SH1 1 OF 1

LEGEND:
 ALTERNATE COMPONENTS USED FOR U3
 THE CONNECTING TRACE BETWEEN U2-28 AND U2-23 IS REMOVED/OUT ON MEMORY UNITS ONLY.
 U1 - U9 AND U11 - U16 ARE LOCATED UNDER D51 - D56.
 COMPONENTS USED ON MEMORY UNITS ONLY.
 COMPONENTS INSIDE DASHED AREA ARE ONLY USED ON ALARM UNITS.

DATE: 10/15/98
 DRAWN: TDA
 CHECKED: DG
 DATE: 01/23/91
 CHECKED: PI
 DATE: 01/23/91
 CHECKED: VH

NONIN MEDICAL, INC. CONFIDENTIAL
 THE DESIGN AND CONSTRUCTION OF THIS DEVICE IS THE PROPERTY OF NONIN MEDICAL, INC. AND IS NOT TO BE REPRODUCED OR TRANSMITTED IN ANY FORM OR BY ANY MEANS, ELECTRONIC OR MECHANICAL, INCLUDING PHOTOCOPYING, RECORDING, OR BY ANY INFORMATION STORAGE AND RETRIEVAL SYSTEM, WITHOUT THE EXPRESS WRITTEN PERMISSION OF NONIN MEDICAL, INC.
 TITLE 8500 HANDHELD WITH MEMORY
 REF 85500
 PART NUMBER 0825-000
 REV L
 SH1 1 OF 1

